



VistA Blood Establishment Computer Software (VBECS) Version 1.5.0.0

Technical Manual-Security Guide

July 2010

Department of Veterans Affairs
Office of Enterprise Development

This page intentionally left blank.

Revision History

Date	Revision	Description	Author
2-3-09	1.0	<p>Updates made to VBECS 1.4.0.0 v 3.0 Technical Manual-Security Guide for VBECS 1.5.0.0:</p> <p>Server Hardware and System Configuration, Printers: Added text stating that printer names and drivers must be consistent.</p> <p>Troubleshooting: Added a section called Printing fails to Report Printer.</p> <p>Implementation and Maintenance: Added instructions for the shut down of VBECS systems (System Shut down Instructions).</p> <p>Global: Made caption titles consistent. Checked sentences for spacing. Removed letter from appendix references to make them more general. Changed "disc" to "disk." Changed "shutdown" to "shut down."</p> <p>Fixed page numbers so that the Introduction appears on page 1.</p> <p>Table 1: Added reference to 9th, spare disk. Added row to account for Integrated Lights Out.</p> <p>Table 3: Changed screen resolution to 1024x768.</p> <p>Periodic Maintenance Checks Table: Included additional information related to checking the database integrity log. Made sure the correct name of the log file was in the document throughout.</p> <p>External Interfaces: Added text stating that services are cluster aware.</p> <p>Configure VistALink Parameters, Process Flow, step 3: Add instructions for viewing error message when connection fails.</p> <p>System Shut Down and Restart Instructions: Adjusted wording throughout for clarity.</p> <p>Troubleshooting: Removed "VBECS Services" section.</p> <p>Installation Time Tasks: Removed "Create Server Accounts" section.</p> <p>Adjusted "Reinstall the System" and "Reconfiguring the VBECS HL7 Multi Listener Service" sections to account for cluster awareness.</p>	BBM team
5-4-09	2.0	<p>Incorporated the changes made in VBECS 1.4.0.0 Technical Manual-Security Guide Version 4.0.</p> <p>Added a section called "Installing a Printer" to account for adding a new laser printer.</p> <p>External Interfaces, Computerized Patient Record System, VistA Patient Updates and VistA Patient Merges: added sentence to describe the supported HL7 versions.</p> <p>Updated Table 7. – Added rows for VBECS: Order Alerts and Pending Order List, VBECS: Patient Update Alerts, VBECS: Patient Merge Alerts under Possible Cause and Solution, describing handling of invalid patient names in HL7 messages to VBECS.</p> <p>Configure Interfaces, Configure CPRS HL7 Interface Parameters: updated step 2 and 3.</p> <p>Configure Interfaces, Configure Patient Update HL7 Interface Parameters: updated step 3.</p> <p>Configure Interfaces, Configure Patient Merge HL7 Interface Parameters: updated step 3.</p> <p>External Interfaces section: Updated the HL7 Service section to include VistALink listener updates.</p> <p>Configure Users, Caution Box: A note was added stating that the user's Windows login ID must not be changed after being configured in VBECS. Added Appendix F: Database Conversion Updates.</p>	BBM team

Date	Revision	Description	Author
		Introduction section: Changed the first caution box so the statement is consistent with other parts of the documentation.	
8-14-09	3.0	<p>Server Configuration section:</p> <p>Added second warning box containing information about adding VBECS Servers to sites exclusion lists.</p> <p>Added third warning box containing information about network requirements for the VBECS servers</p> <p>Under Implementation and Maintenance:</p> <p>Renamed "Periodic Maintenance Checks" to "Periodic System Maintenance".</p> <p>Renamed column headers for Table 5.</p> <p>Added Monitor MOM alerts action to Table 5.</p> <p>Updated Description for Windows Updates, Firmware Updates and VBECS Updates actions in Table 5.</p> <p>Added SQL Maintenance Jobs section.</p> <p>Added SQL Database Job Alerts.</p> <p>Added Figures 28, 29 and 30.</p> <p>Configure Users, Assumptions section:</p> <p>Added a 4th bullet with information to verify application configuration settings.</p> <p>Added a 7th bullet with information about assigning VBECS VISTALINK CONTEXT as a secondary option for all users of the Blood Bank medical device software.</p> <p>Transmit Workload Data: Additional Information section: Added a 3rd bullet with information for Workload multipliers.</p> <p>Updated Figure 89.</p> <p>Troubleshooting Section:</p> <p>Added section for Restarting VBECS Services.</p> <p>Added Figures 94 and 95.</p>	BBM Team

Date	Revision	Description	Author
3-30-10	4.0	<p>Introduction section: Reworded the 2nd sentence of the 1st caution box so the statement is consistent with other VBECS documentation and satisfy security auditors.</p> <p>Firmware Updates section: Added a sentence about the hardware occasionally requiring firmware updates.</p> <p>Added clarifications about Cluster Administrator to the VBECS Windows Services section and to the Reconfiguring the VBECS HL7 Multi Listener and VistALink Services section.</p> <p>Changed the startup type to manual for services listed in Table 8: Windows Service Manager.</p> <p>Configure Interfaces section: Added a statement that if the Facility ID is not supplied, messaging to VBECS will fail.</p> <p>Troubleshooting section: Added a Performance Improvements section that covers stopping and starting of the test services and verification of NIC Card Configuration.</p> <p>Troubleshooting section: Added two new parts: Zebra Printer Problems and Scanner Problems.</p> <p>Integrated Lights Out section: Added a new sub-section for installing iLO.</p> <p>Configure Patient Update HL7 Interface Parameters section, Row 2 of table under Notes, 2nd paragraph: Added last sentence about messaging to VBECS failing if Facility ID is not supplied. Removed the sentences about the field only being validated when using an interface engine to assist with routing HL7 messages, and about the HL7 interface not requiring the use of an interface engine.</p> <p>Configure Patient Merge HL7 Interface Parameters section, Row 2 of table under Notes, 2nd paragraph: Added last sentence about messaging to VBECS failing if Facility ID is not supplied. Removed the sentences about the field only being validated when using an interface engine to assist with routing HL7 messages, and about the HL7 interface not requiring the use of an interface engine.</p> <p>Appendix E: Updated VLAN instructions with new ePolicy servers.</p> <p>Added Appendix G to describe which services are allowed to run on VBECS servers.</p> <p>Added Appendix H to describe what is audited on VBECS servers.</p> <p>Remote Desktop Configuration section, Changed step 6 from (Click, hold, and slide the pointer to a screen resolution of 1024 by 768 pixels) to (Click, hold, and slide the pointer to a screen resolution of Full Screen.) Updatd Figure 3.</p>	BBM Team
7-12-10	5.0	<p>Modified VistA Blood Establishment Computer Software (VBECS) Technical Manual-Security Guide for VBECS 1.5.0.0, Version 4.0:</p> <p>Replaced "Date software turned over from VHIT to VA Product Support" with "July 2010" on the title page.</p> <p>Replaced "March 2010" with "July 2010" in the footer.</p> <p>Replaced "4.0" with "5.0" in the footer.</p>	BBM Team

This page intentionally left blank.

Table of Contents

REVISION HISTORY	I
INTRODUCTION.....	1
RELATED MANUALS AND REFERENCE MATERIALS	1
HOW THIS TECHNICAL MANUAL-SECURITY GUIDE IS ORGANIZED	3
Terms.....	3
Figures and Tables	3
Screen Shots	3
Appendices	3
REMOTE DESKTOP CONFIGURATION.....	5
SCREEN RESOLUTION	5
SOUND	7
CONNECTION SPEED	8
SAVE SETTINGS	9
CREATE A REMOTE DESKTOP CONNECTION SHORTCUT FOR VBECS	10
SERVER HARDWARE AND SYSTEM CONFIGURATION	11
SERVER AND SHARED ARRAY DISKS	11
Server Disk Configuration.....	11
Shared Array Configuration	11
Replacing a Disk	12
PRINTERS	12
Laser Printer	12
Label Printer.....	19
SCANNERS.....	20
SERVER CONFIGURATION	23
REQUIRED HARDWARE.....	24
WORKSTATION CONFIGURATION	24
OFF-THE-SHELF SOFTWARE REQUIREMENTS	24
IMPLEMENTATION AND MAINTENANCE	25
PERIODIC SYSTEM MAINTENANCE	25
SQL MAINTENANCE JOBS	26
SQL Database Job Alerts	26
WINDOWS UPDATES	28
EPOLICY AND VIRUS DEFINITIONS	29
COMMONLY USED SYSTEM RULES.....	29
FIRMWARE UPDATES	30
HARDWARE UTILITIES AND BACKUP EXEC ALERTS	30
HP Event Notifier.....	30
HP System Utilities	33
Backup Exec Alerts	36

INTEGRATED LIGHTS OUT	40
To install iLO	40
To access iLO.....	45
SYSTEM SHUT DOWN AND RESTART INSTRUCTIONS	50
To shut down the system.....	50
To start the system.....	51
MAINTENANCE OPERATIONS.....	53
CONFIGURE INTERFACES	56
CONFIGURE DIVISIONS	65
CONFIGURE SYSTEM ADMINISTRATORS.....	73
CONFIGURE USERS	76
TRANSMIT WORKLOAD DATA	85
NOTIFY VBECS CENTRAL ADMINISTRATOR	86
EXTERNAL INTERFACES.....	87
HEALTH LEVEL SEVEN INTERFACES	87
Client-Server	87
Transport Layers and Lower Layer Protocols	88
TCP Client (Sender)	88
TCP Server (Listener)	89
Computerized Patient Record System.....	89
VistA Patient Updates	89
VistA Patient Merges	89
VISTALINK REMOTE PROCEDURE CALLS	89
VBECS WINDOWS SERVICES	91
RECONFIGURING THE VBECS HL7 MULTI LISTENER AND VISTALINK SERVICES	92
VBECS HL7 Multi Listener Service (Test)	92
VBECS VistALink Service (Test)	93
TROUBLESHOOTING	97
Performance Improvements	97
Stopping and Starting VBECS Test Services.....	97
Verify NIC Card Configuration	100
VistA Query Timeout.....	108
VBECS Exception Logging	110
VBECS Exception Workarounds	111
Restarting VBECS Services	112
VBECS Application Interfaces	114
VBECS Build Version Numbers.....	117
Cluster Connectivity Lost.....	117
Printing Fails to Report Printer	117
Zebra Printer Problems.....	119
Scanner Problems	121
ARCHIVING AND RECOVERY	125
VBECS BACKUP.....	125
VBECS RECOVERY	125
Reinstall the System	126

Inventory the Tape	128
Catalog the Tape.....	129
Restore Files.....	130
Restore the Databases.....	133
FAILOVER	135
PERFORMANCE.....	137
LOCKING.....	137
SECURITY	139
ACTIVE DIRECTORY.....	139
GROUP POLICY	139
VIRTUAL LOCAL AREA NETWORK	139
MICROSOFT OPERATIONS MANAGER	139
APPLICATION-WIDE EXCEPTIONS	140
GLOSSARY	141
APPENDICES.....	143
APPENDIX A: INSTRUCTIONS FOR CAPTURING SCREEN SHOTS.....	143
APPENDIX B: WORKLOAD PROCESS MAPPING TO APPLICATION OPTION TABLE.....	145
APPENDIX C: KNOWN DEFECTS AND ANOMALIES	153
APPENDIX D: ACTIVE DIRECTORY REQUEST FORM.....	155
APPENDIX E: DATA CENTER INSTRUCTIONS	157
Purpose.....	157
Initial Setup Tasks.....	157
Ongoing Tasks.....	161
Installation Time Tasks	162
APPENDIX F: DATABASE CONVERSION UPDATES	163
Warnings and Notifications Displayed by the DTS Package.....	164
APPENDIX G: SERVICES ALLOWED TO RUN ON VBECS SERVERS.....	167
APPENDIX H: AUDITING ON VBECS SERVERS.....	169
INDEX.....	171

This page intentionally left blank.

Introduction

The main purpose of the VistA Blood Establishment Computer Software (VBECS) is to automate the daily processing of blood inventory and patient transfusions in a hospital transfusion service.



Unauthorized access or misuse of this system and/or its data is a federal crime. Use of all data, printed or electronic, must be in accordance with VA policy on security and privacy.



Do not change the system! The U.S. Food and Drug Administration classifies this software as a medical device. Unauthorized modifications will render this device an adulterated medical device under Section 501 of the Medical Device Amendments to the Federal Food, Drug, and Cosmetic Act. Acquiring and implementing this software through the Freedom of Information Act require the implementer to assume total responsibility for the software and become a registered manufacturer of a medical device, subject to FDA regulations. Adding to or updating VBECS software without permission is prohibited.



Changes to the system configuration must be documented with screen captures and kept with the installation record.

Related Manuals and Reference Materials

- *Health Level Seven Implementation Support Guide for HL7 Standard Version 2.3.1, Message & Interface Services (M&IS), VHA OI - Health Systems Design & Development Web site (©1999).*
- *Kernel Systems Manual Version 8.0, Chapter 1: Sign-On Security/User Interface, pp. 13–20.*
- “Locking Down Windows Server 2003 Terminal Server Sessions,” Microsoft Web site (October 29, 2003).
- *National Software Package Distribution, SOP 196-5.*
- *Release of Patches, SOP 196-8.*
- *VBECS Application Interfacing Support Software Installation and User Configuration Guide.*
- *VistA Blood Establishment Computer Software (VBECS) Installation Guide.*
- *VistA Blood Establishment Computer Software (VBECS) User Guide.*
- *VistALink Version 1.0 Developer-System Manager Manual, Chapter 6: Security Management, pp. 34–35.*
- *Windows Server 2003 Security Guide 2.1, Microsoft Corporation (May 8, 2006).*

This page intentionally left blank.

How This Technical Manual-Security Guide Is Organized

Outlined text is used throughout this guide to highlight warnings, limitations, and cautions:



Warnings, limitations, cautions

Terms

For consistency and space considerations, the pronouns “he,” “him,” and “his” are used as pronouns of indeterminate gender equally applicable to males and females.

In many instances, a user may scan a barcode or enter data manually (by typing). The term “enter” is used throughout this guide to mean “enter manually.”

See the Glossary for definitions of other terms and acronyms used in this guide.

Figures and Tables

If you refer to figures and tables from the technical manual-security guide in your local policy and procedure documents, you may wish to use their titles only, without figure or table numbers: as the technical manual-security guide is updated, those numbers may change.

Screen Shots

Because VBECS is a medical device, screen shots must be captured at various points throughout the technical manual-security guide to meet FDA requirements for objective evidence and documentation. A



(camera) at the beginning of each step that requires a screen capture will identify these points. For more information, see Appendix A: Instructions for Capturing Screen Shots.

Appendices

The appendices contain truth tables and other materials for reference.

While pressing the Ctrl button, left click on a section name or page number in the table of contents to move to that section or page. The index does not incorporate this feature.

This page intentionally left blank.

Remote Desktop Configuration


Configure the screen resolution, sound, and connection speed, and create a Remote Desktop Connection shortcut on each VBECS workstation.

Screen Resolution

To set the screen resolution:

- 1) Double click  (the **Remote Desktop Connection** icon).
- 2) Click **Options** (Figure 1).

Figure 1: Remote Desktop Connection Options




- 3) Click the **General** tab (Figure 2).
- 4) Enter the VBECS server cluster name or cluster IP address in the Computer field. Enter **VHAMASTER** in the Domain field. Do not enter a user name or password.

Figure 2: General Tab: Computer and Domain



- 5) Click the **Display** tab (Figure 3).
- 6) Click, hold, and slide the pointer to a screen resolution of Full Screen.

Figure 3: Display Tab



- 7) Click on the **General** tab.
- 8) Click **SAVE** to save the setting.

Sound

To enable sound:

- 1) Click the **Local Resources** tab (Figure 4).
- 2) Select **Bring to this computer** from the Remote computer sound drop-down list.



Failure to properly configure the sound disables audible alerts throughout VBECS.

Figure 4: Remote Computer Sound




Connection Speed

To set the connection speed:

- 1) Click the **Experience** tab (Figure 5).
- 2) Select **LAN (10 Mbps or higher)** from the **Choose your connection speed to optimize performance** drop-down list.

Figure 5: Connection Speed



Save Settings

To save the settings:

- 1) Click the **General** tab (Figure 6).
- 2) Click **Save As**.


Figure 6: General Tab: Save As



Create a Remote Desktop Connection Shortcut for VBECS

- 1) To create a Remote Desktop Connection shortcut for VBECS (Figure 7), save the file as VBECS.rdp in the All Users, Desktop folder.

Figure 7: Remote Desktop Connection Shortcut for VBECS



- 2) Double click the shortcut to launch the remote desktop connection to VBECS.
- 3) The Windows start-up sound confirms that the sound functions.

Server Hardware and System Configuration

The VBECS application requires that hardware and system software serve five users in a standard configuration and up to 25 users in an integrated Veterans Integrated Service Network (VISN) environment.

The System Schematic diagram (Figure 27) describes the major system components: a Windows 2003 Server system (the execution environment for the VBECS application) and Windows XP workstations, with which the user will access the VBECS application using Windows Terminal Services [Remote Desktop Protocol (RDP)]. The VBECS server will also communicate with and exchange information with VistA applications through messages formatted using Extensible Markup Language (XML) and Health Level Seven (HL7) over Transmission Control Protocol/Internet Protocol (TCP/IP) networking.

Server and Shared Array Disks

Server Disk Configuration

Each VBECS server has two disks in a RAID 1 (mirroring) configuration (Figure 8). This means that if one disk fails, the server will continue to run normally.

Figure 8: Server Disks



Shared Array Configuration

The shared disk array consists of nine disks (Figure 9).

- The first four disks are used to store VBECS specific data. These disks are configured as RAID 5.
- The fifth disk is a hot spare. It can be used if one of the other disks fails. Note that the LED on it will be off.
- Disks 6 and 7 are for log storage. These disks are configured as RAID 1.
- Disks 8 and 9 are for cluster support. These disks are configured as RAID 1.

Figure 9: Shared Array



Replacing a Disk

All disks in the system, both server and array, are hot swappable. This means that if a disk should fail, it can be replaced without powering down the system or disrupting users. Simply remove the failing disk and replace it with a new one. It will take a couple of minutes to rebuild. For more information on monitoring and viewing disk health, please see the HP Array Diagnostic Utility section.

Printers

Laser Printer

A laser printer capable of printing 8.5" x 11" sheets may be used. Printer naming and drivers must be consistent across both servers.

Installing a Printer

To install or reinstall a printer, execute the following instructions on each server node:

- 1) Log into the first server with your Windows ID.
- 2) Click **Start, Control Panel, Printers and Faxes, Add Printer.**


- 3) In the Add Printer Wizard screen, click **Next** (Figure 10).

Figure 10: Add Printer Wizard




- 4) Make sure the **Local printer attached to this computer** radio button is selected.
- 5) Uncheck the **Automatically detect and install my Plug and Play printer** check box.
- 6) Click **Next** (Figure 11).

Figure 11: Add Printer Wizard



- 7) Select the **Create a new port** radio button.
- 8) Select **Standard TCP/IP Port** from the drop-down menu. Click **Next** (Figure 12).

Figure 12: Add Printer Wizard




- 9) In the Add Standard TCP/IP Printer Port Wizard screen, click **Next** (Figure 13).

Figure 13: Add Standard TCP/IP Printer Port Wizard



10) Enter the IP address of the printer in the “Printer Name or IP Address” field (the Port Name field will populate automatically). Click **Next** (Figure 14).

Figure 14: Example of TCP/IP Settings



The screenshot shows a Windows-style dialog box titled "Add Standard TCP/IP Printer Port Wizard". The main heading is "Add Port" with a sub-question: "For which device do you want to add a port?". Below this, there is a printer icon. The instruction reads: "Enter the Printer Name or IP address, and a port name for the desired device." There are two text input fields: "Printer Name or IP Address:" containing "10.3.21.150" and "Port Name:" containing "IP_10.3.21.150". At the bottom, there are three buttons: "< Back", "Next >", and "Cancel".

11) Click **Finish** (Figure 15).


Figure 15: Example of Review Settings



The screenshot shows the same dialog box, now at the "Completing the Add Standard TCP/IP Printer Port Wizard" step. It features a large printer icon on the left. The text says: "You have selected a port with the following characteristics." Below this is a list of settings: "SNMP: Yes", "Protocol: RAW, Port 9100", "Device: 10.3.21.150", "Port Name: IP_10.3.21.150", and "Adapter Type: Hewlett Packard Jet Direct". At the bottom, there is a message: "To complete this wizard, click Finish." and three buttons: "< Back", "Finish", and "Cancel".


12) To select a driver, click **Have Disk** (Figure 16). Note: If your site has chosen to use their own printer, you must point to your own driver at this point.

Figure 16: Add Printer Wizard




13) Enter \\10.3.21.77\HP\english\Win32_2000_XP_S2003\PCL6\. Select **hpc9050c.inf**. Click **Open** (Figure 17).

Figure 17: Navigate to the Driver



14) Click **OK** (Figure 18).


Figure 18: Install From Disk



15) Select **HP LaserJet 9040 PCL 6**. Click **Next** (Figure 19).


Make sure that the 9040 driver is selected.

Figure 19: Add Printer Driver Wizard



16) For a single-division site, enter **VBECS Printer** as the printer name. For a multi-divisional site, enter **VBECS Printer** and the site name (e.g., VBECS Printer Hines). Click **Next** (Figure 20).


Figure 20: Example of Add Printer Wizard



The screenshot shows the 'Add Printer Wizard' dialog box with the title bar 'Add Printer Wizard'. The main heading is 'Name Your Printer' with a printer icon to the right. Below the heading is the instruction 'You must assign a name to this printer.' The main area contains a paragraph: 'Type a name for this printer. Because some programs do not support printer and server name combinations of more than 31 characters, it is best to keep the name as short as possible.' Below this is a text input field labeled 'Printer name:' containing the text 'VBECS Printer'. At the bottom are three buttons: '< Back', 'Next >', and 'Cancel'.

17) Click the **Do not share this printer** radio button. Click **Next** (Figure 21).


Figure 21: Add Printer Wizard



The screenshot shows the 'Add Printer Wizard' dialog box with the title bar 'Add Printer Wizard'. The main heading is 'Printer Sharing' with a printer icon to the right. Below the heading is the instruction 'You can share this printer with other network users.' The main area contains a paragraph: 'If you want to share this printer, you must provide a share name. You can use the suggested name or type a new one. The share name will be visible to other network users.' Below this are two radio buttons: 'Do not share this printer' (which is selected) and 'Share name:'. The 'Share name:' radio button is followed by a text input field containing 'VBECSPri'. At the bottom are three buttons: '< Back', 'Next >', and 'Cancel'.

18) Click **Next** (Figure 22).

Figure 22: Add Printer Wizard



19) Repeat these instructions on the other server node.

Label Printer

VBECS is configured to work only with Zebra printers: VBECS uses Zebra printing language to communicate with the printer. Other requirements:

- Ethernet connectivity: the label printer must have an Ethernet card.
- Must print on 4" x 4" label stock
- Must print at 300DPI

Prior to configuring the label printer, load the ribbon and label stock and ensure that the printer is on. If the printer does not display **PRINTER READY**, there is a problem that must be resolved before proceeding. Refer to the Zebra user guide or printer CD for more information.

Set the IP Address on the Printer

- 1) Press **SETUP/EXIT** to access the configuration menus.
- 2) Press + or – to scroll through the configuration menu options. Stop when **IP PROTOCOL** is displayed and press **SELECT**. If there is a prompt for a password, press – to change positions and + to change numbers. Enter **1234**. Press **SELECT**.
- 3) Press + to select **PERMANENT**. Press **SELECT**. The IP address is configured to be static.
- 4) Press + to navigate to the **IP ADDRESS** menu option. Press **SELECT**.
- 5) Press + or – to change numbers (as in Step 2) to enter the IP address specified in the Configuration Checklist. Press **SELECT**.
- 6) Press **SETUP/EXIT** to save the new configuration. **PERMANENT** is displayed. Press **SETUP/EXIT** to save the changes.

Test the Printer

To print a label, press and hold the Network Configuration button (on the back of the printer just above the Ethernet socket) until the DATA LED on the front of the printer blinks. Retain the test label for validation records. If the printer configuration on the label print is blank or faint or it is printing off center, adjust the settings.

Adjust Label Darkness

If the printer configuration on the label print is blank or faint, adjust the darkness:

- 1) Press **SETUP/EXIT**. Press + or – until DARKNESS is displayed. Press **SELECT**.
- 2) Press + to adjust the darkness to a higher number. Press **SELECT**. Move up in small increments: setting the printer to a setting that is too dark may compromise the quality of the labels.
- 3) Repeat these steps to retest the printer.
- 4) If parts of the label are cut off, adjust the X and Y offsets.
- 5) Press **SETUP/EXIT** twice to permanently change the setting.

Adjust Label Offsets

If the printer is printing off center, adjust the X and Y offsets:

- 1) Press **SETUP/EXIT**. Press + or – until LABEL TOP (if vertical alignment is not correct) or LEFT POSITION (if horizontal alignment is not correct) is displayed. Press **SELECT**.
- 2) Press + or – to adjust the alignment to a higher number. Press + in the LABEL TOP menu to move the printing down on the label. Press + in the LEFT POSITION menu to move the printing to the right on the label.
- 3) Press **SELECT**. Adjust in small increments until the label is centered on the label stock.
- 4) Press **SETUP/EXIT** twice to permanently change the setting.

Scanners

Scanners used with VBECS must be able to scan Codabar, ISBT 128, and PDF-417 barcodes. To configure a scanner:

- 1) Connect the scanner to the workstation.
- 2) To configure a Hand Held 4600 barcode scanner, scan the barcode in Figure 23. Repeat for all scanners.

Figure 23: Configure a Barcode Scanner



- 3) To test the scanner, open Notepad. Print and scan the barcodes in Figure 24, Figure 25, and Figure 26. The Codabar and ISBT barcodes must scan as “~123456789”; the PDF 417 must scan as “~Testing.”
- 4) Save and print the Notepad file for validation records.

Figure 24: Codabar



Figure 25: ISBT 128



Figure 26: PDF 417





Figure 27: System Schematic



Server Configuration



The U.S. Food and Drug Administration classifies this software as a medical device. Unauthorized modifications will render this device an adulterated medical device under Section 501 of the Medical Device Amendments to the Federal Food, Drug, and Cosmetic Act. Acquiring and implementing this software through the Freedom of Information Act require the implementer to assume total responsibility for the software and become a registered manufacturer of a medical device, subject to FDA regulations.



VBECS is a medical device; all updates and changes to it must be tested and documented. This will be centrally managed. The VBECS servers must be added to site exclusion lists so they are not part of local update mechanisms. Ensure that login scripts do not run on VBECS servers as they may attempt to install unauthorized software. Do not install the ePolicy agent on the VBECS systems: exclude them from Systems Management Server (SMS) updates. Install Windows updates only after approval is granted.

Table 1: Server Configuration

Hardware	Clustered Database Server (two identical systems)
Processor	Multiple processors (2–4 processors) Pentium 4 Xeon 2.0 GHz processors (or greater) with 512kb level 1 cache
Memory	2-gigabyte (or greater) main storage (RAM)
Storage	Shared Storage Controller Unit. Disk configuration: 8 hot swappable SCSI hard drives (minimum 10,000 RPM). The system drives require 18 gigabytes (or greater) storage capacity. The application data drives require 36 gigabytes; log volume and historical data drives require 72 gigabytes (or greater) storage capacity. A ninth disk has been included that serves as a hot spare. If a live disk should fail, it can be replaced with this one.
Operating System	Microsoft Windows 2003 Server Enterprise with Microsoft Clustering Services providing failover data-device sharing
Network Controller	Multiple 10/100 network cards configured to provide fallback in event of failure.
Power Supply	Primary and secondary (redundant) power supply to server chassis and an uninterruptible power source (UPS)
Backup	Internal tape backup with software
Integrated Lights Out (iLO)	A hardware device attached to the servers that allows for remote management

This configuration is designed to promote 24/7 availability and use of the application. A clustered database server configuration will provide near immediate failover if one node of the server fails. Multiple processors will provide for more efficient processing of database access requests and operating system processes.

Dual power supply and UPS will ensure that the machine will not lose operating power. The disk storage configuration will allow the server disks to be shadowed; if a main disk fails, the shadow disk will automatically continue system operation until the primary disk is replaced. Hot swappable disk drives can be replaced without shutting down the server. Internal tape backup on the application data disk will allow an image of the application data to be restored to another machine if the server is damaged.



VBECS is a Terminal Server application. The VBECS lab workstations are configured to run at LAN speeds (10 Mbps). If your network cannot support this, please file a Remedy ticket.

Required Hardware

Table 2: Required Hardware

Hardware	Description
Zebra Printer	Zebra printer capable of producing barcode labels (network capable)
Barcode Scanner	Symbol Model LS4006i barcode scanner for each workstation
Report Printer	Laser printer or comparable with sufficient speed to handle high-volume reports (network capable)

Workstation Configuration

Table 3: Workstation Configuration

Hardware	Description
Processor	Suitable for Windows XP
Memory	256 megabytes (or greater) main storage (RAM)
Monitor	17" monitor or greater
Video	Video card capable of displaying minimum of 16-bit color at 1024 x 768 resolution
Disk Storage	9 gigabytes (minimum)
Operating System	Microsoft Windows XP Professional with Microsoft Terminal Services Client
Network Controller	10/100 network card
Input Devices	U.S. 101-key keyboard, mouse
Audio	Sound card and speakers (may be internal)

Off-the-Shelf Software Requirements



Do not upgrade, change, or add software to the VBECS server as this may compromise the integrity of VBECS.

Table 4: Off-the-Shelf Software Requirements

Software	Description
.NET Framework	Version 1.1
SQL Server	SQL Server 2000 Enterprise Edition
Crystal Reports	Crystal Reports .NET
Backup software	VERITAS Backup Exec Version 10.0
McAfee VirusScan	Version 8.0

Implementation and Maintenance



The U.S. Food and Drug Administration classifies this software as a medical device. Unauthorized modifications will render this device an adulterated medical device under Section 501 of the Medical Device Amendments to the Federal Food, Drug, and Cosmetic Act. Acquiring and implementing this software through the Freedom of Information Act require the implementer to assume total responsibility for the software and become a registered manufacturer of a medical device, subject to FDA regulations.

Periodic System Maintenance

The system will fail to function as intended when maintenance checks are not performed or are not performed correctly. Follow all instructions in the *VistA Blood Establishment Computer Software (VBECS) Installation Guide* for configuration.

Table 5: Periodic System Maintenance

Action	Frequency	Description
Backup tape rotation	Daily	If using Backup Exec, backups automatically occur every morning per the time specified in the VBECS Installation Guide. Refer to local policy for data retention and offsite storage requirements.
Monitor Microsoft Operations Manager (MOM) Alerts	Daily	MOM emails alert messages to the VBECS Administrators mail group, which is defined in the Installation Guide, as problems occur on the clustered servers. Investigate all alerts to completion.
Review Database Integrity Reports	Weekly (Saturday)	Every Saturday morning, 6 emails are sent with the results of the Database Integrity check jobs. Each email will contain a report that must be manually reviewed for successful completion. See the SQL Maintenance Jobs section for more details.
Windows Updates	2nd Tuesday of the month	A VistA Informational patch is released when the updates have been tested and approved for installation.
Firmware Updates	As needed	A VistA Informational patch is released when the updates have been tested and approved for installation.
VBECS Updates	As needed	A VistA Informational patch is released when the updates have been tested and approved for installation.

SQL Maintenance Jobs

The VBECS databases are contained within Microsoft SQL Server and require regular maintenance jobs to backup, validate integrity, and improve performance. The jobs are automated and configured to run according to the specifications shown in Table 6. The following is a list of the SQL Server databases needed by the VBECS application:

- msdb (contains information relating to the SQL Server jobs)
- master (required for SQL Server and all databases within to operate)
- VBECS_V1_PROD (VBECS production account database)
- VBECS_V1_PROD_MIRROR (VBECS production account audit database)
- VBECS_V1_TEST (VBECS test account database)
- VBECS_V1_TEST_MIRROR (VBECS test account audit database)

Table 6: VBECS SQL Maintenance Jobs


Database Affected	Job Name	Frequency and Time (local time) Job Runs	Description
N/A	ResetServerLogFile	Daily at 12:00:10 am	Truncates and starts a new error log file for SQL server.
All databases	WeeklyIntegrityCheck	Weekly (Saturday) at 12:11:50 am	Checks the physical integrity of the database and generates a report for manual verification.
vbecs_v1_prod vbecs_v1_test	ExpireComponentsOrder	Daily at 1:00:00 am	Expires component orders when the associated specimen expires or other specific criteria are met.
	MarkedPresumedTransfused	Daily at 1:10:00 am	Marks units as presumed transfused if transfusion or bedside verification information was not returned to the blood bank within 48 hours.
All VBECS databases	ShrinkLog	Daily at 1:50:00 am	Removes free space at the end of the database log file.
All databases	DailyBackup	Daily at 2:00:00 am	Full database backup
All VBECS databases	UpdateStats	Daily at 2:20:00 am	Updates statistics on all user defined tables to improve performance.
	TruncateDataFiles	Daily at 2:30:00 am	Removes unused space from database files
All databases	Copy VBECS DB Backups to L Drive	Daily at 2:40:00 am	Copies the latest database backup files to L:\Program Files\Microsoft SQL Server\MSSQL\Backup\ <database> folder and renames the file to include the current date time.
	L Drive Delete old Backup files	Daily at 2:50:00 am	Deletes database backup files that are more than 7 days old.
All VBECS databases	VBECS_V1_PROD_ReIndexTables	Daily at 3:00:00 am	Re-Indexes the database tables to improve performance.

SQL Database Job Alerts

If any of the SQL database jobs should fail, an email alert is sent via SMTP to expedite intervention. Please file a remedy ticket to report the job failure. The email recipients for these alerts are determined by a configuration setting in the VBECS Service Monitor established when VBECS was initially installed.


See Figure 28. The EmailRecipients attribute value listed in the VbecsServiceMonitor.exe.config file (C:\program files\vista\vbecs service monitor) will receive the database job alerts. This is the same recipient list where VBECS service alerts are sent. Any changes made to this file will need to be duplicated on both VBECS servers.

Figure 28: Example VBECS Service Monitor file’s EmailRecipients Setting




The job failure alerts will be formatted similar to the one shown in Figure 29.

Figure 29: Example of a VBECS Job Failure Email



The Weekly Integrity job will send an alert regardless if the job succeeds or fails. This email will contain an attachment (DatabaseIntegrityReport.log) which will need to be manually reviewed to determine if the job completed successfully. To validate the jobs successful completion, open the log file attachment, and verify that the second line from the bottom contains “CHECKDB found 0 allocation errors and 0 consistency errors.” See Figure 30 for an example of a successful integrity report log. If six emails are not received Saturday morning from the Weekly Integrity jobs or a report does not indicate a successful completion, please file a remedy ticket. The reports are also physically stored on the VBECS cluster under D:\Program Files\Microsoft SQL Server\Backup\

Figure 30: Example of Database Integrity Report



Windows Updates

If your servers reside at a data center that has its own update distribution system, please refer to Appendix E: Data Center Instructions.


The VBECS development team must test every Microsoft Windows update. Once the development team is satisfied that the update causes no adverse effects, they will notify sites that there are Windows or Firmware updates. A VistA information patch in the VBEC namespace will be created by the VBECS team each time an update is available describing where to obtain the update and how to apply it. The patch will be released to customers by VA Product Support.

Updates are approved with Windows Software Update Service. Approved updates will be downloaded to your servers automatically. However, a server administrator must install the updates locally.

VA Product Support will notify the sites of updates required for installation.

- 1) Since most updates require a reboot, coordinate a time with the blood bank manager to apply the updates.
- 2) At the agreed upon time, log onto the first server as a user with administrative privileges.
- 3) A shield shaped icon will appear in the System Tray (lower right corner of desktop). Click on it (Figure 31).

Figure 31: Windows Software Update Notification



- 4) Leave **Express Install** selected. Click **Install** (Figure 32).

Figure 32: Example of Automatic Updates



- 5) When the update process is complete, you may be prompted to reboot. If so, reboot the server at this time.
- 6) After the server completely reboots, repeat this process on the second server.

ePolicy and Virus Definitions

The VBECS development team must test virus definitions before they are applied to the servers. The VBECS development team will send the virus definitions: do not apply virus definitions locally.

 *Do not change the system! The U.S. Food and Drug Administration classifies this software as a medical device. Unauthorized modifications will render this device an adulterated medical device under Section 501 of the Medical Device Amendments to the Federal Food, Drug, and Cosmetic Act. Acquiring and implementing this software through the Freedom of Information Act require the implementer to assume total responsibility for the software and become a registered manufacturer of a medical device, subject to FDA regulations. Adding to or updating VBECS software without permission is prohibited.*

Commonly Used System Rules

This section includes system rules that apply to several or all options.

- Only one instance of Configure Interfaces may run at a time.
- VBECS captures changes to verified data for inclusion in the Audit Trail Report.
- VBECS protects application data through encapsulation. Encapsulation promotes data security by hiding the implementation details.

 *The dialogs defined in Configure Interfaces and Configure Divisions cannot run when VBECS is operational. VBECS cannot run when a dialog in these options is operational.*

Firmware Updates


Occasionally, hardware including the server components, printers and scanners require firmware updates. Forum informational patch messages are posted when the updates have been tested and approved for installation.

Hardware Utilities and Backup Exec Alerts

HP Event Notifier

Hardware alerts are generated with HP Event Notifier. Event Notifier will generate email alerts whenever a hardware failure occurs. Examples of hardware failures include, but are not limited to; controller, network interface card and fan failures. An example of a network interface card losing connectivity is displayed in Figure 33.

Figure 33: Example of an Email Alert from Event Notifier



When an alert is received, a server administrator should investigate the problem as soon as possible in order to prevent VBECS downtime. If necessary, contact HP support for assistance at 800.633.2600.

Configuring Event Notifier

To add or modify hardware alerts on servers, take the following steps:

- 1) Log into the server with administrative rights.
- 2) Click **Start, HP Management Agents, Event Notifier Config.**

3) Click **Next** (Figure 34).

Figure 34: Welcome Screen




4) Enter the following (Figure 35):

- From address: <servername>@va.gov
- Mail server: <servername>
- Reply address: <servername>@va.gov


Click **Next**.

Figure 35: Example of SMTP Configuration



- 5) Click **Add** (Figure 36). Note that **Remove** or **Edit** can be used for modification and deletion of existing groups respectively.


Figure 36: Recipients



- 6) Enter the following (Figure 37):
 - Display name: Arbitrary name that describes the email group being entered.
 - E-mail address: Email group address of support personnel (Figure 37). Note: Use the support email address that was defined in the *VBECS Installation Guide* (Appendix: Contact Information).


Click **OK**.

Figure 37: Example of Notification Recipient Information



7) Click **Finish** (Figure 38). Repeat these instructions on the other server.

Figure 38: Example of Event Recipients Information




HP System Utilities

There are several pre-installed utilities on the system that are useful when checking hardware health and diagnosing problems. All of these tools are launched from the **Start** menu and all require administrative rights. Please see HP documentation for specific information regarding further use of any of these tools.

HP System Management Homepage

This tool quickly lets the administrator see the status of all major components of the system including the shared array (Figure 39).


Figure 39: System Management Homepage



HP Array Configuration Utility

This tool shows the state of disks, both server and shared array (Figure 40).


Figure 40: Array Configuration Utility



HP Array Diagnostic Utility

This tool generates a report showing the status of disks, both server and shared array (Figure 41). It is useful for diagnosing disk problems.


Figure 41: Array Diagnostic Utility



HP ProLiant Integrated Log Viewer

All hardware related issues are logged here (Figure 42).


Figure 42: HP ProLiant Integrated Management Log Viewer



Backup Exec Alerts

Backup Exec job failure alerts are sent by Backup Exec. Whenever the nightly job fails, an alert will be sent. An example of one of these alerts is displayed in the screen capture below (Figure 43).

Figure 43: Example of an Email Alert from Backup Exec



When an alert is received, a server administrator should investigate the problem as soon as possible in order to ensure proper data backup.


Configure Backup Exec Alerts

To add or modify Backup Exec Alerts on servers, take the following steps:

- Log into the server (not the cluster) that has Backup Exec installed with administrative rights.
- Click **Start, All Programs, Symantec Backup Exec 10d for Windows Servers.**


- Click Alerts (Figure 44).

Figure 44: Backup Exec Main Screen




- Click Configure e-mail and pagers (Figure 45).

Figure 45: Example of Alerts




- Enter the following (Figure 46):
 - Check the **Enable** box
 - SMTP mail server: <server name>
 - Sender name: <server name>
 - Sender address: <server name>
- Click **OK**.

Figure 46: Example of SMTP Configuration




- Click **Configure recipients** on the main Alerts screen. Click **New** (Figure 47). Note that **Remove** or **Properties** is used for deletion and modification of existing groups respectively.

Figure 47: Configure Recipients



- Click **OK** to select Person (Figure 48).


Figure 48: Recipient Type



- Enter the following (Figure 49):
 - Name: Arbitrary name that describes the email group being entered.
 - Check the **Enable** box.
 - Address: Email group address of support personnel (Note: Use the support email address that was defined in the *VBECS Installation Guide* (Appendix: Contact Information).


Click **OK**.

Figure 49: Example of Recipient Properties



- Click **Close** (Figure 50).

Figure 50: Example of Configure Recipients



Integrated Lights Out


Integrated Lights Out (iLO) is a separate hardware component of the server that allows for increased remote administrative capabilities via a separate network connection. For example, the server can be turned on and diagnostic information can be viewed through the iLO console. For instructions on installing iLO and defining users, please see Appendix: Implementing Integrated Lights Out of the *VBECs Installation Guide*. This section assumes you have already executed those instructions.

To install iLO

- 1) Attach the iLO ports on the back of each server to the VA network with an Ethernet cable.
- 2) Record the following information:
 - IP address for iLO port on server #1: _____
 - IP address for iLO port on server #1: _____
 - Default Gateway: _____
 - Subnet Mask: _____
 - DNS: _____
 - WINS (if applicable): _____


- 3) Log into Server 1 with your Windows ID. Reboot and watch the startup sequence. Press **F8** when prompted (Figure 51).

Figure 51: Press F8




- 4) The iLO configuration screen will appear. With the arrow keys, select **Network**, **DNS/DHCP** and click **Enter** (Figure 52).

Figure 52: DNS/DHCP




- 5) The Network Autoconfiguration screen launches. Turn off DHCP by pressing the space bar. Press **F1** to launch Advanced options (Figure 53).

Figure 53: Disable DHCP




- 6) Disable DHCP in all four options in the top panel. Enter WINS from Step 2 (if applicable) addresses. Enter DNS server addresses from Step 2. Press **F10** to save (Figure 54).

Figure 54: Advanced




7) Select **Network, NIC and TCP/IP** and click **Enter** (Figure 55).

Figure 55: TCP/IP




8) Enter a static IP address, subnet mask and default gateway (from Step 2). Press **F10** to save (Figure 56).

Figure 56: Network



9) Select **User, Add** (Figure 57).

Figure 57: Add user



10) Enter the following (Figure 58):

- User name: Administrator's first and last name
- Login name: Network ID of the administrator
- Password: A complex password consisting of letters, number and special characters with a minimum length of eight.

Press **F10** to save.



Note that the iLO ID and password operate independently of the Windows credentials. Changing the Windows password will not affect the iLO password!

Figure 58: Add a user


A screenshot of a terminal window titled 'Add User'. The form contains the following fields:
User name: Joe Administrator
Login name: vhaminadminj (highlighted with a blue bar)
Password: *****
Verify password: *****
Below these fields is a section titled 'Lights-Out Privileges' with the following settings:
Administer User Accounts: Yes
Virtual Power and Reset: Yes
Configure Settings: Yes
Remote Console Access: Yes
Virtual Media: Yes
At the bottom, it says '[F10] = Save [ESC] = Cancel'.

- 11) Repeat Steps 9 and 10 to add additional administrators.
- 12) Press **Escape** to close the iLO configuration.
- 13) Repeat this entire section on Server 2.

To access iLO

- 1) From any computer in the VA wide area network (WAN), launch a web browser and enter the iLO IP address of the server you would like to administer (Figure 59). Press **Enter**.

Figure 59: Internet Explorer




2) Click **Yes** to proceed (Figure 60).

Figure 60: Security Alert




3) Enter your username and password and click **Log In** (Figure 61).

Figure 61: Example of iLO Login



4) The iLO summary page is displayed (Figure 62).

Figure 62: iLO Summary Page



System Status tab (Figure 62)

Brief explanation of iLO menu items:

- Status Summary: Basic iLO configuration
- iLO Status: Indicates current condition of iLO
- Server Status: Server configuration and status
- iLO Event Log: Events related to iLO
- Integrated Management Log: Log showing server events and error conditions
- Server and iLO Diagnostics: Results of automatic diagnostic tests

Remote Console tab


Options in this tab are unavailable at this time.

Virtual Devices tab (Figure 63)

Options in this tab allow you to accomplish tasks remotely that would normally require you to be at the server console.

- Virtual Power: Turn the server on or off
- Power Regulator: Adjust power settings
- Virtual Media: Connect to a drive on a remote machine
- Virtual Indicator: Control Server Unit ID light
- Virtual Serial Port: Virtual serial port status
- Virtual Serial Port Configuration: Virtual Serial Port Configuration


Figure 63: Virtual Devices Tab



Administration tab

The **User Administration** item is used to configure iLO users (Figure 64). The other options are not being used at this time.

Figure 64: Administration Tab




System Shut Down and Restart Instructions

The system may need to be shut down occasionally for maintenance. Because of the clustered nature of VBECS, the system has to be shut down in a specific order. Shutting down the system requires that a user be physically present at the VBECS system.

To shut down the system


1) Log into either of the servers. Click **Start, Shut Down**. Enter a comment in the **Comment** field and click **OK** (Figure 65).

Figure 65: Example of Shut Down Window




2) After the first server is completely shut down, log into the other server and click **Start, Shut Down**. Enter a comment and click **OK** (Figure 66).

Figure 66: Example of Shut Down Window



3) After both servers are completely shut down, the shared storage may be shut down by pressing the power button in the lower right corner (Figure 67).

Figure 67: Shared Disks



To start the system

- 1) Press the power button on the share storage. Wait until both controllers display a message of **Startup Complete** before continuing.
- 2) Start up one of the servers and allow it to come to the log on screen before continuing. This one will become the active node.
- 3) Start up the other server.

This page intentionally left blank.

Maintenance Operations

These maintenance operations are performed, using the VBECS Administrator software, during the initial installation of VBECS and during post-installation maintenance activities.

When VBECS Administrator is used for the first time, Configure Interfaces is the only option available. Completion of Configure Interfaces enables Configure Divisions. Completion of Configure Divisions enables Configure Users.

Configured options will be available at startup to perform maintenance operations.



Do not change the system! The U.S. Food and Drug Administration classifies this software as a medical device. Unauthorized modifications will render this device an adulterated medical device under Section 501 of the Medical Device Amendments to the Federal Food, Drug, and Cosmetic Act. Acquiring and implementing this software through the Freedom of Information Act require the implementer to assume total responsibility for the software and become a registered manufacturer of a medical device, subject to FDA regulations. Adding to or updating VBECS software without permission is prohibited.

- VistALink is installed and running on the associated VistA system.
- The user is defined in VistA, and has a DUZ and Access and Verify Codes necessary to establish a VistA connection.
- The user has a valid Windows account and is defined as a member of the Active Directory (AD) domain group (see Add and Maintain Users in Active Directory).
- The user is defined as a member of the Windows Administrator group on the Active Directory domain group.
- The VBECS database is installed and operational.
- The VBECS BUNDLE 1.0 KIDS build is installed and configured in VistA.
- The VistA data conversion patch LR*5.2*335 is installed in VistA.
- The VistA data conversion is complete.

Outcome

- Parameters necessary to establish the connection to VistA through VistALink are available to the main VBECS application, as defined in the Configure Interfaces option.
- VBECS-VistA HL7 interface parameters are defined in the Configure Interfaces option.
- One or more divisions are defined for use in VBECS in the Configure Divisions option.
- One or more divisions are activated as local facilities in VBECS in the Configure Divisions option.
- The System Administrator has VBECS login¹ access to all active divisions.
- VBECS users are defined and able to use VBECS in the Configure Users option.

¹ There is a slight difference in terminology between VistA and VBECS: VistA uses “log on” and “logon,” and VBECS uses “log in” and “login.” Therefore, both terms are used throughout this manual. “Log in” and “login” are used generically when referring to both systems at one time.

Limitations and Restrictions



When the division changes from full service to transfusion only or from transfusion only to full service, information must be in a final state.

- The VBECS Administrator performing the initial installation and setup must have the XOBV VISTALINK TESTER option defined as a secondary option in VistA.

Additional Information

- Refer to the completed Appendix: Configuration Worksheet in *VBECS Application Interfacing Support Software Installation and User Configuration Guide* for required information when performing maintenance operations.

User Roles with Access to This Application

VBECS Administrator

Log into VBECS Administrator

User Action	VBECS Administrator
<p>1. To log into VBECS Administrator, double click  (the Remote Desktop Connection icon).</p> <p>Enter your password.</p>	<p>Displays the user and server names.</p>
<p>2. Double click the VBECS Administrator icon.</p>	<p>Opens VBECS Administrator.</p> <p>NOTES _____</p> <p>When the user logs into VBECS Administrator for the first time to set VistALink parameters, the system does not display the VistA Logon – Authorization screen. Continue at Step 6.</p>
<p>3. Continue to the VistA logon screen (Figure 68).</p>	<p>Opens the VistA Logon – Authorization screen. The user may log onto VistA or continue and log on as needed.</p> <p>NOTES _____</p> <p>The VistA logon screen is displayed only after initial setup of VistALink parameters.</p>
<p>4. Log onto VistA when VBECS Administrator starts up or at the invocation of any option that uses VistALink when VistALink is not connected.</p>	<p>Allows a user to log on by entering VistA Access and Verify Codes, separated by a semicolon (;), in the Access Code data entry field. When a user accesses an option that requires a VistALink connection and the connection becomes unavailable, allows the user to restore the connection. When a reconnection attempt is successful, VBECS closes the connection status window and returns to the desktop. The VistALink Connected icon in the status bar indicates a successful connection. When a reconnection attempt is unsuccessful, attempts to reconnect to VistALink until the user cancels.</p> <p>NOTES _____</p> <p>When a user logs into VBECS Administrator, the connection to</p>

User Action	VBECS Administrator
	VistA is established through VistALink. When the VistALink connection is not restorable, VBECS Administrator displays a message that the requested use cannot be executed because VistALink is unavailable.
5. Enter the VistA Access and Verify Codes.	Verifies that user credentials for the VBECS Administrator and VistA Access and Verify Codes belong to the same user.
6. Continue working in VBECS Administrator (Figure 69).	Displays the main menu.

Figure 68: Example of VistA Logon





Figure 69: Example of VBECS Administrator




Configure Interfaces

The System Administrator sets parameters for the connection to VistA to enable retrieval of VistA data and to configure HL7 interfaces between VBECS and VistA.

User Action	VBECS Administrator
1. To configure VBECS VistALink and HL7 interface parameters, click File on the main menu of the VBECS Administrator software.	Displays the menu options used to configure VBECS.
2. Click Configure Interfaces (Figure 70).	Displays the VBECS Configure Interfaces dialog for data entry.

Figure 70: Configure Interfaces




Configure VistALink Parameters

User Action	VBECS Administrator
1. To configure VistALink parameters, select VistALink from the Select Interface list box (Figure 71).	Displays the Configure VistALink group and allows data entry of the IP address (or domain name) and port number of the VistA system VistALink listener. Allows the user to test the VistALink connection parameters. NOTES _____ The user may modify the IP address (or domain name) and port number, as required.
2. Enter a valid IP address (or domain name) and port number of the VistA system VistALink listener in the M Server group box fields.	Validates that the IP address is in the standard four-octet notation (e.g., 127.0.0.1) or that the Domain field was filled in. Validates that the port number is a whole number from 1024 to 65535.

User Action	VBECS Administrator
	<p>NOTES</p> <p>The IP Address field represents the VistALink IP address to which VBECS will direct messages. Refer to the Hardware Information section of Appendix B, row 6 for test, and row 7 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The Port Number field represents the VistALink port number to which VBECS will direct messages. Refer to the Hardware Information section of Appendix B, row 8 for test, and row 9 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p>
<p>3. Click Test Connection.</p> <p> Capture a screen shot.</p>	<p>NOTES</p> <p>The Test Connection button is enabled only when valid entries exist in the IP Address (or Domain) and Port Number fields.</p> <p>If connection to the VistA system is successful, the VistA Logon – Authorization dialog is displayed and the user is required to enter valid Access and Verify Codes.</p> <p>If connection to the VistA system is unsuccessful, hover over the red square and a detailed error message will display.</p>
4. Click Save to save changes.	Displays a confirmation dialog.
5. Click Yes to commit changes to the database.	

Figure 71: Configure Interfaces: VistALink



Configure CPRS HL7 Interface Parameters

User Action	VBECS Administrator
<p>1. To configure CPRS HL7 Interface Parameters, select CPRS from the Select Interface list box in the VBECS – Configure Interfaces dialog (Figure 72).</p>	<p>Displays the Configure Interface group and allows data entry of HL7 interface-related parameters.</p>
<p>2. To configure Interfaced Application group parameters, enter a valid IP address, port number, and facility ID in the related data fields.</p>	<p>Validates that the IP address is in the standard four-octet notation (e.g., 127.0.0.1) or that the Domain field was filled in. Validates that the port number is a whole number from 1024 to 65535.</p> <p>NOTES _____</p> <p>The IP Address field represents the VistA CPRS IP address to which VBECS will direct messages. The Domain name field represents the fully qualified domain name to which VBECS will direct messages. Refer to the Hardware Information section of Appendix B, row 6 for test, and row 7 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The Port Number field represents the VistA CPRS port number to which VBECS will direct messages. Refer to the Hardware Information section of Appendix B, row 10 for test, and row 11 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The Facility ID is used in the MSH segment of the HL7 interface to help identify the system. This free-text field is usually set to the primary site's station number. Messaging to VBECS will fail if this Facility ID is not supplied.</p>
<p>3. To configure VBECS Application group parameters, enter a valid IP address, port number, and facility ID in the related data fields.</p>	<p>Validates that the IP address is in the standard four-octet notation (e.g., 127.0.0.1). Validates that the port number is a whole number from 1024 to 65535.</p> <p>NOTES _____</p> <p>The IP Address field represents the VBECS cluster server IP address to which CPRS will direct messages. Refer to the Hardware Information section of Appendix: Configuration Worksheet, row 1 in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The Port Number field represents the VBECS cluster server port number to which CPRS will direct messages. Refer to the Hardware Information section of Appendix B, row 4 for test, and row 5 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The VBECS Facility ID must be different from the VistA. The Facility ID is used in the MSH segment of the HL7 interface to help identify the system. This is a free-text field set to the primary site's station number. Messaging to VBECS will fail if this Facility ID is not supplied.</p>

User Action	VBECS Administrator
	<p>The data entered in this group is used by the VBECS CPRS HL7 Listener Service when using a single listener interface. This service is installed as disabled and the VBECS HL7 Multi Listener is enabled. In this configuration, the Port Number field must be set to a port that is not currently used by any other services on the Cluster Server. Refer to the Hardware Information section of Appendix B: Configuration Worksheet, rows 4 and 5 in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p>
<p>(This step is optional.)</p> <p>4. To configure Message Options group parameters, enter an ACK timeout period and a number of retransmission attempts in the related data fields.</p>	<p>Validates that the ACK timeout period is a whole number from 1 to 999 (seconds) (default: 10). Validates that the number of retransmission attempts for failed messages is a whole number from 1 to 99 (default: 5).</p>
<p>(This step is optional.)</p> <p>5. To configure Purge Criteria group parameters, enter the number of days after which completed messages and messages in error are to be purged from the database in the related data fields.</p>	<p>Validates that purge periods are whole numbers from 1 to 30 (days) (default: 7).</p>
<p>6. To configure the Interface Failure Alert Recipient group parameter, enter a valid email address in the related data field.</p>	<p>Validates that the interface administrator's email address is entered and conforms to Internet message format RFC 2822.</p> <p>NOTES _____</p> <p>VBECS Windows Services uses this email address to notify local IRM support or the Blood Bank ADPAC when interface errors occur.</p>
<p>7. To configure the Logging Configuration group parameter, click or clear the Log Events and HL7 Messages to Event Log check box.</p> <p> Capture a screen shot.</p>	<p>NOTES _____</p> <p>This check box indicates whether to record incoming and outgoing HL7 messages in the Application Event Log on the VBECS Cluster Server. (This is the only way to view VBECS HL7 messages on the VBECS server.)</p>
<p>8. Click Save and Yes to confirm the save.</p>	
<p>9. To close the VBECS – Configure Interfaces dialog, click  in the upper right corner.</p>	<p>Validates that the data was saved.</p>

Figure 72: Example of Configure Interfaces: CPRS

The screenshot shows a software window titled "VBECS - Configure Interfaces". On the left is a "Select Interface" list with "CPRS" selected. The main area is titled "Configure Interface" and contains several sections:

- Interfaced Application:** "Connection Method*" is set to "IP Address" with a value of "10.2.2.21". "Port Number*" is "2222". "Facility ID" is empty.
- VBECs Application:** "IP Address*" is "10.1.1.5". "Port Number*" is "5555". "Facility ID" is "VBECs".
- Message Options:** "ACK Timeout*" is "10" secs. "Re-Transmit Attempts*" is "5".
- Purge Criteria:** "Completed Messages*" is "7" days. "Messages in Error*" is "7" days.
- Interface Failure Alert Recipient:** "E-mail Address*" is "foo@foo.com".
- Logging Configuration:** "Log Events and HL7 Messages to Event Log" is unchecked.


At the bottom right are "Clear" and "Save" buttons. A legend at the bottom left indicates "* Required Field".

Configure Patient Update HL7 Interface Parameters

User Action	VBECS Administrator
1. To configure Patient Update HL7 Interface Parameters, select PatientUpdate from the Select Interface list box in the VBECS – Configure Interfaces dialog (Figure 73).	Displays the Configure Interface group and allows data entry of HL7 interface-related parameters.
2. To configure Interfaced Application group parameters, enter a facility ID in the related data fields.	<p>NOTES _____</p> <p>The IP Address and Port Number fields are disabled: no outbound messages are sent to VistA for this interface.</p> <p>The facility ID is used in the MSH segment of the HL7 interface to help identify the system. This is a free-text field set to the primary site's station number. Messaging to VBECS will fail if this Facility ID is not supplied.</p>
3. To configure VBECS Application group parameters, enter a valid IP address, port number, and facility ID in the related data fields.	<p>Validates that the IP address is in the standard four-octet notation (e.g., 127.0.0.1).</p> <p>Validates that the port number is a whole number from 1024 to 65535.</p> <p>NOTES _____</p> <p>The IP Address field represents the VBECS cluster server IP address to which VistA will direct messages. Refer to the Hardware Information section of Appendix: Configuration Worksheet, row 1 in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The Port Number field represents the VBECS cluster server port number to which VistA will direct messages. Refer to the Hardware Information section of Appendix B, row 4 for test, and row 5 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p>
(This step is optional.)	Validates that the ACK timeout period is a whole number from 1 to 999 (seconds) (default: 10).
4. To configure Message Options group parameters, enter an ACK Timeout period and number of retransmission attempts in the related data fields.	Validates that the number of retransmission attempts for failed messages is a whole number from 1 to 99 (default: 5).
(This step is optional.)	Validates that the purge periods are whole numbers from 1 to 30 (days) (default: 7).
5. To configure Purge Criteria group parameters, enter the number of days after which completed messages and messages in error are to be purged from the database in the related data fields.	
6. To configure the Interface Failure Alert Recipient group parameter, enter a valid email address in the related data field.	<p>Validates that the interface administrator's email address is entered and conforms to Internet message format RFC 2822.</p> <p>NOTES _____</p> <p>VBECS Windows Services uses this email address to notify local</p>

User Action	VBECS Administrator
	IRM support or the Blood Bank ADPAC when interface errors occur.
<p>7. To configure the Logging Configuration group parameter, click or clear the Log Events and HL7 Messages to Event Log check box.</p> <p> Capture a screen shot.</p>	<p>NOTES _____</p> <p>This check box indicates whether to record incoming and outgoing HL7 messages in the Application Event Log on the VBECS Cluster Server. (This is the only way to view VBECS HL7 messages on the VBECS server.)</p>
8. Click Save and Yes to confirm the save.	
9. To close the VBECS – Configure Interfaces dialog, click  in the upper right corner.	Validates that the data was previously saved.

Figure 73: Example of Configure Interfaces: PatientUpdate




Configure Patient Merge HL7 Interface Parameters

User Action	VBECS Administrator
1. To configure Patient Merge HL7 Interface Parameters, select PatientMerge from the Select Interface list box in the VBECS – Configure Interfaces dialog (Figure 74).	Displays the Configure Interfaces group and allows data entry of HL7 interface-related parameters.
2. To configure Interfaced Application group parameters, enter a facility ID in the related data field.	<p>NOTES _____</p> <p>The IP Address and Port Number fields are disabled: no outbound messages are sent to VistA for this interface.</p> <p>The facility ID is used in the MSH segment of the HL7 interface to help identify the system. This is a free-text field set to the primary site's station number. Messaging to VBECS will fail if this Facility ID is not supplied.</p>
3. To configure VBECS Application group parameters, enter a valid IP address, port number, and facility ID in the related data fields.	<p>Validates that the IP address is in the standard four-octet notation (e.g., 127.0.0.1).</p> <p>Validates that the port number is a whole number from 1024 to 65535.</p> <p>NOTES _____</p> <p>The IP Address field represents the VBECS cluster server IP address to which VistA will direct messages. Refer to the Hardware Information section of Appendix B, row 1: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p> <p>The Port Number field represents the VBECS cluster server port number to which VistA will direct messages. Refer to the Hardware Information section of Appendix B, row 4 for test, and row 5 for production: Configuration Worksheet in <i>VBECS Application Interfacing Support Software Installation and User Configuration Guide</i>.</p>
(This step is optional.)	Validates that the ACK Timeout period is a whole number from 1 to 999 (seconds) (default: 10).
4. To configure Message Options group parameters, enter an ACK Timeout period and number of retransmission attempts in the related data fields.	Validates that the number of retransmission attempts for failed messages is a whole number from 1 to 99 (default: 5).
(This step is optional.)	Validates that the purge periods are whole numbers from 1 to 30 (days) (default: 7).
5. To configure Purge Criteria group parameters, enter the number of days after which completed messages and messages in error are to be purged from the database in the related data fields.	
6. To configure the Interface Failure Alert Recipient group parameter, enter a valid email address in the related data field.	<p>Validates that the interface administrator's email address is entered and conforms to Internet message format RFC 2822.</p> <p>NOTES _____</p>

User Action	VBECS Administrator
	VBECS Windows Services uses this email address to notify local IRM support or the Blood Bank ADPAC when interface errors occur.
<p>7. To configure the Logging Configuration group parameter, click or clear the Log Events and HL7 Messages to Event Log check box.</p> <p> Capture a screen shot.</p>	<p>NOTES _____</p> <p>This check box indicates whether to record incoming and outgoing HL7 messages in the Application Event Log on the VBECS Cluster Server. (This is the only way to view VBECS HL7 messages on the VBECS server.)</p>
8. Click Save and Yes to confirm the save.	
9. To close the VBECS – Configure Interfaces dialog, click  in the upper right corner.	Validates that the data was previously saved.

Figure 74: Example of Configure Interfaces: PatientMerge



Configure Divisions

The System Administrator configures VBECS as a single division or as multidivisional.

Assumptions

- The VistA data conversion is complete.
- VBECS-VistA connection parameters are set.
- VistALink is installed and running on the associated VistA system.
- The user is defined in VistA, and has a DUZ and Access and Verify Codes necessary to establish a VistA connection.
- The user has a valid Windows account and is defined as a member of the Active Directory domain group (see Add and Maintain Users in Active Directory).
- The IP address of the label printer is known.
- The name of the division report printer is known (if multi-divisional).
- The VBECS database is installed and operational.

Outcome

- One or more divisions are defined in VBECS.
- One or more divisions are activated as local facilities in VBECS.
- The System Administrator has VBECS login² access to all active divisions.

Limitations and Restrictions

- All units in a division must be in a final status to allow the division to change from full service to transfusion only or from transfusion only to full service.

Additional Information

- A VBECS Administrator/Supervisor may further configure:
 - VBECS users in Update User Roles.
 - VBECS division parameters in Configure Division, Product Modifications, and Configure Testing.
- The user must log onto VistA using Access and Verify Codes.

User Roles with Access to This Option

System Administrator

Add and Maintain Divisions

The user defines and maintains division attributes.

 *Changes made in the VBECS Administrator option mapping orders to another VBECS division do not affect delivered orders. Orders delivered to a VBECS division must be completed, rejected, or canceled in that division. Resubmit orders after mapping is completed to send an order to another VBECS division.*

² There is a slight difference in terminology between VistA and VBECS: VistA uses “log on” and “logon,” and VBECS uses “log in” and “login.” Therefore, both terms are used throughout this manual. “Log in” and “login” are used generically when referring to both systems at one time.

User Action	VBECS Administrator
1. To add and maintain divisions in VBECS, click File on the main menu of the VBECS Administrator software.	<ul style="list-style-type: none"> • Displays the menu options used to configure VBECS.
2. Select Configure Divisions (Figure 75).	<ul style="list-style-type: none"> • Displays the Configure Division dialog and allows entry of division parameters.
3. To edit a defined division, click the Division Identification tab (Figure 76). Select a division code or name from the drop-down menu or, to configure a new division, click the ellipsis button. Select a division from the list (Figure 77).	<p>NOTES _____</p> <p>The user may not edit the division code or name.</p> <p>A division may be full service (default) or transfusion only. When a unit not in a final status exists, a user may not change the type of transfusion service.</p> <p>When a division is transfusion only, VBECS disables electronic crossmatch.</p> <p>When a division changes from full service to transfusion only, units already in inventory are not restricted to patients and must be returned to the blood center.</p> <p>When a division changes from transfusion only to full service, inventory units are restricted to patients without ABO/Rh confirmation. The facility must decide how to handle this existing inventory.</p> <p>VBECS prevents the user from changing a division from full service to transfusion only or from transfusion only to full service when there are open or partially completed worksheets or processes in the division.</p> <p>The Division Name and Division Code are identified in the VistA INSTITUTION file (#4). The Division Name stored in VBECS is the INSTITUTION file NAME field (#.01); the Division Code stored in VBECS is the STATION NUMBER field (#99). When either value change in VistA, rerun these steps to update the VBECS database with the current values from VistA.</p>

User Action	VBECS Administrator
<p>4. To receive orders from VistA Institutions to the selected Division, check the Map orders from VistA institutions check box. Click the Active checkbox for each institution that applies.</p>	<p>NOTES</p> <p>Changes made to institution mappings require a restart of the VBECS HL7 Multi Listener service. For more information, see Table 8 in the VBECS Windows Services section.</p> <p>One or more VistA institutions from the list of valid institutions retrieved from VistA may be associated with the selected VBECS division from the list of valid institutions retrieved from VistA.</p> <p>A VistA institution may be associated with only one VBECS division.</p> <p>A VistA institution defined as a VBECS division is not eligible for selection as an associated institution to a different VBECS division.</p> <p>To associate additional institutions, enable an optional VistALink query to retrieve a list of all institutions associated with the VistA site that are currently defined within the VistA database but not in the selected VBECS division. VBECS displays the list to the user for selection.</p>
<p>5. Select the FDA Registered Facility associated with the division or, to search for the facility by name or FDA Registration Number, click the ellipsis button (Figure 76).</p>	<ul style="list-style-type: none"> • Allows the user to associate a division with a facility from the National Facility Table. <p>NOTES</p> <p>The user must associate a division with a facility from the National Facility Table. If there is no matching facility, VBECS Administrator asks the user to contact the VA Service Desk.</p> <p>When this occurs, wait for customer support to respond or, to continue establishing a division, select and configure any facility from the National Facility Table. When the configuration is complete, use the Local Facilities option in VBECS to define the local facility that matches the information missing from the National Facility Table.</p> <p>Return to Configure Divisions to re-associate your division with the newly entered local facility.</p> <p>When a division is configured, VBECS displays, "I certify that the blood products listed were properly maintained, in accordance with the Code of Federal Regulations, while in storage at this institution. Components were inspected when packed for shipment and found to be satisfactory in color and appearance."</p>
<p>6. Select the VistA Lab Blood Bank Accession Area associated with the selected division from the drop-down menu (Figure 76).</p>	<p>NOTES</p> <p>The Lab package uses the Accession Area to track blood bank-related workload for the division.</p>
<p>7. Enter the desired number of minutes in the Lock Inactivity Timeout field.</p>	<ul style="list-style-type: none"> • Allows the user to set the lock inactivity timeout period (5 to 15 minutes) (default: 5 minutes). <p>NOTES</p> <p>The lock inactivity timeout period specifies how long a user can be idle and in control of data being edited. VBECS warns the</p>

User Action	VBECS Administrator
	<p>user 60 seconds before the lock inactivity period expires that he will lose priority for the data. When he responds within 60 seconds, VBECS clears the warning and resets the lock activity timer. Otherwise, VBECS informs him that his lock was released and he must reenter his changes.</p> <p>VBECS uses optimistic and pessimistic locking to prevent data corruption. If a user attempts to edit data locked by another user, VBECS alerts him that the record is in use and prevents access (pessimistic locking).</p> <p>If more than one user attempts to change data simultaneously, VBECS accepts only the first update and warns the other users that the record changed (optimistic locking, which is non-configurable and a fail-safe to pessimistic locking).</p>
<p>8. To activate or inactivate the division, click or clear the Active VBECS Division? check box (Figure 76).</p> <p> Capture a screen shot.</p>	<ul style="list-style-type: none"> When the user saves a previously active division as inactive, inactivates user roles for that division. <p>NOTES _____</p> <p>The system will not allow the user to activate a division that has orders mapped to another VBECS division. VBECS displays, "Unable to activate. The VBECS division currently has orders mapped to another VBECS division."</p> <p>The system will not allow the user to inactivate a division that has orders mapped to it. VBECS displays, "Unable to inactivate. This VBECS division currently has orders mapped to it. Release this mapping prior to inactivation,"</p>
<p>9. Click the Service Type tab. Click the Full-Service Facility or Transfusion-Only Facility radio button (Figure 79).</p> <p> Capture a screen shot.</p>	<ul style="list-style-type: none"> Allows the user to identify the facility as full service or transfusion only. <p>NOTES _____</p> <p> When the division changes from full service to transfusion only or from transfusion only to full service, information must be in a final state. VBECS does not check for pending orders or active units in inventory, so there is a risk of corrupting information. There is a risk of having unconfirmed units available for transfusion if any are issued.</p>
<p>10. Click the Printers tab.</p> <p>Clear or click the Division Uses Label Printer check box.</p> <p>Edit the COM port number and/or the TCP port number.</p> <p>Enter the IP address (Figure 80).</p> <p> Capture a screen shot.</p>	<ul style="list-style-type: none"> Allows the user to enter the COM and TCP port numbers and the IP address for the label printer. Allows the user to select the default printer for the division when more than one printer is installed on the system. <p>NOTES _____</p> <p>Standard values for COM and TCP ports:</p> <ul style="list-style-type: none"> COM = 2 TCP = 9100
<p>11. Click the Time Zone tab.</p> <p>Select a time zone.</p> <p>In the Daylight Savings field, select</p>	<ul style="list-style-type: none"> Allows the user to set the time zone and daylight saving parameters.

User Action	VBECS Administrator
<p>US Standard DST, Do not observe DST, or Custom DST.</p> <p>Enter start and end dates for custom DST (Figure 81).</p> <p> Capture a screen shot.</p> <p>Click Save.</p>	
<p>12. Click Save and OK to commit the changes or add the new division to the VBECS database.</p>	<ul style="list-style-type: none"> Commits changes and additions to the database. <p>NOTES _____</p> <p>Multidivisional sites must repeat Steps 3–11 for each division.</p> <p>The VBECS Administrator/Supervisor who configured the divisions must add himself as a user to all divisions to enable the functionality of canned comments in the VBECS system.</p>
<p>13. To close the VBECS – Configure Divisions dialog, click  in the upper right corner.</p>	

Figure 75: Configure Divisions




Figure 76: Example of Configure Division: Division Identification




Figure 77: Example of Select VistA Divisions



Figure 78: Example of Facility Search




Figure 79: Example of Configure Division: Service Type




Figure 80: Example of Configure Division: Label Printing





Figure 81: Example of Configure Division: Time Zone



Configure System Administrators

Each non-data center site must assign an onsite system administrator to perform regular maintenance tasks such as applying a Windows update and troubleshooting. If your servers reside at a data center, personnel at that location will be administering the servers and you may skip this section.

Assumptions

- The user has a valid Windows login and was given permission to manage the Active Directory administrator group (set up at installation).
- Users to be configured have a valid Windows account.

Outcome

- Administrators are defined and able to administer the VBECS servers from the client.

Limitations and Restrictions



Each VBECS user must have a unique Windows login ID. If a Windows login ID becomes inactive and is eligible for re-use in Active Directory, do not re-use it for VBECS: it may result in corrupted data in VBECS.

Additional Information

- None

Add or Remove System Administrators

The user adds and inactivates VBECS users.

User Action	Active Directory Users and Computers
1. Install Active Directory tools (on the Administrator's computer only) from the Windows Server 2003 Enterprise Edition installation CD or as a free download from Microsoft.	
2. Open the Control Panel. Double click Administrative Tools . Double click Active Directory Users and Computers (Figure 82).	<ul style="list-style-type: none"> • Allows the user to view and add users in Active Directory for VBECS.
3. Navigate to the Organizational Unit (OU) in which your VBECS local groups reside. Double click the name of the user group (on the right) to which you wish to add the user (Figure 83).	<ul style="list-style-type: none"> • Displays administrator group in the right panel. • Displays the properties window. <p>NOTES _____</p> <ul style="list-style-type: none"> • Add a user to the administrator group to allow administrative access to the server through Remote Desktop Connection.
4. Click the Members tab (Figure 84). Click Add to add a user.	<p>NOTES _____</p> <p>If the Add button is disabled, you do not have access to this</p>

User Action	Active Directory Users and Computers
To remove a user, select the user name and click Remove .	group. File a Remedy ticket to gain access.
5. If the From this location field does not display the location of the user to be added, click Locations and enter the correct domain (Figure 85).	<ul style="list-style-type: none"> Allows the user to enter the domain.
6. In the Enter the object names to select field, enter the Windows login ID for the user to be added. Click OK .	<p>NOTES _____</p> <p>Click Check Names to verify that the login ID is valid.</p>
7. Click OK .	<ul style="list-style-type: none"> Closes the Properties window.
8. Exit.	

Figure 82: Example of Active Directory User and Computers Console




Figure 83: Example of Administrator User Group




Figure 84: Example of Group Properties




Figure 85: Example of Select Users



Configure Users

The System Administrator matches VistA users to VBECS users and sets user security levels. If this is a data center site, use the form (Appendix D: Active Directory Request Form) to submit Active Directory modifications and skip the “Add and Maintain Users in Active Directory” section (proceed to the “Configure VBECS Users” section after the data center has completed your request).

Assumptions

- The VistA data conversion is complete.
- VBECS-VistA connection parameters are set.
- VistALink is installed and running on the associated VistA system.
- VBECS application configuration files have the correct values for Domain and user group fields.
- At least one division in VBECS is configured.
- The user is defined in VistA, and has a DUZ and Access and Verify Codes necessary to establish a VistA connection.
- All users of the Blood Bank medical device software are assigned the VBECS VISTALINK CONTEXT option as a secondary option. VistALink uses the VBECS VISTALINK CONTEXT option to provide user context sign-on security to VistA.
- The user has a valid Windows login and is defined as a member of the Active Directory domain group.
- The System Administrator created Active Directory local groups, as directed in Appendix: Blood Bank Configuration Checklist, Create Local Groups, in *VistA Blood Establishment Computer Software (VBECS) Installation Guide*.
- The VBECS database is installed and operational.

Outcome

- VBECS users are defined and able to use VBECS.

Limitations and Restrictions



Each VBECS user must have a unique Windows login ID. If a Windows login ID becomes inactive and is eligible for re-use in Active Directory, do not re-use it for VBECS: it may result in corrupted data in VBECS.

A user must not change their Windows login ID after being configured in VBECS. If the user's name changes, the name fields in Active Directory can be modified without changing the login ID.

Additional Information

- A VBECS Administrator/Supervisor may further configure VBECS users in Update User Roles.
- The user must log onto VistA using Access and Verify Codes.

User Roles with Access to This Option

System Administrator

Add and Maintain Users in Active Directory

The user adds and inactivates VBECS users.

User Action	Active Directory Users and Computers
1. Install Active Directory tools (on the Administrator's computer only) from the Windows Server 2003 Enterprise Edition installation CD or as a free download from Microsoft.	
2. Open the Control Panel. Double click Administrative Tools . Double click Active Directory Users and Computers (Figure 86).	<ul style="list-style-type: none"> • Allows the user to view and add users in Active Directory for VBECS.
3. Navigate to the OU in which your VBECS local groups reside. Double click the name of the user group (on the right) to which you wish to add the user (Figure 87).	<ul style="list-style-type: none"> • Displays two user groups in the right panel, one for VBECS Administrator and one for VBECS. • Displays the properties window. <p>NOTES</p> <p>The VBECS local groups (<i>VnnxxxVbecsUsers</i> and <i>VnnxxxVbecsAdministrators</i>, where <i>nn</i> is your VISN number and <i>xxx</i> is your site identifier) were created in Appendix : Blood Bank Configuration Checklist, Create Local Groups, in <i>VistA Blood Establishment Computer Software (VBECS) Installation Guide</i>.</p> <p>The VBECS Administrator/Supervisor who configured the divisions must add himself as a user to all divisions to enable the functionality of canned comments in the VBECS system. He may inactivate himself later without affecting canned comments.</p> <ul style="list-style-type: none"> • Add a user to either group to allow access to the server through Remote Desktop Connection and to VBECS Administrator or

User Action	Active Directory Users and Computers
	VBECS (depending on the group).
<p>4. Click the Members tab (Figure 88).</p> <p>Click Add to add a user.</p> <p>To remove a user, select the user name and click Remove.</p>	<p>NOTES _____</p> <p>If the Add button is disabled, you do not have access to this group. File a Remedy ticket to gain access.</p>
<p>5. If the From this location field does not display the location of the user to be added, click Locations and enter the correct domain (Figure 89).</p>	<ul style="list-style-type: none"> Allows the user to enter the domain.
<p>6. In the Enter the object names to select field, enter the Windows login ID for the user to be added.</p> <p>Click OK.</p>	<p>NOTES _____</p> <p>Click Check Names to verify that the login ID is valid.</p>
<p>7. Click OK.</p>	<ul style="list-style-type: none"> Closes the Properties window.
<p>8. Exit.</p>	

Figure 86: Example of Active Directory Users and Computers




Figure 87: Example of Active Directory Users




Figure 88: Example of Group Properties





Figure 89: Example of Select Users



Configure VBECS Users

The Active Directory setup must be completed prior to configuring users in VBECS.

User Action	VBECS Administrator
1. To add and maintain users in VBECS, click File on the main menu of the VBECS Administrator software.	<ul style="list-style-type: none"> Displays the menu options used to configure VBECS.
2. Select Configure Users (Figure 90).	<ul style="list-style-type: none"> Allows the user to enter or edit user information.
3. To edit an existing user, select a user ID from the drop-down list (Figure 91) or, to search for a new user ID to add to VBECS, click the ellipsis button to the right of the drop-down list (Figure 92). Enter user parameters. For each user, VBECS stores: <ul style="list-style-type: none"> VistA DUZ Windows Login ID Windows Username Email Address (optional) User Initials Active Status Division Code User Role Division Active Status 	<ul style="list-style-type: none"> Displays the Windows user ID and name. <p>NOTES</p> <p>VistALink lists active VistA Blood Bank users. VistA Blood Bank users are identified by the LRBLOODBANK and LRBLSUPER security keys.</p> <p>When VBECS finds users that are inactive in VistA, it asks whether the user wishes to inactivate them in VBECS. Yes inactivates the VBECS users. No allows the user to continue without inactivating the users (Figure 95).</p> <p>The user may not edit the VistA DUZ or user name, the Windows login ID or user name, or the division code or name.</p> <p>There is a one-to-one correspondence between Windows and VistA users. A VistA DUZ may be associated with only one Windows login ID and vice versa.</p> <p>The user may:</p> <ul style="list-style-type: none"> Activate or inactivate but not delete a defined user from VBECS. Rescind a defined user's access privileges at one or more divisions but not delete his record or ID from the database. <p>The user ID stored in VBECS is the user's Windows Logon ID. VBECS displays the data that a user enters in a session. The</p>

User Action	VBECS Administrator
	<p>user may edit and save the data. When a user cancels, VBECS warns that it will not save the data. VBECS closes the form and returns the user to the main menu screen that may include unrelated open windows.</p> <p>VBECS associates the technologist ID, date, time, and division with each process for retrieval by division.</p>
<p>4. To search for a VistA user, click the ellipsis button to the right of the VistA DUZ field (Figure 93).</p>	<ul style="list-style-type: none"> Allows the user to search for VistA Blood Bank users by name or DUZ. <p>NOTES _____</p> <p>The user may not edit the VistA DUZ or user name, the Windows login ID or user name, or the division code or name.</p>
<p>5. Enter the email address of the user in the E-mail field in the Additional Info group. VistA provides the initials, if available. If not, enter them.</p>	<ul style="list-style-type: none"> Allows the user to enter Additional Information about the user for identification. <p>NOTES _____</p> <p>User initials may be loaded from VistA. VBECS requires unique user initials for use as the technologist ID.</p>
<p>6. To select a VistA division to associate with the user, click the ellipsis button to the right of the Division Code drop-down menu (Figure 94).</p>	<ul style="list-style-type: none"> Allows the user to select a division to associate with the user <p>NOTES _____</p> <p>A single user may be associated with multiple divisions.</p>
<p>7. Select a user role from the User Role drop-down menu. Click or clear the Active Role? check box to activate or inactivate the role.</p>	<ul style="list-style-type: none"> Allows the user to assign security roles to the Blood Bank user. If a user was removed from the role of Administrator/Supervisor and was the only Administrator/Supervisor user left for a division, displays "You are trying to remove the last Administrator/Supervisor for your division, which would disallow system configuration in the future. You may not proceed." If all entered data is satisfactory, saves user details and access changes to the file and adds or updates the user information in the list view. <p>NOTES _____</p> <p>One role at a time may be assigned to a user at a division. A user may have only one active user role per division.</p> <p>VBECS allows the assignment of a security level to one or more users at a time. VBECS warns that there must be at least one level 6 VBECS Administrator/Supervisor in the division and does not allow the user to change the last Administrator/Supervisor.</p>
<p>8. Click Update and Save.</p>	<ul style="list-style-type: none"> Displays a confirmation dialog.
<p>9. Click Yes to commit changes to the database.</p>	<ul style="list-style-type: none"> Click Yes to commit changes to the database.
<p>10. To close the Edit Users dialog box, click  in the upper right corner.</p>	

Figure 90: Configure Users




Figure 91: Example of Edit User




Figure 92: Example of Windows Users




Figure 93: Example of VistA Users




Figure 94: Example of VistA Divisions




Figure 95: Example of Inactive Users



Transmit Workload Data

VBECS workload data is recorded in VBECS when records that qualify as Workload Events are saved in VBECS. This data is transmitted to the VistA Laboratory workload recording system for national and local workload reporting.

Assumptions

- Workload codes were assigned to VBECS processes using Workload Codes.
- Healthcare Common Procedure Coding System (HCPCS) codes were assigned to blood products using Blood Products.
- A record was saved or inactivated immediately preceding workload data collection.
- The connection to VistA is active.

Outcome

- Information was transmitted to VistA for inclusion in appropriate reports.

Limitations and Restrictions

- None

Additional Information

- Workload Event data must include information required for Decision Support System (DSS), Patient Care Encounter (PCE), and Billing Awareness. Once in VistA, existing VistA functionality will handle required reporting.
- The system accumulates and periodically transmits workload information to the VistA Lab workload recording process. The data is transmitted from VBECS to VistA by the VBECS Workload Capture Remote Procedure called by a nightly Lab background process.
- Workload multipliers for all Blood Bank activities in VistA File #64 must be set to one (1) to avoid excessive LMIP counts. This allows the workload multiplier set in VBECS to be correctly reflected on VistA reports.

User Roles with Access to This Option

All users

Transmit Workload Data

These steps are associated with the “Save” function within any class that performs a Workload Event such as recording a blood test result or interpretation for a unit or a patient, modifying a unit, and pooling units. VBECS must know which classes perform Workload Events and how to classify the work accomplished for reporting. When the database is updated, the VistA technologist ID of the updater, the division, and the date and time of the update are recorded. In some instances, a mechanism to capture Laboratory Management Index Program (LMIP) workload information exists. In addition, for certain events that involve patient processing, the patient location, treating specialty, service, etc., are captured to satisfy PCE or DSS reporting requirements.

These steps address the initial recording of these events.

User Action	VBECS
1. Click Save to save a record from an option.	<p>Creates a Workload Event for every process record saved. Recognizes the activity as a new Workload Event. Checks for required reporting properties based on the type of record being saved. Determines the proper workload codes and other related information to be included.</p> <p>NOTES _____</p> <p>One or more workload codes can be collected with each Workload Event saved. A workload code may be multiplied for certain Workload Events.</p>
2. Exit.	

Inactivate a Workload Event

VBECS updates VistA to inactivate the associated workload information (for a patient or a unit) so that PCE and Billing Awareness can be updated to reflect that the transaction is not valid.

User Action	VBECS
1. Inactivate a saved record.	<p>Recognizes the activity performed as an inactivation of an existing Workload Event record.</p> <p>NOTES _____</p> <p>See Appendix B: Workload Process Mapping to Application Option Table.</p>
2. Complete the update and choose to save.	<p>Prompts to confirm the save. Saves workload data.</p> <p>NOTES _____</p> <p>When a previously saved workload-generating event is invalidated (such as in Remove Final Status, Invalidate Test Results, or invalidating previously logged-in units through Edit Unit Information or Invalidate Shipment), VBECS must create and transmit the same Workload Event information to VistA as a negative number.</p>
3. Confirm the save.	<p>Saves workload data.</p> <p>NOTES _____</p> <p>When a saved Workload Event is associated with a patient, VBECS needs to link the Workload Event to the patient for future reports.</p>
4. The option ends when the record is saved.	

Notify VBECS Central Administrator

When maintenance operations are configured, the Implementation Manager notifies the VBECS Central Administrator to install ePolicy and MOM.

External Interfaces

VBECS uses VistALink Remote Procedure Calls (RPCs) and HL7 messaging with Microsoft Windows Services for data exchange using a client-server mode interfacing architecture. These services are cluster aware and continue to function in the event of a server failover.

Data exchange between the VBECS medical device software and other VistA applications is maintained by private Database Integration Agreements (DBIAs) with the VBECS Application Interfacing Support Software (VAISS) M software and HL7 messaging specifications with other VistA applications. The VAISS M software in the VistA environment is not classified as a medical device and is; therefore, exempt from the VBECS Blood Bank software FDA 510(k) submission. The purpose of this software is to provide data exchange with other VistA applications through a controlled environment.

When communication failures occur in the VistA environment between VBECS and other VistA applications, MailMan sends an email message to the G.VBECS INTERFACE ADMIN mail group. The message includes details of the error to assist with troubleshooting. Refer to Table 9 in the Troubleshooting section for a list of potential error messages and their solutions.

VBECS is not Clinical Context Management compliant. VBECS utilizes Remote Desktop Connection to connect to its dedicated server. If VBECS were to implement Clinical Context Management, the context would be with the VBECS server environment and require other software such as CPRS to be installed on the VBECS server. This is not compatible with the basic design of the encapsulated medical device.

Health Level Seven Interfaces

The VBECS Health Level Seven (HL7) software is a set of Microsoft .NET libraries written in C sharp (C#) that provide HL7 messaging support for VBECS.

The C# software is invoked by Microsoft Windows Services that run outside the VBECS application on the VBECS Cluster server to allow messaging transactions to occur without user intervention or the need for the VBECS application to be running. Some of the key common functionality provided by the software includes:

- Client-Server Transport Layer with HL7 Lower Layer Protocol support
- Message Queuing
- Message parsing and building libraries


Client-Server

The C# software provides a transport layer with HL7 Lower Layer Protocol support that uses a client-server architecture to allow bidirectional HL7 message exchange between VBECS and other VistA HL7 enabled applications. The software includes a common communications driver that allows VBECS to send and receive HL7 messages to and from multiple VistA applications. The software was designed to support multiple interfaces running concurrently without the operations of one interface interfering with another.

Each interface requires two separate roles of the client and server (Figure 96).

- Sending System = TCP Client (initiates connection to the Receiving System)
- Receiving System = TCP Server (listens for connections)

Figure 96: Client-Server Over TCP/IP Channels



Transport Layers and Lower Layer Protocols

A transport layer defines the physical connections between VBECS and other systems. Examples include TCP/IP networks and serially cabled connections.

The VBECS HL7 software supports multiple HL7 interfaces developed for VBECS and configured through VBECS Administrator by an authorized user. Some of the information, such as TCP/IP addresses and port numbers, are required by the transport layer and lower layer protocols to provide network connectivity and data exchange with an interfaced system.

An HL7 Lower Layer Protocol (LLP) defines how the systems communicate and exchange HL7 messages across a transport layer. While not defined within the HL7 standard itself, several LLPs are defined in *Health Level Seven Implementation Support Guide*.

LLPs provide the lower layer communication functionalities to exchange messages between systems, such as flow control and error recovery. “Lower layer” refers to a portion of the Open Systems Interconnect (OSI) model, which is divided into seven layers. The lower layers (1 through 4) include the physical connection between the systems and the communications protocol used. The HL7 standard itself defines the seventh and highest application layer.

The VBECS HL7 software supports only the Minimal Lower Layer Protocol (MLLP) over the VA TCP/IP transport layer. More information regarding the MLLP can be found in Section C.4: Minimal Lower Layer Protocol, Appendix: Lower Layer Protocols, of *Health Level Seven Implementation Support Guide*.

TCP Client (Sender)

The VBECS HL7 software allows VBECS to send outbound HL7 messages to a TCP/IP listener that supports the MLLP and receive an HL7 acknowledgement message over the same connection. The software provides the transport layer used to deliver the messages and receive the acknowledgement to the message.

To provide guaranteed message delivery of outbound messages from VBECS, all outbound messages will be created when certain events occur and are queued in the VBECS message log. A client monitor service polls the message log periodically to check for new outbound messages and sends them to the receiving system associated with the message type.

TCP Server (Listener)

All VBECS HL7 Listeners are implemented as Windows Services to provide minimal downtime with minimal user interaction. The default services are configured to start automatically on system reboot by default, but can be changed. HL7 interfaces operate using a single or multi listener Windows Service. The multi listener windows service is the default HL7 listener and can accept and process HL7 messages for all VBECS HL7 interfaces.

Computerized Patient Record System

Computerized Patient Record System (CPRS) is used to create requests for blood products and diagnostic tests performed in the blood bank with VBECS. An HL7 interface exists between CPRS and VBECS to transmit requests and provide updates regarding the requests to both sides of the interface. VBECS and CPRS exchange data using OMG-O19 General Order Messages and ORG-O20 Response to General Order Message (Acknowledgement) messages.

Orders in VBECS are directed to a VBECS division based on the division associated with the patient location (hospital location) selected in CPRS during the order entry process. If a patient order is associated with a hospital location for a division other than one defined in VBECS, the order will be returned to CPRS and canceled immediately. MailMan will send an email message to the ordering physician in VistA indicating that the order was canceled. The error text associated with the order will indicate that the division is not supported in VBECS. A new order must be created for a hospital location with a valid blood bank division. The CPRS interface supports HL7 version 2.4.

VistA Patient Updates

VBECS maintains a separate patient table for blood bank patients with a limited subset of patient-specific data, provided by the VistA system, for blood bank patient orders created through CPRS. VBECS must maintain updates on patient-specific data when changes are made in the VistA system. The patient-specific data that VBECS maintains includes the patient name, date of birth, date of death, gender, social security number, Integration Control Number (ICN), and the VistA internal entry number from the VistA Patient file. The Registration HL7 interface allows VBECS to receive ADT-A08 HL7 messages for all VistA patient data update events. The Patient Update interface supports HL7 version 2.3.

VistA Patient Merges

Occasionally, two entries in the VistA patient file are identified as duplicate records for the same patient and the two records must be merged into one. The duplicate records are validated through existing processes in VistA and are merged into a single record. When this occurs, VBECS must receive notification of the merge event and determine whether either of the two patient records exists in the VBECS Patient table. When matching records are identified, VBECS alerts the user. The user must update the patient record manually to match the VistA record. The MPI Patient Merge HL7 interface allows VBECS to receive ADT-A40 HL7 messages when two VistA patient records are merged into one. The Patient Merge interface supports HL7 version 2.4.

VistALink Remote Procedure Calls

Remote Procedure Calls (RPCs) provide a method of data exchange through VistALink for VBECS. The VBECS software provides data to or receives data from the VAISS located in the VistA M environment through RPCs. This data exchange is controlled through DBIAs between the blood bank medical device software and the VAISS VistA M software.

The VAISS software provides a set of M Application Programmer Interfaces (APIs) that call VBECS RPCs through the VBECS VistALink RPC XML Listener Windows Service and return blood bank data to other VistA applications. The VAISS software also provides a set of VistA RPCs under the VBECS namespace in the Remote Procedure File (#8994) that are called by the VistA VistALink Listener client-server software. These calls are not public utilities and may be subject to change.

Table 7: Remote Procedure Calls

RPC Name	Database Integration Agreement (DBIA)	This RPC:
VBECS Order Entry	4619	Supports order entry of Blood Bank requests from the Blood Bank order entry dialog in CPRS.
VBECS Patient Available Units	4620	Provides a list of assigned, crossmatched, autologous and directed blood units that are available for a patient.
VBECS Patient Transfusion History	4621	Provides a list of past transfusions performed for a patient.
VBECS Blood Products	4622	Provides a list of orderable blood products, or component classes, to the VistA Surgery package.
VBECS Patient Report	4623	Provides patient specimen testing results, component requests, and available blood units for a patient to be displayed in CPRS.
VBECS Patient ABO_RH	4624	Provides the most current ABO Group and Rh Type identified for a patient.
VBECS Patient ABID	4625	Provides a list of antibodies identified for a patient.
VBECS Patient TRRX	4626	Provides a list of transfusion reactions for a patient.
VBECS Workload Capture	4627	Provides Blood Bank workload data to the VistA Laboratory Service package for workload reporting to national and local entities.
VBECS Workload Update Event	4628	Inserts completed workload-related data into the VBECS database after the VistA Laboratory Services package has completed workload-reporting transactions. Upon completion of the insert, the RPC returns an XML response to the VBECS Application Interfacing Support Software that initiated the communication indicating a successful or unsuccessful transaction.
VBECS Accession Area Lookup	4607	Provides a list of all Laboratory Blood Bank Accession Areas in VistA and their associated divisions to VBECS for workload reporting purposes.
VBECS Blood Bank User Lookup	4608	Returns a list of all Blood Bank users identified in the VistA system to VBECS. Blood Bank users are identified by the Security Keys of either LRBLOODBANK or LRBLSUPER.
VBECS Division Lookup	4609	Returns a list of all VAMC divisions associated with a VistA system.
VBECS HCPCS Codes Lookup	4610	Returns a list of Blood Bank related HCPCS codes to be associated with processes, or procedures, performed in VBECS.
VBECS Laboratory Test Lookup	4611	Returns a list of VistA Laboratory tests to be associated with blood components in VBECS.
VBECS Lab Test Results Lookup	4612	Returns a list of VistA Laboratory test results for a patient.
VBECS Medication Profile Lookup	4613	Returns a list of medications for a patient from the VistA Pharmacy package.
VBECS Lab Accession UID Lookup	4614	Returns data from the VistA Laboratory Services package based on a Lab order number. The data is used to validate a VBECS specimen test request for a patient and specimen received in the Blood Bank for that test.
VBECS Workload Codes	4615	Returns a list of Blood Bank related workload related data that

RPC Name	Database Integration Agreement (DBIA)	This RPC:
Lookup		is associated with processes in VBECS.
VBECS Patient Lookup	4616	Provides a patient lookup function using standard VistA patient lookup criteria. A list of matching patients found in the lookup is returned to VBECS along with required patient identifiers and demographics.
VBECS Provider Lookup	4617	Provides a lookup of VistA users that hold the PROVIDER security key.
VBECS Hospital Location Lookup	4618	Returns a list of hospital locations associated with a division in VistA.
VBECS Lab Order Lookup by UID	4633	Returns a list of Laboratory Services data related to an order based on a specimen UID.
VBECS Dss Extract	4956	Provides BloodBank post-transfusion related data to the VistA DSS Blood Bank Extract application for DSS reporting.

VBECS Windows Services



Changes made to individual HL7 listeners must be validated in the test account before using in production.




*The VBECS Service Monitor must be stopped before stopping another VBECS service: the VBECS Service Monitor will attempt to restart any VBECS service that was stopped. This service needs to be stopped within the Cluster Administrator. The Cluster Administrator utility can be accessed from the cluster server by clicking **Start, Administrative Tools, Cluster Administrator** (Figure 97).*

Stopping the service through the services window in the control panel will not stop the service.

VBECS uses Microsoft Windows Services (services) to provide minimal downtime and minimal user interaction. These services are installed on each physical server of the VBECS cluster server group. The Cluster Administrator controls the state and operation of the VBECS services. See Table 8 for a complete listing of VBECS services. The Install VBECS Services and the VBECS Application section of the VistA Blood Establishment Computer Software (VBECS) Installation Guide describe how these services are installed. For details on stopping and starting VBECS services see the Restarting VBECS Services section.

Figure 97: Example of VBECS Services in Cluster Administrator



Reconfiguring the VBECS HL7 Multi Listener and VistALink Services

VBECS HL7 Multi Listener Service

If changes need to be made to the configuration of the VBECS HL7 Multi Listener service due to a change in IP address or port number, first take the VBECS Service Monitor resource and VBECS HL7 Multi Listener resources offline. Navigate to the C:\Program Files\Vista\VBECS\WinServices\VBECS HL7 Multi Listener, and locate the file named VbecsHL7ListenerService.exe.config. The file contents will look similar to the following example:

```
<?xml version="1.0" encoding="utf-8"?>
<configuration>
  <appSettings>0
    <add key="PrimaryDbConnectionString" value="Connection Timeout=90;Data
Source=VHAISHSQLZ1;Initial Catalog=VBECS_V1_PROD;persist security info=False;packet
size=8192;integrated security=SSPI;Application Name=VBECS HL7 Multi Listener" />
    <add key="serviceName" value="VBECS HL7 Multi Listener" />
    <add key="allowPing" value="true" />
    <add key="listenerIpAddress" value="10.3.21.82" />
    <add key="listenerPortNumber" value="21994" />
    <add key="monitorService" value="true" />
    <add key="monitorInterval" value="5000" />
    <add key="monitorMaxRetries" value="3" />
    <add key="monitorServiceStartTimeout" value="5" />
    <add key="BuildNumber" value="1.0.6.2" />
  </appSettings>
</configuration>
```

Modify the value for the key named listenerIpAddress and the value for the key named listenerPortNumber. Save the file, close it and bring the VBECS HL7 Multi Listener and the VBECS

Service Monitor resources online. Repeat the update of the configuration file on the other server. There is no need to bring any more resources online; the Cluster Administrator handles both nodes at the same time.

Test account: The test account listener (VBECS Test HL7 Multi Listener) is changed in the same manner. It is located at C:\Program Files\Vista\VBECS Test\WinServices\VBECS Test HL7 Multi Listener.

If troubleshooting requires use of the other listener services, take the VBECS Service Monitor and VBECS HL7 Multi Listener resources offline. Bring the single listeners online as required. Once they are configured properly in the Configure Interfaces section of this guide, then bring the VBECS Service Monitor resource online.

All of the services communicate directly with the VBECS database. Therefore, prior to restoring the database, all of the VBECS service must be stopped and restarted accordingly.

VBECS VistALink Service

If changes need to be made to the configuration of the VBECS VistALink RPC XML Listener service due to a change in IP address or port number, first stop the VBECS Service Monitor service, then stop the VBECS VistALink RPC XML Listener service. Navigate to the c:\Program Files\Vista\VBECS\WinServices\VBECS VistALink RPC XML Listener, and locate the file named VistALink.Listener.WinService.exe.config. The file contents will look similar to the following example:

```
<?xml version="1.0" encoding="utf-8"?>
<configuration>
  <configSections>
    <sectionGroup name="VistALink">
      <section name="RpcList"
type="gov.va.med.vbecs.DAL.VistALink.Listener.Core.RpcListConfigSectionHandler,VistALink.Listen
er.Core" />
    </sectionGroup>
  </configSections>
  <appSettings>
    <add key="PrimaryDbConnectionString" value="Connection Timeout=90;Data
Source=vhaishsqlz1;Initial Catalog=VBECS_V1_PROD;persist security info=False;packet
size=8192;integrated security=SSPI;Application Name=VBECS VistALink RPC XML Listener" />
    <add key="serviceName" value="VBECS VistALink RPC XML Listener" />
    <add key="serverName" value="vhaishsqlz1" />
    <add key="databaseName" value="VBECS_V1_PROD" />
    <add key="listenerPortNumber" value="21992" />
    <add key="allowPing" value="true" />
    <add key="listenerIpAddress" value="10.3.21.81" />
    <add key="monitorService" value="true" />
    <add key="monitorInterval" value="3000" />
    <add key="monitorMaxRetries" value="3" />
    <add key="monitorServiceStartTimeout" value="5" />
    <add key="BuildNumber" value="1.0.6.2" />
  </appSettings>
```

<VistALink>

Modify the value for the key named listenerIpAddress and the value for the key named listenerPortNumber. Save the file, close it and restart the VBECS VistALink RPC XML Listener service and the VBECS Service Monitor service.

Test account: The test listener (VBECS Test VistALink RPC XML Listener) is changed in the same manner. It is located at C:\Program Files\Vista\VBECS Test\WinServices\VBECS Test VistALink RPC XML Listener.

All VBECS services start with the VBECS namespace prefix. There are duplicate services for production and test accounts that provide functionality for their respective databases.

Table 8: Windows Service Manager

Windows Service Name	This Service:
VBECS CPRS HL7 Client Monitor	The startup type is set to manual. The cluster administrator will manage the starting of this service. It polls the VBECS Production database for HL7 update messages to be sent to CPRS in the Vista Production account.
VBECS CPRS HL7 Listener	Is initially installed as disabled. It is a single listener HL7 service for the Production CPRS HL7 interface. It should be used only as a backup for the VBECS HL7 Multi Listener service or for troubleshooting HL7 interface problems so that other HL7 interfaces using the multi listener are not adversely affected.
VBECS HL7 Multi Listener	The startup type is set to manual. The cluster administrator will manage the starting of this service. This is the default HL7 listener service for all Production HL7 interfaces.
VBECS Patient Merge HL7 Listener	Is installed as disabled. It is a single listener HL7 service for the Production Patient Merge HL7 interface. It should be used only as a backup for the VBECS HL7 Multi Listener service or for troubleshooting HL7 interface problems so that other HL7 interfaces using the multi listener are not adversely affected.
VBECS Patient Update HL7 Listener	Is installed as disabled. It is a single listener HL7 service for the Production Patient Update HL7 interface. It should be used only as a backup for the VBECS HL7 Multi Listener service or for troubleshooting HL7 interface problems so that other HL7 interfaces using the multi listener are not adversely affected.
VBECS Scheduled Report Runner	The startup type is set to manual. The cluster administrator will manage the starting of this service. It runs scheduled VBECS reports for the Production database.
VBECS VistALink RPC XML Listener	The startup type is set to manual. The cluster administrator will manage the starting of this service. It provides a client-server TCP/IP listener service for VistALink RPC XML messages from the VAISS APIs. It calls VBECS RPCs to provide Blood Bank data from the VBECS Production database to Vista Production account applications.
VBECS Test CPRS HL7 Client Monitor	The startup type is set to manual. The cluster administrator will manage the starting of this service. It polls the VBECS Test database for HL7 update messages to be sent to CPRS in the Vista Test account.
VBECS Test CPRS HL7 Listener	Is installed as disabled. It is a single listener HL7 service for the Test CPRS HL7 interface. It should be used only as a backup for the VBECS Test HL7 Multi Listener service or for troubleshooting HL7 interface problems so that other HL7 interfaces using the multi listener are not adversely affected.

Windows Service Name	This Service:
VBECs Test HL7 Multi Listener	The startup type is set to manual. The cluster administrator will manage the starting of this service. This is the default HL7 listener service for all Test HL7 interfaces.
VBECs Test Patient Merge HL7 Listener	Is installed as disabled. It is a single listener HL7 service for the Test Patient Merge HL7 interface. It should be used only as a backup for the VBECs Test HL7 Multi Listener service or for troubleshooting HL7 interface problems so that other HL7 interfaces using the multi listener are not adversely affected.
VBECs Test Patient Update HL7 Listener	Is installed as disabled. It is a single listener HL7 service for the Test Patient Update HL7 interface. It should be used only as a backup for the VBECs Test HL7 Multi Listener service or for troubleshooting HL7 interface problems so that other HL7 interfaces using the multi listener are not adversely affected.
VBECs Test Scheduled Report Runner	The startup type is set to manual. The cluster administrator will manage the starting of this service. It runs scheduled VBECs reports for the Test database.
VBECs Test VistALink RPC XML Listener	The startup type is set to manual. The cluster administrator will manage the starting of this service. It provides a client-server TCP/IP listener service for VistALink RPC XML messages from the VAISS APIs. It calls VBECs RPCs to provide Blood Bank data from the VBECs Test database to VistA Test account applications.
VBECs Service Monitor	The startup type is set to manual. The cluster administrator will manage the starting of this service. It monitors all VBECs Production and Test services to ensure that they are running and accepting incoming requests, where appropriate.

This page intentionally left blank.

Troubleshooting

Performance Improvements


Stopping and Starting VBECS Test Services

The VBECS Test Services needs to be stopped when a newly released patch is completely installed in the Production environment and no further testing is required in the Test Environment. Stopping the services will increase the overall performance of the system because only a few services will be running.

Stopping VBECS Test Services


- 1) Click **Start, Administrative Tools, Cluster Administrator**.
- 2) If Open Connection to Cluster window does not appear, click **File, Open Connection**.
- 3) Type **<CLUSTER_NAME>** in the **Cluster or server name** field and click **OK** (Figure 98).

Figure 98: Example of Open Connection to Cluster




- 4) The Cluster Administrator window populates. Expand the Groups folder, and verify that **VBECS Group PROD**, **VBECS Group SM** and **VBECS Group TEST** exists as shown in (Figure 99).

Figure 99: Example of All VBECS Groups Services Online



5) Right-click on **VBECS Group TEST** and select **Take Offline** (Figure 100).

Figure 100: Example of VBECS Group Test Services Offline




The VBECS Test Services needs to be started when installing a new patch in the Test Environment and during the testing phase.

Starting VBECS Test Services


- 1) Click **Start, Administrative Tools, Cluster Administrator**.
- 2) If Open Connection to Cluster window does not appear, click **File, Open Connection**.
- 3) Type **<CLUSTER_NAME>** in the **Cluster or server name** field and click **OK** (Figure 101).

Figure 101: Example of Open Connection to Cluster




- 4) The Cluster Administrator window populates. Expand the Groups folder, and verify that **VBECS Group PROD**, **VBECS Group SM** and **VBECS Group TEST** exists as shown in Figure 102.

Figure 102: Example of VBECS Group TEST Services Offline



- 5) Right-click on **VBECs Group Test** and select **Bring Online** (Figure 103).

Figure 103: Example VBECs Group TEST Services Online




Verify NIC Card Configuration

If the VBECs application experiences network latency issues, such as problems when scanning barcodes, check the NIC card configuration settings.

- 1) Log into Server #1.
- 2) Click Start, Control Panel, Network Connections, Public. Click Properties.
- 3) Click **Configure** (Figure 104).

Figure 104: Example of Public Properties




- 4) Click on the **Advanced** tab (Figure 105).

Figure 105: Example of NIC Properties



- 5) Click **Speed and Duplex** (Figure 106) (e.g. 100Mb Full).


Figure 106: Example of HP NC7782 Gigabit Server Adapter Properties



- 6) Verify NIC Card Configuration section with the Switch Port Speed.
- 7) If both values are the same, click **Cancel** and continue to Step 11 for Server #2.

8) If the values are different, make the values match. (Figure 107).

Figure 107: Example of Updated HP NC7782 Gigabit Server Adapter Properties



9) Click **OK**.

10) The remote desktop reconnection message popup will be received (Figure 108).

Figure 108: Example of Reconnecting Message




11) Log off Server #1 when the remote session is restored. **Log into Server #2.**

12) Click Start, Control Panel, Network Connections, Public. Click Properties.

13) Click **Configure** (Figure 109).

Figure 109: Example of Public Properties




14) Click on the **Advanced** tab (Figure 110).

Figure 110: Example of NIC properties




15) Click **Speed and Duplex** (Figure 111) (e.g., 100Mb Full).

Figure 111: Example of HP NC7782 Gigabit Server Adapter Properties



- 16) Verify NIC Card Configuration section with the Switch Port Speed (Figure 112).
- 17) If both values are the same, click **Cancel** and do not proceed with these remaining steps.
- 18) If the values are different, make the values match.

Figure 112: Example of Update HP NC7782 Gigabit Server Adapter Properties



- 19) Click **OK**.
- 20) The remote desktop reconnection message popup will be received (Figure 113).


Figure 113: Example of Reconnecting Message



- 21) After the remote session is restored, click **Start, Administrative Tools, Cluster Administrator**.


22) If the Passive Cluster Node (Server #2) is marked  (Figure 114).

Figure 114: Example of Passive Cluster Node Offline




23) Wait a few minutes for the Passive Cluster Node (Server #2) to come back online  (Figure 115).

Figure 115: Example of Passive Cluster Node Online



24) Verify that all Active Resources of the Active Cluster Node have State marked Online (Figure 116).

Figure 116: Example of Active Cluster Node Resources Online



 If resource(s) state remains offline, please file a Remedy ticket immediately.


25) Log off Server #2.

VistA Query Timeout

The VistA cache refresh interval is the time (in seconds) that VBECS waits before it attempts to copy new VistA data to the VBECS database (to cache it). VistA data is cached for Workload Codes, CPT Codes, HCPCS Codes, and Hospital Locations.


- 1) To update the refresh interval, locate the VBECS.exe.config file in the installation directory for VBECS: C:\Program Files\VistA\VBECS (Figure 117).

Figure 117: Example of a Directory Structure




2) To open the file, right click it. Select **Notepad** (Figure 118). Click **OK**.

Figure 118: Example of the Open With Dialog



- 3) In the VBECS.exe.config file, find the entry for “VistACacheRefreshIntervalInMinutes” (Figure 119).

Figure 119: Example of a Configuration File



- 4) Edit the value to whatever is required. Save the file. This value is in minutes, so the current value of 1440 minutes is equivalent to 24 hours (to convert minutes to hours, divide by 60).

VBECS Exception Logging

VBECS logs all errors that occur in the system in the application event viewer on the cluster. A user defined as an administrator on the cluster can connect to the cluster through Remote Desktop Connection to view these errors.

- 1) Click **Start, Control Panel, Administrative Tools**.
- 2) Open the Event Viewer and see the application section to view the errors that VBECS logs.
- 3) Double click the application icon on the right side of the screen list view.
- 4) In the list view on the right side of the screen, click the date column header to sort the errors by date.
- 5) Evaluate “Error” and warning errors and submit a Remedy ticket if the error was logged at the same time a VBECS user reported an error. Ignore informational messages. The VBECS development and maintenance team will investigate the ticket.

VBECS Exception Workarounds

When an exception occurs in VBECS, click **Details**. Copy the details to the clipboard. Include all details of the exception in the Remedy ticket. A common exception that occurs within VBECS was traced to a Microsoft .NET 2003 problem that will not be resolved until VBECS is upgraded with the implementation of Microsoft .NET 2005. The exception shows in the details:

1) Exception Information

Exception Type: System.NullReferenceException

Message: Object reference not set to an instance of an object.

TargetSite: IntPtr CallWindowProc(IntPtr, IntPtr, Int32, IntPtr, IntPtr)

HelpLink: NULL

Source: System.Windows.Forms

StackTrace Information

at System.Windows.Forms.UnsafeNativeMethods.CallWindowProc(IntPtr wndProc, IntPtr hWnd, Int32 msg, IntPtr wParam, IntPtr lParam)

at System.Windows.Forms.NativeWindow.DefWndProc(Message& m)

at System.Windows.Forms.Control.DefWndProc(Message& m)

at System.Windows.Forms.Control.WmUpdateUIState(Message& m)

at System.Windows.Forms.Control.WndProc(Message& m)

at System.Windows.Forms.ScrollableControl.WndProc(Message& m)

at System.Windows.Forms.ContainerControl.WndProc(Message& m)

at System.Windows.Forms.ParkingWindow.WndProc(Message& m)

at System.Windows.Forms.ControlNativeWindow.OnMessage(Message& m)


at System.Windows.Forms.ControlNativeWindow.WndProc(Message& m)

at System.Windows.Forms.NativeWindow.Callback(IntPtr hWnd, Int32 msg, IntPtr wparam, IntPtr lparam)

This exception occurs randomly when a screen is loading. When this occurs, the user must click **Shut down** on the exception message and try the option again.

When the user prints a report that accepts a given date range, a Crystal Report Windows Forms Viewer window may appear (Figure 120).

Figure 120: Crystal Reports Message



The user may change the date range given (alter the start or end date by plus or minus one day) to resolve this problem. (This documented Crystal problem will be fixed in a future version of VBECS when Crystal Reports is upgraded.)

Restarting VBECS Services


When troubleshooting VBECS application interfaces, it may be necessary to stop and restart the VBECS services. To do so, use the Cluster Administrator utility. (Do not use the Services utility found under the Administrative Tools.) VBECS services are organized into three groups:

- **VBECS Group PROD** contains the services for the VBECS production environment.
- **VBECS Group TEST** contains the services for the VBECS test environment.
- **VBECS Group SM** contains the monitoring services used by both VBECS environments.

To manipulate VBECS services using Cluster Administrator:


- 1) Click **Start, Administrative Tools, Cluster Administrator**.
- 2) If prompted, enter the cluster alias or IP address in the **Cluster or Server name field** and click **OK** (Figure 121).

Figure 121: Opening a Connection in Cluster Administrator



- 3) Navigation within Cluster Administrator (See Figure 122):
 - a) Click the **Resources** folder in the left panel to populate the right panel with a list of the active resources.
 - b) To stop a service, right-click the service Name and Group combination and select **Take Offline**.
 - c) To start a service, right-click the service Name and Group combination and select **Bring Online**.

Figure 122: Troubleshooting VBECS Services with Cluster Administrator



VBECs Application Interfaces

Table 9: Troubleshooting VBECs Application Interfaces

Source	Description of Problem	Possible Cause	Solution
VBECs: Order Alerts and Pending Order List	New orders or cancellations of existing orders in CPRS are not showing up in VBECs.	The OERR-VBECs Logical Link is not running on the VistA system.	Start the OERR-VBECs Logical Link.
		The VBECs HL7 Multi Listener Windows Service is not running or is locked on the VBECs Cluster server.	Start or restart the VBECs HL7 Multi Listener Windows Service.
		Network connectivity issue	Contact local system support.
		The HL7 message is missing patient last or first name or one or more name components length(s) exceed(s) the VBECs maximum supported value.	VBECs responds to the new order request with an application reject (AR) acknowledgement message indicating Patient Name(s) not found in HL7 Message or Patient's Name(s) field size(s) exceed(s) VBECs maximum supported value. Rejected patient order messages due to invalid patient name message content are recorded on the Windows Event Log and an email message is sent to the interface failure alert recipient set in VBECs Administrator for immediate action.
VBECs Admin: Configure Division	New orders are not showing up in VBECs.	Order mappings to institutions within a division's configuration were changed.	Stop and restart the VBECs HL7 Multi Listener Service.
VBECs: Patient Update Alerts	VistA patient updates are not showing up in VBECs.	The patient being updated in VistA is not in the VBECs Patient table and is, therefore, not a Blood Bank patient.	No action is required.
		The fields that were updated in VistA are not stored in VBECs, therefore, no data will be updated.	No action is required.
		The Taskman scheduled option VAFC BATCH UPDATE is not scheduled to run or has not reached the time limit in the schedule.	Schedule the VAFC BATCH UPDATE option to run at the desired increment or use the option "One-time Option Queue" in the Taskman Management Options to start the task.
		The VBECsPTU Logical Link is not running on the VistA system.	Start the VBECsPTU Logical Link.
		The VBECs HL7 Multi Listener Windows Service is not running or is locked on the VBECs Cluster server.	Start or restart the VBECs HL7 Multi Listener Windows Service.
		Network connectivity issue	Contact local system support.

Source	Description of Problem	Possible Cause	Solution
		The HL7 message is missing patient last or first name or one or more name components length(s) exceed(s) the VBECS maximum supported value.	VBECS responds to the patient update request with an application reject (AR) acknowledgement message indicating Patient Name(s) not found in HL7 Message or Patient's Name(s) field size(s) exceed(s) VBECS maximum supported value. Rejected patient update messages due to invalid patient name message content are recorded on the Windows Event Log and an email message is sent to the interface failure alert recipient set in VBECS Administrator for immediate action.
VBECS: Patient Merge Alerts	VistA Patient Merge events are not showing up in VBECS.	The two patient identifiers in the merge do not exist in VBECS and, therefore, cannot be merged.	No action is required.
		The VBECPTM Logical Link is not running on the VistA system.	Start the VBECSPTM Logical Link.
		The VBECS HL7 Multi Listener Windows Service is not running or is locked on the VBECS Cluster server.	Start or restart the VBECS HL7 Multi Listener Windows Service.
		Network connectivity issue	Contact local system support.
		The HL7 message is missing patient last or first name or one or more name components length(s) exceed(s) the VBECS maximum supported value.	Failed patient merge messages due to invalid patient name message content are recorded on the Windows Event Log and an email message is sent to the interface failure alert recipient set in VBECS Administrator for immediate action.
VistA: HL7 System Link Monitor	The VistA HL7 System Link Monitor shows more MESSAGES TO SEND than MESSAGES SENT for the OERR-VBECS Logical Link and is hung in an "Open" state.	The VBECS HL7 Multi Listener Windows Service is not running or is locked on the VBECS Cluster server.	Start or restart the VBECS HL7 Multi Listener Windows Service.
		Network connectivity issue	Contact local system support.
	The VistA HL7 System Link Monitor shows more MESSAGES TO SEND than MESSAGES SENT for the VBECSPTU Logical Link and is hung in an "Open" state.	The VBECS HL7 Multi Listener Windows Service is not running or is locked on the VBECS Cluster server.	Start or restart the VBECS HL7 Multi Listener Windows Service.
		Network connectivity issue.	Contact local system support.

Source	Description of Problem	Possible Cause	Solution
	The VistA HL7 System Link Monitor shows more MESSAGES TO SEND than MESSAGES SENT for the VBECSPTM Logical Link and is hung in an "Open" state.	The VBECs HL7 Multi Listener Windows Service is not running or is locked on the VBECs Cluster server.	Start or restart the VBECs HL7 Multi Listener Windows Service.
		Network connectivity issue.	Contact local system support.
CPRS: Orders Tab	CPRS does not display the correct status of a Blood Bank order after it was updated in VBECs.	The VBECs CPRS Client Monitor Windows Service is not running or is locked on the VBECs Cluster server.	Start or restart the VBECs CPRS Client Monitor Windows Service.
		The VBECs-OERR Logical Link is not running.	Start the VBECs-OERR Logical Link in Background mode.
		Network connectivity issue	Contact local system support.
CPRS: Blood Bank Order Dialog	CPRS displays "Not able to open port" message in Patient Information screen in Blood Bank Order Dialog.	The VBECs VistALink XML RPC Listener Service is not running or is locked on the VBECs Cluster server.	Start or restart the VBECs VistALink XML RPC Listener Service.
		Network connectivity issue	Contact local system support.
CPRS: Reports Tab, Blood Bank Report	CPRS displays "---- BLOOD BANK REPORT IS UNAVAILABLE----"	The VBECs VistALink XML RPC Listener Service is not running or is locked on the VBECs Cluster server.	Start or restart the VBECs VistALink XML RPC Listener Service.
		Network connectivity issue.	Contact local system support.
CPRS: Blood Bank Order Dialog: Signing an Order	CPRS displays an "Error Saving Order" dialog screen with the text "The error, One or more orders to the VBECs system failed and are queued for later delivery."	An error occurred in the VBECs HL7 Multi Listener Windows Service, which caused a failure to respond to CPRS with acceptance.	Log onto the VBECs Cluster Server and review the System Application Event Log for error details. Click Start, Administrative Tools, Event Viewer. Select Application.
		Network connectivity issue.	Contact local system support.
VBECs Cluster Server Application Event Log: Source is VBECs SimpleListener	An application error has been logged to the Event Log where the Message under Exception Information is "Could not access 'CDO.Message' object."	The HL7 Multi Listener Windows Service has encountered an error trying to send an email message to the Interface Administrator.	Disable port 25 blocking in McAfee. Open the VirusScan Console and select Access Protection. Click the Task menu option, the Properties. Uncheck Prevent mass mailing worms from sending mail, port 25 under Ports to block.
		If the IP address is associated with the local VistA system, the HL7 Application Parameters in VistA were not set up correctly for the supported protocols.	Refer to the VBECs Application Interfacing Support Software Installation and User Configuration Guide for HL7 setup procedures in VistA.
		If the IP address is not from the local VistA system, a rogue HL7 system is sending messages to the VBECs server.	Contact IRM to identify the location of the server with which the IP address is associated. Notify the site that the message is coming from the problem so that the messages can be routed to the correct location.
	The IP address in the description of the error will indicate where the message is coming from.		

Source	Description of Problem	Possible Cause	Solution
VBECs Cluster Server Application Event Log: Source is VBECs HL7 MailServer	An application error was logged in the Event Log with the source of VBECs HL7 MailServer where the Message under Exception Information is, "Could not access 'CDO.Message' object."	The HL7 Multi Listener Windows Service encountered an error trying to send an email message to the Interface Administrator.	Disable port 25 blocking in McAfee. Open the VirusScan Console and select Access Protection. Click the Task menu option, Properties. Uncheck Prevent mass mailing worms from sending mail, port 25 under Ports to block.
VBECs Cluster Server Application Event Log: Source is CPRS HL7 Parser	An HL7 message sent from CPRS to VBECs was rejected. The description in the Event Log is "Exception message: Division [<i>division</i>] is not supported by this instance of VBECs." An HL7 message sent from CPRS to VBECs was rejected. The description in the Event Log is "Exception message: Division [<i>division</i>] is not active in this instance of VBECs."	An invalid or unsupported division associated with the Patient Location was selected in CPRS when the order was created. The division associated with the Patient Location that was selected in CPRS when the order was created is not active in VBECs.	The order must be created in CPRS again with a valid Patient Location associated with a VBECs-supported division. The order must be created in CPRS again with a valid Patient Location associated with a VBECs-active division.

VBECs Build Version Numbers

VBECs builds are numbered as "Major.Minor.Patch.Build." "Major" is the version of the product. The "Minor" number is incremented for minor system changes. The "Patch" number is incremented for minor bug fixes. The "Build" number is incremented with each build but is not displayed publically to customers. For example, "1.2.1.0" represents the first version of VBECs with two minor system changes and one patch. VA Product Support requires the full four digits of the VBECs version number.

Cluster Connectivity Lost

Problem: Connections to the cluster are lost. The cluster is not pingable by name or IP address, but individual nodes are still up.

Probable Cause: A network outage that affects both nodes simultaneously will cause the cluster to fail.

Solution:

1. Log into one of the cluster nodes and restart. Wait 1 minute.
2. Restart the other cluster node.
3. After the node in #1 has finished rebooting, verify that the cluster is back up.
4. When both nodes have restarted, stop and start services per the instructions in the previous section.

Printing Fails to Report Printer

Problem: The printer fails to print.

Probable Cause: A printer name is not consistent with what is configured in VBECs or a driver is incorrect.

Solution:

Verify Printer Name

1. Log into VBECS Administrator and note the default printer in Configure Division.
2. Verify that the printer name on the server is consistent with the name noted in step 1.
3. If still broken, verify printer drivers are consistent.

Verify Printer Drivers

1. Log into one of the servers with administrator rights.
2. Open **Control Panel, Printers and Faxes**.
3. Double click the printer noted in step 1 under **Verify Printer Name**.
4. Select **Printer, Properties** and click the **Advanced** tab.
5. Note the driver name in the **Driver** field.
6. Repeat Steps 1 through 5 on the other server. If drivers are inconsistent, update the server that is not working with the correct driver.

Zebra Printer Problems

Problem: The printer prints, but there is not text on the label or text is too light.

Probable Cause: The printer is out of ribbon or the DARKNESS setting is too light (Figure 123).

Solution: Increase the DARKNESS setting after verifying printer has ribbon.

Figure 123: Example Zebra Printer Settings


View Printer Configuration	
VA 060876.06 GY090205.34901-010.E.VT	
+10	DARKNESS
2 IPS	PRINT SPEED
+000	TEAR OFF
TEAR OFF	PRINT MODE
NON-CONTINUOUS	MEDIA TYPE
WEB	SENSOR TYPE
AUTO SELECT	SENSOR SELECT
THERMAL-TRANS.	PRINT METHOD
105 08/12 MM	PRINT WIDTH
1221	LABEL LENGTH
39.0IN 988MM	MAXIMUM LENGTH
BIDIRECTIONAL	PARALLEL COMM.
RS232	SERIAL COMM.
9600	BAUD
8 BITS	DATA BITS
NONE	PARITY
XON/XOFF	HOST HANDSHAKE
NONE	PROTOCOL
000	NETWORK ID
NORMAL MODE	COMMUNICATIONS
<~> 7EH	CONTROL PREFIX
<^> 5EH	FORMAT PREFIX
<,> 2CH	DELIMITER CHAR
ZPL II	ZPL MODE
CALIBRATION	MEDIA POWER UP
CALIBRATION	HEAD CLOSE

Problem: The printer doesn't print. It also cannot be pinged or be seen in a web browser (Figure 124).

Probable Cause: Network settings are not correct on the printer

Solution: Correct the printer's network settings (see the section titled "Set the IP Address on the Printer").

Figure 124: Zebra Printer Web Console



Problem: The printer doesn't print and network settings have been verified (see previous).

Probable Cause: One or more settings are incorrect.

Solution: Verify that the PRINT METHOD, CONTROL PREFIX, FORMAT PREFIX, DELIMITER CHAR and ZPL MODE match the settings in Figure 123.

Scanner Problems

Problem: When scanning, characters appear in the field that do not match the label being scanned. Often, the bad characters are not alphanumeric.

Probable Cause: Network latency causes data to become corrupted.

Solution: The lab supervisor will program an inter-character delay into the scanner to fix the issue. This puts a small time delay between each character as it is sent over the network, which results in slightly slower scan speeds.

Figure 125 through Figure 132 are configuration barcodes arranged from a 10 millisecond inter-character delay all the way up to an 80 millisecond delay respectively. We suggest that you start with the 10 millisecond delay. If that does not resolve the problem, proceed with larger delays until the problem is corrected.

Note that these barcodes include all of the configuration information for the scanners. There is no need to scan any additional barcodes to configure the scanner.

Figure 125: 10 milliseconds



Figure 126: 20 milliseconds



Figure 127: 30 milliseconds



Figure 128: 40 milliseconds



Figure 129: 50 milliseconds



Figure 130: 60 milliseconds



Figure 131: 70 milliseconds



Figure 132: 80 milliseconds



This page intentionally left blank.

Archiving and Recovery

The VBECS database will be backed up once daily at an established time to a tape drive. If a disaster occurs, the data in VBECS can be recovered from the backup media.

Assumptions

- The SQL Server job that backs up the database is running correctly.
- Replacement hardware will have a tape drive that is compatible with the one lost in the disaster.

Outcome

- VBECS data is successfully recovered.

Limitations and Restrictions

- Only the VBECS data is backed up. The operating system is not backed up. In the event of a disaster, the operating system will have to be reinstalled and configured.

Additional Information

- None

VBECS Backup

If your servers are maintained at a data center, ignore this section since data center personnel will perform this task.

To preserve VBECS data in case of database corruption or destruction of hardware, the VBECS databases are copied over to shared storage via a scheduled job configured with the VBECS installation. VBECS is comprised of the following SQL databases: VBECS_V1_PROD and VBECS_V1_PROD_MIRROR (production) VBECS_V1_TEST and VBECS_V1_TEST_MIRROR (test VBECS account). Both production and test share the use of the msdb and master SQL databases. It is critical that every VBECS database is backed up nightly to tape. Remove the tape and take it to another location in accordance with local policy. For more technical details on backups, see *VistA Blood Establishment Computer Software (VBECS) Installation Guide*. For details on tape storage and backup frequency, refer to local policy.

VBECS Recovery



Unauthorized modifications will render this device an adulterated medical device under Section 501 of the Medical Device Amendments to the Federal Food, Drug, and Cosmetic Act. Acquiring and implementing this software through the Freedom of Information Act require the implementer to assume total responsibility for the software and become a registered manufacturer of a medical device, subject to FDA regulation.

If your servers are maintained at a data center, ignore this section since data center personnel will perform this task.

File a remedy ticket in the event of a disaster that destroys or damages the VBECS system. The VBECS team and VA Product Support will work to recover or rebuild the system.


Reinstall the System

If your servers are maintained at a data center, ignore this section since data center personnel will perform this task.

 *This section should not be followed once application data has been entered. Following these steps will cause all VBECS application data to be lost.*

- 1) Install the image on the server hard drive.
- 2) Reinstall VBECS using *VistA Blood Establishment Computer Software (VBECS) Installation Guide*.
- 3) Make sure all VBECS Services are stopped on both servers. All VBECS service names begin with “VBECS” (Figure 133). To stop a service, open Cluster Administrator and take all VBECS Services offline.


Figure 133: Example of VBECS Services



- 4) Log onto the server that is connected to the tape drive and has Backup Exec installed on it. Log in as an Administrator.


- In Cluster Administrator (Figure 134), make sure this node is the active node in the cluster. If not, drag Cluster Group and SQL Server to the Active Groups folder of this node to make it the active node.

Figure 134: Example of Cluster Administrator



- Click **Start, All Programs, Symantec Backup Exec 10d for Windows Servers**. The main Backup Exec console is displayed (Figure 135).


Figure 135: Example of Backup Exec Console



Inventory the Tape


- 1) Place the tape that reflects the most recent system backup in the tape drive.
- 2) Click the **Devices** button (Figure 136).
- 3) Right click **HP 1** under the server node (not the drive pool).
- 4) Select **Inventory**. The Inventory Job Properties window appears (Figure 137).

Figure 136: Example of Devices




- 5) Click **Run Now**. Click **OK** to close information messages that appear.

Figure 137: Example of Inventory



- Click **Job Monitor** (Figure 138) and make sure the job completed successfully.


Figure 138: Example of Successful Inventory



Catalog the Tape


- Click **Devices** again. Right click **HP 1** under the server node.
- Select **Catalog**. The Catalog Job Properties window appears (Figure 139). Click **Run Now**. Click **OK** to close information messages that appear.

Figure 139: Example of Catalog



- 3) Click **Job Monitor** (Figure 140) and make sure the job completed successfully.


Figure 140: Example of Successful Catalog



Restore Files


- 1) Click **Restore**.
- 2) Select all four folders under temp\Backup (Figure 141).

Figure 141: Example of Restore Properties




- 3) Create the “temp\yyyymmddBackup” directory on the D: drive (Figure 142).

Figure 142: Example of Backup Directory




- 4) Click **File Redirection** on the left (Figure 143). Click the **Redirect file sets** check box.
- 5) In the Restore to drive field, enter **D:** (Backup Exec automatically populates the field with the server name).
- 6) In the Restore to path field, enter **D:\temp\yyyymmddBackup** (yyyymmdd represents the current date).
- 7) Click **Run Now**.
- 8) Click **OK** on information messages that appear.

Figure 143: Example of Restore Properties



9) Click **Job Monitor** (Figure 144) and make sure the job completed successfully.

Figure 144: Example of Successful Restore



Restore the Databases



If you find the need to perform a database restore, contact customer support to have qualified personnel assist you with the database restore.

VA Service Desk Primary Contact

For Information Technology (IT) support, call the VA Service Desk (VASD), 888-596-HELP (4357) toll free, 24 hours per day, 7 days per week. [Users with access to the VASD-supported request tool (e.g., Remedy) may file a ticket in lieu of calling the VASD.]

VA Service Desk Alternate Contacts

- During business hours: As an alternate to the toll-free number, call 205-554-4710 through 205-554-4725, Monday through Friday (excluding holidays), 8:00 a.m. to 7:30 p.m. (Eastern Time).
- Outside business hours: Call 205-554-3459 through 205-554-3465, 205-554-3472, 205-554-3475, or 205-554-3482 through 205-554-3485).
- Web site: http://vawww.va.gov/emc/index.asp?s=6&p=nhd_home (VHA Enterprise Management Center).
- Email: vhacionhd@va.gov.

This page intentionally left blank.

Failover

VBECS does not have a seamless failover mechanism. If one server fails, the user will receive a message that the remote connection was lost. VBECS will lose information entered since the last save. The user must reopen a Remote Desktop Connection session. It may take 30 to 60 seconds for the Windows cluster and SQL Server running on it to fail over, which will open on the secondary server (without the user being aware of it). The user will have to reenter all information that was lost since the last save.

The connection between VBECS and VistA can be lost for a number of reasons:

- A server can fail in the VBECS cluster or the VistA server can fail. When this connection is lost, no messages can be exchanged. When the connection between VBECS and VistA is lost due to a failure of VBECS, the messages are queued on the VistA side. Orders placed during this downtime will remain in the queue. Once the VBECS system fails over and a connection is reestablished with VistA, the messages come across. The order alerts icon located in the VBECS status bar will display the orders that were in the queue at the time of failure.
- VBECS can fail because of a power outage. The UPS device will sound an alarm to alert the staff that the power is out. The IRM staff will inform the VBECS users to save their work and exit the system before the battery runs out.
- A server may fail because of a subcomponent failure such as a network interface card failure. MOM will monitor the servers for subcomponent failures. If a failure occurs, MOM will alert the IRM.

If only one node in a cluster is damaged, failover will occur. The IRM must check the MOM alerts for notification that the act occurred and fix the other node immediately to restore it to operation. When only one node is operating, no further failover can occur.

If a user's client workstation fails in the middle of a VBECS session, the session remains active on the server for a period set by the server administrator. The standard session time out is 15 minutes. If the user resolves the issues with the client workstation and reconnects to the VBECS server through Remote Desktop Connection before the session times out, the session will remain as it was when the client failed.

If a server fails due to a hardware issue, such as a network interface card failure, a Remedy ticket must be entered. If this failure occurs on only one node, users may continue to use the software after the system successfully fails over. The failover process will occur in 90 seconds. If both nodes in the cluster fail, file a Remedy ticket and refer VBECS users to Downtime Forms and Instructions in the *VistA Blood Establishment Computer Software (VBECS) User Guide*.

This page intentionally left blank.

Performance

VBECS may delay a critical function such as patient transfusion if the network suffers latency issues. File a Remedy ticket when latency issues arise.

VBECS was re-factored after performance testing results showed latency issues for VistA queries. As a result, many queries are cached in the VBECS database. Due to the criticality of having correct and current patient data, patient lookups cannot be cached.

Locking

VBECS is designed with pessimistic locking controlled within the application code: if one user selects a record for edit, the record is locked by that user. If another user tries to edit that record, a message will tell him that the record is locked and who has the record. The second user is not granted access to the record.

Locks have a timeout period defined in the configure division portion of the VBECS Administrator application. When a lock times out or is released by a user completing his edit, another user can edit that record.

If the application code fails due to a logic bug, optimistic locking is in place to prevent data corruption. When a record is retrieved, a row version is also retrieved. When a record is saved, the row in the database gets an updated row version; before the save takes place, the save routine checks that the row version supplied matches the row version in the table. If it does not match, the routine notifies the caller that another user changed the data. The save does not complete; the user must retrieve the updated record and start his edits again.

This page intentionally left blank.

Security

VBECs contains sensitive data and performs a critical function, so it is critical to secure the system. It is important to secure the server from both users and malicious attacks from an individual who is trying to gain access to the system. This information section describes the measures taken to secure VBECs.

Active Directory

Access to the VBECs servers is controlled through AD. Each VBECs site will have two groups set up in AD, one for normal VBECs users and one for VBECs Administrators (this is not a server administrator). Unless the user is a system administrator, he must be a member of one of these two groups to gain access to the server. Users will use their normal Windows user names to log in.

These groups also play a role in application level security. Even if a user were able to access the server, he would not be able to access VBECs.

Group Policy

Group policy controls the user experience (what the user sees and has access to on the VBECs server). To configure this correctly, the recommendations in “Locking Down Windows Server 2003 Terminal Server Sessions” and “Windows Server 2003 Security Guide” (Microsoft Web site) were followed to establish a baseline for group policy.

Group policy can be applied to user accounts or to the servers directly. In the case of VBECs, group policy is applied to the servers (it is easier to manage). It is also undesirable to have group policy associated with the user, which may inhibit his use of other systems. Enabling loopback processing applies the policy to any user that logs into the server.

Virtual Local Area Network

As a medical device, VBECs must exist in a segregated part of the LAN [Virtual Local Area Network (VLAN)]. The VLAN is configured to only allow necessary communication in and out of the VBECs system. Unneeded ports are blocked.

Microsoft Operations Manager

Microsoft Operations Manager (MOM) is a proactive monitoring tool. MOM will constantly monitor each server for system abnormalities. If MOM detects a problem, an email will be sent to the system administrator defined during the installation process. MOM will monitor these high-level categories:

- Windows Server 2003 Operating System
- CPU health and usage
- Network interface cards
- SQL Server
- Clustering
- Memory usage
- Hard disk health and usage
- VBECs executables and services
- Windows Services

Application-Wide Exceptions

Table 10 explains system exceptions to aid VA Product Support in determining the cause and resolving system issues.

Table 10: Application-Wide Exceptions

System Exceptions	Description
ArgumentException	Base class for all argument exceptions.
ArgumentNullException	Thrown by methods that do not allow an argument to be null.
ArgumentOutOfRangeException	Thrown by methods that verify that arguments are in a given range.
ComException	Exception encapsulating COM HRESULT information.
Exception	Base class for all exceptions.
ExternalException	Base class for exceptions that occur or are targeted at environments outside the runtime.
IndexOutOfRangeException	Thrown by the runtime only when an array is indexed improperly.
InvalidOperationException	Thrown by methods when in an invalid state.
NullReferenceException	Thrown by the runtime only when a null object is referenced.
SEHException	Exception encapsulating Win32 structured exception handling information.
System.ArithmeticException	A base class for exceptions that occur during arithmetic operations, such as System.DivideByZeroException and System.OverflowException.
System.ArrayTypeMismatchException	Thrown when a store into an array fails because the actual type of the stored element is incompatible with the actual type of the array.
System.DivideByZeroException	Thrown when an attempt to divide an integral value by zero occurs.
System.IndexOutOfRangeException	Thrown when an attempt to index an array via an index that is less than zero or outside the bounds of the array.
System.InvalidCastException	Thrown when an explicit conversion from a base type or interface to a derived type fails at run time.
System.NullReferenceException	Thrown when a null reference is used in a way that causes the referenced object to be required.
System.OutOfMemoryException	Thrown when an attempt to allocate memory (via new) fails.
System.OverflowException	Thrown when an arithmetic operation in a checked context overflows.
System.StackOverflowException	Thrown when the execution stack is exhausted by having too many pending method calls; typically indicative of very deep or unbounded recursion.
System.TypeInitializationException	Thrown when a static constructor throws an exception, and no catch clauses exist to catch it.
SystemException	Base class for all runtime-generated errors.

Glossary

Acronym, Term	Definition
ABO	A group for classifying human blood, based on the presence or absence of specific antigens in the blood, which contains four blood types: A, B, AB, and O. The ABO group is the most critical of the human blood systems. It is used to determine general compatibility of donor units to a recipient.
ABS	Antibody screen, antibody screen test.
Access Code	A field in the VistA New Person file used to uniquely identify a user on the VistA system.
Active Directory	A hierarchical directory service built on the Internet's Domain Naming System (DNS).
API	Application Programmer Interface.
CPRS	Computerized Patient Record System.
DBIA	Database Integration Agreement.
DSS	Decision Support System.
HCPCS	Healthcare Common Procedure Coding System.
HL7	Health Level Seven.
ICN	Integration Control Number.
LLP	Lower Layer Protocol.
LMIP	Laboratory Management Index Program.
MLLP	Minimal Lower Layer Protocol.
MOM	Microsoft Operations Manager.
OSI	Open Systems Interconnect.
OU	Organizational Unit.
PCE	Patient Care Encounter.
RDP	Remote Desktop Protocol.
RPC	Remote procedure call.
TCP/IP	Transmission Control Protocol/Internet Protocol.
UPS	Uninterruptible power source.
VAISS	VBECS Application Interfacing Support Software.
VBECS	VistA Blood Establishment Computer Software.
VDL	VistA Documentation Library.
Verify Code	A field in the VistA New Person file used to verify the identity of a user associated with an Access Code.
VISN	Veterans Integrated Service Network.
VLAN	Virtual Local Area Network.
XML	Extensible Markup Language.

This page intentionally left blank.

Appendices

Appendix A: Instructions for Capturing Screen Shots


Throughout the technical manual-security guide, the Administrator is asked to capture screen shots to document configuration options. To capture a screen shot:

- 1) Open a blank document (for example, in Microsoft Word) and save it as (click **File, Save As**) “mmyydd Technical-Security Validation Record,” or another easily identified file name.

 *If you wish to place a document on both servers for ease of copying and pasting, assign file names similar to “mmyydd Technical-Security Validation Record Server1” and “mmyydd Technical-Security Validation Record Server2.”*


- 2) When the screen you wish to capture is displayed, press the **Print Screen** key.
- 3) In the Technical-Security Validation Record document, place the cursor where you want to insert the picture.
- 4) Click  (the paste icon) or select **Edit, Paste** (Figure 145).

Figure 145: Paste



- 5) Label the screen shot within the document with the technical manual-security guide step, page number, and server on which the picture was taken (Figure 146).

Figure 146: Screen Shot



This page intentionally left blank.

Appendix B: Workload Process Mapping to Application Option Table

Table 11 associates record saves with workload processes. The data fields identified for transmission at the completion of a Workload Event are based on current VistA workload-related files and fields. VBECS will transmit information to a new flat file. There are no donor workload types in VBECS.

Table 11: Workload Process Mapping to Application Option

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Record a Transfusion Reaction Workup	ABO Forward and reverse typing (patient)	P	An ABO/Rh test for “pre” or “post” is enabled and a valid interpretation other than Not Tested is selected. A workload event is accrued separately for “Pre” and “Post” entries.
Record Patient ABO/Rh		P	Accrue workload when a CPRS-ordered ABO/Rh test is performed.
Invalidate Patient Test Results*		P	Accrue workload when a completed ABO/Rh test is invalidated.
Record Patient ABO/Rh	ABO Forward and reverse typing (patient) Repeat Test	M	Accrue workload when a reflex or repeat ABO/Rh test is performed, completed, and saved.
Invalidate Patient Test Results*		M	Accrue workload when a reflex or repeat ABO/Rh test is invalidated.
ABO/Rh Confirmation	ABO forward typing (unit)	U	An ABO confirmation test is performed. When multiple units are selected in a batch, each unit in the batch accrues a workload event. Note: Workload generated during Anti-D testing is not included in the unit’s confirmation test. Workload is not accrued when an ABO or Rh discrepancy override is processed and VBECS releases all patient assignments. Workload is not accrued when VBECS quarantines the unit due to a discrepancy. There is no special handling for workload collection for additional confirmation tests on a unit.
Edit Unit Information*		U	Accrue workload when an ABO confirmation test is invalidated.
ABO/Rh Confirmation	ABO/Rh forward typing (unit)	U	An ABO/Rh confirmation test is performed. When multiple units are selected in a batch, each unit in the batch accrues a workload event. Note: Workload generated during Anti-D is part of the unit’s confirmation test. Workload is not accrued when an ABO or Rh discrepancy override is processed and VBECS releases all patient assignments. Workload is not accrued when VBECS quarantines the unit due to a discrepancy. Any unit successfully confirmed accrues workload. For split modifications: workload is not inherited by split units. A split unit that requires confirmation accrues confirmation workload at the time of testing. There is no special handling for workload collection for additional confirmation tests on a unit.
Edit Unit Information*		U	Accrue workload when an ABO/Rh confirmation test is invalidated.

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Accept Order	Accept Order	M	Accrue workload when an order is accepted. When a multiple orders are selected, each order accrues workload.
Enter Antibody Identification Results	Antibody identification Work-Up	P	User enters additional workload associated with the individual reflex-ordered ABID. The selected VBECS multiplier will multiply against the VistA multiplier and display the (multiplication) product on the Division Workload Report.
Invalidate Patient Test Results*		P	Accrue workload when the ABID is invalidated.
Record a Transfusion Reaction Workup	Antibody Screen (patient)	P	An ABS test for "pre" or "post" is enabled and a valid interpretation other than Not Tested is selected. A workload event is accrued separately for "Pre" and "Post" entries.
Record Patient Antibody Screen		P	Accrue workload when an ordered ABS test is performed.
Invalidate Patient Test Results*		P	Accrue workload when a completed ABS test is invalidated.
Record Patient Antibody Screen	Antibody Screen (patient) Repeat Test	M	Accrue workload when a reflex or repeat ABS test is performed, completed, and saved.
Invalidate Patient Test Results*		M	Accrue workload when a reflex or repeat ABS test is invalidated.
Unit Antigen Typing / Patient Antigen Typing	Antigen phenotyping, Single Test phase (QC)	M	Accrue workload when Antiserum QC in Unit or Patient Antigen Typing includes the testing of both the positive and negative control cells, per specificity by lot number, when only one phase of reactivity is chosen for the test grid (IS or AHG). One workload event is collected per completed tab for regular or repeat antigen tests.
Unit Antigen Typing / Patient Antigen Typing	Antigen phenotyping, Multiple Test phases (QC)	M	Accrue workload when Antiserum QC in Unit or Patient Antigen Typing includes the testing of both the positive and negative control cells, per specificity by lot number, when only multiple phases of reactivity are chosen for the test grid, IS/RT, RT/37, or weak D. One workload event is collected per completed tab for regular or repeat antigen tests. When weak D is the selected test, QC may not be accrued for the rack selection. QC is accrued when positive and negative cells must be tested for the lot number.
Cancel Pending Order	Cancel Order	M	Accrue workload when an order on the pending order list is canceled. When multiple orders are canceled, each order accrues workload.
Cancel Active Order	Cancel Order	M	Accrue workload when an order on the pending task list is canceled. When multiple orders are canceled, each order accrues workload.
Select Units for Crossmatch	Crossmatch unit, electronic	P	This process is invoked when an individual unit is selected for patient assignment and the unit is electronically crossmatched. When multiple units are selected, each unit accrues workload.
Enter Crossmatch Results	Crossmatch unit, serologic immediate spin	P	Accrue workload when an individual unit crossmatch is selected to include only the IS phase, is completed, and is saved. When multiple units are selected, each unit accrues workload.

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Invalidate Patient Test Results*		P	Accrue workload when a completed crossmatch test is invalidated. This applies to the workload originally saved with the serologic immediate spin test.
Record a Transfusion Reaction Workup	Crossmatch unit, serological Coombs	P	A crossmatch test for “pre” or “post” is enabled and a valid interpretation other than Not Tested is selected. A workload event is accrued separately for “Pre” and “Post” entries. When multiple units are selected, each unit accrues workload.
Enter Crossmatch Results		P	Accrue workload when an individual unit crossmatch is selected to include all phases or only the AHG phase, is completed, and is saved. When multiple units are selected, each unit accrues workload.
Invalidate Patient Test Results*		P	Accrue workload when a completed crossmatch test is invalidated. This applies to the workload originally saved with the test, serological Coombs.
Enter Crossmatch Results	Crossmatch, Repeat Test	M	Accrue workload when an individual unit crossmatch is selected to include all phases or IS or only the AHG phase, is completed, and is saved. When multiple units are selected, each unit accrues workload.
Invalidate Patient Test Results*		M	Accrue workload when an individual unit crossmatch is invalidated.
Enter Daily QC Results	Daily Rack Quality Control (QC)	M	Accrue workload when Daily QC rack completed for one individual rack includes all rows in configured QC. When multiple racks are tested, each completed and saved tab accrues a workload event.
Record Patient Direct Antiglobulin Test	DAT (QC)	M	Accrue workload when Reagent QC completed in Patient DAT testing includes the testing of both the positive and negative control cells, per specificity per lot number, when only one phase of reactivity is chosen for the test grid (IS or AHG). One workload event is collected per completed tab for regular or repeat antiglobulin tests (PS, IgG, Comp).
Record a Transfusion Reaction Workup	Direct Antiglobulin Test (DAT)	P	A DAT test for “pre” or “post” is enabled and a valid interpretation other than Not Tested is selected. A workload event is accrued separately for “Pre” and “Post” entries.
Record Patient Direct Antiglobulin Test		P	Accrue workload when a DAT is completed and saved. This count is used for all antiglobulin tests (PS, IgG, Comp) when ordered from CPRS or Reflex testing.
Invalidate Patient Test Results*		P	Accrue workload when a completed DAT, PS, IgG, or Comp is invalidated.
Record Patient Direct Antiglobulin Test	Direct Antiglobulin Test (DAT)	M	Accrue workload when a reflex or repeat DAT test is performed, completed, and saved. This applies to all repeat antiglobulin tests (PS, IgG, Comp).
Invalidate Patient Test Results*	Repeat test	M	Accrue workload when a completed Repeat DAT, PS, IgG, or Comp is invalidated.
Modify Units	Deglycerolize unit	U	Accrue workload when an individual blood unit is processed by the Deglycerolize modification type. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Deglycerolize."
Discard or Quarantine Unit	Discard unit	U	Accrue workload when an individual blood unit's status is invalidated. When a batch of units is selected, each unit accrues workload.
Remove Final Status*		U	Accrue workload when a unit is discarded for waste or credit. When a batch of units is selected, each unit accrues workload.
Modify Units	Freeze unit	U	Accrue workload when an individual blood unit is processed by the Freeze modification type. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Freeze."
Modify Units	Irradiate unit	U	Accrue workload when an individual blood unit is processed by the Irradiate modification type. When a batch of units is irradiated, each unit accrues workload. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Irradiate."
Modify Units	Leukoreduce unit	U	Accrue workload when an individual blood unit is processed by the Leukoreduce modification type. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Leukoreduce."
Split a Unit	Split unit	U	Accrue workload when a unit modification of Split and a single workload event is recorded regardless of the number of units created by the modification. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	A Split Unit has its unit status invalidated. A single workload event is recorded regardless of the number of units originally created by the modification.
Modify Units	Rejuvenate unit	U	Accrue workload when an individual blood unit is processed by the Rejuvenate modification type. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Rejuvenate."

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Modify Units	Thaw	U	Accrue workload when an individual blood unit is processed by the Thaw modification type. When a batch of units is thawed, each unit accrues workload. This applies to Thaw FFP and Thaw Cryo. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Thaw." This modification type is applicable to Thaw FFP and Thaw Cryo.
Modify Units	Wash unit	U	Accrue workload when an individual blood unit is processed by the Wash modification type. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was "Wash."
Modify Units	Volume Reduce	U	Accrue workload when an individual blood unit is processed by the Volume Reduce modification type. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Remove Final Status*		U	An individual blood unit's status is invalidated when the original modification process was Volume Reduce.
Issue Blood Components	Issue unit	P	Accrue workload when a unit is issued to a patient. When a batch of units is processed, each unit invokes one workload process.
Justify Patient ABO/Rh Change	Justification	M	Workload is accrued when a patient's ABO or Rh typing is justified. One workload event is accrued per patient justification.
Login Equipment	Login Equipment	M	Accrue workload when a lot number of any type of equipment is logged into the system. When multiple lot numbers are processed in a batch, each lot number's workload is counted.
Login Reagent	Login Reagent	M	Accrue workload when a lot number of any type of reagent is logged into the system. When multiple lot numbers are processed in a batch, each lot number's workload is counted.
Login Supply	Login Supply	M	Accrue workload when a lot number of any type of supply is logged into the system. When multiple lot numbers are processed in a batch, each lot number's workload is counted.
Maintain Specimen	Maintain Specimen	M	Accrue workload when a specimen is maintained during order acceptance and is required for acceptance of the order. Note: This is collected in addition to the accept order workload accrued by accepting an order. Marking a specimen unacceptable does not create a negative workload event.

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Patient antigen phenotype	Patient antigen phenotype (multiple phases)	P	Accrue workload when a patient antigen phenotype test with IS/RT or IS/37 phases is completed and saved. One workload event is collected per completed tab for repeat or regular antigen tests.
Invalidate Patient Test Results*		P	Accrue workload when a patient antigen phenotype test as defined by the antiserum specificity tested with any phases is invalidated.
Patient antigen phenotype	Patient antigen phenotype (single phase)	P	Accrue workload when a patient antigen phenotype test with AHG or IS phase is completed and saved. One workload event is collected per completed tab for repeat or regular antigen tests.
Invalidate Patient Test Results*		P	Accrue workload when a patient antigen phenotype test as defined by the antiserum specificity tested with a single phases is invalidated.
Pool Units	Pool unit	U	Accrue workload when a pooled unit is created and a single workload event is recorded regardless of the number of units included in the pooled unit. This applies to the Pool modification type. Add/Remove unit from a pool does not accrue any workload. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Edit Unit Information*		U	Accrue workload when a unit is inactivated if the pooled unit was created in VBECS.
Remove Final Status		N/A	No effect on workload accrual when a unit is removed from a modified status that was included in a pool.
Discard or Quarantine Unit	Quarantine unit	U	Accrue workload when a unit is marked for quarantine. When a batch of units is selected, each unit accrues workload.
Free Directed Unit For Crossover	Release directed unit	U	Accrue workload when an individual blood unit with the restriction type of "directed" is released for use as an allogeneic unit.
Release Unit from Patient Assignment	Release unit from patient back to inventory	U	Accrue workload when an individual unit is released from patient assignment. When multiple units are selected, each unit accrues workload.
Discard or Quarantine Unit	Release unit from Quarantine	U	Accrue workload when a unit is released from quarantine. When a batch of units is selected, each unit accrues workload.
Return Issued Unit	Return Issued unit	U	Accrue workload when a unit is returned from issue status.
Modify Units	Thaw/pool Cryo	U	Accrue workload when an individual unit has a modification of Thaw/Pool Cryo. A single workload event is recorded regardless of the number of units included in the pooled unit. Note: Workload is not accrued when a patient assignment is processed and VBECS releases all other patient assignments. Workload is not accrued when VBECS is required to quarantine the unit.
Edit Unit Information*		U	Accrue workload when a unit is inactivated (unit record inactivated) when the pooled unit was created in VBECS.

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Remove Final Status		N/A	There is no effect on workload accrual when a unit is removed from a modified status that was included in a Thaw/pool Cryo pool.
Enter Post-Transfusion Details	Transfuse Unit	U	Accrue workload when an individual blood unit's status is updated to "transfused."
Remove Final Status*		U	An individual blood unit's status is invalidated when the unit was in a status of "transfused."
Record a Transfusion Reaction Workup	Transfusion Reaction Investigation	P	Accrue workload when a transfusion reaction investigation is saved. This does not include workload accrued by the optional TRW serologic testing.
Invalidate Patient Test Results*		P	Accrue workload when a transfusion reaction investigation previously saved is invalidated.
Unit Antigen Typing	Unit Antigen phenotyping, Multiple Test phases	U	Accrue workload when a unit antigen phenotype test with IS/RT or IS/37 phases is selected and completed for an individual blood unit. There is no special handling for workload collection for additional repeat antigen typing tests on a unit.
Edit Unit Information*		U	Accrue workload when a unit antigen phenotype test with Multiple Test phases is invalidated for an individual blood unit.
Unit Antigen Typing	Unit Antigen phenotyping, Single Test phase	U	A unit antigen phenotype test with AHG or IS phase is selected and completed for an individual blood unit. There is no special handling for workload collection for additional repeat antigen typing tests on a unit.
Edit Unit Information*		U	Accrue workload when a unit antigen phenotype test with Single Test phase is invalidated for an individual blood unit.
Incoming Shipment	Unit login	U	An individual unit record is activated as "saved" to an incoming shipment invoice. When multiple units are entered, each unit added to the database accrues workload.
Edit Unit Information*		U	Accrue workload when a unit is inactivated and logged in through incoming shipment or is a pooled unit created in VBECS. When the unit was created by split modification, no workload is invalidated in this option.
Outgoing Shipment	Unit logout	U	An individual unit's status is updated to "transferred" on a confirmed outgoing shipment invoice. When multiple units are selected, each unit accrues workload. Accrue workload on confirmation of the invoice, not the addition of a unit to a temporary outgoing shipment invoice: an invoice may be confirmed only once.
Remove Final Status*		U	An individual unit status is invalidated when the unit had a previous unit status of "transferred."
Update Equipment Record	Update Equipment Record	M	Accrue workload when a lot number of any type of equipment is updated in the system.
Update Reagent Inventory	Update Reagent Inventory	M	Accrue workload when a lot number of any type of reagent is updated in the system. When multiple lot numbers are processed in a batch, each lot number's workload is counted.

Record Save Option	VBECS Process	Transaction Type [P (Patient), U (Unit), M (Miscellaneous)]	Explanation
Update Supply Inventory	Update Supply Inventory	M	Accrue workload when a lot number of any type of supply is updated in the system. When multiple lot numbers are processed in a batch, each lot number's workload is counted.

*Accumulates negative workload when it is associated with inactivation of a unit or removal of a final status.

Appendix C: Known Defects and Anomalies

Copies of *Known Defects and Anomalies* are available at the VDL: VistA Documentation Library (VDL), VHA OI – Health Systems Design & Development Web page.

This page intentionally left blank.

Appendix D: Active Directory Request Form

Fill out this form and email or fax it to your data center contact to have users added or deleted from the VBECS Active Directory groups. Email or fax it to your data center contact for action. Contact the Implementation Team to verify your data center contact, if necessary. The data center administrator facilitating this request will return this form to you when the changes are completed.

Blood bank information

Site Name:	
Site identifier:	VISN number:
Contact name:	Phone number:
Email:	Fax Number:

Data Center information

Technician name:	Phone number:
Email:	Fax number:

VBECS Users (users of normal VBECS): RnnxxxVbecsUsers group (nn is data center identifier and xxx is site identifier)

Specify the action, name and Windows ID of each technician requiring a change in access. The data center administrator will fill in his/her initials in the last column to confirm the change.

Row	Action	Last name, first name	Windows ID	Initials (for data center administrator only)
1	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
2	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
3	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
4	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
5	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
6	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
7	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
8	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
9	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
10	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
11	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
12	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
13	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
14	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
15	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
16	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
17	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
18	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
19	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
20	Add <input type="checkbox"/> Delete <input type="checkbox"/>			

**VBECS Administrators (users of administrative unit of VBECS):
RnnxxxVbecsAdministrators group (nn is data center identifier and xxx is site identifier)**

Specify the action, name and Windows ID of each technician requiring a change in access. The data center administrator will fill in his/her initials in the last column to confirm the change.

Row	Action	Last name, first name	Windows ID	Initials (for data center administrator only)
1	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
2	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
3	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
4	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
5	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
6	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
7	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
8	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
9	Add <input type="checkbox"/> Delete <input type="checkbox"/>			
10	Add <input type="checkbox"/> Delete <input type="checkbox"/>			

Appendix E: Data Center Instructions

Purpose

This appendix describes the tasks that must be completed by the data center for a successful VBECS installation, and is divided into 3 main sections depending on when the activities take place:

- Initial Setup Tasks: These tasks must be completed prior to installation of any VBECS systems.
- Ongoing Tasks: These are continual maintenance tasks.
- Installation Time Tasks: These tasks are to be completed at the time of a VBECS installation.

Initial Setup Tasks

Execute these tasks once, prior to setting up the VBECS systems in the data center.

Active Directory

VBECS User and Server Administrator Requirements

VBECS depends on Active Directory for remote server access for both VBECS and administration.

Set up two groups set up in Active Directory. The groups must have a “Universal” scope and a “Security” type.

- *RmxxxVbecsUsers* (replace *nn* with your two-digit region number and *xxx* with the site location code): These are normal users of the VBECS system. Members of this group will have access to the server and are allowed to launch the VBECS application.
- *RmxxxVbecsAdministrators* (replace *nn* with your two-digit region number and *xxx* with the site location code): These are users who must access the administrative component of VBECS. Members of this group will have access to the server and are allowed to launch the VBECS Administrator application.

Create a server administrator group to be shared across servers. This group must have a “Universal” scope and a “Security” type. This group will have administrative access to the VBECS servers at installation:

- *RxxVbecsServerAdmins* (replace *xx* with your two-digit region number): These are traditional server administrators who need full administrative privileges to the system. For MOM support, add the VA IT Engineering CIS Monitoring Team group to this administrator group.

VBECS Server Requirements

For Group Policy purposes, VBECS servers will reside in their own OU, which will contain only VBECS servers. You may also create OUs under the main OU for organizational purposes. For more information, see the Group Policy section.

Group Policy


Import the VHA VBECS Terminal Server Policy from the VHAMASTER domain. If the VBECS development team changes the policy, import it again.

When importing the policy, clear the VBECS Windows Software Update Services settings (see Computer Configuration/Administrative Templates/Windows Components/Windows Update).

Place the group policy in the top-level server OU. For more information about OUs and server organization, see the Active Directory section.

Configure the policy so that it is not applied to the RxxVbecsServerAdmins Active Directory group. See the example in Figure 147.

Figure 147: Example of a Group Policy Not Applied to VBECSAdministrators Group



Service Accounts

VBECS requires dedicated service accounts for Microsoft Cluster and Microsoft SQL Server. Add these accounts to your RxxVbecsServerAdmins group. Define these service accounts once to be shared across VBECS servers (xx represents the two-digit region number):

- Microsoft Cluster: RxxVBESVCCLU01
- Microsoft SQL Server: RxxVBESVCSQL01

At installation, give the passwords for these accounts to the installer.

Terminal Server License Server

VBECS is a Terminal Server application and requires a license. Ensure that there is at least one Terminal Server License server set up for your domain.

VLAN

Since VBECS is a medical device, VBECS servers and printers must reside in a VLAN. Do not turn on the VLAN until installation is complete. Since this is a data center installation, the servers will reside on a VLAN separate from that of the printers, which reside at the blood bank.


Table 12 details the communication requirements for the VLAN. Figure 148 depicts how VBECS resides in the network.

Table 12: VBECS Communication Requirements

Servers, Workstations, Printers	IP Address	Notes
Data center domain controllers (all), WINS, DNS	See data center network administrator	
Data center WSUS	See data center network administrator	
VHAMASTER WINS	10.3.29.33 10.3.29.34 10.39.129.200	
va.gov domain controllers	10.3.21.197 10.3.30.1 10.2.21.27 10.204.1.10 10.3.21.193	
med.va.gov domain controllers	10.2.21.26 10.4.229.41 10.3.30.2 10.3.21.194 10.30.20.27	
VHAMASTER (vha.med.va.gov) domain controllers	10.4.229.2 10.5.21.192 10.4.21.192 10.6.21.192 10.2.21.192 10.191.10.7 10.6.197.13 10.1.21.192 10.3.27.33 10.3.21.192 10.224.151.15 10.190.40.20 10.40.198.21 10.63.196.55 10.15.32.250 10.189.77.230 10.222.228.3 10.189.110.3 10.61.192.172 10.61.192.139 10.189.1.1 10.224.151.90 10.208.13.14 10.3.30.25 10.189.37.217 10.189.46.203	
VISN WINS	See VISN network administrator	
VISN domain controllers	See VISN network administrator	Due to DNS "round robinning," all local domain controllers must be accessible.

Servers, Workstations, Printers	IP Address	Notes
VistA	See your network administrator	
MOM	10.3.31.51 10.3.31.52	
ePolicy	10.204.9.190 10.254.36.43 10.254.36.45	
SMTP support	10.2.27.92 10.3.27.92 10.208.13.3 10.6.27.92 10.252.92.14 10.252.92.15 10.252.93.14 10.252.93.15 10.252.94.14 10.252.94.15 10.252.95.14 10.252.95.15	
VBECs workstations	See Appendix: Blood Bank Configuration Checklist (installation guide)	
VBECs printers (label and report)	See Appendix: Blood Bank Configuration Checklist (installation guide)	If the printers reside at the same location as the servers, just place them in the same VLAN.

Figure 148: VLAN Schematic



Ongoing Tasks

Execute the tasks in this section continually.

Back Up the VBECS Database

Back up the VBECS database nightly:

- Back up all folders and files in the \\<cluster name>\d\$\Program Files\Microsoft SQL Server\MSSQL\BACKUP directory.
- Maintain backups for at least seven days.

VBECS Updates

When the VBECS development team releases a VBECS patch, install the patch in accordance with instructions supplied by the development team.

Windows Updates

The VBECS development team must test every Microsoft Windows update. Once the development team is satisfied that the update causes no adverse effects, a Vista information patch in the VBEC (yes VBECS) namespace will be created by the VBECS. This patch will describe where to obtain the update and how to apply it. The patch will be released to customers by VA Product Support.

Installation of patches needs to be coordinated with the blood bank manager since most updates require a reboot.

Installation Time Tasks

Complete the Checklists and Password List

Complete these checklists and password list in the *Vista Blood Establishment Computer Software (VBECS) Installation Guide* prior to installation:

- Appendix B: Blood Bank Hardware Checklist: This checklist helps ensure that the correct server hardware is on-site.
- Appendix E: Server Configuration Checklist: This checklist contains server details such as names and IP addresses.
- Appendix H: Password List: This list includes passwords for the cluster and SQL server user IDs.

Update the VBECS Server Administrators Group

Refer to the appendices in the *Vista Blood Establishment Computer Software (VBECS) Installation Guide* to complete the installation of VBECS:

- Add the installers to the VBECS Server Administrators (RxxVbecsServerAdmins) group. See the Windows IDs of VBECS Installers cell in the Contact Information table of the Server Configuration Checklist (Appendix E). Upon successful completion, delete the installers from the group.
- Add the executor of the VBECS data conversion to the VBECS Server Administrators group. See the Data Conversion section of the Blood Bank Configuration Checklist (Appendix).


Appendix F: Database Conversion Updates

Changes to the VBECS 1.5.0.0 database required Data Transformation Services (DTS) package changes. Sites performing database conversion must follow these steps to complete conversion on the VBECS server:

The DTS package confirms the target VBECS database is not in use as configured by the *Maintenance Operation* section of the *VBECS Technical Manual-Security Guide*. The DTS package also checks if a completed conversion exists in the VBECS database prior to inserting the converted records in the VBECS database. If either condition is true, the conversion process terminates.

- 1) To prepare for the execution of the DTS package:
Log onto an account on the VBECS system with Administrator privileges.
Click **Start** and **Run** from the Windows taskbar.
Enter **cmd** in the Run Window. Click **OK**.
- 2) To run the conversion:
Enter **cd c:\dbconv\dts** at the command prompt.
Enter **dtsrun /f conversionpackagemultidb.dts** at the command prompt. Press the **Enter** key.
- 3) After the DTS package verifies the VBECS database is able to be converted a message displays (Figure 149) to the user to confirm the conversion settings where:
 - a) **SQL SERVER** – Displays the value supplied in the **SERVER NAME** field of the **DBCONV.INI** settings file (e.g. VHAxxxSQLZ1). This is the target SQL Server name where conversion will occur.
 - b) **DATABASE** – Displays the value supplied in the **DATABASE NAME** field of the **DBCONV.INI** settings file (e.g. VBECS_V1_PROD). This is the target database name where conversion will occur

Figure 149: Example of DTS Run Message



- 4) Verify the Server and Database names are correct and click **Yes** to proceed with conversion, if the settings are incorrect click **No** to stop the conversion.
Note: If you reply **No**, repeat the *CONV Utilities Used for the Database Conversion* section of *Blood Bank Pre-Implementation Data Validation, Mapping, and Conversion LR*5.2*335 ADPAC Guide* to correct the Server and Database name. FTP the files back to the VBECS server when conversion is complete.

5) The DTS will run and display (Figure 150) when complete.

Figure 150: Example of Congratulation Message




Warnings and Notifications Displayed by the DTS Package

The target VBECs database already configured for use:

After the user presses **Enter** to start the DTS package target database is checked to see if it is in use. An error message displays to the user if the database is in use and the conversion process terminates (Figure 151).

Figure 151: Database in Use



- SQL SERVER – Displays the value supplied in the **SERVER NAME** field of the DBCONV.INI settings file (e.g. VHAxxxSQLZ1). This is the target SQL Server name for conversion.
- DATABASE – Displays the value supplied in the **DATABASE NAME** field of the DBCONV.INI settings file (e.g. VBECs_V1_PROD). This is the target database name that contains existing VBECs data.

Conversion already executed against target database:

After the user presses **Enter** to start the DTS package checks the target database to see if a database conversion has already successfully completed. If a completed conversion is detected, a message (Figure 152) is displayed and the conversion process terminates.

Figure 152: Conversion Already Run Message



- SQL SERVER – Displays the value supplied in the SERVER NAME field of the DBCONV.INI settings file (e.g. VHAxxxSQLZ1). This is the target SQL Server name.
- DATABASE – Displays the value supplied in the DATABASE NAME field of the DBCONV.INI settings file (e.g. VBECS_V1_PROD). This is the target database name where conversion has been completed.


Patient name length error:

When the conversion process is running, patient names supplied from VistA are evaluated. The conversion process terminates and a failure message (Figure 153) appears if any of the following conditions occur:

- The length of PatientLastName and PatientFirstName is greater than 29.
- The length of PatientLastName, PatientMiddleName and PatientFirstName is greater than 28.

The Patient record will need to be updated on VistA and the files retransmitted to the VBECS server, at which point the DTS conversion can be started again.

Figure 153: Example of Patient Name Length Failure




DTS set-up error:

A failure message is displayed and the conversion process terminates (Figure 154) if one of these conditions occur:

- Required conversion text files are not found.
- DBCONV.INI file contains settings pointing to an invalid SQL Server name.
- DBCONV.INI file contains settings pointing to an invalid database name.
- The user executing conversion does not have sufficient privileges for the database.

Figure 154: DTS Failure Message



Appendix G: Services Allowed to run on VBECS Servers

The following services are permitted to run on VBECS servers.

- Application Experience Lookup Service
- Automatic Updates
- Cluster Service
- COM+ Event System
- Computer Browser
- Cryptographic Services
- DCOM Server Process Launcher
- DHCP Client
- Distributed Link Tracking Client
- Distributed Transaction Coordinator
- DNS Client
- Error Reporting Service
- Event Log
- FTP Publishing Service
- HID Input Service
- HP Insight Notifier
- HP Insight Foundation Agents
- HP Insight NIC Agents
- HP Insight Server Agents
- HP Insight Storage Agents
- HP Proliant Remote Monitor Service
- HP Smart Array SAS/SATA Event Notificaion Service
- HP System Management Homepage
- HP Version Control Agent
- HTTP SSL
- IIS Admin Service
- IPSEC Services
- Logical Disk Manager
- McAfee Framework Service
- McAfee McShield

- McAfee Task Manager
- Microsoft Search
- MOM
- MSSQLSERVER
- Net Logon
- Network Connections
- Network Location Awareness (NLA)
- NT LM Security Support Provider
- Plug and Play
- Pml Driver HPZ12
- Print Spooler
- Protected Storage
- Remote Access Connection Manager
- Remote Procedure Call (RPC)
- Remote Registry
- Secondary Logon
- Security Accounts Manager
- Server
- Shell Hardware Detection
- Simple Mail Transfer Protocol (SMTP)
- SQLSERVERAGENT
- System Event Notification
- Task Scheduler
- TCP/IP NetBIOS Helper
- Terminal Services
- VBECS Services (enabled services depends on site configuration)
- Windows Management Instrumentation
- Windows Time
- Workstation
- World Wide Web Publishing Service

Appendix H: Auditing on VBECS Servers

The following events are audited on VBECS servers. These events may be viewed in Event Viewer logs (under **Administrative Tools**).

- Account logon events (Success, Failure)
- Account management (Success, Failure)
- Directory service access (Success, Failure)
- Logon events (Success, Failure)
- Object access (Success, Failure)
- Policy Change (Success, Failure)
- System events (Success, Failure)

This page intentionally left blank.

Index

A

Active Directory	139, 157
Active Directory Request Form	155
Additional Required Hardware	24
Appendices	143
Application-Wide Exceptions	140
Archiving and Recovery	125

B

Back Up the VBECS Database	161
----------------------------------	-----

C

Commonly Used System Rules	29
Complete the Checklists and Password List	162
Configure Interfaces	56
Configure System Administrators	73
Connection Speed	8
Create a Remote Desktop Connection Shortcut for VBECS	10

D

Data Center Instructions	157
Database Conversion Updates	163, 167, 169
Delete Patch Files	
Prod	97

E

ePolicy and Virus Definitions	29
External Interfaces	87

F

Failover	135
Firmware Updates	30

G

Glossary	141
Group Policy	139, 157

H

Hardware and Backup Exec Alerts	30
Hardware and System Configuration	11
Hardware Specifications and Settings	5
Health Level Seven Interfaces	87
How This Technical Manual-Security Guide Is Organized	3

I

Implementation and Maintenance	25
Initial Setup Tasks	157
Installation Time Tasks.....	162
Instructions for Capturing Screen Shots	143
Integrated Lights Out	40
Introduction	1

K

Known Defects and Anomalies	153
-----------------------------------	-----

L

Locking.....	137
--------------	-----

M

Maintenance Operations	53
Microsoft Operations Manager	139

N

Notify VBECS Central Administrator	86
--	----

O

Off-the-Shelf Software Requirements.....	24
Ongoing Tasks	161

P

Performance.....	137
Periodic Maintenance Checks.....	25
Printers.....	12
Purpose	157

R

Reconfiguring the VBECS HL7 Multi Listener Service and VistALink	92
Related Manuals and Reference Materials.....	1
Remote Desktop Configuration	5

S

Save Settings.....	9
Scanners.....	20
Screen Resolution	5
Screen Shots	3
Security.....	139
Server and shared array disks	11
Server Configuration	23
Service Accounts	158
Sound.....	7
System Shut down Instructions.....	50

T

Terminal Server License Server.....	158
Transmit Workload Data	85

U

Update the VBECS Server Administrators Group.....	162
---	-----

V

VBECS Backup	125
VBECS Recovery	125
VBECS Updates	162
VBECS Windows Services.....	iii, 91
Virtual Local Area Network	139
VistALink Remote Procedure Calls.....	89
VLAN	159

W

Warnings and Notifications Displayed by the DTS Package	164
Windows Updates	28, 162
Workload Process Mapping to Application Option Table	145
Workstation Configuration	24

This is the last page of the *VistA Blood Establishment Computer Software (VBECS) 1.5.0.0 Technical Manual-Security Guide*.