

September 2006

This distribution contains change pages for patch MD*1.0*4 of the Clinical Procedures 1.0 User Manual.

The following documentation change pages should be inserted before these replacement pages:

File Name:

MD_1_P2_UM.PDF

MD_1_P10_UM.PDF

Patch:

MD*1.0*2

MD*1.0*10

Patch MD*1.0*4 pages:

Replace Pages:

Revision History

3-7 to 3-16

With Pages:

Revision History

3-7 to 3-16

Revision History

Description	Date	Technical Writer
Originally released.	April 2004	
¹ Patch MD*1.0*2 released.	July 2004	
² Patch MD*1.0*10 released.	March 2005	
³ Patch MD*1.0*4 released.	September 2006	Alfred Bustamante

¹ Patch MD*1.0*2 July 2004 Patch 2 release added.

² Patch MD*1.0*10 March 2005 Patch 10 release added.

³ Patch MD*1.0*4 September 2006 Patch 4 release added.

Vista CPRS in use by: CPUSER, TWO (DHCPSEVER1-9100-DEV-DEV)

File Edit View Action Options Tools Help

CPPATIENT, ONE **4AS 400-2** Primary Care Team Unassigned Remote Postings
 000-66-9999 Dec 25, 1941 (59) Provider: CPUSER, ONE Attending: CPPROVIDER, ONE Data WA

Order Sheet Active Orders (includes Pending & Recent Activity) - ALL SERVICES

Service	Order	Start / Stop	Provider	Nrs	Clk	Chart
Admit...	MRI OF HEAD	Start: 09/10/99	CPPROVIDER, T			
Transfer...	Consults GASTROENTEROLOGY Cons Bedside	Start: 06/21/01 16:12	CPPROVIDER, F			
Discharge	PARACENTESIS Cons Bedside	Start: 04/24/01 14:33	CPPROVIDER, F			
Write Orders	Consult: DENTAL, TEST		CPPROVIDER, F			
Allergies	Procedur PULMONARY FUNCTION TEST	Start: 08/17/01 09:52	CPPROVIDER, T			
Diet	PULMONARY Proc Bedside					
Meds, Inpatient	COLONOSCOPY GASTROENTEROLOGY	Start: 08/17/01 09:52	CPPROVIDER, T			
Meds, Outpatient	Proc Bedside					
IV Fluids	PROCEDURE ENDOSCOPY Proc Bedside	Start: 08/13/01 16:53	CPPROVIDER, G			
Lab Tests	ELECTROCARDIOGRAM CARDIOLOGY	Start: 07/10/01 08:25	CPPROVIDER, G			
Imaging	Proc Bedside					
Consult	ELECTROCARDIOGRAM EKG ISC Proc	Start: 07/09/01 11:08	CPPROVIDER, G			
Procedure	Bedside					
Vitals	CATHETERIZATION CARDIOLOGY	Start: 07/09/01 11:08	CPPROVIDER, G			
Text Only Order	CLINIC Proc Bedside					
	CATHETERIZATION CARDIOLOGY	Start: 06/25/01 10:51	CPPROVIDER, G			
	CLINIC Proc Bedside					
	PROCEDURE ENDOSCOPY Proc Bedside	Start: 06/21/01 16:15	CPPROVIDER, F			
	COLONOSCOPY GASTROENTEROLOGY	Start: 06/06/01 13:43	CPPROVIDER, G			
	Proc Bedside					
	PULMONARY FUNCTION TEST	Start: 06/06/01 13:40	CPPROVIDER, G			
	PULMONARY Proc Bedside					
	BONE MARROW ASPIRATE					

Cover Sheet Problems Meds Orders Notes Consults D/C Summ Labs Reports

Fig. 3-7

Check in a new study

Checking in a new study is the next step in the Clinical Procedures process. You need to check in a new study in CP User after a procedure has been ordered. (Keep in mind that the CP check-in is not related to the Scheduling check-in process.)

¹If you want to link multiple results to one procedure, you can check in multiple studies for the same procedure that you ordered through Consults. A warning screen displays telling you that this consult procedure order has already been checked in (Fig. 3-9). After you ensure that you have the correct consult procedure order, you can continue to check in the study. In this way, you do not have to order multiple procedure requests. In this example, the colonoscopy procedure was ordered and a new study for the colonoscopy procedure is being checked in.

1. To check in a new study, first logon to **CP User** and select the patient. Refer to [Selecting a Patient, 2-2](#).
2. Choose **File > Check in New Study** to check in the patient.

Clinical Procedure	Date Ordered	Urgency/Status
Colonoscopy	3/27/2003 2:23:54 PM	Routine/c
Colonoscopy	3/26/2003 9:22:35 AM	Routine/p
Paracentesis	4/24/2001 2:33:31 PM	Routine/p

Location	Date/Time	Status
GI LAB	3/11/2003 9:00:00 ...	INPATIENT APPOI...

²Fig. 3-8

3. Select a Consult procedure order for the selected patient. See Figure 3-8. The Clinical Procedure column lists the consult procedure orders. Notice that the colonoscopy procedure is selected.

¹ Patch MD*1.0*4 September 2006 Check in multiple consult procedure orders warning screen added.

² Patch MD*1.0*4 September 2006 Clinical procedures and outpatient visits listed in reverse chronological order.

Note: You can only select from Clinical Procedure request orders that are in the Pending (p), Scheduled (s), Partial Results (pr), Complete (c), and Active (a) statuses. Discontinue (d) and Cancel statuses are excluded.

¹If the consult procedure order you selected has already been checked in, a warning screen displays (Fig. 3-9). After you ensure that you have the correct consult procedure order, you can continue to check in the study.




Fig. 3-9

4. Depending on the consult procedure you selected, the appropriate instruments for that procedure are displayed. Click the appropriate instrument if more than one is listed, or click **No Instrument** if no instrument is associated with this procedure. OLYMPUS is the appropriate instrument in this case and is selected.
5. You must associate each CP study with a PCE visit, which is the hospital location where the procedure is performed. Required.

For the majority of TIU notes created through CP, the visit association is completed in the background. If a visit has already been recorded but the note wasn't linked (standalone visits, such as telephone or walk-in visits), you can select a visit from the Clinical Procedures Check In edit screen (Fig. 3-8).

To link the CP study to the visit, select information from the Outpatients Visits tab on Figure 3-8. You can also select the New Visit tab and enter NOW for the date and time.

6. Click **Check-In**. The main CP User window, Figure 3-10, is displayed.

¹ Patch MD*1.0*4 September 2006 Check in multiple consult procedure orders warning screen added.




Fig. 3-10

7. If the study is checked-in for an instrument with a uni-directional interface, the status is Ready to Complete. If the study is checked-in for an instrument supported by a bi-directional interface, the status is Pending Instrument Data. (Notice on Figure 3-10, the colonoscopy status for 3/25 and 3/26 is Pending Instrument Data.)
8. At this point, the clinician performs the procedure on the instrument and transmits the results back to VistA.

If the instrument is bi-directional and the Auto Submit to VistA Imaging checkbox is selected for the procedure in CP Manager, the study status changes from Pending Instrument Data to Complete. This occurs after the result has been transmitted to VistA, matched to the study, and copied over to VistA Imaging successfully. The study is ready for interpretation. At this time, the CP process is complete and attachments cannot be associated with this study. See Clinical Procedures Process, Part 2.

If the instrument is uni-directional or if the instrument is bi-directional and the Auto Submit to VistA Imaging checkbox is not selected, the study status is Ready to Complete. Go to the next step (9) to manually complete the CP process.

9. Open the study (Figure 3-10) and add the instrument results and/or external attachments. You can only open studies that have an Error, Complete, Ready to Complete, or New status. When a study is in the Ready to Complete or New status, you can open the study and finish entering any data that was missed. An example of missed data is an external attachment that was not associated with the study.
 - Open this study and add results and/or external attachments. Click **Open Study** or select **File > Open Study**. Figure 3-11 is displayed.
 - Click **+Results** to select and submit the result to Vista Imaging. Only results for the patient and instrument used for the procedure are displayed. To select multiple results, hold down the CTRL key. To select a range of results, highlight the initial result, hold down the Shift key, and then click the last result, Figure 3-12.
 - You can also click **+Files** (Figure 3-11) to add additional attachments from the External Attachment Directory. If the External Attachment Directory has not been defined for this procedure, the last directory that was accessed may be displayed. You can browse for other attachments to link to the study.

Note: If the system parameter **Allow Non-Instrument Attachments** was not selected in CP Manager, **+Files** does not appear on the Clinical Procedures Study screen, you are not permitted to associate additional attachments with the procedure.

10. Submit the study. The images are copied to the RAID and the TIU document is created and associated with the procedure order.

Clinical Procedures Study

File View Attachments

CPPATIENT, ONE Ward: 4AS Rm: 400-2
000-66-9999 MALE DEC 25, 1941 (61)

Allergies

-CP Procedure:
COLONOSCOPY

-IU Note:
TIU Note Title:
COLONOSCOPY

-Consult:
No. 1395

-Instrument:
OLYMPUS

-Status:
Ready to Complete

Attachments:
\\was-gi3a142wd\iman\export\000000bi.bmp

+ Files
+ Results
- Remove

Submit OK Apply Cancel

Fig. 3-11

Instrument Results

Study	Date/Time Performed	Attachment
58665-47309	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\58665-47309-3.txt
58665-47309	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\000000ci.bmp
58665-47309	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\000000bi.bmp
58665-47096	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\58665-47096-3.txt
58665-47096	8/14/2001 11:05:17 AM	\\was-gi3a142wd\iman\export\000000ci.bmp
58665-47096	8/14/2001 11:05:17 AM	\\was-gi3a142wd\iman\export\000000bi.bmp
58665-40711	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\58665-40711-3.txt
58665-40711	8/14/2001 11:05:17 AM	\\was-gi3a142wd\iman\export\000000ci.bmp
58665-40711	8/14/2001 11:05:17 AM	\\was-gi3a142wd\iman\export\000000bi.bmp
58665-38135	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\58665-38135-1.txt
58665-38134	8/14/2001 11:05:17 AM	\\ish-imaging\uploads\58665-38134-1.txt
58664-60699	8/13/2001 4:05:18 PM	\\ish-imaging\uploads\58664-60699-2.txt
58664-60699	8/13/2001 4:05:18 PM	\\olympus_iman\vol1\iman\export\0000004i.bmp

OK Cancel

Fig. 3-12

11. From Figure 3-11, click the magnifying glass under TIU Note to view the TIU Note for that study if it is available. The magnifying glass for the TIU document is unavailable if the result has not been submitted to Vista imaging. Once the result is copied to Vista Imaging, you can view the TIU document of the study before or after the interpretation has been entered, Figure 3-13.




Fig. 3-13

12. From Figure 3-11, you can also click the magnifying glass under Consult to view the Consult report for that study.

```

CONSULT DISPLAY
DATE OF NOTE: AUG 17, 2001@09:59:23  ENTRY DATE: AUG 17, 2001@09:59:23
AUTHOR: CPUSER,TWO                      EXP COSIGNER:
URGENCY:                                STATUS: COMPLETED

PROCEDURE SUMMARY CODE: Normal
DATE/TIME PERFORMED: AUG 14, 2001@11:05

Normal colonoscopy.

/es/ CPUSER,TWO
TECHNICAL WRITER
Signed: 08/17/2001 10:07

=====
AUTHOR & TITLE:                                |
                                                    |DATE:
-----
ID #: _____|ORGANIZATION:HINES DEVELOPMENT |REG #: _____|LOC: 4AS
                |                               |               |RM/BD: 400-
-----
CPPATIENT,ONE   SC VETERAN                      CONSULTATION SHEET
000-66-9999      12/25/1941                      Standard Form 513 (Rev 9-
445 ANY STREET
METROPOLIS      ILLINOIS      62407      Phone: 1-800-PHONE-HOME

```

Fig. 3-14

Updating study status to correct errors

If you open a study in the Error status and have the MD MANAGER key, the Update Study Status window is displayed. You must have the MD Manager key to access the Update Study Status menu option. See your clinical application coordinator or IRM for access to Update Study Status.

You can use Update Study Status to change the status of any study. Be careful when changing the status of a study. With Update Study Status, you can force a status change of a study if a problem occurs that you cannot fix with the Open a Study option.

1. Select **File > Update Study Status**, Fig. 3-15.
2. Select the status you want to change and click **OK**.
3. After you change the status, choose **File > Open a Study** and click **Submit** to resubmit the study.




Fig. 3-15

