

# Setup for the Dental Card Reader

The following is a Card Reader Device setup:

## Device File setup:

NAME: **C103 200/C129D DENTAL CARD READER** \$I: **\_LTA1188:**  
ASK DEVICE: **YES** ASK PARAMETERS: **YES**  
LOCATION OF TERMINAL: **200/C129D DENTAL CARD READER (D/C)**  
USE PARAMETERS: **NOPACK** SUBTYPE: **C-HEI CARD READER**  
TYPE: **TERMINAL** LAT SERVER NODE: **DSV324**  
LAT SERVER PORT: **PORT\_33** VMS DEVICE TYPE: **NOT SPOOLED**  
LAT PORT SPEED: **48** TIMED READ (# OF SECONDS): **999**

## Terminal Type File setup:

NAME: **C-HEI CARD READER** SELECTABLE AT SIGN-ON: **NO**  
RIGHT MARGIN: **80** FORM FEED: **#, \$C(27,91,50,74,27,91,72)**  
PAGE LENGTH: **24** BACK SPACE: **\$C(8)**  
OPEN EXECUTE: **W \*27,"C21",\*13,\*27,"C30",\*13,\*27,"D31",\*13**  
OPEN EXECUTE DESCRIPTION: **SETUP HEI CARD READER FOR DSR/DTR**