

Choose Your Own Device (Don't use STANDARD CAPTIONED OUTPUT)

If you enter YES at the "STANDARD CAPTIONED OUTPUT?" prompt, output is sent to the HOME device, which ordinarily is the terminal you are using. To send output to a printer, you need to answer NO at the "STANDARD CAPTIONED OUTPUT?" prompt. This means you must choose your own print fields, as described above. Doing this allows you to select the output device, however.

If you answer NO to the "STANDARD CAPTIONED OUTPUT?" prompt, you can *still* get CAPTIONED output, however. To do this, choose the CAPTIONED PRINT template at the first "PRINT FIELD:" prompt, as shown below.

Enter a left bracket, "[", and the word CAPTIONED. CAPTIONED OUTPUT from the template works exactly as if you had answered YES at the "STANDARD CAPTIONED OUTPUT?" to request CAPTIONED OUTPUT.

```

STANDARD CAPTION OUTPUT? YES// N <RET> (NO)
FIRST PRINT FIELD: [CAPTIONED

Include COMPUTED fields: (N/Y/R/B): NO// <RET> - No record number (IEN),
no Computed Fields

*****
Heading (S/C): PATIENT LIST// <RET>
DEVICE:
    
```

Even though you don't choose Captioned Output, you can still use the CAPTIONED Print Template for the same effect.

By not choosing "Standard Captioned Output," you can choose your output device.

Figure 7: Choosing a Device and Getting STANDARD CAPTIONED OUTPUT