

Release Notes

USER RELEASE NOTES

The DHCP Patient Data Exchange (PDX) package is a set of software modules designed to manage the transfer of patient information (demographics, episode of care, medication, and diagnostic evaluations) between VA Medical Centers using the MailMan electronic mail utility. Once transferred, this information may be combined with pertinent local information and assembled into a coherent composite record.

PDX V. 1.5 still uses the five main functions for sending, receiving, and displaying patient data. PDX V. 1.5 has been totally redesigned using the List Manager for the user interface. A new file design was implemented to ease the move to an HL7 messaging architecture.

Functionality Changes

- A. The biggest change is in the user interface. PDX V. 1.5 uses the List Manager extensively. The List Manager is a tool designed to display a list of items to the user. It allows the user to browse back and forth through the items one at a time or by screen. It also allows the user to select items from the list to perform some kind of action. List Manager was initially released as part of MAS V. 5.2 for the scheduling application. It will be released in its own name space with the PIMS V. 5.3 release. This is the version PDX uses.
- B. The other big change is the ability to send/request Health Summary components. Health Summary components are being added to provide clinical information from other facilities. This information will be "display only" and cannot be uploaded. The following Health Summary components have been provided.

Adverse Reactions/Allergies	Clinical Warnings	Crisis Notes
Dietetics	Lab Blood Availability	Lab Blood Transfusions
Lab Chemistry & Hematology	Lab Cumulative-Selected	Lab Cumulative-Selected 1
Lab Cumulative-Selected 2	Lab Cumulative-Selected 3	Lab Cumulative-Selected 4
Lab Cytopathology	Lab Microbiology	Lab Microbiology Brief
Lab Orders	Lab Orders Brief	Lab Surgical Pathology
Lab Tests Selected	MAS ADT History	MAS Admissions/Discharges
MAS Clinic Visits Future	MAS Clinic Visits Past	MAS Demographics
MAS Demographics Brief	MAS Disabilities	MAS Discharge Diagnosis
MAS Discharges	MAS Procedures ICD Codes	MAS Surgeries ICD
MAS Transfers	MAS Treating Specialty	Medicine Summary
Orders Current	Pharmacy Intravenous	Pharmacy Outpatient
Pharmacy Unit Dose	Progress Notes	Progress Notes Brief
Radiology Impression	Radiology Profile	Radiology Status
Surgery Reports	Surgery Reports Brief	Vital Signs
Vital Signs Selected		

- C. The interface for requesting/sending PDX information provides the ability to place time and occurrence limits on appropriate data segments. When defining the data segment to use, the user will be prompted for time and/or occurrence limits if the data segment is a Health Summary component which requires limits. A default value for each limit will be taken from the VAQ - PARAMETER file which the user may accept or change. Time and occurrence limits will be checked against the maximum values of information the facility is willing to release without user intervention.
- D. Other components added as "display only"
Means Test
Copay
- E. Other segment changes
Pharmacy is now two separate components (long and short)
A PDX prefix was attached to all the V. 1.0 components
- F. The ability to send/request a predefined group of segments. PDX V. 1.0 provided the ability to transfer patient demographics and pharmacy data within a single transmission. V. 1.5 will enable the user to request only those segments of data that are of interest to the user. With this capability, users and application programs can limit their request to only Pharmacy, PIMS, Copay, or certain Health Summary data.
- G. The ability to send/request data for a patient from multiple sites and multiple segments. In V. 1.0 of PDX, you did not get the option to select the data you wanted. It was predetermined for you. In V. 1.5, you will be able to select any combination of the components.

- H. To expedite user selection of data segments, it will be possible for users to predefine groups of commonly used data segments. These segment groups will work much like mail groups in that the user may enter a group name to select all segments contained in the group.
- I. To expedite user selection of facilities, it will be possible to predefine groups of commonly used facilities. These facility groups will work much like mail groups.
- J. The ability to request data on behalf of another user. This is a new feature for PDX V. 1.5. It allows multiple people to be notified via a mail message when data is returned. This could be done for people who do not have access to PDX.
- K. The ability to include requested information in a mail message. This is a new feature of PDX V. 1.5. It allows the requested data to be included in the notification message. The prompt will allow the input of local mail groups.
- L. The ability to authenticate users via electronic signature. This was added as a security measure to ensure the identity of the requester. All users will have to set up an electronic signature.
- M. The ability to encrypt data using existing data encryption methods. This was added as an optional security measure. Sites wishing to encrypt certain fields for security reasons will have to turn the flag on and add those fields to the VAQ - ENCRYPTED FIELDS file (#394.73).
- N. An application program interface (API) will be provided so that PDX functions may be directly called by other DHCP applications.
- O. Sensitive Patient
When data pertaining to a sensitive patient is transferred, the data will be treated as sensitive at the receiving facility.
- P. A status screen was added to send/request PDX information. The status screen will indicate current information already on file or pending requests for the selected patient. If desired, the user may display information already on file.
- Q. The comments on "process external" have been expanded from 80 characters to a word processing field. This will allow the user to enter as much information as needed.
- R. Load/Edit was enhanced by running all patients through the PIMS duplicate check.

- S. **Display PDX by user**
This option will display all current transactions originated by the user.
- T. **System reports enhanced by providing more predefined and customizable sort.**
New reports were added showing transactions which require manual processing and transactions which are currently on file.
- U. **Security keys have been placed on the following options.**

	Mnemonic	Option	Key	Sub Option key
1.	LD	Load/Edit PDX Data	VAQ LOAD	none
2.	PRGE	Purging	VAQ PURGE	none
3.	RPT	System Reports	VAQ RPT	VAQ RPT USER
4.	EDT	PDX Edit Files	none	VAQ EDIT FILE

TECHNICAL RELEASE NOTES

The technical release notes outline changes made to Patient Data Exchange (PDX) files, routines, templates, options, mail groups, and bulletins.

I. Files

V. 1.5 of PDX required a full redesign of the file structure. Because of this, the files used in V. 1.0 are being marked for deletion. The following table correlates V. 1.0 files to their V. 1.5 counterparts.

Version 1.0		Version 1.5	
PDX TRANSACTION	394	VAQ - TRANSACTION	394.61
PDX DATA	394.1	VAQ - DATA	394.62
PDX PARAMETER	394.2*	VAQ - PARAMETER	394.81
		VAQ - RELEASE GROUP	394.82
PDX STATUS	394.3	VAQ - STATUS	394.85
PDX STATISTICS	394.4	VAQ - WORKLOAD	394.87

*Information contained in File 394.2 was broken into two files

When V. 1.5 of PDX is installed, information contained in V. 1.0 files will be transferred to the new files. After doing this, V. 1.0 files will be marked for deletion. Deletion of V. 1.0 files will be done during the installation of PDX V. 2.0.

New Files

VAQ - DATA SEGMENT	394.71	Defines each segment currently supported by PDX
VAQ - ENCRYPTION METHOD	394.72	Defines each encryption method currently supported by PDX
VAQ - ENCRYPTED FIELDS	394.73	Contains all fields that should be encrypted in PDX requests and unsolicited PDXs transmitted by the facility. This file is only active when encryption has been turned on.
VAQ - OUTGOING GROUP	394.83	Contains groups of facilities commonly accessed by PDX
VAQ - SEGMENT GROUP	394.84	Contains groups of data segments commonly referenced by the facility
VAQ - AUTO-NUMBERING	394.86	Used to implement auto-numbering in PDX. Fields with auto-numbering capability will have an entry in this file.
VAQ - WORK	394.88	Contains the type of work being tracked by the VAQ - WORKLOAD file

II. Routines

V. 1.5 of PDX required a full redesign of all routines. Because of this, the routines used in V. 1.0 are being deleted during the install process. The list of V. 1.0 routines being deleted and a complete list of V. 1.5 routines can be found in the installation guide portion of this manual. The following table provides the routine naming conventions followed for this version.

Functional Area	Namespace
Bulletins	VAQBULxx
Cross-references	VAQXRFxx
Data Display	VAQDISxx
Database Interface	VAQDBLxx
Encryption/Decryption	VAQHSHxx
Exported Routines	VAQPSExx
Inits	VAQIxxxx
Load/Edit PDX Data	VAQLEDxx
Message Administration	VAQADMxx
Message Construction	VAQCONxx
Message Filing	VAQFILxx
Message Parsing	VAQPARxx
Onits	VAQOxxxx
PDX Data Lookup	VAQUPDxx
Post Inits	VAQPSTxx
Process External PDX Request	VAQEXTxx
Purging	VAQPURxx
Request Patient Information	VAQREQxx
Send Unsolicited PDX	VAQUNSxx
System Administration	VAQADSxx
User Authentication	VAQAUTxx
User Interface	VAQUINxx
Utility Routines	VAQUTLxx

III. Templates

A. Print Templates

VAQ PRINT STATS from V. 1.0 of PDX is replaced by VAQ WORK-LOAD REPORT. This report is associated with VAQ - WORKLOAD file (#394.87).

The following new print templates have been added.

VAQ CUR. TRANSACTIONS REPORT	associated with VAQ - TRANSACTION file (#394.61)
VAQ REQUIRES PROCESSING REPORT	associated with VAQ - TRANSACTION file (#394.61)
VAQ DATA SEGMENT LIST	associated with VAQ - DATA SEGMENT file (#394.71)

B. Sort Templates

The following templates, all of which are associated with VAQ -WORKLOAD file (#394.87), replace VAQ SORT STATS from V. 1.0 of PDX.

VAQ WORK-LOAD BY DATE
 VAQ WORK-LOAD BY FACILITY
 VAQ WORK-LOAD BY PATIENT
 VAQ WORK-LOAD BY WORK

The following new sort templates (all associated with VAQ - TRANSACTION file (#394.61)) have been added.

VAQ REQUIRES PROCESSING
 VAQ TRANSACTIONS BY AUTHORIZER
 VAQ TRANSACTIONS BY DATE SENT
 VAQ TRANSACTIONS BY FACILITY
 VAQ TRANSACTIONS BY PATIENT
 VAQ TRANSACTIONS BY RECEIPT
 VAQ TRANSACTIONS BY REQUESTOR
 VAQ TRANSACTIONS BY STATUS

C. Input Templates

The following new input templates have been added.

VAQ EDIT FILE (#394.71)	allow for modification to time and occurrence limits
VAQ EDIT FILE (#394.73)	allows encrypted fields to be added/edited/deleted
VAQ EDIT FILE (#394.81)	allows modification to the site's VAQ - PARAMETER file (#394.81)
VAQ EDIT FILE (#394.82)	allows release groups to be added/edited/deleted
VAQ EDIT FILE (#394.83)	allows outgoing groups to be added/edited/deleted
VAQ EDIT FILE (#394.84)	allows all segment groups to be added/edited/deleted
VAQ EDIT FILE (PRIVATE) (#394.84)	allows private segment groups to be added/edited/deleted
VAQ EDIT FILE (PUBLIC) (#394.84)	allows public segment groups to be added/edited/deleted

IV. Options

All V. 1.0 options have been deleted except VAQ-PDX-SERVER and VAQ PDX PURGE. Deleted options are listed in the installation guide portion of this manual.

The new options are as follows.

VAQ (EDIT) ENCRY FIELD
 VAQ (EDIT) MAX LIMITS
 VAQ (EDIT) OUTGOING GROUP
 VAQ (EDIT) PARAMETER
 VAQ (EDIT) RELEASE GROUP
 VAQ (EDIT) SEGMENT GRP - ALL
 VAQ (EDIT) SEGMENT GRP - PRIV
 VAQ (EDIT) SEGMENT GRP - PUBL
 VAQ (MENU) CUR TRANSACTIONS
 VAQ (MENU) DISPLAY PDX
 VAQ (MENU) EDIT FILES
 VAQ (MENU) MAIN
 VAQ (MENU) PURGING
 VAQ (MENU) SYSTEM REPORTS
 VAQ (MENU) WORK-LOAD
 VAQ PDX DISPLAY (TRN)
 VAQ PDX DISPLAY (USER)
 VAQ PDX LOAD/EDIT
 VAQ PDX PROCESS EXTERNAL
 VAQ PDX PURGE
 VAQ PDX REQUEST
 VAQ PDX UNSOLICITED
 VAQ PURGE BY ENTERED DATE
 VAQ PURGE BY ENTERED LIFE
 VAQ REQUIRES PROCESSING
 VAQ TRANSACTIONS BY AUTHORIZER

VAQ TRANSACTIONS BY DATE SENT
VAQ TRANSACTIONS BY FACILITY
VAQ TRANSACTIONS BY PATIENT
VAQ TRANSACTIONS BY RECEIPT
VAQ TRANSACTIONS BY REQUESTOR
VAQ TRANSACTIONS BY STATUS
VAQ TRANSACTIONS USER DEFINED
VAQ WORKLOAD BY DATE
VAQ WORKLOAD BY FACILITY
VAQ WORKLOAD BY PATIENT
VAQ WORKLOAD BY WORK
VAQ WORKLOAD USER DEFINED
VAQ-PDX-SERVER

V. Mail Groups

The V. 1.0 mail group VAQ PDX USER has been deleted and replaced by the following.

VAQ PDX ERRORS	all PDX errors
VAQ MANUAL PROCESSING	those transactions that require manual processing
VAQ UNSOLICITED RECEIVED	receipt of an unsolicited request

In V. 1.5 of PDX, specialized mail groups were set up to allow PDX support personnel to only subscribe to those groups which meet their needs and eliminate unwanted messages.

VI. Bulletins

The following V. 1.0 bulletins have been deleted.

VAQ PDX SERVER - PROCESS
VAQ PDX SERVER - RETURNED
VAQ PDX SERVER - UNSOLICITED

NOTE: The V. 1.5 messages are not in the BULLETIN file and are produced by the PDX routines. These routines are contained in the subname space VAQBUL*.

Installation Guide

I. Introduction

PDX V. 1.5 uses the newest version of the List Manager which was released with V. 5.3 of PIMS. PIMS V. 5.3 was mandated for October 1, 1993.

Please instruct PDX users to set up an electronic signature code. This can be accomplished through the Edit Electronic Signature Code option of the User's Toolbox menu.

Note: User responses are shown in boldface type.

It is recommended you slave print the installation process. The data printed out during the post-init should be reviewed.

II. Package Integration

A. The following package versions (or higher) must be installed prior to loading this version of PDX.

Kernel	7.0
VA FileMan	19.0
VA MailMan	7.0
PIMS	5.3
Pharmacy	5.6
Integrating Billing	1.5
Health Summary	1.2
OE/RR	1.96
Allergy Tracking	2.2

B. Integration between PDX V. 1.0 and PDX V. 1.5.

PDX V. 1.0 is not forward compatible with PDX V. 1.5 unless patch VAQ*1*11 has been applied by V. 1.0 sites.

The message structure for PDX transmissions has been changed with this version. V. 1.5 sites are able to accept messages in V. 1.0 format. If the V. 1.5 site must respond to the message (e.g., return requested data), it sends the response in the appropriate V. 1.0 format.

However, there is a problem when a V. 1.5 site originates communication with a site using V. 1.0 (e.g., sends a PDX Request or an Unsolicited PDX). The V. 1.0 site cannot interpret the format and the message is lost. This patch will make it possible for V. 1.0 sites to receive these PDX messages.

After the patch is installed, the V. 1.0 site will automatically request a retransmission for each PDX Request or Unsolicited PDX from the V. 1.5 site. This is accomplished without user intervention. The retransmission will be sent in the proper format.

Please note that the extractions for PDX V. 1.0 and PDX V. 1.5 are different and some fields may be in different formats.

C. The following IB routines are distributed with this version of PDX to allow copay to be an extractable component: IBAPDX IBAPDX0 IBAPDX1.

D. GMTSPDX is distributed with this version of PDX to allow for the extraction of the Health Summary components. This routine uses the spooler. If a spooler is not set up at your site, the site will not be able to transmit Health Summary Data. This is significant to note, as it effects the major enhancement to PDX V. 1.5.

III. Package Initialization/Conversion Time Estimates

The following information was derived from our alpha and beta sites. The information presented are averages. Actual times may vary depending on system utilization. The total install process not including the conversion takes on average one hour or less.

<u>Step</u>	<u>MSM</u>	<u>VAX</u>
VAQINIT	2 min	3 min
CONVERSION	4 min per transaction (avg)*	5 min per transaction (avg)*

* These numbers were calculated during peek usage. Conversion time for non peek was about 1 minute per transaction.

Boston MA	conversion ran for 37 hrs, converting 2673 transactions (VAX)
Bedford MA	conversion ran for 4 hrs, converting 180 transactions (VAX)

Users can be brought back on the system before the conversion has finished. Converted transactions will become available for use as they complete the conversion process. Transactions in the process of conversion will not be available for use.

IV. Resource Requirements

Sites will see growth in the following three PDX files.

VAQ - TRANSACTION file (394.61)	200 bytes per entry*
VAQ - DATA file (394.62)	50 bytes per entry*
VAQ - WORKLOAD file (394.87)	110 bytes per entry*

*Since variable fields are in these records, average record size was determined by sampling test sites.

It was determined that a single transaction on average will generate 185 data records. This was determined by taking the total number of entries in the VAQ - DATA file and dividing by the total number of entries in the VAQ - TRANSACTION file.

Disk storage requirements are estimated at approximately 10K per transaction. This was determined by taking the 185 data records at 50 bytes each plus a single transaction record at 200 bytes plus a variable number of workload records at 110 bytes each and dividing the total by 1024 bytes.

Global Growth can be controlled by the purge features of PDX.

V. Pre-Installation

<u>Step</u>	<u>Description</u>
-------------	--------------------

1. Insure that the patches which affect the PDX purge routine, VAQPRG, have been installed. These are VAQ*1*8 and VAQ*1*15.
2. Request PDX staff to process all V. 1.0 requests. This will cut down on conversion time.
3. Run V. 1.0 purger.
This can be accomplished by rescheduling the start time of the background job (VAQ PDX PURGE) to current date and time plus two minutes.

Go into the Task Manager and select Schedule/Unschedule Options. At the "Select OPTION to schedule or reschedule" prompt, enter **VAQ PDX PURGE**. Enter the **current date/time** at the "Queued to run at what time:" prompt.

4. Determine how PDX V. 1.0 has been distributed (PDX main menu or custom menus) by running Option Access By User. This report is under the Menu Management menu.
5. Backup System(s). [optional]
6. Remove VAQ namespaced routines from the routine map. V. 1.0 had recommended mapping of the following routines.

VAQDSP3 VAQDSP6 VAQUTL VAQUTL1 VAQUTL2

7. A spooling device must be installed in order to transmit Health Summary data. Typically VAX sites will already have one installed, while 486 sites will not. 486 sites should refer to the 486 Cookbook for assistance in setting up the spooler. The Health Summary Developers recommend allocating 2 Mb of disk space for the spooler.
8. Determine Spool device name.
9. The POSTMASTER must be given access to the spool device used by your facility. This can be accomplished by using the Edit an Existing User option [XUSEREDIT] and entering YES to the ALLOWED TO USE SPOOLER field.

VI. Installation

<u>Step</u>	<u>Description</u>
-------------	--------------------

1. Request PDX users be off the system.
2. Sign into UCI where package will be initially installed.
(use 40K partition for MSM)
3. Disable Journaling.
4. Delete VAQ namespaced routines.

D ^%RDELETE (VAX) or **D ^%RDEL** (MSM)

Routine(s) ? > **VAQ***

Routine(s) ? > **-VAQWK** (do not delete this routine - used by PIMS)

(See Attachment A for deleted routines.)

5. Load PDX V. 1.5 tape/disk using routine restore.
Note: Sites using MSM will find the name of the file containing the routines on the disk label.

D ^%RR

(See Attachment B for distributed routine list.)

6. Verify integrity values.

D ^VAQNTEG

7. Setup system variables.

D ^XUP

Verify DUZ, DT, DTIME and U are defined and DUZ(0)="@". DUZ variable must be defined as an active user number and DUZ(0) variable must equal "@" in order to initialize.

<u>Step</u>	<u>Description</u>
-------------	--------------------

8.	D ^VAQINIT
----	-------------------

(Follow Attachment C for VAQINIT responses.)

Please answer all initialization questions carefully.

It is strongly recommended you take the bolded responses from the guide.

A file conversion will be queued at the end of the initialization process. This utility has been included for your convenience. If you do not want to run the conversion, you must enter an up-arrow (^) when prompted for the conversion's device. Please refer to Attachment E if you have chosen to skip the conversion or the file conversion does not run to completion.

9.	Update Health Summary Site Parameter option.
----	--

The option GMTS IRM/ADPAC PARAMETER EDIT must be modified to prompt for the SPOOL DEVICE NAME. This is done by adding ";.04" to the end of the existing DR {DIE} field. A sample of this is shown below.

----- SAMPLE -----

D ^XUP

Setting up programmer environment

Terminal Type set to: C-VT100

Select OPTION NAME: **EVE**

Select Systems Manager Menu Option: **MENU MANAGEMENT**

Select Menu Management Option: **EDIT OPTIONS**

Select OPTION to edit: **GMTS IRM/ADPAC PARAMETER EDIT**

NAME: GMTS IRM/ADPAC PARAMETER EDIT Replace **^DR {DIE}**

DR {DIE}: .02;.03// **.02;.03;.04**

NO UP-ARROW: ^ (*user to enter up-arrow*)

Select OPTION NAME: **GMTS IRM/ADPAC MAINT MENU**

Select Health Summary IRM Maintenance Menu Option: **7 Edit Health Summary Site Parameters**

Select HEALTH SUMMARY PARAMETERS: **HOSPITAL**

PROMPT FOR ACTION PROFILE:// **<RET>** (*take the value in field*)

INCLUDE COMMENTS FOR LABS:// **<RET>** (*take the value in field*)

SPOOL DEVICE NAME: < **enter device name** >

- | <u>Step</u> | <u>Description</u> |
|-------------|--|
| 10. | Using the Mail Group Edit option [XMEDITMG], add members to the newly created mail groups. The mail groups are as follows. |

VAQ PDX ERRORS

It is recommended that the System Manager be added to this group in order to monitor incoming messages which cannot be served (errors). These messages are saved in S.VAQ-PDX-SERVER mail basket.

VAQ MANUAL PROCESSING

VAQ UNSOLICITED RECEIVED

- | | |
|-----|--|
| 11. | The PDX server (VAQ-PDX-SERVER) needs to be edited in order to associate a mail group with the server. The mail group to add is VAQ PDX ERRORS. The installer will also have to change the server action from QUEUE SERVER ROUTINE to RUN IMMEDIATELY. |
|-----|--|

----- SAMPLE -----

>**D ^XUP**

Setting up programmer environment

Terminal Type set to: C-VT100

Select OPTION NAME: **EVE** Systems Manager Menu

Select Systems Manager Menu Option: **MENU Management**

Select Menu Management Option: **EDIT options**

Select OPTION to edit: **VAQ-PDX-SERVER** PDX Server

NAME: VAQ-PDX-SERVER// **^SERVER ACTION**

SERVER ACTION: QUEUE SERVER ROUTINE// **RUN IMMEDIATELY**

SERVER MAIL GROUP: VAQ PDX ERRORS// **^**

- | | |
|-----|---|
| 12. | Enter Task Manager and select Schedule/Unschedule Options. At the "Select OPTION to schedule or reschedule" prompt, enter VAQ PDX PURGE and enter the values for the following fields. Fields not listed need to be reviewed for site specific values. |
|-----|---|

QUEUED TO RUN AT WHAT TIME: **T+15@2400**

RESCHEDULING FREQUENCY: **15D**

- | | |
|-----|--|
| 13. | Verify that the global attributes of ^VAT are as follows.
SYSTEM = RWD WORLD = RWD GROUP = RWD UCI = RWD |
|-----|--|

- | <u>Step</u> | <u>Description</u> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---------------------------|--|-----------|----------|-------------------|------------------|-----------------------|--------------------------------|---------------------------|--------------------------------|-----------------------|--------------------------------|----------------|---|---------------------|--|-------------------------|--|---------------------|--|--------------------|------------------------------|---------------------|-------------------------------------|-------------------------|---|---------------------|---|---------------------|----------------------------------|---------------------|---|------------------|-------------------------|
| 14. | For all other systems move PDX routines, and translate ^VAT (if necessary). | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 15. | Delete obsolete options. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | <table border="0"> <tr> <td>1 VAQ PDX</td> <td>PDX Menu</td> </tr> <tr> <td>2 VAQ PDX DISPLAY</td> <td>Display PDX Data</td> </tr> <tr> <td>3 VAQ PDX DISPLAY MAS</td> <td>Extended - View A Specific PDX</td> </tr> <tr> <td>4 VAQ PDX DISPLAY MAS/PHA</td> <td>Extended - View A Specific PDX</td> </tr> <tr> <td>5 VAQ PDX DISPLAY PHA</td> <td>Extended - View A Specific PDX</td> </tr> <tr> <td>6 VAQ PDX LOAD</td> <td>Load/Edit Patient File with PDX Information</td> </tr> <tr> <td>7 VAQ PDX PRINT MAS</td> <td>Brief - View Most Recent PDX From A Site</td> </tr> <tr> <td>8 VAQ PDX PRINT MAS/PHA</td> <td>Brief - View Most Recent PDX From A Site</td> </tr> <tr> <td>9 VAQ PDX PRINT PHA</td> <td>Brief - View Most Recent PDX From A Site</td> </tr> <tr> <td>10 VAQ PDX PROCESS</td> <td>Process External PDX Request</td> </tr> <tr> <td>11 VAQ PDX SHOW MAS</td> <td>MAS - Display PDX's MAS Information</td> </tr> <tr> <td>12 VAQ PDX SHOW MAS/PHA</td> <td>Both - Display PDX's MAS & Pharmacy Information</td> </tr> <tr> <td>13 VAQ PDX SHOW PHA</td> <td>Pharmacy - Display PDX's Pharmacy Information</td> </tr> <tr> <td>14 VAQ PDX STAT DAY</td> <td>Today - Print Entries Made Today</td> </tr> <tr> <td>15 VAQ PDX STAT GEN</td> <td>General - Print Entries Made During A Given Range</td> </tr> <tr> <td>16 VAQ PDX STATS</td> <td>Statistics File Toolbox</td> </tr> </table> | 1 VAQ PDX | PDX Menu | 2 VAQ PDX DISPLAY | Display PDX Data | 3 VAQ PDX DISPLAY MAS | Extended - View A Specific PDX | 4 VAQ PDX DISPLAY MAS/PHA | Extended - View A Specific PDX | 5 VAQ PDX DISPLAY PHA | Extended - View A Specific PDX | 6 VAQ PDX LOAD | Load/Edit Patient File with PDX Information | 7 VAQ PDX PRINT MAS | Brief - View Most Recent PDX From A Site | 8 VAQ PDX PRINT MAS/PHA | Brief - View Most Recent PDX From A Site | 9 VAQ PDX PRINT PHA | Brief - View Most Recent PDX From A Site | 10 VAQ PDX PROCESS | Process External PDX Request | 11 VAQ PDX SHOW MAS | MAS - Display PDX's MAS Information | 12 VAQ PDX SHOW MAS/PHA | Both - Display PDX's MAS & Pharmacy Information | 13 VAQ PDX SHOW PHA | Pharmacy - Display PDX's Pharmacy Information | 14 VAQ PDX STAT DAY | Today - Print Entries Made Today | 15 VAQ PDX STAT GEN | General - Print Entries Made During A Given Range | 16 VAQ PDX STATS | Statistics File Toolbox |
| 1 VAQ PDX | PDX Menu | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 VAQ PDX DISPLAY | Display PDX Data | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 VAQ PDX DISPLAY MAS | Extended - View A Specific PDX | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 VAQ PDX DISPLAY MAS/PHA | Extended - View A Specific PDX | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 VAQ PDX DISPLAY PHA | Extended - View A Specific PDX | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 VAQ PDX LOAD | Load/Edit Patient File with PDX Information | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 7 VAQ PDX PRINT MAS | Brief - View Most Recent PDX From A Site | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 8 VAQ PDX PRINT MAS/PHA | Brief - View Most Recent PDX From A Site | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 9 VAQ PDX PRINT PHA | Brief - View Most Recent PDX From A Site | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 10 VAQ PDX PROCESS | Process External PDX Request | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 11 VAQ PDX SHOW MAS | MAS - Display PDX's MAS Information | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 12 VAQ PDX SHOW MAS/PHA | Both - Display PDX's MAS & Pharmacy Information | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 13 VAQ PDX SHOW PHA | Pharmacy - Display PDX's Pharmacy Information | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 14 VAQ PDX STAT DAY | Today - Print Entries Made Today | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 15 VAQ PDX STAT GEN | General - Print Entries Made During A Given Range | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 16 VAQ PDX STATS | Statistics File Toolbox | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 16. | Delete obsolete templates. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | Input Templates - None
Sort Templates - VAQ SORT STATS
Print Templates- VAQ PRINT STATS | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 17. | Delete obsolete mail group. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | VAQ PDX USERS | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 18. | Delete obsolete bulletins. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | VAQ PDX SERVER-PROCESS
VAQ PDX SERVER-RETURNED
VAQ PDX SERVER-UNSOLICITED | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 19. | Delete the initialization routines VAQIN*, VAQON*, and VAQPSL*. Do not delete VAQINITY. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

<u>Step</u>	<u>Description</u>
20.	At the end of the process, a mail message is sent to the development ISC indicating that the install is complete.
21.	Enable Journaling.

Version 1.0 Routines - Attachment A

VAQADD	VAQADD1	VAQADD2	VAQADD3	VAQADD4	VAQADD5	VAQADD7	VAQDEBUG
VAQDEBUG1	VAQDEBUG2	VAQDSP1	VAQDSP10	VAQDSP11	VAQDSP12	VAQDSP13	VAQDSP14
VAQDSP15	VAQDSP2	VAQDSP3	VAQDSP3A	VAQDSP4	VAQDSP5	VAQDSP6	VAQDSP7
VAQDSP8	VAQDSP9	VAQIN001	VAQIN002	VAQIN003	VAQIN004	VAQIN005	VAQIN006
VAQIN007	VAQIN008	VAQIN009	VAQIN010	VAQIN011	VAQIN012	VAQIN013	VAQIN014
VAQIN015	VAQINIT	VAQINIT0	VAQINIT1	VAQINIT2	VAQINIT3	VAQMAS	VAQMAS1
VAQMAS1A	VAQMAS2	VAQMAS2A	VAQMAS3	VAQMAS4	VAQMAS5	VAQMAS6	VAQMAS6A
VAQNTEG	VAQPHA1	VAQPRG	VAQPRT	VAQPRT1	VAQRQP	VAQRQP1	VAQRQR
VAQRQR1	VAQRQU	VAQST1	VAQST2	VAQST3	VAQUTL	VAQUTL1	VAQUTL2

72 routines

Version 1.5 Routines - Attachment B

VAQADM2	VAQADM21	VAQADM22	VAQADM23	VAQADM5	VAQADM50	VAQADM51	VAQADS01
VAQAUT	VAQBUL	VAQBUL01	VAQBUL02	VAQBUL03	VAQBUL04	VAQBUL05	VAQBUL06
VAQBUL07	VAQCON	VAQCON0	VAQCON1	VAQCON2	VAQCON3	VAQCON4	VAQCON5
VAQCON6	VAQCON7	VAQCON8	VAQCON93	VAQCON94	VAQCON95	VAQCON96	VAQCON97
VAQCON98	VAQCON99	VAQDBI	VAQDBIH1	VAQDBIH2	VAQDBIH3	VAQDBIH4	VAQDBIM
VAQDBIM0	VAQDBIM1	VAQDBIM2	VAQDBIM3	VAQDBIM4	VAQDBIP	VAQDBIP1	VAQDBIP2
VAQDBIP3	VAQDBIP4	VAQDBIP5	VAQDBIP6	VAQDBIP7	VAQDBIP8	VAQDIS01	VAQDIS10
VAQDIS11	VAQDIS12	VAQDIS15	VAQDIS16	VAQDIS17	VAQDIS20	VAQDIS21	VAQDIS22
VAQDIS23	VAQDIS24	VAQDIS25	VAQDIS26	VAQDIS27	VAQDIS28	VAQDIS29	VAQDIS30
VAQDIS31	VAQDIS32	VAQDIS33	VAQDIS40	VAQDIS41	VAQDIS42	VAQDIS43	VAQEXT01
VAQEXT02	VAQEXT03	VAQEXT04	VAQEXT05	VAQEXT06	VAQFIL10	VAQFIL11	VAQFIL12
VAQFIL13	VAQFIL14	VAQFIL15	VAQFIL16	VAQFIL17	VAQFIL18	VAQFILE	VAQFILE1
VAQSHH	VAQSHH1	VAQIN001	VAQIN002	VAQIN003	VAQIN004	VAQIN005	VAQIN006
VAQIN007	VAQIN008	VAQIN009	VAQIN010	VAQIN011	VAQIN012	VAQIN013	VAQIN014
VAQIN015	VAQIN016	VAQIN017	VAQIN018	VAQIN019	VAQIN020	VAQIN021	VAQIN022
VAQIN023	VAQIN024	VAQIN025	VAQIN026	VAQIN027	VAQIN028	VAQIN029	VAQIN030
VAQIN031	VAQIN032	VAQIN033	VAQIN034	VAQIN035	VAQIN036	VAQIN037	VAQIN038
VAQIN039	VAQIN040	VAQIN041	VAQIN042	VAQIN043	VAQIN044	VAQIN045	VAQIN046
VAQIN047	VAQIN048	VAQIN049	VAQIN050	VAQIN051	VAQIN052	VAQIN053	VAQIN054
VAQIN055	VAQIN056	VAQIN057	VAQIN058	VAQIN059	VAQIN060	VAQIN061	VAQIN062
VAQIN063	VAQIN064	VAQIN065	VAQIN066	VAQIN067	VAQIN068	VAQIN069	VAQIN070
VAQINIT	VAQINIT1	VAQINIT2	VAQINIT3	VAQINIT4	VAQINIT5	VAQINITY	VAQLED01
VAQLED02	VAQLED03	VAQLED04	VAQLED05	VAQLED07	VAQLED09	VAQLED10	VAQNTEG
VAQNTEG0	VAQON001	VAQON002	VAQON003	VAQON004	VAQON005	VAQONIT	VAQONIT1
VAQONIT2	VAQONIT3	VAQPAR1	VAQPAR10	VAQPAR11	VAQPAR6	VAQPAR60	VAQPAR61
VAQPSE01	VAQPSE02	VAQPSE03	VAQPSE04	VAQPSL	VAQPSL1	VAQPSL2	VAQPSL3
VAQPST01	VAQPST02	VAQPST03	VAQPST04	VAQPST05	VAQPST10	VAQPST20	VAQPST21
VAQPST22	VAQPST23	VAQPST24	VAQPST25	VAQPST30	VAQPST31	VAQPST40	VAQPUR
VAQPUR10	VAQPUR11	VAQREQ01	VAQREQ02	VAQREQ03	VAQREQ04	VAQREQ05	VAQREQ06
VAQREQ07	VAQREQ08	VAQREQ09	VAQREQ10	VAQREQ11	VAQUIN01	VAQUPD1	VAQUPD2
VAQUPD25	VAQUTL1	VAQUTL2	VAQUTL3	VAQUTL4	VAQUTL92	VAQUTL93	VAQUTL94
VAQUTL95	VAQUTL96	VAQUTL97	VAQUTL98	VAQUTL99	VAQXRF1	VAQXRF2	VAQXRF3

256 routines

PDX V. 1.5 Sample Installation - Attachment C

>D ^VAQINIT

This version (#1.5) of 'VAQINIT' was created on 12-NOV-1993
(at PDX DEVELOPMENT (PRIMARY), by VA FileMan V.20.0)

I AM GOING TO SET UP THE FOLLOWING FILES:

142.99 HEALTH SUMMARY PARAMETERS (Partial Definition)
Note: You already have the 'HEALTH SUMMARY PARAMETERS' File.

394 *PDX TRANSACTION
*BUT YOU ALREADY HAVE 'PDX TRANSACTION' AS FILE #394!
Shall I change the NAME of the file to *PDX TRANSACTION? NO// **YES**

394 *PDX TRANSACTION
Note: You already have the '*PDX TRANSACTION' File.

394.1 *PDX DATA
*BUT YOU ALREADY HAVE 'PDX DATA' AS FILE #394.1!
Shall I change the NAME of the file to *PDX DATA? NO// **YES**

394.1 *PDX DATA
Note: You already have the '*PDX DATA' File.

394.2 *PDX PARAMETER
*BUT YOU ALREADY HAVE 'PDX PARAMETER' AS FILE #394.2!
Shall I change the NAME of the file to *PDX PARAMETER? NO// **YES**

394.2 *PDX PARAMETER
Note: You already have the '*PDX PARAMETER' File.

394.3 *PDX STATUS (including data)
*BUT YOU ALREADY HAVE 'PDX STATUS' AS FILE #394.3!
Shall I change the NAME of the file to *PDX STATUS? NO// **YES**

394.3 *PDX STATUS (including data)
Note: You already have the '*PDX STATUS' File.
I will OVERWRITE your data with mine.

394.4 *PDX STATISTICS
*BUT YOU ALREADY HAVE 'PDX STATISTICS' AS FILE #394.4!
Shall I change the NAME of the file to *PDX STATISTICS? NO// **YES**

394.61 VAQ - TRANSACTION

394.62 VAQ - DATA

394.71 VAQ - DATA SEGMENT (including data)
I will OVERWRITE your data with mine.

394.72 VAQ - ENCRYPTION METHOD (including data)
I will OVERWRITE your data with mine.

394.73 VAQ - ENCRYPTED FIELDS (including data)
I will MERGE your data with mine.

- 394.81 VAQ - PARAMETER
- 394.82 VAQ - RELEASE GROUP
- 394.83 VAQ - OUTGOING GROUP
- 394.84 VAQ - SEGMENT GROUP
- 394.85 VAQ - STATUS (including data)
I will OVERWRITE your data with mine.
- 394.86 VAQ - AUTO-NUMBERING
- 394.87 VAQ - WORKLOAD
- 394.88 VAQ - WORK (including data)
I will OVERWRITE your data with mine.

SHALL I WRITE OVER FILE SECURITY CODES? NO// **<YES>**
 NOTE: This package also contains SORT TEMPLATES
 SHALL I WRITE OVER EXISTING SORT TEMPLATES OF THE SAME NAME? YES// **<RET>**
 NOTE: This package also contains INPUT TEMPLATES
 SHALL I WRITE OVER EXISTING INPUT TEMPLATES OF THE SAME NAME? YES// **<RET>**
 NOTE: This package also contains PRINT TEMPLATES
 SHALL I WRITE OVER EXISTING PRINT TEMPLATES OF THE SAME NAME? YES// **<RET>**
 NOTE: This package also contains FUNCTIONS
 SHALL I WRITE OVER EXISTING FUNCTIONS OF THE SAME NAME? YES// **<RET>**
 NOTE: This package also contains SECURITY KEYS
 SHALL I WRITE OVER EXISTING SECURITY KEYS OF THE SAME NAME? YES// **<RET>**
 NOTE: This package also contains OPTIONS
 SHALL I WRITE OVER EXISTING OPTIONS OF THE SAME NAME? YES// **<RET>**

ARE YOU SURE EVERYTHING'S OK? NO// **YES**

...SORRY, LET ME THINK ABOUT THAT A MOMENT.....
 'VAQ (EDIT) ENCRY FIELDS' Option Filed
 'VAQ (EDIT) MAX LIMITS' Option Filed
 'VAQ (EDIT) OUTGOING GROUP' Option Filed
 'VAQ (EDIT) PARAMETER' Option Filed
 'VAQ (EDIT) RELEASE GROUP' Option Filed
 'VAQ (EDIT) SEGMENT GRP - ALL' Option Filed
 'VAQ (EDIT) SEGMENT GRP - PRIV' Option Filed
 'VAQ (EDIT) SEGMENT GRP - PUBL' Option Filed
 'VAQ (MENU) CUR TRANSACTIONS' Option Filed
 'VAQ (MENU) DISPLAY PDX' Option Filed
 'VAQ (MENU) EDIT FILES' Option Filed
 'VAQ (MENU) MAIN' Option Filed
 'VAQ (MENU) PURGING' Option Filed
 'VAQ (MENU) SYSTEM REPORTS' Option Filed
 'VAQ (MENU) WORK-LOAD' Option Filed
 'VAQ PDX DISPLAY (TRN)' Option Filed
 'VAQ PDX DISPLAY (USER)' Option Filed
 'VAQ PDX LOAD/EDIT' Option Filed
 'VAQ PDX PROCESS EXTERNAL' Option Filed

```
'VAQ PDX PURGE' Option Filed
'VAQ PDX REQUEST' Option Filed
'VAQ PDX UNSOLICITED' Option Filed
'VAQ PURGE BY ENTERED DATE' Option Filed
'VAQ PURGE BY ENTERED LIFE' Option Filed
'VAQ REQUIRES PROCESSING' Option Filed
'VAQ TRANSACTIONS BY AUTHORIZER' Option Filed
'VAQ TRANSACTIONS BY DATE SENT' Option Filed
'VAQ TRANSACTIONS BY FACILITY' Option Filed
'VAQ TRANSACTIONS BY PATIENT' Option Filed
'VAQ TRANSACTIONS BY RECEIPT' Option Filed
'VAQ TRANSACTIONS BY REQUESTOR' Option Filed
'VAQ TRANSACTIONS BY STATUS' Option Filed
'VAQ TRANSACTIONS USER DEFINED' Option Filed
'VAQ WORKLOAD BY DATE' Option Filed
'VAQ WORKLOAD BY FACILITY' Option Filed
'VAQ WORKLOAD BY PATIENT' Option Filed
'VAQ WORKLOAD BY WORK' Option Filed
'VAQ WORKLOAD USER DEFINED' Option Filed
'VAQ-PDX-SERVER' Option Filed.....
NOTE THAT FILE SECURITY-CODE PROTECTION HAS BEEN MADE
```

Begin of Post Init Process

Installing protocols for use by the list processor

This version of 'VAQONIT' was created on 21-JUL-1993
 (at PDX DEVELOPMENT (PRIMARY), by OE/RR V.2.5)

PROTOCOL INSTALLATION

```
...OK, this may take a while, hold on PLEASE.....
'VAQ ADD PATIENT' Protocol Filed
'VAQ ADD/EDIT REQUEST' Protocol Filed
'VAQ CHANGE PATIENT' Protocol Filed
'VAQ COPY REQUEST' Protocol Filed
'VAQ CREATE REQUEST' Protocol Filed
'VAQ DIS ALL SEGMENT' Protocol Filed
'VAQ DIS SELECTED SEGMENT' Protocol Filed
'VAQ DIS1 (MENU)' Protocol Filed
'VAQ DISPLAY BY REQUESTOR' Protocol Filed
'VAQ DISPLAY PDX' Protocol Filed
'VAQ DISPLAY SELECT' Protocol Filed
'VAQ DUPLICATE' Protocol Filed
'VAQ LOAD DATA' Protocol Filed
'VAQ LOAD EDIT' Protocol Filed
'VAQ LOAD FIELD' Protocol Filed
'VAQ NEW PATIENT' Protocol Filed
  VAQ DUPLICATE added as item to VAQ NEW PATIENT.
'VAQ PDX1 (MENU)' Protocol Filed
  VAQ DISPLAY PDX added as item to VAQ PDX1 (MENU).
  VAQ CREATE REQUEST added as item to VAQ PDX1 (MENU).
'VAQ PDX10 (MENU)' Protocol Filed
  VAQ DISPLAY BY REQUESTOR added as item to VAQ PDX10 (MENU).
'VAQ PDX11 (MENU)' Protocol filed
  VAQ DIS SELECTED SEGMENT added as item to VAQ PDX11 (MENU).
  VAQ DIS ALL SEGMENT added as item to VAQ PDX11 (MENU).
'VAQ PDX12 (MENU)' Protocol filed
```

```
'VAQ PDX2 (MENU)' Protocol Filed
  VAQ CHANGE PATIENT added as item to VAQ PDX2 (MENU).
  VAQ ADD/EDIT REQUEST added as item to VAQ PDX2 (MENU).
  VAQ TRANSMIT REQUEST added as item to VAQ PDX2 (MENU).
  VAQ COPY REQUEST added as item to VAQ PDX2 (MENU).
'VAQ PDX3 (MENU)' Protocol Filed
  VAQ PROCESS MANUAL added as item to VAQ PDX3 (MENU).
'VAQ PDX4 (MENU)' Protocol Filed
  VAQ PROCESS RELEASE added as item to VAQ PDX4 (MENU).
  VAQ PROCESS REJECT added as item to VAQ PDX4 (MENU).
'VAQ PDX5 (MENU)' Protocol Filed
  VAQ LOAD EDIT added as item to VAQ PDX5 (MENU).
'VAQ PDX6 (MENU)' Protocol Filed
  VAQ LOAD DATA added as item to VAQ PDX6 (MENU).
  VAQ LOAD FIELD added as item to VAQ PDX6 (MENU).
'VAQ PDX7 (MENU)' Protocol Filed
  VAQ ADD PATIENT added as item to VAQ PDX7 (MENU).
'VAQ PDX8 (MENU)' Protocol Filed
  VAQ DUPLICATE added as item to VAQ PDX8 (MENU).
  VAQ NEW PATIENT added as item to VAQ PDX8 (MENU).
'VAQ PDX9 (MENU)' Protocol Filed
  VAQ DISPLAY SELECT added as item to VAQ PDX9 (MENU).
'VAQ PROCESS MANUAL' Protocol Filed
'VAQ PROCESS REJECT' Protocol Filed
'VAQ PROCESS RELEASE' Protocol Filed
'VAQ TRANSMIT REQUEST' Protocol Filed
OK, Protocol Installation is Complete.
```

Protocol install completed

Installing list templates for use by list processor

```
'VAQ DIS MIN NUPD' List Template...Filed.
'VAQ DIS MIN UPD' List Template...Filed.
'VAQ DIS PATIENT PDX9' List Template...Filed.
'VAQ DIS REQUESTOR PDX10' List Template...Filed.
'VAQ DISPLAY DATA PDX12' List Template...Filed.
'VAQ DISPLAY MINIMUM' List Template...Filed.
'VAQ DISPLAY SEGMENT PDX11' List Template...Filed.
'VAQ DUPLICATE PDX8' List Template...Filed.
'VAQ LED DIFFERENCES PDX6' List Template...Filed.
'VAQ LED STATUS PDX5' List Template...Filed.
'VAQ MATCHES PDX8' List Template...Filed.
'VAQ PROCESS PDX3' List Template...Filed.
'VAQ PROCESS PDX4' List Template...Filed.
'VAQ REQUEST PDX2' List Template...Filed.
'VAQ STATUS PDX1' List Template...Filed.
```

** List Template install completed

>>> Exported routines will now be loaded

```
    Copying of VAQPSE01 into GMTSPDX      Done
    Copying of VAQPSE02 into IBAPDX       Done
    Copying of VAQPSE03 into IBAPDX0     Done
    Copying of VAQPSE04 into IBAPDX1     Done
```

Done

Installation Guide

Initialization of VAQ Parameter file *(SEE ATTACHMENT D #1)*

Enter missing field(s)

Default Time Limit:

Default Occurrence Limit:

** Missing fields added, initialization complete

Updating pointers to Health Summary components and initializing
maximum time and occurrence limits (when appropriate).....

Updating completed

Initialization of VAQ - Auto-numbering file

** Initialization of VAQ - Auto-numbering file complete

Initialization of VAQ - Encrypted Fields File... (add/edit/delete)

Select VAQ - ENCRYPTED FIELDS Encrypt Field: *(SEE ATTACHMENT D #2)*

** Initialization of VAQ - Encrypted Fields File complete

Initialization of VAQ - Release Group File...

Updating VAQ - Release Group file from version 1.0

Update from version 1 completed

Add/Edit/Delete entries in VAQ - Release Group

Select VAQ - RELEASE GROUP Remote Facility: *(SEE ATTACHMENT D #3)*

** Initialization of VAQ - Release Group File complete

Initialization of VAQ - Outgoing Group File... (add/edit/delete)

Select VAQ - OUTGOING GROUP Group Name: *(SEE ATTACHMENT D #4)*

** Initialization of VAQ - Outgoing Group File complete

Initialization of VAQ - Segment Group File...

Creating a segment group called "ALL"
This group will contain all data segments
Adverse Reactions/Allergies - added
Clinical Warnings - added
Crisis Notes - added
Dietetics - added
Integrated Billing - added
Lab Blood Availability - added
Lab Blood Transfusions - added
Lab Chemistry & Hematology - added

Lab Cumulative-Selected - added
Lab Cumulative-Selected 1 - added
Lab Cumulative-Selected 2 - added
Lab Cumulative-Selected 3 - added
Lab Cumulative-Selected 4 - added
Lab Cytopathology - added
Lab Microbiology - added
Lab Microbiology Brief - added
Lab Orders - added
Lab Orders Brief - added
Lab Surgical Pathology - added
Lab Tests Selected - added
MAS ADT History - added
MAS Admissions/Discharges - added
MAS Clinic Visits Future - added
MAS Clinic Visits Past - added
MAS Demographics - added
MAS Demographics Brief - added
MAS Disabilities - added
MAS Discharge Diagnosis - added
MAS Discharges - added
MAS Minimum Patient Information - added
MAS Procedures ICD Codes - added
MAS Registration Information - added
MAS Surgeries ICD Codes - added
MAS Transfers - added
MAS Treating Specialty - added
Means Test Information - added
Medication Profile - Long Form - added
Medication Profile - Short For - added
Medicine Summary - added
Orders Current - added
Pharmacy Intravenous - added
Pharmacy Outpatient - added
Pharmacy Unit Dose - added
Progress Notes - added
Progress Notes Brief - added
Radiology Impression - added
Radiology Profile - added
Radiology Status - added
Surgery Reports - added
Surgery Reports Brief - added
Vital Signs - added
Vital Signs Selected - added

Done

Create entries in Segment Groups for Health Summary Type File? NO//

Add/Edit/Delete entries in VAQ - Segment Group File

Select VAQ - SEGMENT GROUP Group Name: **(SEE ATTACHMENT D #5)**

** Initialization of VAQ - Segment Group File complete

Creating Mail Groups for PDX

'VAQ PDX ERRORS'	mail group created
'VAQ MANUAL PROCESSING'	mail group created
'VAQ UNSOLICITED RECEIVED'	mail group created

Mail Groups created

- Conversion of version 1.0 files will now be tasked -

Entering 'HOME' as the device for output will cause conversion to be run without an output device. It is recommended that a device be chosen so that errors during the conversion can be reported.

Entering '^' as the device for output will skip the conversion process. Please refer to the INSTALLATION GUIDE if you choose to do this.

Enter device to use during conversion: HOME// *{conditional response}*

{If you have chosen to convert V. 1.0 data, respond as follows:}

Enter device to use during conversion: HOME// **<DEVICE>** *{user should enter a device on the system}*

Requested Start Time: NOW// **<RET>**

Conversion tasked (TASK #)

{If you have chosen not to convert V. 1.0 data, respond as follows:}

Enter device to use during conversion: HOME// ^

- Conversion will not be done at this time -

To run conversion at a later date the entry point TASK^VAQPST20 should be used.

If you have chosen to skip the conversion the entry point DELETE^VAQPST24(1) must be used in order to delete entries contained in the 1.0 files.

{The following information will be printed after answering the device prompt shown above.}

```
*****
* IMPORTANT * The following things need to be done.
*
* - Members need to be added to the newly created mail groups.
*   The option is XMEDITMG.
*
* - The PDX Server (VAQ-PDX-SERVER) needs to be edited in order
*   to associate a mail group with the server. The mail group
*   to add is 'VAQ PDX ERRORS'. The installer will also have
*   to change the server action from 'QUEUE SERVER ROUTINE' to
*   'RUN IMMEDIATELY'
*****
```

In order to effectively use PDX, the following terminal type attributes must be defined for each terminal type used. Please verify these attributes against the Terminal Type file at your facility.

Attribute	Value for a VT series terminal
-----	-----
Form Feed	#, \$C(27,91,50,74,27,91,72)
XY CRT	W \$C(27,91)_(DY+1)_\$C(59)_(DX+1)_\$C(72)
Erase to End of Page	\$C(27,91,74)
Insert Line	\$C(27)_"[1L"
Underline On	\$C(27,91,52,109)
Underline Off	\$C(27,91,109)
High Intensity	\$C(27,91,49,109)
Normal Intensity	\$C(27,91,109)
Save Cursor Position	\$C(27,55)
Restore Cursor Position	\$C(27,56)
Set Top/Bottom Margin	\$C(27,91)_(+IOTM)_\$C(59)_(+IOBM)_\$C(114)

Post init process completed
>

File Descriptions - Attachment D

- 1) **VAQ - PARAMETER file (#394.81)**
This file contains site specific information concerning the use of PDX. This file is initialized from V. 1.0 data
- 2) **VAQ - ENCRYPTED FIELDS file (#394.73)**
This file contains all fields that should be encrypted when PDX Requests and Unsolicited PDXs are transmitted by the facility. This file is only relevant when encryption has been turned on. The installer has the option to populate this file at this time. Entries are not required.
- 3) **VAQ - RELEASE GROUP file (#394.82)**
This file contains the facilities that have been granted "automatic processing". In order for a request to be automatically processed, the requesting facility must have an entry in this file. This file has been populated from the entries found in V. 1.0. The init also allows the installer to make changes to this file at this time.
- 4) **VAQ - OUTGOING GROUP file (#394.83)**
This file contains groups of facilities commonly accessed by PDX. This will work much like mail groups. When requesting a PDX, you could access an outgoing group called G.BETA. This would be the equivalent of entering in the beta sites individually.
- 5) **VAQ - SEGMENT GROUP file (#394.84)**
This file contains groups of data segments commonly referenced by the facility. Groups marked as PUBLIC may be referenced by all users of PDX. Groups marked as PRIVATE may only be referenced by the individual that created the group. This works the same as mail groups and outgoing groups (i.e., G.CLINICAL - this could contain any combination of the Health Summary components).

File Conversion - Attachment E

If the conversion does not run to completion, the entry point `TASK^VAQPST20` should be used to requeue the conversion. This is the same entry point used by the post-init, and all prompts and responses will be consistent with those previously mentioned.

If you would prefer to run the conversion in the foreground (i.e., from your terminal), the entry point `INTER^VAQPST20` should be used. If you do not want information concerning the conversion printed to the current device, the entry point `NONINTER^VAQPST20` should be used.

If you choose not to run the conversion, the data contained in the V. 1.0 files should still be deleted. To do this, the entry point `DELETE^VAQPST24` should be used. If you want information concerning the deletion printed to the current device, a parameter of "1" should be passed.

Sample Calls From Programmer Mode

D TASK^VAQPST20	- prompts for device and queues conversion (choose the HOME device if silent mode is desired)
D NONINTER^VAQPST20	- runs conversion in silent mode (nothing printed)
D INTER^VAQPST20	- runs conversion using current device
D DELETE^VAQPST24()	- deletes entries in V. 1.0 files in silent mode (nothing printed)
D DELETE^VAQPST24(1)	- deletes entries in V. 1.0 files using current device

Department of Veterans Affairs
Decentralized Hospital Computer Program

PDX
PATIENT DATA EXCHANGE
RELEASE NOTES &
INSTALLATION GUIDE

Version 1.5
November 1993

Information Systems Center
Albany, New York

Table of Contents

Release Notes.....	1
User Release Notes.....	1
Functionality Changes.....	1
Technical Release Notes.....	5
Files	5
Routines.....	6
Templates.....	7
Options	8
Mail Groups.....	9
Bulletins	9
Installation Guide.....	10
Introduction.....	10
Package Integration.....	10
Package Initialization/Conversion Time Estimates	11
Resource Requirements.....	12
Pre-Installation	13
Installation	14
Attachment A - Version 1.0 Routines.....	19
Attachment B - Version 1.5 Routines.....	20
Attachment C - PDX V. 1.5 Sample Installation	21
Attachment D - File Descriptions	29
Attachment E - File Conversion	30