

# CPRS PATCH OR\*3.0\*296 (CPRS GUI V.27.90) AND ASSOCIATED PATCHES INSTALLATION GUIDE

**AUGUST 2009** 

# **Table of Contents**

Table of Contents	. 2
CPRS GUI v27.90 Installation Requirements and Notes	. 3
Required Patches	
Related Patches	
Internet Explorer Requirements	. 3
Installation Compliance Date	
Local Modifications to OR*3.0*296 Source Code	
Software Retrieval	
OR*3.0*296 (CPRS GUI v27.90) and Associated Patches Software Installation	
Instructions	. 4
OR*3.0*296 (CPRS GUI v27.90) and Associated Patches Software Installation	
Sequence	. 5
OR*3.0*296 Changed Help text of Existing Parameter – ORPF ACTIVE ORDERS	
CONTEXT HRS	. 6
OR*3.0*296 (CPRS GUI v27.90) Associated Software and Documentation	
Required Patches:	
Related Patches:	
OR 30 296.ZIP:	
OR*3*296 Installation Capture	13
GMRC*3*65 Installation Capture	15

# OR\*3.0\*296 (CPRS GUI v27.90) Installation Requirements and Notes

OR\*3.0\*296 includes one host file that exports routines, two options, three Remote Procedure Calls (RPCs), and parameters. No data dictionary changes are included.

Required Patches
Before you can install OR*3.0*296 (CPRS GUI v.27.90), you must verify that the
following required patches are properly installed on your system:
□ OR*3.0*283
□ OR*3.0*292
□ OR*3.0*301
□ OR*3.0*303
□ OR*3.0*304
□ YS*5.01*98
Related Patches
To fully implement the new features in OR*3.0*296, the following related patches must
be installed on your system:
☐ GMRC*3.0*65
☐ GMRC*3.0*69
Internet Explorer Requirements
Internet Explorer v4.0 (IE4) or later is REQUIRED in order for OR*3.0*296 to run.
However, IE5.5 or later is required for PKI functionality.
Installation Compliance Date
The compliance period for installation of OR*3.0*296 is <b>30</b> days.
The compliance period for installation of oil 3.0 270 is ev days.
Local Modifications to OR*3.0*296 Source Code
OR*3.0*296 is compiled with Delphi 2006. Modification of OR*3.0*296 Source Code requires Delphi 2006.

#### **Software Retrieval**

Download the necessary files from the FTP anonymous directories: download.vista.med.va.gov.

- YS 501 98.KID
- YS MHA A.DLL and YS 501 98 SETUP.EXE
- OR 30 296.KID
- OR 30 296.ZIP

Software distributed via the National Patch Module (FORUM) in a packman message.

- GMRC\*3.0\*65
- GMRC\*3.0\*69

# OR\*3.0\*296 (CPRS GUI v27.90) and Associated Patches Software Installation Instructions

This patch should be installed during non-peak hours to minimize disruption. No users should be accessing CPRS during the install as this prevents you from replacing CPRSChart.exe. Installation will take approximately 15 minutes or more, depending on the menu structure at your site.

Listed below are general installation instructions for installing the OR\*3.0\*296 KIDS build and other patches. For specific installation details, refer to the software installation captures below. Additionally, review the National Patch Module messages for each patch for patch-specific information.

- 1. From the Kernel Installation and Distribution System (KIDS) Menu, select the Installation menu.
- 2. Use Load a Distribution. You may need to precede the host file name with a directory name.
- 3. If given the option to run any Environment Check Routine(s), answer "YES."
- 4. From this menu, you may then elect to use the following options:
  - Backup a Transport Global
  - o Compare Transport Global to Current System
  - Verify Checksums in Transport Global
- 5. When ready, select the Install Packages option.
- 6. When prompted "Want KIDS to Rebuild Menu Trees Upon Completion of Install? Yes//", respond "YES."
- 7. When prompted 'Want KIDS to INHIBIT LOGONs during the install? YES//', respond "YES."

8. When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? YES//', respond YES. When prompted to select the options you would like to place out of order, enter the following:

OR OE/RR MENU CLINICIAN CPRS Clinician Menu
OR OE/RR MENU NURSE CPRS Nurse Menu

OR OE/RR MENU WARD CLERK CPRS Ward Clerk Menu

OR CPRS GUI CHART

9. When prompted 'Delay Install (Minutes): (0-60): 0//; respond '0.'

# OR\*3.0\*296 (CPRS GUI v27.90) and Associated Patches Software Installation Sequence

- 1. Make sure all Nationally Released Required patches are installed
  - o OR\*3\*283
  - o OR\*3\*292
  - o OR\*3\*301
  - o OR\*3\*303
  - o OR\*3\*304
- 2. Install Mental Health required patch and file:
  - YS\*5.01\*98.KID
  - YS\_MHA\_A.DLL (using YS\_501\_98\_SETUP.EXE if appropriate, see step c below)
  - a.) Install YS\_501\_98.KID
  - b.) Ensure that the option YS BROKER1 [YS BROKER1] is on the Mental Health users' **secondary menus**, as well as the **secondary menus** of all users who currently have the OR CPRS GUI CHART option.
  - c.) Distribute new DLL named YS\_MHA\_A.DLL in one of the following ways:
    - Install manually on the client workstation: Run YS\_501\_98\_Setup.exe to place the YS MHA A.dll in the Common Files directory.

#### OR

Pushing files to the workstation over the network: If your site uses application distribution software, you can distribute the YS\_MHA\_A.dll to client workstations. Running the YS\_501\_98\_Setup.exe requires two interactions with the user (to start and to signify finished). To avoid the user/install interaction, the file YS\_MHA\_A.dll can be copied directly to the <drive>\Programs Files\Vista\Common Files directory. This may be the more desirable way to push the dll to many workstations.

Please refer to the YS\*5.01\*98 Installation Capture below in this document for installation details.

3. Install Patch OR\*3.0\*296: OR\_30\_296.KID.

Please refer to the OR\*3.0\*296 Installation Capture for installation details.

4. Distribute CPRSCHART.EXE, CPRS.HLP, CPRS.CNT located in OR 30 296.ZIP.

These three files should be distributed to the same directory as the borlndmm.dll. Usually, these files are located in the CPRS directory.

- 5. Place the User and Technical Manual files in a location that can be accessed by CPRS users.
- 6. Install Patch GMRC\*3.0\*65.

Please refer to the GMRC\*3.0\*65 Installation Capture for installation details.

7. Install Patch GMRC\*3.0\*69.

Please refer to the GMRC\*3.0\*69 Installation Capture for installation details.

# OR\*3.0\*296 Changed Help text of Existing Parameter – ORPF ACTIVE ORDERS CONTEXT HRS

• This parameter defines the number of hours that orders remain in the "Active/Current Orders" context after they have been completed.

Updated Help Values:

Before:

VALUE HELP: Number of hours to include completed orders in Active Orders display.

After:

VALUE HELP: Enter the number of hours to include terminated orders in the Active Orders view

# OR\*3.0\*296 (CPRS GUI v27.90) Associated Software and **Documentation**

**Required Patches:** 

PATCH	SUPPORTED FUNCTIONALITY		
OR*3.0*283	Event Delayed Orders		
	Releases Two Copies		
OR*3.0*292	Locking and Parameter Help Text Problem		
OR*3*301	Changes the ORB PROCESSING FLAG Parameter Updates CPRS v27 Release Notes and Technical Manual		
	Adds new option, Convert IV Inpatient QO to Infusion QO [OR CONVERT INP TO QO]		
OR*3*303	Corrects Display of DCd Med OrdersRemoves CR from Reason for Study fieldCorrects "Invalid Rate" error in Renewed OrdersFixes <subscript>GUI+13^ORCDLR2</subscript>		
YS*5.01*98	Distributes new MH DLL		
Exports three files:	Stores Temporary file in VistA not on Client		
YS_501_98.KID YS_501_98_SETUP.EXE YS_MHA_A.DLL			

#### **Related Patches:**

РАТСН	SUPPORTED FUNCTIONALITY
GMRC*3*65	GWOT COMBAT VET INDICATORS
GMRC*3*69	REFORMATS LONG COMMENTS

### OR\_30\_296.ZIP:

FILE	VERSION	CONTENTS
CPRSChart.exe	1.0.27.90	CPRS GUI v27.90
CPRSGUIUM.doc	N/A	CPRS GUI User Manual
CPRSLMTM.doc	N/A	CPRS ListManager
		Technical Manual
CPRSGUITM.doc	N/A	CPRS GUI Technical Manual
OR_30_296IG.doc	N/A	CPRS Patch OR*3*296 (CPRS
		GUI V.27.90) and Associated
		Patches Installation Guide
OR_30_296RN.doc	N/A	OR*3.0*296 Release Notes
CPRS.cnt	N/A	CPRS help contents file
CPRS.hlp	N/A	CPRS help file
CPRSGUIUM.pdf	N/A	CPRS GUI User Manual
CPRSLMTM.pdf	N/A	CPRS ListManager
		Technical Manual
CPRSGUITM.pdf	N/A	CPRS GUI Technical Manual
OR_30_296ig.pdf	N/A	CPRS Patch OR*3*296 (CPRS GUI
		V.27.90) and Associated Patches
		Installation Guide
OR_30_296RN.pdf	N/A	OR*3.0*296 Release Notes

## YS\*5.01\*98 Installation Capture

1. Install the YS\*5.01\*98.KID patch. See example capture below.

Example

```
1
          Load a Distribution
   2
          Verify Checksums in Transport Global
          Print Transport Global
  3
   4
          Compare Transport Global to Current System
         Backup a Transport Global
   5
         Install Package(s)
   6
          Restart Install of Package(s)
         Unload a Distribution
Select Installation Option: 6 Install Package(s)
Select INSTALL NAME: YS*5.01*98
Select Installation Option: 1 Load a Distribution
Enter a Host File: [LEXPLT]YS 501 98.KID;1
KIDS Distribution saved on Aug 28, 2008@11:10:33
Comment: CPRS V27 SUPPORT
This Distribution contains Transport Globals for the following Package(s):
  YS*5.01*98
Distribution OK!
Want to Continue with Load? YES//
Loading Distribution...
  YS*5.01*98
Use INSTALL NAME: YS*5.01*98 to install this Distribution.
  1
         Load a Distribution
  2
         Verify Checksums in Transport Global
         Print Transport Global
         Compare Transport Global to Current System
  5
         Backup a Transport Global
          Install Package(s)
   6
         Restart Install of Package(s)
          Unload a Distribution
Select Installation Option: 6 Install Package(s)
Select INSTALL NAME: YS*5.01*98
                                  Loaded from Distribution 6/22/09@11:52:35
     => CPRS V27 SUPPORT ;Created on Aug 28, 2008@11:10:33
This Distribution was loaded on Jun 22, 2009@11:52:35 with header of
  CPRS V27 SUPPORT ;Created on Aug 28, 2008@11:10:33
   It consisted of the following Install(s):
    YS*5.01*98
Checking Install for Package YS*5.01*98
Install Questions for YS*5.01*98
Incoming Files:
   601.94 MH CR SCRATCH
Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO// NO
Want KIDS to INHIBIT LOGONs during the install? NO// NO
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO// NO
Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.
DEVICE: HOME// TELNET TERMINAL
Complete
```

- 2. Ensure that the option YS BROKER1 [YS BROKER1] is on the Mental Health users' **secondary menus**, as well as the **secondary menus** of all users who currently have the OR CPRS GUI CHART option.
- 3. After installation of the KIDS build, download the YS\_501\_98\_Setup.exe and YS MHA A.dll file located in the Anonymous Software directory.
- 4. Distribute new DLL named YS MHA A.DLL in one of the following ways:
  - o Install manually on the client workstation: Run YS\_501\_98\_Setup.exe to place the YS\_MHA\_A.dll in the Common Files directory.

**NOTE:** Do NOT queue the installation, because this installation contains questions requiring a response, and queuing will stop the installation.

#### OR

Pushing files to the workstation over the network: If your site uses application distribution software, you can distribute the YS\_MHA\_A.dll to client workstations. Running the YS\_501\_98\_Setup.exe requires two interactions with the user (to start and to signify finished). To avoid the user/install interaction, the file YS\_MHA\_A.dll can be copied directly to the <drive>\Programs Files\Vista\Common Files directory. This may be the more desirable way to push the dll to many workstations.

Please refer to the YS\*5.01\*98 Installation Capture below in this document and the *MENTAL HEALTH ASSISTANT DLL used by CPRS Clinical Reminders INSTALLATION GUIDE: PATCH YS\*5.01\*98* for installation details.

#### OR\*3.0\*296 Checksum Value Verification

```
Load a Distribution
          Verify Checksums in Transport Global
   3
         Print Transport Global
          Compare Transport Global to Current System
   4
         Backup a Transport Global
   5
         Install Package(s)
   6
          Restart Install of Package(s)
          Unload a Distribution
Select Installation Option: 2 Verify Checksums in Transport Global
Select INSTALL NAME: OR*3.0*296
PACKAGE: OR*3.0*296 Jun 22, 2009 11:11 am
______
OR3CONV Calculated 57296998
OR3CONV1 Calculated 41773160
ORB3 Calculated 93373249
ORCDLR1 Calculated 81262557
ORCDPSIV Calculated 96653807
ORCMEDTO Calculated 17670670
ORCMEDT1 Calculated 45036674
ORCMEDT5 Calculated 35418320 ORCSEND Calculated 65578810
ORCSEND Calculated 65578810
ORINQIV Calculated 192421660
ORMARKER Calculated 978608
ORMRA
          Calculated 62138554
Calculated 64231033
ORPRPM
ORWDX Calculated 67409395
ORWDX1 Calculated 52053311
ORWDXM3 Calculated 72094070
ORWDXM4 Calculated 32166098
ORWOD Calculated 72302805
ORWOD1 Calculated 38695937
ORWORB Calculated 69468448
          Calculated
ORWTPUA Calculated 2236253
ORY296 Calculated 1180350
ORYCHKE Calculated 2490043
   23 Routines checked, 0 failed.
```

## **OR\*3.0\*296 Transport Global Print**

```
Load a Distribution
   2
          Verify Checksums in Transport Global
   3
           Print Transport Global
   4
           Compare Transport Global to Current System
   5
          Backup a Transport Global
   6
          Install Package(s)
          Restart Install of Package(s)
           Unload a Distribution
Select Installation Option: 3 Print Transport Global
Select INSTALL NAME: OR*3.0*296
PACKAGE: OR*3.0*296 Jun 22, 2009 11:13 am
                                                                        PAGE 1
SINGLE PACKAGE
                              TRACK NATIONALLY: YES
NATIONAL PACKAGE: ORDER ENTRY/RESULTS REPORTING ALPHA/BETA TESTING: NO
DESCRIPTION:
ENVIRONMENT CHECK: DELETE ENV ROUTINE: PRE-INIT ROUTINE: PRE-ORY296 DELETE PRE-INIT ROUTINE: NO POST-INIT ROUTINE: POST-ORY296 DELETE POST-INIT ROUTINE: NO
PRE-TRANSPORT RTN:
                                          UP
                                                SEND DATA
                                          DATE SEC. COMES SITE RSLV OVER
FILE # FILE NAME
                                         DD CODE W/FILE DATA PTRS RIDE
                                          NO NO YES OVER YES
101.41 ORDER DIALOG
 DATA SCREEN: I Y=130
ROUTINE:
                                              ACTION:
  OR3CONV
                                                 SEND TO SITE
  OR3CONV1
                                                  SEND TO SITE
  ORB3
                                                 SEND TO SITE
                                                 SEND TO SITE
  ORCDLR1
  ORCDESTV
                                                  SEND TO SITE
  ORCMEDT0
                                                  SEND TO SITE
                                                  SEND TO SITE
  ORCMEDT1
  ORCMEDT5
                                                  SEND TO SITE
  ORCSEND
                                                  SEND TO SITE
  ORINOIV
                                                  SEND TO SITE
  ORMARKER
                                                  SEND TO SITE
  ORMRA
                                                  SEND TO SITE
  ORPRPM
                                                  SEND TO SITE
                                                  SEND TO SITE
  ORWDX
  ORWDX1
                                                  SEND TO SITE
  ORWDXM3
                                                  SEND TO SITE
  ORWDXM4
                                                  SEND TO SITE
                                                  SEND TO SITE
  ORWOD1
                                                  SEND TO SITE
  ORWORB
                                                  SEND TO SITE
                                                  SEND TO SITE
  ORWTPUA
  ORY243
                                                  DELETE AT SITE
  ORY2430
                                                  DELETE AT SITE
  ORY24301
                                                  DELETE AT SITE
  ORY24302
                                                 DELETE AT SITE
  ORY24303
                                                  DELETE AT SITE
  ORY24304
                                                  DELETE AT SITE
  ORY24305
                                                 DELETE AT SITE
  ORY24306
                                                  DELETE AT SITE
  ORY24307
                                                  DELETE AT SITE
  ORY24308
                                                 DELETE AT SITE
  ORY24309
                                                 DELETE AT SITE
  ORY2431
                                                  DELETE AT SITE
  ORY2432
                                                  DELETE AT SITE
```

```
ORY2433
                                                     DELETE AT SITE
  ORY2434
                                                     DELETE AT SITE
   ORY243A
                                                     DELETE AT SITE
                                                     DELETE AT SITE
  ORY243ES
  ORY243R
                                                     DELETE AT SITE
  ORY296
                                                     SEND TO SITE
  ORYCHKE
                                                     SEND TO SITE
                                                 ACTION:
OPTION:
  OR CPRS GUI CHART
                                                    SEND TO SITE
                                                    SEND TO SITE
  OR MEDICATION QO CHECKER
PARAMETER DEFINITION:
                                                 ACTION:
  ORPF ACTIVE ORDERS CONTEXT HRS
                                                    SEND TO SITE
  ORWRP TIME/OCC LIMITS INDV
                                                    SEND TO SITE
                                                 ACTION:
REMOTE PROCEDURE:
  GMV CLOSEST READING
                                                    SEND TO SITE
  OR GET COMBAT VET
                                                     SEND TO SITE
  ORWDX UNLKOTH
                                                    SEND TO SITE
INSTALL QUESTIONS:
 Default Rebuild Menu Trees Upon Completion of Install: NO
Default INHIBIT LOGONs during the install: NO
Default DISABLE Scheduled Options, Menu Options, and Protocols: NO
REQUIRED BUILDS:
                                                 ACTION:
  YS*5.01*98
                                                     Don't install, leave global
                                                     Don't install, leave global
Don't install, leave global
Don't install, leave global
  OR*3.0*292
   OR*3.0*283
  OR*3.0*301
  OR*3.0*303
                                                     Don't install, leave global
  OR*3.0*304
                                                     Don't install, leave global
```

## **OR\*3\*296 Installation Capture**

```
Select OPTION NAME: XPD MAIN
                                  Kernel Installation & Distribution System
          Edits and Distribution ...
          Utilities ...
          Installation ...
          Patch Monitor Main Menu ...
Select Kernel Installation & Distribution System Option: INSTALLation
         Load a Distribution
          Verify Checksums in Transport Global
  2
          Print Transport Global
          Compare Transport Global to Current System
          Backup a Transport Global
  5
          Install Package(s)
          Restart Install of Package(s)
          Unload a Distribution
Select Installation Option: 6 Install Package(s)
Select INSTALL NAME: or*3.0*296 Loaded from Distribution 6/22/09@10:50:58
     => OR*3.0*296 KIDS BUILD ;Created on May 07, 2009@11:54:07
This Distribution was loaded on Jun 22, 2009@10:50:58 with header of
  OR*3.0*296 KIDS BUILD ;Created on May 07, 2009@11:54:07
  It consisted of the following Install(s):
     OR*3.0*296
Checking Install for Package OR*3.0*296
Install Questions for OR*3.0*296
Incoming Files:
  101.41
            ORDER DIALOG (including data)
Note: You already have the 'ORDER DIALOG' File.
I will OVERWRITE your data with mine.
Want KIDS to Rebuild Menu Trees Upon Completion of Install? YES// YES
Want KIDS to INHIBIT LOGONs during the install? YES//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? YES//
Enter options you wish to mark as 'Out Of Order': OR OE/RR MENU CLINICIAN
CPRS Clinician Menu
Enter options you wish to mark as 'Out Of Order': OR OE/RR MENU NURSE
                                                                            CPRS
Nurse Menu
Enter options you wish to mark as 'Out Of Order': OR OE/RR MENU WARD CLERK
CPRS Ward Clerk Menu
Enter options you wish to mark as 'Out Of Order': OR CPRS GUI CHART
                                                                          CPRSCh
art version 1.0.27.83
Enter options you wish to mark as 'Out Of Order':
Enter protocols you wish to mark as 'Out Of Order':
Delay Install (Minutes): (0-60): 0//
Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.
DEVICE: HOME// TELNET TERMINAL
                                   OR*3.0*296
```

```
Installing OPTION
 Installing PARAMETER DEFINITION
                 Jun 22, 2009@14:05:48
 Running Post-Install Routine: POST^ORY296
 Updating Routine file...
 Updating KIDS files...
Menu Rebuild Complete: Jun 22, 2009@11:53:46
 OR*3.0*296 Installed.
                Jun 22, 2009@14:05:48
 Not a production UCI
 NO Install Message sent
                                    Rebuilding Menus

        WVMENU
        Women's Health Menu
        1
        03/21/06
        05/27/09

        LRLOINC
        LOINC Main Menu
        1
        02/27/07
        05/27/09

        RMPO-MENU-MAIN
        Home Oxygen Main Menu
        2
        03/05/07
        05/27/09

        ZPSJU DENVER RPH
        Inpatient Medication Menu...
        29
        03/05/07
        05/27/09

        452 CONTRACT TRANSCRIPTION

452 CONTRACT TRANSCRIPTION
Building secondary menu trees....
Merging.... done.
           R_{1},...,T_{1} . 25 50 75 .
  100%
Install Completed
```

Page 14

# **GMRC\*3\*65 Installation Capture**

```
Load a Distribution
  2
       Verify Checksums in Transport Global
  3
       Print Transport Global
       Compare Transport Global to Current System
  4
      Backup a Transport Global
       Install Package(s)
  6
        Restart Install of Package(s)
       Unload a Distribution
You have 6 new messages. (Last arrival: 08/06/08@12:46)
Select Installation Option: 6 Install Package(s)
Select INSTALL NAME: GMRC*3.0*65
                             Loaded from Distribution 6/23/09@10:01:38
   => GMRC*3*65
This Distribution was loaded on Jun 23, 2009@10:01:38 with header of
  GMRC*3*65
  It consisted of the following Install(s):
   GMRC*3.0*65
Checking Install for Package GMRC*3.0*65
Install Questions for GMRC*3.0*65
Want KIDS to INHIBIT LOGONs during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//
Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.
DEVICE: HOME// TELNET TERMINAL
                           GMRC*3.0*65
......
Install Started for GMRC*3.0*65 :
            Jun 23, 2009@10:02:34
Build Distribution Date: Sep 15, 2008
Installing Routines:
            Jun 23, 2009@10:02:34
Updating Routine file...
Updating KIDS files...
GMRC*3.0*65 Installed.
            Jun 23, 2009@10:02:35
Not a production UCI
NO Install Message sent
R,,,,,,T
. 25 50 75 .
Install Completed
```

## **GMRC\*3\*69 Installation Capture**

```
Load a Distribution
   2
          Verify Checksums in Transport Global
   3
          Print Transport Global
          Compare Transport Global to Current System
  5
         Backup a Transport Global
   6
          Install Package(s)
         Restart Install of Package(s)
         Unload a Distribution
Select Installation Option: 6 Install Package(s)
Select INSTALL NAME: GMRC*3.0*69
                                      Loaded from Distribution 3/11/09@15:24:5
     => GMRC*3*69
This Distribution was loaded on Mar 11, 2009@15:24:52 with header of
  GMRC*3*69
   It consisted of the following Install(s):
   GMRC*3.0*69
Checking Install for Package GMRC*3.0*69
Install Questions for GMRC*3.0*69
Want KIDS to INHIBIT LOGONs during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//
Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.
DEVICE: HOME//
               HOME
Install Started for GMRC*3.0*69 :
              Mar 11, 2009@15:25:45
Build Distribution Date: Mar 11, 2009
Installing Routines:
              Mar 11, 2009@15:25:45
Running Pre-Install Routine: PRE^GMRCYP69
Running Post-Install Routine: POST^GMRCYP69
Post-install queued as task #322367
If any records are locked, task will re-run 5 minutes after completion.
 Updating Routine file...
Updating KIDS files...
 GMRC*3.0*69 Installed.
              Mar 11, 2009@15:25:45
Not a production UCI
NO Install Message sent
Install Completed
```