Computerized Patient Record System (CPRS) Version 31b Follow-Up Build

Setup and Configuration Guide



January 2022

Department of Veterans Affairs Office of Information & Technology (OI&T) Enterprise Program Management Office (EPMO)

Revision History

Date	Version	Description	Author		
1/13/2022	0.17	OR*3*498: Redacted for the VDL: removed names and links to software builds	CPRS Development Team		
1/11/2022	0.16	OR*3*498: Updated Section 5.7 (Resolve Any Object Conflicts).	CPRS Development Team		
1/10/2022	0.15	OR*3*498: Added Section 5.7 (Resolve Any Object Conflicts).	CPRS Development Team		
12/9/2021	0.14	OR*3.0*498: More updates to Step #2 in Section 5.2 (Gather SMART Dialog Information).	CPRS Development Team		
12/8/2021	0.13	OR*3.0*498: Updated Step #1 in Section 5.2 (Gather SMART Dialog Information).	CPRS Development Team		
11/10/2021	0.12	OR*3.0*498: Updated the following sections:	CPRS Development Team		
		5.2 (Gather SMART Dialog Information), Step #1; 6.1.1 (Install SMART UPDATES Content), Step #18; 6.1.2 (Install High Risk Medications Content, Step #5			
11/10/2021	0.11	OR*3.0*498: Added Section 5.1 (Verify that Sites Have Installed Clinical Reminder Update_2_0_210 WH Mammogram Screening Reminder). Fixed a typo in Section 5.2 (Gather SMART Term Information). Updated Section 5.2 (Install SMART	CPRS Development Team		
		Updates Content): Removed the note at the beginning of the section; Updated instruction #18.			
10/06/2021	0.10	OR*3.0*498: Updated instruction 5.b ("select the Merge action") in Section 6.1.2 (Install High Risk Medications Content).	CPRS Development Team		
09/01/2021	0.09	OR*3.0*498: Updated the following sections: Section 6 - Post Installation Steps;	CPRS Development Team		

		Section 6.1.1 - Install SMART UPDATES Content;	
		Section 6.1.2 - Install High Risk Medications Content	
07/13/2021	0.08	OR*3.0*498: In Sections 6.1.1 and 6.1.2, changed the "Date Packed" date to 07/07/2021@09:59:53. In section 6.1.1, added a note about the Pittsburgh and Northport sites.	TW: J. Callahan Dev: A. Puleo
07/01/2021	0.07	OR*3.0*498: Fixed a typo in Section 5.2 ('Gather SMART Term Information').	TW: J. Callahan Dev: R. Ruff
		In Section 6.1 ('Install Reminder Content'), fixed instructions #12 and #18: Replaced a reference to section 5.3 with a reference to section 5.2.	
		Removed Section 5.6 ('Resolve TIU Object Naming Conflicts').	
05/20/2021	0.06	OR*3.0*498: In the Install SMART UPDATES Content section, updated the instructions that pertained to tomosynthesis.	TW: J. Callahan Dev: R. Ruff
05/11/2021	0.05	OR*3.0*498: Updated the following sections: Gather SMART Team Information, Backup the SMART Content, Back-Up High-Risk Medications for Women Content, Install SMART UPDATES Content, Install High Risk Medications Content	TW: J. Callahan Dev: R. Ruff
04/13/2021	0.04	OR*3.0*498: Updated the Pre- and Post- Installation Steps sections.	TW: J. Callahan Dev: R. Ruff
03/09/2021	0.03	OR*3.0*498: Updated the Overview, Reporting Issues, Pre- and Post- Installation Checklist, Pre-Installation Steps and Post-Installation Steps sections.	TW: J. Callahan Dev: R. Ruff
01/11/2021	0.02	OR*3.0_498: Fixed a typo in Section 5: Pre-Installation Steps. Changed "VA- WH SMART BREAST IMAGING DIALOG" to "VA-WH SMART BREAST IMAGING FOLLOW-UP".	TW: J. Callahan Dev: R. Ruff
12/20/2020	0.01	OR*3.0*498: Initial Draft	TW: J. Callahan Dev: R. Ruff

Table of Contents

1	CF	PRS v31b Follow-Up Build1
	1.1	Overview1
	1.2	Recommended Audience1
	1.3	About this Guide1
	1.4	Document Conventions1
	1.5	Related Documents2
2	Pr	e-requisites2
	2.1	Pre-requisite Patches2
3	Re	porting Issues
4	Pr	e- and Post-Installation Checklist
5	Pr	e-Installation Steps
	5.1 Mam	Verify that Sites Have Installed Clinical Reminder Update_2_0_210 WH mogram Screening Reminder
	5.2	Gather SMART Dialog Information3
	5.3	Gather SMART Term Information4
	5.4	Gather Pregnancy Dialog Information4
	5.5	Backup the SMART Content5
	5.6	Back-Up the High-Risk Medications for Women Content
	5.7	Resolve Any Object Conflicts18
6	Ро	st-Installation Steps
	6.1	Install Reminder Content20
	6.1	.1 Install SMART UPDATES Content
	6.1	.2 Install High Risk Medications Content
	6.2	Enable Clinical Reminder Order Check Rules42
	6.3	Link Reminder Terms to Women's Health Procedure Type Entries 47
7	Ve	rify Successful Installation 48

1 CPRS v31b Follow-Up Build

1.1 Overview

CPRS v31b Follow-Up Build is a multi-package build that addresses several defects identified during the deployment of the CPRS v31b series of patches.

The CPRS v31b Follow-Up build consists of the following patches:

- OR*3.0*498
- WV*1.0*26
- PXRM*2.0*71
- PSO*7.0*622
- TIU*1.0*341

1.2 Recommended Audience

This guide provides information specifically for those personnel who need to perform set up and configuration steps before and after the CPRS v31b Follow-Up Build installation. These groups include Information Technology Operations and Support (ITOPS) staff, Clinical Application Coordinator (CAC) personnel, the site's Women's Health group, and others who will be needed so that the CPRS v31b Follow-Up Build will work correctly at sites.

1.3 About this Guide

This set up/configuration guide provides instructions for:

- Pre-installation steps, which must be performed before the CPRS v31b Follow-Up Build installation can proceed.
- Post-installation tasks—including configuration tasks—that require knowledge of the underlying VistA system.

1.4 Document Conventions

Examples of VistA "Roll and Scroll" interface actions will be shown in a box such as this:

Select OPTION NAME: XPAR EDIT PARAMETER Edit Parameter Values

Emphasis of important points may be displayed in this manner:

NOTE: This is an important point and must not be omitted.

Call-outs may be used to draw attention to part of a block of text or a table without disrupting the flow of the block or table. For example:

-	°
Sample text, Sample text, Sample text, Samp	informational
Sample text, Sample text	│ text about the │ point of note
Sample text, Sample text, Sample text, Samp Sample text	ole text
Sample text, Sample text, Sample text, SOM	METHING OF NOTE !!
Sample text, Sample text, Sample text, Samp	ole text, Sample text, Sample text

1.5 Related Documents

The following documents, in addition to this document, will be available on the VA Software Document Library (VDL) when the patch is released:

CPRS on the VDL

- CPRS User Guide: GUI Version
- CPRS Technical Manual
- CPRS Technical Manual: GUI Version
- CPRS Release Notes: v31b Follow-Up Build
- CPRS v31b Follow-Up Build Deployment, Installation, Back Out and Rollback Guide

2 Pre-requisites

N/A

2.1 Pre-requisite Patches

N/A

3 Reporting Issues

To report issues with CPRS v31b Follow-Up build, please enter a ticket with the National Help Desk.

4 Pre- and Post-Installation Checklist

The following activities should be completed in order. Use this checklist and the following sections for both your test/mirror system as well as your production system.

No.	Task	Done
1.	Complete the pre-installation steps. (See Section 5.)	
2.	Notify your ITOPS regional installer that your site is ready for installation.	
3.	Make sure that the installation was completed. (See the "CPRS v31b Follow-Up Build Deployment, Installation, Back Out and Rollback Guide")	
4.	Complete the post-installation steps. (See Section 6.)	
5.	Verify the installation was successful. (See Section 7.)	

5 Pre-Installation Steps

These pre-installation steps must be performed by the CAC.

5.1 Verify that Sites Have Installed Clinical Reminder Update_2_0_210 WH Mammogram Screening Reminder

Sites need to install Clinical Reminder Update_2_0_210 WH Mammogram Screening Reminder.

The Install Guide and PowerPoint presentations pertaining to the Clinical Reminder Update for both the clinical staff and clinical informatics staff are available in REDACTED.

5.2 Gather SMART Dialog Information

- 1. Use option Inquire about Reminder Taxonomy Management [PXRM TAXONOMY MANAGEMENT] to capture an inquire of the terms listed below.
 - a. VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - b. VA-WH MAMMOGRAM SCREENING CODES
 - c. VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - d. VA-WH MRI BREAST PROCEDURE CODES
 - e. VA-WH ULTRASOUND BREAST PROCEDURE CODES
- 2. Use option Reminder Dialogs [PXRM DIALOG/COMPONENT EDIT] to record the orders mapped to the VA-WH SMART BREAST IMAGING FOLLOW-UP reminder dialog. The Site will need to re-map these orders to the reminder dialog after the install.

5.3 Gather SMART Term Information

- 1. Use option Inquire about Reminder Term [PXRM TERM INQUIRY] to capture a copy of the following terms:
 - a. VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - b. VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - c. VA-WH ULTRASOUND OF THE BREAST CODES
 - d. VA-WH MRI OF THE BREASTS CODES
 - e. VA-WH MAMMOGRAM SCREENING CODES
- 2. Use option Inquire about Reminder Term [PXRM TERM INQUIRY] to review the VA-WH NEXT BREAST PROCEDURE reminder term. From the output, note the condition statement for the VA-WH NEXT BREAST PROCEDURE finding.

5.4 Gather Pregnancy Dialog Information

- 1. Use option Reminder Dialogs [PXRM DIALOG/COMPONENT EDIT] to record the findings mapped to the elements listed below. You are prompted during the content installation for these findings.
 - a. VA-WH TD PREGNANCY ORDER EMERGENCY CONTRACEPTIVE
 - b. VA-WH TD PREGNANCY ORDER PREGNANCY TEST

5.5 Backup the SMART Content

Pack the existing SMART content into a new clinical reminders exchange entry in case a software rollback is needed:

1. On the Clinical Reminders Exchange screen, at the "Select Action" prompt, enter "CFE" for Create Exchange File Entry.

Clin	ica	al Reminder Exchange May	, 11,	2021009	9:30:12	Page:	1 of	61
Exch	ang	ge File Entries.						_
Ite	m	Entry		Source	9	Date	Packed	
	1	BDI II RESULT GROUP		CPRSCI	RM20@SALT LAKE CI	I 04/13	3/2004@1	5:53
	2	DEPRESSION/PTSD REMINDER 7	TERM	CPRSCH	RM21@NORTHERN CAI	L 02/20	6/201001	2:07
		UPDATES - PATCH 17						
	3	ECOE REMINDER DIALOGS		CPRSCI	RM22@SALT LAKE CI	I 08/28	8/201300	9:38
	4	GMTS FOR HRMH		CPRSCE	RM230SALT LAKE CI	I 12/01	1/201101	1:10
	5	GMTS SKIN RISK HS OBJECTS		CPRSCE	RM20@SALT LAKE CI	I 07/09	9/200701	3:21
	6	GMTS SKIN RISK HS TYPES		CPRSCH	RM20@SALT LAKE CI	E 07/09	9/200701	3:20
	7	GMTSMHV		CPRSCE	RM20@SALT LAKE CI	E 07/0	6/200401	5:06
	8	NATIONAL BLOOD PRESSURE		CPRSCE	RM24@TUSCALOOSA	04/05	5/201101	3:20
		CHANGES						
+		+ Next Screen - Prev	Scree	n ??	More Actions			>>>
CFE	C	reate Exchange File Entry		LHF	Load Host File			
CHF	C	reate Host File		LMM	Load MailMan Mes	ssage		
CMM	C	reate MailMan Message		LR	List Reminder De	efinition	ns	
DFE	De	elete Exchange File Entry		LWH	Load Web Host Fi	ile		
IFE	Ιı	stall Exchange File Entry		RI	Reminder Definit	zion Inqu	uiry	
ΙH	II	nstallation History		RP	Repack	-	-	
Sele	ct	Action: Next Screen// CFE	Cre	ate Exc	change File Entry	Y		

Example – Create Exchange File Entry Action

- 2. Add the reminder dialog VA-WH SMART BREAST IMAGING FOLLOW-UP:
 - a. At the Select a File prompt, type "4" for reminder dialog.
 - b. At the Select REMINDER DIALOG NAME prompt, type "VA-WH SMART BREAST IMAGING FOLLOW-UP".
 - c. At the ... OK? prompt, press the Enter key.
 - d. At the Select REMINDER DIALOG NAME prompt, press the Enter key.

Example – Adding a Reminder Dialog

Select from the following reminder files:	
1 REMINDER COMPUTED FINDINGS	
2 REMINDER COUNTING GROUP	
3 REMINDER DEFINITION	
4 REMINDER DIALOG	
5 REMINDER EXTRACT COUNTING RULE	
6 REMINDER EXTRACT DEFINITION	
7 REMINDER LIST RULE	
8 REMINDER LOCATION LIST	
9 REMINDER SPONSOR	
10 REMINDER TAXONOMY	
11 REMINDER TERM	
12 REMINDER ORDER CHECK ITEMS GROUP	
13 REMINDER ORDER CHECK RULES	
Select a file: (1-13): 4	
Select REMINDER DIALOG NAME: VA-WH SMART BREAST IMAGING FOLLOW-UP	reminder
dialog NATIONAL	
OK? Yes// (Yes)	
Enter another one or just press enter to go back to file selection.	
Select REMINDER DIALOG NAME:	
SELECT REMINDER DIALOG NAME.	

- 3. Add the following reminder taxonomies:
 - VA-WH ULTRASOUND BREAST PROCEDURE CODES
 - VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - VA-WH MRI BREAST PROCEDURE CODES
 - VA-WH MAMMOGRAM SCREENING CODES
 - a. At the Select a File prompt, type "10" for reminder taxonomy.
 - b. At the Select REMINDER TAXONOMY NAME prompt, enter the name of the taxonomy from the list above.
 - c. At the ... OK? prompt, press the Enter key.
 - d. Repeat steps b and c for the remaining taxonomies in the list above.
 - e. At the Select REMINDER TAXONOMY NAME prompt, press the Enter key.

Example – Adding a Reminder Taxonomy

Select from the following reminder files:	
1 REMINDER COMPUTED FINDINGS	
2 REMINDER COUNTING GROUP	
3 REMINDER DEFINITION	
4 REMINDER DIALOG	
5 REMINDER EXTRACT COUNTING RULE	
6 REMINDER EXTRACT DEFINITION	
7 REMINDER LIST RULE	
8 REMINDER LOCATION LIST	
9 REMINDER SPONSOR	
10 REMINDER TAXONOMY	
11 REMINDER TERM	
12 REMINDER ORDER CHECK ITEMS GROUP	
13 REMINDER ORDER CHECK RULES	
Select a file: (1-13): 10	
Select REMINDER TAXONOMY NAME: VA-WH ULTRASOUND BREAST PROCEDURE CODES	LOC
AL	
OK? Yes// (Yes)	
Enter another one or just press enter to go back to file selection.	
Select REMINDER TAXONOMY NAME:	

- 4. Add the following reminder terms:
 - VA-WH MAMMOGRAM SCREENING CODES
 - VA-WH ULTRASOUND OF THE BREAST CODES
 - VA-WH MRI OF THE BREASTS CODES
 - VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - VA-WH NEXT BREAST PROCEDURE
 - a. At the Select a File prompt, type "11" for reminder term.
 - b. At the Select REMINDER TERM NAME prompt, enter the name of the term from the list above.
 - c. At the ... OK? prompt, press the Enter key.
 - d. Repeat steps b and c for the remaining terms in the list above.
 - e. At the Select REMINDER TERM NAME prompt, press the Enter key.

Example – Adding a Reminder Term

```
Select from the following reminder files:
  1 REMINDER COMPUTED FINDINGS
  2 REMINDER COUNTING GROUP
  3 REMINDER DEFINITION
  4 REMINDER DIALOG
  5 REMINDER EXTRACT COUNTING RULE
  6 REMINDER EXTRACT DEFINITION
  7 REMINDER LIST RULE
  8 REMINDER LOCATION LIST
  9 REMINDER SPONSOR
 10 REMINDER TAXONOMY
 11 REMINDER TERM
 12 REMINDER ORDER CHECK ITEMS GROUP
 13 REMINDER ORDER CHECK RULES
Select a file: (1-13): 11
Select REMINDER TERM NAME: VA-WH MAMMOGRAM SCREENING CODES
                                                                NATIONAL
        ...OK? Yes// (Yes)
Enter another one or just press enter to go back to file selection.
Select REMINDER TERM NAME:
```

- 5. At the Select a file prompt, press the Enter key.
- 6. At the Press ENTER to continue or ' $^{\prime}$ to exit prompt, press the Enter key.
- 7. At the Enter the Exchange File entry name prompt, type "SMART CONTENT BACKUP".

Example – Naming the Exchange Entry

Select from the following reminder files: 1 REMINDER COMPUTED FINDINGS 2 REMINDER COUNTING GROUP 3 REMINDER DEFINITION 4 REMINDER DIALOG 5 REMINDER EXTRACT COUNTING RULE 6 REMINDER EXTRACT DEFINITION 7 REMINDER LIST RULE 8 REMINDER LOCATION LIST 9 REMINDER SPONSOR 10 REMINDER TAXONOMY 11 REMINDER TERM 12 REMINDER ORDER CHECK ITEMS GROUP 13 REMINDER ORDER CHECK RULES Select a file: (1-13): Checking reminder term(s) for errors. Checking reminder term VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES No fatal term errors were found. Press ENTER to continue or '^' to exit: Checking reminder term VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES No fatal term errors were found. Checking reminder term VA-WH MRI OF THE BREASTS CODES No fatal term errors were found. Checking reminder term VA-WH ULTRASOUND OF THE BREAST CODES No fatal term errors were found. Checking reminder term VA-WH MAMMOGRAM SCREENING CODES No fatal term errors were found. Checking reminder term VA-WH NEXT BREAST PROCEDURE No fatal term errors were found. No fatal reminder term problems were found, packing will continue. Checking reminder dialog(s) for errors.. **WARNING** VA-WH SMART BREAST IMAGING FOLLOW-UP contains the following errors. The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 3 is using all possible evaluation statuses for an evaluation item. The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 2 is using all possible evaluation statuses for an evaluation item. The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 1 is using all possible evaluation statuses for an evaluation item. The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 0 is using all possible evaluation statuses for an evaluation item. The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 6 is using all possible evaluation statuses for an evaluation item. The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 5 is using all possible evaluation statuses for an evaluation item.

The dialog group VA-WH GP BR BI-RAD SELECTOR: BIRAD 4 is using all possible evaluation statuses for an evaluation item. The dialog element VA-WH SMART BR MALE ALERT RESULT is using all possible evaluation statuses for an evaluation item.

- Enter the Exchange File entry name: SMART CONTENT BACKUP
- 8. You will now see the Text Editor. A description for the exchange entry is auto generated and appears in the editor. You are not required to edit this description. Press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.

Example – Adding a Description

```
Enter a description of the Exchange File entry you are packing.....
. . . . . . . . . . . . . . . . . . .
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
 VA-WH SMART BREAST IMAGING FOLLOW-UP
REMINDER TAXONOMY
 VA-WH ULTRASOUND BREAST PROCEDURE CODES
 VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 VA-WH MRI BREAST PROCEDURE CODES
 VA-WH MAMMOGRAM SCREENING CODES
REMINDER TERM
 VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 VA-WH MRI OF THE BREASTS CODES
 VA-WH ULTRASOUND OF THE BREAST CODES
 VA-WH MAMMOGRAM SCREENING CODES
 VA-WH NEXT BREAST PROCEDURE
```

9. You will now see the Text Editor a second time. You are not required to enter any keywords or phrases. Press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.

Example – Adding Keywords

- 10. The exchange entry is now packed.
- 11. You are returned to the Clinical Reminders Exchange screen. Note at the top of the screen that the entry you created was saved.

	Clinical Reminder Exchange May 11, 2021@09:34:52 Page: 1 of 61							
-	I CONTENT BACKUP was saved in the E		1000. 101 01					
Iter	n Entry	Source	Date Packed					
	1 BDI II RESULT GROUP	CPRSCRM20@SALT LAKE	CI 04/13/2004@15:53					
2	2 DEPRESSION/PTSD REMINDER TERM	CPRSCRM21@NORTHERN C	AL 02/26/2010@12:07					
	UPDATES - PATCH 17							
	3 ECOE REMINDER DIALOGS	CPRSCRM22@SALT LAKE	CI 08/28/2013009:38					
4	4 GMTS FOR HRMH	CPRSCRM23@SALT LAKE	CI 12/01/2011@11:10					
		CPRSCRM20@SALT LAKE						
		CPRSCRM20@SALT LAKE						
		CPRSCRM20@SALT LAKE						
8		CPRSCRM24@TUSCALOOSA	04/05/2011013:20					
	CHANGES							
+	+ Next Screen - Prev Screen		>>>					
CFE	Create Exchange File Entry	LHF Load Host File						
CHF	Create Host File	LMM Load MailMan M						
CMM	Create MailMan Message	LR List Reminder						
DFE	Delete Exchange File Entry	LWH Load Web Host						
IFE	Install Exchange File Entry		ition Inquiry					
ΙH	Installation History	RP Repack						
Seled	ct Action: Next Screen//							

5.6 Back-Up the High-Risk Medications for Women Content

Pack the existing high-risk medications for women content into a new clinical reminders exchange entry in case a software rollback is needed:

1. On the Clinical Reminders Exchange screen, at the "Select Action" prompt, enter "CFE" for Create Exchange File Entry.

Clin	ica	al Rer	ninder Exchange	May 11,	202101	0:10:07	Pa	age:	1 of	61
Exch	anç	ge Fil	le Entries.							
Ite	m	Entry	7		Source	9			Packed	
	1	BDI 1	II RESULT GROUP		CPRSCI	RM20@SALT	LAKE CI	04/1	3/200401	5:53
	2	DEPRE	ESSION/PTSD REMINI	DER TERM	CPRSCI	RM21@NORTH	HERN CAL	02/2	6/201001	2:07
		UPDA?	TES - PATCH 17							
i i	3	ECOE	REMINDER DIALOGS		CPRSCI	RM22@SALT	LAKE CI	08/2	8/201300	9:38
	4	GMTS	FOR HRMH		CPRSCI	RM23@SALT	LAKE CI	12/0	1/201101	1:10
	5	GMTS	SKIN RISK HS OBJE	ECTS	CPRSCI	RM20@SALT	LAKE CI	07/0	9/2007@1	3:21
	6	GMTS	SKIN RISK HS TYPE	IS	CPRSCI	RM20@SALT	LAKE CI	07/0	9/2007@1	3:20
	7	GMTSN	1HV		CPRSCI	RM20@SALT	LAKE CI	07/0	6/2004@1	5:06
	8	NATIC	NAL BLOOD PRESSUF	₹E	CPRSCI	RM24@TUSCA	ALOOSA	04/0	5/201101	3:20
		CHANC	GES							
+		+	Next Screen - H	Prev Scree	en ??	More Act:	ions			>>>
CFE	Cı	reate	Exchange File Ent	ry	LHF	Load Host	t File			
CHF	Cı	reate	Host File		LMM	Load Mail	lMan Messa	age		
CMM	Cı	reate	MailMan Message		LR	List Rem:	inder Defi	Initio	ns	
DFE	De	elete	Exchange File Ent	ry	LWH	Load Web	Host File	9		
IFE	Ir	nstall	L Exchange File Er	ntry	RI	Reminder	Definitio	on Inq	uiry	
ΙH	Ir	nstall	lation History		RP	Repack				
Sele	ct	Actio	on: Next Screen//	CFE Cre	eate Exc	change Fil	le Entry			

Example – Create Exchange File Entry Action

- 2. Add the following reminder definitions:
 - VA-WH CHANGE IN PREGNANCY STATUS
 - VA-WH NO PREGNANCY END DOC RPT
 - VA-WH POTENTIALLY UNSAFE MEDICATIONS REPORT COHORT
 - VA-WH PREGNANCY AND LACTATION DATA ENTRY ALLOWED
 - VA-WH UPDATE LACTATION STATUS
 - VA-WH UPDATE PREGNANCY STATUS
 - a. At the Select a File prompt, type "3" for reminder definition.
 - b. At the Select REMINDER DEFINITION NAME prompt, enter the name of the definition from the list above.
 - c. Repeat step b for the remaining definitions in the list above.
 - d. At the Select REMINDER DEFINITION NAME prompt, press the Enter key.

Example – Adding a Reminder Definition

Select from the following reminder files:	
1 REMINDER COMPUTED FINDINGS	
2 REMINDER COUNTING GROUP	
3 REMINDER DEFINITION	
4 REMINDER DIALOG	
5 REMINDER EXTRACT COUNTING RULE	
6 REMINDER EXTRACT DEFINITION	
7 REMINDER LIST RULE	
8 REMINDER LOCATION LIST	
9 REMINDER SPONSOR	
10 REMINDER TAXONOMY	
11 REMINDER TERM	
12 REMINDER ORDER CHECK ITEMS GROUP	
13 REMINDER ORDER CHECK RULES	
Select a file: (1-13): 3	
Select REMINDER DEFINITION NAME: VA-WH CHANGE IN PREGNANCY STATUS	NATIONAL
Enter another one or just press enter to go back to file selection.	
Select REMINDER DEFINITION NAME:	

3. Add the reminder dialog PXRM PATCH 45 TIU/HS OBJECTS:

- a. At the Select a File prompt, type "4" for reminder dialog.
- b. At the Select REMINDER DIALOG NAME prompt, type "PXRM PATCH 45 TIU/HS OBJECTS".
- c. At the ... OK? prompt, press the Enter key.
- d. At the Select REMINDER DIALOG NAME prompt, press the Enter key.

Example – Adding a Reminder Dialog

Onland from the following reminder files.
Select from the following reminder files:
1 REMINDER COMPUTED FINDINGS
2 REMINDER COUNTING GROUP
3 REMINDER DEFINITION
4 REMINDER DIALOG
5 REMINDER EXTRACT COUNTING RULE
6 REMINDER EXTRACT DEFINITION
7 REMINDER LIST RULE
8 REMINDER LOCATION LIST
9 REMINDER SPONSOR
10 REMINDER TAXONOMY
11 REMINDER TERM
12 REMINDER ORDER CHECK ITEMS GROUP
13 REMINDER ORDER CHECK RULES
Select a file: (1-13): 4
Select REMINDER DIALOG NAME: PXRM PATCH 45 TIU/HS OBJECTS dialog element
LOCAL
OK? Yes// (Yes)
Enter another one or just press enter to go back to file selection.
Select REMINDER DIALOG NAME:

- 4. Add the following reminder taxonomies:
 - VA-WH CURRENTLY LACTATING
 - VA-WH POSSIBLE PREGNANCY
 - a. At the Select a File prompt, type "10" for reminder taxonomy.
 - b. At the Select REMINDER TAXONOMY NAME prompt, enter the name of the taxonomy from the list above.
 - c. At the ... OK? prompt, press the Enter key.
 - d. Repeat steps b and c for the remaining taxonomies in the list above.
 - e. At the Select REMINDER TAXONOMY NAME prompt, press the Enter key.

Example – Adding a Reminder Taxonomy

```
Select from the following reminder files:
   1 REMINDER COMPUTED FINDINGS
   2 REMINDER COUNTING GROUP
   3 REMINDER DEFINITION
   4 REMINDER DIALOG
   5 REMINDER EXTRACT COUNTING RULE
   6 REMINDER EXTRACT DEFINITION
   7 REMINDER LIST RULE
   8 REMINDER LOCATION LIST
  9 REMINDER SPONSOR
  10 REMINDER TAXONOMY
 11 REMINDER TERM
 12 REMINDER ORDER CHECK ITEMS GROUP
 13 REMINDER ORDER CHECK RULES
Select a file: (1-13): 10
Select REMINDER TAXONOMY NAME: VA-WH CURRENTLY LACTATING
                                                             NATTONAL
         ...OK? Yes// (Yes)
Enter another one or just press enter to go back to file selection.
Select REMINDER TAXONOMY NAME:
```

- 5. Add the following reminder order check items groups:
 - VA-WH HIRISK CONTRACEPTIVES GROUP
 - VA-WH HIRISK IMAGING AGENTS GROUP
 - VA-WH HIRISK IMAGING PROCEDURES (MRI) GROUP
 - VA-WH HIRISK IMAGING PROCEDURES (NON MRI) GROUP
 - VA-WH HIRISK MEDICATIONS (EXTREME RISK) GROUP
 - VA-WH HIRISK MEDICATIONS (LACTATION LEVEL 1) GROUP
 - VA-WH HIRISK MEDICATIONS (LACTATION LEVEL 2) GROUP
 - VA-WH HIRISK MEDICATIONS (MOD/HIGH RISK DURING PREGNANCY) GROUP
 - VA-WH HIRISK MEDICATIONS (MODERATE/HIGH RISK) GROUP
 - a. At the Select a File prompt, type "12" for reminder order check items group.
 - b. At the Select REMINDER ORDER CHECK ITEMS GROUP GROUP NAME prompt, enter the name of the order check items group from the list above.
 - c. Repeat step b for the remaining order check items groups in the list above.
 - d. At the Select REMINDER ORDER CHECK ITEMS GROUP GROUP NAME prompt, press the Enter key.

Example – Adding a Reminder Taxonomy

Selec	t from the following reminder files:
1	REMINDER COMPUTED FINDINGS
2	REMINDER COUNTING GROUP
3	REMINDER DEFINITION
4	REMINDER DIALOG
5	REMINDER EXTRACT COUNTING RULE
6	REMINDER EXTRACT DEFINITION
7	REMINDER LIST RULE
8	REMINDER LOCATION LIST
9	REMINDER SPONSOR
10	REMINDER TAXONOMY
11	REMINDER TERM
12	REMINDER ORDER CHECK ITEMS GROUP
13	REMINDER ORDER CHECK RULES
Selec	t a file: (1-13): 12
Selec	REMINDER ORDER CHECK ITEMS GROUP GROUP NAME: VA-WH HIRISK CONTRACEPTIVESP
Enter	another one or just press enter to go back to file selection.
Selec	ET REMINDER ORDER CHECK ITEMS GROUP GROUP NAME:

- 6. At the Select a file prompt, press the Enter key.
- 7. At the Press ENTER to continue or '^' to exit prompt, press the Enter key.
- 8. At the Enter the Exchange File entry name prompt, type "HI RISK MEDICATIONS CONTENT BACKUP".

Example – Naming the Exchange Entry

Checking reminder definition(s) for errors. Checking reminder definition VA-WH PREGNANCY AND LACTATION DATA ENTRY ALLOWED WARNING: There is no Resolution logic. Press ENTER to continue or '^' to exit: No fatal reminder definition errors were found. Checking reminder definition VA-WH UPDATE LACTATION STATUS No fatal reminder definition errors were found. Checking reminder definition VA-WH UPDATE PREGNANCY STATUS No fatal reminder definition errors were found. Checking reminder definition VA-WH POTENTIALLY UNSAFE MEDICATIONS REPORT - COHOR WARNING: There is no Resolution logic. No fatal reminder definition errors were found. Checking reminder definition VA-WH CHANGE IN PREGNANCY STATUS WARNING: There is no Resolution logic. No fatal reminder definition errors were found. Checking reminder definition VA-WH NO PREGNANCY END DOC RPT WARNING: There is no Resolution logic. No fatal reminder definition errors were found. No fatal reminder definition problems were found, packing will continue. Checking reminder dialog(s) for errors.... No fatal dialog problems were found, packing will continue. Enter the Exchange File entry name: HI RISK MEDICATIONS CONTENT BACKUP

9. You will now see the Text Editor. A description for the exchange entry is auto generated and appears in the editor. You are not required to edit this description. Press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.

```
Example – Adding a Description
```

```
Enter a description of the Exchange File entry you are packing ......
. . . . . . . . . . . . . . . . . . .
The following Clinical Reminder items were selected for packing:
REMINDER ORDER CHECK ITEMS GROUP
 VA-WH HIRISK IMAGING AGENTS GROUP
 VA-WH HIRISK MEDICATIONS (LACTATION LEVEL 1) GROUP
 VA-WH HIRISK IMAGING PROCEDURES (NON MRI) GROUP
 VA-WH HIRISK IMAGING PROCEDURES (MRI) GROUP
 VA-WH HIRISK MEDICATIONS (MOD/HIGH RISK DURING PREGNANCY) GROUP
 VA-WH HIRISK MEDICATIONS (EXTREME RISK) GROUP
 VA-WH HIRISK MEDICATIONS (LACTATION LEVEL 2) GROUP
 VA-WH HIRISK MEDICATIONS (MODERATE/HIGH RISK) GROUP
 VA-WH HIRISK CONTRACEPTIVES GROUP
REMINDER DIALOG
 PXRM PATCH 45 TIU/HS OBJECTS
REMINDER TAXONOMY
 VA-WH CURRENTLY LACTATING
 VA-WH POSSIBLE PREGNANCY
```

10. You will now see the Text Editor a second time. You are not required to enter any keywords or phrases. Press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.

Example – Adding Keywords

- 11. The exchange entry is now packed.
- 12. You are returned to the Clinical Reminders Exchange screen. Note at the top of the screen that the entry you created was saved.
- 13. At the Select Action prompt, type "Q" to exit Reminder Exchange.

Example – Exiting Reminder Exchange

Clin	ic	al Rep	minder Exchange May 11,	202101	0:14:57	Page:	1 of	61
HI R	IS	K MED	ICATIONS CONTENT BACKUP was	s saved	in the Exchange	File.		
Ite	m	Entr	V	Sourc	e	Date	Packed	
	1	BDI	II RESULT GROUP	CPRSC	RM20@SALT LAKE C	04/13	/2004@1	5:53
	2	00110	ESSION/PTSD REMINDER TERM FES - PATCH 17	CPRSC	RM21@NORTHERN CAI	L 02/26	/201001	L2:07
	3		REMINDER DIALOGS	CPRSC	RM22@SALT LAKE CI	I 08/28	/201300	9:38
	4		FOR HRMH		RM23@SALT LAKE CI		/201101	
	5	GMTS	SKIN RISK HS OBJECTS		RM20@SALT LAKE CI		/2007@1	
	6	GMTS	SKIN RISK HS TYPES	CPRSC	RM200SALT LAKE CI	I 07/09	/2007@1	L3:20
	7	GMTSI	VHN	CPRSC	RM20@SALT LAKE CI	I 07/06	/2004@1	L5:06
	8	HI R	ISK MEDICATIONS CONTENT	CPRSC	RM25@WEST PALM BE	E 05/11	/202101	L0:10
		BACK	JP					
+		+	Next Screen - Prev Scree	en ??	More Actions			>>>
CFE	С	reate	Exchange File Entry	LHF	Load Host File			
CHF	C	reate	Host File	LMM	Load MailMan Mes	ssage		
CMM	C	reate	MailMan Message	LR	List Reminder De	efinition	S	
DFE			Exchange File Entry			ile		
IFE	I	nstal	l Exchange File Entry	RI	Reminder Definit	tion Inqu	iry	
ΙH	I	nstal	lation History	RP	Repack			
Sele	ct	Acti	on: Next Screen// Q					

5.7 Resolve Any Object Conflicts

Two new national TIU objects are being released in TIU*1.0*341. As part of the installation, there is a check to ensure that an object currently on the system doesn't duplicate the print name or abbreviation.

WARNING: If the object that duplicates the print name or abbreviation is found during the OR*3*498 installation process, the installation will abort.

To ensure there is no issue with the installation, you should identify any problematic object(s) and modify the print name or abbreviation, as appropriate. You need to use FileMan to search for the following:

- ABBREVIATION equal to POX
- ABBREVIATION equal to EDD
- PRINT NAME equal to PULSE OXIMETRY
- PRINT NAME equal to EXPECTED DUE DATE

The example on the next page shows the search for all of the above-listed TIU Objects.

After you have searched for the TIU objects, if you need to change the abbreviation and/or print name, you can either use FileMan or go to the following menu:

• TIU IRM Maintenance Menu--> Document Definitions (Manager) --> Edit Document Definitions.

Example – Search for TIU Objects

Select VA FileMan Option: Search File Entries Output from what File: TIU DOCUMENT DEFINITION// -A- SEARCH FOR TIU DOCUMENT DEFINITION FIELD: TYPE -A- CONDITION: EQUALS -A- EQUALS: OBJECT -B- SEARCH FOR TIU DOCUMENT DEFINITION FIELD: PRINT NAME -B- CONDITION: EQUALS -B- EQUALS: PULSE OXIMETRY -C- SEARCH FOR TIU DOCUMENT DEFINITION FIELD: PRINT NAME -C- CONDITION: EQUALS -C- EQUALS: EXPECTED DUE DATE -D- SEARCH FOR TIU DOCUMENT DEFINITION FIELD: ABBREVIATION -D- CONDITION: EQUALS -D- EOUALS: POX -E- SEARCH FOR TIU DOCUMENT DEFINITION FIELD: ABBREVIATION -E- CONDITION: EQUALS -E- EQUALS: EDD -F- SEARCH FOR TIU DOCUMENT DEFINITION FIELD: TYPE EQUALS "O" (OBJECT) IF: AB and PRINT NAME EQUALS (case-insensitive) "PULSE OXIMETRY" Or TYPE EQUALS "O" (OBJECT) OR: AC and PRINT NAME EQUALS (case-insensitive) "EXPECTED DUE DATE" OR: AD Or TYPE EQUALS "O" (OBJECT) and ABBREVIATION EQUALS (case-insensitive) "POX" Or TYPE EQUALS "O" (OBJECT) OR: AE and ABBREVIATION EQUALS (case-insensitive) "EDD" OR: STORE RESULTS OF SEARCH IN TEMPLATE: TIU 341 Are you adding 'TIU 341' as a new SORT TEMPLATE? No// Y (Yes) DESCRIPTION: Edit? NO// Sort by: NAME// Start with NAME: FIRST// First Print FIELD: [CAPTIONED Include COMPUTED fields: (N/Y/R/B): NO// BOTH Computed Fields and Record Number (IEN) Heading (S/C): TIU DOCUMENT DEFINITION Search Replace DEVICE: 0;0;999 VIRTUAL TELNET Right Margin: 80// TIU DOCUMENT DEFINITION Search JAN 10, 2022@11:10 PAGE 1 _____

6 Post-Installation Steps

These post-installation steps must be performed by the CAC.

Post-Installation steps are:

- Install reminder content.
- Enable clinical reminder order check rules.
- Link reminder terms to women's health procedure type entries.

This section will take approximately 45 minutes to complete.

6.1 Install Reminder Content

Install the following exchange file entries:

- PXRM*2.0*71 SMART UPDATES
- PXRM*2.0*71 HI RISK MEDS CONTENT

NOTE: These steps need to be performed by the site's Clinical Reminders Manager.

6.1.1 Install SMART UPDATES Content

- 1. In the Reminder Exchange, search for the text PXRM*2.0*71 SMART UPDATES and then, locate an entry titled PXRM*2.0*71 SMART UPDATES.
- 2. At the "Select Action" prompt, enter "IFE" for Install Exchange File Entry.
- 3. At the "Enter a list or range of numbers" prompt, enter the number that corresponds with your entry, PXRM*2.0*71 SMART UPDATES. (In the "Example-Reminder Exchange" below, it is 38.) It will vary by site. Make sure you select the Item whose "Date Packed" is set to 07/07/2021@09:59:53.

Example -	Reminder	Exchange

+Item Entry	Source	Date Packed
38 PXRM*2.0*71 SMART UPDATES	CPRSCRM1@CAMP MASTER	07/07/2021009:59
39 RT ALCOHOL NONE PAST 1YR	CPRSCRM2@NORTHERN CAL	02/03/2009008:32
40 TIU TEMPLATE URL FIX	CPRSCRM3@SALT LAKE CI	06/07/2016006:18
41 TIU*1*112 20040325	CPRSCRM4@ALBANY	03/25/2004@14:20
42 TIU*1*242 20080814	CPRSCRM5@ZZ ALBANY-PR	08/14/2008@08:20
43 TIU*1*246 20090128 DC/TL	CPRSCRM5@SALT LAKE CI	01/28/2009017:05
44 TIU*1*246 20090128 TL	CPRSCRM6@SALT LAKE CI	01/28/2009017:09
45 TX HIGH RISK FLU AND PNEUMONIA	CPRSCRM7@NORTHERN CAL	02/09/2009@09:22
46 UPDATE VA-DIABETES	CRRSCRM7@NORTHERN CAL	02/18/2009@10:04
47 UPDATE_2_0_1 BENEFICIARY	CPRSCRM8@SALT LAKE CI	02/18/2016013:34
+ + Next Screen - Prev Scre	en ?? More Actions	>>>
CFE Create Exchange File Entry	LHF Load Host File	
CHF Create Host File	LMM Load MailMan Messa	ge
CMM Create MailMan Message	LR List Reminder Defin	itions
DFE Delete Exchange File Entry	LWH Load Web Host File	
IFE Install Exchange File Entry	RI Reminder Definition	Inquiry
IH Installation History	RP Repack	

```
Select Action: Next Screen// IFE Install Exchange File Entry Enter a list or range of numbers (1-334): 38
```

4. At the "Select Action" prompt in the Exchange File Entry, enter "IA" for Install All Components.

```
Example – Exchange File Entry
      Component
                                                              Category
                                                                           Exists
Source:
             CPRSCRM1, ONE at CAMP MASTER
Date Packed: 07/07/2021@09:59:53
Package Version: 2.0P42
Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-WH SMART BREAST IMAGING FOLLOW-UP
REMINDER TAXONOMY
  VA-WH ULTRASOUND BREAST PROCEDURE CODES
  VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
  VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
  VA-WH MRI BREAST PROCEDURE CODES
  VA-WH MAMMOGRAM SCREENING CODES
          Enter ?? for more actions
                                                                               >>>
IA Install all Components
                                         IS
                                            Install Selected Component
Select Action: Next Screen// IA
```

- 5. The general rules for component installation are as follows:
 - a. If a component does not exist (it is NEW), select the Install action.
 - b. If these reminder terms exist (they are DIFFERENT), select the "Skip, do not install this entry" action:
 - VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - VA-WH ULTRASOUND BREAST CODES
 - VA-WH MRI BREASTS CODES
 - VA-WH MAMMOGRAM SCREENING CODES
 - c. For all other reminder terms and components, if a component does exist (it is DIFFERENT), select the "Overwrite" action.
- 6. At this point, the dialog installation begins. Install the FIRST reminder dialog:
 - a. At the "Select Action" prompt, enter "IA" for Install All. It will install the dialog, VA-WH SMART BREAST IMAGING FOLLOW-UP.

```
      Example - VA-WH SMART BREAST IMAGING FOLLOW-UP dialog install

      Packed reminder dialog: VA-WH SMART BREAST IMAGING FOLLOW-UP

      Item
      Seq. Dialog Findings
      Type
      Exists

      1
      VA-WH SMART BREAST IMAGING FOLLOW-UP
      dialog
      X

      2
      10
      VA-WH GP BR OPEN PROCEDURES*
      group
      X

      Finding:
      *NONE*
      State
      State
```

3	10.10 VA-WH BR NO OPEN PROCEDURES* e Finding: *NONE*	lement	Х
4		lement	Х
	Finding: VIEW EPISODE (REMINDER GENERAL FINDING)		Х
5	PXRM GF VIEW BUTTON	prompt	Х
6	10.20 VA-WH BR LAST 3 TRMTS e	lement	Х
	Finding: VIEW LAST THREE BR TRMTS (REMINDER GENERAL	FINDING)	Х
7	PXRM GF VIEW BUTTON	prompt	Х
8	10.25 VA-WH BLANK LINE e	lement	Х
	Finding: *NONE*		
9	10.30 VA-WH GP BR BI-RAD SELECTOR*	group	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DT Dialog Text IS Install	Selected	
DF	Dialog Findings DU Dialog Usage QU Quit		
DS	Dialog Summary IA Install All		
Sele	ct Action: Next Screen// IA		

7. At the "Install reminder dialog and all components with no further changes" prompt, enter "Yes".

- 8. Respond to the prompt on how you want to handle the non-existent finding, Q. BREAST TREATMENT MENU:
 - a. At the "Enter response" prompt, enter "P" for Replace with an existing entry.
 - b. At the "Select ORDER DIALOG NAME" prompt, enter the name of the first finding you recorded in Section 5.2, Step $\frac{#2}{2}$.
 - c. After the dialog installation completes, you will see the Dialog Components screen. At the "Select Action" prompt, enter "QU" for Quit.

Example – Dialog Components screen

	- Dialog Components Sere					
Packed	reminder dialog: VA-W	NH SMART BREAST IMA	GING FOLL	JOW-UP		
VA-WH S	MART BREAST IMAGING 1	OLLOW-UP (reminder	dialog)	installed	from exc	hange f
Item	Seq. Dialog Findings				Туре	Exists
1	VA-WH SMART BREAST :	MAGING FOLLOW-UP		(dialog	Х
2	10 VA-WH GP BR OPEN I	PROCEDURES*			group	Х
	Finding: *NONE*					
3	10.10 VA-WH BR NO OPP	IN PROCEDURES*		e	lement	Х
	Finding: *NONE*					
4	10.15 VA-WH BR EPISON				lement	
	5	PISODE (REMINDER GE	NERAL FIN			Х
5	PXRM GF VIEW BUTTON			-	prompt	
6	10.20 VA-WH BR LAST 3				lement	
_	2	AST THREE BR TRMTS	(REMINDER			
	PXRM GF VIEW BUTTON			-	prompt	
8	10.25 VA-WH BLANK LI	IE		e	lement	Х
	Finding: *NONE*					
9	10.30 VA-WH GP BR BI-				group	Х
+		Prev Screen ?? M				
	alog Details D	-			Selected	
	alog Findings DU		QU	Quit		
	alog Summary II					
Select	Action: Next Screen/,	Ý QU				

9. You will now see the Exchange File Components screen. At the "Select Action" prompt, enter "Q" for Quit.

Example – Exchange File Components screen

```
Component
                                                           Category
                                                                       Exists
            CPRSCAC1, ONE at CAMP MASTER
Source:
Date Packed: 07/07/2021@09:59:53
Package Version: 2.0P42
Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
 VA-WH SMART BREAST IMAGING FOLLOW-UP
REMINDER TAXONOMY
 VA-WH ULTRASOUND BREAST PROCEDURE CODES
 VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 VA-WH MRI BREAST PROCEDURE CODES
 VA-WH MAMMOGRAM SCREENING CODES
         + Next Screen - Prev Screen ?? More Actions
                                                                           >>>
                           IS Install Selected Component
IA Install all Components
Select Action: Next Screen// Q
```

- 10. You will now see the main exchange screen. At the "Select Action" prompt, enter "Q" for Quit.
- 11. Use option Inquire about Reminder Term [PXRM TERM INQUIRY] to capture a copy of the following terms:
 - VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - VA-WH ULTRASOUND BREAST CODES
 - VA-WH MRI BREASTS CODES
 - VA-WH MAMMOGRAM SCREENING CODES
- 12. Compare the capture for each term in the previous step with the capture you created in Section 5.2, step <u>#1</u>. Identify the findings in the pre-install capture that do not appear in the post-install capture. If there are any pre-install findings that are not included in the post-install findings, use option Add/Edit Reminder Term [PXRM TERM EDIT] to add them:

NOTE: During installation, the VA-WH ULTRASOUND OF THE BREAST CODES term listed in Section 5.2, step <u>#1</u> was renamed to VA-WH ULTRASOUND BREAST CODES and the VA-WH MRI OF THE BREASTS CODES term was renamed to VA-WH MRI BREASTS CODES. Make sure you use the new name for these terms in step a below.

- a. At the "Select Reminder Term" prompt, enter the name of the term listed above in the previous step.
- b. At the "Select Finding" prompt, enter the letters found to the right of the equals sign on the line below the name of the finding you are adding (highlighted in gray in the example on the next page), then enter a period and then finally enter the name of the finding.

Example – Pre-Install Capture Finding Not Found in Post-Install Capture

```
Finding Item: LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES
(FI(2)=TX(129))
Finding Type: REMINDER TAXONOMY
```

- c. At the "...OK?" prompt, press the Enter key.
- d. At the "Are you adding '...' as a new FINDINGS" prompt, enter the word YES.
- e. At the "FINDING ITEM" prompt, press the Enter key.
- f. At the "BEGINNING DATE/TIME" prompt, enter a caret (^).
- g. Repeat steps b through f for the remaining pre-install findings that were not found in the post-install findings.

Example – Adding a Finding to a Term

```
Select Reminder Term: VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
                                                                   NATIONAL
        ...OK? Yes// (Yes)
Choose from:
TX VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
                                                              Finding # 1
Select Finding: TX.LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES
    Searching for a REMINDER TAXONOMY, (pointed-to by FINDING ITEM)
    Searching for a REMINDER TAXONOMY
 LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES LOCAL
        ...OK? Yes// (Yes)
 Are you adding 'LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES' as
   a new FINDINGS (the 2ND for this REMINDER TERM)? No// YES (Yes)
Editing Finding Number: 2
FINDING ITEM: LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES
        - / /
BEGINNING DATE/TIME: ^
Choose from:
                                                               Finding # 2
TX LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES
TX VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
Select Finding:
```

- h. At the "Select Finding" prompt, press the Enter key.
- i. At the "Press ENTER to continue or '^' to exit" prompt, press the Enter key.
- j. At the "Edit?" prompt, enter the word YES.

Example – Closing the Edit Session

```
Select Finding:
Checking integrity of the term ...
No fatal term errors were found.
Press ENTER to continue or '^' to exit:
Input your edit comments.
Edit? NO// YES
```

k. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Added local findings removed during installation per CPRS v31b Follow-up Setup and Configuration Guide."

Example – Text Editor

```
==[WRAP ]==[INSERT ]===========< >=====(Press <PF1>H for help]====
Added local findings removed during installation per CPRS v31b Follow-up
Setup and Configuration Guide.
```

- 1. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.
- m. Repeat steps a through 1 for the remaining reminder terms listed in step #11.
- n. At the "Select Reminder Term" prompt, press the ENTER key.
- 13. Use option Inquire about Reminder Term [PXRM TERM INQUIRY] to capture a copy of the following terms:
 - VA-WH MAMMOGRAM SCREENING CODES
 - VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
- 14. Review the term captures from the previous step and identify the findings related to tomosynthesis as well as the term those findings belong to.

15. For each term with tomosynthesis-related findings you identified in the previous step, use option Add/Edit Reminder Term [PXRM TERM EDIT] to add those findings to the corresponding term listed in the "Add Tomosynthesis Findings to Term" column in the table below:

Term with Tomosynthesis Findings	Add Tomosynthesis Findings to Term
VA-WH MAMMOGRAM	VA-WH BREAST TOMOSYNTHESIS
SCREENING CODES	SCREENING
VA-WH MAMMOGRAM BILAT	VA-WH BREAST TOMOSYNTHESIS
DIAGNOSTIC CODES	BILAT
VA-WH MAMMOGRAM UNILAT	VA-WH BREAST TOMOSYNTHESIS
DIAGNOSTIC CODES	UNILAT

- a. At the "Select Reminder Term" prompt, enter the name of the term listed in the Add Tomosynthesis Findings to Term column in the table above.
- b. At the "Select Finding" prompt, enter the letters found to the right of the equals sign on the line below the name of the finding you are adding (highlighted in gray in the example below), then enter a period and finally, enter the name of the finding.

Example – Step #13 Capture Finding

```
Finding Item: LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES
(FI(2)=TX(393))
Finding Type: REMINDER TAXONOMY
```

- c. At the "...OK?" prompt, press the Enter key.
- d. At the "Are you adding '...' as a new FINDINGS" prompt, enter the word YES.
- e. At the "FINDING ITEM" prompt, press the Enter key.
- f. At the "BEGINNING DATE/TIME" prompt, enter T-90.
- g. At the "ENDING DATE/TIME" prompt, enter a caret (^).
- h. Repeat steps b through g for the remaining tomosynthesis-related findings you are adding.

Example – Adding a Finding to a Term

Select Reminder Term: VA-WH BREAST TOMOSYNTHESIS SCREENING NATIONAL ...OK? Yes// (Yes) Choose from: Finding # 1 TX VA-WH BREAST TOMOSYNTHESIS SCREENING Select Finding: TX.LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES Searching for a REMINDER TAXONOMY, (pointed-to by FINDING ITEM) Searching for a REMINDER TAXONOMY LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES LOCAL ...OK? Yes// (Yes) Are you adding 'LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES' as a new FINDINGS (the 2ND for this REMINDER TERM)? No// YES (Yes) Editing Finding Number: 2 FINDING ITEM: LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES 11 BEGINNING DATE/TIME: T-90 ENDING DATE/TIME: ^ Choose from: TX LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES Finding # 2 TX VA-WH BREAST TOMOSYNTHESIS SCREENING Finding # 1 Select Finding:

- i. At the "Select Finding" prompt, press the Enter key.
- j. At the "Press ENTER to continue or '^' to exit" prompt, press the Enter key.
- k. At the "Edit?" prompt, enter the word YES.

Example – Closing the Edit Session

```
Select Finding:
Checking integrity of the term ...
No fatal term errors were found.
Press ENTER to continue or '^' to exit:
Input your edit comments.
Edit? NO// YES
```

1. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Added local tomosynthesis-related findings per the CPRS v31b Follow-up Setup and Configuration Guide."

```
Example – Text Editor
```

- m. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.
- n. Repeat steps a through 1 for the remaining reminder terms listed in the Add Tomosynthesis Findings to Term column in the table above.
- o. At the "Select Reminder Term" prompt, press the ENTER key.
- 16. For each term in the table below, use option Add/Edit Reminder Term [PXRM TERM EDIT] to remove all tomosynthesis-related findings.

Term with Tomosynthesis Findings
VA-WH MAMMOGRAM SCREENING CODES
VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES

- a. At the "Select Reminder Term" prompt, enter the name of the desired term from the table above.
- b. At the "Select Finding" prompt, enter an accent (`) followed by the number of the tomosynthesis-related finding you need to remove.
- c. At the "FINDING ITEM" prompt, enter an at-sign (@).

- d. At the "SURE YOU WANT TO DELETE THE ENTIRE FINDING ITEM" prompt, enter the word YES.
- e. Repeat steps b through d for the remaining tomosynthesis-related findings you need to remove.

```
Example – Removing a Finding from a Term
```

```
Select Reminder Term: VA-WH MAMMOGRAM SCREENING CODES
                                                           NATIONAL
         ...OK? Yes// (Yes)
Choose from:
                                                               Finding # 2
TX LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES
TX VA-WH MAMMOGRAM SCREENING CODES
                                                               Finding # 1
Select Finding: `2 LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES
Editing Finding Number: 2
FINDING ITEM: LOCAL DIGITAL BREAST TOMOSYNTHESIS CODES
        // @
   SURE YOU WANT TO DELETE THE ENTIRE FINDING ITEM? YES (Yes)
Choose from:
TX VA-WH MAMMOGRAM SCREENING CODES
                                                               Finding # 1
Select Finding:
```

- f. At the "Select Finding" prompt, press the Enter key.
- g. At the "Press ENTER to continue or '^' to exit" prompt, press the Enter key.
- h. At the "Edit?" prompt, enter the word YES.

Example - Closing the Edit Session Select Finding:

```
Checking integrity of the term ...
No fatal term errors were found.
Press ENTER to continue or '^' to exit:
Input your edit comments.
Edit? NO// YES
```

i. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Removed local tomosynthesis-related findings per the CPRS v31b Follow-up Setup and Configuration Guide."

```
Example – Text Editor
```

- j. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.
- k. Repeat steps a through j for the remaining reminder terms listed in the Term with Tomosynthesis Findings column in the table above.
- 1. At the "Select Reminder Term" prompt, press the ENTER key.
- 17. Use option Add/Edit Reminder Term [PXRM TERM EDIT] to add a BEGINNING DATE/TIME value to every finding in the following list of terms:
 - VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES
 - VA-WH MAMMOGRAM BILAT DIAGNOSTIC CODES
 - VA-WH ULTRASOUND BREAST CODES
 - VA-WH MRI BREASTS CODES
 - VA-WH MAMMOGRAM SCREENING CODES
 - a. At the "Select Reminder Term" prompt, enter the name of the term listed above.
 - b. At the "Select Finding" prompt, enter an accent (`) followed by the number of the desired finding.
 - c. At the "FINDING ITEM" prompt, press the ENTER key.
 - d. At the "BEGINNING DATE/TIME" prompt, enter T-90.
 - e. At the "ENDING DATE/TIME" prompt, enter a caret (^).
 - f. Repeat steps b through e for every remaining finding.

Example – Adding a BEGINNING DATE/TIME Value to All Findings Select Reminder Term Management <TEST ACCOUNT> Option: TE Add/Edit Reminder Ter m Select Reminder Term: VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES NATIONAL ...OK? Yes// (Yes) Choose from: TX LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES Finding # 2 TX VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES Finding # 1 Select Finding: `1 VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES Editing Finding Number: 1 FINDING ITEM: VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES 11 BEGINNING DATE/TIME: T-90 ENDING DATE/TIME: ^ Choose from: Finding # TX LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES 2 TX VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES Finding # 1 Select Finding: `2 LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES Editing Finding Number: 2 FINDING ITEM: LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES 11 BEGINNING DATE/TIME: T-90 ENDING DATE/TIME: ^ Choose from: TX LOCAL MAMMOGRAM UNILAT DIAGNOSTIC CODES Finding # 2 TX VA-WH MAMMOGRAM UNILAT DIAGNOSTIC CODES Finding # 1

g. At the "Select Finding" prompt, press the ENTER key.

h. At the "Edit?" prompt, enter YES.

Example – Closing the Edit Session

Select Finding: Checking integrity of the term ... Input your edit comments. Edit? NO// YES i. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Added a beginning date of T-90 to all findings per CPRS v31b Follow-up Setup and Configuration Guide."

Example – Text Editor

- j. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.
- k. Repeat steps a through j for the remaining reminder terms.
- 1. At the "Select Reminder Term" prompt, press the ENTER key.
- 18. Use option Inquire about Reminder Term [PXRM TERM INQUIRY] to review the VA-WH NEXT BREAST PROCEDURE reminder term. From the output, compare the condition statement for the VA-WH NEXT BREAST PROCEDURE finding to the value you noted in Section 5.3, step <u>#2</u>. If the condition statements do not match, use option Add/Edit Reminder Term [PXRM TERM EDIT] to change the condition statement:
 - a. At the "Select Reminder Term" prompt, enter VA-WH NEXT BREAST PROCEDURE.
 - b. At the "Select Finding" prompt, enter an accent (`) followed by the number one (1).
 - c. At the "FINDING ITEM" prompt, press the ENTER key.
 - d. At the "BEGINNING DATE/TIME" prompt, press the ENTER key.
 - e. At the "ENDING DATE/TIME" prompt, press the ENTER key.
 - f. At the "OCCURRENCE COUNT" prompt, press the ENTER key.
 - g. At the "COMPUTED FINDING PARAMETER" prompt, press the ENTER key.
 - h. At the "CONDITION: I V("Procedure")["Mammogram" Replace" prompt, enter three periods (...).

- i. At the "With" prompt, please contact REDACTED for directions on the appropriate data to enter.
- j. At the "Replace" prompt, press the ENTER key.
- k. At the "CONDITION CASE SENSITIVE" prompt, enter YES.
- 1. At the "USE STATUS/COND IN SEARCH:" prompt, enter a caret (^).

```
Example – Changing the Condition Statement
```

```
Select Reminder Term: VA-WH NEXT BREAST PROCEDURE
                                                         NATIONAL
         ...OK? Yes//
                       (Yes)
Choose from:
CF VA-WH NEXT PROCEDURE
                                                                 Finding # 1
Select Finding: `1 VA-WH NEXT PROCEDURE
Display help for CF.VA-WH NEXT PROCEDURE? N// O
Editing Finding Number: 1
FINDING ITEM: VA-WH NEXT PROCEDURE//
BEGINNING DATE/TIME: T//
ENDING DATE/TIME:
OCCURRENCE COUNT:
COMPUTED FINDING PARAMETER: BR//
CONDITION: I V("Procedure") ["Mammogram" Replace ... With I
((V("Procedure")["Mammogram")!(V("Procedure")["MAMMOGRAM")!
(V("Procedure") ["mammogram"))
CONDITION CASE SENSITIVE: YES
USE STATUS/COND IN SEARCH: ^
Choose from:
CF VA-WH NEXT PROCEDURE
                                                                Finding #
                                                                             1
```

- m. At the "Select Finding" prompt, press the ENTER key.
- n. At the "Press ENTER to continue or '^' to exit" prompt, press the ENTER key.
- o. At the "Edit" prompt, enter YES.

Example – Closing the Edit Session

```
Select Finding:
Checking integrity of the term ...
No fatal term errors were found.
Press ENTER to continue or '^' to exit:
Input your edit comments.
Edit? NO// YES
```

p. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Reverted the condition statement for the VA-WH NEXT PROCEDURE finding per the CPRS v31b Follow-up Setup and Configuration Guide."

Examp	le –	Text	Editor

==[WRAP]==[INSERT]==========< >======[Press <PF1>H for help]==== Reverted the condition statement for the VA-WH NEXT PROCEDURE finding per the CPRS v31b Follow-up Setup and Configuration Guide.

- q. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.
- r. At the "Select Reminder Term" prompt, press the ENTER key.
- 19. Use option Reminder Dialogs [PXRM DIALOG/COMPONENT EDIT] to delete the finding item from the VA-WH PT CONTACT:PACT TEAM reminder dialog element:
 - a. At the "Select Item" prompt, enter CV.
 - b. At the "TYPE OF VIEW" prompt, enter E.
 - c. At the "Select Item" prompt, enter SL.
 - d. At the "Search for" prompt, enter VA-WH PT CONTACT: PACT TEAM.
 - e. At the "Find Next 'VA-WH PT CONTACT: PACT TEAM" prompt, enter NO.

f. At the "Select Item" prompt, enter the number that corresponds with your entry, VA-WH PT CONTACT:PACT TEAM. (In the "Example-Dialog Elements" below, it is 18560.) It will vary by site.

Example – Dialog Elements

+Item	Dialog	Name I	Dialog type	Status
18560	VA-WH	PT CONTACT: PACT TEAM	Dialog Element	
18561	VA-WH	PT CONTACT: SECURE MESSAGING	Dialog Element	
18562	VA-WH	REFER BIOPSY	Dialog Element	
18563	VA-WH	REFER SURGEON	Dialog Element	
18564	VA-WH	SMART BR MALE ALERT RESULT	Dialog Element	
18565	VA-WH	SMART BR MALE BI-RAD SELECTOR: BI	Dialog Element	
18566	VA-WH	SMART BR MALE BI-RAD SELECTOR: BI	Dialog Element	
18567	VA-WH	SMART BR MALE BIOPSY ALREADY OBTA	Dialog Element	
18568	VA-WH	SMART BR MALE CONSULT	Dialog Element	
18569	VA-WH	SMART BR MALE CURRENTLY UNDER TRE	Dialog Element	
18570	VA-WH	SMART BR MALE FU MAM 12M	Dialog Element	
18571	VA-WH	SMART BR MALE FU MAM 1M	Dialog Element	
18572	VA-WH	SMART BR MALE FU MAM 2M	Dialog Element	
18573	VA-WH	SMART BR MALE FU MAM 3M	Dialog Element	
+		Next Screen - Prev Screen ?? Mo	ore Actions	>>>
AD A	Add	PT List/Print All	QU Quit	
		alog INQ Inquiry/Print		
CV C	Change N	/iew TE Dialog Taxonomy	Edit	
Select	: Item:	Next Screen// 18560		

- g. At the "FINDING ITEM" prompt, enter @.
- h. At the "SURE YOU WANT TO DELETE" prompt, enter YES.
- i. At the "Select ADDITIONAL FINDINGS" prompt, enter ^.
- j. At the "Edit" prompt, enter YES.

```
Example – Removing the Finding Item
```

```
Dialog Name: VA-WH PT CONTACT:PACT TEAM
Current dialog element/group name: VA-WH PT CONTACT:PACT TEAM
Used by: VA-WH GP BR PT CONTACT METHOD (Dialog Group)
FINDING ITEM: PACT/PC TO NOTIFY// @
SURE YOU WANT TO DELETE? YES (Yes)
Select ADDITIONAL FINDINGS: ^
Checking reminder dialog for errors..
NO ERRORS FOUND
Input your edit comments.
Edit? NO// YES
```

k. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Removed the PACT/PC TO NOTIFY finding per the CPRS v31b Follow-up Setup and Configuration Guide."

```
Example – Text Editor
```

1. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.

m. At the "Select Item" prompt, enter QU.

Installation of the SMART UPDATES reminder content is now complete.

6.1.2 Install High Risk Medications Content

- 1. In the reminder exchange, search for the text, PXRM*2.0*71 HI RISK MEDS and then, locate an entry titled PXRM*2.0*71 HI RISK MEDS CONTENT.
- 2. At the "Select Action" prompt, enter "IFE" for Install Exchange File Entry.
- 3. Enter the number that corresponds with your entry, **PXRM*2.0*71 HI RISK MEDS CONTENT**. (In the "Example-Reminder Exchange" below, it is entry 37.) It will vary by site. Make sure you select the item whose "Date Packed" is set to 07/07/2021@16:35:08.

Example –Reminder Exchange

+Item	Entry	Source	Date Packed
37	PXRM*2.0*71 HI RISK MEDS	CPRSCRM9@CAMP MASTER	07/07/2021016:35
	CONTENT		
38	PXRM*2.0*71 SMART UPDATES	CPRSCRM1@CAMP MASTER	07/07/2021009:59
39	RT VA-ALCOHOL NONE PAST 1YR	CPRSCRM2@NORTHERN CAL	02/03/2009@08:32
40	TIU TEMPLATE URL FIX	CPRSCRM3@SALT LAKE CI	06/07/2016006:18
41	TIU*1*112 20040325	CPRSCRM4@ALBANY	03/25/2004@14:20
42	TIU*1*242 20080814	CPRSCRM5@ZZ ALBANY-PR	08/14/2008@08:20
43	TIU*1*246 20090128 DC/TL	CPRSCRM5@SALT LAKE CI	01/28/2009017:05
44	TIU*1*246 20090128 TL	CPRSCRM6@SALT LAKE CI	01/28/2009017:09
45	TX HIGH RISK FLU AND PNEUMONIA	CPRSCRM7@NORTHERN CAL	02/09/2009@09:22
+	+ Next Screen - Prev Screen	n ?? More Actions	>>>
CFE Cr	reate Exchange File Entry	LHF Load Host File	
CHF Cı	reate Host File	LMM Load MailMan Message	e
CMM Ci	reate MailMan Message	LR List Reminder Defin:	itions
	elete Exchange File Entry	LWH Load Web Host File	
IFE Ir	nstall Exchange File Entry	RI Reminder Definition	Inquiry
IH Ir	nstallation History	RP Repack	
Select	Action: Next Screen// IFE Inst	tall Exchange File Entry	
Enter	a list or range of numbers (1-334	4): 37	

4. At the "Select Action" prompt, enter "IA" for Install All Components.

Example –Install All Components

Component	Category	Exists
Source: CPRSCRM9,NINE at CAMP MASTER		
Date Packed: 07/07/2021016:35:08		
Package Version: 2.0P42		
Description:		
The following Clinical Reminder items were selected for	packing:	
REMINDER ORDER CHECK ITEMS GROUP		
VA-WH HIRISK MEDICATIONS (LACTATION LEVEL 1) GROUP		
VA-WH HIRISK MEDICATIONS (MOD/HIGH RISK DURING PREGNAN	ICY) GROUP	
VA-WH HIRISK MEDICATIONS (EXTREME RISK) GROUP		
VA-WH HIRISK MEDICATIONS (LACTATION LEVEL 2) GROUP		
VA-WH HIRISK MEDICATIONS (MODERATE/HIGH RISK) GROUP		
REMINDER DIALOG		
VA-WH UPDATE LACTATION STATUS		
VA-WH UPDATE PREGNANCY STATUS		
+ Enter ?? for more actions		>>>
	ected Component	
Select Action: Next Screen// IA	eccea component	
Select Action: Next Screen// IA		

- 5. The general rules for component installation are as follows:
 - a. If a component does not exist (it is NEW), select the Install action.
 - b. If a component does exist (it is DIFFERENT), select the Overwrite the current entry action.

For the following reminder terms, select the Skip action:

- i. VA-WH PREGNANCY TEST ORDERED
- ii. VA-WH POSITIVE LAB PREGNANCY TEST
- iii. VA-WH NEGATIVE LAB PREGNANCY TEST

For all reminder order check items group entries, select the Update action. If a finding does not exist, select the Delete action to continue the installation.

- **NOTE:** For test systems, you may safely ignore these missing findings. For production systems, after this section is complete, contact the Enterprise Service Desk for assistance in determining why each finding is missing and how to populate it. At this point, the dialog installations begin. Install the FIRST reminder dialog:
- c. At the "Select Action" prompt, enter "IA" for Install All. It will install the dialog VA-WH UPDATE PREGNANCY STATUS.
- d. At the "Install reminder dialog and all components with no further changes" prompt, enter "Yes".

LVallihie	-Install VA-WIT OFDATE FREGNANCE STATUS dialog	
Packed	reminder dialog: VA-WH UPDATE PREGNANCY STATUS	
Item	Seq. Dialog Findings	Туре
1	VA-WH UPDATE PREGNANCY STATUS	dialog
2	2 VA-WH TD PREGNANCY DOCUMENTATION NEEDED* Finding: *NONE*	group
3	2.2 VA-WH TD PREGNANCY CONTRA RESOURCES Finding: *NONE*	element
4	2.5 VA-WH TD PREGNANCY STATUS ACTION PROMPT* Finding: *NONE*	group
5	2.5.5 VA-WH TD PREGNANCY STATUS EDIT RECORD* Finding: WH-TD PREGNANCY EDIT RECORD (REMINDER FINDING)	group GENERAL

2.5.5.5 VA-WH TD PREGNANCY ENDED HISTORICAL

Example –Install VA-WH UPDATE PREGNANCY STATUS dialog

FINDING)

Select Action: Next Screen// IA

Finding: *NONE*

6. Respond to the prompts on how you want to handle the first non-existent finding, Q.EMERGENCY CONTRACEPTIVE QUICK ORDER:

+ Next Screen - Prev Screen ?? More Actions

a. At the "Enter response" prompt, enter "P" for Replace with an existing entry.

Add. Finding: WH-PREGNANCY DATA SOURCE (REMINDER GENERAL

Dialog DetailsDTDialog TextISInstall SelectedDialog FindingsDUDialog UsageQUQuitDialog SummaryIAInstall All

6

DD DF DS Exists X

Х

Х

Х

X X

Х

Χ

qroup

- b. At the "Select ORDER DIALOG NAME" prompt, enter the name of the first finding you recorded in Section 5.4, Step #1a.
- 7. Respond to the prompts on how you want to handle the second non-existent finding Q.PREGNANCY TEST QUICK ORDER:
 - a. At the "Enter response" prompt, enter "P" for Replace with an existing entry.
 - b. At the "Select ORDER DIALOG NAME" prompt, enter the name of the first finding you recorded in Section 5.4, Step $\frac{\#1b}{\#1b}$.
- 8. After the dialog installation completes, you will see the Dialog Components screen. At the "Select Action" prompt, enter "QU" for Quit.

Packed	reminder dialog: VA-WH UPDATE PREGNANCY STATUS		
Item	Seq. Dialog Findings	Туре	Exists
1	VA-WH UPDATE PREGNANCY STATUS	dialog	Х
2	2 VA-WH TD PREGNANCY DOCUMENTATION NEEDED* Finding: *NONE*	group	Х
3	2.2 VA-WH TD PREGNANCY CONTRA RESOURCES Finding: *NONE*	element	Х
4	2.5 VA-WH TD PREGNANCY STATUS ACTION PROMPT* Finding: *NONE*	group	Х
5	2.5.5 VA-WH TD PREGNANCY STATUS EDIT RECORD*	group	Х
	Finding: WH-TD PREGNANCY EDIT RECORD (REMINDER GEI FINDING)	JERAL	Х
	Add. Finding: WH-PREGNANCY DATA SOURCE (REMINDER (FINDING)	GENERAL	Х
6	2.5.5.5 VA-WH TD PREGNANCY ENDED HISTORICAL Finding: *NONE*	group	Х
F	+ Next Screen - Prev Screen ?? More Actions		
	alog Details DT Dialog Text IS Insta	all Selecte	d
	alog Findings DU Dialog Usage QU Quit		
	alog Summary IA Install All Action: Next Screen// QU		

- 9. Install the SECOND reminder dialog:
 - a. At the "Select Action" prompt, enter "IA" for Install All. It will install the dialog VA-WH UPDATE LACTATION STATUS.
 - b. At the "Install reminder dialog and all components with no further changes" prompt, enter "Yes".

Example –Install VA-WH UPDATE LACTATION STATUS dialog

Packed	reminder dialog: VA-WH UPDATE LACTATION STATUS			
Item	Seq. Dialog Findings	Туре	Exists	
1	VA-WH UPDATE LACTATION STATUS	dialog	Х	
2	2 VA-WH TD LACTATION DOCUMENTATION NEEDED* Finding: *NONE*	group	Х	
3	2.5 VA-WH TD LACTATION STATUS ACTION PROMPT* Finding: *NONE*	group	Х	
4	4 2.5.5 VA-WH TD LACTATION STATUS EDIT RECORD* group X			
	Finding: WH-TD LACTATION EDIT RECORD (REMINDER GENE FINDING)	RAL	Х	
Add. Finding: WH-LACTATION DATA SOURCE (REMINDER GENERAL X FINDING)				
5	2.5.5.5 VA-WH TD LACTATION STATUS EDIT Finding: *NONE*	group	Х	
6	2.5.5.5.1 VA-WH TD LACTATION STATUS YES	group	Х	
	Finding: PATIENT IS LACTATING (REMINDER GENERAL FINDING) X		Х	
+	+ Next Screen - Prev Screen ?? More Actions			
DD D	ialog Details DT Dialog Text IS Instal	l Selected	1	
DF D	ialog Findings DU Dialog Usage QU Quit			
DS D	ialog Summary IA Install All			
Select	Action: Next Screen// IA			

10. After the dialog installation completes, you will see the Dialog Components screen. At the "Select Action" prompt, enter "QU" for Quit.

Example – Dialog Components screen

Packed	reminder dialog: VA-WH UPDATE LACTATION STATUS			
Item	Seq. Dialog Findings	Туре	Exists	
1	VA-WH UPDATE LACTATION STATUS	dialog	Х	
2	2 VA-WH TD LACTATION DOCUMENTATION NEEDED* Finding: *NONE*	group	Х	
3	2.5 VA-WH TD LACTATION STATUS ACTION PROMPT* Finding: *NONE*	group	Х	
4				
	Finding: WH-TD LACTATION EDIT RECORD (REMINDER GENE FINDING)			
	Add. Finding: WH-LACTATION DATA SOURCE (REMINDER GE FINDING)	NERAL	Х	
5	2.5.5.5 VA-WH TD LACTATION STATUS EDIT Finding: *NONE*	group	Х	
6	2.5.5.5.1 VA-WH TD LACTATION STATUS YES	group	Х	
	Finding: PATIENT IS LACTATING (REMINDER GENERAL FIN			
+	+ Next Screen - Prev Screen ?? More Actions			
	ialog Details DT Dialog Text IS Instal	l Selecte	d	
DF D	ialog Findings DU Dialog Usage QU Quit			
	ialog Summary IA Install All			
Select	Action: Next Screen// QU			

11. You will now see the Exchange File Components screen. At the "Select Action" prompt, enter "Q" for Quit.

Example – Dialog Components screen

```
Component
                                                          Category
                                                                      Exists
          CPRSCRM9,NINE at CAMP MASTER
Source:
Date Packed: 07/07/2021@16:35:08
Package Version: 2.0P42
Description:
The following Clinical Reminder items were selected for packing:
REMINDER ORDER CHECK ITEMS GROUP
 VA-WH HIRISK MEDICATIONS (MOD/HIGH RISK DURING PREGNANCY) GROUP
 VA-WH HIRISK MEDICATIONS (EXTREME RISK) GROUP
 VA-WH HIRISK MEDICATIONS (MODERATE/HIGH RISK) GROUP
REMINDER DIALOG
 VA-WH TD LACTATION DO NOTHING SRN
 VA-WH TD PREGNANCY PLACE ORDERS GROUP
 VA-WH TD PREGNANCY DO NOTHING SRN
 VA-WH TD PREGNANCY CONTRA CURRENT PARTNER
     + Next Screen - Prev Screen ?? More Actions
                                                                          >>>
IA Install all Components IS Install Selected Component
Select Action: Next Screen// Q
```

Installation of the Women's Health High Risk Medications reminder content is now complete.

6.2 Enable Clinical Reminder Order Check Rules

NOTE: This step must be performed by the site's CAC.

After CPRS v31b was installed, Pharmacy Benefits Management (PBM) and Women's Health (WH) recommended that sites disable the reminder order check rules for medications classified as Lactation Level 2.

But that recommendation has been rescinded because the CPRS v31b Follow-up Build is now installed. Thus, sites must now activate the following order check rules:

- VA-WH HIRISK MEDS (LACT 2) EXPIRED DOC RULE
- VA-WH HIRISK MEDS (LACT 2) RULE

To activate the order check rules, do the following:

1. On the Reminder Managers Menu screen, enter "ROC" for Reminder Order Check Menu.

```
Example – Reminder Managers Menu
Select Core Applications <TEST ACCOUNT> Option: REMINDER MANagers Menu
   CF
          Reminder Computed Finding Management ...
   RM
          Reminder Definition Management ...
        Reminder Sponsor Management ...
   SM
   ΤΧΜ
          Reminder Taxonomy Management
   TRM
          Reminder Term Management ...
         Reminder Location List Management ...
   T.M
         Reminder Exchange
   RX
   RТ
         Reminder Test
   OS
        Other Supporting Menus ...
   INFO Reminder Information Only Menu ...
        Reminder Dialog Management ...
   DM
   CP
        CPRS Reminder Configuration ...
   RP
         Reminder Reports ...
   MST Reminders MST Synchronization Management ...
          Reminder Patient List Menu ...
   ΡL
   PAR
          Reminder Parameters ...
   VS
          NLM Value Set Menu
          Reminder Order Check Menu ...
   ROC
   СQМ
          NLM Clinical Quality Measures Menu
         Reminder Extract Menu ...
   ХМ
   GEC
         GEC Referral Report
Select Reminder Managers Menu <TEST ACCOUNT> Option: ROC Reminder Order Check
Menu
```

2. On the Reminder Order Check Menu screen, enter "RE" for Add/Edit Reminder Order Check Rule.

Example – Reminder Order Check Menu screen

```
GE Add/Edit Reminder Order Check Items Group
GI Reminder Order Check Items Inquiry
RE Add/Edit Reminder Order Check Rule
RI Reminder Order Check Rule Inquiry
TEST Reminder Order Check Test

        Select Reminder Order Check Menu <TEST ACCOUNT> Option: RE Add/Edit Reminder
Order Check Rule
```

3. When prompted with "Select Reminder Order Check Rule By", enter "N" for Order Check Rule Name.

4. When prompted again with "Select Reminder Order Check Rule By", enter "VA-WH HIRISK MEDS (LACT 2) EXPIRED DOC RULE".

Example – "Select Reminder Order Check Rule By" prompts

```
Select Reminder Order Check Rule by one of the following:
    N: ORDER CHECK RULE NAME
    R: REMINDER DEFINITION
    T: REMINDER TERM
    Q: QUIT
Select Reminder Order Check Rule by: (N/R/T/Q): N// <Enter> ORDER CHECK RULE NAME
Select Reminder Order Check Rule: VA-WH HIRISK MEDS (LACT 2) EXPIRED DOC RULE
```

5. You will now see a screen with "Rule Name: VA-WH HIRISK MEDS (LACT 2) EXPIRED DOC RULE". When prompted with "Status", enter "P" for Prod and hit the Tab key.

Example – "Status" prompt

```
RULE NAME: VA-WH HIRISK MEDS (LACT 2) EXPIRED DOC RULE
DISPLAY NAME: Known or Potential Unsafe Medication (FDB Lactation Level 2)
STATUS: P
CLASS: NATIONAL
SPONSOR: Women Veterans Health Program
REVIEW DATE:
I=INACTIVE, P=PRODUCTION, T=TESTING
Press <PF1>H for help Insert
```

6. When prompted with "Command", enter "S" for Save.

Example – Entering "S" for Save

7. When prompted again with "Command", enter "E" for Exit.

Example – Entering "E" for Exit

8. You will now see the Edit History screen. When prompted with "Edit Comments", press the <Enter> key.

Example – Edit History screen

```
Edit History
Edit by: ADPAC,ONE on 01/13/2020@09:40:03
EDIT COMMENTS:
Press enter to add a description of the changes made.
Press <PF1>H for help Insert
```

9. You will now see the Text Editor. Enter a comment. In the example below, the comment is "Changed the rule's status to production."

Example – Text Editor

- 10. After you enter a comment, press the keyboard key that is mapped to PF1 (either the NUM LOCK key on the keypad or the F1 key) and then, immediately press the E key to exit the text editor.
- 11. On the Edit History screen, when prompted with "Command", enter "E" for Exit.

```
Example - Edit History screen after a comment has been added

Edit History

Edit by: ADPAC,ONE on 01/13/2020@09:40:03

EDIT COMMENTS: Changed the rule's status to active.

EXIT COMMENTS: Changed the rule's status to active.

Exit Save Refresh Quit

Enter a COMMAND, or "^" followed by the CAPTION of a FIELD to jump to.

COMMAND: E Press <PFI>H for help Insert
```

12. To activate the VA-WH HIRISK MEDS (LACT 2) RULE, repeat steps #3 through #11.

6.3 Link Reminder Terms to Women's Health Procedure Type Entries

NOTE: This step must be performed either by the site's CAC or the site's ITOPS VistA installer, depending on who has FileMan read/write access to the WV PROCEDURE TYPE file (#790.2).

To link Reminder Terms to Women's Health Procedure Type Entries:

- 1. In FileMan, select the "ENTER OR EDIT FILE ENTRIES [DIEDIT]" option.
- 2. At the Input to what File: prompt, select "790.2".
- 3. At the EDIT WHICH FIELD: prompt, select "Reminder Term".
- 4. At the THEN EDIT FIELD: prompt, press the <enter> key.

5. Enter the Reminder Term value for the corresponding Women's Health Procedure entries.

WV Procedure Type	Reminder Term
BREAST TOMOSYNTHESIS BILAT	VA-WH BREAST TOMOSYNTHESIS BILAT
BREAST TOMOSYNTHESIS SCREENING	VA-WH BREAST TOMOSYNTHESIS SCREENING
BREAST TOMOSYNTHESIS UNILAT	VA-WH BREAST TOMOSYNTHESIS UNILAT

Example

VA F <u>ILEMAN</u> 22.2
S <u>ELECT</u> OPTION: ENTER OR EDIT FILE ENTRIES
I <u>NPUT TO WHAT</u> F <u>ILE</u> : IMAGING TYPE// 790.2 WV PROCEDURE TYPE (31 <u>ENTRIES</u>)
EDIT WHICH FIELD: ALL// <u>REMINDER</u> <u>TERM</u> THEN EDIT FIELD:
SELECT WV PROCEDURE TYPE: BREAST TOMOSYNTHESIS BILAT
REMINDER TERM: VA-WH BREAST TOMOSYNTHESIS BILAT NATIONAL OK? YES// (YES)

7 Verify Successful Installation

At this point, sites may use a locally developed smoke test script to verify that the installation and configuration were successful.