Mental Health YS*5.01*204 Deployment, Installation, Back-Out, and Rollback Guide



March 2023 Version 1.0

Department of Veterans Affairs
Office of Information and Technology (OIT)

Revision History

Date	Version	Description	Author
3/6/2023	1.0	Initial Version	Redacted

Table of Contents

1	Int	troduction1		
	1.1	Purpose	1	
	1.2	Dependencies	1	
	1.3	Constraints	1	
2	Ro	les and Responsibilities	1	
3		ployment		
	3.1	Timeline	2	
	3.2	Site Readiness Assessment	2	
	3.2	.1 Deployment Topology (Targeted Architecture)	2	
	3.2	.2 Site Information (Locations, Deployment Recipients)	2	
	3.2	.3 Site Preparation	2	
	3.3	Resources	3	
	3.3	.1 Facility Specifics (optional)	3	
	3.3			
	3.3			
	3.3			
4	Ins	stallation	3	
	4.1	Pre-installation and System Requirements		
	4.2	Platform Installation and Preparation	3	
	4.3	Access Requirements and Skills Needed for the Installation	4	
	4.4	Installation Procedure	4	
	4.5	Post-installation		
	4.6	Installation Verification Procedure		
	4.7	Database Tuning		
5	Ba	ck-Out Procedure	5	
	5.1	Back-Out Strategy		
	5.2	Back-Out Considerations	5	
	5.3	Back-Out Criteria	5	
	5.4	Back-Out Risks		
	5.5	Authority for Back-Out		
	5.6	Back-Out Procedure		
	5.7	Back-out Verification Procedure		
6	Ro	Ilback Procedure	6	
	6.1	Rollback Considerations	6	

	6.2	Rollback Criteria6	
	6.3	Rollback Risks7	
	6.4	Authority for Rollback7	
	6.5	Rollback Procedure7	
	6.6	Rollback Verification Procedure7	
7	Ar	pendix A – Acronyms 8	
		List of Tables	
Ta	ble 1:	Deployment, Installation, Back-out, and Rollback Roles and Responsibilities	1

1 Introduction

This document describes how to deploy and install the patch YS*5.01*204 of the Mental Health package, as well as how to back-out the product and rollback to a previous version or data set.

This document is a companion to the project charter and management plan for this effort.

1.1 Purpose

The purpose of this plan is to provide a single, common document that describes how, when, where, and to whom Mental Health patch YS*5.01*204 will be deployed and installed, as well as how it is to be backed out and rolled back, if necessary. The plan also identifies resources, communications plan, and rollout schedule. Specific instructions for installation, back-out, and rollback are included in this document.

1.2 Dependencies

It is assumed that this patch is being installed into a fully patched Veterans Health Information System and Technology Architecture (VistA) system. Patch YS*5.01*217 must be installed prior to this patch.

1.3 Constraints

For all other VistA sites, there are no constraints beyond the installation into an up-to-date VistA system.

2 Roles and Responsibilities

The following describes the roles and responsibilities associated with the testing and release of YS*5.01*204. This application requires both a VistA installation and an update to the web application. The Azure application manager will install the web application part of the patch. The VistA patch will be deployed via the normal PackMan route.

Table 1: Deployment, Installation, Back-out, and Rollback Roles and Responsibilities

Team	Phase / Role	Tasks	Project Phase (See Schedule)
Project Manager	Deployment	Determine and document the roles and responsibilities of those involved in the deployment	Design
Software Quality Assurance (SQA), Test Sites	Deployment	Test for operational readiness	Test
Project Manager, Release Manager	Deployment	Execute deployment	Release
Individual VistA Sites	Installation	Plan and schedule installation	Release

Team	Phase / Role	Tasks	Project Phase (See Schedule)
Azure Manager	Installation	Plan and schedule installation	Release
Release Manager	Back-out	Confirm availability of back-out instructions and back-out strategy (what are the criteria that trigger a back-out)	Release
Sustainment Team	Post Deployment	Hardware, Software and System Support	Sustain

3 Deployment

The deployment is planned as a simultaneous (National Release) rollout. Once approval has been given to nationally release, YS*5.01*204 will be available for installation and deployment at all sites.

Scheduling of test installs, testing, and production deployment will be at the site's discretion. It is anticipated there will be a 30-day compliance period.

3.1 Timeline

The deployment and installation are scheduled to run during February 2023.

3.2 Site Readiness Assessment

This section discusses the locations that will receive the YS*5.01*204 deployment.

3.2.1 Deployment Topology (Targeted Architecture)

The web part of the application for YS*5.01*204 will be deployed to the Azure application server and will be available at each site once the VistA patch is installed. Local sites, as well as regional data centers, will need to execute the VistA installation steps during the required installation period to stay synchronized with the updates to the web application.

3.2.2 Site Information (Locations, Deployment Recipients)

The initial deployment will be to Initial Operating Capability (IOC) sites for verification of functionality. Once testing is completed and approval is given for national release, YS*5.01*204 will be deployed to all VistA systems.

The Production IOC testing sites are:

- Clement J. Zablocki VAMC (Milwaukee, WI)
- Orlando VAMC (Orlando, FL)

3.2.3 Site Preparation

YS*5.01*204 requires a fully patched VistA system. Patch YS*5.01*217 must be installed prior to the installation of YS*5.01*204.

3.3 Resources

3.3.1 Facility Specifics (optional)

No specific facility instructions needed.

3.3.2 Hardware

No hardware instructions needed.

3.3.3 Software

No software instructions needed.

3.3.4 Communications

When YS*5.01*204 is released, the released-patch notification will be sent from the National Patch Module to all personnel who have subscribed to notifications for the Mental Health package patch.

4 Installation

4.1 Pre-installation and System Requirements

There are no pre-installation requirements.

4.2 Platform Installation and Preparation

This patch can be loaded with users in the system. Installation time will be less than 5 minutes.

To ensure the integrity of the transport global, use the "Verify Checksums in Transport Global" to compare the checksums with the list that follows:

```
The checksums below are new checksums, and can be checked with CHECK1^XTSUMBLD.
```

```
Select BUILD NAME: YS*5.01*204
                                     MENTAL HEALTH
YS204PST value = 15281137
          value = 35139367
YTQREST
YTQRQAD
          value = 77970444
         value = 123502753
YTQRQAD1
YTQRQAD2
         value = 52539210
YTQRQAD3
         value = 86650740
         value = 209919666
YTQRQAD4
YTQRQAD5
         value = 55478874
YTQRQAD6
         value = 58612471
         value = 254756664
YTQRQAD7
YTQRQAD8 value = 9735324
done
```

4.3 Access Requirements and Skills Needed for the Installation

Installation of YS*5.01*204 requires access to Kernel Installation and Distribution System (KIDS) options to be able to load and install the KIDS build.

4.4 Installation Procedure

This patch may be installed with users on the system although it is recommended that it be installed during non-peak hours to minimize potential disruption to users. This patch should take less than 5 minutes to install.

- 1. Choose the PackMan message containing this build. Then select the INSTALL/CHECK MESSAGE PackMan option to load the build.
- 2. From the Kernel Installation and Distribution System Menu, select the Installation Menu. From this menu:
 - A. Select the Verify Checksums in Transport Global option to confirm the integrity of the routines that are in the transport global. When prompted for the INSTALL NAME enter the patch or build name (YS*5.01*204).
 - B. Select the Backup a Transport Global option to create a backup message. You must use this option for each patch contained in the Host File. For each patch you can specify what to backup, the entire Build or just Routines. The backup message can be used to restore just the routines or everything that will restore your system to prepatch condition.
 - C. You may also elect to use the following options:
 - i. Print Transport Global This option will allow you to view the components of the KIDS build.
 - ii. Compare Transport Global to Current System This option will allow you to view all changes that will be made when this patch is installed. It compares all of the components of this patch, such as routines, DDs, templates, etc.
 - D. Select the Install Package(s) option and choose the patch to install (YS*5.01*204).
 - i. If prompted 'Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//', answer NO.
 - ii. When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//', answer NO.
 - iii. When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//', answer NO.

4.5 Post-installation

A post-install routine will run to change the name of the SCL9R instrument to SCL90R and update the launch URL for MHA Web.

4.6 Installation Verification Procedure

Open CPRS and launch the MHA Web application. Create a new assignment on a test patient and verify that the SCL9R instrument is now SCL90R. Press the CANCEL button to cancel the assignment.

4.7 Database Tuning

No database tuning required.

5 Back-Out Procedure

5.1 Back-Out Strategy

This patch updates the Mental Health Assistant – Web (MHA Web) application. If MHA Web does not perform as desired, it is possible to back out to the previous implementation.

5.2 Back-Out Considerations

If the YS*5.01*204 patch is backed out, there will be minimal impact to users.

5.3 Back-Out Criteria

A back-out should only be considered if there is a patient safety issue, if MHA Web no longer functions, or if there is some other catastrophic failure.

5.4 Back-Out Risks

The risks vary depending on what is causing the failure of the system. The main risk is that the MHA Web will be unavailable.

5.5 Authority for Back-Out

The VistA system manager determines if a back-out of YS*5.01*204 should be considered.

5.6 Back-Out Procedure

The following routines need to be restored to their previous versions:

- YTOREST
- YTQRQAD
- YTQRQAD1
- YTQRQAD2
- YTQRQAD3
- YTQRQAD4
- YTQRQAD5
- YTQRQAD6
- YTQRQAD7
- YTQRQAD8

Use the KIDS utility restore the routines backed up in section 4.4, 2B.

The SCL90R instrument name should not be changed back to SCL9R. If that is required, a new patch will be needed.

The link in CPRS will need to be reverted back to its original format. Log into VistA and go to the GUI TOOL MENU. Select SYS for system level. Identify the Sequence number for the MHA Web option. Replace /b/ in the URL with /a/. In the example screenshot below, the Sequence is 14 and the station is 965. At the Replace prompt, enter /b/. At the With prompt, enter /a/ and press [Enter]. *NOTE: In the example above 965 is a test station. Your actual station number should be substituted.

```
Select GUI Parameters <**SPPSQA**> Option: GUI tool Menu Items
CPRS GUI Tools Menu may be set for the following:
                       USR
                              [choose from NEW PERSON]
        User
       Location
                              [choose from HOSPITAL LOCATION]
    2
                      LOC
    2.5 Service
                      SRV
                             [choose from SERVICE/SECTION]
                      DIV
                             [SPPSQA]
    3
       Division
                            [SPPSQA.DEVSLC.FO-SLC.MED.VA.GOV]
        Svstem
                       SYS
        Package
                       PKG
                             [ORDER ENTRY/RESULTS REPORTING]
Enter selection: 4 System
                           SPPSQA.DEVSLC.FO-SLC.MED.VA.GOV
-- Setting CPRS GUI Tools Menu for System: SPPSQA.DEVSLC.FO-SLC.MED.VA.GOV --
Select Sequence: ?
Sequence Value
1
         VistA Imaging Display="C:\Program Files\VistA\Imaging\MagImageDisplay
2
         VistA Imaging Capture="C:\Program Files\VistA\Imaging\MagImageCapture
3
         Mental Health Assistant="C:\Program files (x86)\Vista\YS\MHA3\YS_MHA.
         MHA Web=cmd /c start msedge.exe -new-window "https://mha.med.va.gov/a
Select Sequence: 14
Sequence: 14//
Name=Command: MHA Web=cmd /c start msedge.exe -new-window "https://mha.med.va.gov
/app/home/b/?poi=%DFN&station=965" Replace /b/ With /a/
  MHA Web=cmd /c start msedge.exe -new-window "https://mha.med.va.gov/app/home/a/?p
oi=%DFN&station=965"
Select Sequence:
```

5.7 Back-out Verification Procedure

Open CPRS and launch MHA Web. After logging in with your PIV card, look at the URL. Verify that it now has /a/ in it.

6 Rollback Procedure

6.1 Rollback Considerations

There are no rollback considerations.

6.2 Rollback Criteria

There are no rollback criteria.

6.3 Rollback Risks

There are no rollback risks.

6.4 Authority for Rollback

If the backout of the patch is authorized, then that same authorization is required for rollback.

6.5 Rollback Procedure

There are no rollback procedures.

6.6 Rollback Verification Procedure

There are no rollback verification procedures.

7 Appendix A – Acronyms

Table 2: Acronyms

Acronym	Definition	
CAG	Citrix Access Gateway	
CPRS	Computerized Patient Record System	
DIBRG	Deployment, Installation, Back-out, and Rollback Guide	
IOC	Initial Operating Capability	
KIDS	Kernel Installation and Distribution System	
МНА	Mental Health Assistant	
OIT	Office of Information and Technology	
PIN	Personal Identification Number	
PIV	Personal Identity Verification	
SPP	Suicide Prevention Package	
SQA	Software Quality Assurance	
SSOi	Single Sign-On Integration	
VA	Department of Veterans Affairs	
VAMC	Veterans Affairs Medical Center	
VIP	Veteran-focused Integration Process	
VistA	Veterans Health Information System and Technology Architecture	