Mental Health YS*5.01*218

Deployment, Installation, Back-Out, and Rollback Guide



March 2023 Version 1.0

Department of Veterans Affairs Office of Information and Technology (OIT)

Revision History

Date	Version	Description	Author	
2/1/2023	1.0	Initial Version	Booz Allen Hamilton	

Table of Contents

1	Introduction1				
	1.1	1.1 Purpose			
	1.2	De	pendencies1		
	1.3	Со	nstraints1		
2	Ro	oles	and Responsibilities1		
3	De	eplo	yment2)	
	3.1	Tir	neline	2	
	3.2	Sit	e Readiness Assessment2	2	
	3.2	2.1	Deployment Topology (Targeted Architecture)2	2	
	3.2	2.2	Site Information (Locations, Deployment Recipients)	2	
	3.2	2.3	Site Preparation	2	
	3.3	Re	sources	3	
	3.3	8.1	Facility Specifics (optional)	3	
	3.3	8.2	Hardware	3	
	3.3	3.3	Software	3	
	3.3	8.4	Communications	3	
4	Ins	stal	lation3	}	
	4.1	Pre	e-installation and System Requirements	3	
	4.2	Pla	atform Installation and Preparation	3	
	4.3	Ac	cess Requirements and Skills Needed for the Installation	3	
	4.4	Ins	stallation Procedure4	ŀ	
	4.5	Ро	st-installation4	ŀ	
	4.6	Ins	stallation Verification Procedure4	ŀ	
	4.7	Da	tabase Tuning5	5	
5	Ba	ack-	Out Procedure)	
	5.1	Ва	ck-Out Strategy5	5	
	5.2	Ва	ck-Out Considerations5	5	
	5.3	Ва	ck-Out Criteria5	5	
	5.4	Ва	ck-Out Risks	5	
	5.5	Au	thority for Back-Out5	5	
	5.6	Ва	ck-Out Procedure5	5	
	5.7	Ва	ck-out Verification Procedure5	5	
6	Ro	ollba	ack Procedure6)	
	6.1	Ro	Ilback Considerations6	3	
	6.2	Ro	IIback Criteria6	3	

	6.3	Rollback Risks	6		
	6.4	Authority for Rollback	6		
		Rollback Procedure			
	6.6	Rollback Verification Procedure	8		
7	Appendix A – Acronyms9				

List of Tables

Table 1: Deployment, Installation,	Back-out, and Rollback Roles and Responsibilities1	
Table 2: Acronyms	9)

1 Introduction

This document describes how to deploy and install the patch YS*5.01*218 of the Mental Health package, as well as how to back-out the product and rollback to a previous version or data set.

This document is a companion to the project charter and management plan for this effort.

1.1 Purpose

The purpose of this plan is to provide a single, common document that describes how, when, where, and to whom Mental Health patch YS*5.01*218 will be deployed and installed, as well as how it is to be backed out and rolled back, if necessary. The plan also identifies resources, communications plan, and rollout schedule. Specific instructions for installation, back-out, and rollback are included in this document.

1.2 Dependencies

It is assumed that this patch is being installed into a fully patched Veterans Health Information System and Technology Architecture (VistA) system. Patch YS*5.01*202 must be installed prior to this patch.

1.3 Constraints

For all other VistA sites, there are no constraints beyond the installation into an up-to-date VistA system.

2 Roles and Responsibilities

The following describes the roles and responsibilities associated with the testing and release of YS*5.01*218. This application requires both a VistA installation and an update to the web application. The Azure application manager will install the web application part of the patch. The VistA patch will be deployed via the normal PackMan route.

Team	Phase / Role Tasks		Project Phase (See Schedule)	
Project Manager	Deployment	Determine and document the roles and responsibilities of those involved in the deployment	Design	
Software Quality Assurance (SQA), Test Sites	Deployment	Test for operational readiness	Test	
Project Manager, Release Manager	Deployment	Execute deployment	Release	
Individual VistA Sites	Installation	Plan and schedule installation	Release	
Azure Manager	Installation	Plan and schedule installation	Release	

Table 1: Deployment, Installation, Back-out, and Rollback Roles and Responsibilities

Team Phase / Role		Tasks	Project Phase (See Schedule)	
Release Manager	Back-out	Confirm availability of back-out instructions and back-out strategy (what are the criteria that trigger a back-out)	Release	
Sustainment Team	Post Deployment	Hardware, Software and System Support	Sustain	

3 Deployment

The deployment is planned as a simultaneous (National Release) rollout. Once approval has been given to nationally release, YS*5.01*218 will be available for installation and deployment at all sites.

Scheduling of test installs, testing, and production deployment will be at the site's discretion. It is anticipated there will be a 30-day compliance period.

3.1 Timeline

The deployment and installation are scheduled to run during March 2023.

3.2 Site Readiness Assessment

This section discusses the locations that will receive the YS*5.01*218 deployment.

3.2.1 Deployment Topology (Targeted Architecture)

The web part of the application for YS*5.01*218 will be deployed to the Azure application server and will be available at each site once the VistA patch is installed. Local sites, as well as regional data centers, will need to execute the VistA installation steps during the required installation period to stay synchronized with the updates to the web application.

3.2.2 Site Information (Locations, Deployment Recipients)

The initial deployment will be to Initial Operating Capability (IOC) sites for verification of functionality. Once testing is completed and approval is given for national release, YS*5.01*218 will be deployed to all VistA systems.

The Production IOC testing sites are:

- Clement J. Zablocki VAMC (Milwaukee, WI)
- Orlando VAMC (Orlando, FL)

3.2.3 Site Preparation

YS*5.01*218 requires a fully patched VistA system. Patch YS*5.01*202 must be installed prior to the installation of YS*5.01*218.

3.3 Resources

3.3.1 Facility Specifics (optional)

No specific facility instructions needed.

3.3.2 Hardware

No hardware instructions needed.

3.3.3 Software

No software instructions needed.

3.3.4 Communications

When YS*5.01*218 is released, the released-patch notification will be sent from the National Patch Module to all personnel who have subscribed to notifications for the Mental Health package patch.

4 Installation

4.1 Pre-installation and System Requirements

There are no pre-installation requirements.

4.2 Platform Installation and Preparation

This patch can be loaded with users in the system. Installation time will be less than 5 minutes.

To ensure the integrity of the transport global, use the "Verify Checksums in Transport Global" to compare the checksums with the list that follows:

```
The checksums below are new checksums, and can be checked with CHECK1^XTSUMBLD.
```

Select BUI	ILD NAME:	YS*5.01*2	218	MENTAL	HEALTH
YS218PST	value =	13216877			
YTQRCAT	value =	47845218			
YTQRRPT	value =	58943245			
YTRPWRP	value =	17144134			
YTSMIDAS	value =	5453302			
YTSMIOSB	value =	24768449			
YTSPRM10	value =	18948086			
YTSSBAF	value =	21134576			
YTSSBAFP	value =	8480703			
YTXCHG	value =	69915923			
YTXCHGL	value =	39988739			
done					

4.3 Access Requirements and Skills Needed for the Installation

Installation of YS*5.01*218 requires access to Kernel Installation and Distribution System (KIDS) options to be able to load and install the KIDS build.

4.4 Installation Procedure

This patch may be installed with users on the system although it is recommended that it be installed during non-peak hours to minimize potential disruption to users. This patch should take less than 5 minutes to install.

- 1. Choose the PackMan message containing this build. Then select the INSTALL/CHECK MESSAGE PackMan option to load the build.
- 2. From the Kernel Installation and Distribution System Menu, select the Installation Menu. From this menu:
 - A. Select the Verify Checksums in Transport Global option to confirm the integrity of the routines that are in the transport global. When prompted for the INSTALL NAME enter the patch or build name (YS*5.01*218).
 - B. Select the Backup a Transport Global option to create a backup message. You must use this option for each patch contained in the Host File. For each patch you can specify what to backup, the entire Build or just Routines. The backup message can be used to restore just the routines or everything that will restore your system to prepatch condition.
 - C. You may also elect to use the following options:
 - i. Print Transport Global This option will allow you to view the components of the KIDS build.
 - ii. Compare Transport Global to Current System This option will allow you to view all changes that will be made when this patch is installed. It compares all of the components of this patch, such as routines, DDs, templates, etc.
 - D. Select the Install Package(s) option and choose the patch to install (YS*5.01*218).
 - i. If prompted 'Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//', answer NO.
 - ii. When prompted 'Want KIDS to INHIBIT LOGONs during the install? NO//', answer NO.
 - iii. When prompted 'Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//', answer NO.

4.5 Post-installation

A post-install routine will run to update instrument categories, update the PCL-5, and change the Case Mix, PSS-3 2ND, and ADS instruments to be Staff Entry Only, and the COPD to be eligible for Patient Entry.

4.6 Installation Verification Procedure

Open CPRS and launch the MHA Web application. Create a new assignment on a test patient and verify that the PROMIS10, MIOS+B-IPF, SBAF, SBAF-PTSD, HIT-6, and MIDAS instruments now appear in the instrument list. Press the CANCEL button to cancel the assignment.

4.7 Database Tuning

No database tuning required.

5 Back-Out Procedure

5.1 Back-Out Strategy

This patch updates the Mental Health Assistant – Web (MHA Web) application. If MHA Web does not perform as desired, it is possible to back out to the previous implementation.

5.2 Back-Out Considerations

If the YS*5.01*218 patch is backed out, there will be minimal impact to users.

5.3 Back-Out Criteria

A back-out should only be considered if there is a patient safety issue, if MHA Web no longer functions, or if there is some other catastrophic failure.

5.4 Back-Out Risks

The risks vary depending on what is causing the failure of the system. The main risk is that the MHA Web will be unavailable.

5.5 Authority for Back-Out

The VistA system manager determines if a back-out of YS*5.01*218 should be considered.

5.6 Back-Out Procedure

The following routines need to be restored to their previous versions:

- YTQRCAT
- YTQRRPT
- YTRPWRP
- YTXCHG
- YTXCHGL

Use the KIDS utility restore the routines backed up in section 4.4, 2B.

```
*** NOTE ***
```

Follow the instrument Deactivation procedure in section 6.5 Rollback Procedure to Drop the PROMIS10, MIOS+B-IPF, SBAF, SBAF-PTSD, HIT-6, and MIDAS instruments.

5.7 Back-out Verification Procedure

Open CPRS and launch MHA Web. Click on the + sign to create a new Assignment. Verify that the PROMIS10, MIOS+B-IPF, SBAF, SBAF-PTSD, HIT-6 and MIDAS instruments are no longer visible.

6 Rollback Procedure

6.1 Rollback Considerations

To rollback this patch, six entries in the MH TESTS AND SURVEYS file need to be updated. The six instruments that were added in this patch need to be set to DROPPED.

6.2 Rollback Criteria

If the patch is backed out, the MH TESTS AND SURVEYS file must be rolled back.

6.3 Rollback Risks

By only changing the OPERATIONAL flag, LAST EDIT DATE, and LAST EDITED BY fields in the MH TESTS AND SURVEYS file, the risks are minimal as any existing MH ADMINISTRATION that used these instruments will still be available.

6.4 Authority for Rollback

If the backout of the patch is authorized, then that same authorization is required for rollback.


6.5 Rollback Procedure

Deactivate:

- PROMIS10
- MIOS+B-IPF
- SBAF
- SBAF-PTSD
- HIT-6
- MIDAS


The following example is for the PROMIS10 but will be the same for each instrument.

Log in to VistA. Go to the VA FileMan menu Go to **Enter or Edit File Entries**




At the **Input to what File:** prompt enter MH TESTS AND SURVEYS At the **EDIT WHICH FIELD** prompt enter:

OPERATIONAL LAST EDIT DATE LAST EDITED BY



At the Select MH TESTS AND SURVEYS NAME: prompt enter PROMIS10 At the OPERATIONAL prompt enter DROPPED At the LAST EDIT DATE prompt enter N for NOW At the LAST EDITED BY prompt enter your name.



Repeat this process for the remaining instruments.

6.6 Rollback Verification Procedure

Open CPRS and launch the MHA Web application. Create a new assignment on a test patient and verify that the PROMIS10, MIOS+B-IPF, SBAF, SBAF-PTSD, HIT-6, and MIDAS instruments **do not** appear in the instrument list. Press the CANCEL button to cancel the assignment.

7 Appendix A – Acronyms

Acronym Definition CAG Citrix Access Gateway **CPRS Computerized Patient Record System** DIBRG Deployment, Installation, Back-out, and Rollback Guide IOC Initial Operating Capability **KIDS** Kernel Installation and Distribution System MHA Mental Health Assistant MIDAS Maintaining Implementation through Dynamic Adaptations OIT Office of Information and Technology PIN Personal Identification Number PIV Personal Identity Verification PROMIS10 Patient-Reported Outcomes Measurement Information System SPP Suicide Prevention Package SQA Software Quality Assurance SSOi Single Sign-On Integration VA **Department of Veterans Affairs** VAMC Veterans Affairs Medical Center VIP Veteran-focused Integration Process VistA Veterans Health Information System and Technology Architecture

Table 2: Acronyms