



**PROSTHETICS Delayed Order Report
(DOR)
PATCH RMPR*3*59
INSTALLATION GUIDE**

Version 3.0

June 2003

Department of Veterans Affairs
VistA Health System Design and Development (HSD&D)

Table of Contents

Patch RMPR*3*59 Installation Guide	1
Overview	1
Client Installation.....	3
Extract the file.....	6
Install the DOR File	8

Patch RMPR*3*59 Installation Guide

Overview

Introduction

With Patch RMPR*3*59, Prosthetics offers the **Delayed Order Report (DOR)** feature, which is a GUI windows-based functionality available to users. This ***DOR Installation Guide*** will help IRM and users to download the GUI portion of the patch. See the patch description for the instructions for the Vista portion of this patch. There are two user manuals that pertain to this patch:

- Automated Delayed Order Report (DOR) User Manual
- Prosthetics Main Menu User Manual

The new **Prosthetics Main Menu** window allows you to now access the NPPD Detail Display functionality as well as the Delayed Order Report.

IMPORTANT: **Purging** **Suspense files** **BEFORE** **Installation**

Each local site has the capability to purge obsolete data; typically the default setting is set at 365 days. When the patch for the ***Delayed Orders Report*** is released, it may take several hours to load due to a large suspense file. Thus purging is probably best done on the weekend or evenings.

Purging will eliminate the amount of time to run the installation. The loading time can be reduced if everyone identifies ahead of time the suspense data, which is no longer necessary and purges it prior to installation. **Recommendation:** It is recommended to keep no more than two years of suspense data.

VPRs will need to work with their IRM staff to make the default adjustment. When data is purged, there is an automatic print queue. You will need to over ride or send the data to a null device (non-printing source). This will save reams of paper. **Do NOT send this data to the printer** or it will print every record that is purged! Also, if a job is halted, make sure the records you want to purge actually were purged.

Important: **Remove the** **NPPD Detail** **Display**

You **must** uninstall the **NPPD Detail Display** application. To do this, go to the **Control Panel → Add/Remove Programs** option and select the **NPPD Detail Display** item. Click the **REMOVE** button and the program will automatically uninstall.

RPC Broker

The assumption is made that the user's PC has the required RPC Broker Client Workstation set up. If this is not the case, the user should first install this software. The installation is beyond the scope of this installation guide, but the RPC Broker files and installation instructions can be found on this VA's RPC Broker web page: <http://vista.med.va.gov/broker/index.asp>. These files must be installed before running the **DOR** GUI application.

Continued on next page

Overview, Continued

HOSTS file

A link must be established between the client and the server via the PC's Hosts file and the PC's Windows Registry. The Hosts file (no extension) is in the following:

C:\Windows	for Windows 95 or 98
C:\WinNt\System32\Drivers\Etc	for Windows NT 4

If the Hosts file does not exist, it must be created. (There may be a hosts.sam sample file already in place.) Note that if you create the file with Notepad, Notepad automatically adds a .txt extension. To remove the extension, use Explorer to rename it, leaving off the .txt extension.


Client Installation

Hosts File	If the Hosts file already exists, do not delete anything.
-------------------	---

Steps	Follow these instructions for updating the Hosts file.
--------------	--

Step	Action
1	If the user needs to access only one DOR GUI server, and the ServerPort is 9200, then add a line to the Hosts file, as shown in the example (on the next page) using the IP address provided by the server administrator. (There must be at least one space between the address and the name. There can be no spaces in the name itself.) Continue to the next page.

Notepad – Sample HOSTS file



```
File Edit Format Help
# This is a sample HOSTS file used by Microsoft TCP/IP for Windows.
#
# This file contains the mappings of IP addresses to host names. Each
# entry should be kept on an individual line. The IP address should
# be placed in the first column followed by the corresponding host name
# The IP address and the host name should be separated by at least one
# space.
#
# Additionally, comments (such as these) may be inserted on individual
# lines or following the machine name denoted by a '#' symbol.
#
# For example:
#
#       102.54.94.97       rhino.acme.com       # source server
#       38.25.63.10       x.acme.com          # x client host
127.0.0.1       localhost
10.3.29.201     brokerserver
```

Continued on next page

Client Installation, Continued

Steps (continued)

Continue to follow these instructions for updating the Hosts file (if the Hosts file already exists.)

Step	Action
2	If the user needs to access only one DOR GUI Server and the ServerPort is not 9200, then add one line to the Hosts file as shown in scenario 1 on the previous page, and run the ServerList utility described in scenario 3 below, or use the command line parameters shown in the next section.
3	<p>If the user needs to access multiple DOR GUI Servers and/or the ServerPort is NOT 9200, and the DOR GUI Server names to the Hosts file:</p> <ul style="list-style-type: none">• 152.111.222.333 BrokerServer1 (The names can be anything you like.)• 152.111.222.444 BrokerServer2

Continued on next page

Client Installation, Continued

Multiple servers

Then run the ServerList Utility, installed in the **C:\Program Files\Vista\PROS** file of the GUI install package.

- The left windowpane of the ServerList utility shows those servers that have been set up as Host Servers in the HOSTS file.
- The right windowpane shows those servers that have been set up as BrokerServers.

Steps


Continue to follow these instructions for updating the Hosts file.

Step	Action
4	<p>To enable a BrokerServer, double-click on the name in the left window to move it to the right window, and then enter the corresponding ServerPort number.</p> <p>Note: The order of the names in the right list is the order they will be displayed in the drop-down box when starting the DOR GUI feature. If only one server is listed on the right, then you won't see the drop down box – it will automatically select that server as the default).</p>

“Del” column

A left-click in the “**Del**” column marks and un-marks server for deletion. See the following ServerList screen print for an example.

ServerList



Extract the file

Anonymous directory

The file for this patch can be obtained from the ANONYMOUS.SOFTWARE directory at one of the OI Field Offices. The preferred method is to FTP the file from DOWNLOAD.VISTA.MED.VA.GOV, which will transmit the file from the first available server.

Alternatively, you may elect to retrieve the file from a specific OI Field Office. The documentation can be read via the Adobe Acrobat Reader browser program. The .EXE file is binary file and must be transferred using binary file transfer tools.

Steps

To begin the install process and unzip the files, follow these steps:

Step	Action
1	From the Anonymous Directory, download the file RMPR_3_59.EXE .
2	Double click the RMPR_3_59.exe icon.

Warning: Wrong filename

Warning: If the filename is **RMPR_3_59.exe;1** when you try to download it, then you **must** rename it to **RMPR_3_59.exe** so that the filename does not have an “;1” in it.

Rename the file

To change the filename, right click on **RMPR_3_59.exe;1**, (executable application) and select the **Rename** option. Then delete the “;1” from the filename.

Continued on next page


Extract the file, Continued

Unzipping files

The file RMPR_3_59.EXE is a zip file that contains the “setup” files needed for the automated installation. The files will unzip to the C:\PROS folder for setup purposes. As an alternative to extracting these files to every client computer, IRM staff can unzip this file to a shared folder (directory) on a network server and then run the installation “SETUP.EXE” remotely.

Recommendation: Do not change the file folder, as future patches will need to overwrite these files.

WinZip dialog box sample
(Your dialog box will contain the most current version number for the files.)



Steps (continued)

To continue to unzip the files, follow these steps:

Step	Action
3	In the WinZip Self-Extractor dialog box, click the Unzip button (as shown above).
4	A successful message box displays when the unzip process has completed. (See below.)
5	Click the OK button on the WinZip Self-Extractor confirmation box.
6	The WinZip Self-Extractor dialog box continues to display (as shown above).
7	Click the Close button.

Successful message box



Install the DOR File

Automatic Wizard


Once you have extracted the files onto your c-drive, you can perform an installation using an automatic wizard that is provided for you. This feature has combined several steps to make installation easy.

Steps

To begin the installation, follow these steps:

Step	Action
1	Open the My Computer icon on your desktop by double clicking it.
2	Double click the C-drive in your Windows Explorer view.
3	Double click the folder that contains the files: PROS .

Windows Explorer



Continued on next page


Install the DOR File, Continued

Steps (continued)

To run the installation program, follow these steps:

Step	Action
4	Double click the Setup.exe file in the C:\PROS folder shown on your Windows Explorer view.

Windows Explorer



Continued on next page


Install the DOR File, Continued

Steps (continued)


To run the wizard for the **DOR** installation, follow these steps:

Step	Action
5	The Install Shield Wizard begins and displays the download process.
6	Click the Next button on the first window of the Install Shield Wizard for the DOR .

Install Shield Wizard



NPPD Detail Display Setup



Continued on next page


Install the DOR File, Continued

Steps (continued)

To continue the installation, follow these steps:

Step	Action
7	Click the Yes button on the License Agreement if you agree to the terms of the license. Note: If you choose No , the setup will automatically close.

License Agreement



Continued on next page


Install the DOR File, Continued

Recommendation The destination folder “**C:\Program Files\Vista\PROS**” has been established by VA policy so as not to conflict with other VA programs. It is recommended that you do NOT change this portion of the path.


Steps To continue the installation, follow these steps:

Step	Action
8	Click the Next button on the Destination Location dialog box.
9	Click OK on the confirmation box. Notice that it states that the PROS file folder was created.

File Destination Location



Create Directory message



Continued on next page


Install the DOR File, Continued

Steps (continued)

To continue the installation, follow these steps:

Step	Action
10	Click the Next button on the Start Copy Files dialog box.

Start Copying Files



Continued on next page


Install the DOR File, Continued

Steps (continued)

To continue to install the Prosthetics feature onto your desktop, follow these steps:

Step	Action
11	The Install Shield Wizard will be complete with this step.
12	Click the Finish button.

Install Shield Wizard



Installation Complete

The application is now installed on C:\Program Files\Vista\PROS. An icon has been placed on your desktop to run the DOR application. The application can also be run from **Start → Programs → Prosthetics Vista Suite**.

Uninstall procedures

To uninstall the application, go to the **Control Panel → Add/Remove Programs** option and select the **Prosthetics Vista Suite** item. Click the **REMOVE** button and the program will automatically uninstall. To complete the full removal, the setup files will need to be deleted from the file location you selected.