



**SURGERY
PATCH SR*3.0*41
INSTALLATION GUIDE**

June 1998

Department of Veterans Affairs
Veterans Health Administration
Office of Chief Information Officer

Table of Contents

Introduction.....	1
Pre-Installation Tasks	1
Installation Procedures	1
Example Installation.....	3
Post-Installation Procedures.....	7
Post-Installation Procedures (HL7 V. 1.5)	7
Post-Installation Procedures (HL7 V. 1.6)	8
1. Confirm SR AAIS and SR SURGERY Application Entries	9
2. Confirm SR AAIS Server Protocol Entries	10
3. Confirm SR AAIS Client Protocol Entries	10
4. Confirm SR SURGERY Server Protocol Entries	12
5. Confirm SR SURGERY Client Protocol Entries.....	14
6. Confirm Logical Links for HLLP.....	15
7. Confirm Logical Links for TCP/IP.....	15
8. Edit SR AAIS and SR SURGERY Logical Link Device Field.....	15
9. Start LLP(s).....	16
10. Flag Interface Fields.....	17
Additional Technical Information.....	19
New Protocol File (# 101) Entries	19
New HL7 NON-DHCP APPLICATION PARAMETER (# 770) Entry.....	19
New HL7 APPLICATION PARAMETER (# 771) Entries.....	19
New MAIL GROUP (#3.8) Entry.....	20
New/Modified Routines.....	20
New Field in SURGERY SITE PARAMETERS File #133.....	20
New File	20
New Options.....	21
Modified Options.....	22
Example VISTA Device File Entries.....	23

Table of Contents

Introduction

This patch introduces new functionality to the Surgery package, enabling the Surgery package to send Health Level 7 (HL7) messages to Commercial Off-The-Shelf (COTS) systems. Sites that do not need this new functionality now should still install the patch to keep the Veterans Health Information Systems and Technology Architecture (**VISTA**) Surgery software up to date, but do not have to perform the post-installation procedures until needed.

Pre-Installation Tasks

1. Interfaces over RS-232 ports will need the Information Resource Management (IRM) staff to add new device entries. Some sites have identified a need to add two, incoming and outgoing, devices to prevent message collision. Message collision can be an issue if the system cannot be adjusted for intermittent delay (read/acknowledgment) times. The **VISTA** HL7 package can be adjusted to prevent message collision by simply changing the default Read Timeout (3) and ACK Timeout (10). These times can be adjusted by using the Interface Workbench option from within the HL7 package. However, if adjusting the delay times does not fix the problem then it might be necessary to add a second device to prevent message collision. To add these entries, IRM should follow the sample in the Technical Information section of this document.
2. Interfaces over TCP/IP ports will need to identify the system's IP address and two available port numbers. The range of recommended port numbers is 5130 to 5139.

Installation Procedures

The routines and files contained in this release are exported as Surgery patch SR*3.0*41.

<u>Step</u>	<u>Description</u>
-------------	--------------------

- | | |
|----|---|
| 1. | This patch may be installed at any time with users remaining on the system. Installation should take less than one minute. (See Example Installation.) |
| 2. | Download the KIDS file SR_3_P41.KID from the anonymous.software directory at the CIOFO at Hines, Albany, or Salt Lake City to the appropriate directory on your system. |

3. From the *Kernel Installation and Distribution System (KIDS) Menu* [XPD MAIN], select the *Installation* menu. Use the *Load a Distribution* [XPD LOAD DISTRIBUTION] option to load the distribution from the downloaded distribution file.
4. Use *Verify Checksums in Transport Global* [XPD PRINT CHECKSUM] option to confirm the patch for correctness.
5. Review your mapped set. If any of the following routines are mapped, remove them from the mapped set.

SRHLDW	SRHLDW1	SRHLENV	SRHLMFN	SRHLOORU	SRHLORU	SRHLPOST
SRHLPRE	SRHLQRY	SRHLSCRN	SRHLU	SRHLUI	SRHLUO	SRHLUO1
SRHLUO2	SRHLUO3	SRHLUO4	SRHLUO4C	SRHLVOOR	SRHLVORU	SRHLVQRY
SRHLVU	SRHLVUI	SRHLVUI2	SRHLVUO	SRHLVUO1	SRHLVUO2	SRHLVUO4
SRHLVZIU	SRHLVZQR	SRHLVZSQ	SRHLXTMP	SRHLZIU	SRHLZQR	SROCODE
SROERR	SROERR0	SROFILE	SRONITE	SRSCHD2		

6. If you are a test site for this patch, use the *Stop LLP* [HL STOP] option under the *HL7 Communications Server Option* [HL COMMUNICATIONS SERVER] menu to stop any SR* lower level protocols (LLP). If you are not a test site for this patch, skip this step. If you are not sure, use the *Systems Link Monitor* option under the *HL7 Communications Server Option* [HL COMMUNICATIONS SERVER] to see if any SR* lower level protocols (LLP) that are actively running.
7. From the *Kernel Installation and Distribution System (KIDS) Menu* [XPD MAIN] menu, select the *Installation* [XPD INSTALLTION MENU] menu. Use the *Install Package(s)* [XPD INSTALL BUILD] option and select the package SR*3.0*41 and proceed with install. (See Example Installation.)
8. If any routines were unmapped as part of step 5, they should be returned to the mapped set once the installation has run to completion.
9. The checksum routine SR41NTEG may be run from programmer mode to confirm the integrity of the routines distributed with this patch. Afterwards, SR41NTEG may be deleted from the routine directory.

Example Installation

KIDS Step #1

Select Kernel Installation & Distribution System Option: Installation

Select Installation Option: Load a Distribution

Enter a Host File: SR_3_P41.KID

KIDS Distribution saved on Feb 24, 1998@09:49:54

Comment: SR*3*41 - Surgery HL7 Interface Patch

This Distribution contains Transport Globals for the following Package(s):
SR*3.0*41

Want to Continue with Load? YES// <RET>

Loading Distribution...

Want to RUN the Environment Check Routine? YES// <RET>

SR*3.0*41

Will first run the Environment Check Routine, SRHLENV

Checking the SURGERY HL7 Environment

Environment check was successful!!

Use INSTALL NAME: SR*3.0*41 to install this Distribution.

KIDS Step #2

Select Installation Option: Verify Checksums in Transport Global

Select INSTALL NAME: SR*3.0*41 Loaded from Distribution

2/24/98@10:07:53

=> SR*3*41 - Surgery HL7 Interface Patch ;Created on Feb 24, 1998@09:49:

DEVICE: HOME// <RET> DECSERVER

PACKAGE: SR*3.0*41 Feb 24, 1998 10:08 am

PAGE 1

40 Routine checked, 0 failed.

KIDS Step #3

Select Installation Option: Install Package(s)

Select INSTALL NAME: SR*3.0*41 Loaded from Distribution 2/24/98@10:07:

53

=> SR*3*41 - Surgery HL7 Interface Patch ;Created on Feb 24, 1998@09:49:

This Distribution was loaded on Feb 24, 1998@10:07:53 with header of

Surgery Patch SR*3.0*41

SR*3*41 - Surgery HL7 Interface Patch ;Created on Feb 24, 1998@09:49:54
It consisted of the following Install(s):
SR*3.0*41

SR*3.0*41
Will first run the Environment Check Routine, SRHLENV

Checking the SURGERY HL7 Environment

Environment check was successful!!

Install Questions for SR*3.0*41

Incoming Files:

133 SURGERY SITE PARAMETERS (Partial Definition)
Note: You already have the 'SURGERY SITE PARAMETERS' File.

133.2 SURGERY INTERFACE PARAMETER (including data)

Want to DISABLE Scheduled Options, Menu Options, and Protocols? YES// **NO**

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME// **;P-DEC DECSERVER**

Install Started for SR*3.0*41 :
Feb 24, 1998@10:09:25

Installing Routines:.....
Feb 24, 1998@10:09:27

Running Pre-Install Routine: ^SRHLPRE.

Mail Group SRHL DISCREPANCY installed

Installing Data Dictionaries: ...
Feb 24, 1998@10:09:27

Installing Data: .
Feb 24, 1998@10:09:28

Installing PACKAGE COMPONENTS:

Installing HL LOWER LEVEL PROTOCOL PARAMETER.....

Installing HL7 APPLICATION PARAMETER...

Installing PROTOCOL.....
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.


```
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.
Located in the SR (SURGERY) namespace.

Installing OPTION.....
      Feb 24, 1998@10:09:30

Running Post-Install Routine: ^SRHLPOST.

Updating the HL LOWER LEVEL PROTOCOL PARAMETER file (#869.2)

      Adding SRAAISTSURGERY-TCP

Updating HL LOGICAL LINK file (#870).

      Adding SR AAIS

      Adding SR SURGERY

Updating the PROTOCOL file (#101) with the new TCP/IP connect

SURGERY HL7 SETUP IS COMPLETE.

Updating Routine file.....

Updating KIDS files.....

SR*3.0*41 Installed.
      Feb 24, 1998@10:09:31

Install Completed
```


Post-Installation Procedures

For applications that are not compatible with **VISTA** HL7 V. 1.6, this patch includes a form of the Surgery interface that is compatible with **VISTA** HL7 V.1.5. (Most applications by far should be compatible with **VISTA** HL7 V. 1.6.) Determine which form of the Surgery interface is appropriate and edit the site parameter that determines which Surgery HL7 interface will be used. Use the *Update Interface Parameter Field* [SRHL PARAMETER] option on the *Surgery Interface Management Menu* [SRHL INTERFACE] to update this parameter to YES if, and only if, the application interfacing with Surgery is NOT compatible with **VISTA** HL7 V. 1.6. Otherwise, enter NO or leave the field blank.

Select Surgery Package Management Menu Option: SI Surgery Interface Management Menu

Select Surgery Interface Management Menu Option: P Update Interface Parameter Field

This option may be used to edit the site parameter that determines which Surgery HL7 interface will be used, the interface compatible with VISTA HL7 v1.6 or the older one compatible with VISTA HL7 v1.5.

If applications communicating with the Surgery HL7 interface must use the interface designed for HL7 v1.5, enter YES. Otherwise, enter NO or leave this field blank.

Use Surgery Interface Compatible with VISTA HL7 v1.5 (Y/N): NO

If the Surgery interface that is compatible with **VISTA** HL7 V. 1.5 is being used, continue with the post-installation procedures for HL7 V.1.5. Otherwise, continue with the post-installation procedures for HL7 V. 1.6.

Post-Installation Procedures (HL7 V. 1.5)

For applications that are compatible only with **VISTA** HL7 V. 1.5, the following post-installation instructions apply.

1. Use the *Non-DHCP Application Parameter Enter/Edit* [HL EDIT SITE PARAM] option (on the *V. 1.5 OPTIONS* [HL MENU 1.5] menu on the *HL7 Main Menu* [HL7 MAIN MENU]) to define an HL7 device for SR AAIS.
2. Use the *Initiate Background Task* [HL TASK] option (on the *V. 1.5 OPTIONS* [HL MENU 1.5] menu on the *HL7 Main Menu* [HL MAIN MENU]) to create a background task to start up the lower level protocol routine for SR AAIS.
3. Using the option *Flag Interface Fields* [SRHL INTERFACE FIELDS], flag identifiers as appropriate to control the flow of data between the **VISTA** Surgery package and the ancillary system.

Post-Installation Procedures (HL7 V. 1.6)

For applications that are compatible with **VISTA** HL7 V. 1.6, the following post-installation instructions apply.

Using the *HL7 Interface Workbench* [HL INTERFACE WORKBENCH] option identify the Lower Level Protocol for your interface. Interfaces using a Hybrid Lower Level Protocol (HLLP) or an RS-232 port connection with one port for incoming and outgoing messages need to match the SR AAIS logical link with the SR HLLP LLP Parameter (option 1). Interfaces using an HLLP or an RS-232 port connection with two ports, incoming and outgoing, need to match the SR AAIS (shown in option 1) as well as SR SURGERY logical link with the SR HLLP2 LLP Parameter (option 2). The logical links that are used will also need to be linked to the devices setup in the Pre-Installation Instructions. Interfaces using the Minimal Lower Level Protocol (MLLP) or TCP/IP connections will need to match the SR AAIS logical link with the SRAAISTSURGERY-TCP LLP Parameter (option 3) as well as SR SURGERY logical link with the SRSURGERY-TCP-RECV LLP Parameter (option 4). Once the SR AAIS and SR SURGERY logical links are matched up with the appropriate LLP Parameters, use VA FileMan to add the TCP/IP addresses and ports to the HL LOWER LEVEL PROTOCOL PARAMETER file (#869.2) entries. See option 5 for the appropriate FileMan setups for SRAAISTSURGERY-TCP and SRSURGERY-TCP-RECV entries.

OPTION 1

(10) SR AAIS

LLP Parameter: SR HLLP
 LLP Type: HLLP (H)
 Device:
 Version ID: 2.2
 Block Size: <DEFAULT>
 Queue Size: <DEFAULT>

Read Timeout: 3
 ACK Timeout: <DEFAULT>
 Re-transmission Attempts: 3
 Start Block Character: <DEFAULT>
 End Block Character: <DEFAULT>

OPTION 2

(10) SR SURGERY

LLP Parameter: SR HLLP2
 LLP Type: HLLP (H)
 Device:
 Version ID: 2.2
 Block Size: <DEFAULT>
 Queue Size: <DEFAULT>

Read Timeout: 3
 ACK Timeout: <DEFAULT>
 Re-transmission Attempts: 3
 Start Block Character: <DEFAULT>
 End Block Character: <DEFAULT>

OPTION 3

(9) SR AAIS

LLP Parameter: SRAAISTSURGERY-TCP
 LLP Type: TCP
 Queue Size: <DEFAULT>

LLP is not supported by the workbench. FileMan must be used to fill in the LLP specific info.

OPTION 4

(10) SR SURGERY

LLP Parameter: SRSURGERY-TCP-RECV

LLP is not supported by the workbench. FileMan

LLP Type: TCP
Queue Size: <DEFAULT>

must be used to fill in
the LLP specific info.

OPTION 5

NUMBER: 16
LLP TYPE: TCP
TCP/IP PORT: [REMOTE PORT]

NAME: SRAAISTSURGERY-TCP
TCP/IP ADDRESS: [REMOTE IP ADDRESS]
CLIENT/SERVER: CLIENT

NUMBER: 17
LLP TYPE: TCP
TCP/IP PORT: [LOCAL PORT]

NAME: SRSURGERY-TCP-RECV
TCP/IP ADDRESS: [LOCAL IP ADDRESS]
CLIENT/SERVER: SERVER

1. Confirm SR AAIS and SR SURGERY Application Entries

Confirm that the following Application entries were created by the KIDS install and add facility name and Mail Group (if needed).

HL-7 Interface Workbench

Jun 19, 1997 15:05:48
Currently Defined
Applications

Page: 15 of 16

(32) SR AAIS

Facility Name:
Country Code: USA
Mail Group:

Active/Inactive: ACTIVE
HL7 Field Separator: <DEFAULT>
HL7 Encoding Characters: <DEFAULT>

(33) SR SURGERY

Facility Name:
Country Code: USA
Mail Group:

Active/Inactive: ACTIVE
HL7 Field Separator: <DEFAULT>
HL7 Encoding Characters: <DEFAULT>

For the above two (2) Application entries created by the KIDS install, edit and set the Facility Name and Mail Group fields. To do this, at the "Select Tool:Jump To Next App//" prompt, choose "EA Edit Application":

2. Confirm SR AAIS Server Protocol Entries

Confirm that the following Server Protocol entries were created by the KIDS install for the SR AAIS application:

```

HL-7 Interface Workbench      Jun 19, 1997 15:22:26      Page:      1 of      2
                               Server Protocols For
                               SR AAIS

(1) SR Unsolicited transmission of AAIS Requested Observation
    Text: Unsolicited Transmission of Placer Surgery Observation
    Package: SURGERY
        Message Type Rcvd: ORU                      Accept Acknowledgement: <DEFAULT>
        Event Type Rcvd: R01                      Application Acknowledgement: <DEFAULT>
            Priority: <DEFAULT>
        Processing ID: PRODUCTION
        Version ID: 2.1
        Entry Action: <NONE>
        Exit Action: <NONE>
    Process ACK Routine:
    Subscribers: SR Receiver of Unsolicited Requested Observation

(2) SR Query for Scheduling Information
    Text: Query for Surgery Case Scheduling Information
    Package: SURGERY
        Message Type Rcvd: QRY                      Accept Acknowledgement: <DEFAULT>
        Event Type Rcvd: S25                      Application Acknowledgement: <DEFAULT>
            Priority: <DEFAULT>
        Processing ID: PRODUCTION
        Version ID: 2.1
        Entry Action: <NONE>
        Exit Action: <NONE>
    Process ACK Routine:
    Subscribers: SR Receiver of Scheduling Query

```

3. Confirm SR AAIS Client Protocol Entries

Confirm that the following Client Protocol entries were created by the KIDS install for the SR AAIS application:

```

HL-7 Interface Workbench      Jun 19, 1997 15:22:46      Page:      1 of      9
                               Client Protocols For
                               SR AAIS

(1) SR Receiver of Observation Unsolicited
    Text: Receiver of ORU
    Package: SURGERY
        Message Type Rcvd: ORU                      Sending Facility Required: YES
        Message Type Sent: ACK                      Receiving Facility Required: YES
    Event Type Rcvd/Sent: R01                      Date/Time of Message Required: YES
            Priority: IMMEDIATE                      Security Required: <DEFAULT>
        Processing ID: PRODUCTION                      Logical Link: SR AAIS
        Version ID: 2.1
        Entry Action: <NONE>
        Exit Action: <NONE>
    Process Routine: D ^SRHLORU

```

(2) SR Receiver of New Appointment Booking

Text: Receiver of New Surgery Appointment Booking

Package: SURGERY

Message Type Rcvd: ACK	Sending Facility Required: YES
Message Type Sent: ZIU	Receiving Facility Required: YES
Event Type Rcvd/Sent: S12	Date/Time of Message Required: YES
Priority: IMMEDIATE	Security Required: <DEFAULT>
Processing ID: PRODUCTION	Logical Link: SR AAIS
Version ID: 2.1	
Entry Action: <NONE>	
Exit Action: <NONE>	
Process Routine: Q	

(3) SR Receiver of Appointment Rescheduling

Text: Receiver of Surgery Appointment Rescheduling

Package: SURGERY

Message Type Rcvd: ACK	Sending Facility Required: YES
Message Type Sent: ZIU	Receiving Facility Required: YES
Event Type Rcvd/Sent: S13	Date/Time of Message Required: YES
Priority: IMMEDIATE	Security Required: <DEFAULT>
Processing ID: PRODUCTION	Logical Link: SR AAIS
Version ID: 2.1	
Entry Action: <NONE>	
Exit Action: <NONE>	
Process Routine: Q	

(4) SR Receiver of Appointment Modification

Text: Receiver of Surgery Appointment Modification

Package: SURGERY

Message Type Rcvd: ACK	Sending Facility Required: YES
Message Type Sent: ZIU	Receiving Facility Required: YES
Event Type Rcvd/Sent: S14	Date/Time of Message Required: YES
Priority: IMMEDIATE	Security Required: <DEFAULT>
Processing ID: PRODUCTION	Logical Link: SR AAIS
Version ID: 2.1	
Entry Action: <NONE>	
Exit Action: <NONE>	
Process Routine: Q	

(5) SR Receiver of Appointment Cancellation

Text: Receiver of Surgery Appointment Cancellation

Package: SURGERY

Message Type Rcvd: ACK	Sending Facility Required: YES
Message Type Sent: ZIU	Receiving Facility Required: YES
Event Type Rcvd/Sent: S15	Date/Time of Message Required: YES
Priority: IMMEDIATE	Security Required: <DEFAULT>
Processing ID: PRODUCTION	Logical Link: SR AAIS
Version ID: 2.1	
Entry Action: <NONE>	
Exit Action: <NONE>	
Process Routine: Q	

(6) SR Receiver of Appointment Deletion

Text: Receiver of Surgery Appointment Deletion

Package: SURGERY

Message Type Rcvd: ACK	Sending Facility Required: YES
Message Type Sent: ZIU	Receiving Facility Required: YES
Event Type Rcvd/Sent: S17	Date/Time of Message Required: YES
Priority: IMMEDIATE	Security Required: <DEFAULT>
Processing ID: PRODUCTION	Logical Link: SR AAIS
Version ID: 2.1	
Entry Action: <NONE>	
Exit Action: <NONE>	
Process Routine: Q	

(7) SR Receiver of Staff Master File Notification

Text: Receiver of Surgery Staff Master File Notification

Package: SURGERY
Message Type Rcvd: ACK
Message Type Sent: MFN
Event Type Rcvd/Sent: M02
Priority: IMMEDIATE
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process Routine: Q
Sending Facility Required: YES
Receiving Facility Required: YES
Date/Time of Message Required: YES
Security Required: <DEFAULT>
Logical Link: SR AAIS

(8) SR Receiver of Master File Notification
Text: Receiver of Master File Notification

Package: SURGERY
Message Type Rcvd: ACK
Message Type Sent: MFN
Event Type Rcvd/Sent: M01
Priority: IMMEDIATE
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process Routine: Q
Sending Facility Required: YES
Receiving Facility Required: YES
Date/Time of Message Required: YES
Security Required: <DEFAULT>
Logical Link: SR AAIS

4. Confirm SR SURGERY Server Protocol Entries

Confirm that the following Server Protocol entries were created by the KIDS install for the SR SURGERY application:

HL-7 Interface Workbench Jun 19, 1997 15:31:34 Page: 1 of 9
Server Protocols For
SR SURGERY

(1) SR Unsolicited transmission of DHCP Requested Observation
Text: Unsolicited Transmission of Filler Surgery Observation

Package: SURGERY
Message Type Rcvd: ORU
Event Type Rcvd: R01
Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine:
Subscribers: SR Receiver of Observation Unsolicited
Accept Acknowledgement: <DEFAULT>
Application Acknowledgement: <DEFAULT>

(2) SR Notification of Appointment Booking
Text: Notification of New Surgery Appointment Booking

Package: SURGERY
Message Type Rcvd: ZIU
Event Type Rcvd: S12
Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of New Appointment Booking
Accept Acknowledgement: <DEFAULT>
Application Acknowledgement: <DEFAULT>

(3) SR Notification of Appointment Rescheduling
Text: Notification of Surgery Appointment Rescheduling
Package: SURGERY


```

Message Type Rcvd: ZIU                               Accept Acknowledgement: <DEFAULT>
Event Type Rcvd: S13                                 Application Acknowledgement: <DEFAULT>
Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of Appointment Rescheduling

```

(4) SR Notification of Appointment Modification
Text: Notification of Surgery Appointment Modification
Package: SURGERY

```

Message Type Rcvd: ZIU                                Accept Acknowledgement: <DEFAULT>
Event Type Rcvd: S14                                Application Acknowledgement: <DEFAULT>
Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of Appointment Modification

```

(5) SR Notification of Appointment Cancellation
Text: Notification of Surgery Appointment Cancellation
Package: SURGERY

```

Message Type Rcvd: ZIU
Event Type Rcvd: S15
Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of Appointment Cancellation

```

(6) SR Notification of Appointment Deletion
Text: Notification of Surgery Appointment Deletion
Package: SURGERY

```

Message Type Rcvd: ZIU                               Accept Acknowledgement: <DEFAULT>
Event Type Rcvd: S17                                Application Acknowledgement: <DEFAULT>
Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of Appointment Deletion

```

(7) SR Other Master File Notification
Text: Master File Update

```

Package: SURGERY
  Message Type Rcvd: MFN
    Event Type Rcvd: M01
      Priority: <DEFAULT>
        Processing ID: PRODUCTION
          Version ID: 2.1
            Entry Action: <NONE>
              Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of Master File Notification

```

(8) SR Staff Master File Notification
Text: Surgery Staff Master File Notification
Package: SURGERY

```

Message Type Rcvd: MFN
Event Type Rcvd: M02
Accept Acknowledgement: <DEFAULT>
Application Acknowledgement: <DEFAULT>

```

Priority: <DEFAULT>
Processing ID: PRODUCTION
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process ACK Routine: Q
Subscribers: SR Receiver of Staff Master File Notification

5. Confirm SR SURGERY Client Protocol Entries

Confirm that the following Client Protocol entries were created by the KIDS install for the SR SURGERY application:

Note: The SR SURGERY client protocols logical links for a HLLP, two device, setup will need to be changed to SR SURGERY all others will remain as distributed by the KIDS build.

HL-7 Interface Workbench Jun 19, 1997 15:45:05 Page: 1 of 2
Client Protocols For
SR SURGERY

(1) SR Receiver of Unsolicited Requested Observation
Text: Receiver of Unsolicited Requested Observation
Package: SURGERY
Message Type Rcvd: ORU Sending Facility Required: YES
Message Type Sent: ACK Receiving Facility Required: YES
Event Type Rcvd/Sent: R01 Date/Time of Message Required: YES
Priority: IMMEDIATE Security Required: <DEFAULT>
Processing ID: PRODUCTION Logical Link: SR AAIS
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process Routine: D ^SRHLORU

(2) SR Receiver of Scheduling Query
Text: Receiver of Surgery Case Query
Package: SURGERY
Message Type Rcvd: QRY Sending Facility Required: NO
Message Type Sent: ZSQ Receiving Facility Required: NO
Event Type Rcvd/Sent: S25 Date/Time of Message Required: YES
Priority: IMMEDIATE Security Required: <DEFAULT>
Processing ID: PRODUCTION Logical Link: SR AAIS
Version ID: 2.1
Entry Action: <NONE>
Exit Action: <NONE>
Process Routine: D ^SRHLQRY

6. Confirm Logical Links for HLLP

Confirm that the following Logical Link entries were created by the KIDS install:

```

HL-7 Interface Workbench      Jun 19, 1997 15:58:43      Page:      7 of      7
                               Currently Defined
                               Logical Links

(11) SR AAIS
    LLP Parameter: SR HLLP
        LLP Type: HLLP (H)
        Device:
        Version ID: 2.2
        Block Size: <DEFAULT>
        Queue Size: <DEFAULT>
        Read Timeout: 3
        ACK Timeout: <DEFAULT>
        Re-transmission Attempts: 3
        Start Block Character: <DEFAULT>
        End Block Character: <DEFAULT>

(12) SR SURGERY      {NOTE : SR SURGERY is only used in a two device HLLP setup}
    LLP Parameter: SR HLLP2
        LLP Type: HLLP (H)
        Device:
        Version ID: 2.2
        Block Size: <DEFAULT>
        Queue Size: <DEFAULT>
        Read Timeout: 3
        ACK Timeout: <DEFAULT>
        Re-transmission Attempts: 3
        Start Block Character: <DEFAULT>
        End Block Character: <DEFAULT>

```

7. Confirm Logical Links for TCP/IP

```

(9) SR SURGERY
    LLP Parameter: SRSURGERY-TCP-RECV
    LLP Type: TCP
    Queue Size: <DEFAULT>
    LLP is not supported by
    the workbench. FileMan
    must be used to fill in
    the LLP specific info.

(10) SR AAIS
    LLP Parameter: SRAAISTSURGERY-TCP
    LLP Type: TCP
    Queue Size: <DEFAULT>
    LLP is not supported by
    the workbench. FileMan
    must be used to fill in
    the LLP specific info.

```

8. Edit SR AAIS and SR SURGERY Logical Link Device Field

At the "Select Tool:Jump To Next Link//" prompt choose the "EL Edit Logical Link". Add the device that your IRM created in the Pre-Installation Tasks section of the installation to the device field of the SR AAIS logical link.

```

Select Tool:Jump To Next Link// E1      Edit Logical Link
Select Logical Link (1-12): 11

Logical Link information
-----
Queue Size: <RET>
LLP Parameter: SR HLLP// <RET>

LLP Parameters
-----
LLP Type: HLLP// <RET>
Device:      < Enter the device that associates with the port that HL7 used
              to connect to the COTS system. For more information about
              setting up and activating this device, see Additional Technical

```

```

                                Information section.
Version ID: 2.2// <RET>
Block Size: <RET>
Read Timeout: 3// <RET>
ACK Timeout: <RET>
Re-transmission Attempts: 3// <RET>
Start Block Character: <RET>
End Block Character: <RET>.....

```

At the "Select Tool:Jump To Next Link/" prompt choose the "EL Edit Logical Link". Add the device that your IRM created in the Pre-Installation Tasks section of the installation to the device field of the SR SURGERY logical link.

```

Select Tool:Jump To Next Link// E1      Edit Logical Link
Select Logical Link (1-12): 12

Logical Link information
-----
Queue Size: <RET>
LLP Parameter: SR HLLP2// <RET>

LLP Parameters
-----
LLP Type: HLLP// <RET>
Device:      ← Enter the device that associates with the port that HL7 used
               to connect to the COTS system. For more information about
               setting up and activating this device, see Additional Technical
               Information section.

Version ID: 2.2// <RET>
Block Size: <RET>
Read Timeout: 3// <RET>
ACK Timeout: <RET>
Re-transmission Attempts: 3// <RET>
Start Block Character: <RET>
End Block Character: <RET>.....

```

9. Start LLP(s)

Use the HL7 Main Menu [HL MAIN MENU] option to start the LLP for entry SR AAIS and SR SURGERY:

For a one port/device HLLP connection, start SR AAIS.

```

Select Systems Manager Menu Option: HL7 Main Menu

Select HL7 Main Menu Option: 2  V1.6 OPTIONS

Select V1.6 OPTIONS Option: 1  Communications Server

Select Communications Server Option: 4  Start LLP
This option is used to launch the lower level protocol for the
Appropriate device. Please select the node with which you want
to communicate
Select HL LOGICAL LINK NODE: SR AAIS
Select one of the following:

```

F	FOREGROUND
B	BACKGROUND
Q	QUIT

Method for running the receiver: B// <RET> ACKGROUND
 Job was queued as 16655.

For a two port/device HLLP or MLLP (TCP/IP) connection, start SR SURGERY.

Select Communications Server Option: **4** Start LLP
 This option is used to launch the lower level protocol for the
 appropriate device. Please select the node with which you want
 to communicate
 Select HL LOGICAL LINK NODE: **SR SURGERY**
 Select one of the following:

F	FOREGROUND
B	BACKGROUND
Q	QUIT

Method for running the receiver: B// <RET> ACKGROUND
 Job was queued as 16658.

10. Flag Interface Fields

Using the option *Flag Interface Fields* [SRHL INTERFACE FIELDS], flag
 identifiers as appropriate to control the flow of data between the **VISTA** Surgery
 package and the ancillary system.

Additional Technical Information

New Protocol File (# 101) Entries

SR Unsolicited transmission of AAIS Requested Observation
SR Query for Scheduling Information
SR Receiver of Observation Unsolicited
SR Receiver of New Appointment Booking
SR Receiver of Appointment Rescheduling
SR Receiver of Appointment Modification
SR Receiver of Appointment Cancellation
SR Receiver of Appointment Deletion
SR Receiver of Master File Notification
SR Receiver of Staff Master File Notification
SR Unsolicited transmission of VISTA Requested Observation
SR Notification of Appointment Booking
SR Notification of Appointment Rescheduling
SR Notification of Appointment Modification
SR Notification of Appointment Cancellation
SR Notification of Appointment Deletion
SR Staff Master File Notification
SR Other Master File Notification
SR Receiver of Unsolicited Requested Observation
SR Receiver of Scheduling Query

New HL7 NON-DHCP APPLICATION PARAMETER (# 770) Entry

SR AAIS

New HL7 APPLICATION PARAMETER (# 771) Entries

SR AAIS
SR SURGERY

New LOGICAL LINKS Entries

SR AAIS
SR SURGERY

New MAIL GROUP (#3.8) Entry

SRHL DISCREPANCY

New/Modified Routines

SRHLDW	SRHLDW1	SRHLENV	SRHLMFN	SRHLOORU	SRHLORU	SRHLPOST
SRHLPRE	SRHLQRY	SRHLSCRN	SRHLU	SRHLUI	SRHLUO	SRHLUO1
SRHLUO2	SRHLUO3	SRHLUO4	SRHLUO4C	SRHLVOOR	SRHLVORU	SRHLVQRY
SRHLVU	SRHLVUI	SRHLVUI2	SRHLVUO	SRHLVUO1	SRHLVUO2	SRHLVUO4
SRHLVZIU	SRHLVZQR	SRHLVZSQ	SRHLXTMP	SRHLZIU	SRHLZQR	SROCODE
SROERR	SROERR0	SROFILE	SRONITE	SRSCHD2		

New Field in SURGERY SITE PARAMETERS File #133

HL7 V. 1.5 COMPATIBLE (34)

If applications communicating with the Surgery HL7 interface must use the interface designed for use with HL7 V. 1.5, enter YES. Otherwise, enter NO or leave this field blank.

New File

SURGERY INTERFACE PARAMETER file (#133.2)

The Surgery package uses this file for its HL7 interface with **VISTA** and Non-**VISTA** packages or systems. The file acts as a mapping and processing tool for receiving information from other applications and transmitting information from the Surgery package. This file is exported with data. The global root is “^SRO(133.2,”.

New Options

File Download [SRHL DOWNLOAD INTERFACE FILES]

This option is used to download Surgery interface files to the Automated Anesthesia Information System (AAIS). The process is currently being done by a screen capture to a file. In the future, this will be changed to a background task that can be queued to send HL7 master file updates.

Table Download [SRHL DOWNLOAD SET OF CODES]

This option downloads the SURGERY file (#130) set of codes to the Automated Anesthesia Information System (AAIS). This process is currently being done by a screen capture to a file. In the future, this will be changed to a background task that can be queued to send HL7 master file updates.

Surgery Interface Management Menu [SRHL INTERFACE]

This menu contains options that allow the user to set up certain interface parameters that control the processing of HL7 messages. The interface adheres to the Health Level 7 (HL7) protocol and forms the basis for the exchange of health care information between the **VISTA** Surgery package and any ancillary system.

Flag Interface Fields [SRHL INTERFACE FLDS]

This option allows the package coordinator to set the INTERFACE field in the SURGERY INTERFACE file (#133.2). The categories listed on the first screen correspond to entries in File 133.2. These categories are listed in the Surgery Interface Specifications document as being the OBR (OBServation Request) text identifiers. Each identifier corresponds to several fields in the **VISTA** Surgery package. This allows the user to control the flow of data between the **VISTA** Surgery package and the ancillary system on a field by field basis.

The option lists each identifier and its current setting. To receive the data coming from the ancillary system for a category, set the flag to R for receive. To ignore the data, set the flag to N for not receive. To see a second underlying layer of OBX text identifiers (SURGERY file #130 fields) and their settings, set the OBR text identifier to R for receive. The option will allow the user to toggle the settings for a range or for individual items.

Update Interface Parameter Field [SRHL PARAMETER]

This option may be used to edit the site parameter that determines which Surgery HL7 interface will be used, the interface compatible with **VISTA** HL7 V. 1.6 or the older one compatible with **VISTA** HL7 V. 1.5.

If applications communicating with the Surgery HL7 interface must use the interface designed for use with HL7 V. 1.5, enter YES. Otherwise, enter NO or leave this field blank.

Modified Options

Make Operation Requests [SROOPREQ]
Schedule Requested Operations [SRSCHD1]
Schedule Unrequested Operations [SROSRES]
Schedule Unrequested Concurrent Cases [SRSCHDC]
Cancel Scheduled Operation [SRSCAN]
Operation Information [SROMEN-OPINFO]
Surgical Staff [SROMEN-STAFF]
Operation Startup [SROMEN-START]
Operation [SROMEN-OP]
Post Operation [SROMEN-POST]
Operation (Short Screen) [SROMEN-OUT]
Surgeon's Verification of Diagnosis & Procedures [SROVER]
Anesthesia Information (Enter/Edit) [SROMEN-ANES]
Edit Non-O.R. Procedure [SRONOP-EDIT]
Flag Drugs for Use as Anesthesia Agents [SROCODE]
Update Site Configurable Files [SR UPDATE FILES]
Operation Information (Enter/Edit) [SROA OPERATION DATA]

Example VISTA Device File Entries

NAME: HL7 PORT	\$I: _LTA9048:
LOCATION OF TERMINAL: VISTA HL7 device	MARGIN WIDTH: 255
FORM FEED: #	PAGE LENGTH: 256
BACK SPACE: \$C(8)	SUBTYPE: P-OTHER
TYPE: TERMINAL	LAT SERVER PORT: LC-3-16
LAT SERVER NODE: DSV10	
NAME: NULL DEVICE	\$I: _NLA0:
LOCATION OF TERMINAL: null device	MARGIN WIDTH: 255
FORM FEED: #	PAGE LENGTH: 256
BACK SPACE: \$C(8)	SUBTYPE: P-OTHER
TYPE: TERMINAL	