

**Change pages for patch GMRV\*4\*7 of the Vitals/Measurements 4.0 user manual.**

Change pages for Nursing 4.0 user manual are contained in NURS\_4\_P23\_UM.PDF.

The following documentation change pages should be inserted before these replacement pages:

File Name:  
VITL\_4\_P1\_UM.PDF

Patch:  
GMRV\*4\*1

Patch GMRV\*4\*7 pages:

Replace Pages:  
iii-iv  
4.1-4.20  
5.1-5.18  
IN.1-IN.4

With Pages:  
iii-iv  
4.1-4.22  
5.1-5.14  
IN.1-IN.4



# Table of Contents

## Section 1 Package Management

Introduction .....	1
Chapter 1 Implementation and Maintenance .....	1.1
Description.....	1.1
Virgin Installation of software.....	1.1
Setting up the software environment.....	1.1
Name spacing and file listing .....	1.1
Editing site configurable files .....	1.2
Queueing TaskMan jobs.....	1.2
Accessing menus.....	1.2
Assigning menus.....	1.2
Printer issues.....	1.3
Non-Virgin Installation of software .....	1.4
Implementation Considerations .....	1.5
Resource Requirements.....	1.5
Chapter 2 Maintenance of Site Files.....	2.1
Vitals/Measurement .....	2.1
Vitals/Measurements Site Files Menu .....	2.2
Edit Vitals Site Parameter File .....	2.3
Change Default Qualifiers for Temp./Pulse.....	2.6
Enter/Edit Vitals Qualifiers.....	2.8
Display Vitals Category/Qualifier/Synonym Table .....	2.11
Display Administration Schedule File .....	2.14
Edit Administration Schedules File .....	2.15
Create Vital Measurement Quick Order Protocol.....	2.18

## Section 2 Package Operation

Chapter 3 Package Operation .....	3.1
Chapter 4 Enter/Edit Vitals/Measurements.....	4.1
Vitals/Measurement Data Entry .....	4.1
Edit a Vital/Measurement Entered in Error.....	4.18
Chapter 5 Vitals/Measurements Reports.....	5.1
Vitals/Measurements Results Reporting.....	5.1
V/M Graphic Reports.....	5.2
Vital Signs Record .....	5.4
B/P Plotting Chart .....	5.4
Weight Chart.....	5.4
Pulse Oximetry/Respiration Graph.....	5.4
Pain Chart.....	5.4
Latest Vitals Display for a Patient .....	5.5
Latest Vitals by Location .....	5.7
Cumulative Vitals Report .....	5.9
Print Vitals Entered in Error for a Patient.....	5.12
Glossary .....	GL.1
Index.....	IN.1

## Table of Contents

# Chapter 4 Enter/Edit Vitals/Measurements

## GMRV V/M ENTRY MENU

### Vitals/Measurement Data Entry

#### Description:

This option permits users to enter different types of patient vitals signs and selected measurements into the clinical patient record system. Vitals/Measurements data is stored in the GMRV Vital Measurement (#120.5) file.

#### Additional Information:

Various vitals and other measurement data can be retrieved through any of the Vitals/Measurements Results Reporting menu options. If data was entered in error, the user can correct the misinformation through the Edit a Vitals/Measurements Entered in Error option.

Patient vital signs can be entered by patient, selected room and bed and by ward location (i.e., Hospital Location (#44) file). When vitals are entered, both a date and a time are required as well as hospital location. Vital signs and other patient measurements are documented in the specific order indicated by the vital type (temperature, weight, pulse oximetry). Data representing each vital type entry is also separated by a hyphen: Temp-Pulse-Resp-B/P-Height-Weight. To omit entering a vital/measurement reading:

Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the reason for omission.

The software supports the documentation of measurement omissions by displaying three reasons (i.e., unavailable, pass, and refused) when an '\*' is entered for a measurement value.

Pain tracking is included in the software, where a user can document the level of pain that the patient is experiencing by using the following pain scale:

- 0 No pain
- 1 (Minimal discomfort/pain) through 10 (worst imaginable pain)
- 99 Unable to respond

Qualifiers (alpha characters appended to the measurement's numeric value) are used to provide additional information about the vital sign(s) and measurements taken by the clinician. Qualifiers are grouped by categories such as: location (e.g., right arm, left leg), position (e.g., lying, sitting, standing), method (e.g., cuff, Doppler, auscultate, assisted ventilator, T-piece), site (e.g., right, left), quality (e.g., actual, estimated), and cuff size (e.g., adult, small adult, pediatric). Vital types can be printed by category with their associated qualifiers and synonyms through the Display Vitals Category/Qualifier/Synonym Table option. When documenting patient vital signs and measurements, the user may see a listing of valid qualifiers by entering a '?' or '??'. The application coordinator has the ability to add qualifiers and change other software parameters through the Vitals/Measurements Site Files Menu option.

Note: After entering vitals/measurements information for a patient, the software displays this data and asks you, "Is this correct? YES/". The user must answer YES to store the data after it is verified correct. If the user enters an '^' at the prompt, the application does NOT store the data and displays 'DATA DELETED' on the screen. Of course, if NO is entered after the prompt, patient information is NOT stored and all vitals/measurements data must be re-entered.

Defaults: The software has defaults for each vital/measurement except for blood pressure. These defaults are displayed when the clinician validates the measurements entered. The following list summarizes the software's defaults:

- Temperature - Oral
- Pulse - Radial
- Respirations - Spontaneous
- Height - Actual
- Weight - Actual

An application coordinator can change the temperature and pulse defaults through the Change Default Qualifiers for Temp./Pulse option.

Asterisks have two different meanings in the software.

1. When entering vitals/measurements information, the user may enter the \* to prompt the software to display reasons for not documenting a measurement. The user then selects the most appropriate entry.
2. When an \* appears in a report appended to a vitals/ measurement entry, the value of the measurement is either higher or lower than the acceptable value defined in the GMRV Vitals Parameters (#120.57) file. These file values can be edited by the package coordinator through the Edit Vitals Site Parameter File option.

All patient measurements are stored in the U.S. Customary System in File #120.5 and displayed in both U.S. customary units and metric equivalents. Small discrepancies between the entering and displaying of data is due to the conversion between scales.

Temperature scales are automatically defined by the software when a user does not append a 'F' or a 'C' to the numeric entry. A number that is larger than or equal to 45 will be considered as a Fahrenheit temperature, and any number less than 45 will be considered as a Centigrade temperature.

The scales for weight and height cannot be automatically determined. Therefore, the user must enter the scale (U.S. customary unit or metric equivalent) associated with the weight and height.

### Menu Display:

```
Select OPTION NAME:  GMRVMGR                      Vitals/Measurement

1      Vitals/Measurement Data Entry ...
2      Vitals/Measurements Results Reporting ...
3      Edit a Vital/Measurement Entered in Error
4      Vitals/Measurements Site Files Menu ...
```

### Screen Prints:

```
Select Vitals/Measurement Option: 1  Vitals/Measurement Data Entry
ENTER DATE (TIME Required) VITALS WERE TAKEN:  T@8  (APR 04, 1993@08:00)
```

The first prompt is a system parameter used to identify the date and time the vitals/ measurements were taken; time is required. The date/time is stamped on every vital/measurement that is entered into the system. The Change Date/Time Taken option allows the user to identify a new date/time the measurements were taken, a different vital type, or a new selection of patients without leaving the Vitals/ Measurement Data Entry menu. This Change Date/Time Taken option is further explained later in this chapter.

```
1      TPR
2      TPR B/P
3      TPR B/P, Ht and Wt.
4      TPR, B/P and Wt.
5      Temp, Detailed PR and B/P
6      Detailed B/P and Associated Pulse
7      Pulse
8      Weight
9      Circumference/Girth
10     Pulse Oximetry
11     CVP (Central Venous Pressure)
12     User Configurable Combination
13     Change Date/Time Taken
14     Pain
```

```
Select Vitals/Measurement Data Entry Option: 1  TPR
Vitals by (A)ll patients on a unit, (S)electd Rooms on unit, or
(P)atient? P
```

**Example: Process for entering vitals for an individual patient:**

Enter patient name.

## Enter/Edit Vitals/Measurements

Select PATIENT NAME: **PARKER, PETER P.**                      --                      500-1                      YES  
SC VETERAN

Enter appropriate hospital location.

Select Hospital Location: 4E// **<RET>**

To omit entering a vital/measurement reading:  
Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the reason for omission.

The temp-pulse-resp are entered at this prompt.

Temp-Pulse-Resp: **102.1-70-24**

The following is representative of the software's on-line help:

**\*\* Temperature:** a number + a location qualifier(optional), e.g., 99.6A or 99.6AX (F AXILLARY) or 37R (C RECTAL).

Default: ORAL

AXILLARY, CORE, ORAL, RECTAL, SKIN, TYMPANIC

**\*\* Pulse:** a number + a location qualifier(optional), e.g., 70A or 70AP (70 APICAL).

Default: RADIAL

APICAL, BILATERAL PERIPHERALS, BRACHIAL, CAROTID, DORSALIS PEDIS, FEMORAL, OTHER, PERIPHERAL, POPLITEAL, POSTERIOR TIBIAL, RADIAL, ULNAR

**\*\* Respiration:** a number(0-99 with no decimal places) + method qualifier (optional) e.g., 40C or 40c.

Default: SPONTANEOUS

ASSISTED VENTILATOR, CONTROLLED VENTILATOR, SPONTANEOUS

The entry should be in the following format:

TTT.T-PPP-RR

To omit entering a vital/measurement reading:

Enter 'N' or 'n' for the value when NOT documenting a reason for omission.

Enter an \* for the specific value when documenting the reason for omission.

Enter a single \* to document that all measurements were omitted and the reason for omission.

The data entry is one single string of alphanumeric characters with each measurement separated by a hyphen (-). The software redisplay the data as:

Temp.: 102.1 F (38.9 C)\* ORAL  
Pulse: 70 RADIAL  
Resp.: 24 SPONTANEOUS

Note the asterisk (\*) following the value for temperature. An asterisk (\*) following data indicates that the value is higher or lower than the acceptable values which are defined in the site configurable file GMRV Vitals Parameters (#120.57).

Is this correct? YES// **<RET>**



At this point the user can reject or accept the entry. This is the last chance to change any possible errors during data entry. The user is able to correct errors at a later time using the error correction option (Edit a Vital/Measurement Entered in Error) which will be discussed later in this chapter.

Do you want to enter other V/M data for this patient? No// <RET> (No)

Select PATIENT NAME: <RET>

This option also supports entering patient vital signs and measurements for all patients on a ward or unit and by selected hospital rooms on a unit. Navigating through the software to document measurements for all patients on a unit is described next.

- 1 TPR
- 2 TPR B/P
- 3 TPR B/P, Ht and Wt.
- 4 TPR, B/P and Wt.
- 5 Temp, Detailed PR and B/P
- 6 Detailed B/P and Associated Pulse
- 7 Pulse
- 8 Weight
- 9 Circumference/Girth
- 10 Pulse Oximetry
- 11 CVP (Central Venous Pressure)
- 12 User Configurable Combination
- 13 Change Date/Time Taken
- 14 Pain

Select Vitals/Measurement Data Entry Option: 1 TPR

Select the measurement(s) you want to enter.

**Example: Entering vitals for all of the patients on a particular unit.**

Vitals by (A)ll patients on a unit, (S)electd Rooms on unit, or (P)atient? **A**

Enter the appropriate ward/unit name.

Select Unit: **SICU**

Begin entering patient vitals.

YGAFO,PAUL E 550-1 OK? YES// <RET>

The clinician may enter patient information in the format described for the individual patient or bypass an entry by entering <RET>. The software indicates that no data was entered. The name of the next patient is displayed.

Temp-Pulse-Resp: <RET>

NO DATA ENTERED

Do you wish to stop looping through names? YES// **N** (NO)

When the clinician enters a NO behind the patient's name, he may continue to the next patient by answering NO to the question 'Do you wish to stop looping through names?', or exit by entering YES.

## Enter/Edit Vitals/Measurements

KJER,JAMES E                    550-2    OK? YES// **N**    (NO)

KCKK,GEORGE R                560-1    OK? YES// **<RET>**    (YES)

Temp-Pulse-Resp: **100.2R-60A-34C**

Temp.: 100.2 F    (37.9 C) RECTAL

Pulse: 60\* APICAL

Resp.: 34\* CONTROLLED VENTILATOR

Is this correct? YES// **<RET>**

GEVFCC,MICHAEL E            560-10    OK? YES// **N**    (NO)

Do you wish to stop looping through names? YES// **<RET>**    (YES)

Enter return to continue **<RET>**

- 1        TPR
- 2        TPR B/P
- 3        TPR B/P, Ht and Wt.
- 4        TPR, B/P and Wt.
- 5        Temp, Detailed PR and B/P
- 6        Detailed B/P and Associated Pulse
- 7        Pulse
- 8        Weight
- 9        Circumference/Girth
- 10       Pulse Oximetry
- 11       CVP (Central Venous Pressure)
- 12       User Configurable Combination
- 13       Change Date/Time Taken
- 14       Pain

Select Vitals/Measurement Data Entry Option: **3**    TPR B/P, Ht and Wt.

**Example: Entering vital signs by selecting specific beds. The software prompts for entering information is identical with the previous two pathways, i.e., entering data by patient or by unit.**

Vitals by (A)ll patients on a unit, (S)elected Rooms on unit, or (P)atient? **S**

Select Unit: **SICU**

**Enter appropriate ward location.**

Ward SICU has the following rooms:

- 1.       550-1
- 2.       550-2
- 3.       550-3
- 4.       550-4
- 5.       560-1
- 6.       560-10
- 7.       560-2
- 8.       570-1

Select the NUMBER(S) of the rooms: **5-6**

**Enter room or range of rooms to enter patient vitals.**

Begin entering patient vitals.

KCKK,GEORGE R

560-1 OK? YES// &lt;RET&gt;

The software loops through the patients located in Rooms 560-1 and 560-10 and the user enters the vitals/measurements values.

- 1 TPR
- 2 TPR B/P
- 3 TPR B/P, Ht and Wt.
- 4 TPR, B/P and Wt.
- 5 Temp, Detailed PR and B/P
- 6 Detailed B/P and Associated Pulse
- 7 Pulse
- 8 Weight
- 9 Circumference/Girth
- 10 Pulse Oximetry
- 11 CVP (Central Venous Pressure)
- 12 User Configurable Combination
- 13 Change Date/Time Taken
- 14 Pain

### Example: Entering detailed B/P and associated pulse.

Select Vitals/Measurement Data Entry Option: **6** Detailed B/P and Associated Pulse

Vitals by (A)ll patients on a unit, (S)elected Rooms on unit, or (P)atient? **P**

Select PATIENT NAME: **PARKER, PETER P.** -- 500-1 YES  
SC VETERAN

Select Hospital Location: 4E// <RET> 4E

To omit entering a vital/measurement reading:

Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the reason for omission.

BP-Pulse: **120/70-60**

The software displays selections for documenting detailed blood pressures. Each qualifier must be separated by a comma. The user is not required to enter qualifiers to proceed to the next prompt.

Qualifiers for BLOOD PRESSURE:

LOCATION	POSITION	METHOD	CUFF SIZE
1 L ARM	5 LYING	8 CUFF	12 ADULT
2 L LEG	6 SITTING	9 DOPPLER	13 LG ADULT
3 R ARM	7 STANDING	10 NON-INVASIVE	14 PEDIATRIC
4 R LEG		11 PALPATED	15 SM ADULT
			16 THIGH

Select a number under each category (optional).

Separate the numbers with ',': **1,6,8**

L ARM SITTING CUFF

## Enter/Edit Vitals/Measurements

The software displays selections for documenting detailed pulse.

Qualifiers for PULSE:

SITE	METHOD	POSITION
1 LEFT	3 AUSCULTATE	6 LYING
2 RIGHT	4 DOPPLER	7 SITTING
	5 PALPATED	8 STANDING

Select a number under each category (optional).

Separate the numbers with ',': **1**

LEFT

B/P: 120/70 L ARM SITTING CUFF

Pulse: 60\* RADIAL LEFT

Is this correct? YES// **<RET>**

After reviewing the entry, a user can reject or accept the data. This is the last chance to change any possible errors during data entry. The user is able to correct errors at a later time using the error correction option, Edit a Vital/Measurement Entered in Error, which will be discussed later in this chapter.

Enter another B/P? NO// **<RET>**

If YES is entered, the user can enter another B/P and pulse under the same time for a different position, location, etc..

1	TPR
2	TPR B/P
3	TPR B/P, Ht and Wt.
4	TPR, B/P and Wt.
5	Temp, Detailed PR and B/P
6	Detailed B/P and Associated Pulse
7	Pulse
8	Weight
9	Circumference/Girth
10	Pulse Oximetry
11	CVP (Central Venous Pressure)
12	User Configurable Combination
13	Change Date/Time Taken
14	Pain

### Example: Girth measurement.

Select Vitals/Measurement Data Entry Option: **9** Circumference/Girth

Vitals by (A)ll patients on a unit, (S)electd Rooms on unit, or

(P)atient? **P**

Select PATIENT NAME: **PARKER, PETER P.** -- 500-1 YES  
SC VETERAN

Select Hospital Location: 4E// **<RET>**

To omit entering a vital/measurement reading:  
 Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
 Enter an \* for the specific value when documenting the reason for omission.  
 Enter a single \* to document that all measurements were omitted and the reason for omission.

The following is representative of the software's on-line help:

**\*\* Circumference: a number + 'I' or 'C' (2 decimals allowed)**

For example: 72.25I (inches) 147C (centimeters)

The entry should be in the following format:

NNN.NN

Circumference/Girth: **32.25I**

Qualifiers for CIRCUMFERENCE/GIRTH:

LOCATION	SITE
1 ABDOMINAL	10 LEFT
2 ANKLE	11 RIGHT
3 CALF	
4 HEAD	
5 LOWER ARM	
6 OTHER	
7 THIGH	
8 UPPER ARM	
9 WRIST	

The user identifies the type of circumference.

Enter a number under each category, separate numbers with a ','  
 DO NOT select SITE if this is a HEAD/ABDOMINAL girth measurement: **1**  
 ABDOMINAL

Circumference/Girth: 32.25 IN (81.92 CM) ABDOMINAL  
 Is this correct? YES// **Y**

**Example: Circumference of left lower arm.**

To omit entering a vital/measurement reading:  
 Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
 Enter an \* for the specific value when documenting the reason for omission.  
 Enter a single \* to document that all measurements were omitted and the reason for omission.

Circumference/Girth: **12.25I**

Qualifiers for CIRCUMFERENCE/GIRTH:

LOCATION	SITE
1 ABDOMINAL	10 LEFT
2 ANKLE	11 RIGHT
3 CALF	

## Enter/Edit Vitals/Measurements

4 HEAD  
5 LOWER ARM  
6 OTHER  
7 THIGH  
8 UPPER ARM  
9 WRIST

Enter a number under each category, separate numbers with a ','  
DO NOT select SITE if this is a HEAD/ABDOMINAL girth measurement: 5,10  
LOWER ARM LEFT

Circumference/Girth: 12.25 IN (31.12 CM) LOWER ARM LEFT  
Is this correct? YES// <RET>

Do you want to enter other V/M data for this patient? No// <RET> (No)

1 TPR  
2 TPR B/P  
3 TPR B/P, Ht and Wt.  
4 TPR, B/P and Wt.  
5 Temp, Detailed PR and B/P  
6 Detailed B/P and Associated Pulse  
7 Pulse  
8 Weight  
9 Circumference/Girth  
10 Pulse Oximetry  
11 CVP (Central Venous Pressure)  
12 User Configurable Combination  
13 Change Date/Time Taken  
14 Pain

### Example: Pulse Oximetry.

Select Vitals/Measurement Data Entry Option: 10 Pulse Oximetry  
Vitals by (A)ll patients on a unit, (S)elected Rooms on unit, or  
(P)atient? P  
Select PATIENT NAME: PARKER, PETER P. -- 500-1 YES  
SC VETERAN

### Answer with patient name.

Select Hospital Location: 4E// <RET> 4E

### Enter appropriate hospital location.

To omit entering a vital/measurement reading:  
Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the  
reason for omission.

The following is representative of the software's on-line help:

**\*\* Pulse Oximetry:** Enter the numeric value of the patients O2 saturation obtained via pulse oximetry. The value will be interpreted as a percentage.

The entry should be in the following format:

NNN

PO2: 89

At this prompt the user may indicate that patient was receiving oxygen at the time the pulse oximetry reading was taken. If YES is entered, additional information on the flow rate, O2 concentration, and method of O2 administration can be entered.

Is the patient on supplemental oxygen? Yes// <RET> (Yes)

Enter a number between 0.5-20 for liters/minute, and/or 21-100 for percent of oxygen concentration. If you wish to enter both rates, separate the values with a ','.

Enter the numeric value(s) for amount of supplemental oxygen provided (Separate values with a ','): 3.5,40

Oxygen is supplied by

- 1 AEROSOL/HUMIDIFIED MASK
- 2 FACE TENT
- 3 MASK
- 4 NASAL CANNULA
- 5 NON RE-BREATHING
- 6 PARTIAL RE-BREATHING
- 7 T-PIECE
- 8 TRACHEOSTOMY COLLAR
- 9 VENTILATOR
- 10 VENTURI MASK

Enter a number: 1

Pulse Ox.: 89 % with supplemental O2 40% 3.5 l/min  
via AEROSOL/HUMIDIFIED MASK

Is this correct? YES// <RET>

Do you want to enter other V/M data for this patient? No// <RET> (No)

- 1 TPR
- 2 TPR B/P
- 3 TPR B/P, Ht and Wt.
- 4 TPR, B/P and Wt.
- 5 Temp, Detailed PR and B/P
- 6 Detailed B/P and Associated Pulse
- 7 Pulse
- 8 Weight
- 9 Circumference/Girth
- 10 Pulse Oximetry
- 11 CVP (Central Venous Pressure)
- 12 User Configurable Combination
- 13 Change Date/Time Taken
- 14 Pain

### Example: Central Venous Pressure (CVP).

Select Vitals/Measurement Data Entry Option: **11** CVP (Central Venous Pressure)  
Vitals by (A)ll patients on a unit, (S)electd Rooms on unit, or (P)atient? **P**  
Select PATIENT NAME: **PARKER, PETER P.** -- 500-1 YES  
SC VETERAN

Select Hospital Location: 4E// **<RET>** 4E

To omit entering a vital/measurement reading:  
Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the reason for omission.

The following is representative of the software's on-line help:

**\*\* Central venous pressure: a number for cmH2O measurement or a number + 'G' for mmHg measurement (1 decimal allowed).**

**A negative number can be entered up to and including -13 cmH2O or -9.6 mmHg. The entry should be in the following format:**

**NN**

**CVP: 14G**

CVP: 19.04 cmH2O (14.0 mmHg)  
Is this correct? YES// **<RET>**

Do you want to enter other V/M data for this patient? No// **<RET>** (No)

- 1 TPR
- 2 TPR B/P
- 3 TPR B/P, Ht and Wt.
- 4 TPR, B/P and Wt.
- 5 Temp, Detailed PR and B/P
- 6 Detailed B/P and Associated Pulse
- 7 Pulse
- 8 Weight
- 9 Circumference/Girth
- 10 Pulse Oximetry
- 11 CVP (Central Venous Pressure)
- 12 User Configurable Combination
- 13 Change Date/Time Taken
- 14 Pain

### Example: User Configurable Combination.

Select Vitals/Measurement Data Entry Option: **12** User Configurable Combination

This option allow the users to create their own vital signs and measurements combination(s). Although these selections display all prompts that support documentation of detailed measurements, the user can bypass these prompts.



Select the combination of Vitals/Measurements you want to enter.

- 1 T
- 2 P
- 3 R
- 4 B/P
- 5 Wt
- 6 Ht
- 7 Circumference/Girth
- 8 Pulse Oximetry

Data on weight and height is documented on a patient.

Select from 1 to 9 (enter 1,3-5 etc.) or ^ to exit: **5,6**

Vitals by (A)ll patients on a unit, (S)electd Rooms on unit, or  
(P)atient? **P**

Select PATIENT NAME: **PARKER, PETER P.** -- 500-1 YES  
SC VETERAN

Select Hospital Location: 4E// **<RET>**

To omit entering a vital/measurement reading:

Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the  
reason for omission.

The following is representative of the software's on-line help:

**\*\* Weight: a number+L(l) or a number+K(k) + 2nd quality (optional) 2  
decimals allowed, e.g., 120.25L or 120l (pounds),  
45.25K or 45k (kilograms).**

**120L, 120LA, or 120K for actual weight.**

**Default: ACTUAL**

**ACTUAL, DRY, ESTIMATED**

**\*\* Height: a number + ',' + qualifier (2 decimals allowed)**

**72I (inches) 147C or 147CM (centimeters)**

**5F10I;E or 5FT10IN;E or 5'10";E (for 5 feet 10 inches ESTIMATED)**

**Default: ACTUAL**

**ACTUAL, ESTIMATED**

The entry should be in the following format:

**WWW.WW-HH**

**Ht.-Wt.: 5'6";E-165.25LE**

Qualifiers for WEIGHT:

METHOD

- 1 BED
- 2 CHAIR
- 3 STANDING

## Enter/Edit Vitals/Measurements

The user may enter a method indicating how the weight was taken or enter <RET> to display, and then validate data.

Select a number under each category (optional).

Separate the numbers with ',': 1

BED

Wt.: 165.25 LB (75.11 KG) ESTIMATED BED

Ht.: 5 FT 6 IN (167.64 CM) ESTIMATED

Is this correct? YES// <RET>

Do you want to enter other V/M data for this patient? No// <RET> (No)

- 1 TPR
- 2 TPR B/P
- 3 TPR B/P, Ht and Wt.
- 4 TPR, B/P and Wt.
- 5 Temp, Detailed PR and B/P
- 6 Detailed B/P and Associated Pulse
- 7 Pulse
- 8 Weight
- 9 Circumference/Girth
- 10 Pulse Oximetry
- 11 CVP (Central Venous Pressure)
- 12 User Configurable Combination
- 13 Change Date/Time Taken
- 14 Pain

### Example: Change Date/Time Taken.

Select Vitals/Measurement Data Entry Option: 7 Pulse

When information is entered using the All patients on a unit or the Selected Rooms on a unit pathway, the user may enter new data without exiting and re-entering the option by selecting (13) Change Date/Time Taken. This option prompts for a new date/time measurements were taken, a different vital type, and a new selection of patients. As an example, the user documents 6 am pulse readings for 2AS patients, then enters the Change Date/Time Taken option and documents 8 am weights.

Vitals by (A)ll patients on a unit, (S)electd Rooms on unit, or (P)atient? A

Select Unit: 2AS

Begin entering patient vitals.

ANESMJS, ESTELL H

OK? YES// <RET>

To omit entering a vital/measurement reading:

Enter 'N' or 'n' for the value when NOT documenting a reason for omission.

Enter an \* for the specific value when documenting the reason for omission.

Enter a single \* to document that all measurements were omitted and the reason for omission.

Pulse: 86

Qualifiers for PULSE:

SITE	METHOD	POSITION
1 LEFT	3 AUSCULTATE	6 LYING
2 RIGHT	4 DOPPLER	7 SITTING
	5 PALPATED	8 STANDING

Select a number under each category (optional).

Separate the numbers with ',': **1,3,7**

LEFT AUSCULTATE SITTING

Pulse: 86 RADIAL LEFT AUSCULTATE SITTING

Is this correct? YES// **<RET>**

.  
.  
.

The user continues to enter 6 am pulse readings, until he selects the Change Date/Time Taken option to enter 8 am weights.

Do you wish to stop looping through names? YES// **<RET>**

Enter return to continue **<RET>**

1	TPR
2	TPR B/P
3	TPR B/P, Ht and Wt.
4	TPR, B/P and Wt.
5	Temp, Detailed PR and B/P
6	Detailed B/P and Associated Pulse
7	Pulse
8	Weight
9	Circumference/Girth
10	Pulse Oximetry
11	CVP (Central Venous Pressure)
12	User Configurable Combination
13	Change Date/Time Taken
14	Pain

At the next prompt, the user enters, 13, to enter another date/time, new vital types and select patient(s).

Select Vitals/Measurement Data Entry Option: **13** Change Date/Time Taken

ENTER DATE (TIME Required) VITALS WERE TAKEN: **T@8:00** (MAR 03, 1997@08:00)

1	TPR
2	TPR B/P
3	TPR B/P, Ht and Wt.
4	TPR, B/P and Wt.
5	Temp, Detailed PR and B/P
6	Detailed B/P and Associated Pulse
7	Pulse
8	Weight
9	Circumference/Girth
10	Pulse Oximetry
11	CVP (Central Venous Pressure)
12	User Configurable Combination
13	Change Date/Time Taken
14	Pain

## Enter/Edit Vitals/Measurements

Select Vitals/Measurement Data Entry Option: **8** Weight  
Vitals by (A)ll patients on a unit, (S)elected Rooms on unit, or (P)atient? **A**  
Select Unit: **2AS**  
Begin entering patient vitals.

ANESMJS,ESTELL H OK? YES// **<RET>**

To omit entering a vital/measurement reading:  
Enter 'N' or 'n' for the value when NOT documenting a reason for omission.  
Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the reason for omission.

Wt.: **400L**

Qualifiers for WEIGHT:

METHOD

- 1 BED
- 2 CHAIR
- 3 STANDING

Select a number under each category (optional).  
Separate the numbers with ',': **2**  
CHAIR

Wt.: 400 LB (181.82 KG) ACTUAL CHAIR  
Is this correct? YES// **<RET>**

The user continues to select the patients that have 8 am weight readings and exits when documentation is completed.

- 1 TPR
- 2 TPR B/P
- 3 TPR B/P, Ht and Wt.
- 4 TPR, B/P and Wt.
- 5 Temp, Detailed PR and B/P
- 6 Detailed B/P and Associated Pulse
- 7 Pulse
- 8 Weight
- 9 Circumference/Girth
- 10 Pulse Oximetry
- 11 CVP (Central Venous Pressure)
- 12 User Configurable Combination
- 13 Change Date/Time Taken
- 14 Pain

### Example: Pain.

Select Vitals/Measurement Data Entry Option: **14** Pain  
Vitals by (A)ll patients on a unit, (S)elected Rooms on unit, or (P)atient? **P**  
Select PATIENT NAME: **PARKER, PETER P.** -- 500-1 YES  
SC VETERAN

Select Hospital Location: 4E// **<RET>** 4E

Enter an \* for the specific value when documenting the reason for omission.  
Enter a single \* to document that all measurements were omitted and the reason for omission.

The following is representative of the software's on-line help:

**\*\* Pain:** Enter a numeric value. 0 for no pain. 1-10 for pain scale  
10 is worst imaginable pain or enter 99 if the patient is unable to respond  
For example: 99 or 5

The entry should be in the following format:

NN

Pain: 3

Pain: 3  
Is this correct? YES// <RET>

Do you want to enter other V/M data for this patient? No// <RET> (No)

### Menu Access:

The Vitals/Measurement Data Entry option is accessed through the Vitals/Measurement menu.

## **GMRV ERROR EDIT**

### **Edit a Vital/Measurement Entered in Error**

#### **Description:**

This option allows users to correct errors in vitals/measurements. When data is edited, a new record is created, and the old record is marked entered in error. Vitals/measurements data is stored in the GMRV Vital Measurement (#120.5) file.

#### **Additional Information:**

Vitals/measurements data can be retrieved through the Print Vitals Entered in Error for a Patient option. The information entered through this option is also reflected in the Cumulative Vitals Report.

#### **Menu Display:**

```
Select OPTION NAME:  GMRVMGR                      Vitals/Measurement

1      Vitals/Measurement Data Entry ...
2      Vitals/Measurements Results Reporting ...
3      Edit a Vital/Measurement Entered in Error
4      Vitals/Measurements Site Files Menu ...
```

#### **Screen Prints:**

```
Select Vitals/Measurement Option: 3  Edit a Vital/Measurement Entered in Error

Select PATIENT NAME:  PARKER,PETER P.           08-04-18       931848807       SC
VETERAN
```

Answer with patient name.

```
Select DATE (TIME optional) of this error: T  (APR 04, 1993)
```

Enter an appropriate date. When the exact time is not remembered, the user may enter the date, and all entries for that date display.

CHOOSE FROM:

- 1 APR 4, 1993@08:00  
BP: 120/80 (R ARM/LYING)  
T: 102.1 F (38.9 C)\* (ORAL)  
R: 24 (SPONTANEOUS)  
P: 70 (RADIAL)
- 2 APR 4, 1993@08:01  
BP: 130/85 (L ARM/LYING)

Select VITAL/MEASUREMENT from the above list: **1**

Enter the number of the measurement or group of measurements you want to enter in error.

A detailed summary of the data is provided.

```
*** DATA TO BE ENTERED IN ERROR ***
APR  4, 1993@08:00
BP:  120/80 R ARM/LYING
T:   102.1 F (38.9 C)* ORAL
R:   24 (SPONTANEOUS)
P:   70 RADIAL
```

Select a VITAL type record, or press return to correct DATE/TIME, PATIENT or set INVALID for the record: **RESPIRATION**

In this example, the respirations were incorrectly documented. The user must select an appropriate reason(s) from 1-4.

In this case, the invalid record selection (4) cannot be used in combination with 1, 2, or 3. The user receives the following prompt.

- 1 INCORRECT DATE/TIME
- 2 INCORRECT READING
- 3 INCORRECT PATIENT
- 4 INVALID VITAL/RECORD

The INVALID RECORD reason cannot be used in combination with any other selections.

The following are examples of identifying incorrect patient data using the above four selections:

**INCORRECT DATE/TIME:** changes the date/time that a vital measurement was taken for a patient.

```
Select the reason(s) for entering this vital/measurement in error: 1
Enter new date/time vitals were taken: 9A (APR 04, 1993@09:00)
```

**INCORRECT READING:** changes the data for a vital/measurement that was taken for a patient. The following prompt displays only if a single/specific vital has been identified for editing.

```
Select the reason(s) for entering this vital/measurement in error: 2
```

```
NEW RESPIRATION: 20A
```

..
The corrected information is entered.

**INCORRECT PATIENT:** changes the patient for a particular vital measurement (e.g., the data was entered for the wrong patient).

## Enter/Edit Vitals/Measurements

Select the reason(s) for entering this vital/measurement in error: **3**  
Select the NEW Patient's name: **STARK, ANTHONY L.** 01-19-25 589632111  
NSC VETERAN

**INVALID VITAL/RECORD:** deletes the record without adding any new data. In actuality, the bad record is still in the computer, but is not a part of the current patient data.

Select the reason(s) for entering this vital/measurement in error: **4**

If a combination of the above applies, choices can be selected at one time as described below.

Select PATIENT NAME: **GEVFCC, MICHAEL E** 03-22-22 541969342 NSC  
VETERAN

Select DATE (TIME optional) of this error: **T** (APR 04, 1993)  
CHOOSE FROM:

1 APR 4, 1993@08:00  
BP: 120/80  
T: 100.4 F (38.0 C) (ORAL)  
R: 20 (SPONTANEOUS)  
P: 80 (RADIAL)  
HT: 5 FT 8 IN (172.72 CM)  
WT: 200 LB (90.91 KG) (ACTUAL)

Select VITAL/MEASUREMENT from the above list: **1**

\*\*\* DATA TO BE ENTERED IN ERROR \*\*\*

APR 4, 1993@08:00  
BP: 120/80  
T: 100.4 F (38.0 C) ORAL  
R: 20 (SPONTANEOUS)  
P: 80 RADIAL  
HT: 5 FT 8 IN (172.72 CM)  
WT: 200 LB (90.91 KG) ACTUAL

Select a VITAL type record, or press return to correct DATE/TIME, PATIENT  
or set INVALID for the record: **RESPIRATION**

- 1 INCORRECT DATE/TIME
- 2 INCORRECT READING
- 3 INCORRECT PATIENT
- 4 INVALID VITAL/RECORD

Select the reason(s) for entering this vital/measurement in error: **1-3**

Enter new date/time vitals were taken: **9A** (APR 04, 1993@09:00)

NEW RESPIRATION: **24A**

Select the NEW Patient's name: **STARK, ANTHONY L.** 01-19-25 589632111  
NSC VETERAN



Select DATE (TIME optional) of this error: **T** (APR 05, 1993)  
CHOOSE FROM:

1 APR 5, 1993@08:00  
BP: 130/80 (L ARM/LYING)  
T: 101.2 F (38.4 C) (RECTAL)  
R: 20 (ASSISTED VENTILATOR)  
P: 80 (APICAL)

Select VITAL/MEASUREMENT from the above list: **1**

\*\*\* DATA TO BE ENTERED IN ERROR \*\*\*

APR 5, 1993@08:00  
BP: 130/80 L ARM/LYING  
T: 101.2 F (38.4 C) RECTAL  
R: 20 ASSISTED VENTILATOR  
P: 80 APICAL

Select a VITAL type record, or press return to correct DATE/TIME, PATIENT  
or set INVALID for the record: **<RET>**

1 INCORRECT DATE/TIME  
2 INCORRECT PATIENT  
3 INVALID VITAL/RECORD

Select the reason(s) for entering these vital/measurements in error:

### Menu Access:

The Edit a Vital/Measurement Entered in Error option is accessed through the Vitals/Measurement menu.



# Chapter 5 Vitals/Measurements Reports

## GMRV PRINT MENU

### Vitals/Measurements Results Reporting

#### Description:

This menu contains options which allow the users to print various types of vitals/measurements reports. Data printed on these reports is stored in the GMRV Vital Measurement (#120.5) file.

#### Additional Information:

Data printed in these reports was entered through the Vitals/Measurement Data Entry option, or the Edit a Vital/Measurement Entered in Error option. All abnormal patient values are identified on the vitals/measurements reports by an asterisk (\*). Patient vital signs and measurements can be printed by patient, selected room and bed and by ward location (i.e., entry in the Hospital Location (#44) file). All vitals/measurements data for a patient is reflected on the Cumulative Vitals Report and the Vital Signs Record (SF 511) graphic report. The division name is obtained from the Hospital Location (#44) file.

#### Menu Display:

```
Select OPTION NAME:  GMRVMGR                      Vitals/Measurement

1      Vitals/Measurement Data Entry ...
2      Vitals/Measurements Results Reporting ...
3      Edit a Vital/Measurement Entered in Error
4      Vitals/Measurements Site Files Menu ...

Select Vitals/Measurement Option:  2  Vitals/Measurements Results Reporting

1      V/M Graphic Reports
2      Latest Vitals Display for a Patient
3      Latest Vitals by Location
4      Cumulative Vitals Report
5      Print Vitals Entered in Error for a Patient
```

#### Menu Access:

The Vitals/Measurements Results Reporting option is accessed through the Vitals/Measurement menu.

## **GMRV SF511**

### **V/M Graphic Reports**

#### **Description:**

This option provides the user with 5 different reports, the first is a report simulating the Standard Form 511 - Vital Signs Record (Vital Flow Sheet), second is a B/P Plotting Chart, third is a Weight Chart, fourth is a Pulse Oximetry/Respiratory Graph, and last is a Pain Chart. Patient information contained in this report is stored in the GMRV Vital Measurement (#120.5) file, and the GMRY Patient I/O (#126) file.

#### **Additional Information:**

The information printed in this report was entered through the Vitals/Measurements and the GMRY Intake and Output applications. All 24 hour intake and output totals are in cubic centimeters (i.e., milliliters). Reports can be printed by selected patient, for all patients on a ward, and for patients in a specific hospital room/bed. Refer to Printer Issues in the Package Management section for additional information on setting up Kyocera F-800A and Hewlett-Packard LaserJet printers for printing linear graphic reports. Users may print the graphic reports using a dot matrix printer, data plotted is not connected by lines.

#### **Menu Display:**

Select Vitals/Measurement Option: **2** Vitals/Measurements Results Reporting

- 1 V/M Graphic Reports
- 2 Latest Vitals Display for a Patient
- 3 Latest Vitals by Location
- 4 Cumulative Vitals Report
- 5 Print Vitals Entered in Error for a Patient

#### **Screen Prints:**

Select Vitals/Measurements Results Reporting Option: **1** V/M Graphic Reports

- 1 Vital Signs Record
- 2 B/P Plotting Chart
- 3 Weight Chart
- 4 Pulse Oximetry/Respiratory Graph
- 5 Pain Chart

Select a number between 1 and 5: **1** Vital Signs Record// **<RET>**

Vitals by (A)ll patients on a ward, (S)electd Rooms on ward, or (P)atient? **P**

Enter the patient's name.

Select PATIENT NAME: **STARK, ANTHONY L.** 01-19-25 589632111 NSC  
VETERAN

Enter appropriate start date.

Start DATE (TIME optional): T-7// **4/1/97** (APR 01, 1997)

Enter appropriate end date.

Go to DATE (TIME optional): NOW// **4/22/97** (APR 22, 1997)

DEVICE: Enter appropriate response

This report must be queued to a line printer or sent to a slave printer configured to print a 132 column report.

## Vitals/Measurements Reports

The following is an example of an SF511 (Vital Signs Record) graph sent to a dot matrix printer:

Date/Time	04/01/97 08:00	04/02/97 10:19	04/02/97 10:19	04/02/97 11:03	04/02/97 11:49	04/03/97 07:06	04/03/97 10:54	04/03/97 10:58	04/17/97 14:56	04/22/97 10:53
Pulse Temp/F/C										
160 105/40.6-					T*					
150 104/40.0-										
140 103/39.4-										
130 102/38.9-				T*				T*		
120 101/38.3-					P*					
110 100/37.8-	T								T	
100 99/37.2- 98.6/37.0-						T			P	
90 98/36.7-							T			P
80 97/36.1-		P		P		P		P		P
70 96/35.6-							P			
60 95/35.0-										
50 94/34.4-			P**							
Temperature		100.0 A		102.0*T	105.0*R	99.0 T	98.0 T	102.0*T	100.5 T	
Pulse		77	45*	76	120*	80 Lt	72	82 Lt	99	88
Respiration		A Aus	Rad Pal	Rad	A Pal	Rad Aus Ly	Rad	Rad Aus Ly	Rad	Lt Rad Si
Pulse Ox.	0*S	22 S Si		30*S	45*AV Ly	20 S Ly	14 S	12 S Ly	33*S	
L/Min	80				70*			92		
%	2.5				2			4		
Method	30				NC			AM		
BLOOD		99/88	120/80	144/90	200/100	120/88	130/90	132/92		167/97
PRESSURE		RA	Si Dop		LA	LA		LA		LA SmA
Weight (lb)		Si Cu Ad			Ly N-I LgA	Ly Cu Ad		Ly Cu Ad		
(kg)					166.00 D C		140.00 A	152.00 A		
Body Mass Index					75.45		63.64	69.09		
Height (in)					23		19	20		
(cm)	73.00 A				72.00		73.00 E	73.00		
C/G (in)	185.42				182.88		185.42	185.42		
(cm)					12.00 Ank			52.00 Abd		
CVP (cm H2O)					30.48			132.08		
(mm Hg)					22.0			29.9*		
Intake(24 Hr)(cc)					16.2			22.0*		
Output(24 Hr)(cc)										
Pain			5		4			0		
T: Temperature P: Pulse C/G: Circumference/Girth * - Abnormal value ** - Abnormal value off of graph										
TEMP - A: Axillary R: Rectal T: Tympanic PULSE - A: Apical Aus: Auscultate Lt: Left Ly: Lying Pal: Palpated Rad: Radial Si: Sitting RESP - AV: Assisted ventilator Ly: Lying S: Spontaneous Si: Sitting BP - Ad: Adult Cu: Cuff Dop: Doppler LA: L arm LgA: Lg adult Ly: Lying N-I: Non-invasive RA: R arm Si: Sitting SmA: Sm adult HT - A: Actual E: Estimated WT - A: Actual C: Chair D: Dry PULSE OX - AM: Aerosol/humidified mask NC: Nasal cannula CIRC/GIRTH - Abd: Abdominal Ank: Ankle STARK, ANTHONY L. 112-33-4443 NOV 4, 1985 11 YRS MALE Unit: ICU Room: 4M1-A Division: HINES APR 1, 1997 - APR 22, 1997@24:00										
					Page 1	MEDICAL RECORD VITAL FLOW SHEET VAF 10-7987 VICE SF 511				

The B/P Plotting Chart, Weight Chart, and Pulse Oximetry/Respiratory Graph have basically the same appearance with different data. When printing these graphic reports to a Kyocera F-800A or a Hewlett-Packard LaserJet printer, the data on the graph will be connected with lines. The Pain Chart can only be printed on a Hewlett-Packard LaserJet printer.

## Menu Access:

The V/M Graphic Reports option is accessed through the Vitals/Measurements Results Reporting option of the Vitals/Measurement menu.

**GMRV DISPLAY V/M****Latest Vitals Display for a Patient****Description:**

This option prints the latest vitals/measurements for an individual patient. Patient data is retrieved from the GMRV Vital Measurement (#120.5) file.

**Additional Information:**

Data displayed in this report is entered through the Vitals/Measurement Data Entry option, or the Edit a Vital/Measurement Entered in Error option. This report includes vitals/measurements obtained as an outpatient, inpatient, or both.

**Menu Display:**

Select Vitals/Measurement Option: **2** Vitals/Measurements Results Reporting

- 1 V/M Graphic Reports
- 2 Latest Vitals Display for a Patient
- 3 Latest Vitals by Location
- 4 Cumulative Vitals Report
- 5 Print Vitals Entered in Error for a Patient

**Screen Prints:**

Select Vitals/Measurements Results Reporting Option: **2** Latest Vitals Display for a Patient

Select PATIENT NAME: **STARK, ANTHONY L.** 01-19-25 589632111 NSC  
VETERAN

DEVICE: Enter appropriate response

MAR 25, 1997 (11:21) LATEST VITALS REPORT PAGE 1

```

4M1-A      STARK, ANTHONY L.      112-33-4443

Temp.:      (03/21/97@09:48) 99.4 F (37.4 C)(ORAL)
Pulse:      (03/21/97@09:48) 55*(RADIAL)
Resp.:      (03/21/97@09:48) 20
Pulse Ox:   (03/21/97@09:48) 99% with supplemental O2 30% 3L/min
              - VENTURI MASK
B/P:        (03/13/97@09:28) 120/80
Ht.:        (03/20/97@10:39) 2 ft 10 in (86.36 cm)
Wt.:        (03/20/97@10:39) 500 lb (227.27 kg)(ACTUAL, BED)
BMI:        305*
CVP:        (03/21/97@09:48) 24 cmH2O (17.6 mmHg)
Circ/Girth: (01/23/97@08) 32.1 in (81.53 cm)(ABDOMINAL)
Pain:       (03/22/97@09:22) 5
  
```

## Vitals/Measurements Reports

Press return to continue <RET>

### **Menu Access:**

The Latest Vitals Display for a Patient option is accessed through the Vitals/Measurements Results Reporting option of the Vitals/Measurement menu.



## GMRV V/M BY LOCATION

### Latest Vitals by Location

#### Description:

This option prints the latest vitals/measurements for all patients on a specific ward. Vitals data for a patient is stored in the GMRV Vital Measurement (#120.5) file.

#### Additional Information:

Data printed in this report is entered through the Vitals/Measurement Data Entry option, or the Edit a Vital/Measurement Entered in Error option.

#### Menu Display:

Select Vitals/Measurement Option: **2** Vitals/Measurements Results Reporting

- 1 V/M Graphic Reports
- 2 Latest Vitals Display for a Patient
- 3 Latest Vitals by Location
- 4 Cumulative Vitals Report
- 5 Print Vitals Entered in Error for a Patient

#### Screen Prints:

Select Vitals/Measurements Results Reporting Option: **3** Latest Vitals by Location

Select WARD LOCATION NAME: **MICU**

Enter appropriate ward name.

DEVICE: Enter appropriate response

MAR 25,1997 (11:22) VITALS REPORT FOR UNIT: MICU - HINES, IL PAGE 1

-----

	TRAMG, HENRY	793-13-3012
Temp.:	(03/20/97@1)	99.6 F (37.6 C) (ORAL)
Pulse:	(03/20/97@1)	88 (RADIAL)
Resp.:	(03/20/97@1)	24
B/P:	(07/17/96@14:58)	120/80

Press return to continue or "^" to quit **<RET>**

## Vitals/Measurements Reports

MAR 25,1997 (11:22) VITALS REPORT FOR UNIT: MICU - HINES, IL

PAGE 2

500-1 PARKER,PETER P.

Temp.: (03/20/97@1) 100 F (37.8 C)(ORAL)  
Pulse: (03/20/97@1) 55\*(RADIAL)  
Resp.: (03/20/97@1) 20  
Pulse Ox: (02/20/97@07) 89% with supplemental O2 40% 3.5L/min  
- AEROSOL/HUMIDIFIED MASK  
B/P: (02/21/97@14:17) 120/70(L ARM, LYING, CUFF)  
Ht.: (02/20/97@08) 6 ft (182.88 cm)(ESTIMATED)  
Wt.: (02/20/97@08) 200 lb (90.91 kg)  
BMI: 27  
CVP: (02/21/97@14:28) 19.04 cmH2O (14.0 mmHg)  
Circ/Girth: (02/21/97@14:28) 32.25 in (81.92 cm)(THIGH)  
| Pain: (03/22/97@09:22) 5

Press return to continue <RET>

### Menu Access:

The Latest Vitals by Location option is accessed through the Vitals/Measurements Results Reporting option of the Vitals/Measurement menu.

## **GMRV CUMULATIVE V/M**

### **Cumulative Vitals Report**

#### **Description:**

This option prints vitals/measurement information for patients over a user specified period of time. Data comes from the GMRV Vital Measurement (#120.5) file.

#### **Additional Information:**

Data printed in the reports is entered into the database through the Vitals/Measurement Data Entry option, or the Edit a Vital/Measurement Entered in Error option. Reports can be printed by individual patient, for all patients on a ward, or for patients in a specific hospital room/bed.

#### **Restrictions:**

Data prints in chronological order.

#### **Menu Display:**

Select Vitals/Measurement Option: **2** Vitals/Measurements Results Reporting

- 1 V/M Graphic Reports
- 2 Latest Vitals Display for a Patient
- 3 Latest Vitals by Location
- 4 Cumulative Vitals Report
- 5 Print Vitals Entered in Error for a Patient

#### **Screen Prints:**

Select Vitals/Measurements Results Reporting Option: **4** Cumulative Vitals Report

#### **Print by patient.**

Vitals by (A)ll patients on a ward, (S)elected Rooms on ward, or (P)atient? **P**

Select PATIENT NAME: **STARK, ANTHONY L.** 01-19-25 589632111 NSC  
VETERAN

#### **Enter patient name.**

Start with DATE (TIME optional): JAN 3,1993@10:00// **1/1/95** (JAN 01, 1995)

## Vitals/Measurements Reports

Enter an appropriate beginning date for this report. The default date is the date of admission for inpatients or T-7 for outpatients.

Go to DATE (TIME optional): NOW// <RET> (FEB 25, 1997@09:46)

Enter appropriate ending date for this report.

DEVICE: Enter appropriate response

FEB 25, 1997 (10:52)

Cumulative Vitals/Measurements Report

Page 1

-----  
07/05/95

14:00

T: 97.7 F (36.5 C) (ORAL)  
P: 89 (RADIAL)  
R: 20  
B/P: 145/90 (R ARM, SITTING)  
Ht: 72.00 in (182.88 cm)  
Wt: 225.00 lb (102.27 kg)  
Body Mass Index: 31\*

06/18/96

08:00

P: 62 (RADIAL)  
R: 20  
B/P: 122/68 (R ARM, SITTING)

08:01

B/P: 102/60\* (L ARM, SITTING)

08:02

B/P: 100/60\* (R ARM, STANDING)

08:03

B/P: 120/70 (L ARM, STANDING)

08:04

B/P: 110/60\* (R ARM, LYING)

08:05

B/P: 110/70 (L ARM, LYING)

02/17/97

08:00

T: 99.0 F (37.2 C) (ORAL)  
P: 75 (LEFT, RADIAL, PALPATED, SITTING)  
R: 14 (SITTING)  
B/P: 120/70 (L ARM)  
Ht: 72.00 in (182.88 cm)  
Wt: 200.00 lb (90.91 kg) (ACTUAL, STANDING)  
Body Mass Index: 27  
Pulse Ox.: 92 % with supplemental O2 30% 3.5L/min  
via MASK

02/19/97

14:00

T: 98.7 F (37.1 C) (ORAL)  
P: 70 (LEFT, RADIAL)  
R: 20 (STANDING)  
B/P: 120/60\* (L ARM, STANDING)  
Ht: 72.00 in (182.88 cm)  
Wt: 205.00 lb (93.18 kg) (ACTUAL, STANDING)  
Body Mass Index: 28\*  
Pulse Ox.: 98 % with supplemental O2 25% 4L/min  
via MASK

Pain: 5

The following legend is printed on the bottom of the report.

\*\*\* (E) - Error entry

STARK,ANTHONY L 112-33-4443 NOV 4,1985 11 YRS MALE VAF 10-7987j  
Unit: 4E Room: 500-3  
Division: Hines

### **Menu Access:**

The Cumulative Vitals Report option is accessed through the Vitals/Measurements Results Reporting option of the Vitals/Measurement menu.

## **GMRV ERROR REPORT**

### **Print Vitals Entered in Error for a Patient**

#### **Description:**

This option prints a report of all vitals/measurements entered in error for a particular patient for a specific date range. Data comes from the GMRV Vital Measurement (#120.5) file.

#### **Additional Information:**

Patient information printed in this report was entered into the database through the Edit a Vital/Measurement Entered in Error option.

#### **Menu Display:**

Select Vitals/Measurement Option: **2** Vitals/Measurements Results Reporting

- 1 V/M Graphic Reports
- 2 Latest Vitals Display for a Patient
- 3 Latest Vitals by Location
- 4 Cumulative Vitals Report
- 5 Print Vitals Entered in Error for a Patient

#### **Screen Prints:**

Select Vitals/Measurements Results Reporting Option: **5** Print Vitals Entered in Error for a Patient

Select PATIENT NAME: **STARK,ANTHONY L.** 01-19-25 589632111

Enter patient name.

Start with DATE (TIME optional): T-7//**1/1/97** (JAN 01, 1997)

Enter appropriate start date.

go to DATE (TIME optional): NOW// **<RET>** (FEB 25, 1997@12:37)

Enter appropriate end date.

DEVICE: HOME// Enter appropriate response

FEB 25,1997 (11:00) ENTERED IN ERROR VITAL/MEASUREMENT REPORT PAGE 1  
 Patient: STARK,ANTHONY L. 112-33-4443

Date Vit./Meas. taken	User who made error
-----	

FEB 25, 1997@06:30 TEMPERATURE Reason: incorrect reading (Revised) 102.1 F (38.9 C)* (ORAL) (Bad data) 99 F (37.2 C) (ORAL)	WYMAN, SANDRA
--	---------------

FEB 25, 1997@06:30 TEMPERATURE Reason: incorrect reading (Revised) 102.1 F (38.9 C)* (ORAL) (Bad data) 102 F (38.9 C)* (ORAL)	
--	--

FEB 25, 1997@06:30 PAIN Reason: incorrect reading (Revised) 3 (Bad data) 5	WYMAN, SANDRA
---	---------------

Data might not appear in the Revised Data column under certain conditions, for example, when original data was associated with another patient.

### Menu Access:

The Print Vitals Entered in Error for a Patient option is accessed through the Vitals/Measurements Results Reporting option of the Vitals/Measurement menu.





# Index

- Accessing menus, 1.2
- Assigning menus, 1.2
- B/P (Blood Pressure):
  - Diastolic, 2.4-2.5
  - Options containing, 2.3, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20
- Systolic, 2.4-2.5
- BMI (Body Mass Index):
  - Definition, 1
  - Options containing, 4.1
  - Reports containing, 5.2, 5.13, 5.15, 5.17, 5.20
- Central Venous Pressure (example), 4.12
- Change date/time taken (example), 4.14
- Change Default Qualifiers for Temp./Pulse, 2.6
- Circumference/Girth:
  - Options containing, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20
- Circumference measurement (example), 4.9
- CPRS (Computerized Patient Record System), 2
- Create Vital Measurement Quick Order Protocol, 2.18
- Cumulative Vitals Report, 5.17
- CVP (Central Venous Pressure):
  - Options containing, 2.3, 2.18, 4.1, 4.17
  - Reports containing, 5.2, 5.13, 5.15, 5.17, 5.20
- Default qualifiers for temp./pulse, 2.2, 2.6
- Default stop date, 2.3
- Display Administration Schedule File, 2.14
- Display Vitals Category/Qualifier/Synonym Table, 2.11
- Edit a Vitals/Measurement Entered in Error, 4.18
- Edit Administration Schedules File, 2.15
- Edit Vitals Site Parameter File, 2.3
- Editing site configurable files, 1.2
- Enter/Edit Vitals Qualifiers, 2.8
- Enter/Edit Vitals/Measurements, 4.1
- Entering detailed B/P and associated pulse (example), 4.7
- Entering vital signs by selecting specific beds (example), 4.6
- Entering vitals for all of the patients on a particular unit (example), 4.5
- Entering vitals for an individual patient (example), 4.3
- Files:
  - Administration Schedule, 1.1, 1.2, 2.2, 2.14, 2.15
  - GMRV Vital Category, 1.2, 2.2, 2.6, 2.11
  - GMRV Vital Measurement, 2.1, 4.1, 4.17, 5.1, 5.2, 5.13, 5.15, 5.17, 5.20
  - GMRV Vital Qualifier, 1.2, 2.2, 2.8, 2.11
  - GMRV Vital Type, 2.2
  - GMRV Vitals Parameters, 1.2, 2.3, 4.2, 4.4
  - Protocol, 1.2, 2.2, 2.18
- Functionality, 1

- Girth measurement (example), 4.8
- GMRVMGR, 1.3, 2.1, 4.3, 4.17, 5.1
- Health Summary, 1, 2, 1.1
- Height:
  - Options containing, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20
- Implementation and Maintenance, 1.1
- Implementation Considerations, 1.5
- Incorrect date/time, 4.18
- Incorrect patient, 4.18
- Incorrect reading, 4.18
- Interfaces with other packages, 2
- Invalid record, 4.19
- IRMS (Information Resource Management Service), 1.1, 1.4
- Latest Vitals by Location, 5.15
- Latest Vitals Display for a Patient, 5.13
- Main menu for Vitals/Measurements (GMRVMGR), 1.3, 2.1, 4.3, 4.17, 5.1
- Maintenance of Site Files, 2.1
- Name spacing and file listing, 1.1
- Non-Virgin Installation of Software, 1.4
- OE/RR (Order Entry/Results Reporting), 1.1, 1.2, 2.2-2.4, 2.15, 2.18
- Package Operation, 3.1
- PCE (Patient Care Encounter), 2
- Print Vitals Entered in Error for a Patient, 5.20
- Printer issues, 1.3
- Printers:
  - HP LaserJet III, 1.3-1.4, 5.2-5.11
  - HP LaserJet 4, 1.3-1.4, 5.2-5.11
  - HP LaserJet 5, 1.3-1.4, 5.2-5.11
  - Kyocera, 1.3-1.4, 5.2, 5.2-5.11
- Pulse:
  - Options containing, 2.3, 2.6, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20
- Pulse Oximetry:
  - Options containing, 2.3, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20
- Pulse Oximetry (example), 4.10
- Queueing TaskMan jobs, 1.2
- Quick Order, 2, 2.18-2.21
- Resource Requirements, 1.5
- Respiration:
  - Options containing, 2.3, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20
- Scales (Metric and US Customary), 4.2
- Setting up the software environment, 1.1
- SF 511 (Vital Signs Record), 5.4
- Synonym, 2, 2.2, 2.8, 2.9, 2.11, 4.1
- Temperature:
  - Options containing, 2.3, 2.6, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20

- User configurable combination (example), 4.12
- V/M Graphic Reports:, 5.2
  - B/P Plotting Chart, 5.4
  - Pulse Oximetry/Respiration Graph, 5.4
  - Vital Signs Record (SF511), 5.4
  - Weight Chart, 5.4
- Virgin Installation of Software, 1.1
- VISTA, 1, 1.1, 2.14
- Vitals/Measurement, 2.1
- Vitals/Measurement Data Entry, 4.1
- Vitals/Measurements Reports, 5.1
- Vitals/Measurements Results Reporting, 5.1
- Vitals/Measurements Site Files Menu, 2.2
- Weight:
  - Options containing, 2.8, 2.18, 4.1, 4.17
  - Reports containing, 2.11, 5.2, 5.13, 5.15, 5.17, 5.20

