

VA Research Core Waiver for Samples Run by the VA Core Personnel:

All samples run by the VA Research Core will be run to the best of our ability in a timely manner.

We strongly suggest some steps researchers should perform before utilizing some of the instruments here at the VA:

For using the Nanostring Sprint and the NextSeq550

-We strongly suggest running your samples on the Agilent TapeStation to verify the quality and quantity of my sample before running on these instruments.

I have read the above statement and have decided to quantify or check the quality of my sample

I have read the above statement and have decided **NOT** to quantify or check the quality of my sample.

For Flow Core LSR2 and Cytoflex

-We suggest doing a trypan stain before performing your staining or include a live dead stain in your panel. All samples must be fixed before running on any cytometer at the VA.

The VA Research core is not responsible for replacing kits that are faulty for the Nanostring or the NextSeq550.

---

Signature of researcher

---

Date