



Update_2_0_240A

Clinical Reminders

VA-PAVE TEMPLATE UPDATE

Install Guide

March 2024

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION 1
PRE-INSTALL..... 3
INSTALL DETAILS 4
INSTALL EXAMPLE 4
POST INSTALLATION 7

Introduction

The updates to this template are to improve usability based on feedback from the field and meet EPRP (Dmg 61) VVC performance measures criteria.

Changes to PAVE Foot Check

PAVE Foot Check v.1.5

++ Instructions ++

++ Patients with HIGH-RISK level are required to be offered a referral to podiatry. ++

Foot Inspection:

- Document complete NORMAL face-to-face foot check:
- Document complete face-to-face foot check with one or more abnormal findings:
- Document NORMAL VA Video Connect (VVC) foot check:
- Document complete VA Video Connect (VVC) foot check with one or more abnormal findings:
- Patient indicates foot exam (including monofilament test for sensation) was performed in the past year in the private sector:
- Patient refused limb care exam:
- Not applicable for this patient due to:
- Not able to perform complete foot exam today due to:

Veteran Specific Education Material:

- [Field Manual for Foot Health](#)
- [Keeping Feet Healthy](#)
- [Nutrition and Food Diabetes Education](#)

PAVE Foot Check v.1.5

++ Instructions ++

++ Patients with HIGH-RISK level are required to be offered a referral to podiatry. ++

Foot Inspection:

- Document complete NORMAL face-to-face foot check:
- Document complete face-to-face foot check with one or more abnormal findings:
- Document NORMAL VA Video Connect (VVC) foot check:
- Document complete VA Video Connect (VVC) foot check with one or more abnormal findings:
- Patient indicates foot exam (including monofilament test for sensation) was performed in the past year in the private sector:
- Patient refused limb care exam:
- Not applicable for this patient due to:
- Not able to perform complete foot exam today due to:

Veteran Specific Education Material:

- [Field Manual for Foot Health](#)
- [Keeping Feet Healthy](#)
- [Nutrition and Food Diabetes Education](#)

There are no changes to the reminder definition or logic.

UPDATE_2_0_240 contains 1 Reminder Exchange entry:

UPDATE_2_0_240 VA-PAVE TEMPLATE UPDATE

The exchange file contains the following components:

REMINDER GENERAL FINDINGS

VIEW PROGRESS NOTE TEXT

TIU TEMPLATE FIELD

VA-DIABETES NUTRITION AND FOOD
VA-DIABETES:KEEPING FEET HEALTHY
VA-FIELD MANUAL FOR FOOT HEALTH
VA-WORD
OTHER NON-VA
OTHER TEXT
VA-PAVE RESOURCES
VA-PAVE LOW RISK DESCRIPTION
VA-PAVE HIGH RISK DESCRIPTION

EDUCATION TOPICS

VA-DIABETIC FOOT CARE COUNSELING

EXAM

VA-DIABETIC FOOT CHECK

HEALTH FACTORS

VA-REMINDER UPDATES [C]
VA-UPDATE_2_0_240A
VA-PAVE PACT FOOT CHECK [C]
VA-PAVE PACT VVC SENSORY REVIEWED
VA-PAVE PACT VVC PERIPHERAL ARTERY DISEASE REVIEWED
VA-PAVE PODIATRY FOOT EXAM [C]
VA-PAVE FOOT LEVEL 1 EDUCATION
VA-PAVE PACT FOOT LEVEL LOW RISK
VA-PAVE PACT VISUAL FOOT CHECK NORMAL
VA-PAVE PACT VVC FOOT CHECK
VA-PAVE FOOT LEVEL 3 EDUCATION
VA-PAVE PACT PATIENT ESTABLISHED WITH POD/VASC
VA-PAVE PACT PT SEEN BY OUTSIDE FOOT SPECIALIST
VA-PAVE PACT VVC PT HAS ULCER
VA-PAVE PACT VVC PT HAS NO ULCER
VA-PAVE PACT FOOT LEVEL HIGH RISK
VA-PAVE PACT VISUAL FOOT CHECK ABNORMAL
VA-PAVE PACT VVC UNABLE TO PERFORM VISUAL EXAM
VA-PAVE PACT VVC PT DECLINES VISUAL FOOT CHECK
VA-PAVE PACT SENSORY FOOT CHECK NORMAL
VA-PAVE PACT PEDAL PULSES FOOT CHECK NORMAL

VA-PAVE RECEIVES BEHAVIORAL HEALTH AND REFUSE PRIMARY CARE
VA-PAVE BILATERAL FOOT AMPUTATION
VA-PAVE INCORRECT DIAGNOSIS [C]
VA-PAVE INCORRECT NEUROPATHY DIAGNOSIS
VA-PAVE INCORRECT VASCULAR DIAGNOSIS
VA-PAVE INCORRECT RENAL DIAGNOSIS
VA-PAVE INCORRECT DIABETES DIAGNOSIS
REMINDER FACTORS [C]
LIFE EXPECTANCY <6 MONTHS
VA-PAVE PACT PATIENT REFUSED FOOT CHECK
VA-PAVE PACT FOOT CHECK DONE ELSEWHERE
VA-PAVE PACT SENSORY FOOT CHECK ABNORMAL
VA-PAVE PACT PEDAL PULSES FOOT CHECK ABNORMAL

REMINDER TERM

VA-REMINDER UPDATE_2_0_240

REMINDER DIALOG

VA-PACT PAVE FOOT CHECK

Pre-Install

1. **Review the dialog text of the following element and identify what local DATA OBJECT is used for PATIENT ADDRESS, you will need to add it back post install.**
 - a. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 - b. Inquire/Print information for the following element: VAL-TH VVC TELEPHONE TRIAGE PATIENT LOCATION HOME

2. **Review the dialog text of the following element and identify what local DATA OBJECT is used for PATIENT PHONE NUMBER, you will need to add it back post install.**
 - a. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 - b. Inquire/Print information for the following element: VAL-TH VVC TELEPHONE TRIAGE CONTACT PHONE OBJECT

3. **Review the following elements and if you previously added a quick order as an additional finding, you will need to add them back post install:**
 - a. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to group
 - b. Inquire/Print information for the following group:
 - i. VAL-PACT PAVE FOOT REFERRAL TO PODIATRY LEVEL 1 REVISED

- c. Facilities also had the option with the initial install to add additional orders or enable the below elements.
- d. Change your view to element view and Inquire/Print information for the following elements:
 - i. VAL-PACT PAVE PODIATRY CONSULT 2
 - ii. VAL-PACT PAVE PODIATRY CONSULT 3
 - iii. VAL-PACT PAVE PODIATRY CONSULT 4

You can also create an exchange file to back up all of the items listed and install that exchange file backup after this install if you would like.

Install Details

This update is being distributed as a web host file. The address for the host file is:
https://REDACTED/UPDATE_2_0_240A.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 5 minutes.

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: https://REDACTED/UPDATE_2_0_240A.PRD

```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

https://REDACTED/UPDATE_2_0_240A.PRD

You should see a message at the top of your screen that the file successfully loaded.

Search for and locate the reminder exchange entry titled UPDATE_2_0_240A VA-PAVE TEMPLATE UPDATE.

```

+Item  Entry                               Source                               Date Packed
-----
353    UPDATE_2_0_240A VA-PAVE
      TEMPLATE UPDATE
354    UPDATE_2_0_241 VA-HCOE
      TEMPLATES
355    UPDATE_2_0_242 VA-CSP CEAT
      UPDATE
356    UPDATE_2_0_242A VA-CSP CEAT
      UPDATE
357    UPDATE_2_0_243 VA-TBI
      SCREENING UPDATE

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE   Create Exchange File Entry        LHF   Load Host File
CHF   Create Host File                   LMM   Load MailMan Message
CMM   Create MailMan Message             LR    List Reminder Definitions
DFE   Delete Exchange File Entry        LWH   Load Web Host File
IFE   Install Exchange File Entry       RI    Reminder Definition Inquiry
IH    Installation History               RP    Repack
Select Action: Next Screen// IFE   Install Exchange File Entry
Enter a list or range of numbers (1-840): 353

```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
 Enter the number that corresponds with your entry titled **UPDATE_2_0_240A VA-PAVE REMINDER AND TEMPLATE** (in this example it is entry 353, it will vary by site) The date on the exchange file should be 02/21/2024.

```

Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-PACT PAVE FOOT CHECK

REMINDER TERM
  VA-REMINDER UPDATE_2_0_240A

Keywords:

Components:

+      Enter ?? for more actions  >>>
IA   Install all Components        IS   Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
 Select Action: Next Screen// **IA Install all Components**

During the installation, you will see several prompts:

Install all new components

Overwrite any existing components that you are prompted to overwrite

You will be prompted to install the reminder dialog component:

```
Packed reminder dialog: VA-PACT PAVE FOOT CHECK [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PACT PAVE FOOT CHECK	dialog	X
2	1	VA-PACT PAVE SCREENING MAIN Finding: *NONE*	group	X
3	1.2	VA-PACT PAVE FOOT CHECK VERSION Finding: *NONE*	element	X
4	1.3	VA-PACT PAVE FOOT CARE INSTRUCTIONS BTTN Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
5		PXRM GF VIEW BUTTON	prompt	X
6	1.5	VA-PACT PAVE FOOT CARE INSTRUCTIONS Finding: *NONE*	element	X
7	1.10	VA-PACT PAVE FOOT ASSESSMENT PRIMARY CARE OPTIONS Finding: *NONE*	group	X
8	1.10.2	VA-PACT PAVE FOOT CHECK REVIEW NORMAL	group	
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details		DT Dialog Text	IS Install Selected
DF	Dialog Findings		DU Dialog Usage	QU Quit
DS	Dialog Summary		IA Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog – **VA-PODIATRY PAVE FOOT EXAM**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```
Packed reminder dialog: VA-PACT PAVE FOOT CHECK [NATIONAL DIALOG]
```

VA-PACT PAVE FOOT CHECK (reminder dialog) installed from exchange file.

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PACT PAVE FOOT CHECK	dialog	X
2	1	VA-PACT PAVE SCREENING MAIN Finding: *NONE*	group	X
3	1.2	VA-PACT PAVE FOOT CHECK VERSION Finding: *NONE*	element	X
4	1.3	VA-PACT PAVE FOOT CARE INSTRUCTIONS BTTN Finding: VIEW PROGRESS NOTE TEXT (REMINDER GENERAL FINDING)	element	X
5		PXRM GF VIEW BUTTON	prompt	X
6	1.5	VA-PACT PAVE FOOT CARE INSTRUCTIONS Finding: *NONE*	element	X
7	1.10	VA-PACT PAVE FOOT ASSESSMENT PRIMARY CARE OPTIONS Finding: *NONE*	group	X
8	1.10.2	VA-PACT PAVE FOOT CHECK REVIEW NORMAL	group	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details		DT Dialog Text	IS Install Selected
DF	Dialog Findings		DU Dialog Usage	QU Quit
DS	Dialog Summary		IA Install All	
Select Action: Next Screen// Q				

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```
Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-PACT PAVE FOOT CHECK

REMINDER TERM
  VA-REMINDER UPDATE_2_0_240A

Keywords:

Components:

+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components          IS  Install Selected Component
Select Action: Next Screen// Q
```

You will be returned to this screen. Type Q to quit.
Install Completed

Post Installation

1. **In CPRS, Open reminder PAVE Foot Check (reminder dialog VA-PACT PAVE FOOT CHECK) confirm display is version 2.2.**
2. **Add Order to the reminder dialog – review your pre-install inquiry and add back any order as an additional finding.**
 - a. Add your local quick order for a referral to Podiatry for a PAVE evaluation as an additional finding to the reminder dialog group: VAL-PACT PAVE FOOT REFERRAL TO PODIATRY LEVEL 1 REVISED
 - b. If your site would like to have more than one consult type for PAVE, there are three additional elements under the refer to podiatry section that are installed as disabled. You can enable one or more of these elements and add orders as additional findings if needed, you can also edit the dialog text to provide any additional information if needed. .
 - i. VAL-PACT PAVE PODIATRY CONSULT 2
 - ii. VAL-PACT PAVE PODIATRY CONSULT 3
 - iii. VAL-PACT PAVE PODIATRY CONSULT 4
3. **Add your local DATA OBJECT FOR PATIENT ADDRESS** to the following element: VAL-TH VVC TELEPHONE TRIAGE PATIENT LOCATION HOME
 - a. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 - b. Search for and edit the following element: VAL-TH VVC TELEPHONE TRIAGE PATIENT LOCATION HOME

4. **Add your local DATA OBJECT FOR PATIENT PHONE NUMBER** to the following element: VAL-TH VVC TELEPHONE TRIAGE CONTACT PHONE OBJECT.
 - a. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 - b. Search for and edit the following element: VAL-TH VVC TELEPHONE TRIAGE CONTACT PHONE OBJECT

5. (Optional) **Add a local DATA OBJECT** for Podiatry Clinic Visit to the following element: VAL-PACT PAVE FOOT PODIATRY/VASCULAR LAST SEEN.
 - a. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 - b. Search for and edit the following element: VAL-PACT PAVE FOOT PODIATRY/VASCULAR LAST SEEN
 - i. Instructions for how to create this object are located here:REDACTED