



Update _2_0_265

Clinical Reminders

**VA-COVID-19 TEMPLATE UPDATES
Install Guide**

March 2022

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION 1
PRE-INSTALL 5
INSTALL DETAILS 6
INSTALL EXAMPLE 7
POST INSTALLATION 14

Introduction

DESCRIPTION:

Based on updated CDC guidance, feedback from the field and working with clinical leadership updates have been made to the following COVID-19 templates.

VA-COVID-19 DISCHARGE TEMPLATE

Updated isolation criteria for positive options and guidance for quarantine for negative options

VA-COVID-19 TELEPHONE FOLLOW-UP

Removed monoclonal antibodies

Added branching logic if patient had been screened that day, does not require symptom screen unless the user would like to complete it

updated isolation criteria for positive options and guidance for quarantine for negative options

VA-COVID-19 PROVIDER FOLLOW-UP

Re-organized medical history to be alphabetical

Updated social determinants section

Re-organized so that labs and imaging fall before the assessment

Positive assessment options of mild, moderate and severe disease no longer expand

In discharge home options - updated isolation criteria for positive options and guidance for quarantine for negative options

VA-COVID-19 RETURN TO WORK

Simplified options.

VA-COVID-19 STATUS ADJUSTMENT

Updated isolation/quarantine guidance under clinical status section

VA-SCREENING COVID

If positive screen, added two new fields, one for date symptoms began and one for date of exposure at least one is required

This screen is embedded in the ED triage note and the SCI template

VA-COVID-19 CLC SCREEN

If positive screen, added two new fields, one for date symptoms began and one for date of exposure at least one is required

UPDATE_2_0_265 contains 1 Reminder Exchange entry:

UPDATE_2_0_265 VA-COVID-19 TEMPLATE UPDATES

The exchange file contains the following components:

REMINDER GENERAL FINDINGS

VIEW PROGRESS NOTE TEXT

TIU TEMPLATE FIELD

VA-BLANK

VA-CDC UP TO DATE VACCINES

VA-CDC QUARANTINE ISOLATION
HYPERLINK COVID-19 PREVENT SPREAD
HYPERLINK COVID-19 PATIENT EDU
HYPERLINK COVID-19 CARE AT HOME
HYPERLINK COVID-19 DISP NONHOSP PT
HYPERLINK COVID-19 HOW TO PROTECT SELF
VA-COVID TEST RESULTS PHONE INFORMATION
VA-CRC SOCIAL
VA-WORD
VA-PPE TYPE PROVIDER
TEXT (1-10 CHAR) REQ
VA-COVID 19 CDC SYMPTOMS
VA-COVID SYMPTOM TEXT
VA-COVID IM STATUS
OTHER NON-VA
OTHER TEXT
VA-TIME SPENT
DATE
WORD PROC 70X5
ADDITIONAL INFO (DISP ONLY)
VA-CDC IF YOU ARE SICK

HEALTH FACTORS

VA-REMINDER UPDATES [C]
VA-UPDATE_2_0_265
VA-COVID-19 [C]
VA-COVID-19 DC TEST NEGATIVE
VA-COVID-19 DC SUSPECTED FALSE NEGATIVE
VA-COVID-19 DC TEST POSITIVE
VA-COVID-19 DC TEST PENDING RESULTS
VA-COVID-19 DC SUSPECTED, NO TEST
VA-COVID-19 RESIDENTIAL LIVING
VA-COVID-19 HOUSEHOLD CONTACT
VA-COVID-19 HEALTH CARE WORKER
VA-COVID-19 PLAN ADMISSION
VA-COVID-19 PLAN TRANSFER
VA-COVID-19 PLAN REMAIN HOME
VA-COVID-19 PLAN REFERRAL
VA-COVID-19 OUTSIDE PCR SPEC COLLECTED
VA-COVID-19 NON PCR AG LAB OUTSIDE NEG
VA-COVID-19 UNK OUTSIDE LAB TEST NEG
VA-COVID-19 PCR LAB OUTSIDE NEGATIVE
VA-COVID-19 NON PCR AG LAB OUTSIDE POS
VA-COVID-19 UNK OUTSIDE LAB TEST POS
VA-COVID-19 PCR LAB OUTSIDE POSITIVE
VA-COVID-19 SUSPECTED SEVERE

VA-COVID-19 SUSPECTED MILD
VA-COVID-19 SUSPECTED MODERATE
VA-COVID-19 SUSPECTED
VA-COVID-19 SCREENED LOW SUSPICION
VA-COVID-19 FE COUGH DURING INTERVIEW
VA-COVID-19 FE LOOKS ILL
VA COVID-19 FE NO RESPIRATORY DIFFICULTY
VA-COVID-19 FE CYANOSIS
VA-COVID-19 FE ACCESSORY MUSC RESP
VA-COVID-19 FE RESP DISTRESS
VA-COVID-19 FE ALERT
VA-COVID-19 ROS OTHER
VA-COVID-19 ROS GI
VA-COVID-19 ROS CV
VA-COVID-19 ROS RESP
VA-COVID-19 ROS ENT
VA-COVID-19 ROS GENERAL
VA-COVID-19 MED HX PREGNANCY
VA-COVID-19 MED HX HOMELESSNESS
VA-COVID-19 MED HX OBESITY
VA-COVID-19 MED HX BREASTFEEDING
VA-COVID-19 MED HX OTHER
VA-COVID-19 MED HX CANCER
VA-COVID-19 MED HX HTN
VA-COVID-19 MED HX IMMUNOCOMPROMISED
VA-COVID-19 MED HX LIVER DISEASE
VA-COVID-19 MED HX ESRD
VA-COVID-19 MED HX DIABETES
VA-COVID-19 MED HX HEART DISEASE
VA-COVID-19 MED HX LUNG DISEASE
VA-COVID-19 MED HX AGE
VA-COVID-19 VISIT TYPE PHONE
VA-COVID-19 VISIT TYPE F2F
VA-TRAVEL SCREEN [C]
VA-EXPOSURE TO COVID NO
VA-EXPOSURE TO COVID YES
VA-HOME MONITORING BLOOD PRESSURE
VA-HOME MONITORING PULSE OX
VA-HOME MONITORING TEMP
VA-COVID-19 FOLLOW UP PC
VA-COVID-19 PLAN OTHER FOLLOW UP
VA-COVID-19 PLAN SECURE MESSAGE
VA-COVID-19 PLAN ANNIE AP FOLLOW UP
VA-COVID-19 PLAN CONTINUE CALLS
VA-COVID-19 PLAN PROVIDE PREVENTION INFO
VA-COVID-19 PLAN CONTACT VA NURSE LINE

VA-COVID-19 UNABLE TO CARE FOR SELF HOME
VA-COVID-19 CAN CARE FOR SELF AT HOME
VA-COVID-19 NO SYMPTOMS
VA-SYMPTOM VOMITING
VA-SYMPTOM SORE THROAT
VA-SYMPTOM SHORTNESS OF BREATH
VA-SYMPTOM RUNNY NOSE
VA-SYMPTOM NAUSEA
VA-SYMPTOM MUSCLE PAIN
VA-SYMPTOM LOSS OF TASTE OR SMELL
VA-SYMPTOM HEADACHE
VA-SYMPTOM FEVER
VA-SYMPTOM FATIGUE
VA-SYMPTOM DIARRHEA
VA-SYMPTOM COUGH
VA-SYMPTOM CHILLS
VA-EXPOSURE TO FEVER/RESP ILL 14D
VA-EXPOSURE TO FEVER/RESP ILL 14D NO
VA-COVID-19 UNABLE TO LEAVE MESSAGE
VA-COVID-19 UNABLE TO REACH PATIENT
VA-COVID-19 EXPOSURE DATE
VA-COVID-19 SYMPTOM ONSET DATE
VA-COVID SCREEN POSITIVE
VA-COVID-19 RESOLVED
VA-COVID-19 NO LONGER SUSPECTED

REMINDER SPONSOR

OFFICE OF EMERGENCY MEDICINE
OFFICE OF NURSING SERVICE

REMINDER TERM

VA-REMINDER UPDATE_2_0_265
VA-SEX MALE
VA-BL COVID SCREEN COMPLETE TODAY

HEALTH SUMMARY TYPE

VA-COVID-19 EXPOSURE
COVID19 LAB RESULTS

HEALTH SUMMARY OBJECTS

VA-COVID-19 EXPOSURE (TIU)
COVID19 LAB RESULTS (TIU)

TIU DOCUMENT DEFINITION

VA-COVID-19 EXPOSURE
COVID19 LAB RESULTS

REMINDER DIALOG

VA-COVID-19 DISCHARGE INSTRUCTIONS
VA-COVID-19 PROVIDER FOLLOW-UP
VA-COVID-19 TELEPHONE FOLLOW-UP
VA-COVID-19 FEB 2022 UPDATE

Pre-Install

Using the Vista menu Reminder Dialog Manager, Reminder Dialog, change your view to GROUP and choose inquire/print for the following reminder dialog groups:

GROUPS:

VAL-COVID-19 PENDING LOCAL INFORMATION - this dialog group was sent out disabled, you may have enabled it and added additional local content.

VAL-COVID 19 RESULTS NEGATIVE LOCAL INFORMATION – this dialog group was sent out disabled, you may have enabled it and added additional local content.

VAL-COVID-19 PREVENTION LOCAL INFORMATION - this dialog group was sent out disabled, you may have enabled it and added additional local content.

VAL-COVID-19 POSITIVE LOCAL INFORMATION - this dialog group was sent out disabled, you may have enabled it and added additional local content.

VAL-GP COVID-19 FOLLOW-UP PLAN EDUCATION- this dialog group was sent with standard components, your site may have added additional items.

VAL-GP COVID ASSESSMENT LAB TESTS - this dialog group was sent with standard components, your site may have added additional items.

ELEMENTS

VAL-TH VVC TELEPHONE TRIAGE CONTACT PHONE OBJECT - This element should have your local object for the patient phone number

VAL-TH VVC TELEPHONE TRIAGE PATIENT LOCATION HOME - This element should have your local object for patient address

VAL-COVID-19 FOLLOW-UP PLAN ANNIE - Your site could have disabled this element

VAL-COVID-19 FOLLOW-UP PLAN ED NURSE HOTLINE - This element should contain your local contact information

VAL-COVID PATIENT INFORMATION - This element contains data objects for patient name and patient age, you may have replaced them locally.

VAL-COVID ALLERGY AND MEDICATIONS - This element should contain your local object for medication

VAL-COVID FOCUSED EXAM TEMP - This element has a data object for temperature, you may have replaced it with a local object

VAL-COVID PLAN DISCHARGE HOME FOLLOW UP - This element has information on virtual care, you may have added additional text to this element.

VAL-COVID-19 DIRECTIONS TO PACT - This element contains instructions to add PACT team as an additional signer, you may have added additional text to this element.

VAL-COVID PROVIDER NOTE - This is an element you may have used to add additional local information to the VA-COVID-19 PROVIDER FOLLOW-UP note.

This TIU template field contains information on how patients can access lab results. Review and save for reference post install:

VA-COVID TEST RESULTS PHONE INFORMATION

Save this info for the post install, it will make adding information back post -install easier if your site has made local edits. Not all sites will have made local edits to every one of these dialog groups/elements.

Install Details

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 30 minutes. You will need to allow additional time for the pre and post install steps.

This update is being distributed as a web host file. The address for the host file is:

https://REDACTED /UPDATE_2_0_265.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

If you install in your test account, because your test account does not have the COVID immunizations in the immunization file, you may get prompted to replace the immunization finding or you may get an error.

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LWH  Load Web Host File
DFE  Delete Exchange File Entry       RI   Reminder Definition Inquiry
IFE  Install Exchange File Entry       RP   Repack
IH   Installation History
Select Action: Next Screen// LWH   Load Web Host File
Input the URL for the .prd file: https://[REDACTED]/UPDATE_2_0_265.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

https://REDACTED/UPDATE_2_0_265.PRD

You should see a message at the top of your screen that the file loaded successfully

```
https://[REDACTED]/UPDATE_2_0_265.PRD successfully loaded.
```

Search and locate an entry titled **UPDATE_2_0_265 VA-COVID-19 TEMPLATE UPDATES** in reminder exchange.

```
+Item  Entry                               Source                               Date Packed
383    UPDATE_2_0_265 VA-COVID-19
      TEMPLATE UPDATES
384    UPDATE_2_0_27 VA-CHOICE
      DIALOGS
385    UPDATE_2_0_27 VA-CHOICE
      DIALOGS
386    UPDATE_2_0_28
      VA-MENINGOCOCCAL IMMUNIZATIONS
387    UPDATE_2_0_29 VA-HOMELESSNESS
      FOOD INSECURITY
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE    Create Exchange File Entry       LHF    Load Host File
CHF    Create Host File                  LMM    Load MailMan Message
CMM    Create MailMan Message            LR     List Reminder Definitions
DFE    Delete Exchange File Entry        LWH    Load Web Host File
IFE    Install Exchange File Entry       RI     Reminder Definition Inquiry
IH     Installation History              RP     Repack
Select Action: Next Screen// IFE    Install Exchange File Entry
Enter a list or range of numbers (1-624): 383
```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
Enter the number that corresponds with your entry **UPDATE_2_0_265 VA-COVID-19 TEMPLATE UPDATES** (in this example it is entry 383 it will vary by site). The date of the exchange file should be 03/03/2022.

```

Component                                     Category   Exists
Source:
Date Packed:
Package Version: 2.0P71

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-COVID-19 PROVIDER FOLLOW-UP
  VA-COVID-19 TELEPHONE FOLLOW-UP
  VA-COVID-19 DISCHARGE INSTRUCTIONS
  VA-COVID-19 FEB 2022 UPDATE

REMINDER TERM
  VA-REMINDER UPDATE_2_0_265
Non-exchangeable TIU object(s):
-----
+      Enter ?? for more actions                                     >>>
IA   Install all Components           IS   Install Selected Component
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
 Select Action: Next Screen// **IA Install all Components**
 You will see several prompts, for all new entries you will choose **I to Install**
*For components that already exists but the packed component is different, you will choose **S to Skip** with one exception, if you are prompted to overwrite the reminder term VA-BL COVID SCREEN COMPLETE TODAY you should **overwrite**.*
 REMINDER TERM entry named VA-BL COVID SCREEN COMPLETE TODAY already exists but the packed component is different, what do you want to do?
 Select one of the following:
 C Create a new entry by copying to a new name
 M Merge findings
 O Overwrite the current entry
 U Update
 Q Quit the install
 S Skip, do not install this entry
 Enter response: O// overwrite the current entry
 Are you sure you want to overwrite? N// y YES

You will be prompted to install the FIRST reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

```
Packed reminder dialog: VA-COVID-19 FEB 2022 UPDATE [NATIONAL DIALOG]
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-COVID-19 FEB 2022 UPDATE                             group
      Finding: *NONE*
  2    5 VA-COVID-19 STATUS VERSION                             element  X
      Finding: *NONE*
  3    10 VA-GP COVID-19 STATUS CLINICAL                       group    X
      Finding: *NONE*
  4    10.5 VA-COVID-19 STATUS CLINICAL                        element  X
      Finding: *NONE*
  5    10.10 VA-GP COVID-19 STATUS CLINICAL MAIN              group    X
      Finding: *NONE*
  6    10.10.5 VA-COVID-19 STATUS SUSPECTED                   element  X
      Finding: VA-COVID-19 SUSPECTED (HEALTH FACTOR)
  7    VISIT DATE TODAY                                        forced   X
+      + Next Screen    - Prev Screen    ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-COVID-19 FEB 2022 UPDATE**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**

```
Packed reminder dialog: VA-COVID-19 FEB 2022 UPDATE [NATIONAL DIALOG]
VA-COVID-19 FEB 2022 UPDATE (group) installed from exchange file.
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-COVID-19 FEB 2022 UPDATE                             group  X
      Finding: *NONE*
  2    5 VA-COVID-19 STATUS VERSION                             element  X
      Finding: *NONE*
  3    10 VA-GP COVID-19 STATUS CLINICAL                       group    X
      Finding: *NONE*
  4    10.5 VA-COVID-19 STATUS CLINICAL                        element  X
      Finding: *NONE*
  5    10.10 VA-GP COVID-19 STATUS CLINICAL MAIN              group    X
      Finding: *NONE*
  6    10.10.5 VA-COVID-19 STATUS SUSPECTED                   element  X
      Finding: VA-COVID-19 SUSPECTED (HEALTH FACTOR)
  7    VISIT DATE TODAY                                        forced   X
+      + Next Screen    - Prev Screen    ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

After completing this dialog install, you will type Q.

You will be prompted to install the next reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

```
Packed reminder dialog: VA-COVID-19 DISCHARGE INSTRUCTIONS [NATIONAL DIALOG]
Item  Seq.  Dialog Findings                                     Type  Exists
 1     VA-COVID-19 DISCHARGE INSTRUCTIONS                dialog  X
 2     5 VA-COVID-19 DISCHARGE INSTRUCTIONS VERSION      element  X
      Finding: *NONE*
 3     10 VA-COVID-19 DC INSTRUCTIONS HEADER            group   X
      Finding: *NONE*
 4     10.10 VA-COVID 19 RESULTS NEGATIVE                group   X
      Finding: VA-COVID-19 DC TEST NEGATIVE (HEALTH FACTOR)
 5     10.10.5 VAL-COVID 19 RESULTS NEGATIVE            group   X
      RECOMMENDATIONS
      Finding: *NONE*
 6     10.10.5.5 VA-COVID 19 RESULTS NEGATIVE INFORMATION group   X
      Finding: *NONE*
 7     10.10.5.10 VA-COVID-19 PREVENTION STEPS NEG     group   X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog –**VA-COVID-19 DISCHARGE INSTRUCTIONS**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**

```
Packed reminder dialog: VA-COVID-19 DISCHARGE INSTRUCTIONS [NATIONAL DIALOG]
VA-COVID-19 DISCHARGE INSTRUCTIONS (reminder dialog) installed from exchange fil
Item  Seq.  Dialog Findings                                     Type  Exists
 1     VA-COVID-19 DISCHARGE INSTRUCTIONS                dialog  X
 2     5 VA-COVID-19 DISCHARGE INSTRUCTIONS VERSION      element  X
      Finding: *NONE*
 3     10 VA-COVID-19 DC INSTRUCTIONS HEADER            group   X
      Finding: *NONE*
 4     10.10 VA-COVID 19 RESULTS NEGATIVE                group   X
      Finding: VA-COVID-19 DC TEST NEGATIVE (HEALTH FACTOR)
 5     10.10.5 VAL-COVID 19 RESULTS NEGATIVE            group   X
      RECOMMENDATIONS
      Finding: *NONE*
 6     10.10.5.5 VA-COVID 19 RESULTS NEGATIVE INFORMATION group   X
      Finding: *NONE*
 7     10.10.5.10 VA-COVID-19 PREVENTION STEPS NEG     group   X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

After completing this dialog install, you will type Q.

You will be prompted to install the next reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

```
Packed reminder dialog: VA-COVID-19 TELEPHONE FOLLOW-UP [NATIONAL DIALOG]
Item  Seq.  Dialog Findings                                     Type  Exists
  1    VA-COVID-19 TELEPHONE FOLLOW-UP                    dialog  X
  2    5 VA-COVID-19 FOLLOW-UP CALL VERSION                element  X
      Finding: *NONE*
  3    8 VA-GP COVID-19 FOLLOW-UP CONTACT                 group   X
      Finding: *NONE*
  4    8.5 VA-GP COVID-19 FOLLOW-UP CALL                  group   X
      Finding: *NONE*
  5    8.5.5 VA-GP COVID-19 SYMPTOMS FOLLOW UP*          group   X
      Finding: *NONE*
  6    8.5.5.10 VA-GP COVID UPDATED SCREENING QUESTIONS group   X
      Finding: *NONE*
  7    8.5.5.10.5 VA-COVID SYMPTOM CHILLS                element  X
      Finding: VA-SYMPTOM CHILLS (HEALTH FACTOR)          X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-COVID-19 TELEPHONE FOLLOW-UP**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**

```
Packed reminder dialog: VA-COVID-19 TELEPHONE FOLLOW-UP [NATIONAL DIALOG]
VA-COVID-19 TELEPHONE FOLLOW-UP (reminder dialog) installed from exchange file.
Item  Seq.  Dialog Findings                                     Type  Exists
  1    VA-COVID-19 TELEPHONE FOLLOW-UP                    dialog  X
  2    5 VA-COVID-19 FOLLOW-UP CALL VERSION                element  X
      Finding: *NONE*
  3    8 VA-GP COVID-19 FOLLOW-UP CONTACT                 group   X
      Finding: *NONE*
  4    8.5 VA-GP COVID-19 FOLLOW-UP CALL                  group   X
      Finding: *NONE*
  5    8.5.5 VA-GP COVID-19 SYMPTOMS FOLLOW UP*          group   X
      Finding: *NONE*
  6    8.5.5.10 VA-GP COVID UPDATED SCREENING QUESTIONS group   X
      Finding: *NONE*
  7    8.5.5.10.5 VA-COVID SYMPTOM CHILLS                element  X
      Finding: VA-SYMPTOM CHILLS (HEALTH FACTOR)          X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// Q
```

After completing this dialog install, you will type Q.

You will be prompted to install the LAST reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

```
Packed reminder dialog: VA-COVID-19 PROVIDER FOLLOW-UP [NATIONAL DIALOG]
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-COVID-19 PROVIDER FOLLOW-UP                          dialog  X
  2    5 VA-COVID-19 PROVIDER FOLLOW-UP VERSION                 element  X
      Finding: *NONE*
  3    10 VA-GP COVID 19 SCREEN METHOD                          group   X
      Finding: *NONE*
  4    10.5 VA-GP COVID-19 PROVIDER MODALITY CVT               group   X
      Finding: *NONE*
  5    10.5.3 VA-GP TH VVC MODALITY COVID                      group   X
      Finding: *NONE*
  6    10.5.3.5 VA-GP TH VVC TELEPHONE TRIAGE PT ID           group   X
      Finding: *NONE*
  7    10.5.3.5.10 VA-GP TH VVC TELEPHONE TRIAGE PT ID       group   X
      VERIFIED
+      + Next Screen    - Prev Screen    ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// IA
```

At the Select Action prompt, type **IA** to install the dialog – **VA-COVID-19 PROVIDER FOLLOW-UP**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**

```
Packed reminder dialog: VA-COVID-19 PROVIDER FOLLOW-UP [NATIONAL DIALOG]
VA-COVID-19 PROVIDER FOLLOW-UP (reminder dialog) installed from exchange file.
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-COVID-19 PROVIDER FOLLOW-UP                          dialog  X
  2    5 VA-COVID-19 PROVIDER FOLLOW-UP VERSION                 element  X
      Finding: *NONE*
  3    10 VA-GP COVID 19 SCREEN METHOD                          group   X
      Finding: *NONE*
  4    10.5 VA-GP COVID-19 PROVIDER MODALITY CVT               group   X
      Finding: *NONE*
  5    10.5.3 VA-GP TH VVC MODALITY COVID                      group   X
      Finding: *NONE*
  6    10.5.3.5 VA-GP TH VVC TELEPHONE TRIAGE PT ID           group   X
      Finding: *NONE*
  7    10.5.3.5.10 VA-GP TH VVC TELEPHONE TRIAGE PT ID       group   X
      VERIFIED
+      + Next Screen    - Prev Screen    ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings     DU  Dialog Usage     QU  Quit
DS  Dialog Summary      IA  Install All
Select Action: Next Screen// Q
```

After completing this dialog install, you will type Q.

```

Component                                     Category   Exists
Source:
Date Packed:
Package Version: 2.0P71

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-COVID-19 PROVIDER FOLLOW-UP
  VA-COVID-19 TELEPHONE FOLLOW-UP
  VA-COVID-19 DISCHARGE INSTRUCTIONS
  VA-COVID-19 FEB 2022 UPDATE

REMINDER TERM
  VA-REMINDER UPDATE_2_0_265
Non-exchangeable TIU object(s):
-----
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components          IS  Install Selected Component
Select Action: Next Screen// Q

```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

Post Installation

1. Open the COVID-19 Screen Template (VA-SCREEN COVID) and confirm the version is 9.1
2. Open the COVID-19 Screen Template (VA-COVID-19 CLC SCREEN) and confirm the version is 5.1
3. Open the COVID-19 Provider Screen Template (VA-COVID-19 PROVIDER FOLLOW-UP) and confirm the version is 5.1
4. Open the COVID-19 Telephone Follow-Up (VA-COVID-19 TELEPHONE FOLLOW-UP) and confirm the version is 8.1
5. Open the COVID-19 Discharge Template (VA-COVID-19 DISCHARGE INSTRUCTIONS) and confirm the version is 5.1
6. Open the COVID-19 Status Adjustment Template (VA-COVID-19 STATUS ADJUSTMENTS) and confirm the version is 5.1
7. Open the COVID-19 Return to Work Template (VA-COVID-19 RETURN TO WORK) and confirm the version is 3.1

8. Add local content to reminder dialog groups and elements if needed.

Review the pre-install for these items and edit to make any local changes.

- If you had made changes to the below items before, you will need to re-do any local modifications to the following groups.
- If your pre-install inquiry matches the post install review, no changes are needed.
- Some items will require that you enter your local data objects, those are specified below.

GROUPS:

VAL-COVID-19 PENDING LOCAL INFORMATION - if in the pre-install you had enabled and added local information, add that information back and re-enable this group.

VAL-COVID 19 RESULTS NEGATIVE LOCAL INFORMATION - if in the pre-install you had enabled and added local information, add that information back and re-enable this group.

VAL-COVID-19 PREVENTION LOCAL INFORMATION - if in the pre-install you had enabled and added local information, add that information back and re-enable this group.

VAL-COVID-19 POSITIVE LOCAL INFORMATION - if in the pre-install you had enabled and added local information, add that information back and re-enable this group.

VAL-GP COVID ASSESSMENT LAB TESTS – we have added a new sequence to this group but you can add back any additional local items you previously had, do not remove any sequences

VAL-GP COVID-19 FOLLOW-UP PLAN EDUCATION - this dialog group was sent with standard components; your site may have added additional items. If in the pre-install you have additional items, you may add them back.

ELEMENTS

VAL-TH VVC TELEPHONE TRIAGE CONTACT PHONE OBJECT – replace the dialog text *PHONE OBJECT HERE* with your local object that represents the patient phone number

VAL-TH VVC TELEPHONE TRIAGE PATIENT LOCATION HOME – replace the dialog text *PATIENT ADDRESS OBJECT HERE* with your local object that represents the patient phone number

VAL-COVID-19 FOLLOW-UP PLAN ANNIE - Your site could have disabled this element, compare with your pre-install and make changes if needed.

VAL-COVID-19 FOLLOW-UP PLAN ED NURSE HOTLINE -This element should contain your local contact information compare with your pre-install and make changes if needed.

VAL-COVID PATIENT INFORMATION- This element has data objects for patient name and patient age, you may have replaced them locally, compare with your pre-install and make changes if needed.

VAL-COVID ALLERGY AND MEDICATIONS - replace the dialog text *****local medication object to be added here for display only***** with your local object for current patient medications.

VAL-COVID FOCUSED EXAM TEMP- This element has data objects for patient temperature, you may have replaced it locally, compare with your pre-install and make changes if needed.

VAL-COVID PLAN DISCHARGE HOME FOLLOW UP- replace the dialog text ****SITE SPECIFIC INFO CAN BE ADDED HERE**** with any local information you had added.

VAL-COVID-19 DIRECTIONS TO PACT- This element contains instructions to add PACT team as an additional signer, you may have added additional text to this element, compare with your pre-install and make changes if needed.

VAL-COVID PROVIDER NOTE - This is an element you may have used to add additional local information to the VA-COVID-19 PROVIDER FOLLOW-UP note, compare with your pre-install and make changes if needed.

In addition to the items identified in the pre-install the following items also should be edited:

VAL-COVID-19 PREVENTION INFORMATION NOT TESTED PN – replace the alternate progress note text ****ADD LOCAL INFO FOR TESTING**** with your site local information for how Veterans can get tested for COVID-19.

VAL-COVID RACE ETHNICITY CONFIRM NO - **replace the dialog text *****local info here******* with information on who the Veteran should contact at your site to get the race and ethnicity changed in the medical record if it is incorrectly listed.

This TIU template field contains information on how patients can access lab results. Review the pre install and make any needed changes:

VA-COVID TEST RESULTS PHONE INFORMATION

This TIU template field should contain information about your site, please make sure it is updated appropriately.

VA COVID SITE NAME