



**Update \_2\_0\_295**

**Clinical Reminders**

**VA-DIAGNOSTIC IMAGING ORDER/EXAMINATION  
REVIEW**

**Install Guide**

**July 2022**

Product Development  
Office of Information Technology  
Department of Veterans Affairs

Table of Contents

**INTRODUCTION** ..... 1

**INSTALL DETAILS** ..... 1

**INSTALL EXAMPLE** ..... 2

**POST INSTALLATION** ..... 5

## Introduction

---

### Description:

The template will be used by clinical radiology staff to document a diagnostic imaging order or examination review. This review will include determining if an imaging order/examination will be placed in the hold status, cancelled, or changed to another imaging modality. An example of when this note would be used is the current national shortage of iodinated contrast.

UPDATE\_2\_0\_295 contains 1 Reminder Exchange entry:

**UPDATE\_2\_0\_295 VA-DIAGNOSTIC IMAGING ORDER REVIEW**

**The exchange file contains the following components:**

### TIU TEMPLATE FIELD

TEXT (1-60 CHAR)REQUIRED  
VA-DIAGNOSTIC IMAGING TXT (1-50)R  
VA-BLANK

### HEALTH FACTORS

VA-REMINDER UPDATES [C]  
VA-UPDATE\_2\_0\_295  
VA-DIAGNOSTIC IMAGE ORDER REVIEW [C]  
VA-RAD DX IMAG ORDER REV ASSESS POT HARM NO  
VA-RAD DX IMAG ORDER REV ASSESS POT HARM YES REC INSTIT DISCLOS  
VA-RAD DX IMAG ORDER REV ASSESS POT HARM YES REC CLIN DISCLOS  
VA-RAD DX IMAG ORDER REV REASON OTHER  
VA-RAD DX IMAG ORDER REV REASON GLOBAL CONTRAST SHORTAGE  
VA-RAD DX IMAG ORDER REV REASON COMPANION CT SAFETY LOOKBACK  
VA-RAD DX IMAG ORDER REV REASON SERVICE CURTAILMENT  
VA-RAD DX IMAG ORDER REV ACTION REC ALT IMAG MODALITY  
VA-RAD DX IMAG ORDER REV ACTION OTHER  
VA-RAD DX IMAG ORDER REV ACTION CANCELLED  
VA-RAD DX IMAG ORDER REV ACTION PLACED ON HOLD

### REMINDER SPONSOR

RADIOLOGY

### REMINDER TERM

VA-REMINDER UPDATE\_2\_0\_295

### REMINDER DIALOG

VA-DIAGNOSTIC IMAGING ORDER REVIEW

## Install Details

---

This update is being distributed as a web host file. The address for the host file is:  
[https://REDACTED/UPDATE\\_2\\_0\\_295.PRD](https://REDACTED/UPDATE_2_0_295.PRD)

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 5 minutes.

## Install Example

To Load the exchange file, use LWH. The URL is

[https://REDACTED/UPDATE\\_2\\_0\\_295.PRD](https://REDACTED/UPDATE_2_0_295.PRD)

```
+      + Next Screen  - Prev Screen  ?? More Actions
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                LHM  Load MailMan Message
CMM  Create MailMan Message          LR   List Reminder Definitios
DFE  Delete Exchange File Entry      LWH  Load Web Host File
IFE  Install Exchange File Entry     RI   Reminder Definition Inqy
IH   Installation History            RP   Repack
Select Action: Next Screen// LWH   Load Web Host File
Input the URL for the .prd file: https://[REDACTED]/UPDATE_2_0_295.PRD
```

If you use LWH and it works, you will see something like this message. Make sure the file was loaded successfully.

[https://\[REDACTED\]/UPDATE\\_2\\_0\\_295.PRD](https://[REDACTED]/UPDATE_2_0_295.PRD) successfully loaded.

Search and locate an entry titled UPDATE\_2\_0\_295 VA-DIAGNOSTIC IMAGING ORDER REVIEW in reminder exchange.

```
+Item  Entry                               Source          Date Packed
317    UPDATE_2_0_295 VA-DIAGNOSTIC          G*              07/05/2022
      IMAGING ORDER REVIEW
318    UPDATE_2_0_30 VA-IRAQ &             G*              07/05/2022
      AFGHAN POST-DEPLOY
319    UPDATE_2_0_31 VA-REACH              G*              07/05/2022
320    UPDATE_2_0_31 VA-REACH              G*              07/05/2022
321    UPDATE_2_0_32 VA-MH                 G*              07/05/2022
      CONTINGENCY MANAGEMENT
322    UPDATE_2_0_33 MISC UPDATE 4         G*              07/05/2022
323    UPDATE_2_0_34 VA-WH PAP SHEAR       G*              07/05/2022
+      + Next Screen  - Prev Screen  ?? More Actions
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                LHM  Load MailMan Message
CMM  Create MailMan Message          LR   List Reminder Definitios
DFE  Delete Exchange File Entry      LWH  Load Web Host File
IFE  Install Exchange File Entry     RI   Reminder Definition Inqy
IH   Installation History            RP   Repack
Select Action: Next Screen// IFE   Install Exchange File Entry
Enter a list or range of numbers (1-519): 317
```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry  
Enter the number that corresponds with your entry titled UPDATE\_2\_0\_295 VA-DIAGNOSTIC IMAGING ORDER REVIEW (in this example it is entry 317, it will vary by site) that is dated 07/05/2022.

```

Component                                     Category   Exists
Source:
Date Packed:
Package Version: 2.0P71

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-DIAGNOSTIC IMAGING ORDER REVIEW

REMINDER TERM
  VA-REMINDER UPDATE_2_0_295

Keywords:

Components:

+      Enter ?? for more actions
IA  Install all Components                      IS  Install S
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.  
**For all new entries in any category you will choose I to Install**  
**Skip any items** in any category you are prompted to overwrite

You will then be prompted to install the dialog component.

```

Packed reminder dialog: VA-DIAGNOSTIC IMAGING ORDER REVIEW [NATIONAL DIALOG]

```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-DIAGNOSTIC IMAGING ORDER REVIEW	dialog	
2	5	VA-TEXT DIAGNOSTIC IMAGING ORDER REVIEW VERSION Finding: *NONE*	element	
3	7	VA-BLANK LINE (EXCLUDE FROM PN) *ELE Finding: *NONE*	element	X
4	10	VA-GP DIAGNOSTIC IMAGE REVIEW ORDER Finding: *NONE*	group	
5	10.5	VA-EL DIAGNOSTIC IMAGE REVIEW ORDER SPECIFY Finding: *NONE*	element	
6		PXRM COMMENT	prompt	X
7	12	VA-BLANK LINE (EXCLUDE FROM PN) *ELE	element	X

```

+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details          DT  Dialog Text
DF  Dialog Findings        DU  Dialog Usage
DS  Dialog Summary         IA  Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-DIAGNOSTIC IMAGING ORDER REVIEW**

Install reminder dialog and all components with no further changes: Y//ES

```

Packed reminder dialog: VA-DIAGNOSTIC IMAGING ORDER REVIEW [NATIONAL DIALOG]
VA-DIAGNOSTIC IMAGING ORDER REVIEW (reminder dialog) installed from exchange file
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-DIAGNOSTIC IMAGING ORDER REVIEW                     dialog  X
  2    5 VA-TEXT DIAGNOSTIC IMAGING ORDER REVIEW VERSION      element  X
      Finding: *NONE*
  3    7 VA-BLANK LINE (EXCLUDE FROM PN) *ELE                 element  X
      Finding: *NONE*
  4    10 VA-GP DIAGNOSTIC IMAGE REVIEW ORDER                  group    X
      Finding: *NONE*
  5    10.5 VA-EL DIAGNOSTIC IMAGE REVIEW ORDER SPECIFY      element  X
      Finding: *NONE*
  6    PXR COMMENT                                           prompt   X
  7    12 VA-BLANK LINE (EXCLUDE FROM PN) *ELE                 element  X
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text
DF  Dialog Findings    DU  Dialog Usage
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```

Component                                     Category  Exists
Source:
Date Packed:
Package Version: 2.0P71

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-DIAGNOSTIC IMAGING ORDER REVIEW

REMINDER TERM
  VA-REMINDER UPDATE_2_0_295

Keywords:

Components:

+      + Next Screen  - Prev Screen  ?? More Actions
IA  Install all Components                                     IS  Install S
Select Action: Next Screen// Q

```

You will then be returned to this screen. At the **Select Action** prompt, type **Q**.

Install complete.

## Post Installation

---

### 1. Setup of Note Title.

The note title to create to associate with this template is:

DIAGNOSTIC IMAGING ORDER EXAMINATION REVIEW

VHA Enterprise Standard Title: RADIOLOGY NOTE

**You can add these to your document class where IMAGING progress notes are located.**

*Note: Comprehensive information on Creating new Document Classes and Titles can be found beginning on page 54 of the TIU/ASU Implementation Guide at the following link: [http://www.va.gov/vdl/documents/Clinical/CPRS-Text Integration Utility \(TIU\)/tiuim.pdf](http://www.va.gov/vdl/documents/Clinical/CPRS-Text%20Integration%20Utility%20(TIU)/tiuim.pdf)*

### 2. Make the Dialog templates available to be attached to a Personal or Shared templates or to be able to attach the dialog to a progress note title.

Select Reminder Managers Menu <TEST ACCOUNT> Option: **CP CPRS**

#### Reminder Configuration

CA Add/Edit Reminder Categories  
CL CPRS Lookup Categories  
CS CPRS Cover Sheet Reminder List  
MH Mental Health Dialogs Active  
PN Progress Note Headers  
RA Reminder GUI Resolution Active  
**TIU TIU Template Reminder Dialog Parameter**  
DEVL Evaluate Coversheet List on Dialog Finish  
DL Default Outside Location  
PT Position Reminder Text at Cursor  
LINK Link Reminder Dialog to Template  
TEST CPRS Coversheet Time Test  
NP New Reminder Parameters  
GEC GEC Status Check Active  
WH WH Print Now Active

Select CPRS Reminder Configuration <TEST ACCOUNT> Option: **TIU TIU**

#### Template Reminder Dialog Parameter

Reminder Dialogs allowed as Templates may be set for the following:

1 User        USR [choose from NEW PERSON]  
3 Service     SRV [choose from SERVICE/SECTION]  
4 Division    DIV [SALT LAKE CITY]  
5 System      SYS [NATREM.FO-SLC.MED.VA.GOV]

Enter selection: **5** *Choose the appropriate number for your site. Your site may do this by System or other levels. The example below uses SYSTEM level.*

Setting Reminder Dialogs allowed as Templates for System: NATREM.FO-SLC.MED.VA.GOV

Select Display Sequence: **?**

388 VA-COVID-19 TELEPHONE FOLLOW-UP

389 VA-SRA-CONSULT RESPONSE (D)

*When you type a question mark above, you will see the list of #'s (with dialogs) that are already taken. Choose a number **NOT** on this list. For this example, 390 is not present so I will use 390.*

Select Display Sequence: **390**

Are you adding 390 as a new Display Sequence? Yes// **YES**


Display Sequence: 390// **<Enter>** 390

Clinical Reminder Dialog: **VA-DIAGNOSTIC IMAGING ORDER REVIEW**  
reminder dialog NATIONAL

...OK? Yes// **<Enter>** (Yes)



3. **Associate the reminder dialogs with the note title in CPRS**
- a. In CPRS, open the Template Editor, expand Document Titles, select the appropriate folder (if applicable), then click New Template (top right-hand corner).



Type in the following information:

Name: **DIAGNOSTIC IMAGING ORDER EXAMINATION REVIEW**

Template Type: **Reminder Dialog**

Reminder Dialog: **VA-DIAGNOSTIC IMAGING ORDER REVIEW**

Next, type in the Associated Title name **DIAGNOSTIC IMAGING ORDER EXAMINATION REVIEW**

— then click Apply.

#### 4. UPDATE TIU TEMPLATE FIELD:

Using the Template field editor in CPRS, edit the template field names VA-DIAGNOSTIC IMAGING TXT (1-50)R

Replace the characters “(xxx)xxx-xxxx Ext xxxx” with the phone number for your facility’s Imaging Department number.

The screenshot shows the 'Template Field Editor' window. The left pane lists various template fields, with 'VA-DIAGNOSTIC IMAGING TXT (1-50)R' selected. The right pane shows the details for this field: Name: VA-DIAGNOSTIC IMAGING TXT (1-50)R, Type: Text, and Items: (xxx)xxx-xxxx Ext xxxx. The 'Miscellaneous' section includes checkboxes for Inactive, Separate Lines, Required, and Exclude From Note. The 'Indent' section includes spinners for Indent Field and Indent Text. At the bottom, there are buttons for Preview, OK, Cancel, and Apply.

If your site has an integrated Vista database or you want to provide multiple phone number options, you can change the type to a combo, check box or radio buttons and add multiple lines of phone numbers.