

Update _2_0_300

Clinical Reminders

VA-TH-VVC TEMPLATE UPDATES Install Guide

July 2022

Product Development Office of Information Technology Department of Veterans Affairs

Table of Contents	
INTRODUCTION	1
PRE-INSTALL	5
INSTALL DETAILS	5
INSTALL EXAMPLE	6
POST INSTALLATION	9

Introduction

Description:

Guidance from the telehealth program office was recently released updating the required synchronous telehealth disclosure statement:

VVC Disclosure Statement				
 Old – Visit conducted by synchronous telehealth. Patient verbal consent obtained. Location/emergency number confirmed. 	 New - Visit conducted by synchronous telehealth. Veteran verbal consent obtained. Location/emergency number confirmed. Environment surveyed and all participants identified. Virtual conference room locked 			

Templates that included one of the following statements:

Before:

EXAMPLE 1 ******* Telehealth Disclosure: Visit conducted by clinical video telehealth. Patient verbal consent obtained. Location/emergency number confirmed. The virtual conference room was locked after all participants joined.

EXAMPLE 2

Patient confirms the location is private and safe for visit?

🖸 *Yes

Visit conducted by clinical video telehealth.

*Patient verbal consent for VA Video Connect visit obtained.

*Location/emergency number confirmed.

After:

```
Telehealth Disclosure:
Visit conducted by synchronous telehealth.
Patient verbal consent obtained.
Location/emergency number confirmed.
Environment surveyed and all participants identified.
Virtual conference room locked.
```

These are the dialogs that are impacted by this change:

VA-TH VVC HYPERTENSION VA-COVID-19 PROVIDER FOLLOW-UP VA-CSP VETERAN WELLNESS CONTACT **VA-PGCSS VETERAN INTAKE** VA-SAH GENERAL NOTE VA-SAH GROUP NOTE VA-TH VVC EM CASE MANAGEMENT VA-TH VVC EM CLINICAL PHARMACY VA-TH VVC EM MENTAL HEALTH VA-TH VVC EM PROVIDER VA-CSP CAREGIVERS FIRST NOTE VA-CSP PCAFC APP INTAKE SECONDARY CG VA-CSP PGCSS INTERMITTENT NOTE VA-CSP CAREGIVER ASSESSMENT REV VA-CSP CAREGIVER WELLNESS CONTACT VA-CSP PCAFC INTAKE VA-CSP PCAFC RE-ASSESSMENT VA-CSP PGCSS CAREGIVER ASSESSMENT

In the reminder dialog VA-TELE URGENT CARE NOTE: Before:



🖸 Video Visit: Patient verbal consent obtained. Location/emergency number confirmed.

After:

```
Type of Visit:

Video Visit:

Telehealth Disclosure:

Visit conducted by synchronous telehealth.

Patient verbal consent obtained.

Location/emergency number confirmed.

Environment surveyed and all participants identified.

Virtual conference room locked.
```

UPDATE_2_0_300 contains 1 Reminder Exchange entry: UPDATE_2_0_300 VA-TH-VVC TEMPLATE UPDATES

The exchange file contains the following components: REMINDER GENERAL FINDINGS

VIEW PROGRESS NOTE TEXT

TIU TEMPLATE FIELD

VA-BLANK SP VA-EDIT 12+ VA-BLANK IN PN TEXT VA-EDIT 7+ VA-HTN GUIDANCE VA-BLANK VA-OTHER TEXT VA-WORD PROCESS 65 VA-EDIT 20 VA-PROVIDER NOTIFIED VA-HTN FALL PREVENTION VA-WORD

EDUCATION TOPICS

VA-HTN MEDICATION ADHERENCE VA-HTN NUTRITION EDUCATION VA-HTN EXERCISE VA-BP CUFF USE VA-SAFETY/HOME/FALLS

HEALTH FACTORS

VA-REMINDER UPDATES [C] VA-UPDATE 2 0 300 HYPERTENSION [C] HTN LIFESTYLE MODIFICATIONS RECOMMENDED VA-HTN SELF-RECORDED DIASTOLIC BLOOD PRESSURE VA-HTN SELF-RECORDED SYSTOLIC BLOOD PRESSURE VA-HTN REFERRED OTHER VA-HTN REFERRED TO TELEHEALTH **HTN REFERRED TO MOVE** HTN REFERRED TO DIETICIAN HTN REFERRED TO PACT PHARMACIST VA-VIDEO CONNECT [C] VA-VVC BP FOLLOW-UP APPOINTMENT OTHER VA-VVC BP FOLLOW-UP APPOINTMENT NA VA-VVC BP FOLLOW-UP APPOINTMENT RTC VA-VVC BP FOLLOW-UP APPOINTMENT

VA-PERSONAL HEALTH INVENTORY [C] VA-PHI DAY TO DAY LIFE WHY VA-PHI DAY TO DAY LIFE VA-PHI PHYSICAL WELL BEING WHY VA-PHI PHYSICAL WELL BEING VA-PHI MENTAL/EMOTIONAL WELL BEING WHY VA-PHI MENTAL/EMOTIONAL WELL BEING VA-PHP [C] VA-PHP MAP HTN REFUSAL OF MED INTERVENTION HTN NO MED CHANGE - LIMITED LIFE EXP HTN NO MED CHANGE - SIDE EFFECTS FROM RX **HTN NO MED CHANGE - NONCOMPLIANCE HTN NO MED CHANGE - COMORBID ILLNESS** HTN NO MED CHANGE - BP CONTROLLED HTN MEDICATIONS ADJUSTED FALL [C] FALL INJURY DID NOT CHANGE FUNCTIONALITY FALL INJURY CHANGED FUNCTIONALITY FALL CAUSED AN INJURY FALL DID NOT CAUSE INJURY FALL-MORE THAN 1 IN PAST 12 MONTHS FALL-NONE IN PAST 12 MONTHS FALL-ONE IN LAST 12 MONTHS FALL IN LAST YEAR **TELE URGENT CARE [C]** TELE URGENT TYPE VIDEO HF

REMINDER SPONSOR

VHA OFFICE OF CONNECTED CARE (10P8) OFFICE OF QUALITY & PERFORMANCE

REMINDER TAXONOMY

VA-HT ENCOUNTER PHONE 21 VA-HT ENCOUNTER PHONE 11 VA-HT ENCOUNTER PHONE 5

REMINDER TERM

VA-REMINDER UPDATE_2_0_300

HEALTH SUMMARY TYPE

VA-FALLS

HEALTH SUMMARY OBJECTS

VA-FALLS (TIU)

TIU DOCUMENT DEFINITION VA-FALLS

VA-FALLS

REMINDER DIALOG VA-GP TH VVC UPDATE JUL 22

Pre-Install

You should have installed Update 249 VA-STRENGTH AT HOME TEMPLATES and UPDATE 226 VA-VVC BLOOD PRESSURE TEMPLATE UPDATE prior to installing this update.

Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element, search for and save the inquiry for the following reminder dialog elements. You will need to update these in the post install.

VAL-TH VVC HYPERTENSION DISCLOSURE LAB OBJECT VAL-TH VVC HYPERTENSION DISCLOSURE MED AND ALLERGY OBJECT VAL-TH VVC HYPERTENSION PHONE NUMBER OBJECT VAL-TH VVC HYPERTENSION LOCATION HOME

Identify any orders that you have as additional findings on the following elements, you will add any orders back in the post install. In the post-install they will be named VAL-.

VA-TH VVC HTN REFERRAL PHARM VA-TH VVC HTN REFERRAL DIETITIAN VA-TH VVC HTN REFERRAL MOVE VA-TH VVC HTN REFERRAL OTHER VA-TH VVC HTN REFERRAL HT VA-TH VVC HYPERTENSION BP CUFF SIZE

Install Details

Installation:

=================

This update can be loaded with users on the system. Installation will take less than 5 minutes. This update is being distributed as a web host file.

The address for the host file for sites that have not upgraded to CPRS VERSION 32B is:

https://Redacted /UPDATE_2_0_300.PRD

The address for sites that HAVE upgraded to CPRS VERSION 32B is:

https://Redacted /UPDATE_2_0_30032B.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

You will only need to install ONE of these two files. Depending on your 32B install status.

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

At the **<u>Select Action</u>**: prompt, enter **<u>LWH</u>** for Load Web Host File

FOR SITES THAT HAVE NOT INSTALLED 32B:

At the Input the url for the .prd file: prompt, type the following web address

https://vaww.va.gov/reminders/docs/UPDATE_2_0_300.PRD

+	+	Next	Scree	n -	Prev	Screen	??	More	Actions	>>>
CFE	Create	Excha	ange F	ile E	Intry		LHF	Load	Host File	
CHF	Create	Host	File				LMM	Load	MailMan Messa	ige
CMM	Create	Mail	Man Me	ssage)		LR	List	Reminder Defi	nitions
DFE	Delete	Excha	ange F	ile E	Intry		LWH	Load	Web Host File	•
IFE	Instal	L Excl	hange	File	Entry		RI	Remin	nder Definitio	on Inquiry
IH	Instal	latio	n Hist	ory			RP	Repa	ck	
Selec	t Actio	on: Ne	ext Sc	reen/	/ LWH	Load	Web	Host H	File	
Input	t the UF	RL for	r the	.prd	file:	https:/			-	/UPDATE_2_0_300.PR

FOR SITES THAT HAVE INSTALLED 32B:

At the **Input the url for the .prd file:** prompt, type the following web address **https://vaww.va.gov/reminders/docs/UPDATE_2_0_30032B.PRD**

+	+ Next Screen	- Prev S	creen ?	? More	Actions	>>>
CFE	Create Exchange File E	Entry	LHF	Load	Host File	
CHF	Create Host File		LMM	Load	MailMan Message	3
CMM	Create MailMan Message	e	LR	List	Reminder Defin:	itions
DFE	Delete Exchange File E	Entry	LWH	Load	Web Host File	
IFE	Install Exchange File	Entry	RI	Remir	der Definition	Inquiry
IH	Installation History		RP	Repac	:k	
Selec	t Action: Next Screen,	// LWH	Load Web	Host F	ile	
Input	the URL for the .prd	Tile: h	ittps://			/UPDATE_2_0_30032B.PRD

You should see a message at the top of your screen that the file loaded successfully

FOR SITES THAT HAVE NOT INSTALLED 32B

Search and locate an entry titled **UPDATE_2_0_300 VA-TH-VVC TEMPLATE UPDATES** in reminder exchange.

+Item	Entry	Source	Date Packed
301	UPDATE_2_0_300 VA-TH VVC		07/26/2022@13:10
	UPDATE		
302	UPDATE_2_0_31 VA-REACH		07/05/2017@07:55
303	UPDATE_2_0_32 VA-MH		01/26/2018@04:12
	CONTINGENCY MANAGEMENT		
304	UPDATE_2_0_33 MISC UPDATE 4		07/07/2017@12:39
305	UPDATE_2_0_34 VA-WH PAP SMEAR		09/11/2019@11:39
	SCREENING		
306	UPDATE_2_0_35 VA-BENEFICIARY		12/22/2017@05:23
	TRAVEL UPDATES		
+	+ Next Screen - Prev Screen	?? More Actions	>>>
CFE C	reate Exchange File Entry	LHF Load Host File	
CHF C	reate Host File	LMM Load MailMan Messag	ge
CMM C	reate MailMan Message	LR List Reminder Defir	nitions
DFE D	elete Exchange File Entry	LWH Load Web Host File	
IFE I	nstall Exchange File Entry	RI Reminder Definition	n Inquiry
IH I	nstallation History	RP Repack	
Select	Action: Next Screen// ife Inst	all Exchange File Entry	
Enter	a list or range of numbers (1-522	2): 301	

At the <u>Select Action</u> prompt, enter <u>IFE</u> for Install Exchange File Entry Enter the number that corresponds with your entry **UPDATE_2_0_300 VA-TH-VVC TEMPLATE UPDATES** (*in this example it is entry 301 it will vary by site*). The date of the exchange file should be 07/27/2022.

FOR SITES THAT HAVE INSTALLED 32B

Search and locate an entry titled **UPDATE_2_0_30032B VA-TH-VVC TEMPLATE UPDATES** in reminder exchange.

+Item	Entry	Source	Date Packed
434	UPDATE_2_0_30032B VA-TH VVC		-31
	UPDATE		
435	UPDATE_2_0_31 VA-REACH		
436	UPDATE_2_0_32 VA-MH		
	CONTINGENCY MANAGEMENT		
437	UPDATE_2_0_33 MISC UPDATE 4		
438	UPDATE 2 0 34 VA-WH PAP SMEAR		
	SCREENING		
439	UPDATE_2_0_35 VA-BENEFICIARY		4
	TRAVEL UPDATES		
+	+ Next Screen - Prev Scree	n ?? More Actions	>>>
CFE	Create Exchange File Entry	LHF Load Host Fil	e
CHF	Create Host File	LMM Load MailMan	Message
CMM	Create MailMan Message	LR List Reminder	Definitions
DFE	Delete Exchange File Entry	LWH Load Web Host	File
IFE	Install Exchange File Entry	RI Reminder Defi	nition Inquiry
IH	Installation History	RP Repack	
Selec	t Action: Next Screen// IFE Ins	tall Exchange File E	ntry
Enter	a list or range of numbers (1-66	8): 434	

At the <u>Select Action</u> prompt, enter <u>IFE</u> for Install Exchange File Entry Enter the number that corresponds with your entry **UPDATE_2_0_30032B VA-TH-VVC TEMPLATE UPDATES** (*in this example it is entry 434 it will vary by site*). The date of the exchange file should be 07/27/2022. The following instructions apply to all sites (including CPRS 32B sites)

0.00000000				0		Eviete
Compone	nt			Ca	ategory	Exists
Source:						
Date Packed:						
Package Versi	on: 2.0P71					
Description:						
The following	Clinical Remind	er items were se	elected	for packin	ng:	
REMINDER DIAL	.OG				Ť	
VA-GP TH VV	C UPDATE JUL 22					
REMINDER TERM	1					
VA-REMINDER	UPDATE 2 0 300					
Non-exchangea	ble order dialog	(s):				
	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·				
	Name:	GMRAOR ALLERGY	ENTER/E	DIT		
	Type:	Action				
	Display Text:	Allergy/Adverse	e Reacti	on		
+ Ent	er ?? for more a	ctions				>>>
IA Install	all Components	IS	Install	Selected	Component	
Select Action	: Next Screen//	IA				

At the **<u>Select Action</u>** prompt, type <u>IA</u> for Install all Components and hit enter.

Select Action: Next Screen// IA Install all Components

You will see several prompts, for all new entries you will choose **I to Install** For components that already exists but the packed component is different, you will choose **S to SKIP**.

You will be prompted to install the reminder dialog component – You will choose to install all even though some components may have an X to indicate that they exist.

Packe	d reminder dialog: VA-GP TH VVC UPDATE JUL 22 [NATIONAL D	IALOG]	
Item	Seq. Dialog Findings	Туре	Exists
1	VA-GP TH VVC UPDATE JUL 22	group	
	Finding: *NONE*		
2	5 VA-TELE URGENT TYPE VIDEO	element	X
	Finding: TELE URGENT TYPE VIDEO HF (HEALTH FACTOR)		Х
3	10 VA-GP TH VVC HYPERTENSION DISCLOSURE	group	Х
	Finding: *NONE*		
4	10.5 VA-GP TH VVC HYPERTENSION LOCATION	group	Х
5	10 5 5 VAL-TH VVC HYPERTENSION LOCATION HOME	element	x
Ŭ	Finding: *NONE*	Crement	X
6	10.5.10 VA-TH VVC HYPERTENSION LOCATION OTHER	element	х
	Finding: *NONE*		
7	PXRM COMMENT	prompt	Х
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DT Dialog Text IS Insta	ll Selecte	d
DF	Dialog Findings DU Dialog Usage QU Quit		
DS	Dialog Summary IA Install All		
Selec	t Action: Next Screen// IA		

At the <u>Select Action</u> prompt, type <u>IA</u> to install the dialog – VA-GP TH VVC UPDATE JUL 22

Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// Yes

Packed reminder dialog: VA-GP TH VVC UPDATE JUL 22 [NATIONAL DI	ALOG]	
VA-GP TH VVC UPDATE JUL 22 (group) installed from exchange file	<mark>e.</mark> –	
Item Seq. Dialog Findings	Туре	Exists
1 VA-GP TH VVC UPDATE JUL 22	group	Х
Finding: *NONE*		
2 5 VA-TELE URGENT TYPE VIDEO	element	Х
Finding: TELE URGENT TYPE VIDEO HF (HEALTH FACTOR)		Х
3 10 VA-GP TH VVC HYPERTENSION DISCLOSURE	group	X
Finding: *NONE*		
4 10.5 VA-GP TH VVC HYPERTENSION LOCATION	group	X
Finding: *NUNE*	alamant	V
5 10.5.5 VAL-TH VVC HYPERTENSION LOCATION HOME	erement	X
FINGING: "NONE"	alamant	v
5 10.5.10 VA-TH VVC HIPERTENSION LOCATION OTHER Finding: *NONE*	erement	^
7 DVDM COMMENT	ppompt	V
+ + Next Screen - Prov Screen 22 Mere Actions	prompt	^
DD Dialog Dotails DT Dialog Toxt IS Instal	1 Soloctor	4
DE Dialog Eindings DU Dialog leage OU Ouit	T OFTECTED	4
DS Dialog Summary IA Install All		
Select Action: Next Screen// O		

After completing this dialog install, you will type Q.

Component	Category	Exists
Source: Date Packed: Package Version: 2.0P/1		
Description: The following Clinical Reminder items were selected for REMINDER DIALOG	packing:	
VA-GP TH VVC UPDATE JUL 22 REMINDER TERM		
VA-REMINDER UPDATE_2_0_300 Non-exchangeable order dialog(s):		
Name: GMRAOR ALLERGY ENTER/EDIT Type: Action		
Ulsplay Text: Allergy/Adverse Reaction	10	
IA Install all Components IS Install Sel Select Action: Next Screen// O	ected Component	22

You will then be returned to this screen. At the Select Action prompt, type Q.

Install complete.

Post Installation

 Add your local DATA OBJECT FOR PATIENT ADDRESS to the following element: VAL-TH VVC HYPERTENSION LOCATION HOME Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element. Search for and edit the following element: VAL-TH VVC HYPERTENSION LOCATION HOME.

- Add your local DATA OBJECT FOR PATIENT PHONE NUMBER to the following element: VAL-TH VVC HYPERTENSION PHONE NUMBER OBJECT Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 Search for and edit the following element: VAL-TH VVC HYPERTENSION PHONE NUMBER OBJECT
- Add your local DATA OBJECTS FOR MEDICATIONS AND ALLERGIES to the following element: VAL-TH VVC HYPERTENSION DISCLOSURE MED AND ALLERGY OBJECT Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element. Search for and edit the following element: VAL-TH VVC HYPERTENSION DISCLOSURE MED AND ALLERGY OBJECT
- Add a local DATA OBJECT FOR BMP WITHIN THE LAST 2 YEARS to the following element: VAL-TH VVC HYPERTENSION DISCLOSURE LAB OBJECT Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.
 Search for and edit the following element: VAL-TH VVC HYPERTENSION DISCLOSURE LAB OBJECT
- 5. Add referral orders identified in the pre-install to the following reminder dialog elements as an additional finding. Using the menu Reminder Dialog Management, Reminder Dialogs, change your view to Element.

VAL-TH VVC HTN REFERRAL PHARM VAL-TH VVC HTN REFERRAL DIETITIAN VAL-TH VVC HTN REFERRAL MOVE VAL-TH VVC HTN REFERRAL OTHER VAL-TH VVC HTN REFERRAL HT VA-TH VVC HYPERTENSION BP CUFF SIZE