

Update_2_0_331

Clinical Reminders

VA-Community Care COMPACT Act Update Install Guide

January 2023

Product Development Office of Information Technology Department of Veterans Affairs

Table of Contents	
INTRODUCTION	1
INSTALL DETAILS	
PRE-INSTALL	4
INSTALL EXAMPLE	5
POST INSTALLATION	7

Introduction

The Community Care-Care Coordination Plan Note and Community Care-Emer Self Presenting Care Coord Plan Notes are being updated to include documentation for COMPACT ACT.

The COMPACT ACT will be considered for approval beginning on 1-16-2023. Each of the above notes provides text indicating that the care provided is covered by the COMPACT ACT as per legislation.

Additional health factors were also added to these notes.

VA-COMMUNITY CARE-CARE COORDINATION PLAN NOTE:



VA-COMMUNITY CARE-EMER SP CARE COORD PLAN NOTE:



The Office of Community Care will provide training for staff using the updated templates.

UPDATE_2_0_331 contains 1 Reminder Exchange entry: UPDATE_2_0_331 VA-COMMUNITY CARE COMPACT ACT UPDATE

The exchange file contains the following components:

TIU TEMPLATE FIELD WORD 2 LINES1 REQ DATE & TIME ENTER DATE/TIME OTHER

TEXT (1-30 CHAR) REQ TEXT (1-50 CHAR) TEXT (1-60 CHAR)R SELECT ONE VA CC COORD NOTE FINISH WP 1 LINE REQ DATE (*) TEXT CC TRANSFER FORMS INFO WP 2/74 TEXT NOTES: **CC-TEXT TRANSFER INSTRUCTIONS 2 RB YES/NO** VA-CC TEXT REFUSAL FOR TRANFER INFO **CC-TEXT TRANSFER INSTRUCTIONS** TEXT (1-20 CHAR) TEXT 30 CHAR (REQ) TEXT (1-30 CHAR) **TEXT (1-60 CHAR)REQUIRED** WP 1 LINE **OPTIONAL COMMENT-DISPLAY ONLY** WP 2/74REQ TEXT (1-40 CHAR) REQ TEXT (1-20 CHAR) REQ TEXT (1-40 CHAR) TEXT (1-60 CHAR) **GEN TEXT BOX 50** DATE WP 1 LINE REQ1 TEXT RISK INFO **GEN TEXT BOX 20 SPACES** GEN BED TYPE **GEN TEXT BOX 25** CC CARE COORD PLAN INSTRUCTIONS TEXT (1-20 CHAR REQ) OPTIONAL TEXT TEXT (1-10 CHAR)

HEALTH FACTORS

VA-REMINDER UPDATES [C] VA-UPDATE_2_0_331 COMMUNITY CARE [C] VA-CCET SELF PRESENTING EOC COMPLETE VA-CCET SELF PRESENTING VETERAN HANDOFF CCET SELF PRESENTING DISCH DISP DEATH VA-CCET SELF PRESENTING DISCH DISP OTHER

CCET SELF PRESENTING DISCH DISP SKILLED NRS FACILITY CCET SELF PRESENTING DISCH DISP LTAC CCET SELF PRESENTED DISCH DISP ACUTE REHAB CCET SELF PRESENTING DISCH DISP CNH CCET SELF PRESENTING DISCH DISP FAMILY/OTHER PERSON CCET SELF PRESENTING DISCH DISP HOME CCET TRANS SELF PRESENTED OTHER CMF VA-COVID-19 [C] VA-COVID-19 OUTSIDE PCR SPEC COLLECTED VA-COVID-19 PCR LAB OUTSIDE NEGATIVE VA-COVID-19 PCR LAB OUTSIDE POSITIVE CCET TRANS SELF PRESENTED TO VA CCET TRANSFER SELF PRESENTED CCPN CARE COORD EPISODE OF CARE COMPLETE CCPN CARE COORD DISCH DISP HOME CCPN CARE COORD TRANSFER OTHER CMF CCPN CARE COORD TRANSFER TO VA CCPN CARE COORD TRANSFER CCET SUICIDE PREVENTION COORD NOTIFIED COM CARE HOSPITAL NOTE [C] ADMIN NOT ELIGIBLE FOR 38 U.S.C. 1703 CCET COMPACT ACT CCET SUICIDAL BEH SELF PRESENTED CCPN COMPACT ACT CCPN CAREGIVER AGENCY CCPN CAREGIVER INDIVIDUAL CCPN REFERRAL HOSPITAL INITIAL NOTE CCPN CARE COORD ADMIT DATE VA REF ADMIT YES CCPN LOC URGENT CCPN LOC COMPLEX/CHRONIC CCPN LOC MODERATE CCPN LOC BASIC CCPN CARE COORD INITIAL NOTE

REMINDER SPONSOR

VHA OFFICE OF COMMUNITY CARE

REMINDER TERM

VA-REMINDER UPDATE_2_0_331

REMINDER DIALOG

VA-GP COMMUNITY CARE COMPACT UPDATE DEC22

Install Details

This update is being distributed as a web host file. The address for the host file is: https://REDACTED/UPDATE_2_0_331.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

==============

This update can be loaded with users on the system. Installation will take less than 5 minutes.

Pre-Install

 Clinical Reminder Update_2_0_161 VA-Community Care Coordination Plan -Emergency Care Update should be installed before installing this update. Review the install history in reminder exchange to confirm that you have installed.

The install guide is available here: <u>EMERGENCY PLAN (REDACTED)</u>

2. Prior to installing this update, use the action INQ (Inquiry/Print) to print and save the existing reminder dialog for the following items. This will ensure the local modifications to these items are easier to complete when you complete the Post-Install steps. *Local modifications are optional, your facility may have no modifications*

ELEMENT	VAL-EL OPTIONAL LOCAL INFORMATION
ELEMENT	VAL-EL NON VA FACILITY INTAKE S/I SUICIDE COOR NAME
GROUP	VAL-GP CC REFERRAL HOSPITALS/PROVIDER CCPN
GROUP	VAL-GP CC REFERRAL HOSPITAL CCPN
GROUP	VAL-GP CC REFERRAL HOSPITAL INFO CCPN
GROUP	VAL-GP CC REFERRAL PROVIDER CCPN
GROUP	VAL-GP CC REFERRAL PROVIDER INFO CCPN
GROUP	VAL-GP COMMUNITY CARE COORD ADMIT SAME
GROUP	VAL-GP COMMUNITY CARE COORD PLAN VA FACILITY
	FOLLOW UP
GROUP	VAL-GP COMMUNITY CARE COORD TRAN FACILITY

Install Example

T	D Load the Web Host File. Navigate to	ס Reminder exchange in Vista
+	+ Next Screen - Prev Scree	n ?? More Actions >>>
С	FE Create Exchange File Entry	LHF Load Host File
C	IF Create Host File	LMM Load MailMan Message
Ç	4M Create MailMan Message	LR List Reminder Definitions
D	E Delete Exchange File Entry	LWH Load Web Host File
I	'E Install Exchange File Entry	RI Reminder Definition Inquiry
I	Installation History	RP Repack
8	elect Action: Next Screen// LWH Loa	d Web Host File
I	put the URL for the .prd file: https	:// /UPDATE_2_0_331.PRD
Δ.		

At the <u>Select Action:</u> prompt, enter <u>LWH</u> for Load Web Host File At the <u>Input the url for the .prd file:</u> prompt, type the following web address: https://REDACTED/UPDATE 2 0 331.PRD

You should see a message at the top of your screen that the file successfully loaded.



At the <u>Select Action</u> prompt, enter <u>IFE</u> for Install Exchange File Entry Enter the number that corresponds with your entry titled UPDATE_2_0_331 VA-COMMUNITY CARE COMPACT ACT UPDATE(*in this example it is entry 481, it will*

vary by site, the date packed should be 12/14/2022) Package Version: 2.0P65 Description: The following Clinical Reminder items were selected for packing: REMINDER DIALOG VA-GP COMMUNITY CARE COMPACT UPDATE DEC22 REMINDER TERM VA-REMINDER UPDATE_2_0_331 Keywords: Components: >>> IS IA Install Selected Component Select Action: Next Screen// IA

At the <u>Select Action</u> prompt, type **IA** for Install all Components and hit enter.

Select Action: Next Screen// IA Install all Components

During the installation, you will see several prompts:

Install all new components

If you are prompted to OVERWRITE any existing components select S for Skip. You will be promped to install the dialog component:

Packe	d reminder dialog: VA-GP COMMUNITY CARE COMPACT UPDATE	DEC22 [NATIO	NAL DIAL
Item	Seq. Dialog Findings	Туре	Exists
1	VA-GP COMMUNITY CARE COMPACT UPDATE DEC22	group	
	Finding: *NONE*		
2	5 VA-CC EMER SP CARE COORD NOTE VERSION	element	×
	Finding: *NONE*		
3	10 VA-CC CARE COORD PLAN NOTE VERSION	element	×
	Finding: *NONE*		
4	15 VA-GP NVCC NON VA FACILITY INTAKE S/I QUESTION	group	X
	Finding: CCET SUICIDAL BEH SELF PRESENTED (HEALT	H FACTOR)	X
5	15.5 VA-GP NVCC NON VA COMPACT APPROVED	group	×
	Finding: *NONE*		
6	15.5.10 VA-GP NVCC NON VA COMPACT APPROVED YES	group	X
	Finding: CCET COMPACT ACT (HEALTH FACTOR)		X
+	+ Next Screen - Prev Screen ?? More Actions		
DD	Dialog Details DI Dialog lext IS Ins	tall Selecte	d
DF	Dialog Findings DU Dialog Usage QU Qui	t	
bs	Dialog Summary IA Install All		
Selec	t Action: Next Screen// IA		

At the **<u>Select Action</u>** prompt, type <u>IA</u> to install the dialog – <u>VA-GP COMMUNITY CARE</u> COMPACT UPDATE DEC22 Select Action: Next Screen// IA Install All

Install reminder dialog and all components with no further changes: Y// ES

Packe	d reminder dialog: V/	A-GP CO	MMUNITY CA	ARE COMP	ACT UPDA	ATE DEC22	[NATIO	NAL DIAL
VA-GP	COMMUNITY CARE COMP	ACT UPD	ATE DEC22	(group)	instal	led from	exchang	e file.
Item	Seq. Dialog Finding	gs					Туре	Exists
1	VA-GP COMMUNITY C/ Finding: *NON	ARE COM E*	IPACT UPDAT	TE DEC22			group	×
2	5 VA-CC EMER SP CAN Finding: *NON	RE COOR E*	D NOTE VEF	RSION		el	ement	x
З	10 VA-CC CARE COORI Finding: *NONI	D PLAN E*	NOTE VERSI	ION		el	ement	x
4	15 VA-GP NVCC NON	VA FACI	LITY INTAK	KE S/I QI	UESTION		group	×
	Finding: CCET	SUICID	AL BEH SEL	F PRESE	NTED (HE	EALTH FAC	TOR)	×
5	15.5 VA-GP NVCC NO Finding: *NON	N VA CO E*	MPACT APP	ROVED			group	×
6	15.5.10 VA-GP NVCC	NON VA	COMPACT A	APPROVED	YES		group	×
	Finding: CCET	COMPAC	T ACT (HEA	ALTH FAC	TOR)			X
+	+ Next Screen	- Prev	Screen	?? More	Actions	5		
DD	Dialog Details	DT D	ialog Text	È	IS	Install	Selecte	d
DF	Dialog Findings	DU D	ialog Usag	ge	QU	Quit		
DS	Dialog Summary	IA I	nstall All	L				
Selec	t Action: Next Scree	n// Q						

When the dialog has completed installation, you will then be returned to this screen. At the <u>Select Action</u> prompt, type Q.



You will be returned to this screen. Type Q to quit. Install Completed

Post Installation

- 1. Open the note titled COMMUNITY CARE-CARE COORDINATION PLAN NOTE reminder dialog VA-COMMUNITY CARE-CARE COORD PLAN NOTE and confirm the version is 2.1.
- Open the note titled COMMUNITY CARE-EMER SELF PRESENTING CARE COORD PLAN NOTE reminder dialog VA-COMMUNITY CARE-EMER SP CARE and confirm the version is 2.1.
- **3.** Refer to your saved Inquiry from the pre-installation steps, and add back any local modifications to the following items. *Local modifications are optional, your facility may have no modifications*

ELEMENT	VAL-EL OPTIONAL LOCAL INFORMATION	
ELEMENT	VAL-EL NON VA FACILITY INTAKE S/I SUICIDE COOR NAME	
GROUP	VAL-GP CC REFERRAL HOSPITALS/PROVIDER CCPN	
GROUP	VAL-GP CC REFERRAL HOSPITAL CCPN	
GROUP	VAL-GP CC REFERRAL HOSPITAL INFO CCPN	
GROUP	VAL-GP CC REFERRAL PROVIDER CCPN	
GROUP	VAL-GP CC REFERRAL PROVIDER INFO CCPN	
GROUP	VAL-GP COMMUNITY CARE COORD ADMIT SAME	
GROUP	VAL-GP COMMUNITY CARE COORD PLAN VA FACILITY	
	FOLLOW UP	
GROUP	VAL-GP COMMUNITY CARE COORD TRAN FACILITY	