



**Update\_2\_0\_455**

**Clinical Reminders**

**VA-Asynchronous Telehealth Templates**

**Install Guide**

---

**July 2024**

Product Development

Office of Information Technology  
Department of Veterans Affairs

# Table of Contents

<b>INTRODUCTION .....</b>	<b>1</b>
<b>INSTALL DETAILS .....</b>	<b>20</b>
<b>PRE-INSTALL .....</b>	<b>20</b>
<b>POST INSTALLATION .....</b>	<b>27</b>



## **Introduction**

This reminder update contains changes to 4 existing reminder dialog templates.

1. TIME SENSITIVE EMERGENT URGENT DISCLOSURE DIALOG (Note: this template is embedded into the two existing **Telehealth Reader notes** for Eye and Dermatology)
  - Directive 1088 requires closed looped communication from interpreting provider to referring provider only for emergent/imminently life-threatening findings.
  - While closed loop communication is not required for urgent findings, it does require communication of that, which happens through the consult closure and view alerts packaging process which is in place currently.
  - This template supports the documentation of time sensitive, urgent/emergent results and communication of these results to the Referring Provider or designee.
  
2. STORE & FORWARD (SFT) RESULTS COMMUNICATION DIALOG (Note: this template is embedded into the two existing **Telehealth Imager notes** for Eye and Dermatology)
  - Per the directive, all results requiring action require documentation of communication to the patient of the results and care plan BY the REFERRING provider or designee.
  - This template supports awareness of time sensitive, urgent/emergent result communication with an actionable process and enables reporting capabilities.

Additional changes were made to the tele-eye templates based on feedback from the field.

### **UPDATE\_2\_0\_455 contains 1 Reminder Exchange entry:**

UPDATE\_2\_0\_455 VA-ASYNCHRONOUS TELEHEALTH TEMPLATES

### **The exchange file contains the following components:**

#### **REMINDER GENERAL FINDINGS**

VIEW PROGRESS NOTE TEXT

#### **TIU TEMPLATE FIELD**

VA-TDR ATTENDING NAMES  
VA\*WP 65/6 INDENT2  
ALLFT30C  
VA-TELE METHOD OF COMMUNICATION  
DATE & TIME (REQ)  
VA-TELE REFERRING PROVIDER  
VHA DIRECTIVE 1088  
VA\*TDI TREATMENT INFO  
DATE  
VA-EB 65/240 I3\*

VA-TELE VET METHOD OF COMMUNICATION  
VA-TDI WP 65/2 REQ  
VA\*TDI WP 74/20 CONSULT REQ  
VA\*TDI CONSULT REQUEST INFO  
VA-TRR DATE 9M  
VA-TRR DATE 2Y  
VA-TRR DATE 1Y  
VA-TRR DATE 6M  
VA-TRR DATE 4M  
VA-TRR DATE 3M  
VA-TRR DATE 2M  
VA-TRR DATE 1M  
NUM REQ  
TEXT (1-40 CHAR)  
VA-TRI IOP DROPS  
BLANK LINE 1  
DATE (\*)  
VA-TRI EDUC GIVEN  
TEXT (1-50 CHAR) REQ  
DATE (MONTH/YEAR ONLY REQ)  
VA-TOBACCO STATUS  
BLANK SPACE1  
VA-YES NO REQUIRED  
WP FREE TEXT FL55 I3  
COMMENT OPTIONAL  
VA-TECS OD/OS/OU  
MSCOE BLANK LINE FOR DIALOG USE  
TEXT (1-60 CHAR)  
VA-TEXT (1-75 CHAR)

## **EDUCATION TOPICS**

VA-TDI EDUCATION

## **EXAM**

TRR EXAM OTH ASSESS-U/A  
TRR EXAM OTH ASSESS-ABN  
TRR EXAM OTH ASSESS-NORM  
TRR EXAM OPTC NRV-NORMAL  
TRR EXAM OPTC NRV-ABN  
TRR EXAM OPTC NRV-U/A  
TRR EXAM MACULA-U/A  
TRR EXAM MACULA-ABN  
TRR EXAM MACULA-NORMAL  
VA-TRR DM EXAM RTNOPTY-U/A

VA-TRR DM EXAM RTNOPTY-ABN  
VA-TRR DM EXAM RTNOPTY-NORMAL

## HEALTH FACTORS

VA-REMINDER UPDATES [C]  
VA-UPDATE\_2\_0\_455  
VA-TELE SFT RESULTS FINDING DISCLOSURE [C]  
VA-TELE SFT RESULTS TIME-SENSITIVE EMERGENT  
VA-TELE SFT RESULTS TIME-SENSITIVE URGENT  
VA-TELE SFT RESULTS NOT TIME-SENSITIVE  
VA-TDR CATEGORY-PROBLEM J [C]  
VA-TDR READER TX CRYOTHERAPY PROB J  
VA-TDR CATEGORY-EXAM RECOMMENDATIONS [C]  
VA-TDR READER TX RECOMMEND NONE  
VA-TDR READER TX CULTURES PROBLEM J  
VA-TDR READER TX MEDS ADD PROBLEM J  
VA-TDR READER TX SKIN CARE PROBLEM J  
VA-TDR READER TX OTHER PROBLEM J  
VA-TDR READER TX BIOPSY PROBLEM J  
VA-TDR READER TX MEDS PROB J  
VA-TDR CATEGORY-PROBLEM I [C]  
VA-TDR READER TX CRYOTHERAPY PROB I  
VA-TDR READER TX CULTURES PROBLEM I  
VA-TDR READER TX MEDS ADD PROBLEM I  
VA-TDR READER TX SKIN CARE PROBLEM I  
VA-TDR READER TX OTHER PROBLEM I  
VA-TDR READER TX BIOPSY PROBLEM I  
VA-TDR READER TX MEDS PROB I  
VA-TDR CATEGORY-PROBLEM H [C]  
VA-TDR READER TX CRYOTHERAPY PROB H  
VA-TDR READER TX CULTURES PROBLEM H  
VA-TDR READER TX MEDS ADD PROBLEM H  
VA-TDR READER TX SKIN CARE PROBLEM H  
VA-TDR READER TX OTHER PROBLEM H  
VA-TDR READER TX BIOPSY PROBLEM H  
VA-TDR READER TX MEDS PROB H  
VA-TDR CATEGORY-PROBLEM G [C]  
VA-TDR READER TX CRYOTHERAPY PROB G  
VA-TDR READER TX CULTURES PROBLEM G  
VA-TDR READER TX MEDS ADD PROBLEM G  
VA-TDR READER TX SKIN CARE PROBLEM G  
VA-TDR READER TX OTHER PROBLEM G  
VA-TDR READER TX BIOPSY PROBLEM G  
VA-TDR READER TX MEDS PROB G

VA-TDR CATEGORY-PROBLEM F [C]  
VA-TDR READER TX CRYOTHERAPY PROB F  
VA-TDR READER TX CULTURES PROBLEM F  
VA-TDR READER TX MEDS ADD PROBLEM F  
VA-TDR READER TX SKIN CARE PROBLEM F  
VA-TDR READER TX OTHER PROBLEM F  
VA-TDR READER TX BIOPSY PROBLEM F  
VA-TDR READER TX MEDS PROB F  
VA-TDR CATEGORY-PROBLEM E [C]  
VA-TDR READER TX CRYOTHERAPY PROB E  
VA-TDR READER TX CULTURES PROBLEM E  
VA-TDR READER TX MEDS ADD PROBLEM E  
VA-TDR READER TX SKIN CARE PROBLEM E  
VA-TDR READER TX OTHER PROBLEM E  
VA-TDR READER TX BIOPSY PROBLEM E  
VA-TDR READER TX MEDS PROB E  
VA-TDR CATEGORY-PROBLEM D [C]  
VA-TDR READER TX CRYOTHERAPY PROB D  
VA-TDR READER TX CULTURES PROBLEM D  
VA-TDR READER TX MEDS ADD PROBLEM D  
VA-TDR READER TX SKIN CARE PROBLEM D  
VA-TDR READER TX OTHER PROBLEM D  
VA-TDR READER TX BIOPSY PROBLEM D  
VA-TDR READER TX MEDS PROB D  
VA-TDR CATEGORY-PROBLEM C [C]  
VA-TDR READER TX CRYOTHERAPY PROB C  
VA-TDR READER TX CULTURES PROBLEM C  
VA-TDR READER TX MEDS ADD PROBLEM C  
VA-TDR READER TX SKIN CARE PROBLEM C  
VA-TDR READER TX OTHER PROBLEM C  
VA-TDR READER TX BIOPSY PROBLEM C  
VA-TDR READER TX MEDS PROB C  
VA-TDR CATEGORY-PROBLEM B [C]  
VA-TDR READER TX CRYOTHERAPY PROB B  
VA-TDR READER TX CULTURES PROBLEM B  
VA-TDR READER TX MEDS ADD PROBLEM B  
VA-TDR READER TX SKIN CARE PROBLEM B  
VA-TDR READER TX OTHER PROBLEM B  
VA-TDR READER TX BIOPSY PROBLEM B  
VA-TDR READER TX MEDS PROB B  
VA-TDR CATEGORY-PROBLEM A [C]  
VA-TDR READER TX CRYOTHERAPY PROB A  
VA-TDR READER TX CULTURES PROBLEM A  
VA-TDR READER TX MEDS ADD PROBLEM A



VA-TDR READER TX SKIN CARE PROBLEM A  
VA-TDR READER TX OTHER PROBLEM A  
VA-TDR READER TX BIOPSY PROBLEM A  
VA-TDR READER TX MEDS PROB A  
VA-TDR CATEGORY-EXAM F/U [C]  
VA-TDR REFER COMMUNITY CARE  
VA-TDR RECONSULT TELEDERM  
VA-TDR RETURN F/U NOT REQ  
VA-TDR REFER DERM  
VA-TDR RETURN PCC F/U  
VA-TDR READER UNABLE TO ASSESS PROB J  
VA-TDR DX OTHER J  
VA-TDR DX LESION J-SEB KERATOSIS  
VA-TDR DX LESION J-NEVUS  
VA-TDR DX LESION J-DYSPLASTIC NEVUS  
VA-TDR DX LESION J-OTHER MALIGNANT  
VA-TDR DX LESION J-SCC  
VA-TDR DX LESION J-BCC  
VA-TDR DX LESION J-NEOPLASM UNC BEHAV  
VA-TDR DX LESION J-ACTINIC KERATOSIS  
VA-TDR DX LESION J-NON MALIG TUMOR  
VA-TDR DX LESION J-MELANOMA  
VA-TDR DX RASH J  
VA-TDR READER UNABLE TO ASSESS PROB I  
VA-TDR DX OTHER I  
VA-TDR DX LESION I-SEB KERATOSIS  
VA-TDR DX LESION I-NEVUS  
VA-TDR DX LESION I-DYSPLASTIC NEVUS  
VA-TDR DX LESION I-OTHER MALIGNANT  
VA-TDR DX LESION I-SCC  
VA-TDR DX LESION I-BCC  
VA-TDR DX LESION I-NEOPLASM UNC BEHAV  
VA-TDR DX LESION I-ACTINIC KERATOSIS  
VA-TDR DX LESION I-NON MALIG TUMOR  
VA-TDR DX LESION I-MELANOMA  
VA-TDR DX RASH I  
VA-TDR READER UNABLE TO ASSESS PROB G  
VA-TDR DX OTHER G  
VA-TDR DX LESION G-SEB KERATOSIS  
VA-TDR DX LESION G-NEVUS  
VA-TDR DX LESION G-DYSPLASTIC NEVUS  
VA-TDR DX LESION G-OTHER MALIGNANT  
VA-TDR DX LESION G-SCC  
VA-TDR DX LESION G-BCC

VA-TDR DX LESION G-NEOPLASM UNC BEHAV  
VA-TDR DX LESION G-ACTINIC KERATOSIS  
VA-TDR DX LESION G-NON MALIG TUMOR  
VA-TDR DX LESION G-MELANOMA  
VA-TDR DX RASH G  
VA-TDR READER UNABLE TO ASSESS PROB F  
VA-TDR DX OTHER F  
VA-TDR DX LESION F-SEB KERATOSIS  
VA-TDR DX LESION F-NEVUS  
VA-TDR DX LESION F-DYSPLASTIC NEVUS  
VA-TDR DX LESION F-OTHER MALIGNANT  
VA-TDR DX LESION F-SCC  
VA-TDR DX LESION F-BCC  
VA-TDR DX LESION F-NEOPLASM UNC BEHAV  
VA-TDR DX LESION F-ACTINIC KERATOSIS  
VA-TDR DX LESION F-NON MALIG TUMOR  
VA-TDR DX LESION F-MELANOMA  
VA-TDR DX RASH F  
VA-TDR READER UNABLE TO ASSESS PROB H  
VA-TDR DX OTHER H  
VA-TDR DX LESION H-SEB KERATOSIS  
VA-TDR DX LESION H-NEVUS  
VA-TDR DX LESION H-DYSPLASTIC NEVUS  
VA-TDR DX LESION H-OTHER MALIGNANT  
VA-TDR DX LESION H-SCC  
VA-TDR DX LESION H-BCC  
VA-TDR DX LESION H-NEOPLASM UNC BEHAV  
VA-TDR DX LESION H-ACTINIC KERATOSIS  
VA-TDR DX LESION H-NON MALIG TUMOR  
VA-TDR DX LESION H-MELANOMA  
VA-TDR DX RASH H  
VA-TDR READER UNABLE TO ASSESS PROB E  
VA-TDR DX OTHER E  
VA-TDR DX LESION E-SEB KERATOSIS  
VA-TDR DX LESION E-NEVUS  
VA-TDR DX LESION E-DYSPLASTIC NEVUS  
VA-TDR DX LESION E-OTHER MALIGNANT  
VA-TDR DX LESION E-SCC  
VA-TDR DX LESION E-BCC  
VA-TDR DX LESION E-NEOPLASM UNC BEHAV  
VA-TDR DX LESION E-ACTINIC KERATOSIS  
VA-TDR DX LESION E-NON MALIG TUMOR  
VA-TDR DX LESION E-MELANOMA  
VA-TDR DX RASH E

VA-TDR READER UNABLE TO ASSESS PROB D  
VA-TDR DX OTHER D  
VA-TDR DX LESION D-SEB KERATOSIS  
VA-TDR DX LESION D-NEVUS  
VA-TDR DX LESION D-DYSPLASTIC NEVUS  
VA-TDR DX LESION D-OTHER MALIGNANT  
VA-TDR DX LESION D-SCC  
VA-TDR DX LESION D-BCC  
VA-TDR DX LESION D-NEOPLASM UNC BEHAV  
VA-TDR DX LESION D-ACTINIC KERATOSIS  
VA-TDR DX LESION D-NON MALIG TUMOR  
VA-TDR DX LESION D-MELANOMA  
VA-TDR DX RASH D  
VA-TDR READER UNABLE TO ASSESS PROB C  
VA-TDR DX OTHER C  
VA-TDR DX LESION C-SEB KERATOSIS  
VA-TDR DX LESION C-NEVUS  
VA-TDR DX LESION C-DYSPLASTIC NEVUS  
VA-TDR DX LESION C-OTHER MALIGNANT  
VA-TDR DX LESION C-SCC  
VA-TDR DX LESION C-BCC  
VA-TDR DX LESION C-NEOPLASM UNC BEHAV  
VA-TDR DX LESION C-ACTINIC KERATOSIS  
VA-TDR DX LESION C-NON MALIG TUMOR  
VA-TDR DX LESION C-MELANOMA  
VA-TDR DX RASH C  
VA-TDR READER UNABLE TO ASSESS PROB B  
VA-TDR DX OTHER B  
VA-TDR DX LESION B-NEVUS  
VA-TDR DX LESION B-NON MALIG TUMOR  
VA-TDR DX LESION B-SEB KERATOSIS  
VA-TDR DX LESION B-DYSPLASTIC NEVUS  
VA-TDR DX LESION B-ACTINIC KERATOSIS  
VA-TDR DX LESION B-OTHER MALIGNANT  
VA-TDR DX LESION B-MELANOMA  
VA-TDR DX LESION B-SCC  
VA-TDR DX LESION B-BCC  
VA-TDR DX LESION B-NEOPLASM UNC BEHAV  
VA-TDR DX RASH B  
VA-TDR READER UNABLE TO ASSESS PROB A  
VA-TDR DX OTHER A  
VA-TDR DX LESION A-NON MALIG TUMOR  
VA-TDR DX LESION A-SEB KERATOSIS  
VA-TDR DX LESION A-NEVUS

VA-TDR DX LESION A-DYSPLASTIC NEVUS  
VA-TDR DX LESION A-ACTINIC KERATOSIS  
VA-TDR DX LESION A-OTHER MALIGNANT  
VA-TDR DX LESION A-MELANOMA  
VA-TDR DX LESION A-SCC  
VA-TDR DX LESION A-BCC  
VA-TDR DX LESION A-NEOPLASM UNC BEHAV  
VA-TDR DX RASH A  
VA-TDR CATEGORY-IMAGE QUALITY [C]  
VA-TDR IMAGE QUALITY IMPROVEMENT  
VA-TDR CONSULT QUALITY INADEQUATE  
VA-TDR IMAGE QUALITY INADEQUATE  
VA-TDR IMAGE QUALITY HIGH  
VA-TELE SFT TIME-SENSITIVE RESULTS COMMUNICATION [C]  
VA-TELE SFT TIME-SENSITIVE DATE COMMUNICATED  
VA-TELE SFT TIME-SENSITIVE RESULTS COMMUNICATION  
VA-TDI CATEGORY-IMAGE QUALITY [C]  
VA-TDI IMAGE ACQUIRED FROM VETERAN  
VA-TDI IMAGE ACQUIRED TRADITIONALLY  
VA-TDI CATEGORY-HISTORY INFO [C]  
VA-TDI HX IMAGER COMMENT  
VA-TDI HX PROBLEM-MELANOMA FAMILY NO  
VA-TDI HX PROBLEM-MELANOMA FAMILY YES  
VA-TDI PCP HX PRIOR SKIN DISORDERS  
VA-TDI HX IMMUNOSUPPRESSION  
VA-TDI PCP SIG MED HX  
VA-TDI HX NEW MED YES  
VA-TDI CATEGORY-PROBLEM J [C]  
VA-TDI HX PROBLEM J BX YES  
VA-TDI HX PROBLEM J SX  
VA-TDI HX PROBLEM J TX YES  
VA-TDI HX PROBLEM J CHANGE  
VA-TDI HX PROBLEM J DURATION  
VA-TDI CATEGORY-PROBLEM I [C]  
VA-TDI HX PROBLEM I BX YES  
VA-TDI HX PROBLEM I SX  
VA-TDI HX PROBLEM I TX YES  
VA-TDI HX PROBLEM I CHANGE  
VA-TDI HX PROBLEM I DURATION  
VA-TDI CATEGORY-PROBLEM H [C]  
VA-TDI HX PROBLEM H BX YES  
VA-TDI HX PROBLEM H SX  
VA-TDI HX PROBLEM H TX YES  
VA-TDI HX PROBLEM H CHANGE

VA-TDI HX PROBLEM H DURATION  
VA-TDI CATEGORY-PROBLEM G [C]  
VA-TDI HX PROBLEM G BX YES  
VA-TDI HX PROBLEM G SX  
VA-TDI HX PROBLEM G TX YES  
VA-TDI HX PROBLEM G CHANGE  
VA-TDI HX PROBLEM G DURATION  
VA-TDI CATEGORY-PROBLEM F [C]  
VA-TDI HX PROBLEM F BX YES  
VA-TDI HX PROBLEM F SX  
VA-TDI HX PROBLEM F TX  
VA-TDI HX PROBLEM F CHANGE  
VA-TDI HX PROBLEM F DURATION  
VA-TDI CATEGORY-PROBLEM D [C]  
VA-TDI HX PROBLEM D BX YES  
VA-TDI HX PROBLEM D SX  
VA-TDI HX PROBLEM D TX YES  
VA-TDI HX PROBLEM D CHANGE  
VA-TDI HX PROBLEM D DURATION  
VA-TDI CATEGORY-PROBLEM E [C]  
VA-TDI HX PROBLEM E BX YES  
VA-TDI HX PROBLEM E SX  
VA-TDI HX PROBLEM E TX YES  
VA-TDI HX PROBLEM E CHANGE  
VA-TDI HX PROBLEM E DURATION  
VA-TDI CATEGORY-PROBLEM C [C]  
VA-TDI HX PROBLEM C BX YES  
VA-TDI HX PROBLEM C SX  
VA-TDI HX PROBLEM C TX YES  
VA-TDI HX PROBLEM C CHANGE  
VA-TDI HX PROBLEM C DURATION  
VA-TDI CATEGORY-PROBLEM B [C]  
VA-TDI HX PROBLEM B BX YES  
VA-TDI HX PROBLEM B SX  
VA-TDI HX PROBLEM B TX YES  
VA-TDI HX PROBLEM B CHANGE  
VA-TDI HX PROBLEM B DURATION  
VA-TDI CATEGORY-PROBLEM A [C]  
VA-TDI HX PROBLEM A BX YES  
VA-TDI HX PROBLEM A TX YES  
VA-TDI HX PROBLEM A SX  
VA-TDI HX PROBLEM A CHANGE  
VA-TDI HX PROBLEM A DURATION  
VA-TDI DERM VISIT PRIOR TDI CONSULT

VA-TDI DERM VISIT PRIOR DERM VISIT  
VA-TDI DERM VISIT NEW CONDITION  
VA-TDI CONSENT NO  
VA-TDI CATEGORY-EDUCATION [C]  
VA-TDI EDUCATION DONE  
VA-TRR CATEGORY-RECOMMENDATIONS [C]  
VA-TRR RECOMMEND FULL EYE EXAM-1Y  
VA-TRR REFER OUTSIDE PROV KEEP SCHEDULED APPT  
VA-TRR REFER OUTSIDE PROV 9 MONTHS  
VA-TRR REFER OUTSIDE PROV 2 YRS  
VA-TRR REFER OUTSIDE PROV 1 YR  
VA-TRR REFER OUTSIDE PROV 6 MONTHS  
VA-TRR REFER OUTSIDE PROV 4 MONTHS  
VA-TRR REFER OUTSIDE PROV 3 MONTHS  
VA-TRR REFER OUTSIDE PROV 2 MONTHS  
VA-TRR REFER OUTSIDE PROV 1 MONTH  
VA-TRR REFER OUTSIDE PROV ASAP  
VA-TRR RECOMMEND VA EYE CLINIC APPT  
VA-TRR REFER OPTOMETRY 1 MONTH  
VA-TRR REFER OPTOMETRY 2 MONTHS  
VA-TRR REFER OPTOMETRY 9 MONTHS  
VA-TRR REFER OPTOMETRY 2 YRS  
VA-TRR REFER OPTOMETRY 1 YR  
VA-TRR REFER OPTOMETRY 6 MONTHS  
VA-TRR REFER OPTOMETRY 4 MONTHS  
VA-TRR REFER OPTOMETRY 3 MONTHS  
VA-TRR REFER OPTOMETRY ASAP  
VA-TRR REFER OPHTHALMOLOGY 9 MONTHS  
VA-TRR REFER OPHTHALMOLOGY 4 MONTHS  
VA-TRR REFER OPHTHALMOLOGY 1 MONTH  
VA-TRR REFER OPHTHALMOLOGY 2 YRS  
VA-TRR REFER OPHTHALMOLOGY 1 YR  
VA-TRR REFER OPHTHALMOLOGY 6 MONTHS  
VA-TRR REFER OPHTHALMOLOGY 3 MONTHS  
VA-TRR REFER OPHTHALMOLOGY 2 MONTHS  
VA-TRR REFER OPHTHALMOLOGY ASAP  
VA-TRR REFER OPHTH/OPTOM 2 MONTHS  
VA-TRR REFER OPHTH/OPTOM 1 MONTH  
VA-TRR REFER OPHTH/OPTOM 9 MONTHS  
VA-TRR REFER OPHTH/OPTOM 2 YRS  
VA-TRR REFER OPHTH/OPTOM 1 YR  
VA-TRR REFER OPHTH/OPTOM 6 MONTHS  
VA-TRR REFER OPHTH/OPTOM 4 MONTHS  
VA-TRR REFER OPHTH/OPTOM 3 MONTHS

VA-TRR REFER OPHTH/OPTOM ASAP  
VA-TRR REPEAT IMAGING-RETURN 6 MONTHS  
VA-TRR REPEAT IMAGING-RETURN OTH  
VA-TRR REPEAT IMAGING-RETURN 2 YRS  
VA-TRR REPEAT IMAGING-RETURN 1 YR  
VA-TRR REPEAT IMAGING-RTN SELECT DATE  
VA-TRR CATEGORY-SCREEN [C]  
VA-TRR NOT CANDIDATE FOR SCRNM  
VA-TRR CATEGORY-EYE SCREENING EXAM FINDINGS [C]  
VA-TRR OTH ASSESS EXAM UNABLE TO ASSESS  
VA-TRI CATEGORY-DM EYE EXAM [C]  
VA-TRI OTH ASSESS UNABLE TO ASSESS  
VA-TRR OTH ASSESS EXAM ABNORMAL  
VA-TRI OTH ASSESS EXAM ABNORMAL  
VA-TRR OTH ASSESS EXAM NORMAL  
VA-TRI OTH ASSESS EXAM NORMAL  
VA-TRR OPTIC NRV EXAM UNABLE TO ASSES  
VA-EYE CARE AT-RISK SCREEN [C]  
VA-TELE EYE GLAUC EXAM UNABLE TO ASSESS  
VA-TRI OPTIC NRV EXAM UNABLE TO ASSESS  
VA-TRR OPTIC NRV EXAM ABNORMAL  
VA-TELE EYE GLAUCOMA EXAM ABNORMAL  
VA-TRI OPTIC NRV EXAM ABNORMAL  
VA-TRR OPTIC NRV EXAM NORMAL  
VA-TELE EYE GLAUCOMA EXAM NORMAL  
VA-TRI OPTIC NRV EXAM NORMAL  
VA-TRR MACULA EXAM UNABLE TO ASSESS  
VA-TRI MACULA EXAM UNABLE TO ASSESS  
VA-TRR MACULA EXAM ABNORMAL  
VA-TRI MACULA EXAM ABNORMAL  
VA-TRR MACULA EXAM NORMAL  
VA-TRI MACULA EXAM NORMAL  
VA-TRR DM RETINOP EXAM UNABLE TO ASSESS  
VA-TRI DM RETINOP EXAM UNABLE TO ASSESS  
VA-TRR DM RETINOP EXAM ABNORMAL  
VA-TRI DM RETINOP EXAM ABNORMAL  
VA-TRR DM RETINOP EXAM NORMAL  
VA-TRI DM RETINOP EXAM NORMAL  
VA-TRI DM EYE CONSULT DONE  
VA-TRI DM IOP LT  
VA-TRI DM IOP RT  
VA-TRI DM LAST RET EVAL-COMP F-T-F  
VA-TRI DM LAST RET EVAL-SELF REPORTED  
VA-TRI LAST RET EVAL-RECORDS NOT REVIEWED

VA-TRI LAST RET EVAL-RECORDS REVIEWED  
VA-TRI DM LAST RET EVAL-DATE UNKNOWN  
VA-TRI DM LAST RETINAL EVAL  
VA-TRI DM EXAM PUPIL DILATED  
VA-TRI DM EXAM PUPIL NOT DILATED  
VA-TRI DM EYECARE PREFERENCE-NON VA  
VA-TRI CATEGORY-DM LAB RESULTS [C]  
VA-TRI DM LAB HGBA1C  
VA-INTRAOCULAR INJECTION HX  
VA-TRI CATEGORY-DM COMPLICATIONS [C]  
VA-TRI DM HX LASER THERAPY RETINOPATHY  
VA-TRI CATEGORY-DM THERAPY [C]  
VA-TRI DM CURRENT THERAPY-NONE  
VA-TRI DM CURRENT THERAPY-UNKNOWN  
VA-TRI DM CURRENT THERAPY-INSULIN  
VA-TRI DM CURRENT THERAPY-ORAL AGENT  
VA-TRI DM CURRENT THERAPY-DIET  
VA-TRI CATEGORY-DM DURATION [C]  
VA-TRI DM DURATION-CHART REVIEW  
VA-TRI DM DURATION-SELF REPORT  
VA-TRI DM DURATION 6 OR MORE YEARS  
VA-TRI DM DURATION 1-5 YEARS  
VA-TRI DM DURATION < 1 YEAR  
VA-TRI CATEGORY-DIABETES TYPE [C]  
VA-TRI DIABETES TYPE 2 PATIENT  
VA-TRI DIABETES TYPE 1 PATIENT  
VA-TRR REFER UNABLE TO ASSESS  
VA-TRR REFER COMP EYE CARE  
VA-TRR REFER VISUALLY SIGNF CATARACT  
VA-TRR REFER OTHER EYE DISEASE  
VA-TRR REFER OPTIC NERVE DISEASE  
VA-TRR REFER MACULA DISEASE  
VA-TRR REFER RETINOPATHY  
VA-TRR CATEGORY-NON DM EYE EXAM [C]  
VA-TRR AT RISK OTHER EYE CONDITION  
VA-TRR AT RISK FOR MAC DEGEN & GLAUCOMA  
VA-TRR AT RISK FOR GLAUCOMA  
VA-TRR AT RISK FOR MAC DEGENERATION  
VA-TRR CATEGORY-NON DM FINDINGS LT EYE [C]  
VA-TRR OTH ASSESS FIND INADQ ASSESS LT  
VA-TRR OTH ASSESS FINDING INTRARET HEM LT  
VA-TRR OTH ASSESS FINDING FAMILY HX GLAUCOMA LT  
VA-TRR OTH ASSESS FINDING PERIPH RET DEGEN LT  
VA-TRR OTH ASSESS FINDING COTTON WOOL SPOT LT



VA-TRR OTH ASSESS FINDING RETINAL VASC CHANGE LT  
VA-TRR OTH ASSESS FINDING CHORIORETINAL SCAR LT  
VA-TRR OTH ASSESS FINDING LENTICULAR OPACIFICATION LT  
VA-TRR OTH ASSESS FINDING LENS IMPLANT LT  
VA-TRR ON FINDING ELEVATED IOP LT  
VA-TRR OTH ASSESS FINDING OTH LT  
VA-TRR OTH ASSESS FINDING CHOR NEVUS LT  
VA-TRR OTH ASSESS FINDING HTN RTNPHTY LT  
VA-TRR OTH ASSESS FINDING OCCLUSION LT  
VA-TRR OTH ASSESS FINDING RET EMBOLUS LT  
VA-TRR OTH ASSESS FINDING CATARACT LT  
VA-TRR OTH ASSESS FINDING NO ABN FIND LT  
VA-TRR ON FINDING PALLOR LT  
VA-TRR ON FINDING PAPILLEDEMA SUSPECTED LT  
VA-TRR ON FINDING DISC EDEMA SUSP LT  
VA-TRR ON FINDING DISC W/ SHARP MARGINS LT  
VA-TRR ON FINDING DISC DRUSEN LT  
VA-TRR ON FINDING KNOWN GLAUCOMA LT  
VA-TRR ON FINDING PHYSIOLOGIC CUPPING SUSPECTED LT  
VA-TRR ON FINDING PERIPAPILLARY ATROPHY LT  
VA-TRR ON FINDING OTHER LT  
VA-TRR ON FINDING DISC NEOVASCULARIZATION LT  
VA-TRR ON FINDING LG CUP/DISC LT  
VA-TRR ON FINDING NO APP ABN LT  
VA-TRR ON FINDING INADQ ASSESS LT  
VA-TRR CATEGORY-NON DM FINDINGS RT EYE [C]  
VA-TRR MAC FINDING MICROANEURYSM LT  
VA-TRR MAC FINDING ATROPHY LT  
VA-TRR MAC FINDING MAC HOLE-LT  
VA-TRR MAC FINDING EXUDATES LT  
VA-TRR MAC FINDING INTRARETINAL HEMORRHAGE LT  
VA-TRR MAC FINDING DRUSEN LT  
VA-TRR MAC FINDING PIGMENTARY CHANGE LT  
VA-TRR MAC FINDING OTHER LT  
VA-TRR MAC FINDING NEOVASC MACULOPHTY LT  
VA-TRR MAC FINDING GEOG ATROPHY LT  
VA-TRR MAC FINDING PUCKERING LT  
VA-TRR MAC FINDING DRUSEN-LG LT  
VA-TRR MAC FINDING DRUSEN-SMALL LT  
VA-TRR MAC FINDING NO APPARENT ABN LT  
VA-TRR MAC FINDING INADEQ ASSESS LT  
VA-TRR OTH ASSESS FIND INADQ ASSESS RT  
VA-TRR OTH ASSESS FINDING INTRARET HEM RT  
VA-TRR OTH ASSESS FINDING FAMILY HX GLAUCOMA RT

VA-TRR OTH ASSESS FINDING PERIPH RET DEGEN RT  
VA-TRR OTH ASSESS FINDING COTTON WOOL SPOT RT  
VA-TRR OTH ASSESS FINDING RETINAL VASC CHANGE RT  
VA-TRR OTH ASSESS FINDING CHORIORETINAL SCAR RT  
VA-TRR OTH ASSESS FINDING LENTICULAR OPACIFICATION RT  
VA-TRR OTH ASSESS FINDING LENS IMPLANT RT  
VA-TRR ON FINDING ELEVATED IOP RT  
VA-TRR OTH ASSESS FINDING OTH RT  
VA-TRR OTH ASSESS FINDING CHOR NEVUS RT  
VA-TRR OTH ASSESS FINDING HTN RTNPTHY RT  
VA-TRR OTH ASSESS FINDING OCCLUSION RT  
VA-TRR OTH ASSESS FINDING RET EMBOLUS RT  
VA-TRR OTH ASSESS FINDING CATARACT RT  
VA-TRR OTH ASSESS FINDING NO ABN FIND RT  
VA-TRR ON FINDING PALLOR RT  
VA-TRR ON FINDING PAPILLEDEMA SUSPECTED RT  
VA-TRR ON FINDING DISC EDEMA SUSP RT  
VA-TRR ON FINDING DISC DRUSEN RT  
VA-TRR ON FINDING DISC W/ SHARP MARGINS RT  
VA-TRR ON FINDING KNOWN GLAUCOMA RT  
VA-TRR ON FINDING PHYSIOLOGIC CUPPING SUSPECTED RT  
VA-TRR ON FINDING PERIPAPILLARY ATROPHY RT  
VA-TRR ON FINDING OTHER RT  
VA-TRR ON FINDING DISC NEOVASCULARIZATION RT  
VA-TRR ON FINDING LG CUP/DISC RT  
VA-TRR ON FINDING NO APP ABN RT  
VA-TRR ON FINDING INADQ ASSESS RT  
VA-TRR MAC FINDING INADEQ ASSESS RT  
VA-TRR MAC FINDING INTRARETINAL HEMORRHAGE RT  
VA-TRR MAC FINDING ATROPHY RT  
VA-TRR MAC FINDING MAC HOLE-RT  
VA-TRR MAC FINDING EXUDATES RT  
VA-TRR MAC FINDING MICROANEURYSM RT  
VA-TRR MAC FINDING DRUSEN RT  
VA-TRR MAC FINDING PIGMENTARY CHANGE RT  
VA-TRR MAC FINDING OTHER RT  
VA-TRR MAC FINDING NEOVASC MACULOPHTHY RT  
VA-TRR MAC FINDING GEOG ATROPHY RT  
VA-TRR MAC FINDING PUCKERING RT  
VA-TRR MAC FINDING DRUSEN-LG RT  
VA-TRR MAC FINDING DRUSEN-SMALL RT  
VA-TRR MAC FINDING NO APPARENT ABN RT  
VA-TRR CATEGORY-DM EXAM FINDINGS RT EYE [C]  
VA-TRR DM RETINOP EXAM-NO DM MAC EDEMA RT

VA-TRR DM RETINOP EXAM-DM MAC EDEMA RT  
VA-TRR DM2 FINDING DR NO RETINOPATHY RT  
VA-TRR DM2 FINDING DR INADQ ASSESS RT  
VA-TRR DM2 FINDING DR OTHER RT  
VA-TRR DM2 FIND DR PREV RET LASER TX RT  
VA-TRR DM2 FINDING DR PDR MAC EDEMA RT  
VA-TRR DM2 FIND DR PDR MAC NO EDEMA RT  
VA-TRR DM2 FINDING DR SEV MAC EDEMA RT  
VA-TRR DM2 FIND DR SEV MAC NO EDEMA RT  
VA-TRR DM2 FINDING DR MOD MAC EDEMA RT  
VA-TRR DM2 FIND DR MOD MAC NO EDEMA RT  
VA-TRR DM2 FIND DR MLD MAC NO EDEMA RT  
VA-TRR DM2 FINDING DR MLD MAC EDEMA RT  
VA-TRR CATEGORY-DM EXAM FINDINGS LT EYE [C]  
VA-TRR DM2 FINDING DR INADQ ASSESS LT  
VA-TRR DM RETINOP EXAM-DM MAC EDEMA LT  
VA-TRR DM RETINOP EXAM-NO DM MAC EDEMA LT  
VA-TRR DM2 FINDING DR NO RETINOPATHY LT  
VA-TRR DM2 FINDING DR OTHER LT  
VA-TRR DM2 FIND DR PREV RET LASER TX LT  
VA-TRR DM2 FINDING DR PDR MAC EDEMA LT  
VA-TRR DM2 FIND DR PDR MAC NO EDEMA LT  
VA-TRR DM2 FINDING DR SEV MAC EDEMA LT  
VA-TRR DM2 FIND DR SEV MAC NO EDEMA LT  
VA-TRR DM2 FINDING DR MOD MAC EDEMA LT  
VA-TRR DM2 FIND DR MOD MAC NO EDEMA LT  
VA-TRR DM2 FIND DR MLD MAC NO EDEMA LT  
VA-TRR DM2 FINDING DR MLD MAC EDEMA LT  
VA-TRR CATEGORY-DM EYE EXAM [C]  
VA-TRR DIABETES TYPE 2  
VA-TRR DM1 FINDING DR INADQ ASSESS LT  
VA-TRR DM1 FINDING DR NO RETINOPATHY LT  
VA-TRR DM1 FINDING DR OTHER RT  
VA-TRR DM1 FIND DR PREV RET LASER TX LT  
VA-TRR DM1 FINDING DR PDR MAC EDEMA LT  
VA-TRR DM1 FIND DR PDR MAC NO EDEMA LT  
VA-TRR DM1 FINDING DR SEV MAC EDEMA LT  
VA-TRR DM1 FIND DR SEV MAC NO EDEMA LT  
VA-TRR DM1 FINDING DR MOD MAC EDEMA LT  
VA-TRR DM1 FIND DR MOD MAC NO EDEMA LT  
VA-TRR DM1 FIND DR MLD MAC NO EDEMA LT  
VA-TRR DM1 FINDING DR MLD MAC EDEMA LT  
VA-TRR DM1 FINDING DR NO RETINOPATHY RT  
VA-TRR DM1 FINDING DR INADQ ASSESS RT

VA-TRR DM1 FIND DR PREV RET LASER TX RT  
VA-TRR DM1 FINDING DR PDR MAC EDEMA RT  
VA-TRR DM1 FIND DR PDR MAC NO EDEMA RT  
VA-TRR DM1 FINDING DR SEV MAC EDEMA RT  
VA-TRR DM1 FIND DR SEV MAC NO EDEMA RT  
VA-TRR DM1 FINDING DR MOD MAC EDEMA RT  
VA-TRR DM1 FIND DR MOD MAC NO EDEMA RT  
VA-TRR DM1 FIND DR MLD MAC NO EDEMA RT  
VA-TRR DM1 FINDING DR MLD MAC EDEMA RT  
VA-TRR DIABETES TYPE 1  
VA-TRR CATEGORY-IMAGE QUALITY [C]  
VA-TRR IMAGE QUAL ADEQUATE  
VA-TRR IMAGE QUAL INADQ-PT POSITION  
VA-TRR IMAGE QUAL INADQ-OTHER REASON  
VA-TRR IMAGE QUAL INADQ-MEDIA  
VA-TRR IMAGE QUAL INADQ-PUPIL SIZE  
VA-TRR IMAGE QUAL INADQ-FIELD MISSING  
VA-TRR IMAGE QUAL INADQ-BLURRED IMAGE  
VA-TRR IMAGE QUAL INADEQUATE

#### **REMINDER SPONSOR**

OFFICE OF TELEHEALTH SERVICES STORE AND FORWARD TELEHEALTH  
OFFICE OF PATIENT CARE SERVICES

#### **REMINDER TAXONOMY**

VA-TELEDERM READER TIME SPENT < 5 MIN  
VA-TELEDERM READER TIME SPENT > 5 MIN  
VA-TRI IMAGER COMP EYE EXAM  
VA-TRI IMAGER CPT CODES  
VA-DIABETES HEDIS  
VA-TRI DM2 NO RETINOPATHY NO MACULAR EDEMA  
VA-TRI DM1 NO RETINOPATHY NO MACULAR EDEMA  
VA-DIABETES  
VA-GLAUCOMA RISK FACTORS  
VA-AMD RISK FACTORS  
VA-ABNORMAL BMI  
VA-CAD  
TOBACCO USE CODES  
VA-TRI READER TIME SPENT < 5 MIN  
VA-TRI READER TIME SPENT > 5 MIN  
VA-TRR MACULAR INTRARETINAL HEMORRHAGE LT  
VA-TRR OTH ASSESS CATARACT-VISUAL SIGN SUSP  
VA-TRR FAMILY HX OF GLAUCOMA  
VA-TRR PERIPHERAL RETINAL DEGENERATION

VA-TRR COTTON WOOL SPOT  
VA-TRR RETINAL VASCULAR CHANGE LT  
VA-TRR CHORIORETINAL SCAR LT  
VA-TRR LENTICULAR OPACIFICATION LT  
VA-TRR POSTERIOR IO LENS IMPLANT  
VA-TRI ELEVATED IOP L  
VA-TRI CHOROIDAL NEVUS L  
VA-TRI HYPERTENSIVE RETINOPATHY L  
VA-TRI BRANCH RETINAL VEIN OCCLUSION  
VA-TRI RETINAL EMBOLUS BRANCH RETINAL ARTERY OCCLUSION L  
VA-TRR ON PALLOR LT  
VA-TRR ON PAPILLEDEMA SUSPECT  
VA-TRR ON DISC EDEMA SUSPECT LT  
VA-TRR ON DISC DRUSEN LT  
VA-TRR ON KNOWN GLAUCOMA  
VA-TRR ON DISC NEOVASCULARIZATION LT  
VA-TRI GLAUCOMA SUSPECT L  
VA-TRR MACULAR MICROANEURYSM LT EYE  
VA-TRR MACULAR ATROPHY LT  
VA-TRR MACULAR HOLE LT  
VA-TRR MACULAR DRUSEN LT  
VA-TRI EVIDENCE OF NEOVASCULAR MACULOPATHY ADV EXUDATIVE AMD L  
VA-TRI GEOGRAPHIC ATROPHY INVOLVING FOVEA ADV AMD L  
VA-TRI MACULAR PUCKER L  
VA-TRI INT DRUSEN DRUSE ATROPHY NO FOVEA INT NONEXUDATIVE AMD L  
VA-TRI SMALL INT DRUSEN EARLY NONEXUDATIVE AMD L  
VA-TRR MACULAR INTRARETINAL HEMORRHAGE RT  
VA-TRR RETINAL VASCULAR CHANGE RT  
VA-TRR CHORIORETINAL SCAR RT  
VA-TRR LENTICULAR OPACIFICATION RT  
VA-TRI ELEVATED IOP R  
VA-TRI CHOROIDAL NEVUS R  
VA-TRI HYPERTENSIVE RETINOPATHY R  
VA-TRI RETINAL EMBOLUS BRANCH RETINAL ARTERY OCCLUSION R  
VA-TRR ON PALLOR RT  
VA-TRR ON DISC EDEMA SUSPECT RT  
VA-TRR ON DISC DRUSEN RT  
VA-TRR ON DISC NEOVASCULARIZATION RT  
VA-TRI GLAUCOMA SUSPECT R  
VA-TRR MACULAR ATROPHY RT  
VA-TRR MACULAR HOLE RT  
VA-TRR MACULAR MICROANEURYSM RT EYE  
VA-TRR MACULAR DRUSEN RT  
VA-TRI EVIDENCE OF NEOVASCULAR MACULOPATHY ADV EXUDATIVE AMD R

VA-TRI GEOGRAPHIC ATROPHYINVOLVING FOVEA ADV AMD R  
VA-TRI MACULAR PUCKER R  
VA-TRI INT DRUSEN DRUSE ATROPHY NO FOVEA INT NONEXUDATIVE AMD R  
VA-TRI SMALL INT DRUSEN EARLY NONEXUDATIVE AMD R  
TRI EX DM RETINOPATHY U/A  
VA-TRI DM1/DM2 NO RETINOPATHY NO MACULAR EDEMA R/L  
VA-TRI DM2 PROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM2 PROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM2 SEVERE NONPROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM2 SEVERE NONPROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM2 MOD NONPROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM2 MOD NONPROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM2 MILD NONPROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM2 MILD NONPROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM2 PROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM2 PROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM2 SEVERE NONPROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM2 SEVERE NONPROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM2 MOD NONPROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM2 MOD NONPROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM2 MILD NONPROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM2 MILD NONPROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM1 PROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM1 PROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM1 SEVERE NONPROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM1 SEVERE NONPROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM1 MOD NONPROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM1 MOD NONPROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM1 MILD NONPROLIF RETINOPATHY NO MACULAR EDEMA L  
VA-TRI DM1 MILD NONPROLIF RETINOPATHY PROB MACULAR EDEMA L  
VA-TRI DM1 PROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM1 PROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM1 SEVERE NONPROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM1 SEVERE NONPROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM1 MOD NONPROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRI DM1 MOD NONPROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM1 MILD NONPROLIF RETINOPATHY NO MACULAR EDEMA R  
VA-TRI DM1 MILD NONPROLIF RETINOPATHY PROB MACULAR EDEMA R  
VA-TRR IMAGE ADEQUATE ICD CODE  
VA-TRR IMAGE INADEQUATE ICD CODE

#### **REMINDER TERM**

VA-REMINDER UPDATE\_2\_0\_455  
VA-BL TELE EYE DIABETIC PT

VA-BL SEMAGLUTIDE RX  
VA-BL TELE EYE GLAUCOMA PT  
VA-AMD RISK FACTORS  
VA-BL TELE EYE AMD PT

#### **REMINDER DEFINITION**

VA-OBJ TELE-EYE MEDICATION-SEMAGLUTIDE  
VA-DIABETES OBJECT  
VA-BL TELE-EYE SCREEN-GLAUCOMA  
VA-OBJ GLAUCOMA RISK FACTORS  
VA-BL TELE-EYE SCREEN-AMD  
VA-OBJ AMD RISK FACTORS

#### **HEALTH SUMMARY TYPE**

TRI HGBA1C  
VA-OBJ SEMAGLUTIDE MED  
VA-DIABETES OBJECT  
VA-GLAUCOMA RISK FACTORS  
VA-AMD RISK FACTORS

#### **HEALTH SUMMARY OBJECTS**

TRI HGBA1C (TIU)  
VA-OBJ SEMAGLUTIDE MED  
VA-DIABETES OBJECT (TIU)  
VA-GLAUCOMA RISK FACTORS (TIU)  
VA-AMD RISK FACTORS (TIU)

#### **TIU DOCUMENT DEFINITION**

TRI HGBA1C  
VA-DM HIGH-RISK MEDS  
VA-DIABETES OBJECT  
VA-GLAUCOMA RISK FACTORS  
VA-AMD RISK FACTORS

#### **REMINDER DIALOG**

VA-TELEDERMATOLOGY READER NOTE  
VA-TELEDERMATOLOGY IMAGER NOTE  
VA-TELERETINOPATHY IMAGER NOTE  
VA-TELERETINOPATHY READER NOTE

## Install Details

This update is being distributed as a web host file. The address for the host file is:

[https://REDACTED/UPDATE\\_2\\_0\\_455\\_3.PRD](https://REDACTED/UPDATE_2_0_455_3.PRD)

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 5 minutes.

## Pre-Install

- Using the VistA menu Reminder Dialog Manager, Reminder Dialog, change your view to Element and do an inquire/print for the following reminder dialog elements. Identify orders mapped to the elements as additional findings. You will use this as your reference for post install. These orders will need to be re-mapped after installation.

**VAL-OI EYE CARE AT-RISK RETINAL IMAGE CONSULT**  
**VAL-OI EYE CARE AT-RISK EYE GLASS CONSULT**  
**VAL-OI EYE CARE AT-RISK COMP EYE EXAM CONSULT**

- Using the VistA menu Reminder Dialog Manager, Reminder Dialog, change your view to Element. and do an inquire/print for the following reminder dialog elements. Review the dialog text for the use of any local data objects, they will need to be added back after the install.

**VA-TXT TRI IMAGING DEMOGRAPHICS**  
**VAL-OBJ MOST RECENT EYE EXAM DATES**  
**VAL-OBJ FUTURE EYE EXAM DATES**



## Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

**https://REDACTED/UPDATE\_2\_0\_455\_3.PRD**

You should see a message at the top of your screen that the file successfully loaded.

Search for and locate the entry **UPDATE\_2\_0\_455 VA-ASYNCHRONOUS TELEHEALTH TEMPLATES.**

```
+Item  Entry                                     Source                                     Date Packed
717    UPDATE_2_0_455
       VA-ASYNCHRONOUS TELEHEALTH
       TEMPLATES
718    UPDATE_2_0_46 VA-MH TOBACCO
       USE SCREEN
719    UPDATE_2_0_47 VA-HEPATITIS B
       IMMUNIZATIONS
720    UPDATE_2_0_48 MISC UPDATE 6
721    UPDATE_2_0_49 VA-SUICIDE
       BEHAVIOR AND OVERDOSE REPORT
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History              RP   Repack
Select Action: Next Screen// IFE  Install Exchange File Entry
Enter a list or range of numbers (1-943): 717
```

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry titled **UPDATE\_2\_0\_455 VA-ASYNCHRONOUS TELEHEALTH TEMPLATES** (in this example it is entry 717

, it will vary by site, the date packed should be 07/12/2024)

Package Version: 2.0P65

Description:

The following Clinical Reminder items were selected for packing:

REMINDER DIALOG

VA-TELERETINOPATHY READER NOTE  
VA-TELERETINOPATHY IMAGER NOTE  
VA-TELEDERMATOLOGY IMAGER NOTE  
VA-TELEDERMATOLOGY READER NOTE

REMINDER TERM

VA-REMINDER UPDATE\_2\_0\_455

Non-exchangeable TIU object(s):

```
-----  
+          Enter ?? for more actions                                >>>  
IA  Install all Components          IS  Install Selected Component  
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.

Select Action: Next Screen// **IA Install all Components**

During the installation, you will see several prompts:

**Install** all new components

If you are prompted to OVERWRITE any existing components select **S for Skip**.

**You will be prompted to install the FIRST dialog component:**

```
Packed reminder dialog: VA-TELEDERMATOLOGY READER NOTE [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-TELEDERMATOLOGY READER NOTE	dialog	X
2	5	VA-TELEDERM READER NOTE VERSION Finding: *NONE*	element	X
3	10	VA-TDR RESIDENT/ATTENDING Finding: *NONE*	group	X
4	10.5	VA-TDR TXT READER RESIDENT/ATTENDING Finding: *NONE*	element	X
5	15	VA-TDR-DG HISTORY Finding: *NONE*	group	X
6	20	VA-TDR CONSULT/IMAGE QUALITY Finding: *NONE*	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-TELEDERMATOLOGY READER NOTE**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```
Packed reminder dialog: VA-TELEDERMATOLOGY READER NOTE [NATIONAL DIALOG]
VA-TELEDERMATOLOGY READER NOTE (reminder dialog) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-TELEDERMATOLOGY READER NOTE	dialog	X
2	5	VA-TELEDERM READER NOTE VERSION Finding: *NONE*	element	X
3	10	VA-TDR RESIDENT/ATTENDING Finding: *NONE*	group	X
4	10.5	VA-TDR TXT READER RESIDENT/ATTENDING Finding: *NONE*	element	X
5	15	VA-TDR-DG HISTORY Finding: *NONE*	group	X
6	20	VA-TDR CONSULT/IMAGE QUALITY Finding: *NONE*	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

**You will be prompted to install the NEXT dialog component:**

```

Packed reminder dialog: VA-TELEDERMATOLOGY IMAGER NOTE [NATIONAL DIALOG]
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-TELEDERMATOLOGY IMAGER NOTE                       dialog  X
  2    1 VA-TELEDERMATOLOGY IMAGER NOTE VERSION             element  X
      Finding: *NONE*
  3    5 VA-GP TELEDERMATOLGY IMAGER NOTE                   group
      Finding: *NONE*
  4    5.5 VA-TDI DG CONSULT REQUEST                        group    X
      Finding: *NONE*
  5    5.5.5 VA-TDI HF CONSENT/EDUCATION                    element  X
      Finding: VA-TDI EDUCATION (EDUCATION TOPIC)
      Add. Finding: VA-TDI EDUCATION DONE (HEALTH FACTOR)
  6    VA-TDI PXRМ COMMENT                                  prompt  X
  7    5.5.10 VA-TDI GP TELEDERM CONSENT                    group    X
      Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details    DT Dialog Text      IS Install Selected
DF Dialog Findings   DU Dialog Usage    QU Quit
DS Dialog Summary    IA Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-TELEDERMATOLOGY IMAGER NOTE**

Select Action: Next Screen// **IA** Install All

Install reminder dialog and all components with no further changes: **Y// ES**

```

Packed reminder dialog: VA-TELEDERMATOLOGY IMAGER NOTE [NATIONAL DIALOG]
VA-TELEDERMATOLOGY IMAGER NOTE (reminder dialog) installed from exchange file.
Item  Seq. Dialog Findings                                     Type  Exists
  1    VA-TELEDERMATOLOGY IMAGER NOTE                       dialog  X
  2    1 VA-TELEDERMATOLOGY IMAGER NOTE VERSION             element  X
      Finding: *NONE*
  3    5 VA-GP TELEDERMATOLGY IMAGER NOTE                   group    X
      Finding: *NONE*
  4    5.5 VA-TDI DG CONSULT REQUEST                        group    X
      Finding: *NONE*
  5    5.5.5 VA-TDI HF CONSENT/EDUCATION                    element  X
      Finding: VA-TDI EDUCATION (EDUCATION TOPIC)
      Add. Finding: VA-TDI EDUCATION DONE (HEALTH FACTOR)
  6    VA-TDI PXRМ COMMENT                                  prompt  X
  7    5.5.10 VA-TDI GP TELEDERM CONSENT                    group    X
      Finding: *NONE*
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details    DT Dialog Text      IS Install Selected
DF Dialog Findings   DU Dialog Usage    QU Quit
DS Dialog Summary    IA Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

**You will be prompted to install the NEXT dialog component:**

```

Packed reminder dialog: VA-TELERETINOPATHY IMAGER NOTE [NATIONAL DIALOG]
+-----+-----+-----+-----+-----+
Item  Seq. Dialog Findings                                     Type  Exists
+-----+-----+-----+-----+-----+
  1    VA-TELERETINOPATHY IMAGER NOTE                       dialog  X
  2    5 VA-TELE-EYE IMAGER NOTE VERSION                    element  X
       Finding: *NONE*
  3    10 VA-GP TRI IMAGER CONSULT                          group   X
       Finding: *NONE*
  4    10.5 VA-GP TRI IMAGING CONSULT INITIAL                group   X
       Finding: *NONE*
  5    10.5.5 VA-TXT TRI IMAGING DEMOGRAPHICS              element  X
       Finding: *NONE*
  6    10.5.10 VA-GP TRI IMAGER'S CONSULT INITIAL RESULT   group   X
       Finding: *NONE*
  7    10.5.10.5 VA-TXT TRI BLANK SPACE                    element  X
       Finding: *NONE*
+-----+-----+-----+-----+-----+
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA

```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-TELERETINOPATHY IMAGER NOTE**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```

Packed reminder dialog: VA-TELERETINOPATHY IMAGER NOTE [NATIONAL DIALOG]
VA-TELERETINOPATHY IMAGER NOTE (reminder dialog) installed from exchange file.
+-----+-----+-----+-----+-----+
Item  Seq. Dialog Findings                                     Type  Exists
+-----+-----+-----+-----+-----+
  1    VA-TELERETINOPATHY IMAGER NOTE                       dialog  X
  2    5 VA-TELE-EYE IMAGER NOTE VERSION                    element  X
       Finding: *NONE*
  3    10 VA-GP TRI IMAGER CONSULT                          group   X
       Finding: *NONE*
  4    10.5 VA-GP TRI IMAGING CONSULT INITIAL                group   X
       Finding: *NONE*
  5    10.5.5 VA-TXT TRI IMAGING DEMOGRAPHICS              element  X
       Finding: *NONE*
  6    10.5.10 VA-GP TRI IMAGER'S CONSULT INITIAL RESULT   group   X
       Finding: *NONE*
  7    10.5.10.5 VA-TXT TRI BLANK SPACE                    element  X
       Finding: *NONE*
+-----+-----+-----+-----+-----+
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

**You will be prompted to install the NEXT dialog component:**

```
Packed reminder dialog: VA-TELERETINOPATHY READER NOTE [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-TELERETINOPATHY READER NOTE	dialog	X
2	5	VA-TELE-EYE READER NOTE VERSION Finding: *NONE*	element	X
3	10	VA-GP TRR TELE-EYE CARE IMAGE QUALITY Finding: *NONE*	group	X
4	10.1	VA-GP TRR TELE-EYE EXAM QUAL ADEQUATE Finding: VA-TRR IMAGE ADEQUATE ICD CODE (REMINDER TAXONOMY) Add. Finding: VA-TRR IMAGE QUAL ADEQUATE (HEALTH FACTOR)	group	X X
5		FORCED PRIMARY DIAGNOSIS	forced	X
6	10.1.10	VA-TXT TRR TELE-EYE EXAM QUAL ADEQUATE-FILTERS Finding: *NONE*	element	X
7		PXRM COMMENT	prompt	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-TELERETINOPATHY READER NOTE**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```
Packed reminder dialog: VA-TELERETINOPATHY READER NOTE [NATIONAL DIALOG]
```

**VA-TELERETINOPATHY READER NOTE (reminder dialog) installed from exchange file.**

Item	Seq.	Dialog Findings	Type	Exists
1		VA-TELERETINOPATHY READER NOTE	dialog	X
2	5	VA-TELE-EYE READER NOTE VERSION Finding: *NONE*	element	X
3	10	VA-GP TRR TELE-EYE CARE IMAGE QUALITY Finding: *NONE*	group	X
4	10.1	VA-GP TRR TELE-EYE EXAM QUAL ADEQUATE Finding: VA-TRR IMAGE ADEQUATE ICD CODE (REMINDER TAXONOMY) Add. Finding: VA-TRR IMAGE QUAL ADEQUATE (HEALTH FACTOR)	group	X X
5		FORCED PRIMARY DIAGNOSIS	forced	X
6	10.1.10	VA-TXT TRR TELE-EYE EXAM QUAL ADEQUATE-FILTERS Finding: *NONE*	element	X
7		PXRM COMMENT	prompt	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```

Package Version: 2.0P65

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-TELERETINOPATHY READER NOTE
  VA-TELERETINOPATHY IMAGER NOTE
  VA-TELEDERMATOLOGY IMAGER NOTE
  VA-TELEDERMATOLOGY READER NOTE

REMINDER TERM
  VA-REMINDER UPDATE_2_0_455
Non-exchangeable TIU object(s):
-----
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
IA  Install all Components          IS  Install Selected Component
Select Action: Next Screen// Q

```

You will be returned to this screen. Type Q to quit.  
 Install Completed

## Post Installation

1. In CPRS, from the Notes tab, locate the progress notes linked to the following reminder dialogs and confirm the version number:
  - a. **VA-TELERETINOPATHY READER NOTE** (i.e., TELE-EYE IMAGING CONSULT), confirm the version is 7.2.
  - b. **VA- TELEDERMATOLOGY READER NOTE** (i.e., TELEDERM READER REPORT), confirm the version is 5.2.
  - c. **VA-TELERETINOPATHY IMAGER NOTE** (i.e., TELE-EYE SCREENING CONSULT), confirm the version is 7.2.
  - d. **VA- TELEDERMATOLOGY IMAGER NOTE** (i.e., TELEDERM IMAGING REQUEST CONSULT), confirm the version is 5.2.
  
2. Using the VistA menu Reminder Dialog Manager, Reminder Dialog, change your view to Element. Search for the following reminder dialog element and edit the dialog text to add local data objects:
  - a. **VA-TXT TRI IMAGING DEMOGRAPHICS**  
 Edit the dialog text. If your site does not have the data objects included in this element, replace the dialog text with the local patient data object(s) that your site uses to display [patient ethnicity, patient race, patient date of birth and patient sex].
  
  - b. **VAL-OBJ MOST RECENT EYE EXAM DATES**

- Edit the dialog text. Add your local data object for most recent eye exam.
- c. **VAL-OBJ FUTURE EYE EXAM DATES**  
Edit the dialog text. Add your local data object for future eye exam appointments.
3. Using the Vista Menu- Reminder Dialog Management, Reminder Dialogs, change view to Element view and search for the following reminder dialog elements to add the quick orders identified in the pre-install as an additional finding:
- a. **VAL-OI EYE CARE AT-RISK RETINAL IMAGE CONSULT**  
Add local quick order for Retinal Image Consult as an additional finding.
  - b. **VAL-OI EYE CARE AT-RISK EYE GLASS CONSULT**  
Add local quick order/process for scheduling an eye clinic appointment for glasses as an additional finding.
  - c. **VAL-OI EYE CARE AT-RISK COMP EYE EXAM CONSULT**  
Add local quick order/process for scheduling an eye clinic appointment as an additional finding.

*NOTE: If you require further technical assistance, if there are any questions or problems during the installation, please notify your local IT support to log a national help desk ticket or contact the VA Service Desk and have them submit a national ticket to:  
Category: Affected Service  
Affected Service: VistA -CPRS: Clinical Reminders*