



Update_2_0_270

Clinical Reminders

VA-MISC Update 17

Install Guide

June 2022

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION 1

PRE-INSTALLATION 4

INSTALL DETAILS 5

INSTALL EXAMPLE 5

POST INSTALLATION 22

Introduction

This update contains minor spelling corrections and other changes/edits to existing clinical reminder content, a full list of the changes can be found in the Document MISC_UPDATE_17_CHANGES in docx or pdf format in the VistA Documentation Library, <https://www.va.gov/vdl>.

UPDATE_2_0_270 contains 1 Reminder Exchange entry:
UPDATE_2_0_270 VA-MISC UPDATE 17

The exchange file contains the following components:

MH TESTS AND SURVEYS

PHQ-2

REMINDER DIALOG LINK TYPE

VA-DEPRESSION POSITIVE WARM HANDOFF CHECK

TIU TEMPLATE FIELD

VA-HCOE HA DURATION

VA-HCOE HA SYMPTOMS

YES / NO

VA-HCOE HA INTENSITY

VA-HCOE HA PAIN TYPE

VA-HCOE HA LOCATION

VA-HCOE HA FREQUENCY

TEXT (1-60 CHAR)

VA-HCOE HA CAUSE

NUM REQ

YES/NO/UNSURE

VA-HCOE CBT MOTOR

TEXT (1-10 CHAR)

NCP HERPES ZOSTER DETAILED GUIDANCE

CDC ZOSTER VACCINE RECOMMENDATIONS

NCP HERPES ZOSTER GUIDANCE

WOUND DRESSINGS

EDIT 20

YES NO NA

HEALTH FACTORS

VA-REMINDER UPDATES [C]

VA-UPDATE_2_0_270

VA-SCI/D BLADDER MGMT-SURG PROCEDURES [C]

VA SCI/D BLADDER SURG HX-OTHER

VA-SCI/D BLADDER SURG HX-SUPRAPUBIC CATH

VA-SCI/D BLADDER SURG HX-SPHINCTEROTOMY

VA-SCI/D BLADDER SURG HX-SACR ANT STIM
VA-SCI/D BLADDER SURG HX-ILEOVESICOSTOMY
VA-SCI/D BLADDER SURG HX-ILEAL CONDUIT
VA-SCI/D BLADDER SURG HX-CONTINENT STOMA
VA-SCI/D BLADDER SURG HX-AUGMENTATION
VA-SCI/D BLADDER SURG HX-ART SPHINCTER
VA-SCI/D BLADDER SURG HX-NONE
VA-SCI/D BLADDER MGMT-COLLECTION DEVICE [C]
VA-SCI/D BLADDER COLLECTION-FOLEY/LEG BAG
VA-SCI/D BLADDER COLLECTION-URINAL
VA-SCI/D BLADDER COLLECTION-OSTOMY BAG
VA-SCI/D BLADDER COLLECTION-DIAPER/PAD
VA-SCI/D BLADDER COLLECTION-CONDOM CATH
VA-SCI/D BLADDER COLLECTION-NONE
VA-SCI/D BLADDER MGMT-EMPTYING METHOD [C]
VA-SCI/D BLADDER-SACRAL ANT ROOT STIM
VA-SCI/D BLADDER EMPTY METHOD-UNKNOWN
VA-SCI/D BLADDER INTRMNT CATH-ATNDT
VA-SCI/D BLADDER INTRMNT CATH-SELF
VA-SCI/D BLADDER EMPTY METHOD-OTHER
VA-SCI/D BLADDER-DIALYSIS
VA-SCI/D BLADDER-UROSTOMY
VA-SCI/D BLADDER INDWELL CATH-SUPRAPUBIC
VA-SCI/D BLADDER INDWELL CATH-TRANSURETH
VA-SCI/D BLADDER REFLEX TRIG-EXT COMP
VA-SCI/D BLADDER REFLEX TRIG-STRAINING
VA-SCI/D BLADDER REFLEX TRIG-INVOLUNTARY
VA-SCI/D BLADDER REFLEX TRIG-VOLUNTARY
VA-SCI/D BLADDER-NORMAL VOIDING
VA-SCI/D REVIEW DATE [C]
VA-SCI/D BLADDER MGMT REVIEW DATE
VA-HEADACHE CENTER OF EXCELLENCE COGNITIVE BEH THERAPY
VA-HCOE CBT MOD-SEV DEP SYMP
VA-VAAES LINE [C]
VA-VAAES MIDLINE INSERT PREP DRAPE
VA-VAAES MIDLINE INSERT PREP MASK
VA-VAAES MIDLINE INSERT PREP GOWN
VA-VAAES MIDLINE INSERT PREP STERILE GLOVES
VA-VAAES MIDLINE INSERT PREP CAP
VA-VAAES MIDLINE INSERTION PREP OTHER
VA-VAAES MIDLINE INSERTION PREP ANTIMIC
VA-VAAES MIDLINE INSERTION PREP CHLOR
VA-VAAES MIDLINE INSERT PREP HAND HYGIENE
VA-VAAES CENTRAL LINE #2 INSERT PREP DRAPE
VA-VAAES CENTRAL LINE #2 INSERT PREP MASK
VA-VAAES CENTRAL LINE #2 INSERT PREP GOWN

VA-VAAES CENTRAL LINE #2 INSERT PREP STERILE GLOVES
VA-VAAES CENTRAL LINE #2 INSERT PREP CAP
VA-VAAES CL#2 INSERTION PREP OTHER
VA-VAAES CL#2 INSERTION PREP ANTIMIC
VA-VAAES CL#2 INSERTION PREP CHLOR
VA-VAAES CENTRAL LINE #2 INSERT PREP HAND HYGIENE
VA-VAAES CENTRAL LINE #1 INSERT PREP DRAPE
VA-VAAES CENTRAL LINE #1 INSERT PREP MASK
VA-VAAES CENTRAL LINE #1 INSERT PREP GOWN
VA-VAAES CENTRAL LINE #1 INSERT PREP STERILE GLOVES
VA-VAAES CENTRAL LINE #1 INSERT PREP CAP
VA-VAAES CL#1 INSERTION PREP OTHER
VA-VAAES CL#1 INSERTION PREP ANTIMIC
VA-VAAES CL#1 INSERTION PREP CHLOR
VA-VAAES CENTRAL LINE #1 INSERT PREP HAND HYGIENE
VA-VAAES PICC INSERT PREP DRAPE
VA-VAAES PICC INSERT PREP MASK
VA-VAAES PICC INSERT PREP GOWN
VA-VAAES PICC INSERT PREP STERILE GLOVES
VA-VAAES PICC INSERT PREP CAP
VA-VAAES PICC INSERTION PREP OTHER
VA-VAAES PICC INSERTION PREP ANTIMIC
VA-VAAES PICC INSERTION PREP CHLOR
VA-VAAES PICC INSERT PREP HAND HYGIENE

REMINDER SPONSOR

VHA NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE
PREVENTION (NCP) AND VHA OFFICE OF WOMEN'S HEALTH
OFFICE OF NURSING SERVICE
NATIONAL CENTER FOR HEALTH PROMOTION AND DISEASE PREVENTION
(NCP)
OFFICE OF NURSING INFORMATICS
MENTAL HEALTH SERVICES
OFFICE OF TELEHEALTH SERVICES STORE AND FORWARD TELEHEALTH

REMINDER TAXONOMY

VA-WH HYSTERECTOMY W/CERVIX REMOVED
VA-WH PAP SMEAR SCREEN CODES

REMINDER TERM

VA-REMINDER UPDATE_2_0_270
VA-SCI/D BLADDER MGMT-SURG HX (OBJ)
VA-SCI/D BLADDER MGMT-COLLECTION (OBJ)
VA-SCI/D BLADDER MGMT-EMPTYING METH (OBJ)
VA-SCI/D BLADDER MANAGEMENT STATUS

REMINDER DEFINITION

VA-OBJ SCI/D BLADDER MGMT-SURG HX
VA-OBJ SCI/D BLADDER MGMT-COLLECTION DEVICE
VA-OBJ SCI/D BLADDER MGMT-EMPTYING METHOD

HEALTH SUMMARY COMPONENT

CLINICAL REMINDER FINDINGS
PCE HEALTH FACTORS SELECTED

HEALTH SUMMARY TYPE

VA-SCI/D BLADDER MGMT-SURG HX
VA-SCI/D BLADDER-COLL DEVICE
VA-SCI/D BLADDER MGT-EMP METH
VA-SCI/D BLADDER MGMT REV DATE

HEALTH SUMMARY OBJECTS

VA-SCI/D BLADDER MGMT-SURG HX (TIU)
VA-SCI/D BLADDER MGMT-COLLECTION DEVICE (TIU)
VA-SCI/D BLADDER MGMT-EMPTYING METHOD (TIU)
VA-SCI/D BLADDER MGMT REV DATE

TIU DOCUMENT DEFINITION

VA-SCI/D BLADDER MGMT-SURG HX
VA-SCI/D BLADDER MGMT-COLLECTION DEVICE
VA-SCI/D BLADDER MGMT-EMPTYING METHOD
VA-SCI/D BLADDER MGMT REV DATE

REMINDER DIALOG

VA-GP MISC UPDATE 17

Pre-Installation

Prior to installing this Update:

1. Check Reminder Exchange Install History (IH) and ensure your facility has installed the following (previously released) Updates:
 - UPDATE_2_0_192 VA-MH UPDATES 8.0
 - UPDATE_2_0_212 VA-AES SHIFT ASSESSMENT BUNDLE
 - UPDATE_2_0_213 VA-SCI/D EMERGENCY DIALOGS
 - UPDATE_2_0_241 VA-HCOE TEMPLATES
 - UPDATE_2_0_232 VA-TELEWOUND IMAGER AND READER TEMPLATES
 - UPDATE_2_0_262 VA-HERPES ZOSTER REMINDER UPDATE

Install Details

This update is being distributed as a web host file. The address for the host file is:

https://REDACTED/UPDATE_2_0_270.PRD

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====
This update can be loaded with users on the system. Installation will take less than 15 minutes.

Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LWH  Load Web Host File
DFE  Delete Exchange File Entry       RI   Reminder Definition Inquiry
IFE  Install Exchange File Entry       RP   Repack
IH   Installation History
Select Action: Next Screen// lwh  Load Web Host File ←
Input the URL for the .prd file: https://vaww.va.gov/reminders/docs/UPDATE_2_0_270.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

https://REDACTED/UPDATE_2_0_270.PRD

You should see a message at the top of your screen that the file successfully loaded.

https://vaww.va.gov/reminders/docs/UPDATE_2_0_270.PRD successfully loaded.

Search for and locate the reminder exchange entry titled UPDATE_2_0_270 VA-MISC UPDATE 17 (see screenshot below).

+Item	Entry	Source	Date Packed
569	UPDATE_2_0_270 VA-MISC UPDATE 17		05/03/2022@09:01
570	UPDATE_2_0_271 VA-PNEUMOCOCCAL UPDATE		03/09/2022@12:28
571	UPDATE_2_0_271 VA-PNEUMOCOCCAL UPDATE		03/21/2022@14:26
572	UPDATE_2_0_272 VA-ATHLETES MED INFO TEMPLATE UPDATE		02/22/2022@08:37
573	UPDATE_2_0_273 VA-PHP PERSONAL HEALTH PLAN UPDATE		03/15/2022@04:12

+	+ Next Screen	- Prev Screen	?? More Actions	>>>
CFE	Create Exchange File Entry	LHF	Load Host File	
CHF	Create Host File	LMM	Load MailMan Message	
CMM	Create MailMan Message	LWH	Load Web Host File	
DFE	Delete Exchange File Entry	RI	Reminder Definition Inquiry	
IFE	Install Exchange File Entry	RP	Repack	
IH	Installation History			

Select Action: Next Screen// IFE Install Exchange File Entry
Enter a list or range of numbers (1-931): 569

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry
Enter the number that corresponds with your entry titled **UPDATE_2_0_270 VA-MISC UPDATE 17** (in this example it is entry 569, it will vary by site). The exchange file date should be 05/03/2022@09:01.

Component	Category	Exists
Source: [REDACTED]		
Date Packed: 05/03/2022@09:01:27		
Package Version: 2.0P71		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-GP MISC UPDATE 17		
REMINDER TAXONOMY		
VA-WH HYSTERECTOMY W/CERVIX REMOVED		
VA-WH PAP SMEAR SCREEN CODES		
REMINDER TERM		
VA-REMINDER UPDATE_2_0_270		

+	Enter ?? for more actions		>>>
IA	Install all Components	IS	Install Selected Component

Select Action: Next Screen// **IA**

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.
Select Action: Next Screen// **IA Install all Components**

You will see several prompts:

1 - For all new components, you will choose **I to Install**

2 - For some components, you will choose **O to Overwrite.**

3 - For some components, you will choose **S to Skip**

Each individual component you may be prompted for is listed below, along with the response you should provide

TIU TEMPLATE FIELD entry named **VA-HCOE HA CAUSE** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **verwrite the current entry**

Are you sure you want to overwrite? N// y YES..

TIU TEMPLATE FIELD entry named **YES/NO/UNSURE** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

TIU TEMPLATE FIELD entry named **VA-HCOE CBT MOTOR** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **verwrite the current entry**
Are you sure you want to overwrite? N// y YES.

TIU TEMPLATE FIELD entry named **NCP HERPES ZOSTER DETAILED GUIDANCE** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **verwrite the current entry**
Are you sure you want to overwrite? N// y YES..

TIU TEMPLATE FIELD entry named **WOUND DRESSINGS** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: **O// verwrite the current entry**
Are you sure you want to overwrite? N// y YES...

TIU TEMPLATE FIELD entry named **YES NO NA** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

HEALTH FACTORS entry **VA-UPDATE_2_0_270** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall.....**

HEALTH FACTORS entry **VA-SCI/D BLADDER COLLECTION-FOLEY/LEG BAG** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-SCI/D BLADDER COLLECTION-URINAL** is NEW,
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall.

HEALTH FACTORS entry named **VA-SCI/D BLADDER COLLECTION-DIAPER/PAD**
already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// verwrite the current entry
Are you sure you want to overwrite? N// y YES.....

HEALTH FACTORS entry **VA-VAAES MIDLINE INSERT PREP DRAPE** is NEW,
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall

HEALTH FACTORS entry **VA-VAAES MIDLINE INSERT PREP MASK** is NEW,
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES MIDLINE INSERT PREP GOWN** is NEW,
what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES MIDLINE INSERT PREP STERILE GLOVES** is
NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES MIDLINE INSERT PREP CAP** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall...**

HEALTH FACTORS entry **VA-VAAES MIDLINE INSERT PREP HAND HYGIENE** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #2 INSERT PREP DRAPE** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #2 INSERT PREP MASK** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #2 INSERT PREP GOWN** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #2 INSERT PREP STERILE GLOVES** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #2 INSERT PREP CAP** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall...**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #2 INSERT PREP HAND HYGIENE** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #1 INSERT PREP DRAPE** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #1 INSERT PREP MASK** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #1 INSERT PREP GOWN** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #1 INSERT PREP STERILE GLOVES** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall**

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #1 INSERT PREP CAP** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall...

HEALTH FACTORS entry **VA-VAAES CENTRAL LINE #1 INSERT PREP HAND HYGIENE** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: I// nstall...

REMINDER TAXONOMY entry named **VA-WH PAP SMEAR SCREEN CODES** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// verwrite the current entry

Are you sure you want to overwrite? N// y YES

REMINDER TAXONOMY entry named **VA-WH HYSTERECTOMY W/CERVIX REMOVED** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **verwrite the current entry**
Are you sure you want to overwrite? N// y YES

REMINDER TERM entry **VA-REMINDER UPDATE_2_0_270** is NEW, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: **I// nstall...**

REMINDER TERM entry **VA-SCI/D BLADDER MGMT-COLLECTION (OBJ)** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- I Install
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

REMINDER DEFINITION entry named **VA-OBJ SCI/D BLADDER MGMT-SURG HX** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

REMINDER DEFINITION entry named **VA-OBJ SCI/D BLADDER MGMT-EMPTYING METHOD** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

HEALTH SUMMARY TYPE entry named **VA-SCI/D BLADDER MGMT-SURG HX** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

HEALTH SUMMARY TYPE entry named **VA-SCI/D BLADDER-COLL DEVICE** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// s Skip, do not install this entry

HEALTH SUMMARY TYPE entry named **VA-SCI/D BLADDER MGT-EMP METH** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// s Skip, do not install this entry

HEALTH SUMMARY TYPE entry named **VA-SCI/D BLADDER MGMT REV DATE** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// s Skip, do not install this entry

HEALTH SUMMARY OBJECTS entry named **VA-SCI/D BLADDER MGMT-SURG HX (TIU)** already exists but the packed component is different, what do you want to do?

Select one of the following:

- C Create a new entry by copying to a new name
- O Overwrite the current entry
- U Update
- Q Quit the install
- S Skip, do not install this entry

Enter response: O// **s Skip, do not install this entry**

You will be prompted to install the dialog component:

```
Packed reminder dialog: VA-GP MISC UPDATE 17 [NATIONAL DIALOG]
+-----+-----+-----+-----+-----+
Item  Seq.  Dialog Findings                                     Type  Exists
-----+-----+-----+-----+-----+
  1    VA-GP MISC UPDATE 17
      Finding: *NONE*                                group
  2    10 VA-TXT TWI NURSE WOUND-DRESSING INFO
      Finding: *NONE*                                element  X
  3    20 VA-HERPES ZOSTER GUIDANCE DETAILS
      Finding: *NONE*                                element  X
  4    30 VA-PHQ-2
      Finding: PHQ-2 (MH TESTS AND SURVEY)           element  X
  5    PXRМ PHQ2 RESULT GROUP                          rs.group  X
  6    PXRМ PHQ2 RESULT ELEMENT 1                     rs.element X
  7    PXRМ PHQ2 RESULT ELEMENT 2                     rs.element X
+-----+-----+-----+-----+-----+
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-GP MISC UPDATE 17**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: **Y// ES**

```

Packed reminder dialog: VA-GP MISC UPDATE 17 [NATIONAL DIALOG]
VA-GP MISC UPDATE 17 (group) installed from exchange file.
Item  Seq. Dialog Findings                                     Type      Exists
-----
  1    VA-GP MISC UPDATE 17                                  group      X
      Finding: *NONE*
  2    10 VA-TXT TWI NURSE WOUND-DRESSING INFO              element     X
      Finding: *NONE*
  3    20 VA-HERPES ZOSTER GUIDANCE DETAILS                 element     X
      Finding: *NONE*
  4    30 VA-PHQ-2                                           element     X
      Finding: PHQ-2 (MH TESTS AND SURVEY)
  5    PXRМ PHQ2 RESULT GROUP                               rs.group    X
  6    PXRМ PHQ2 RESULT ELEMENT 1                          rs.element  X
  7    PXRМ PHQ2 RESULT ELEMENT 2                          rs.element  X

+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage     QU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q

```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

```

Component                                     Category      Exists
-----
Source:
Date Packed: 05/03/2022@09:01:27
Package Version: 2.0P71

Description:
The following Clinical Reminder items were selected for packing:
REMINDER DIALOG
  VA-GP MISC UPDATE 17

REMINDER TAXONOMY
  VA-WH HYSTERECTOMY W/CERVIX REMOVED
  VA-WH PAP SMEAR SCREEN CODES

REMINDER TERM
  VA-REMINDER UPDATE_2_0_270

+      + Next Screen  - Prev Screen  ?? More Actions      >>>
IA  Install all Components      IS  Install Selected Component
Select Action: Next Screen// Q

```

You will be returned to this screen. Type Q to quit.
Install Completed

Post Installation

No post-installation steps are needed.