



**Update \_2\_0\_119A**

**Clinical Reminders**

**VA-PTSD SCREENING UPDATE  
Install Guide**

**January 2021**

Product Development  
Office of Information Technology  
Department of Veterans Affairs

## **Table of Contents**

<b>INTRODUCTION .....</b>	<b>1</b>
<b>INSTALL DETAILS .....</b>	<b>1</b>
<b>INSTALL EXAMPLE.....</b>	<b>3</b>
<b>POST INSTALLATION.....</b>	<b>6</b>

# Introduction

---

## DESCRIPTION:

With the release of reminder update 119, the instrument used in screening for PTSD was changed from PC-PTSD-5+I9 to PC-PTSD-5. Along with this, the result group was also changed from PXRm PC-PTSD-5+I9 RESULT GROUP to PXRm PTSD RESULT GROUP. Within this new group, the positive sequence Result Element had a condition of I SCORE>2 which implies that a score of 3 is positive. This was also transferred to the corresponding progress note with verbiage that a score of 3 was positive. In MHA, the result displayed as the score of 3 being negative which was inconsistent. In Accordance to VA guidelines, a PTSD screening score of 4 or more constitutes a positive screen. In this update, we have created new Result Elements and a Result Group and added this new result group to dialog element VA-MH PC-PTSD+I9 SCREEN. This update corrects the scoring and progress note text. This is an emergency install.

UPDATE\_2\_0\_119A contains 1 Reminder Exchange entry:

### **UPDATE\_2\_0\_119A VA-PTSD SCREENING UPDATE**

**The exchange file contains the following components:**

#### **MH TESTS AND SURVEYS**

PC-PTSD-5

#### **HEALTH FACTORS**

VA-REMINDER UPDATES

VA-UPDATE\_2\_0\_119A

#### **REMINDER SPONSOR**

OFFICE OF MH AND SUICIDE PREVENTION

OFFICE OF MENTAL HEALTH SERVICES

#### **REMINDER TERM**

VA-REMINDER UPDATE\_2\_0\_119A

#### **REMINDER DIALOG**

VA-MH PC-PTSD+I9 SCREEN

# Install Details

---

This update is being distributed as a web host file. The address for the host file is:  
**[https://Redacted /UPDATE\\_2\\_0\\_119A.PRD](https://Redacted /UPDATE_2_0_119A.PRD)**

The file will be installed using Reminder Exchange, programmer access is not required.

## Installation:

=====

This update can be loaded with users on the system. Installation will take less than 10 minutes.

## Install Example

To Load the Web Host File. Navigate to Reminder exchange in Vista

```
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History             RP   Repack
Select Action: Next Screen// LWH  Load Web Host File
Input the URL for the .prd file: https://[REDACTED]'UPDATE_2_0_119A.PRD
```

At the **Select Action:** prompt, enter **LWH** for Load Web Host File

At the **Input the url for the .prd file:** prompt, type the following web address:

[https://Redacted/UPDATE\\_2\\_0\\_119A.PRD](https://Redacted/UPDATE_2_0_119A.PRD)

You should see a message at the top of your screen that the file successfully loaded.

https://[REDACTED]'UPDATE\_2\_0\_119A.PRD successfully loaded.

Search and locate an entry titled **UPDATE\_2\_0\_119A VA-PTSD SCREENING UPDATE** in reminder exchange.

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry UPDATE\_2\_0\_119A VA-PTSD SCREENING UPDATE (*in this example it is entry 143 it will vary by site*). The date of the exchange file should be 01/15/2021.

```
+Item  Entry                                     Source      Date Packed
143  UPDATE_2_0_119A VA-PTSD SCREENING UPDATE  [REDACTED]  01/15/2021@07:46
144  UPDATE_2_0_12 VA-DIFFICULT AIRWAY        [REDACTED]  03/10/2016@14:12
145  UPDATE_2_0_120 VA-OSP SUICIDE RISK MANAGEMENT FOLLOW-UP  [REDACTED]  03/18/2020@13:47
146  UPDATE_2_0_120 VA-OSP SUICIDE RISK MANAGEMENT FOLLOW-UP  [REDACTED]  03/19/2020@12:34
147  UPDATE_2_0_120 VA-OSP SUICIDE RISK MANAGEMENT FOLLOW-UP  [REDACTED]  03/20/2020@08:21
+      + Next Screen  - Prev Screen  ?? More Actions  >>>
CFE  Create Exchange File Entry      LHF  Load Host File
CHF  Create Host File                 LMM  Load MailMan Message
CMM  Create MailMan Message           LR   List Reminder Definitions
DFE  Delete Exchange File Entry       LWH  Load Web Host File
IFE  Install Exchange File Entry      RI   Reminder Definition Inquiry
IH   Installation History             RP   Repack
Select Action: Next Screen// IFE  Install Exchange File Entry
Enter a list or range of numbers (1-489): 143
```

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.  
Select Action: Next Screen// **IA Install all Components**

For items that do not exist, use INSTALL action

## You will be prompted to install the reminder dialog component

Item	Seq.	Dialog Findings	Type	Exists
1		VA-MH PC-PTSD+I9 SCREEN	element	X
		Finding: PC-PTSD-5 (MH TESTS AND SURVEY)		X
2		PXRM PC-PTSD-5 RESULT GROUP	rs.group	
3		PXRM PC-PTSD-5 RESULT ELEMENT 1	rs.element	
4		PXRM PC-PTSD-5 RESULT ELEMENT 2	rs.element	

  

+ Next Screen    - Prev Screen    ?? More Actions					
DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Quit// IA

At the **Select Action** prompt, type **IA** to install the dialog – **VA-MH PC-PTSD+I9 SCREEN**

Select Action: Next Screen// **IA Install All**

Install reminder dialog and all components with no further changes: Y// **Yes**

PC-PTSD-5 (element) installed from exchange file.

Item	Seq.	Dialog Findings	Type	Exists
1		VA-MH PC-PTSD+I9 SCREEN	element	X
		Finding: PC-PTSD-5 (MH TESTS AND SURVEY)		X
2		PXRM PC-PTSD-5 RESULT GROUP	rs.group	X
3		PXRM PC-PTSD-5 RESULT ELEMENT 1	rs.element	X
4		PXRM PC-PTSD-5 RESULT ELEMENT 2	rs.element	X

  

+ Next Screen   - Prev Screen   ?? More Actions					
DD	Dialog Details	DT	Dialog Text	IS	Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU	Quit
DS	Dialog Summary	IA	Install All		

Select Action: Quit// Q

After completing this dialog install, you will type Q.

You will then be returned to this screen

Component	Category	Exists
Source:		
Date Packed: 01/15/2021@07:46:59		
Package Version: 2.0P46		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-MH PC-PTSD+I9 SCREEN		
REMINDER TERM		
VA-REMINDER UPDATE_2_0_119A		
Keywords:		
Components:		
+   + Next Screen   - Prev Screen   ?? More Actions   >>>		
IA	Install all Components	IS   Install Selected Component

Select Action: Next Screen// Q

At the **Select Action** prompt, type **Q**.

Install complete.

## Post Installation

---

1. Navigate to dialog elements. Use INQ option and select dialog element VA-MH PC-PTSD+I9 SCREEN and print this to screen. Ensure that the Result Group entry has changed to PXRM PC-PTSD-5 RESULT GROUP

```
REMINDER DIALOG INQUIRY                               Jan 15, 2021 9:25:17 am Page 2
-----
Exclude from P/N: NO
Resolution type:  DONE AT ENCOUNTER
Finding Item:    MH PC-PTSD-5
Additional Findings: none
Vital Prompt Caption:
Orderable item:  N/A
Suppress Checkbox:
Result Group:
                  : PXRM PC-PTSD-5 RESULT GROUP
Exclude Mental Health
Test Text from P/N: NO
```