



Clinical Reminders

ICD-10 Follow-up Installation Guide

PXRM*2.0*47

GMTS*2.7*113

December 2017

**Department of Veterans Affairs
Office of Information and Technology (OIT)
Product Development**

(this page intentionally left blank)

Table of Contents

Introduction.....	1
Background.....	1
Patches in this project build	1
Related Documentation	3
Related Web Sites.....	3
Pre-Installation.....	4
1. Required Software	4
Installation.....	4
1. Retrieve the host file containing the multi-package build.....	4
2. Install the patch first in a training or test account.....	5
3. Backout preparation; these steps are taken in case the patch needs to be backed out.....	5
4. Load the distribution.	5
5. Backup a Transport Global	6
6. Compare Transport Global to Current System.....	6
7. Verify Checksums in Transport Global.....	7
8. Install the build.....	7
Post-Install Instructions	9
Acronyms	9
Appendix A: Installation Example	11
Appendix B: Post-Install Checksums	19
Appendix C: Install File Print Example.....	22
Appendix D: Build File Print Example	24

Introduction

Background

This patch is the follow-up to patch PXR*2*26. It removes the options, protocols, routines, and taxonomy fields made obsolete by PXR*2*26, adds the taxonomy editor actions that could not be completed in time for the release of PXR*2*26, and fixes some defects.

Testing revealed that the ScreenMan variable DDSCHANG can sometimes falsely indicate that a taxonomy was changed by the user during an editing session. To create a reliable indicator, the SHA-256 hash is computed when a taxonomy is selected for editing and again when the editing session is finished. If there is any difference in the values, then the taxonomy has changed. The SHA-256 hash is computed using the Kernel SHA hash utility distributed in XU*8.0*657; consequently XU*8.0*657 is required by PXR*2.0*47.

PXR*2.0*47 and GMTS*2.7*113 are distributed in a multi-package build named: CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0.

Patches in this project build

PXR*2.0*47

This build is the follow-up to PXR*2.0*26. It deletes the Expanded Taxonomies file #811.3, removes the options, protocols, routines, and taxonomy fields made obsolete by PXR*2.0*26, adds the taxonomy editor actions that could not be completed for the release of PXR*2.0*26, and fixes some defects.

The following items are deleted:

PRINT TEMPLATES:

- PXR TAXONOMY DIALOG FILE #811.2
- PXR TAXONOMY INQUIRY FILE #811.2
- PXR TAXONOMY LIST FILE #811.2
- PXR TAXONOMY LIST HEADER FILE #811.2

ROUTINES:

- PXRMBXTL
- PXRMCS
- PXRMCS
- PXRMC
- PXRMTAXS
- PXRMTDUP
- PXRMT
- PXRMTEDT
- PXRMVAL
- PXRMVALC
- PXRMVALU

OPTIONS:

PXRM TAX SELECTED CODES REPAIR
PXRM TAXONOMY CODE SEARCH
PXRM TAXONOMY COPY
PXRM TAXONOMY EDIT
PXRM TAXONOMY EXPANSION
PXRM TAXONOMY EXPANSION ALL
PXRM TAXONOMY EXPANSION VERIFY
PXRM TAXONOMY INQUIRY
PXRM TAXONOMY LIST
PXRM TAXONOMY MANAGEMENT (OLD)

PROTOCOLS:

PXRM DIALOG ADD
PXRM LEXICON VIEW
PXRM TAXONOMY OLD INQUIRE

The National Library of Medicine has a Value Set Authority Center (VSAC) web site, where value sets can be obtained. From the web site: “Value sets are lists of specific values (terms and their codes) derived from single or multiple standard vocabularies used to define clinical concepts (e.g. patients with diabetes, clinical visit, reportable diseases) used in quality measures and to support effective health information exchange.” These value sets cover many clinical areas of relevance to the VA and because they are very similar to taxonomies, they can be used to automatically generate taxonomies. Some value sets contain coding systems that are not supported in taxonomies; if a taxonomy is generated from one of these value sets, only the codes from supported coding systems will be included.

The Value Set functionality requires three new files: NLM QUALITY MEASURE GROUPS, NLM VALUE SET CODING SYSTEMS, and NLM VALUE SETS. These files are populated with the VSAC’s most recent release: January 2017. The following new options were created and added to the PXRM MANAGERS MENU.

1 NLM CLINICAL QUALITY MEASURES PXRM NLM CQM MENU NLM
Clinical Qual

ity Measures Menu

2 NLM VALUE SET MENU PXRM NLM VALUE SET MENU NLM Value Set
Menu

GMTS*2.7*113

One of the defects that is fixed involves display of the status line when the frequency is in hours. This change adds time to the Date Due and Last Done fields. Because there are a number of Clinical Reminders Health Summary components, a Health Summary patch is required, it is: GMTS*2.7*113.

Related Documentation

Documentation	Documentation File name
Installation Guide	PXRM_2_0_47_IG.PDF
Manager's Manual	PXRM_2_0_MM.PDF

NOTE: In this document, you will see references to both PXRM*2*47 and PXRM*2.0*47. The difference is that PXRM*2*47 is the name of the patch and PXRM*2.0*47 is the name of the build.

Related Web Sites

Site	URL	Description
National Clinical Reminders site	http://vista.med.va.gov/reminders	Contains manuals, PowerPoint presentations, and other information about Clinical Reminders
National Clinical Reminders Committee	http://vaww.portal.va.gov/sites/ncrcpublic/default.aspx	This committee directs the development of new and revised national reminders
Vista Document Library	http://www.va.gov/vdl/	Contains manuals for Clinical Reminders and related applications.
Health Information Management (HIM) ICD-10 Implementation	https://vaww.vha.vaco.portal.va.gov/sites/HDI/HIM/vaco_HIM/SitePages/Home.aspx	Contains general information, training resources, and other tools for VA's transition to ICD-10.

Pre-Installation

Previously, the computed finding VA-REMINDER DEFINITION would return false if the Computed Finding Parameter was set to a reminder definition that did not exist. The most likely reason for the definition not to exist is because it was renamed and the Computed Finding Parameter was not updated accordingly. Returning false was somewhat misleading because the reminder was never actually evaluated. To prevent this, the computed finding was changed so if the reminder definition does not exist it sets a fatal error and sends an error message to the Clinical Reminders mail group.

Before installing, it would be a good idea to run the Finding Usage Report and get a list of definitions and terms that use VA-REMINDER DEFINITION as a finding. For each instance, make sure that the Computed Finding Parameter is set to a valid reminder definition. This will prevent getting fatal errors because of incorrectly configured Computed Finding Parameters.

1. Required Software

*PXRM*2.0*47*

Package/Patch	Namespace	Version	Comments
Clinical Reminders PXRM*2*26	PXRM	2	Fully patched
Kernel XU*8.0*657	XU	8.0	Fully patched
VA FileMan	DI	22.2	Fully patched

*GMTS*2.7*113*

Package/Patch	Namespace	Version	Comments
Health Summary GMTS*2.7*89 PXRM*2.0*47	GMTS	2.7	

Installation

This section describes how to install the multi-package build that includes builds , PXRM*2.0*47 and GMTS*2.7*113.

This patch can be loaded with users on the system, but it should be done during off-hours. Estimated installation time is less than 5 minutes.

The installation needs to be done by a person with DUZ(0) set to "@."

1. Retrieve the host file containing the multi-package build.

Use sftp to access the build (the name of the host file is CR_ICD-10_FOLLOWUP .KID) from one of the following locations (with the ASCII file type):

Hines	ftp.fo-hines.med.va.gov	<sftp ftp.fo-hines.med.va.gov>
Salt Lake City	ftp.fo-slc.med.va.gov	<sftp ftp.fo-slc.med.va.gov>

2. Install the patch first in a training or test account.

Installing in a non-production environment will give you time to get familiar with new functionality.

3. Backout preparation; these steps are taken in case the patch needs to be backed out.

- Backup files ^PXD(811.2), ^PXRMD(811.4), and ^PXD(811.9) using your site's policy for backing up data.
- If the steps are unknown, here is a way it can be done, using ^PXD(811.2) as an example:
 - Go to a command prompt.
 - At the prompt, enter D GOGEN^%ZSPECIAL.
 - At the device prompt, enter the directory and file name where the global backup is to be stored.
 - At the Parameters? Prompt, press <enter>.
 - At the Global prompt, enter ^PXD(811.2).
 - Verify that the file was created and exists in the directory specified.

Example

```
DEV5A4:DVFDDEV>D GOGEN^%ZSPECIAL

Device: VA5$:[Local Directory]811_2_BACKUP.GBL
Parameters? ("WNS") =>
Warning: Use a "V" format to avoid problems with control characters.
Global ^PXD(811.2,
Global ^
```

4. Load the distribution.

In programmer mode, type, D ^XUP, select the Kernel Installation & Distribution System menu (XPD INSTALLATION MENU), and then the option LOAD a Distribution. Enter the directory name where you placed the host file followed by CR_ICD-10_FOLLOWUP.KID at the Host File prompt.

Example

```
Select Installation Option: LOAD a Distribution
```



```
Enter a Host File: VA5$:[DOWNLOADS]CR_ICD-10_FOLLOWUP.KID
KIDS Distribution saved on Apr 24, 2017@08:21:25
Comment: CLINICAL REMINDERS ICD-10 FOLLOWUP
```

From the Installation menu, you may elect to use the following options:

5. Backup a Transport Global

Use the KIDS Installation option, **Backup a Transport Global [XPD BACKUP]**. This option creates a MailMan message that will backup all current routines on your VistA/M system that will be replaced by the builds in this transport global. (If you need to preserve components that are not routines, you must back them up separately.)

At the **Select INSTALL NAME:** prompt, enter **Clinical Reminders ICD-10 Followup 1.0**

Example:

```
Select Installation <TEST ACCOUNT> Option: Backup a Transport Global
Select INSTALL NAME: CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 7/11/16@14:17:39
=> PXRM*2.0*47, GMTS*2.7*113 ;Created on Jul 05, 2016@11:14:
```

This Distribution was loaded on Apr 24, 2017@13:31:45with header of

CLINICAL REMINDERS ICD-10 FOLLOWUP ;Created on Apr 24, 2017@08:21:25

It consisted of the following Install(s):

CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 PXR*2.0*47 GMTS*2.7*113Subject:

Backup of CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 install on Apr

Replace

No routines for CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0

```

Loading Routines for PXRM*2.0*47..

```

```
Routine PXRMBXTL is not on the disk.....
```

```
Routine PXRMCSD is not on the disk..
```

```

Routine PXRMCSSC is not on the disk...

```

```

.....
Routine PXRMCUSU is not on the disk.....

```

```
Routine PXRMPCTX is not on the disk.....
```

Routine PXRMTAXS is not on the disk..

```
Routine PXRMTDUP is not on the disk..
```

```
Routine PXRMTCK is not on the disk.
```

```
Routine PXRMTEDT is not on the disk.....
```

Routine PXRMMVAL is not on the disk..

Routine PXRMCVALC is not on the disk..

```
Routine PXRVALU is not on the disk.....
```

Loading Routines for GMTS*2.7*113.

Send mail to: PXRMPROGRAMMER, ONE

```
Select basket to send to: IN// PATCH BACKUP
```

And Send to:

6. Compare Transport Global to Current System

This option will allow you to view all changes that will be made when the patch is installed. It compares all components of the patch (routines, DDs, templates, etc.).

7. Verify Checksums in Transport Global

This option will allow you to ensure the integrity of the routines that are in the transport global. If there are any discrepancies, do not run the Install Package(s) option. Instead, run the Unload a Distribution option to remove the Transport Global from your system. Retrieve the file again from the anonymous directory (in case there was corruption in SFTPing) and Load the Distribution again. If the problem still exists, log a ticket and/or call the national Help Desk (1-888-596-HELP) to report the problem.

Example

```
CHOOSE 1-2: 2 Option: 2 Verify Checksums in Transport Global
Select INSTALL NAME:
CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 4/24/17@13:31:
45
=> CLINICAL REMINDERS ICD-10 FOLLOWUP ;Created on Apr 24,
2017@08:21:25This Distribution was loaded on Apr 24, 2017@13:31:45 with header
of
CLINICAL REMINDERS ICD-10 FOLLOWUP ;Created on Apr 24, 2017@08:21:25
It consisted of the following Install(s):
CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 PXR*2.0*47 GMTS*2.7*113
Want each Routine Listed with Checksums: Yes//
DEVICE: HOME// ;;999 TELNET PORT
```

8. Install the build

From the Installation menu on the Kernel Installation and Distribution System (KIDS) menu, run the option Install Package(s). Select the build CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 and proceed with the install. If you have problems with the installation, log a Remedy ticket and/or call the National Help Desk to report the problem.

```
Select Installation & Distribution System Option: Installation
Select Installation Option: INSTALL PACKAGE(S)
Select INSTALL NAME: CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0
```

Answer the following install questions as follows:

- Although typically the answer is "No," you can answer "Yes," to the question:

```
Want KIDS to Rebuild Menu Trees Upon Completion of Install?
```

Please remember that rebuilding menu trees will increase patch installation time.

- Answer "No" to the question:

```
Want KIDS to INHIBIT LOGONS during the install?
```

- Answer "No" to the question:

```
Want to DISABLE Scheduled Options, Menu Options, and Protocols?
```


Post-Install Instructions

Post-install: After a successful installation, the init routine PXRMP47I can be deleted.

Acronyms

The OIT Master Glossary is available at
<http://vaww.oed.wss.va.gov/process/Lists/glossary/default.aspx>

National Acronym Directory:

<http://vaww1.va.gov/Acronyms/>

Term	Definition
ASU	Authorization/Subscription Utility
Clin2	National Customer Support team that supports Clinical Reminders
CPRS	Computerized Patient Record System
DBA	Database Administration
DG	Registration and Enrollment Package namespace
ESM	Enterprise Systems Management (ESM)
FIM	Functional Independence Measure
GMPL	Problem List namespace
GMTS	Health Summary namespace (also HSUM)
GUI	Graphic User Interface
HRMH/HRMHP	High Risk Mental Health Patient
IAB	Initial Assessment & Briefing
ICD-10-CM	International Classification of Diseases, 10th Revision Clinical Modification
ICD-9-CM	International Classification of Diseases, Ninth Revision, Clinical Modification
ICR	Internal Control Number
IOC	Initial Operating Capabilities
LSSD	Last Service Separation Date
MH	Mental Health
OHI	Office of Health Information

Term	Definition
OI	Office of Information
OIT/OI&T	Office of Information Technology
OMHS	Office of Mental Health Services
ORR	Operational Readiness Review
PCS	Patient Care Services
PD	Product Development
PIMS	Patient Information Management System
PMAS	Program Management Accountability System
PTM	Patch Tracker Message
PXRM	Clinical Reminder Package namespace
RSD	Requirements Specification Document
SD	Scheduling Package Namespace
SG	Scheduling, Registration, Admission/Discharge/Transfer namespace
SME	Subject Matter Expert
SNOMED CT	Systematic Nomenclature of Medicine of Clinical Terms
SQA	Software Quality Assurance
VA	Department of Veteran Affairs
VHA	Veterans Health Administration
VISN	Veterans Integrated Service Network
VistA	Veterans Health Information System and Technology Architecture

Appendix A: Installation Example

This is a capture of a Clinical Reminders ICD-10 Followup installation that provides details of the install.

Example: First-time Install

```
Select INSTALL NAME:      CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0      Loaded
from D
istribution      4/25/17@09:21:38
=> CLINICAL REMINDERS ICD-10 FOLLOWUP ;Created on Apr 24, 2017@08:21:25

This Distribution was loaded on Apr 25, 2017@09:21:38 with header of
  CLINICAL REMINDERS ICD-10 FOLLOWUP ;Created on Apr 24, 2017@08:21:25
  It consisted of the following Install(s):
CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0      PXR*2.0*47      GMTS*2.7*113
Checking Install for Package CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0

Install Questions for CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0

Checking Install for Package PXR*2.0*47

Install Questions for PXR*2.0*47

Incoming Files:

      802.1      NLM VALUE SET CODING SYSTEMS

      802.2      NLM VALUE SETS

      802.3      NLM QUALITY MEASURE GROUPS

      811.2      REMINDER TAXONOMY
Note:  You already have the 'REMINDER TAXONOMY' File.

      811.4      REMINDER COMPUTED FINDINGS (including data)
Note:  You already have the 'REMINDER COMPUTED FINDINGS' File.
I will OVERWRITE your data with mine.

      811.9      REMINDER DEFINITION (Partial Definition)
Note:  You already have the 'REMINDER DEFINITION' File.

Want KIDS to Rebuild Menu Trees Upon Completion of Install? NO//

Checking Install for Package GMTS*2.7*113

Install Questions for GMTS*2.7*113
```

Want KIDS to INHIBIT LOGONs during the install? NO//
Want to DISABLE Scheduled Options, Menu Options, and Protocols? NO//

Enter the Device you want to print the Install messages.
You can queue the install by enter a 'Q' at the device prompt.
Enter a '^' to abort the install.

DEVICE: HOME// TELNET

Install Started for CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 :
Apr 25, 2017@09:22:59

Build Distribution Date: Apr 24, 2017

Installing Routines:
Apr 25, 2017@09:22:59

Install Started for PXRMP*2.0*47 :
Apr 25, 2017@09:22:59

Build Distribution Date: Apr 24, 2017

Installing Routines:
Apr 25, 2017@09:23

Running Pre-Install Routine: PRE^PXRMP47I

DISABLE options.

DISABLE protocols.

Deleting the DATA DICTIONARY...
Deleting the INPUT TEMPLATES...
Deleting the PRINT TEMPLATES.....
Deleting the SORT TEMPLATES...
Deleting the FORMS...
Deleting the BLOCKS...

Installing Data Dictionaries:
Apr 25, 2017@09:23:08

Installing Data:
Apr 25, 2017@09:23:08

Installing PACKAGE COMPONENTS:

Installing PRINT TEMPLATE

Installing FORM

Installing PROTOCOL
Located in the PXRMP (CLINICAL REMINDERS) namespace.
Located in the PXRMP (CLINICAL REMINDERS) namespace.
Located in the PXRMP (CLINICAL REMINDERS) namespace.

Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.
Located in the PXR (CLINICAL REMINDERS) namespace.

Installing LIST TEMPLATE

Installing OPTION

Apr 25, 2017@09:23:08

Running Post-Install Routine: POST^PXRMP47I

Removing obsolete field data

Working on taxonomy GP IM PNEUMOC PCV13 PREVNAR
Working on taxonomy GP IM PNEUMOC PPSV23 PNEUMOVAX
Working on taxonomy HF LIPID LDL 120-129
Working on taxonomy HF LIPID LDL 71-99
Working on taxonomy HF LIPID LDL >190
Working on taxonomy PALI CONS DYSPNEA MILD (E)
Working on taxonomy PALI CONS DYSPNEA SEVERE (E)
Working on taxonomy PALI CONS INPT 99251 (E)
Working on taxonomy PALI CONS INPT 99252 (E)
Working on taxonomy PALI CONS INPT 99253 (E)
Working on taxonomy PALI CONS INPT 99254 (E)
Working on taxonomy PALI CONS INPT 99255 (E)
Working on taxonomy PALI CONS OUTPT 99241 (E)
Working on taxonomy PALI CONS OUTPT 99242 (E)
Working on taxonomy PALI CONS OUTPT 99243 (E)
Working on taxonomy PALI CONS OUTPT 99244 (E)
Working on taxonomy PALI CONS OUTPT 99245 (E)
Working on taxonomy PALI CONS OUTPT ENCOUNTER GP
Working on taxonomy PALI CONS PAIN NONE (E)
Working on taxonomy PALI CONS PAIN SEVERE (E)
Working on taxonomy PALI CONS PT PREFER DOC NO (E)
Working on taxonomy PALI CONS PT PREFER DOC YES (E)
Working on taxonomy VA-ABD AORTIC ANEURYSM
Working on taxonomy VA-ALCOHOL ABUSE
Working on taxonomy VA-BREAST TUMOR
Working on taxonomy VA-CERVICAL CA/ABNORMAL PAP
Working on taxonomy VA-DEPRESSION DIAGNOSIS
Working on taxonomy VA-DEPRESSION DX OUTPT VISIT
Working on taxonomy VA-DIABETES

Working on taxonomy VA-DIABETES HEDIS
 Working on taxonomy VA-DIABETES HEDIS PROB LIST
 Working on taxonomy VA-DRUG ABUSE FOR HCV TESTING
 Working on taxonomy VA-EBOLA RISK TRIAGE DIAGNOSIS CODE
 Working on taxonomy VA-ECOE DIAGNOSIS CODES
 Working on taxonomy VA-HEPATITIS C INFECTION
 Working on taxonomy VA-HEPATITIS C SEROPOSITIVE
 Working on taxonomy VA-HERPES ZOSTER (SHINGLES) IMMUNIZATION
 Working on taxonomy VA-HIV INFECTION
 Working on taxonomy VA-HYPERTENSION
 Working on taxonomy VA-HYPERTENSION CODES
 Working on taxonomy VA-HYSTERECTOMY
 Working on taxonomy VA-IHD AND ASVD
 Working on taxonomy VA-IM FLU HIGH DOSE
 Working on taxonomy VA-IMAGING FOR AAA (NON-SPECIFIC)
 Working on taxonomy VA-IMMUNOCOMPROMISED
 Working on taxonomy VA-MAMMOGRAM/SCREEN
 Working on taxonomy VA-MASTECTOMY
 Working on taxonomy VA-MHV BILATERAL AMPUTEE
 Working on taxonomy VA-MHV COLONOSCOPY
 Working on taxonomy VA-MHV DIABETIC RETINAL DISEASE
 Working on taxonomy VA-MHV RETINAL EXAMINATION
 Working on taxonomy VA-MHV SIGMOIDOSCOPY
 Working on taxonomy VA-PALLI CONS CANCER/HEMA COND
 Working on taxonomy VA-PALLI CONS CARDIO COND OTHER THAN CA
 Working on taxonomy VA-PALLI CONS CNS COND OTHER THAN CA
 Working on taxonomy VA-PALLI CONS DERM CONDITION DX
 Working on taxonomy VA-PALLI CONS GI OTHER THAN CA
 Working on taxonomy VA-PALLI CONS INFECTIOUS COND
 Working on taxonomy VA-PALLI CONS RENAL OTHER THAN CA
 Working on taxonomy VA-PALLI CONS RHEUM/VASC/THROMB
 Working on taxonomy VA-PNEUMOC DZ RISK - CHEMOTHERAPY
 Working on taxonomy VA-PNEUMOC DZ RISK - HIGH
 Working on taxonomy VA-PNEUMOC DZ RISK - HIGHEST/NOT IMMUNO COMP
 Working on taxonomy VA-PNEUMOC DZ RISK - IMMUNOCOMPROMISED
 Working on taxonomy VA-PNEUMOC PCV13 IMMUNIZATION
 Working on taxonomy VA-PNEUMOC PPSV23 IMMUNIZATION
 Working on taxonomy VA-PSYCHOTHERAPY CPT CODES
 Working on taxonomy VA-PTSD DIAGNOSIS
 Working on taxonomy VA-PTSD DX OUTPT PRIMARY
 Working on taxonomy VA-PTSD DX OUTPT VISIT
 Working on taxonomy VA-SCHIZOPHRENIA
 Working on taxonomy VA-TB/POSITIVE PPD
 Working on taxonomy VA-TD VACCINE ADSORBED GENERIC (RD)
 Working on taxonomy VA-TD VACCINE PRESERV FREE ADSORBED (RD)
 Working on taxonomy VA-TDAP VACCINE (RD)
 Working on taxonomy VA-TERATOGENIC MEDICATIONS ORDER CHECK EXCL (TAXONOMIES)
 Working on taxonomy VA-TERMINAL CANCER PATIENTS
 Working on taxonomy VA-TOBACCO USE
 Working on taxonomy VA-WH BILATERAL MASTECTOMY
 Working on taxonomy VA-WH HYSTERECTOMY W/CERVIX REMOVED
 Working on taxonomy VA-WH PAP SMEAR OBTAINED
 Working on taxonomy VA-WH PAP SMEAR SCREEN CODES
 Working on taxonomy VA-WH TUBAL LIGATION CODES (TAXONOMY)
 Working on taxonomy VA-WH TUBAL REANASTOMOSIS (TAXONOMY)
 Working on taxonomy ZZVA-ALCOHOLISM SCREENING
 Working on taxonomy ZZVA-CERVICAL CANCER SCREEN

Working on taxonomy ZZVA-CHOLESTEROL
 Working on taxonomy ZZVA-COLORECTAL CA
 Working on taxonomy ZZVA-COLORECTAL CANCER SCREEN
 Working on taxonomy ZZVA-EXERCISE COUNSELING
 Working on taxonomy ZZVA-FLEXISIGMOIDOSCOPY
 Working on taxonomy ZZVA-FOBT
 Working on taxonomy ZZVA-HIGH RISK FOR FLU/PNEUMONIA
 Working on taxonomy ZZVA-HIGH RISK FOR INFLUENZA
 Working on taxonomy ZZVA-HIGH RISK FOR TB
 Working on taxonomy ZZVA-HYPERTENSION SCREEN
 Working on taxonomy ZZVA-IM FLU H1N1 (1 DOSE)
 Working on taxonomy ZZVA-INFLUENZA IMMUNIZATION
 Working on taxonomy ZZVA-ISCHEMIC HEART 412 DISEASE
 Working on taxonomy ZZVA-ISCHEMIC HEART DISEASE
 Working on taxonomy ZZVA-MHV IHD AND ATHEROSCLEROSIS
 Working on taxonomy ZZVA-NUTRITION
 Working on taxonomy ZZVA-OBESITY
 Working on taxonomy ZZVA-PNEUMOCOCCAL VACCINE
 Working on taxonomy ZZVA-POLYTRAUMA AMPUTATION
 Working on taxonomy ZZVA-POLYTRAUMA AUDITORY
 Working on taxonomy ZZVA-POLYTRAUMA BRAIN INJURY
 Working on taxonomy ZZVA-POLYTRAUMA BURN
 Working on taxonomy ZZVA-POLYTRAUMA INPT REHAB
 Working on taxonomy ZZVA-POLYTRAUMA ORTHO
 Working on taxonomy ZZVA-POLYTRAUMA PTSD
 Working on taxonomy ZZVA-POLYTRAUMA SCI
 Working on taxonomy ZZVA-POLYTRAUMA VISION
 Working on taxonomy ZZVA-POLYTRAUMA WAR INJURY
 Working on taxonomy ZZVA-PROSTATE CA
 Working on taxonomy ZZVA-PSA
 Working on taxonomy ZZVA-SAFETY COUNSELING
 Working on taxonomy ZZVA-TETANUS DIPHTHERIA
 Working on taxonomy ZZVA-WEIGHT AND NUTRITION SCREEN
 Working on taxonomy ZZVA-WH IUD INSERTION (TAXONOMY)
 Working on taxonomy ZZVA-WH IUD REMOVAL (TAXONOMY)

Building the Selected Codes multiple indexes.

Taxonomy: GP IM PNEUMOC PCV13 PREVNA; IEN=100
 Taxonomy: GP IM PNEUMOC PPSV23 PNEUMOVAX; IEN=71
 Taxonomy: HF LIPID LDL 120-129; IEN=101
 Taxonomy: HF LIPID LDL 71-99; IEN=99
 Taxonomy: HF LIPID LDL >190; IEN=98
 Taxonomy: PALLI CONS DYSPNEA MILD (E); IEN=95
 Taxonomy: PALLI CONS DYSPNEA SEVERE (E); IEN=94
 Taxonomy: PALLI CONS INPT 99251 (E); IEN=91
 Taxonomy: PALLI CONS INPT 99252 (E); IEN=90
 Taxonomy: PALLI CONS INPT 99253 (E); IEN=84
 Taxonomy: PALLI CONS INPT 99254 (E); IEN=83
 Taxonomy: PALLI CONS INPT 99255 (E); IEN=82
 Taxonomy: PALLI CONS OUTPT 99241 (E); IEN=78
 Taxonomy: PALLI CONS OUTPT 99242 (E); IEN=77
 Taxonomy: PALLI CONS OUTPT 99243 (E); IEN=74
 Taxonomy: PALLI CONS OUTPT 99244 (E); IEN=73
 Taxonomy: PALLI CONS OUTPT 99245 (E); IEN=72
 Taxonomy: PALLI CONS OUTPT ENCOUNTER GP; IEN=79
 Taxonomy: PALLI CONS PAIN NONE (E); IEN=97
 Taxonomy: PALLI CONS PAIN SEVERE (E); IEN=96

Taxonomy: PALLI CONS PT PREFER DOC NO (E); IEN=92
 Taxonomy: PALLI CONS PT PREFER DOC YES (E); IEN=93
 Taxonomy: VA-ABD AORTIC ANEURYSM; IEN=13
 Taxonomy: VA-ALCOHOL ABUSE; IEN=17
 Taxonomy: VA-BREAST TUMOR; IEN=18
 Taxonomy: VA-CERVICAL CA/ABNORMAL PAP; IEN=5
 Taxonomy: VA-DEPRESSION DIAGNOSIS; IEN=49
 Taxonomy: VA-DEPRESSION DX OUTPT VISIT; IEN=160
 Taxonomy: VA-DIABETES; IEN=28
 Taxonomy: VA-DIABETES HEDIS; IEN=317
 Taxonomy: VA-DIABETES HEDIS PROB LIST; IEN=316
 Taxonomy: VA-DRUG ABUSE FOR HCV TESTING; IEN=55
 Taxonomy: VA-EBOLA RISK TRIAGE DIAGNOSIS CODE; IEN=153
 Taxonomy: VA-ECOE DIAGNOSIS CODES; IEN=61
 Taxonomy: VA-HEPATITIS C INFECTION; IEN=2
 Taxonomy: VA-HEPATITIS C SEROPOSITIVE; IEN=152
 Taxonomy: VA-HERPES ZOSTER (SHINGLES) IMMUNIZATION; IEN=159
 Taxonomy: VA-HIV INFECTION; IEN=57
 Taxonomy: VA-HYPERTENSION; IEN=1
 Taxonomy: VA-HYPERTENSION CODES; IEN=46
 Taxonomy: VA-HYSTERECTOMY; IEN=6
 Taxonomy: VA-IHD AND ASVD; IEN=315
 Taxonomy: VA-IM FLU HIGH DOSE; IEN=243
 Taxonomy: VA-IMAGING FOR AAA (NON-SPECIFIC); IEN=12
 Taxonomy: VA-IMMUNOCOMPROMISED; IEN=156
 Taxonomy: VA-MAMMOGRAM/SCREEN; IEN=16
 Taxonomy: VA-MASTECTOMY; IEN=19
 Taxonomy: VA-MHV BILATERAL AMPUTEE; IEN=87
 Taxonomy: VA-MHV COLONOSCOPY; IEN=85
 Taxonomy: VA-MHV DIABETIC RETINAL DISEASE; IEN=88
 Taxonomy: VA-MHV RETINAL EXAMINATION; IEN=89
 Taxonomy: VA-MHV SIGMOIDOSCOPY; IEN=86
 Taxonomy: VA-PALLI CONS CANCER/HEMA COND; IEN=70
 Taxonomy: VA-PALLI CONS CARDIO COND OTHER THAN CA; IEN=68
 Taxonomy: VA-PALLI CONS CNS COND OTHER THAN CA; IEN=69
 Taxonomy: VA-PALLI CONS DERM CONDITION DX; IEN=65
 Taxonomy: VA-PALLI CONS GI OTHER THAN CA; IEN=67
 Taxonomy: VA-PALLI CONS INFECTIOUS COND; IEN=63
 Taxonomy: VA-PALLI CONS RENAL OTHER THAN CA; IEN=66
 Taxonomy: VA-PALLI CONS RHEUM/VASC/THROMB; IEN=64
 Taxonomy: VA-PNEUMOC DZ RISK - CHEMOTHERAPY; IEN=80
 Taxonomy: VA-PNEUMOC DZ RISK - HIGH; IEN=3
 Taxonomy: VA-PNEUMOC DZ RISK - HIGHEST/NOT IMMUNO COMP; IEN=344
 Taxonomy: VA-PNEUMOC DZ RISK - IMMUNOCOMPROMISED; IEN=81
 Taxonomy: VA-PNEUMOC PCV13 IMMUNIZATION; IEN=149
 Taxonomy: VA-PNEUMOC PPSV23 IMMUNIZATION; IEN=103
 Taxonomy: VA-PSYCHOTHERAPY CPT CODES; IEN=50
 Taxonomy: VA-PTSD DIAGNOSIS; IEN=161
 Taxonomy: VA-PTSD DX OUTPT PRIMARY; IEN=37
 Taxonomy: VA-PTSD DX OUTPT VISIT; IEN=176
 Taxonomy: VA-SCHIZOPHRENIA; IEN=52
 Taxonomy: VA-TB/POSITIVE PPD; IEN=10
 Taxonomy: VA-TD VACCINE ADSORBED GENERIC (RD); IEN=157
 Taxonomy: VA-TD VACCINE PRESERV FREE ADSORBED (RD); IEN=60
 Taxonomy: VA-TDAP VACCINE (RD); IEN=102
 Taxonomy: VA-TERATOGENIC MEDICATIONS ORDER CHECK EXCL (TAXONOMIES); IEN=108
 Taxonomy: VA-TERMINAL CANCER PATIENTS; IEN=56

Taxonomy: VA-TOBACCO USE; IEN=22
 Taxonomy: VA-WH BILATERAL MASTECTOMY; IEN=62
 Taxonomy: VA-WH HYSTERECTOMY W/CERVIX REMOVED; IEN=58
 Taxonomy: VA-WH PAP SMEAR OBTAINED; IEN=48
 Taxonomy: VA-WH PAP SMEAR SCREEN CODES; IEN=59
 Taxonomy: VA-WH TUBAL LIGATION CODES (TAXONOMY); IEN=75
 Taxonomy: VA-WH TUBAL REANASTOMOSIS (TAXONOMY); IEN=76
 Taxonomy: ZZVA-ALCOHOLISM SCREENING; IEN=34
 Taxonomy: ZZVA-CERVICAL CANCER SCREEN; IEN=30
 Taxonomy: ZZVA-CHOLESTEROL; IEN=24
 Taxonomy: ZZVA-COLORECTAL CA; IEN=4
 Taxonomy: ZZVA-COLORECTAL CANCER SCREEN; IEN=31
 Taxonomy: ZZVA-EXERCISE COUNSELING; IEN=32
 Taxonomy: ZZVA-FLEXISIGMOIDOSCOPY; IEN=15
 Taxonomy: ZZVA-FOBT; IEN=27
 Taxonomy: ZZVA-HIGH RISK FOR FLU/PNEUMONIA; IEN=9
 Taxonomy: ZZVA-HIGH RISK FOR INFLUENZA; IEN=7
 Taxonomy: ZZVA-HIGH RISK FOR TB; IEN=11
 Taxonomy: ZZVA-HYPERTENSION SCREEN; IEN=23
 Taxonomy: ZZVA-IM FLU H1N1 (1 DOSE); IEN=227
 Taxonomy: ZZVA-INFLUENZA IMMUNIZATION; IEN=33
 Taxonomy: ZZVA-ISCHEMIC HEART 412 DISEASE; IEN=44
 Taxonomy: ZZVA-ISCHEMIC HEART DISEASE; IEN=14
 Taxonomy: ZZVA-MHV IHD AND ATHEROSCLEROSIS; IEN=40
 Taxonomy: ZZVA-NUTRITION; IEN=21
 Taxonomy: ZZVA-OBESITY; IEN=20
 Taxonomy: ZZVA-PNEUMOCOCCAL VACCINE; IEN=25
 Taxonomy: ZZVA-POLYTRAUMA AMPUTATION; IEN=54
 Taxonomy: ZZVA-POLYTRAUMA AUDITORY; IEN=53
 Taxonomy: ZZVA-POLYTRAUMA BRAIN INJURY; IEN=51
 Taxonomy: ZZVA-POLYTRAUMA BURN; IEN=47
 Taxonomy: ZZVA-POLYTRAUMA INPT REHAB; IEN=45
 Taxonomy: ZZVA-POLYTRAUMA ORTHO; IEN=43
 Taxonomy: ZZVA-POLYTRAUMA PTSD; IEN=42
 Taxonomy: ZZVA-POLYTRAUMA SCI; IEN=41
 Taxonomy: ZZVA-POLYTRAUMA VISION; IEN=39
 Taxonomy: ZZVA-POLYTRAUMA WAR INJURY; IEN=38
 Taxonomy: ZZVA-PROSTATE CA; IEN=8
 Taxonomy: ZZVA-PSA; IEN=26
 Taxonomy: ZZVA-SAFETY COUNSELING; IEN=35
 Taxonomy: ZZVA-TETANUS DIPHTHERIA; IEN=29
 Taxonomy: ZZVA-WEIGHT AND NUTRITION SCREEN; IEN=36
 Taxonomy: ZZVA-WH IUD INSERTION (TAXONOMY); IEN=109
 Taxonomy: ZZVA-WH IUD REMOVAL (TAXONOMY); IEN=110

Installing Value Set data.

Making sure all DO IN ADVANCE TIME FRAMEs are uppercase.

Making sure all reminder frequencies are uppercase.

Found lower case baseline frequency:

IEN=128 IND=2 Frequency=0y

Changing it to upper case.

ENABLE options.

```
ENABLE protocols.

Updating Routine file...

Updating KIDS files...

PXRM*2.0*47 Installed.
      Apr 25, 2017@09:23:40

Not a production UCI

NO Install Message sent

Install Started for GMTS*2.7*113 :
      Apr 25, 2017@09:23:40

Build Distribution Date: Apr 24, 2017

Installing Routines:
      Apr 25, 2017@09:23:40

Updating Routine file...

Updating KIDS files...

GMTS*2.7*113 Installed.
      Apr 25, 2017@09:23:40

Not a production UCI

NO Install Message sent

Updating Routine file...

Updating KIDS files...

CLINICAL REMINDERS ICD-10 FOLLOWUP 1.0 Installed.
      Apr 25, 2017@09:23:40

No link to PACKAGE file

Install Completed
```

Appendix B: Post-Install Checksums

Installation Checksums: PXRМ Routines

CPRS285A4:CPRS30>D CHECK1^XTSUMBLD

New CheckSum CHECK1^XTSUMBLD:

This option determines the current checksum of selected routine(s).
The Checksum of the routine is determined as follows:

1. Any comment line with a single semi-colon is presumed to be followed by comments and only the line tag will be included.
2. Line 2 will be excluded from the count.
3. The total value of the routine is determined (excluding exceptions noted above) by multiplying the ASCII value of each character by its position on the line and position of the line in the routine being checked.

Select one of the following:

P	Package
B	Build

Build from: Build

This will check the routines from a BUILD file.

Select BUILD NAME: PXRМ*2.0*47 CLINICAL REMINDERS

PXRМ	value = 47307370
PXRМASL	value = 4451968
PXRМCDEF	value = 4852404
PXRМCOND	value = 44611363
PXRМCQIN	value = 24112852
PXRМCQLM	value = 20001822
PXRМCSTX	value = 2714904
PXRМDATE	value = 72754452
PXRМDBLD	value = 54628364
PXRМDEV	value = 53641052
PXRМDIEV	value = 75240155
PXRМDLL	value = 147254046
PXRМDLLB	value = 36248685
PXRМDRUG	value = 62235689
PXRМDTAX	value = 198889216
PXRМERRH	value = 21386013
PXRМETXR	value = 76063275
PXRМEXCC	value = 14643642
PXRМEXCO	value = 26994015
PXRМEXCS	value = 16682539
PXRМEXFI	value = 59766502
PXRМEXIC	value = 87332423
PXRМEXIU	value = 86717543
PXRМEXLM	value = 51322537
PXRМEXMH	value = 10558717
PXRМEXU0	value = 29114606
PXRМEXWB	value = 1606210
PXRМFF	value = 75070206
PXRМFF0	value = 19175056

PXRMFIND	value = 14402705
PXRMFMTO	value = 10636061
PXRMFNFT	value = 58123785
PXRMHVET	value = 1748161
PXRMICHK	value = 238423274
PXRMLDR	value = 18348692
PXRMLEX	value = 8264453
PXRMLEXL	value = 190120643
PXRMLOCF	value = 97441474
PXRMLOCL	value = 28870554
PXRMLOG	value = 66402903
PXRMSG	value = 5151324
PXRMORCH	value = 53534150
PXRMOUTC	value = 38936387
PXRMOUTD	value = 15041011
PXRMOUTM	value = 29787683
PXRMOUTU	value = 18283584
PXRMP47I	value = 30350491
PXRMPDEM	value = 65842415
PXRMPDR	value = 52184470
PXRMPDRP	value = 100715785
PXRMPRF	value = 15497099
PXRMRCUR	value = 13247111
PXRMRDF	value = 84791298
PXRMRDT	value = 70476581
PXRMRPCA	value = 57506116
PXRMRUL1	value = 49947013
PXRMSXRM	value = 98450970
PXRMTAXD	value = 79799584
PXRMTAXL	value = 80891432
PXRMTERM	value = 55250159
PXRMTXT	value = 44027439
PXRMTXCR	value = 23196731
PXRMTXDL	value = 1372560
PXRMTXIM	value = 196801557
PXRMTXIN	value = 62837955
PXRMTXLS	value = 150968479
PXRMTXSM	value = 56360268
PXRMUIDE	value = 17671178
PXRMUTIL	value = 164856611
PXRMCPT	value = 52602127
PXRMVPOV	value = 51972823
PXRMVSCS	value = 6081007
PXRMSIN	value = 48824341
PXRMSLM	value = 26634140
PXRMSVTX	value = 61497746
PXRMXEVL	value = 1513078
PXRMXSE1	value = 30585167
done	

Installation Checksums: GMTS Routines

CPRS285A4:CPRS30>D CHECK1^XTSUMBLD

New CheckSum CHECK1^XTSUMBLD:

This option determines the current checksum of selected routine(s).
The Checksum of the routine is determined as follows:

1. Any comment line with a single semi-colon is presumed to be followed by comments and only the line tag will be included.
2. Line 2 will be excluded from the count.
3. The total value of the routine is determined (excluding exceptions noted above) by multiplying the ASCII value of each character by its position on the line and position of the line in the routine being checked.

Select one of the following:

P	Package
B	Build

Build from: Build

This will check the routines from a BUILD file.

Select BUILD NAME: GMTS*2.7*113 HEALTH SUMMARY

GMTSPXHR value = 21877798

done

Appendix C: Install File Print Example

Use the KIDS Install File Print option to print out the results of the installation process. You can select the multi-package build or any of the individual builds included in the multi-package build.

```
Select OPTION NAME:      XPD PRINT INSTALL FILE      Install File Print
Install File Print

Select INSTALL NAME: PXRM*2.0*47      Install Completed      4/25/17@09:23:40
=> CLINICAL REMINDERS ICD-10 FOLLOWUP ;Created on Apr 24, 2017@08:21:25
DEVICE: HOME// ;;999 TELNET

PACKAGE: PXRM*2.0*47      Apr 25, 2017 9:40 am      PAGE 1
                                COMPLETED      ELAPSED
-----
STATUS: Install Completed      DATE LOADED: APR 25, 2017@09:21:38
INSTALLED BY: TESTMASTER,USER
NATIONAL PACKAGE: CLINICAL REMINDERS

INSTALL STARTED: APR 25, 2017@09:22:59      09:23:40      0:00:41

ROUTINES:      09:23      0:00:01

PRE-INIT CHECK POINTS:
XPD PREINSTALL STARTED      09:23:07      0:00:07
XPD PREINSTALL COMPLETED      09:23:07

FILES:
NLM VALUE SET CODING SYSTEMS      09:23:07
NLM VALUE SETS      09:23:07
NLM QUALITY MEASURE GROUPS      09:23:08      0:00:01
REMINDER TAXONOMY      09:23:08
REMINDER COMPUTED FINDINGS      09:23:08
REMINDER DEFINITION      09:23:08

PRINT TEMPLATE      09:23:08
FORM      09:23:08
PROTOCOL      09:23:08
LIST TEMPLATE      09:23:08
OPTION      09:23:08

POST-INIT CHECK POINTS:
XPD POSTINSTALL STARTED      09:23:40      0:00:32
XPD POSTINSTALL COMPLETED      09:23:40

INSTALL QUESTION PROMPT      ANSWER

XP01      Want KIDS to Rebuild Menu Trees Upon Completion of Install      NO
MESSAGES:

      Install Started for PXRM*2.0*47 :
              Apr 25, 2017@09:22:59

Build Distribution Date: Apr 24, 2017

      Installing Routines:
              Apr 25, 2017@09:23

      Running Pre-Install Routine: PRE^PXRMP47I
```

DISABLE options.

DISABLE protocols.

Installing Data Dictionaries:

Apr 25, 2017@09:23:08

Installing Data:

Apr 25, 2017@09:23:08

Installing PACKAGE COMPONENTS:

Installing PRINT TEMPLATE

Installing FORM

Installing PROTOCOL

Installing LIST TEMPLATE

Installing OPTION

Apr 25, 2017@09:23:08

Running Post-Install Routine: POST^PXRMP47I

Removing obsolete field data

Working on taxonomy GP IM PNEUMOC PCV13 PREVNAR

Working on taxonomy GP IM PNEUMOC PPSV23 PNEUMOVAX

Working on taxonomy HF LIPID LDL 120-129

Etc. See Installation example in [Appendix A](#)

Appendix D: Build File Print Example

Use the KIDS Build File Print option to print out the build components. You can select the multi-package build or any of the individual builds included in the multi-package build.

```
Select OPTION NAME:      XPD PRINT BUILD      Build File Print
Build File Print
Select BUILD NAME: PXRM*2.0*47      CLINICAL REMINDERS
DEVICE: HOME// ;;999 TELNET PORT
PACKAGE: PXRM*2.0*47      Apr 25, 2017 9:45 am      PAGE 1
-----
TYPE: SINGLE PACKAGE      TRACK NATIONALLY: YES
NATIONAL PACKAGE: CLINICAL REMINDERS      ALPHA/BETA TESTING: NO

DESCRIPTION:
Clinical Reminders ICD-10 follow-up build.

ENVIRONMENT CHECK:      DELETE ENV ROUTINE:
  PRE-INIT ROUTINE: PRE^PXRM*2.0*47      DELETE PRE-INIT ROUTINE: No
  POST-INIT ROUTINE: POST^PXRM*2.0*47      DELETE POST-INIT ROUTINE: No
  PRE-TRANSPORT RTN: VSSAVE^PXRM*2.0*47

FILE #      FILE NAME      UP      SEND      DATA      USER
DATE      SEC.      COMES      SITE      RSLV      OVER
DD      CODE      W/FILE      DATA      PTRS      RIDE
-----
802.1      NLM VALUE SET CODING SYSTEMS      YES      YES      NO
802.2      NLM VALUE SETS      YES      YES      NO
802.3      NLM QUALITY MEASURE GROUPS      YES      YES      NO
811.2      REMINDER TAXONOMY      YES      YES      NO
811.4      REMINDER COMPUTED FINDINGS      NO      YES      YES      OVER      NO      NO
      DATA SCREEN: I $$$CFINC^PXRM*2.0*47(Y)
811.9      REMINDER DEFINITION      YES      YES      NO      NO
Partial DD: subDD: 811.9256      fld: .01

PRINT TEMPLATE:      ACTION:
  PXRM TAXONOMY DIALOG      FILE #811.2      DELETE AT SITE
  PXRM TAXONOMY INQUIRY      FILE #811.2      DELETE AT SITE
  PXRM TAXONOMY LIST      FILE #811.2      DELETE AT SITE
  PXRM TAXONOMY LIST HEADER      FILE #811.2      DELETE AT SITE

FORM:      ACTION:
  PXRM DIALOG TAXONOMY EDIT      FILE #811.2      SEND TO SITE
  PXRM TAXONOMY CHANGE LOG      FILE #811.2      SEND TO SITE
  PXRM TAXONOMY EDIT      FILE #811.2      SEND TO SITE

ROUTINE:      ACTION:
  PXRM      SEND TO SITE
  PXRMASL      SEND TO SITE
  PXRMBXTL      DELETE AT SITE
  PXRMCDEF      SEND TO SITE
  PXRMCOND      SEND TO SITE
  PXRMCQIN      SEND TO SITE
  PXRMCQLM      SEND TO SITE
```


PXRMTXDL	SEND TO SITE
PXRMTXIM	SEND TO SITE
PXRMTXIN	SEND TO SITE
PXRMTXLS	SEND TO SITE
PXRMTXSM	SEND TO SITE
PXRMUIDE	SEND TO SITE
PXRMUTIL	SEND TO SITE
PXRMVAL	DELETE AT SITE
PXRMVALC	DELETE AT SITE
PXRMVALU	DELETE AT SITE
PXRMVCPT	SEND TO SITE
PXRMVPOV	SEND TO SITE
PXRMVSCS	SEND TO SITE
PXRMVSIN	SEND TO SITE
PXRMVSLM	SEND TO SITE
PXRMVSTX	SEND TO SITE
PXRMXEVL	SEND TO SITE
PXRMXSE1	SEND TO SITE
OPTION:	
PXRM DISABLE/ENABLE EVALUATION	SEND TO SITE
PXRM MANAGERS MENU	USE AS LINK FOR MENU/ITEM/SUBS
CRIBERS	
PXRM NLM CQM MENU	SEND TO SITE
PXRM NLM VALUE SET MENU	SEND TO SITE
PXRM TAX SELECTED CODES REPAIR	DELETE AT SITE
PXRM TAXONOMY CODE SEARCH	DELETE AT SITE
PXRM TAXONOMY COPY	DELETE AT SITE
PXRM TAXONOMY EDIT	DELETE AT SITE
PXRM TAXONOMY EXPANSION	DELETE AT SITE
PXRM TAXONOMY EXPANSION ALL	DELETE AT SITE
PXRM TAXONOMY EXPANSION VERIFY	DELETE AT SITE
PXRM TAXONOMY INQUIRY	DELETE AT SITE
PXRM TAXONOMY LIST	DELETE AT SITE
PXRM TAXONOMY MANAGEMENT (OLD)	DELETE AT SITE
PROTOCOL:	
PXRM CQM INQUIRE	SEND TO SITE
PXRM CQM MENU	SEND TO SITE
PXRM CQM SELECT ENTRY	SEND TO SITE
PXRM DIALOG ADD	DELETE AT SITE
PXRM DIALOG EXIT	SEND TO SITE
PXRM DIALOG HISTORY	SEND TO SITE
PXRM DIALOG LINK	SEND TO SITE
PXRM LEXICON REMOVE FROM DIALOG	SEND TO SITE
PXRM LEXICON VIEW	DELETE AT SITE
PXRM TAXONOMY ADD	SEND TO SITE
PXRM TAXONOMY ALL SELECTED ADD	SEND TO SITE
PXRM TAXONOMY ALL SELECTED EXIT	SEND TO SITE
PXRM TAXONOMY ALL SELECTED MENU	SEND TO SITE
PXRM TAXONOMY ALL SELECTED RFD	SEND TO SITE
PXRM TAXONOMY ALL SELECTED RFT	SEND TO SITE
PXRM TAXONOMY ALL SELECTED SAVE	SEND TO SITE
PXRM TAXONOMY ALL SELECTED SELECT	SEND TO SITE
PXRM TAXONOMY ALL SELECTED UID	SEND TO SITE
PXRM TAXONOMY ALL SELECTED UID EDIT	SEND TO SITE
PXRM TAXONOMY CHANGE LOG	SEND TO SITE
PXRM TAXONOMY CODE SEARCH	SEND TO SITE
PXRM TAXONOMY COPY	SEND TO SITE
PXRM TAXONOMY EDIT	SEND TO SITE
PXRM TAXONOMY IMPORT	SEND TO SITE
PXRM TAXONOMY INQUIRE	SEND TO SITE
PXRM TAXONOMY MENU	SEND TO SITE

PXRM TAXONOMY OLD INQUIRE	DELETE AT SITE
PXRM TAXONOMY UID EDIT	SEND TO SITE
PXRM TAXONOMY UID REPORT	SEND TO SITE
PXRM TAXONOMY UIDE EXIT	SEND TO SITE
PXRM TAXONOMY UIDE SAVE	SEND TO SITE
PXRM TAXONOMY UIDE SELECT ENTRY	SEND TO SITE
PXRM TAXONOMY VALUE SET COMPARE	SEND TO SITE
PXRM VS CODE SEARCH	SEND TO SITE
PXRM VS CREATE TAXONOMY	SEND TO SITE
PXRM VS INQUIRE	SEND TO SITE
PXRM VS MENU	SEND TO SITE
PXRM VS SELECT ENTRY	SEND TO SITE
LIST TEMPLATE:	ACTION:
PXRM CQM MENU	SEND TO SITE
PXRM EX MAIN HELP	DELETE AT SITE
PXRM TAXONOMY ALL SELECTED	SEND TO SITE
PXRM TAXONOMY UID EDIT	SEND TO SITE
PXRM VS MENU	SEND TO SITE
INSTALL QUESTIONS:	
Default Rebuild Menu Trees Upon Completion of Install: NO	
Default INHIBIT LOGONs during the install: NO	
Default DISABLE Scheduled Options, Menu Options, and Protocols: NO	
REQUIRED BUILDS:	ACTION:
PXRM*2.0*26	Don't install, leave global
PX*1.0*201	Don't install, leave global
XU*8.0*657	Don't install, leave global
PXRM*2.0*53	Don't install, leave global