



Update _2_0_63

Clinical Reminders

MENTAL HEALTH REMINDER UPDATES 6.0

Install Guide

February 2019

Product Development
Office of Information Technology
Department of Veterans Affairs

Table of Contents

INTRODUCTION	1
PRE-INSTALL	6
INSTALL DETAILS	8
INSTALL EXAMPLE	8
POST INSTALLATION	16

Introduction

On October 1, 2018, Depression and PTSD screening changed to use the Mental Health instruments PHQ-2+I9 and PC-PTSD-5+I9 respectively. These were released nationally as part of UPDATE_2_0_53. We have found a couple of issues with the logic and have corrected those.

One of the main objectives of this Update was to modify the reminder logic for Depression and Alcohol screenings done during an inpatient stay to not factor into the resolution and due date of these reminders. Due to software limitations, we were not able to accomplish this at this current time. Instead, this will be resolved with new functionality released with CPRS31B. At that time, we will update the MH reminders to meet the EPRP measure. For now, a screening done during an inpatient stay will resolve the outpatient clinical reminder. EPRP is aware of this and should not be counted against your site.

In addition, the following logic issues were addressed:

When a screen (depression or PTSD) is positive using the new tools, and patient has a prior negative screen using the old tool, and follow-up is not done following the positive screen, the screening reminders were not coming DUE after 36 hours.

Patient has a positive PTSD and the follow-up was completed. Later a Depression screen is positive, and the follow-up is completed. This was causing the PTSD screening to come due again because the follow-up looks to be later than the old positive screen.

Conflicting results of depression or PTSD screen on same day did not allow the follow up to be completed. This was fixed by adding an occurrence count of 4 and setting the USE STATUS/CONDITION IN SEARCH to YES for findings 1 and 2 of reminder definition VA-BL SHOW C-SSRS

When a positive Audit-C and a SUD clinic visit was same date/time, this was causing the alcohol follow up reminder to be DUE NOW instead of N/A. This was fixed by updating FF(3) string from MRD(7) > MRD(1) & (DIFF_DATE(1,7) < 15) to Function String:
MRD(7) < MRD(1) & (DIFF_DATE(1,7) < 15), FF(6) string from
MRD(7) > MRD(1) & (DIFF_DATE(1,7) > 14) to
MRD(7) < MRD(1) & (DIFF_DATE(1,7) > 14), FF(7) string from
MRD(1) > MRD(7) & (DIFF_DATE(1,7) > 90) to
MRD(1) < MRD(7) & (DIFF_DATE(1,7) > 90). Cohort logic was changed from
(FI(1) & FF(9) & (FF(7) ! FF(3))) ! ((FI(1) & (FF(7) ! FF(3))) ! (FI(9) & FI(10))) & (FI(1) & FI(7) & FF(6)) to
(FI(1) & FF(9) & (FF(7) ! FF(3))) ! ((FI(1) & (FF(7) ! FF(3))) ! (FI(9) & FI(10))) & (FI(1) & FI(7)).

Other changes included in this Update:

Reminder Term VA-USER CLASS FOR COGNITIVE IMPAIRMENT needs mapped items to have T as the Beginning date/time because the term may not evaluate correctly.

Reminder dialog result group PXRМ PC-PTSD-5+I9 RESULT ELEMENT 1 SRE needed Progress note text changed from “revealed suicidal ideation” to “revealed no suicidal ideation”

Reminder Term VA-DEPRESSION SCREEN OEF/OIF had PHQ-2+I9 added as a mapped item.

The following reminder dialogs emphasized the word “positive” by making it all caps “POSITIVE”:

PXRМ C-SSRS RESULT ELEMENT 3
PXRМ C-SSRS RESULT ELEMENT 4
PXRМ C-SSRS RESULT ELEMENT 5
PXRМ C-SSRS RESULT ELEMENT 8

In reminder dialog for follow-up for positive depression/PTSD, the results of a previous C_SSRS (done on same day) are now displayed to provider within the template as they are made aware there is no need to perform a second C_SSRS on same day. This will prevent the provider from taking time to search for this information.

Updated reminder taxonomy VA-DEPRESSION DX OUPT VISIT to meet HEDIS standards.

To help with confusion on “next steps” for positive C-SSRS, a text element VA-C-SSRS POS COMPLETE CSRE TEXT was added to reminder dialog VA-FOLLOW-UP POS PTSD/DEP COMBINATION. The dialog text is “Once completed, the C-SSRS screener with any positive answers that require a same-day Comprehensive Suicide Risk Evaluation will be indicated both in a pop-up window and in progress note text”

Reminder Definition VA-BL C-SSRS DONE TODAY POSITIVE was updated to add an occurrence count of 2 to each of the four findings. The USE STATUS/CONDITION IN SEARCH was also set to YES for all four findings. Findings that were updated are:

VA-C-SSRS DONE TODAY POSITIVE QUESTION 3
VA-C-SSRS DONE TODAY POSITIVE QUESTION 4
VA-C-SSRS DONE TODAY POSITIVE QUESTION 5
VA-C-SSRS DONE TODAY POSITIVE QUESTION 8

In the VA-DEPRESSION SCREENING reminder, term VA-LIFE EXPECTANCY < 6 MONTHS did not have a T-6M as a BEGINNING DATE/TIME entry. This was added to look back 6 months for any occurrence of a mapped item.

On the VA-FOLLOW-UP POS PTSD/DEP COMBINATION dialog, if a user selected “Patient to be evaluated by Mental Health” on dialog group VA-PTSD/DEP SCREEN POS MH OPTIONS GP as a disposition and didn’t select one of the radio buttons, there was no text being entered into progress note due to the EXCLUDE FROM PROGRESS NOTE being set to YES. This was changed to the default of NO and, we are now requiring a selection be made for Stat/Emergent or Routine evaluation.

```

----- MENTAL HEALTH SCREEN DISPOSITIONS -----
I have reviewed the results of the Mental Health screens and have evaluated the patient. Based on the evaluation, the following disposition
plan will be implemented. Select as many options from the list below as needed

 No further intervention is needed at this time. Contact information and instructions for accessing emergency services provided.
 Already receiving needed treatment. Contact information and instructions for accessing emergency services provided.
 Patient to obtain further evaluation/treatment outside this VA. Contact information and instructions for accessing emergency
services provided.
 Patient to be evaluated by Mental Health
 Stat/Emergent Mental Health Evaluation needed.
 Routine/Non-emergent Mental Health Evaluation needed.
  
```

Four reminder definitions and their associated dialog are being distributed with this update except for the Alcohol reminder which is being distributed without a dialog attached. Sites will need to attach the stand alone or the combination dialog to the alcohol screening reminder:

- VA-ALCOHOL USE SCREEN (AUDIT-C)
- VA-ALCOHOL AUDIT-C POSITIVE F/U EVAL
- VA-DEPRESSION SCREENING
- VA-PTSD SCREENING

UPDATE_2_0_63 contains 1 Reminder Exchange entry:
UPDATE_2_0_63 MH UPDATES 6.0

The exchange file contains the following components:

REMINDER DEFINITION

- VA-BL C-SSRS DONE TODAY POSITIVE
- VA-BL SHOW C-SSRS
- VA-ALCOHOL USE SCREEN (AUDIT-C)
- VA-BL AUDC 1 TO 4 NO AUD DX
- VA-BL ALCOHOL WITHIN SAFE LIMIT
- VA-BL AUD-C>7/NO AUD DX
- VA-BL AUD-C 5-7/NO AUD DX
- VA-BL AUD-C>0 & AUD DX
- VA-ALCOHOL AUDIT-C POSITIVE F/U EVAL
- VA-BL NURSE/LPN COG IMPAIRMENT
- VA-PTSD SCREENING
- VA-DEPRESSION SCREENING

REMINDER DIALOG

VA-GP POS PTSD/DEP OBJECTS
VA-FOLLOW-UP POS PTSD/DEP
COMBINATION
VA-PTSD SCREENING
VA-DEPRESSION SCREEN

REMINDER TERM

VA-USER CLASS FOR COGNITIVE IMPAIRMENT
VA-DEPRESSION SCREEN OEF/OIF
VA-REMINDER UPDATE_2_0_63
VA-BL C-SSRS DONE TODAY
VA-C-SSRS DONE TODAY POSITIVE QUESTION 8
VA-C-SSRS DONE TODAY POSITIVE QUESTION 5
VA-C-SSRS DONE TODAY POSITIVE QUESTION 4
VA-C-SSRS DONE TODAY POSITIVE QUESTION 3
VA-BL C-SSRS DONE TODAY POSITIVE
VA-BL SHOW C-SSRS
VA-MH ACUTE ILLNESS
VA-POS ALCOHOL SCREEN
VA-COGNITIVE IMPAIRMENT
VA-REFUSED ALCOHOL SCREENING
VA-LIFE EXPECTANCY < 6 MONTHS
VA-ALCOHOL USE SCREEN
VA-ALC USE INTERVENTIONS
VA-BL NURSE/LPN COG IMPAIRMENT
VA-PC-PTSD NEGATIVE
VA-VETERAN PATIENT
VA-PTSD DIAGNOSIS
VA-SUICIDE SCREEN POSITIVE PTSD
VA-SUICIDE SCREEN NEGATIVE PTSD
VA-PTSD/DEPRESSION ASSESSMENT COMPLETE
VA-REFUSED PTSD SCREEN
VA-PTSD SCREEN POSITIVE
VA-PTSD SCREEN NEGATIVE
VA-PTSD SCREEN
VA-PHQ-2 NEGATIVE
VA-SUICIDE SCREEN POSITIVE DEPRESSION
VA-SUICIDE SCREEN NEGATIVE DEPRESSION
VA-REFUSED DEPRESSION SCREENING
VA-DEPRESSION SCREEN POSITIVE
VA-DEPRESSION SCREEN NEGATIVE
VA-DEPRESSION DIAGNOSIS

REMINDER TAXONOMY

VA-PTSD DX OUTPT VISIT
VA-TERMINAL CANCER PATIENTS

VA-DEPRESSION DX OUTPT VISIT

REMINDER COMPUTED FINDINGS

VA-REMINDER DEFINITION
VA-FILEMAN DATE
VA-ASU USER CLASS
VA-VETERAN
VA-SERVICE SEPARATION DATES
VA-WAS INPATIENT

REMINDER SPONSOR

Office of Mental Health Services
Office of Public Health and Environmental Hazards
Mental Health Services
Office of Quality & Performance
Mental Health Services and Women Veterans Health
Program
National Clinical Practice Guideline Council
Office of Patient Care Services
Mental Health and Behavioral Science Strategic
Group
National Center for Health Promotion and Disease
Prevention (NCP)

HEALTH FACTORS

VA-REMINDER UPDATES
VA-UPDATE_2_0_63
PTSD ON GUARD
PTSD SCREEN - ON GUARD
PTSD SCREEN - NO ON GUARD
PTSD NIGHTMARES
PTSD SCREEN - NO NIGHTMARES
PTSD SCREEN - NIGHTMARES
PTSD DETACHMENT
PTSD SCREEN - DETACHED
PTSD SCREEN - NO DETACHMENT
PTSD AVOIDANCE
PTSD SCREEN - NO AVOIDANCE
PTSD SCREEN - AVOIDANCE
MENTAL HEALTH
REFUSED PTSD SCREEN
REFUSED DEPRESSION SCREENING
VA-NURSE/LPN-PATIENT COG IMPAIRMENT
SEVERE CHRONIC COGNITIVE IMPAIRMENT

UNABLE TO SCREEN - ACUTE ILLNESS
REMINDER FACTORS
LIFE EXPECTANCY < 6 MONTHS
DEP/PTSD/SI ASSESSMENT COMPLETE
PTSD ASSESSMENT COMPLETE
DEPRESSION ASSESS COMPLETE

TIU TEMPLATE FIELD

VA-PTSD CPG
BLANK TEXT FIELD FOR DIALOGS
VA-PTSD PCL INFO
BLANK SPACE1
VA-URL PHQ9
VA-DEP CPG

HEALTH SUMMARY COMPONENT

MHA Score
PCE HEALTH FACTORS SELECTED

HEALTH SUMMARY TYPE

PC-PTSD PCL+I9
PHQ2+I9
VA-CSRE DONE TODAY
OB AUDC

HEALTH SUMMARY OBJECTS

PC-PTSD PCL+I9 (TIU)
PHQ2+I9 (TIU)
VA-CSRE DONE TODAY (TIU)
AUDC (TIU)

TIU DOCUMENT DEFINITION

PC-PTSD PCL+I9
PHQ2+I9
VA-CSRE DONE TODAY
AUDICT-C

Pre-Install

1. Sites may want to check for quick order(s) in the VA-FOLLOW-UP POS PTSD/DEP COMBINATION dialog and have those available when installing this dialog. Even though the name of the elements leads one to think this is only for PTSD, this should be a generic consult to MH for either a depression and/or a PTSD issue. Below are the elements you can check for Quick Orders.

VA-PTSD SCREEN POS MH EMERGENT EVAL

VA-PTSD SCREEN POS MH NON-EMERGENT EVAL

2. Sites may want to check for quick order(s) in the VA-AUD-C POS BRIEF ALCOHOL COUNSELING reminder dialog and have available when installing this dialog. Below is the element you can check for Quick Orders.

VA-ALC REFER TO SPECIALTY SUD PROG

3. Reminder update 37 (UPDATE_2_0_37) and update 53 (UPDATE_2_0_53) should be considered a prerequisite to Update 63.

Verify that you have installed the following update PRIOR to this install. Do this by looking for the entry in reminder exchange and check the Installation History (IH) for the entries:

UPDATE_2_0_53 MH UPDATES

If you need the install guide for UPDATE 53, email me at Alan.Montgomery@va.gov

Install Details

This update is being distributed as a web host file. The address for the host file is:

http://vista.med.va.gov/reminders/UPDATE_2_0_63.PRD

If your site is unable to use LWH, please use LMM for load mail message functionality. You should have been mailed the reminder exchange file when the update was released. If you do not have it in your VistA mail, please contact us.

The file will be installed using Reminder Exchange, programmer access is not required.

Installation:

=====

This update can be loaded with users on the system. Installation will take less than 1 hour.

Install Example

To Load the exchange file, use LWH or if your system doesn't allow use of LWH, use LMM to load the mailman message. Navigate to the message with UPDATE_2_0_63 in the subject and choose it. This will load UPDATE_2_0_63 MH UPDATES 6.0 into your reminder exchange entries. If you are using LWH, the URL is

http://vista.med.va.gov/reminders/update_2_0_63.prd

If you use LWH and it works, you will see something similar to this message. Make sure the file was loaded successfully.

http://vista.med.va.gov/reminders/UPDATE_2_0_63.PRD successfully loaded.

Most components will be the same and will be skipped. Here are general rules for prompts to which you will have to respond.

If a component exists, but is different, use **OVERWRITE** action

If a component is new, use **INSTALL** action

Search and locate an entry titled UPDATE_2_0_63 MH UPDATES 6.0 in reminder exchange. If you were a test site, you may have more than 1 entry with this name. Be sure and choose the entry with date of 02/27/2019@05:43:48

At the **Select Action** prompt, enter **IFE** for Install Exchange File Entry

Enter the number that corresponds with your entry titled UPDATE_2_0_63 MH UPDATES 6.0 (in this example it is entry 141, it will vary by site) that is dated 02/27/2019@05:43:48

+Item	Entry	Source	Date Packed
	NOTE		
140	UPDATE_2_0_63 MH UPDATES 6.0	MONTGOMERY@SALT LAKE CI	01/07/2019@08:17
141	UPDATE_2_0_63 MH UPDATES 6.0	MONTGOMERY@SALT LAKE CI	02/27/2019@05:43
142	UPDATE_2_0_65 VA-MH VETERAN COMMUNITY STATUS NOTE	SMALLEY@SALT LAKE CI	12/20/2018@10:44
143	UPDATE_2_0_66 VA-STAR DIALOGS	GRIFFITH@SALT LAKE CI	02/07/2019@10:20
144	UPDATE_2_0_67 VA-POST 9/11 CASE MANAGEMENT SCREENING	GRIFFITH@SALT LAKE CI	12/27/2018@09:37
145	UPDATE_2_0_68 VA-NUR SUR FIRE RISK ASSESSMENT	PLOTT@SALT LAKE CI	02/22/2019@12:22
+ + Next Screen - Prev Screen ?? More Actions >>>			
CFE	Create Exchange File Entry	LHF	Load Host File
CHF	Create Host File	LMM	Load MailMan Message
CMM	Create MailMan Message	LR	List Reminder Definitions
DFE	Delete Exchange File Entry	LWH	Load Web Host File
IFE	Install Exchange File Entry	RI	Reminder Definition Inquiry
IH	Installation History	RP	Repack
Select Action: Next Screen// IFE Install Exchange File Entry			
Enter a list or range of numbers (1-243): 141			

Component	Category	Exists
Source: MONTGOMERY,ALAN at SALT LAKE CITY		
Date Packed: 02/21/2019@10:40:47		
Package Version: 2.0P35		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-FOLLOW-UP POS PTSD/DEP COMBINATION		
VA-GP POS PTSD/DEP OBJECTS		
REMINDER TAXONOMY		
VA-DEPRESSION DX OUTPT VISIT		
REMINDER TERM		
VA-DEPRESSION SCREEN OEF/OIF		
VA-USER CLASS FOR COGNITIVE IMPAIRMENT		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// IA		

At the **Select Action** prompt, type **IA** for Install all Components and hit enter.

You will see several prompts. For all new entries you will choose **I to Install** :

For Items that exist but is different, use the OVERWRITE action **except for Reminder Terms**

VA-USER CLASS FOR COGNITIVE IMPAIRMENT

VA-DEPRESSION SCREEN OEF/OIF

For these two terms, use the MERGE option

Dialog installation may not necessarily appear in the order listed in the install guide

You will be prompted to install the first reminder dialog component:

```
Packed reminder dialog: VA-GP POS PTSD/DEP OBJECTS
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP POS PTSD/DEP OBJECTS Finding: *NONE*	group	X
2	5	VA-POS PTSD/DEP OBJECTS Finding: *NONE*	element	X

```
+ Next Screen - Prev Screen ?? More Actions
DD Dialog Details DT Dialog Text IS Install Selected
DF Dialog Findings DU Dialog Usage QU Quit
DS Dialog Summary IA Install All
Select Action: Quit// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-GP POS PTSD/DEP OBJECTS**

Select Action: Next Screen// **IA Install All**

All dialog components for **VA-GP PTSD/DEP OBJECTS** Already exist.

Install reminder dialog and all components with no further changes: **Y// YES**

```
Packed reminder dialog: VA-GP POS PTSD/DEP OBJECTS
VA-GP POS PTSD/DEP OBJECTS (group) installed from exchange file.
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-GP POS PTSD/DEP OBJECTS Finding: *NONE*	group	X
2	5	VA-POS PTSD/DEP OBJECTS Finding: *NONE*	element	X

```
+ Next Screen - Prev Screen ?? More Actions
DD Dialog Details DT Dialog Text IS Install Selected
DF Dialog Findings DU Dialog Usage QU Quit
DS Dialog Summary IA Install All
Select Action: Quit// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the second reminder dialog component:

```
Packed reminder dialog: VA-FOLLOW-UP POS PTSD/DEP COMBINATION
```

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog	X		
2	2	VA-DEPRESSION/PTSD FOLLOW UP VERSION Finding: *NONE*	element	X		
3	5	VA-PTSD/DEPRESSION COMBINATION MASTER GP Finding: DEP/PTSD/SI ASSESSMENT COMPLETE (HEALTH FACTOR)	group	X		
4	5.5	VA-GP POS PTSD/DEP HEADER Finding: *NONE*	group	X		
5	5.5.5	VA-GP POS PTSD/DEP OBJECTS Finding: *NONE*	group	X		
6	5.5.5.5	VA-POS PTSD/DEP OBJECTS Finding: *NONE*	element	X		
7	5.6	VA-C-SSRS DONE TODAY Finding: *NONE*	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog	Details	DT	Dialog Text	IS	Install Selected
DF	Dialog	Findings	DU	Dialog Usage	QU	Quit
DS	Dialog	Summary	IA	Install All		

Select Action: Next Screen// IA

At the **Select Action** prompt, type **IA** to install the dialog – **VA-FOLLOW-UP POS PTSD/DEP COMBINATION**

Select Action: Next Screen// **IA Install All**

All dialog components for **VA-FOLLOW-UP POS PTSD/DEP COMBINATION** already exist.

Install reminder dialog and all components with no further changes: **Y// YES**

```
Packed reminder dialog: VA-FOLLOW-UP POS PTSD/DEP COMBINATION
```

Item	Seq.	Dialog Findings	Type	Exists		
1		VA-FOLLOW-UP POS PTSD/DEP COMBINATION	dialog	X		
2	2	VA-DEPRESSION/PTSD FOLLOW UP VERSION Finding: *NONE*	element	X		
3	5	VA-PTSD/DEPRESSION COMBINATION MASTER GP Finding: DEP/PTSD/SI ASSESSMENT COMPLETE (HEALTH FACTOR)	group	X		
4	5.5	VA-GP POS PTSD/DEP HEADER Finding: *NONE*	group	X		
5	5.5.5	VA-GP POS PTSD/DEP OBJECTS Finding: *NONE*	group	X		
6	5.5.5.5	VA-POS PTSD/DEP OBJECTS Finding: *NONE*	element	X		
7	5.6	VA-C-SSRS DONE TODAY Finding: *NONE*	group	X		
+ + Next Screen - Prev Screen ?? More Actions						
DD	Dialog	Details	DT	Dialog Text	IS	Install Selected
DF	Dialog	Findings	DU	Dialog Usage	QU	Quit
DS	Dialog	Summary	IA	Install All		

Select Action: Next Screen// Q

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the third reminder dialog component:

```
Packed reminder dialog: VA-AUD-C POS BRIEF ALCOHOL COUNSELING [NATIONAL DIALOG]
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-AUD-C POS BRIEF ALCOHOL COUNSELING	dialog	X
2	1	VA-POS ALCOHOL FOLLOW UP VERSION Finding: *NONE*	element	
3	3	VA-AUDC SCORE=0 Finding: *NONE*	element	X
4	4	VA-GP ALCOHOL FOLLOW-UP 0 Finding: *NONE*	group	X
5	4.5	TEXT AUDC PRIOR RESULTS Finding: *NONE*	element	X
6	4.10	VA-GP ALCOHOL AUDIT C POSITIVE Finding: *NONE*	group	X
7	4.10.10	VA-GP ALCOHOL AUDIT C 5-7*	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    OU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-AUD-C POS BRIEF ALCOHOL COUNSELING**

Select Action: Next Screen// **IA Install All**

All dialog components for **VA- AUD-C POS BRIEF ALCOHOL COUNSELING** already exist.

Install reminder dialog and all components with no further changes: **Y// YES**

```
VA-AUD-C POS BRIEF ALCOHOL COUNSELING (reminder dialog) installed from exchange
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-AUD-C POS BRIEF ALCOHOL COUNSELING	dialog	X
2	1	VA-POS ALCOHOL FOLLOW UP VERSION Finding: *NONE*	element	X
3	3	VA-AUDC SCORE=0 Finding: *NONE*	element	X
4	4	VA-GP ALCOHOL FOLLOW-UP 0 Finding: *NONE*	group	X
5	4.5	TEXT AUDC PRIOR RESULTS Finding: *NONE*	element	X
6	4.10	VA-GP ALCOHOL AUDIT C POSITIVE Finding: *NONE*	group	X
7	4.10.10	VA-GP ALCOHOL AUDIT C 5-7*	group	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD  Dialog Details      DT  Dialog Text      IS  Install Selected
DF  Dialog Findings    DU  Dialog Usage    OU  Quit
DS  Dialog Summary     IA  Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the fourth reminder dialog component:

```
Packed reminder dialog: VA-PTSD SCREENING
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PTSD SCREENING	dialog	X
2	1	VA-PTSD SCREENING VERSION Finding: *NONE*	element	X
3	10	VA-PTSD SCREEN Finding: *NONE*	group	X
4	10.10	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5+I9 (MH TESTS AND SURVEY)	element	X
5	PXRM	PC-PTSD-5+I9 RESULT GROUP	rs.group	X
6	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 1	rs.element	X
7	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 2	rs.element	X
8	PXRM	PC-PTSD-5+I9 RESULT GROUP SRE	rs.group	X
9	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 1 SRE	rs.element	X
10	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 2 SRE	rs.element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details      DT Dialog Text      IS Install Selected
DF Dialog Findings    DU Dialog Usage    QU Quit
DS Dialog Summary     IA Install All
Select Action: Next Screen// IA
```

At the **Select Action** prompt, type **IA** to install the dialog – **VA-PTSD SCREENING**

Select Action: Next Screen// **IA Install All**

All dialog components for **VA-PTSD SCREENING** already exist.

Install reminder dialog and all components with no further changes: **Y// YES**

```
Packed reminder dialog: VA-PTSD SCREENING
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-PTSD SCREENING	dialog	X
2	1	VA-PTSD SCREENING VERSION Finding: *NONE*	element	X
3	10	VA-PTSD SCREEN Finding: *NONE*	group	X
4	10.10	VA-MH PC-PTSD+I9 SCREEN Finding: PC-PTSD-5+I9 (MH TESTS AND SURVEY)	element	X
5	PXRM	PC-PTSD-5+I9 RESULT GROUP	rs.group	X
6	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 1	rs.element	X
7	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 2	rs.element	X
8	PXRM	PC-PTSD-5+I9 RESULT GROUP SRE	rs.group	X
9	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 1 SRE	rs.element	X
10	PXRM	PC-PTSD-5+I9 RESULT ELEMENT 2 SRE	rs.element	X

```
+      + Next Screen  - Prev Screen  ?? More Actions
DD Dialog Details      DT Dialog Text      IS Install Selected
DF Dialog Findings    DU Dialog Usage    QU Quit
DS Dialog Summary     IA Install All
Select Action: Next Screen// Q
```

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.

You will be prompted to install the fifth reminder dialog component:

```
Packed reminder dialog: VA-DEPRESSION SCREEN
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-DEPRESSION SCREEN	dialog	X
2	1	VA-DEPRESSION SCREENING VERSION Finding: *NONE*	element	X
3	5	VA-DEP PHQ2 Finding: *NONE*	group	X
4	5.10	VA-MH PHQ-2+I9 Finding: PHQ-2+I9 (MH TESTS AND SURVEY)	element	X
5		PXRM PHQ-2+I9 RESULT GROUP	rs.group	X
6		PXRM PHQ-2+I9 RESULT ELEMENT 1	rs.element	X
7		PXRM PHQ-2+I9 RESULT ELEMENT 2	rs.element	X
8		PXRM PHQ-2+I9 RESULT GROUP SRE	rs.group	X
9		PXRM PHQ-2+I9 RESULT ELEMENT 1 SRE	rs.element	X
10		PXRM PHQ-2+I9 RESULT ELEMENT 2 SRE	rs.element	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	
Select Action: Next Screen// IA				

At the **Select Action** prompt, type **IA** to install the dialog – **VA-DEPRESSION SCREENING**
 Select Action: Next Screen// IA **Install All**
 All dialog components for **VA-DEPRESSION SCREENING** already exist.
 Install reminder dialog and all components with no further changes: **Y// YES**

```
Packed reminder dialog: VA-DEPRESSION SCREEN
```

Item	Seq.	Dialog Findings	Type	Exists
1		VA-DEPRESSION SCREEN	dialog	X
2	1	VA-DEPRESSION SCREENING VERSION Finding: *NONE*	element	X
3	5	VA-DEP PHQ2 Finding: *NONE*	group	X
4	5.10	VA-MH PHQ-2+I9 Finding: PHQ-2+I9 (MH TESTS AND SURVEY)	element	X
5		PXRM PHQ-2+I9 RESULT GROUP	rs.group	X
6		PXRM PHQ-2+I9 RESULT ELEMENT 1	rs.element	X
7		PXRM PHQ-2+I9 RESULT ELEMENT 2	rs.element	X
8		PXRM PHQ-2+I9 RESULT GROUP SRE	rs.group	X
9		PXRM PHQ-2+I9 RESULT ELEMENT 1 SRE	rs.element	X
10		PXRM PHQ-2+I9 RESULT ELEMENT 2 SRE	rs.element	X
+ + Next Screen - Prev Screen ?? More Actions				
DD	Dialog Details	DT	Dialog Text	IS Install Selected
DF	Dialog Findings	DU	Dialog Usage	QU Quit
DS	Dialog Summary	IA	Install All	
Select Action: Next Screen// Q				

When the dialog has completed installation, you will then be returned to this screen. At the **Select Action** prompt, type **Q**.
 This will complete the Dialog Install section.

Component	Category	Exists
Source: MONTGOMERY,ALAN at SALT LAKE CITY		
Date Packed: 02/27/2019@05:43:48		
Package Version: 2.0P35		
Description:		
The following Clinical Reminder items were selected for packing:		
REMINDER DIALOG		
VA-FOLLOW-UP POS PTSD/DEP COMBINATION		
VA-GP POS PTSD/DEP OBJECTS		
REMINDER TAXONOMY		
VA-DEPRESSION DX OUTPT VISIT		
REMINDER TERM		
VA-DEPRESSION SCREEN OEF/OIF		
VA-USER CLASS FOR COGNITIVE IMPAIRMENT		
+ Enter ?? for more actions >>>		
IA	Install all Components	IS Install Selected Component
Select Action: Next Screen// Q		

Type Q to exit the dialog installation

Post Installation

1. Add any Quick Orders to your dialogs that are needed. See Pre-Install section for names of Dialog components to add Quick Orders as finding items.
2. You will need to re-link the dialog you were using prior to this install to the national alcohol screening reminder. Link either the Alcohol Screening or the Alcohol Screening with follow up to your Alcohol screening reminder

Install complete.

NOTE: If there are any questions or problems during the installation, please notify your local IT support to log a national help desk ticket or contact the VA Service Desk at 1-888-596-4357 and have them submit a national ticket to:

Category: Enterprise Applications

Subcategory: Applications

Enterprise Application: VistA - Clinical Reminders