



# **Dental Record Manager Plus (DRM Plus)**

## **User Manual**

---

**Version 7.7**

**May 2020**



© 2020 Document Storage Systems, Inc. All rights reserved.

*Document Storage Systems (DSS) is a privately held corporation and has been the premier provider of health information and document imaging distribution and storage systems to Veterans Affairs facilities for over twenty-five years. DSS is located at 12575 US Highway One, Suite 200, Juno Beach, Florida 33408.*

*World Wide Web: [www.dssinc.com](http://www.dssinc.com)*

*Dental Record Manager Plus is a trademark of Document Storage Systems, Inc. Outlook, Internet Explorer, and Windows are trademarks of Microsoft Corporation. VistA is a trademark of the Department of Veterans Affairs, Software Services, and Computerized Patient Record System.*

*No portion of this manual or software may be reproduced without the prior written consent of Document Storage Systems, Inc.*

# Table of Contents

<b>Introduction .....</b>	<b>8</b>
Document Storage Systems, Inc.....	8
From the Department of Veterans Affairs .....	8
Introduction.....	8
Quality Improvement/Performance Measures and Benefits .....	9
Customer Support .....	9
<b>Dental Record Manager Plus User and Administrator Requirements.....</b>	<b>10</b>
DRM Plus User Requirements .....	10
Administrator Option .....	11
<b>Accessing DRM Plus .....</b>	<b>12</b>
Access .....	12
<b>Using the DRM Plus Drop-Down Menus .....</b>	<b>16</b>
File.....	16
<i>Refresh Patient Chart.....</i>	<i>16</i>
<i>File Administrative Time .....</i>	<i>16</i>
<i>File Fee Basis .....</i>	<i>17</i>
<i>Print .....</i>	<i>18</i>
<i>Spell Check .....</i>	<i>19</i>
<i>Save Unfiled Data .....</i>	<i>20</i>
<i>Exit.....</i>	<i>21</i>
Edit .....	22
<i>Copy.....</i>	<i>22</i>
<i>Cut .....</i>	<i>22</i>
<i>Paste .....</i>	<i>22</i>
<i>Select All .....</i>	<i>22</i>
Dental Encounter Data .....	23
<i>Create New PCE Visit.....</i>	<i>23</i>
<i>View Scheduled Appointments and Historical Visits .....</i>	<i>24</i>
Treatment & Exam.....	27
<i>Show Configuration .....</i>	<i>28</i>
<i>Add/Edit Personal QuickList.....</i>	<i>38</i>
<i>Add Medical Codes to ADA Table.....</i>	<i>38</i>
<i>Edit Code Information in the ADA Table .....</i>	<i>39</i>
<i>Edit Procedure Costs .....</i>	<i>39</i>
<i>Filter View .....</i>	<i>39</i>
<i>Clean Slate.....</i>	<i>40</i>
<i>Undo Clean Slate .....</i>	<i>40</i>
<i>All Planned Care to Clipboard .....</i>	<i>41</i>
Tools .....	42
<i>Windows Calculator.....</i>	<i>42</i>
<i>Windows Explorer.....</i>	<i>42</i>
<i>Windows Notepad .....</i>	<i>42</i>
<i>User Inquiry.....</i>	<i>43</i>
<i>User Options.....</i>	<i>45</i>
<i>Administrative Toolbox.....</i>	<i>61</i>
<i>Panel Add/Edit.....</i>	<i>61</i>
<i>Provider Add/Edit.....</i>	<i>61</i>
<i>Ancillary Tool Functions – ADA Website.....</i>	<i>61</i>
Reports.....	62

<i>Reports</i> .....	62
<i>Service Reports</i> .....	72
<i>Data Warehouse Reports</i> .....	74
<i>Adverse Events</i> .....	75
<i>Device Tracking</i> .....	75
<i>Extract History File</i> .....	75
<i>Queued Extract History File</i> .....	75
Help .....	75
<i>Contents</i> .....	76
<i>Version Release Notes</i> .....	77
<i>Last Broker Call</i> .....	78
<i>VA Intranet Website</i> .....	78
<i>Find your DRM Plus Administrators</i> .....	78
<i>Have a Question, Comment or Suggestion about DRM?</i> .....	78
<i>About</i> .....	79
<b>DRM Plus Banner</b> .....	<b>81</b>
Patient Information .....	81
Visit Information.....	82
Dental Provider Information .....	83
Vitals Button.....	83
Adverse Events Button.....	84
<i>To add an adverse event:</i> .....	84
<i>To edit a filed Active adverse event:</i> .....	85
<i>To delete a filed adverse event:</i> .....	86
Device Tracking Button .....	88
<i>To add completed procedure device identifiers:</i> .....	89
<i>To add non-transaction device identifiers:</i> .....	91
<i>To remove filed device identifiers for a transaction or non-transaction:</i> .....	93
<i>To delete the filed device identifiers for a transaction or non-transaction:</i> .....	93
Dental Class Information .....	95
Clean Slate .....	95
Icons .....	95
<i>General Coding Standards</i> .....	96
Patient Flags.....	96
<i>Clinical Reminders</i> .....	96
<i>Consult</i> .....	96
<i>Exam Quality Indicator</i> .....	96
<i>Fluoride Quality Indicator</i> .....	97
<i>Alerts</i> .....	98
<b>Cover Page</b> .....	<b>100</b>
Dental Eligibility.....	101
<i>Dental Class</i> .....	101
<i>Service Connected Teeth/Service Trauma</i> .....	101
<i>Adjunctive Medical Condition(s)</i> .....	101
<i>Eligibility Expiration Date</i> .....	101
<i>Anticipated Rehab Date</i> .....	101
Demographics .....	102
Case Management.....	102
<i>Status</i> .....	102
<i>Suggested Recare Date</i> .....	102
Recent Dental Activity .....	103

Fluoride Indicator Prescription Date .....	103
Dental Alerts .....	104
Notes .....	104
Filed Planned Care .....	105
<b>Clinical Record.....</b>	<b>106</b>
Problems .....	108
Consultations .....	109
Notes .....	109
<i>Adding a New TIU Progress Note .....</i>	<i>110</i>
<i>Adding a New TIU Progress Note Addendum .....</i>	<i>112</i>
<b>Dental History .....</b>	<b>113</b>
Viewing Dental Information by Tooth .....	114
Viewing Other Dental History Information.....	114
Viewing All Dental History Information .....	114
Viewing Dental History Information by Episode of Care .....	115
Episode in Date Range .....	115
<b>Chart/Treatment – Treatment &amp; Exam.....</b>	<b>116</b>
Diagnostic Findings .....	118
<i>Editing Diagnostic Finding Descriptions.....</i>	<i>119</i>
<i>Deleting a Diagnostic Finding.....</i>	<i>119</i>
Treatment Plan .....	120
<i>Entering a Treatment Plan.....</i>	<i>120</i>
<i>Editing a Treatment Plan Description.....</i>	<i>122</i>
<i>Deleting a Treatment Plan .....</i>	<i>123</i>
<i>Completing a Treatment Plan .....</i>	<i>123</i>
Completed Care .....	124
<i>Entering Completed Care.....</i>	<i>124</i>
<i>Editing Completed Care Description .....</i>	<i>125</i>
<i>Deleting a Completed Care.....</i>	<i>125</i>
Include “Completed”/Include “Findings and Completed”/Include “Findings” .....	125
Perio Buttons Icon.....	126
Seq Plan/Sequencing Button .....	127
<i>Plan a Treatment Sequence.....</i>	<i>128</i>
<i>Complete a Planned Treatment in the Sequencing Screen .....</i>	<i>128</i>
<i>Deleting a Planned Treatment in the Sequencing Screen.....</i>	<i>128</i>
Chart Hx (History) Button .....	129
Summary Button .....	130
H&N Button.....	131
PSR Button .....	134
OHA (Oral Health Assessment) Button .....	135
TMJ Button.....	137
Occl (Occlusion) Button .....	138
Habits (Parafunctional) Button.....	140
Social Hx (Social History) Button .....	141
<i>Multiple Filings to Same Modal on Same Day.....</i>	<i>142</i>
<b>Chart/Treatment – Periodontal Chart .....</b>	<b>144</b>
History and Compare Buttons .....	145
Summary Button .....	146
H&N Button.....	146
PSR Button .....	146
Stats Button.....	147

OHA (Oral Health Assessment) Button .....	147
Notes Button .....	148
Entering Periodontal Information.....	149
Other Tools .....	150
<b>Exam.....</b>	<b>151</b>
Exam Elements .....	153
<i>Presentation/Chief Complaint Element.....</i>	<i>153</i>
<i>Vitals Elements.....</i>	<i>154</i>
<i>PMH (Past Medical History) and Medications Element .....</i>	<i>155</i>
<i>Social History Element.....</i>	<i>156</i>
<i>H&amp;N (Head and Neck) Findings Element.....</i>	<i>157</i>
<i>Radiographic Findings Element.....</i>	<i>158</i>
<i>Diagnostic Findings Element.....</i>	<i>159</i>
<i>Periodontal Assessment Element.....</i>	<i>160</i>
<i>Parafunctional Habits Element.....</i>	<i>163</i>
<i>TMJ Findings Element .....</i>	<i>164</i>
<i>Occlusal Findings Element .....</i>	<i>165</i>
<i>Salivary Flow Element .....</i>	<i>166</i>
<i>Removable Prostheses Element.....</i>	<i>167</i>
<i>Assessment/Plan Element.....</i>	<i>168</i>
<i>Patient Education Element.....</i>	<i>169</i>
<i>Disposition Element .....</i>	<i>170</i>
<i>Import Previously Filed Data Screen.....</i>	<i>171</i>
<b>Completing the Encounter .....</b>	<b>173</b>
Potential Duplicate Transactions Screen .....	177
File Data Option Screen .....	178
<i>File to PCE/DES with Code.....</i>	<i>178</i>
<i>File to DES-Only Data.....</i>	<i>179</i>
Filing Options Screen .....	179
<i>Filing Options.....</i>	<i>179</i>
<i>Visit Date/Time .....</i>	<i>179</i>
<i>Encounter Dental Class .....</i>	<i>180</i>
<i>Disposition.....</i>	<i>180</i>
<i>Suggested Recare Date .....</i>	<i>180</i>
<i>Primary PCE Diagnosis &amp; Send Dx to CPRS Problem List.....</i>	<i>181</i>
<i>Service Connection .....</i>	<i>181</i>
<i>Additional Providers/Additional Signers.....</i>	<i>182</i>
<i>Station.....</i>	<i>183</i>
Progress Note Screen .....	183
<i>Viewing/Importing DRM Object/Progress Note.....</i>	<i>184</i>
<i>Viewing/Importing CPRS Templates .....</i>	<i>184</i>
<i>Importing VistA Medical Information .....</i>	<i>184</i>
<i>Other Options in the Import Menu .....</i>	<i>185</i>
<i>Accessing Dental CNTs.....</i>	<i>185</i>
<i>Electronic Signature .....</i>	<i>185</i>
<i>Progress Note Addendum.....</i>	<i>186</i>
CNT Navigator .....	187
<i>Navigating Within CNTs .....</i>	<i>188</i>
Consult Notes.....	190
Resident Filing as Cosigners or Distributed Providers .....	192
<b>Appendix A – Glossary of VA Terms .....</b>	<b>194</b>

<b>Appendix B – Common Application Functions .....</b>	<b>198</b>
<b>Appendix C – Hints and Notes.....</b>	<b>200</b>
Save Unfiled Data .....	200
Dental Class Displayed on Banner .....	200
Diagnostic Findings .....	200
Treatment Plan .....	200
Multi-Add Screen .....	201
Ranged Codes .....	201
Speed Codes.....	201
Tx Planning/Sequencing Screen.....	201
Completed Care .....	202
Periodontal Chart .....	202
Completing the Encounter.....	203
Reports – Non-Clinical Time by Provider.....	205
Code Boilerplates.....	205
Last Broker Call.....	205
Recent Dental Activity .....	205
<b>Appendix D – Icon Definitions.....</b>	<b>206</b>
Diagnostic Findings .....	206
Treatment Plan .....	209
Completed Care .....	211
Special Descriptions – Bridge Icon .....	213
Special Descriptions – Conn Bar Icon.....	214
Special Descriptions – Notes Icon.....	214
<b>Appendix E – Using the Keyboard to Enter Periodontal Data .....</b>	<b>215</b>
Overview .....	215
<i>Navigating the Periodontal Screen .....</i>	<i>215</i>
<i>Arch Views .....</i>	<i>215</i>
<i>Cursor Movement.....</i>	<i>215</i>
<i>Entering Data.....</i>	<i>216</i>
<i>Special Buttons.....</i>	<i>216</i>
<i>Other Functions .....</i>	<i>216</i>
<b>Appendix F – Ranged Codes .....</b>	<b>217</b>
<b>Appendix G – Option to Set Dental Patients to Inactive Status .....</b>	<b>218</b>
<b>Appendix H – How to Map Dental CNTs.....</b>	<b>219</b>
<b>Appendix I – Recommendations for Coding of Prosthetic Appliance.....</b>	<b>221</b>
<b>Appendix J – Business Use of DRM Plus.....</b>	<b>223</b>
Local Policy and Practice .....	223
National Policy and Practice Coding Standards .....	223
<b>Appendix K – Data Security .....</b>	<b>224</b>
<b>Appendix L – madExcept Application .....</b>	<b>224</b>
<b>Periodontal Keyboard Shortcuts Tear-Out.....</b>	<b>224</b>

# Introduction

## **Document Storage Systems, Inc.**

DSS, Inc. specializes in the computerization of patient medical charts. Our core specialty within the medical market is building Windows Graphical User Interface (GUI) applications, which insert, update and retrieve patient data held in a MUMPS (M) data repository, or SQL database system. DSS offers an array of GUI products, which allow for the electronic documentation of TIU progress notes and other significant parts of medical records, scanning and viewing of clinical and administrative documents and automated medical record coding through simple points and clicks.

## **From the Department of Veterans Affairs**

Dental Record Manager Plus (DRM Plus) captures specific dentally related information elements not readily available in CPRS. These elements include oral cavity/tooth related diagnostic findings, dental-specific care plans and a superset of completed care information. DRM Plus aids the provider in the entry of dental diagnostic information, coding and crediting dental procedures, completing TIU progress notes, and planning and tracking dental patient care. DRM Plus is adjunctive to CPRS and is NOT designed to replace CPRS for dental users. While some information from CPRS is available, and can be accessed in DRM plus, providers should use all the available tools in the VistA suite of applications. These tools include: VistA Imaging, I-Med Consent, and any clinical system applications specific to the local sites. DRM Plus is a Dental Graphical User Interface front end for data input into the VistA Dental files, as well as the Patient Care Encounter (PCE), Text Integration Utility (TIU) and CPRS Problem List packages.

## **Introduction**

The DRM Plus program is designed to provide dental health care facilities with an intuitive, user-friendly Windows interface for end-users to create encounter information, evaluate patient dental conditions, and develop and maintain the treatment plan. The DRM Plus program is an application that uses RPC Broker technology, which permits the facility users to store and retrieve clinical data within the VistA System.

DRM Plus supports the Veterans Health Administration, Office of Dentistry, continuous quality improvement initiatives by providing added value to the clinical and administrative management of the patient's electronic dental record. The enhanced methods of data capture included in this application continue to eliminate unnecessary paperwork and administrative functions through the automation of clinical dental data.

The use of DRM Plus results in more accurate insurance billing for dental visits, consults and procedures. This application supports the filing of Dental Encounter System (DES) within the guidelines established by the Veterans Health Administration, Office of Dentistry.



Some features of DRM Plus are summarized in the following:

- Entry of dental conditions, plans, and completed procedures using graphic icons with extensive use of color schemes.
- Upper/Lower/Full Views with full color-coded graphics.
- Sequencing of Treatment Plan procedures
- Dental History with date-change capability
- Quadrant or Tooth summaries
- Head/Neck Findings availability
- Periodontal charting
- Full Mouth Plaque Index with definitions
- ADA/Local/Quick Codes
- Creation and maintenance of tooth-specific and general patient notes.

## **Quality Improvement/Performance Measures and Benefits**

DRM Plus supports the VA Administration, Office of Dentistry's continuous quality improvement initiatives by providing added value to the clinical and administrative management of the patient's electronic dental record. The efficient data capture methods included in this product eliminates unnecessary paperwork and administrative functions. Additional quality improvement benefits and sample performance measures include:

- Performances Measures
- Reductions in operating cost and improved services through better integration of VHA resources and data.
- Supports high level job satisfaction by providing clinicians with feedback regarding quality of care and promotes a culture that places a high value on individual and collective accountability through reporting.
- Promotes a VHA culture of ongoing quality improvement that is predicated on providing excellent health care value.
- Accuracy and usefulness of data increases based on the reduction of data entry points and decreased potential for error.
- Enhanced capability to measure quality of care consistent with the VA Dentistry GPRA Performance Plan.

## **Customer Support**

DRM Plus is supported in the same manner as any other nationally supported software product. Problems should be reported to the local site ADPAC and/or Help Desk, who in turn utilize the Computer Associate's Service Desk Manager (CA SDM) system to log and track problems. Help desk support is provided from 8:00 AM to 7:00PM Eastern Standard Time, Monday through Friday. Documenting problems provides a means to find and disseminate solutions to those involved in any area of DRM Plus or VistA.

# Dental Record Manager Plus User and Administrator Requirements

## DRM Plus User Requirements

1. DRM Plus users must have a valid **Person Class** in VistA file 200 (New Person File) to file data in DRM Plus.
  - a. ALL residents and fellows must have one of the following three person classes: **V030300**, **V115500**, or **V115600**, and this requires the resident to select a distributive provider (attending) as the primary provider when filing to DES and PCE.
  - b. DRM Plus users must have a **Person Class** (different than the three listed above) to file data in DRM Plus. DRM Plus users that do NOT have a Person Class will receive an “!”-screen (Informational) stating the user is required to have an active Person Class to file data in DRM Plus. Please contact IT support for assignment of Person Class.
  - c. DRM Plus users without a **Person Class** assigned may file an unsigned encounter for another provider or save unfiled data for another provider.
2. Dental Residents and Fellows using DRM Plus must have a valid **User Class** in VistA TIU if they require a cosigner. Please refer to the Authorization/Subscription Utility (ASU) User Manual to ensure that Dental Residents are prompted for cosignature on all Progress Notes and Consults. This is typically done by a Clinical Coordinator, IT Staff, ADPAC or other manager using the Document Parameter Edit option on the TIU Parameters Menu on the IRM Maintenance Menu. The **USERS REQUIRING COSIGNATURE** field within the Document Parameter Edit option indicates which groups of users (i.e., User Classes) require cosignature for the type of document in question.
3. All DRM Plus users filing TIU progress notes must have an **electronic signature** in VistA file 200 (New Person File).
4. All DRM Plus users must have a **default division (station number)** in VistA file 200 (New Person File) for the station number to appear as the default (preselected) when filing an encounter.
5. All DRM Plus users must have the secondary menu option **DENTV DSS DRM GUI** assigned to access DRM Plus.
6. All DRM Plus users filing encounters should have their **initials defined** in the new person file (200) so that they will appear in the DRM Plus transaction tables.
7. All DRM Plus users (except Cover Page Only users) must have an active **8-Digit Dental Provider ID** in the VistA Dental Provider File (220.5) to open DRM Plus and file data. DRM Plus users that do NOT have an **8-Digit Dental Provider ID** will receive a “**Red X**” screen (Stop) stating that they are required to have an active Dental Provider ID and will be denied access to DRM Plus. Note to DRM Plus application administrators: Use the DRM Plus option **Provider Add/Edit** on the **Tools** menu to enter this information. Select the **New** button if the username does NOT display in the list.
8. Dental Students should NOT be assigned a **Person Class** in VistA; they only need an **8-digit Provider ID** and require a **User Class** in VistA TIU (student). The **User Class** allows them to access CPRS. The dental student may file data for another provider.
9. DRM Plus users that have **Cover Page Only** access require the secondary menu option **DENTV DSS DRM GUI** assigned and no other requirements to access DRM Plus. Cover Page Only access must be granted by a DRM Plus application administrator.

**Note:** DRM Plus users must have permission to write and modify to-and-from the **DOCSTORE** folder.

## **Administrator Option – (Access to the DENTV XPAR EDIT PARAMS option in VistA required)**


All DRM Plus application administrators must be DRM Plus users. Enter VistA in Programmer Mode by typing D ^XUP at the VistA prompt and get to "Select Option".

1. At Select Option Name, type: **DENTV XPAR EDIT PARAMS**.
2. At the Select PARAMETER DEFINITION NAME, type: **DENTV DRM ADMINISTRATOR**.
3. At the Select NEW PERSON NAME, type in the name of the person to be made a dental administrator.
4. Set the value to **YES**.

# Accessing DRM Plus

## Access

To access DRM Plus, first open CPRS and select the desired patient. Open the **Tools** menu in the CPRS tool bar and select the **DRM Plus** submenu from the available submenus.



**Figure 1: Access DRM Plus through CPRS**

DRM Plus opens with the patient information loaded and, unless changed by the user, the **Chart/Treatment** tab as the default opening screen.

**Note:** Users may be required to re-enter Access/Verify codes when opening DRM Plus. The default opening settings of DRM Plus is the **Treatment Plan** screen on the **Chart/Treatment** tab, unless changed by the user.

**Note:** The proper ways to close DRM Plus are listed:

1. Selecting the [X] button in the upper right corner of the DRM Plus screen; or
2. Selecting **File** menu → **Exit** submenu; or
3. Selecting **Task Manager** → **Processes** tab → highlight **DRM Plus** task → right click select the **End Task** submenu.




Figure 2: DRM Plus Chart/Treatment Tab

The “!”- screen (Informational) displays when a DRM Plus user does NOT have an active VistA **Person Class**. This user will need to contact the local Help Desk and request an active VistA **Person Class**. One of the requirements to file data in DRM Plus is to have an active VistA **Person Class**. However, a DRM Plus user with no VistA **Person Class** may file an unsigned encounter for another provider or save unfiled data for another provider.




Figure 3: “!”-screen (Informational)

The “Red X”- screen (Stop) displays when a DRM Plus user does NOT have an active **8-Digit Dental Provider ID**. No access will be allowed to open DRM Plus by the user. A DRM Plus Administrator may create an active **8-Digit Dental Provider ID** from the **Provider Add/Edit** submenu from the **Tools** menu in DRM Plus. One of the requirements to file data in DRM Plus is to have an active **8-Digit Dental Provider ID**.




Figure 4: “Red X” –screen (Stop)

**Note:** No informational warning screen displays to a user with **Cover Page** tab access only.

The following informational screen will display when someone opens a DRM Plus patient’s chart that someone else has open. The informational screen displays the name of the user who first opened the patient’s chart. The informational screen also asks the user to only view the chart until the first user accessing the patient chart closes it, nothing should be entered, edited or deleted. The user receiving the warnings should refresh the patient’s chart or close and reopen DRM Plus after the other user has finished.

There can be issues in DRM Plus when filing entries such as unfiled data or patient’s TIU progress notes when more than one person is accessing the same patient chart.





Figure 5: Warning Informational Screen

In the following pages, the various parts of DRM Plus are highlighted, and the functionality of the program is explored. The main screen is broken into three distinct parts. The drop-down menus allow the user to access various menus throughout the program, regardless of which tab is in use. Some drop-down menu functions are NOT available with every different tab. In this case, the menu function is disabled when the tab is open.




**Figure 6: DRM Plus Drop-Down Menus**

The banner contains patient, visit/location, provider/patient information and limited vitals entry. There are also the adverse events button, device tracking button, dental class button, coding standards and alerts icons on the banner.



**Figure 7: DRM Plus Banner**

The tabs are the heart of DRM Plus. They allow the user to create a new exam template, new treatment plan, view the dental history of a patient, view clinical records, and create a text note or a text note addendum. All providers may perform myriad tasks by simply clicking through each of the tabs and adding the pertinent information that is allowed in the appropriate place.



**Figure 8: DRM Plus Tabs**

The following chapters explore the functionality of each of the areas of the program in detail.

# Using the DRM Plus Drop-Down Menus

The DRM Plus drop-down menus consist of seven menus: **File, Edit, Dental Encounter Data, Treatment & Exam, Tools, Reports and Help.**




Figure 9: DRM Plus Drop-Down Menus

## File

The **File** menu contains seven submenus: **Refresh Patient Chart, File Administrative Time, File Fee Basis, Print, Spell Check, Save Unfiled Data** and **Exit**. The **Spell Check** submenu is only active in the note and note addendum screens.




Figure 10: DRM Plus File Menu

## Refresh Patient Chart

The **Refresh Patient Chart** submenu allows DRM Plus users to refresh the patient's chart while working in DRM Plus.

## File Administrative Time

When the **File Administrative Time** submenu is selected from the **File** menu, the **File Administrative Time** screen displays.



**Figure 11: File Administrative Time Screen**

1. Use the drop-down menu near the top of the screen to select the desired **Station Number**.
2. Click the appropriate **radio button** to select the type of administrative time.
3. Use the **up and down arrows** next to the hours and minutes text boxes to adjust how much time is recorded. Note that the minutes can only be adjusted in 15-minute increments.
4. Click the **OK** button. The screen closes and files that administrative time for report usage.

**Note:** This filing of administrative time is for local use only and does NOT file to the VA-MCA Labor Mapping Access Database Program.

## File Fee Basis

When the **File Fee Basis** submenu is selected, the **Dental Record Manager Fee Basis** screen displays.

**Figure 12: Dental Record Manager Fee Basis Screen**

1. Use the **Report Date** drop-down menu to select a date to edit/delete a previous Fee Basis entry.
2. Choose the station by clicking the appropriate **radio button**.
3. Click the **Dental Category** drop-down menu to choose a **Dental Class**.
4. Click the **Date Authorized for Payment** drop-down menu to display a calendar. The user may toggle through this calendar to choose the authorized date for payment.
5. Enter the **Total Cost** in the text box.
6. Click the **Finish** button.
7. A screen displays stating that a Fee Basis record has been added. Click the **OK** button.

End-user criteria required to allow entering fee basis data within DRM Plus includes:


- Does need to be in the **Dental Provider file** (8-digit provider ID).
- Does NOT need a **Person Class** in VistA.
- Does need access to **CPRS**.
- Does need access to DRM Plus (**DENTV DSS DRM GUI** secondary menu option).
- Does NOT need DRM Plus **administrative access**.

**Note:** DRM Plus Administrators can run all **Fee Basis** reports. Patient care provided by fee should be entered in DRM Plus as **Diagnostic Findings**.

**Note:** Fee basis data entered in DRM Plus is only available for local reports created in DRM Plus.

## Print

Select the **Print** button to view the **Print** screen.



**Figure 13: Print Screen**

Select the check boxes that correspond to what is to be printed.

## Spell Check

Select the **Spell Check** submenu to correct possible spelling errors. This feature is only active in note and note addendum screens.




Figure 14: Spelling Screen

The program goes through the text and highlights words that may have been misspelled and suggests possible correct spellings. Use the buttons to **Ignore**, **Change** or **Add** words. Click the **Options** button to select various options, pick a language/dictionary or add a custom dictionary.




Figure 15: Spelling Options Screen

Click the check boxes beside the desired options located on the **Spelling** tab. Select the language and dictionaries from the **Language** tab and click the **OK** button. The **Spelling Options** screen closes.

## Save Unfiled Data

Select the **Save Unfiled Data** submenu. The **Save DRM Plus Data** screen displays. Click the **Yes** button to save the unfiled data to the listed provider. A screen displays. Click the **Yes** button again to confirm.




Figure 16: Save DRM Plus Data Screen

To change the save unfiled data to another provider, click the **Provider...** button.





Figure 17: Search for Provider Screen

1. Enter the name or partial name of the desired provider in the **search criteria** text box.
2. Press the <Enter> key.
3. Click the needed provider from the list of results.
4. Click the **OK** button to change the provider. The **Save DRM Plus Data** screen displays.
5. Click the **Yes** button to save the unfiled data to the new provider.

When a dental provider is saving unfiled data and assigning it to another dental provider for a selected patient who has previously saved unfiled data that has NOT been filed, the following screen displays.



**Figure 18: Provider Already Has Unfiled Data**


If the user clicks the **Yes** button, previously saved unfiled data originally saved by another dental provider, or this provider, is overwritten. Only one unfiled data entry may be maintained by a single provider, per patient.

## **Exit**

Exit the program by selecting **Exit** from the **File** menu. The CPRS main screen displays.

## **Edit**

The **Edit** menu consists of four submenus: **Copy**, **Cut**, **Paste** and **Select All**.



**Figure 19: Edit Menu**

## **Copy**

To copy, highlight the desired text and choose **Copy**.

## **Cut**

To cut, highlight the desired text and choose **Cut**. The selected text is removed.

## **Paste**

To paste, move the cursor to the area where the copy or cut text is to be replaced. Select the **Paste** submenu to add the text to the chosen area.

## **Select All**

**Select All** submenu highlights all the text visible on the screen which can be copied and/or cut. Use the **Copy** or **Cut** function to complete the desired task.

## Dental Encounter Data

The **Dental Encounter Data** menu has two submenus: **Create New PCE Visit** and **View Scheduled Appointments and Historical Visits**.




Figure 20: Dental Encounter Data Menu

### Create New PCE Visit

Select the **Create New PCE Visit** submenu to display the **Provider and Location for Current Activities** screen.

**Note:** The **Create New PCE Visit** submenu is only available if the DRM Plus Administrator allows new PCE visits to be created in the DRM Plus application. The opening default tab is the **New Visit** tab.

The screenshot shows the 'Provider and Location for Current Activities' window. It has two main sections: 'Encounter Provider' and 'Encounter Location'. The 'Encounter Provider' section has a list box with 'DRMPROVIDER,DENTIST STAFF DENTIST,PROSTHODONTICS' selected. The 'Encounter Location' section has a text field with 'Jun 07, 17@11:00 DENTAL'. Below these are tabs for 'My Clinic Visits', 'Dental Visits', 'All Visits', 'Admissions', and 'New Visit'. The 'New Visit' tab is active, showing a list box with 'DENTAL' selected. To the right of the list box are fields for 'Visit Date' (Jun 07, 2017) and 'Visit Time' (11:00), and an unchecked 'Historical Visit' checkbox.

Figure 21: Provider and Location for Current Activities Screen

The **Encounter Provider** field should default to the correct end-user that is signed into VistA. Select the **Encounter Location** if the **Default Location** parameter is NOT set in advance. The **Default Location** parameter is explained in the Treatment System section in the Using the DRM Plus Drop-Down Menus chapter of this manual.

**Visit Date/Time** defaults to the present date/time for a new visit in the **New Visit** tab. The date and time may be changed if desired.

**Figure 22: New Visit Tab**

To record a new visit other than the present date/time:

1. DRM Plus defaults to the present provider; however, a different provider may be selected using the **Encounter Provider** list.
2. Select the clinic location from the scroll menu if the **Default Location** is NOT set.
3. Use the drop-down arrow to toggle through the calendar screen and select a date.
4. Use the up and down arrows to adjust the **Visit Time**.
5. Check the **Historical Visit** check box, if applicable.
6. Click the **OK** button to create the new PCE visit.

**Note:** Future date appointments may NOT be created in DRM Plus.

**Note:** Creating a new PCE visit in DRM Plus does NOT update VistA **Appointment Manager**.

## View Scheduled Appointments and Historical Visits

The **My Clinic Visits** tab lists the patient visits for the selected clinic. This tab only displays if a default **Dental Location** parameter is selected. When no default **Dental Location** parameter is selected, the **Dental Visits** tab displays.

To record the scheduled appointment for the patient:

1. DRM Plus defaults to the present provider; however, another provider may be selected from the **Encounter Provider** list.
2. If there is only one scheduled visit, it is automatically defaulted.
3. Select the correct **scheduled visit** in the bottom window, if it is NOT defaulted.
4. Click the **OK** button and the provider/location is set for the scheduled visit.



**Provider and Location for Current Activities**

Encounter Provider

DRMPROVIDER,DENTIST STAFF DENTIST,PROSTHODONTICS

BLUYLUI,CXLY  
 CPRSPHYSICIAN,ONE STAFF DENTIST,GENERAL PRA  
 CXYTHYSDYX,LAKHUS  
 DRMPROVIDER,ADMINDENTIST CHIEF, DENTAL SERVICE,ORAL  
 DRMPROVIDER,ASSISTANT DENTAL ASSISTANT  
 DRMPROVIDER,DENTIST THREE STAFF DENTIST,GENERAL PRA  
 DRMPROVIDER,DENTIST TWO STAFF DENTIST,ENDODONTICS  
 DRMPROVIDER,DENTIST STAFF DENTIST,PROSTHODON

OK  
Cancel

Encounter Location

Jun 07, 17@11:00 DENTAL

My Clinic Visits Dental Visits All Visits Admissions New Visit

V	May 31, 2017@10:00	DENTAL
V	May 31, 2017@09:00	DENTAL
V	May 31, 2017@08:00	DENTAL
V	May 26, 2017@09:30	DENTAL
V	May 26, 2017@09:00	DENTAL
V	Apr 28, 2017@15:00	DENTAL
V	Mar 22, 2017@14:00	DENTAL
V	Mar 17, 2017@14:00	DENTAL

**Figure 23: My Clinic Visits Tab**

The **Dental Visits** tab lists all the dental clinic visits.

Encounter Location

< Select a visit from the tabs below...>

My Clinic Visits Dental Visits All Visits Admissions New Visit

V	May 31, 2017@10:00	DENTAL
V	May 31, 2017@09:00	DENTAL
V	May 31, 2017@08:00	DENTAL
V	May 26, 2017@09:30	DENTAL
V	May 26, 2017@09:00	DENTAL
V	Apr 28, 2017@15:00	DENTAL
V	Mar 22, 2017@14:00	DENTAL
V	Mar 17, 2017@14:00	DENTAL

**Figure 24: Dental Visits Tab**

The **All Visits** tab lists all dental visits and admission(s) if the selected patient is an inpatient.

Encounter Location		
< Select a visit from the tabs below...>		
My Clinic Visits   Dental Visits   <b>All Visits</b>   Admissions   New Visit		
A	OCT 18,2005@13:46	7B
V	May 31, 2017@10:00	DENTAL
V	May 31, 2017@09:00	DENTAL
V	May 31, 2017@08:00	DENTAL
V	May 26, 2017@09:30	DENTAL
V	May 26, 2017@09:00	DENTAL
V	Apr 28, 2017@15:00	DENTAL
V	Mar 22, 2017@14:00	DENTAL
V	Mar 17, 2017@14:00	DENTAL

**Figure 25: All Visits Tab**

The **Admissions** tab lists the admissions for the selected patient.

Encounter Location		
< Select a visit from the tabs below...>		
My Clinic Visits   Dental Visits   All Visits   <b>Admissions</b>   New Visit		
A	OCT 18,2005@13:46	7B

**Figure 26: Admissions Tab**

## Treatment & Exam

The **Treatment & Exam** menu has nine submenus: **Show Configuration**, **Add/Edit Personal QuickList**, **Add Medical Codes to ADA Table**, **Edit Code Information in ADA Table**, **Edit Procedure Costs**, **Filter View**, **Clean Slate**, **Undo Clean Slate** and **All Planned Care to Clipboard**.

**Note:** The **Add Medical Codes to ADA Table**, **Edit Code Information in ADA Table**, **Edit Procedure Costs**, **Clean Slate** and **Undo Clean Slate** are DRM Plus administrative submenus. Please contact a DRM Plus Administrator for more information about these submenus.




Figure 27: Treatment & Exam Menu

## Show Configuration

Select the **Show Configuration** submenu to display the **Charting Configuration** screen.

**Charting Configuration**

Speed Codes   Suggestion Links   Statistics   H & N  
Tx & Exam   Periodontal   Report   Voice

Initial Tx Mode: Treatment Plan

☐ Display sequencing upon entry  
☒ Show warning box when adding duplicate transaction on tooth

Display Defaults

When I am viewing: Diagnostic Findings

I would like the following to display:

Graphical Display	Transaction List
<input checked="" type="checkbox"/> Dx Findings	<input checked="" type="checkbox"/> Dx Findings
<input type="checkbox"/> Tx Plan	<input type="checkbox"/> Tx Plan
<input type="checkbox"/> Completed	<input type="checkbox"/> Completed

Reset   Ok

**Figure 28: Charting Configuration Screen**


Use the various tabs to configure the chart. The tabs are: **TX & Exam, Periodontal, Report, Voice, H&N, Statistics, Suggestion Links** and **Speed Codes**. These parameters on each tab are a user specific function; changing it does NOT impact other users. When finished, click the **OK** button.

## ***Tx & Exam***

Use the **Tx & Exam** tab to change the default view screen that displays when DRM Plus is first opened. The original default view screen is the **Treatment Plan** view.

The **Sequencing** screen displayed upon entry is NOT selected as a default parameter; however, showing a warning box when adding duplicate transactions on a tooth for each view chart is a default parameter. Use the check boxes to change these user specific parameter functions.

Use the functions on this tab to fine tune the **Display Defaults**; choose **Graphical Displays** or **Transaction Lists** to display check boxes based on the screen being viewed.



The image shows a 'Charting Configuration' dialog box with the 'Tx & Exam' tab selected. The dialog has a title bar with a close button. Below the title bar are several tabs: 'Speed Codes', 'Suggestion Links', 'Statistics', 'H & N', 'Tx & Exam' (selected), 'Periodontal', 'Report', and 'Voice'. The main area of the dialog contains the following settings:

- Initial Tx Mode:** A dropdown menu set to 'Treatment Plan'.
- ☐ Display sequencing upon entry
- ☒ Show warning box when adding duplicate transaction on tooth
- Display Defaults:** A section containing:
  - When I am viewing:** A dropdown menu set to 'Diagnostic Findings'.
  - I would like the following to display:** A section with two columns:


Graphical Display	Transaction List
<input checked="" type="checkbox"/> Dx Findings	<input checked="" type="checkbox"/> Dx Findings
<input type="checkbox"/> Tx Plan	<input type="checkbox"/> Tx Plan
<input type="checkbox"/> Completed	<input type="checkbox"/> Completed

At the bottom of the dialog are 'Reset' and 'Ok' buttons.

**Figure 29: Tx & Exam Tab**

## *Periodontal*

Choose the **Periodontal** tab to set pocket depth warning and choose the colors that display as pocket warnings and normal pockets on the **Periodontal Chart** screen. Other options on this tab include **Show MGJ Trace** and **Exam Sequence**.



The image shows a software window titled "Charting Configuration" with a close button (X) in the top right corner. The window has a tabbed interface with the following tabs: "Speed Codes", "Suggestion Links", "Statistics", "H & N", "Ix & Exam", "Periodontal" (which is the active tab), "Report", and "Voice". The "Periodontal" tab contains the following settings:


- Pocket Depth Warning:** A numeric input field containing the value "4".
- Normal Pocket Color:** A color selection button showing a black square.
- Warning Pocket Color:** A color selection button showing a red square.
- Enable ProbeOne:** An unchecked checkbox next to a dropdown menu currently showing "cboProbe1".
- Show MGJ Trace:** A checked checkbox.
- Exam Sequence:** A text input field containing "F1, L16, F17, L32" and an "Edit" button to its right.

At the bottom of the window are two buttons: "Reset" and "Ok".

**Figure 30: Periodontal Tab**

To change the exam sequence:

1. Click the **Edit** button.
2. The **Edit Perio Sequence** screen displays.




**Figure 31: Edit Perio Sequence Screen**

3. Click each **Section** in the order in which the perio exam sequence should be performed.
4. Click the **OK** button to save the new exam sequence.

To go back to the original settings, which appeared when this screen was first displayed, click the **Reset** button. Once the exam sequence has been changed and the user has clicked the **OK** button on the **Periodontal** tab, this becomes the permanent default exam sequence.

## ***Report***

Use the functions on the **Report** tab to select certain pieces of information, which appears on individual reports when using the **Print** submenu under the **File** menu. The **Chart** selection prints the graphic chart, displayed on the last view screen of the **Chart/Treatment** tab, prior to the chosen **Print** submenu. The **Transactions** selection prints the transaction table, displayed on the last view screen of the **Chart/Treatment** tab, prior to the chosen **Print** submenu. **Patient Notes** and **Tooth Notes** selections print the entries entered using the **Notes** icon.



**Figure 32: Report Tab**


## ***Voice***

**Voice** is NOT enabled in DRM Plus.



## Speed Codes


Use the **Speed Codes** tab to set/create individual icons in DRM Plus for frequently used procedure codes entered using the **Treatment Plan** or **Completed Care** viewing screens.



**Figure 33: Speed Codes Tab**

To add a speed code:

1. Click the **Add** button. The **Edit Speed Code** screen displays.
2. Add a new **Name**, which cannot exceed 10 characters.
3. Entering a description is optional. Two symbols; semicolon (;) and up-carrot (^) may NOT be added or used in the description text.
4. Use the search function, **ADA Codes**, to look up a procedure code(s) and add it to the new speed code.
5. Entering an icon is optional.
6. Click the **OK** button to begin finalization.



**Edit Speed Code**

Name:  Description:

ADA Codes:


Icon:

**Figure 34: Edit Speed Code Screen**

To edit or delete the speed code, highlight the desired name in the **Speed Codes** tab and click the **Edit** or **Delete** button. Provide appropriate entry in the subsequent screens; otherwise, click the **OK** button to complete this part of the process.

To complete the speed code process:

1. Move to the **Completed Care** or **Treatment Plan** view of the **Treatment & Exam** screen.
2. Click one of the undesignated icon squares. The **Configure Button** screen displays.



**Configure Button**

Speed Code:

Perio Mode: ☐


**Figure 35: Configure Button Screen**

3. Click the **drop-down arrow**, highlight and click the desired **Speed Code** name.
4. Click the **OK** button and the speed code is linked to that icon.

The **Perio Mode** check box on the **Configure Button** screen designates the viewing preference when the **Perio Buttons** icon is clicked. The **Perio Buttons** icon is used as a toggle for displaying another 19 available icon buttons. Clicking the **Perio Buttons** icon displays any 19 **Speed Code** icons that have been designated in the **Perio Mode** (check box clicked) while hiding any non-perio mode **Speed Code** icons from the display. Clicking the **Perio Buttons** icon again reverses the display. This option allows for a total of 38 **Speed Code** icons to be created. The 19 non-perio mode speed codes are the default **Speed Code** icons when DRM Plus is initially opened. Please see the Perio Buttons Icon section, in the Chart/Treatment – Treatment & Exam chapter of this manual, for more information.

## *Suggestion Links*


Use the **Suggestion Links** tab to enter code suggestions, when entering one procedure code which is linked to another procedure code(s), without having to use an icon to find the other code. A screen displays asking if other linked codes should be added providing an opportunity to decline the entry of suggested linked codes.



**Figure 36: Suggestion Links Tab**

To add a suggestion link:

1. Click the **Add** button.
2. A screen displays featuring a list of all DRM Plus procedure codes. Click the desired **primary procedure code** that other procedure codes are linked to, and then click the **OK** button.



**Figure 37: Suggestion Links Code List**

3. A screen requesting the linked codes to the primary procedure code displays.



**Figure 38: Linked Codes Screen**

4. Click the **Add** button to add the first linked code. The list of all DRM Plus procedure codes displays again.
5. Choose the second code to be linked with the primary procedure code and click the **OK** button.
6. Add as many linked codes to the primary procedure code as desired. To finish and return to the tab, click the **OK** button.

**Note:** As many codes as necessary can be linked. Simply continue clicking the **Add** button on the **Linked codes** screen and choosing more codes from the list.


To edit the suggestion link:

1. Select a suggestion link to be edited and click the **Edit** button. The **Linked codes** screen displays.
2. Click the **Add** button for another procedure code, and the list of procedure codes displays. Click the **OK** button.
3. To remove a linked code entry, click the **Delete** button and then the **OK** button.

To delete the suggestion link, select the suggestion link and click the **Delete** button.

## Statistics

Choose the **Statistics** tab to set the warning level for pocket depth, free gingival margin, mucogingival junction, furcation, and mobility found in the **Stats** screen on the **Periodontal Chart**.




The image shows a 'Charting Configuration' dialog box with the 'Statistics' tab selected. The tab contains a table titled 'Statistics Warning Levels:' with five rows. Each row has a 'Statistic' column and a 'Warning Level' column. The warning levels are: Pocket Depth: 4, Free Gingival Margin: 4, Mucogingival Junction: 3, Furcation: 2, and Mobility: 1.5. At the bottom of the dialog are 'Reset' and 'Ok' buttons.

Statistic	Warning Level
Pocket Depth:	4
Free Gingival Margin:	4
Mucogingival Junction:	3
Furcation:	2
Mobility:	1.5

Figure 39: Statistics Tab

To change the warning level:

1. Double-click the box containing the **Warning Level** to be changed.
2. A screen displays. Enter the new warning level in the text box and click the **OK** button.



The image shows a 'DRM' dialog box with the title 'Please enter new value'. It contains a text input field at the bottom and two buttons, 'OK' and 'Cancel', on the right side.

Figure 40: Adjust the Warning Level

3. The **Warning Level** is changed on the tab.

## H&N

The **H&N** tab can be entered by any user but only saved by DRM Plus Administrators. Please contact a DRM Plus Administrator for more information.

### Add/Edit Personal QuickList

Select the **Add/Edit Personal QuickList** submenu to manage a **Quick List** of codes for personal use. For additional convenience, enter frequently used procedure codes that have multi-add functionality associated with the code, into the **Quick List**. The **Manage Personal QuickList** screen displays.




Figure 41: Manage Personal QuickList Screen

To add to the **Quick List**:

1. Type the search criteria into the **Find** text box. Search by words or numbers.
2. A matching list displays on the left side of the screen. Click one of them to select it.
3. Click the **right arrow** button to move the selected code to the **Quick List**.
4. Click the **OK** button to end and close the screen or repeat to add another code to the **Quick List**.

To remove from the **Quick List**:

1. Select an entry from the **Quick List** on the right side of the **Manage Personal QuickList** screen.
2. Click the **left arrow** button to remove it from the list. A screen displays confirming that the entry is to be deleted. Click the **Yes** button to continue.

**Note:** Codes entered in a **Quick List** are accessed through the **Quick Code** icon.

### Add Medical Codes to ADA Table

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Edit Code Information in the ADA Table

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Edit Procedure Costs

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Filter View

Use the **Filter View** submenu to choose which encounters display on the **Chart/Treatment** tab of DRM Plus.




Figure 42: Filter View Submenus

### *Current Episode of Care*


Select this filter to show only those treatments that have been completed for all visits during the current disposition or patient status.

### *All Episodes of Care*

Select this filter to show all treatments completed for all visits. This is the default setting.

### *Select Episode of Care*

Select this filter to see all the treatments completed for all visits during a previous specific disposition or patient status. When this submenu is selected, a screen listing all previous dispositions or patient statuses associated with a given patient displays.



The dialog box titled "Select Episode-of-Care" contains a table with three columns: "Date", "Time", and "Disposition". The table has 10 rows. Below the table are "OK" and "Cancel" buttons.


Date	Time	Disposition

**Figure 43: Select Episode of Care Screen**

To select a previous disposition, click the desired one from the list and click the **OK** button.

### ***Date Range***

Select this filter to show treatments that have been completed within a specified date range. When this filter is selected, the following screen displays. Use this screen to select a date range.



The dialog box titled "Select Date Range" contains two date selection fields. The "From:" field shows "6/ 7/2016" and the "To:" field shows "6/ 7/2017". Both fields have drop-down arrows. Below the fields are "OK" and "Cancel" buttons.

**Figure 44: Select Date Range Screen**

To filter by date range:

1. Use the drop-down menu to select the needed dates. Click the **OK** button.
2. The treatments completed in the entered date range display on the screen. If no entries were made during the selected date range, DRM Plus displays as a clean slate.

### **Clean Slate**

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

### **Undo Clean Slate**

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.



## **All Planned Care to Clipboard**

**All Planned Care to Clipboard** is a submenu located as the last header on the **Treatment & Exam** menu. This submenu when selected will allow the user to copy all planned data; filed and unfilled planned treatment from the **Seq Plan** screen with extra planned details that occur on the cover page. This copy may be pasted on any word document, text document or any application window if allowed.

Select the **All Planned Care to Clipboard** submenu and place your cursor in the word document, text document or any application window where you would like the planned data pasted.

## Tools

The **Tools** menu has 9 submenus: **Windows Calculator**, **Windows Explorer**, **Windows Notepad**, **User Inquiry**, **User Options**, **Administrative Toolbox**, **Panel Add/Edit**, **Provider Add/Edit** and **Vitals**.

The **ADA Website** submenu is an ancillary application that the DRM Plus Administrator may customize for all users. The DRM Plus Administrator may customize up to 10 ancillary applications.

**Note:** **Administrative Toolbox**, **Panel Add/Edit** and **Provider Add/Edit** are DRM Plus administrative functions. Therefore, these submenus are NOT selectable.




Figure 45: Tools Menu

### **Windows Calculator**

Select this submenu to open **Windows Calculator**.

### **Windows Explorer**


Select this submenu to open **Windows Explorer**.

### **Windows Notepad**

Select this submenu to open **Windows Notepad**.

## User Inquiry

Select this submenu to view and change the VistA fields or to view the VistA fields of other users. The **VistA User Inquiry** screen displays.



**Figure 46: VistA User Inquiry Screen**

1. Type the **User Name** into the input text box and press the <Enter> key.
2. The results display on the left side of the screen.
3. Select a **User Name** to view. The user's information displays on the right side of the screen as shown in the next figure.

**VistA User Inquiry...**

View Users or View and change Your Vista Fields  
Click on 'User Tool Box' to change the fields for your name in VistA.

User Name (type and press enter)	User Information
drm	Initials: DEN No NickName on File.
DRM,KDF	No Title on File.
DRM,YHP	Service: DENTAL
DRMPROVIDER,ADMINDENTIST	No Mail Code on File.
DRMPROVIDER,ASSISTANT	No Office Phone on File.
<b>DRMPROVIDER,DENTIST</b>	No Voice Pager on File.
DRMPROVIDER,DENTIST THREE	No Digital Pager on File.
DRMPROVIDER,DENTIST TWO	No Fax on File.
DRMPROVIDER,HYGIENIST	No Location on File.
DRMPROVIDER,RESIDENT	Menu Option: NO MENUS GENERATED
DRMPROVIDER,RESIDENT TWO	Type Ahead Allowed: ALLOWED
	Preferred Editor: SCREEN EDITOR - VA FILEMAN
	Signature Block: Dentist DRMPROVIDER
	No Signature Title on File.

User Tool Box (for DRMPROVIDER,DENTIST)    **OK**    Cancel

**Figure 47: VistA User Entry Screen with User Information Displayed**

Select the **User Tool Box** button to change the personal fields in VistA. Click the **User Tool Box** button at the bottom of the screen and the **User's TBox** screen displays.

**User's TBox...**

User Profile Fields:

- ☒ Initials
- ☐ Nickname
- ☐ Service
- ☐ Mail Code
- ☐ Office Phone
- ☐ Home Phone
- ☐ Voice Pager
- ☐ Digital Pager
- ☐ Fax
- ☐ Location
- ☐ Automenu
- ☐ Type Ahead
- ☐ Preferred Editor
- ☐ Signature Block

Update Initials field.

DEN


**Finished**

**Figure 48: User's TBox Screen**


4. Select the desired **User Profile Fields** by clicking the corresponding radio button.
5. Edit the new text in the text box.
6. Click the **Update Field** button.
7. Click the **Finished** button. The **VistA User Inquiry** screen displays again.

## User Options

Adjust various user settings in this submenu. The screen contains five tabs: **General**, **Printing**, **Progress Note**, **Treatment System** and **Exam Settings**.




**Figure 49: User Settings Screen**

The **Broker Call History** icon  opens the broker calls screen. Please see the Last Broker Call section, in the Using the DRM Plus Drop-Down Menus chapter of this manual, for more information.

## General

The default **General** tab allows the provider to change **Date Range defaults**, **Other Settings** and **File Location** folders.




**Figure 50: General Tab**

The **Delete User Settings** button located on the lower left corner of the screen displays in all the tabs. This button allows the user to delete any new changes in this session before the parameter is saved.

**Note:** The **Delete User Settings** function only applies to the user that is currently logged in. Other users are NOT affected if one chooses to delete user settings.

To change the default date ranges:

1. Click the **Date Range Defaults** button.
2. The **Date Range** screen displays.



**Date Range** [X]

**Visits/Admissions**

Days Back: 365 [up/down] Maximum records: 100 [up/down]

**Progress Notes**

Days Back: 180 [up/down] Maximum records: 100 [up/down]

**Lab Results**

Days Back: 180 [up/down]

**Tooth History**

Days Back: 999 [up/down] Maximum records: 999 [up/down]


OK Cancel

**Figure 51: Date Range Screen**

3. Use the up and down arrows to set the desired date range.
4. Click the **OK** button to return to the **User Setting** screen.

To change other parameter settings:

1. Click the **Other Parameters** button.
2. The **Other Parameters** screen displays.



**Other Parameters** [X]

**Tabs**

Initial tab when DRM starts: Chart/Treatment [down arrow]

Initial chart display in Chart/Treatment: Treatment & Exam [down arrow]

**Note Boilerplate**

☒ Prompt for Note Boilerplate insert


OK Cancel

**Figure 52: Other Parameters Screen**

3. Use the **Tabs** drop-down menu to set the initial tab and chart display in DRM Plus.
4. Use the **Note Boilerplate** check box to indicate whether the program should prompt for the boilerplate insert associated with the VistA TIU progress note title selection.

To change the file folder location:

1. Click the **Set default folder** button.
2. The **Select Default Folder** screen displays.



**Figure 53: Select Default Folder Screen**


3. Navigate to and click the **desired folder**.
4. Click the **OK** button to select it.

**Note:** This option allows the importing of information stored as a .txt file into the TIU progress note.



## *Printing*

Use the **Printing** tab to set print margins, orientations, etc.



**Figure 54: Printing Tab**

To change the page configuration:

1. Click the **Page Setup** button.
2. The **Page Setup** screen displays.





Figure 55: Page Setup Screen

3. Use the up and down arrows to adjust the margins.
4. Use the **Orientation** radio buttons to change the orientation of the printed document.
5. Use the **Page Number** check box to indicate whether page numbers are to be included.
6. Use the **Ellipsis [...]** buttons to choose fonts for the **Page Text**, **Header Text** and **Page Number**.
7. Click the **OK** button to return to the **User Settings** screen.

## *Progress Note*

Use the functions in the **Progress Note** tab to configure progress note data objects, configure note data sequence and configure code boilerplates.



**Figure 56: Progress Note Tab**

To configure progress note data objects:

1. Click the **Progress Note Data** button.
2. The **Progress Note Data Objects** screen displays.

**Progress Note Data Objects**

Include the following data objects

<input checked="" type="checkbox"/> Diagnostic Findings	<input checked="" type="checkbox"/> Modified	<input type="checkbox"/> Deleted
<input checked="" type="checkbox"/> Adverse Events	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> Planned Items	<input checked="" type="checkbox"/> Modified	<input type="checkbox"/> Deleted
<input checked="" type="checkbox"/> Periodontal Exam		
<input checked="" type="checkbox"/> PSR Exam	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> Head/Neck Findings	<input checked="" type="checkbox"/> Modified	<input checked="" type="checkbox"/> Deleted
<input checked="" type="checkbox"/> Tooth Notes	<input checked="" type="checkbox"/> Modified	<input checked="" type="checkbox"/> Deleted
<input checked="" type="checkbox"/> Completed Items	<input type="checkbox"/> Modified	<input checked="" type="checkbox"/> Deleted
<input checked="" type="checkbox"/> Device Tracking	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> Social History	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> OHA Findings	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> Occlusal Findings	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> Parafunctional Habits	<input checked="" type="checkbox"/> Modified	
<input checked="" type="checkbox"/> TMJ Findings	<input checked="" type="checkbox"/> Modified	
<input type="checkbox"/> Code Boilerplate	<p>These items are always included. The user may opt to not include the Note BP when selecting a note title.</p> <input checked="" type="checkbox"/> Note Boilerplate <input checked="" type="checkbox"/> Text Marker	
<input checked="" type="checkbox"/> Next Appointment		
<input type="checkbox"/> Dental Alerts		

OK Cancel


**Figure 57: Progress Note Data Objects Screen**

3. Use the various check boxes to include or exclude desired progress note data objects.
4. Click the **OK** button to return to the **User Settings** screen.

**Note:** The **Code Boilerplate** check box activates the automatic importing into the TIU progress note of any code boilerplate created in DRM Plus.

To configure the note data sequence:

1. Click the **Set Note Sequence** button.
2. The **Note Object Sequence** screen displays.




**Figure 58: Note Object Sequence Screen**

3. Select the note object to be moved in the list.
4. Use the up and down arrows on the right side of the screen to change the sequence of the note object on the list.
5. Click the **OK** button to return to the **User Settings/Progress Note** tab screen.

To configure the code boilerplate:

1. Click the **Configure Code Boilerplate** button.
2. The **Code Boilerplate** screen displays.



**Figure 59: Code Boilerplate Screen**

To add a new code boilerplate:

1. Click the **Add New** button.
2. The **New Boilerplate** screen displays.




Figure 60: New Boilerplate Screen

3. Enter the name in the text box and click the **OK** button.
4. The **Code Boilerplate Text** screen displays.





Figure 61: Code Boilerplate Text Screen

5. Click the **Add Code** button to add a code to the boilerplate.
6. The **Find CPT Code or CPT Description Text** screen displays.




**Figure 62: Find CPT Code Screen**

7. Type in the CPT code. A partial number is acceptable. Press the <Enter> key.
8. The search results display in the screen. Select one and click the **OK** button.
9. The selected code displays on the **Code Boilerplate Text** screen. The provider may add more than one CPT code to this code boilerplate.
10. To delete that code, click the **Delete Code** button in the **Code Boilerplate Text** screen.
11. Type the desired associated text into the right side of the **Code Boilerplate Text** screen.
12. Click **OK**. A confirmation screen displays. Click **OK** to return to the **Code Boilerplate** screen.

To edit a code boilerplate:

1. Select the code boilerplate to be edited from the **Code Boilerplate** screen.
2. Click the **Edit** button.
3. The **Code Boilerplate Text** screen displays.



**Figure 63: Code Boilerplate Text Screen**


4. From here, type in the right side of the screen to add or delete text from the boilerplate. Use the **Add Code** and **Delete Code** buttons to add or delete codes from the boilerplate.
5. Click the **OK** button. An information screen displays. Click the **OK** button to return to the **Code Boilerplate** screen.

To delete a code boilerplate:

1. Select a code boilerplate from the list on the **Code Boilerplate** screen.
2. Click the **Delete** button.
3. A confirmation screen displays. Click the **Yes** button to delete the boilerplate.
4. An information screen displays. Click the **OK** button to return to the **Code Boilerplate** screen.

## *Treatment System*

The **Treatment System** tab allows access to additional options.




**Figure 64: Treatment System Tab**

Use the **check boxes** to choose whether to prompt for a diagnostic code when adding a planned item, or select the default tree view to display DRM Plus note objects on the **Progress Note** screen.



To choose a default location:

1. Click the **Ellipsis (...)** button next to the **Default Location** text box.
2. The **Select Location** screen displays.



**Figure 65: Select Location Screen**

3. Type the location into the text box and press the **<Enter>** key.
4. Search results display on the **Select Location** screen.
5. Choose the desired location and click the **OK** button.
6. Select the **OK** button from the informational screen and the location is saved.
7. Use the **Clear** button if the location should be removed, then the **OK** button on the informational screen.
8. The user may also change the default location by using the **Select Location** screen.
9. Always select the **OK** button from the informational screen to save any changes.
10. Select the **Done** button at the bottom of the screen to close the **User Settings** screen.

The **Delete User Settings** button located on the lower left corner of the screen displays in all the tabs. This button allows the user to delete any new changes in this session before the parameter is saved.

**Note:** The **Delete User Settings** function only applies to the user that is currently logged in. Other users are NOT affected if one chooses to delete user settings.

## Exam Settings

The **Exam Settings** tab provides the user with several options. These include: **Canned Statements**, **Next/Back Button** and **Requirements**.




Figure 66: Exam Settings Tab

The **Canned Statements** parameter allows adding of additional pre-defined statements by the end-user to four elements. All local providers are end-users when utilizing this function from the **User Options** submenu, whether they are DRM Plus administrators. Any changes made from the **User Settings** screen affect only the individual end user.


Pre-defined statements are broken into five categories: **Radiographic**, **Assessment Summary** and **Treatment Plan** (located in the same element), **Patient Education** and **Disposition**. There is a maximum of twelve pre-defined statements allowed per category.

The local DRM Plus Administrator has priority when entering these statements system-wide, utilizing the administrator settings parameter (NOT displayed here).

When any of these element categories are maxed out with pre-defined statements, the DRM Plus Administrator may add another. This can be done by utilizing the administrator settings parameter. This hides the last pre-defined statement entered by any end-user, and only affects those end-users with twelve entered and displayed in the given category.

To add a pre-defined statement (admin or non-admin) from the **User Settings** screen:

1. Select one of the five pre-defined statement buttons, such as **Assessment Summary**.
2. Type or copy/paste a pre-defined statement in the lower text box.
3. Click the green **Add (+)** button.
4. Click the **OK** button to confirm the new pre-defined statement addition.




**Figure 67: Assessment Summary Screen**

The end-user may highlight any of the pre-defined statements that were entered from their **User Settings** screen and either delete that statement or move the statement's position in the list. This deletion or rearranging of the order only affects the end-user's list of pre-defined statements and **NOT** any entered by the DRM Plus Administrator or any national pre-defined comment that was kept by the DRM Plus Administrator; these are listed at the top.

The **Next/Back Button** parameter setting allows end-users, when selecting the **Next** or **Back** buttons located on any **Exam** tab element screen, to go directly to the next proceeding or previous required element screen for that exam code and skip all optional element screens.

**Note:** There is no **Back** button on the **Presentation/Chief Complaint** element screen and there is no **Next** button on the last **Disposition** element screen.



**Figure 68: Next/Back Button Parameters**

Both options are unchecked by default. When unchecked, the **Next** button skips any element that is optional or has been completed from new data entered during this session and opens the next required element. When checked, the **Next** button opens the very next element regardless of whether it is optional or completed during this session.

The **Back** button when unchecked skips any element that is optional but opens all previous required elements that are completed or NOT. When checked, the **Back** button opens the previous element regardless of whether it is optional or required.

The user is required to complete any optional or required element when selecting the **Next** button when trying to move forward. Selecting the **Back** button does NOT require the element to be completed to open a previous element.

**Note:** When this parameter has been formatted in the **User Settings** screen, these selections only affect the end-user's profile and follows that end-user to any computer when loading DRM Plus with their VistA access/verify codes.

The first requirements parameter **Configure Requirements Display**, allows the end-user to keep the requirements display open when selecting any element from the **Exam** tab or the definitions from the **OHA** and **Occlusal** screens. The second requirements parameter, **Configure Radiographs Requirement**, allows the end-user to require a **Radiograph Finding** entry with any exam/consult code entered as completed care and requires data entered using the **Exam** tab.





Figure 69: Requirements Parameters

The **Configure Requirements Display** parameter is checked by default and displays the element's requirements whenever an element screen is open. When unchecked, the end-user must select the **Done** button and then close/reopen DRM Plus to activate. This parameter change requires the end-user to open the **Element Requirements Panel** manually.

The **Requirements** icon  located in the upper right corner of the element screen when selected displays the **Element Requirements Panel**.




Figure 70: Presentation/Chief Complaint Element Screen

The **Configure Radiographs Requirement** is unchecked by default and only requires radiographs for the D0150 and D0180 exams. When checked, the end-user must select the **Done** button and then close/reopen DRM Plus to activate. This parameter change requires the end-user to enter data from the **Radiographic Findings** element with any exam/consult code entered as completed care and requires data entered using the **Exam** tab.

**Note:** The ADA exam codes D0145, D0171, D0190 and D0191 are NOT included with the DRM Plus **Exam** tab functionality when entered in DRM Plus by a provider.

## **Administrative Toolbox**

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

### **Panel Add/Edit**

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

### **Provider Add/Edit**

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.


## **Ancillary Tool Functions – ADA Website**

This American Dental Association website is only available if the DRM Plus Administrator formats this in the administrator's **Ancillary** applications and parameters. Some users may NOT have permission to access the Internet or must enter/re-enter a username/passcode. Please see the DRM Plus Administrator for more details.

## Reports

The **Reports** menu has the seven following submenus: **Reports** (DRM canned reports), **Service Reports** (old DAS reports), **Load Data Warehouse Reports**, **Adverse Events**, **Device Tracking**, **Extract History File** (for small date range extract reports), and **Queued Extract History File** (for large date range extract reports).

**Note:** The following reports: **Adverse Events**, **Device Tracking**, **Extract History File** and **Queued Extract History File** require DRM Plus administrative functionality.



**Figure 71: Reports Menu**

## **Reports**

When the **Reports** submenu is selected, the **Report Selection** screen displays. The **Report Selection** screen has three tabs: **General**, **Patient** and **Planning**.

## General Tab

The screenshot shows the 'Report Selection' dialog box with the 'General' tab selected. The dialog has three tabs: 'General', 'Patient', and 'Planning'. Under the 'General' tab, there are two columns of radio button options for report types: 'Provider Summary' (selected), 'Clinic Summary', 'Visits by Provider', 'Visits by Clinic', 'Non Clinical Time by Provider', 'Fee Basis/Detailed Fee Basis', 'Encounter/Visits by Patient Type', and 'Recare Report'. Below these are three dropdown menus for 'Fiscal Year' (n/a), 'Start Date' (6/ 7/2017), and 'End Date' (6/ 7/2017). There are checkboxes for 'All Stations' (checked), 'All Providers?' (unchecked), 'Use Provider Name on Reports' (checked), and 'Include Distributed Provider Totals' (unchecked). A table titled 'Select Provider(s)' contains one row with the text 'DRMPROVIDER,DENTI...' in the 'Name' column, '03070003' in the 'ID' column, and 'STAFF DENTIST,PROSTHODO' in the 'Type/Specialty' column. To the right of the table are two groups of radio button options: 'Patient Status' with 'Active', 'Inactive', 'Maintenance', and 'All Statuses' (selected), and 'Active/Maint'. Below these are 'Transaction Status' options: 'Complete' (checked), 'Planned', and 'Deleted'. At the bottom left, there is a 'Report Category Type' section with radio buttons for '13 ADA Categories', '131 VA-DSS Prod', and 'ADA/CPT Codes' (selected). At the bottom right, there is a 'Search for records using' section with radio buttons for 'Visit Date' (selected) and 'Create Date'. 'OK' and 'Cancel' buttons are at the bottom right of the dialog.

Name	ID	Type/Specialty
DRMPROVIDER,DENTI...	03070003	STAFF DENTIST,PROSTHODO

Figure 72: General Tab

To create a report:

1. Choose the desired report type.
2. Select the **Fiscal Year** or the **Start Date** and **End Date**.
3. Use the check box to indicate whether the provider name and the distributed provider totals should be included in the report.
4. Choose a **Patient Status**.
5. Indicate what the **Transaction Status** is.
6. Select the **Report Category Type**.
7. Choose the date type that is to be represented on the report.
8. Click the **OK** button to generate the report. The report screen displays.

**DRM Report**

Dental Service Treatment Report - Summary Report by Provider

From Jan 01, 2017 to Jun 07, 2017 Station: All Dental Provider:

ADA/CPT Code	IIC	IV	Total	Tot RVU
D0120 PERIODIC ORAL EVALUATION		3	3	90
D0210 INTRAORAL FULL IMAGE SERII		2	2	36
D0274 DENTAL BITEWING FOUR IMA		1	1	4
D0330 DENTAL PANORAMIC IMAGE		3	3	30
D2160 AMALGAM THREE SURFACES F		6	6	300
D3428 BONE GRAFT W/PRAD SURG-F		1	1	10
D4263 BONE REPLCE GRAFT FIRST S		1	1	10
D4264 BONE REPLCE GRAFT EACH AI	1		1	5
D6010 ENDOSTEAL IMPLANT BODY F		5	5	450
D7140 EXTRACTION ERUPTED TOOT		2	2	60
D7950 MANDIBLE GRAFT		3	3	180
D9215 LOCAL ANESTHESIA		1	1	0
D9930 TREATMENT OF COMPLICATI		2	2	60
<b>TOTAL</b>	<b>1</b>	<b>30</b>	<b>31</b>	<b>1235</b>

Total Visits = 13  
All data is complete ☒ Print Totals Only

**Figure 73: DRM Report Screen**

This screen has options to save an Excel file (**Save to XLS** button) or close. Some of the options may NOT be available with every report type.

Eight report types are accessible through this tab:

- **Provider Summary:** Summary counts of procedures by Station/Provider and Dental Classification.
- **Clinic Summary:** Summary counts of procedures by Station and Dental Classification.
- **Visits by Provider:** Detailed listing of procedures by Station/Provider.
- **Visits by Clinic:** Detailed listing of procedures by Station.
- **Non-Clinical Time by Provider:** Total days by provider for time applied to Education, Administration, Research and Fees.
- **Fee Basis/Detailed Fee Basis:** Total amount authorized and number of cases by Dental Classification (for local dental use only).
- **Encounters/Visits by Pat Type:** Summary counts of encounters/visits by patient type.
- **Recare Report:** List of patients with recare dates.

Click the corresponding radio buttons to select the desired report types. Use the check boxes to customize these reports.

**Note:** Selecting multiple reports from the **General** tab while the **Report Selection** screen is displayed always requires the selection of the report radio button first, followed by the selection of the **Fiscal Year**, even if the same fiscal year is desired for multiple reports.



## ***Patient Tab***

Use the **Patient** tab to run a report on any patient and view a list of visits. The patient is only used for report generation; changing patients in DRM Plus is NOT allowed.


The screenshot shows the 'Report Selection' dialog box with the 'Patient' tab selected. The 'Patient Visit List' radio button is chosen. The 'Patient' field contains 'DRMPATIENT.ONE (D0001)' and the 'Patient Selection' button is visible. Below these are fields for 'Fiscal Year' (n/a), 'Start Date' (6/ 8/2016), and 'End Date' (6/ 8/2017). Checkboxes for 'All Stations' (checked), 'All Providers?' (unchecked), and 'Use Provider Name on Reports' (checked) are present. A 'Select Provider(s)' section contains a table with one row selected.

Name	ID	Type/Specialty
DRMPROVIDER,DENTI...	03070003	STAFF DENTIST,PROSTHODO

At the bottom of the dialog are 'OK' and 'Cancel' buttons.

**Figure 74: Patient Tab**

1. Click the **Patient Selection** button to select a patient. The program automatically defaults to the patient whose record is currently opened in DRM Plus.
2. The **Patient Selection** screen displays.



**Figure 75: Patient Selection Screen**

3. Type a patient name into the text box and press the **<Enter>** key. Partial names are acceptable.
4. Select the desired patient from the results box and click the **OK** button.
5. The selected patient's name now displays on the **Patient** tab.
6. Choose the date and select other information to be included or excluded using the check boxes on the **Patient** tab.
7. Click the **OK** button.
8. The **DRM Report** screen displays. Save to an Excel file, print or close.

## Planning Tab

Use the options in the **Planning** tab to run planning reports, active patients by provider report or unfiled data by provider report.

The screenshot shows the 'Report Selection' dialog box with the 'Planning' tab selected. The dialog has three tabs: 'General', 'Patient', and 'Planning'. The 'Planning' tab contains the following options:

- Radio buttons for report types:
  - ☒ Provider Planning
  - ☐ Planned Items List
  - ☐ Planned Non-VA Care
  - ☐ Active Patients by Provider
  - ☐ Unfiled Data by Provider
- Check boxes for provider information:
  - ☒ Use Primary Provider
  - ☐ Use Secondary Provider
  - ☐ Use Entered By Provider
- Check boxes for station and provider selection:
  - ☒ All Stations
  - ☐ All Providers?
  - ☒ Use Provider Name on Reports
- A table to select a provider(s):

Name	ID	Type/Specialty
DRMPROVIDER,DENTI...	03070003	STAFF DENTIST,PROSTHODO
- Radio buttons for patient status:
  - ☒ Active
  - ☐ Inactive
  - ☐ Maintenance
  - ☐ All Statuses
  - ☐ Active/Maint

At the bottom of the dialog are 'OK' and 'Cancel' buttons.

Figure 76: Planning Tab

1. Select the type of report from the four radio buttons.
2. Use the check boxes to indicate provider information.
3. Choose a **Patient Status** except for the **Unfiled Data by Provider** report.
4. Click the **OK** button.
5. The selected reports screen will display. Print or save the results to Excel.

**Note:** The Primary/Secondary provider option is utilized for these reports: **Provider Planning**, **Planned Items List**, **Planned Non-VA Care** and **Active Patients by Provider**.

## Unfiled Data by Provider Report

The **Unfiled Data by Provider** report displays a list of patients who have unfiled data for providers. Unfiled data is data that resides in a temporary scratch pad-type area and is only visible by the provider the data is saved to.

This data is NOT part of the patient's chart record and should be filed to completion in a timely manner. Unfiled data becomes inactive after eight-calendar days; the saved data is viewable but no longer available to be filed.

Provider	Patient Name & Last 4 SSN	Date saved as "Unfiled"	Inactive
<input type="checkbox"/> DRMPROVIDER,DENTIST	IDOSEPATIENT,ONE M (I6001)	Jun 08, 2017	<a href="#">View data</a>

**Figure 77: Unfiled Data by Provider Screen**


After selecting an **Unfiled Data by Provider** report, the provider needs to select the **View data** button which allows the user to display the data that was saved as unfiled data on that patient.

The **TX Note Preview** screen opens and displays the save unfiled data. This displays the unfiled data saved by this provider or by some other provider who sent it to this provider on a specific patient.

The provider may print this unfiled data, especially if the data was made inactive either by the unfiled data now saved over the 8-day limit, or if a DRM Plus Administrator used the **Clean Slate** submenu option on this patient's chart. An example of inactive unfiled data would have a **Yes** listed in the **Inactive** column. Because DRM Plus is unable to reload this inactive data back into the patient's chart, the provider is required to re-enter the data manually with another encounter.

The provider may delete any unfiled data by selecting the check box under the **Provider** column and then selecting the **Delete Checked** button. The **Check Inactives/Uncheck Inactives** button allows the provider to select/unselect all the inactive unfiled data reports. The **Check All/Uncheck All** button allows the provider to select/unselect all the check boxes in the **Unfiled Data by Provider** report.

The following dialog is an example of unfiled data saved on a patient. The user may print the unfiled data by selecting the **Print** button or close the screen by clicking the **OK** button.



**Figure 78: TX Note Preview**

Non-administrative end-users can delete, view and print the active/inactive unfiled data for all their respective patients when accessing this report. The **Unfiled Data by Provider** report only allows a non-administrative provider to view their own saved unfiled data and NOT of other providers.

The following dialog is the screen that allows the provider to **Load**, **View** (non-load) or **Delete** any unfiled data when opening the DRM Plus chart for a patient. The **Delete** button allows the provider to delete unfiled data before it is loaded into this patient's chart. The provider is NOT able to view the unfiled data when they select the **Delete** button from this screen.

The **Load DRM Plus Data** screen informs the provider the name of who saved the unfiled data to them. The name listed at the beginning of the statement on the screen is the person who saved the unfiled data to the user opening the patient's chart.




Figure 79: Load DRM Plus Data Screen

There are two ways to view the unfiled data before the provider deletes this data. The first is to select the **Load** option and go to the **Unfiled Data by Provider** report located on the **Reports** menu → **Reports** submenu → **Planning** tab → **Unfiled Data by Provider** report.

The second way to view unfiled data when the provider cannot remember exactly what was saved as unfiled data is to select the **View** button. This option directs the user to the **Unfiled Data by Provider** report, but it does **NOT** load the unfiled data. If the user wants the unfiled data loaded and filed, then they must close the report and select the **Refresh Patient Chart** submenu under the **File** menu.

Selecting the **View** button displays a screen giving instructions to the end-user on how to load the data into the patient's chart. Selecting the **Load** button when the patient's chart reopens after the refresh allows the user to file the encounter. The following informational screen displays the steps to **Load** the unfiled data.





Figure 80: Instructional Steps to Load Unfiled Data

If **Load** or **View** were selected upon entry into the patient's chart and the provider wants to delete the unfiled data after viewing it, use the **Delete Checked** button in the **Unfiled Data by Provider** report. If the **Load** button was selected, then the user also needs to select the **Refresh Patient Chart** submenu from the **File** menu.

**Note:** All data, including unfiled data saved or filed for ‘test’ patients configured with the first three digits of the SSN as zeros will NOT appear in any DRM Plus report.

Unfiled data becomes inactive after eight calendar days. The end-user receives a screen message on the ninth day after saving data whenever they enter the patient’s chart. This message provides two button options to either **View** or **Delete** the inactive unfiled data.

The following screen displays when the patient’s chart which has inactive unfiled data by the provider who saved the unfiled data or from another DRM Plus provider with access to the patient.



**Figure 81: Previously Saved Inactive Data**

The **View** button takes the user directly to the **Unfiled Data by Provider** report, where they can view, print or delete the inactive data. There is no way to load inactive unfiled data into the patient’s chart except to re-enter all the data manually.

The **Delete** button deletes the patient’s inactive unfiled data from the VistA scratch pad.

## Service Reports

Use the **Service Reports** submenu from the **Report** menu to select and create a service report.

**Service Report Selection**

Service Reports

☒ All Reports

☐ 1: Observations per Month

☐ 2: Service Profile ADA/CPT Codes

☐ 3: Service by Product Group

☐ 4: Provider Profile ☒ Use Provider Name

☐ 5: Patient Category Profile

☐ 6: Patients by Eligibility Setting

☐ 7: Patients by Eligibility Group

☐ 8: Patients by Dental Patient Category

☐ 9: Outpatient Service by Product Group

☐ 10: Inpatient Service by Product Group

☐ 11: Distribution of Dental Services

Fiscal Year n/a Start Date 6/ 8/2017 End Date 6/ 8/2017

☒ All Stations

☐ All Providers?

Select Provider(s) [Provider data is aggregated on the reports]

Name	ID	Type/Specialty
DRMPROVIDER,DENTI...	03070003	STAFF DENTIST,PROSTHOD

Search for records using

☒ Visit Date ☐ Create Date

OK Cancel

**Figure 82: Service Report Selection Screen**

1. Choose the desired type of service report.
2. Set the **Fiscal Year** or date range, if applicable.
3. Select the date type that is to be represented on the report.
4. Click the **OK** button.
5. The **Service Reports** screen displays with the results.



DRM

Service Reports

7: Patients by Elig Category

8: Patients by Pt Category

9: Outpt Service

10: Inpt Service

11: Distribution of Dental Services

1: Obs/Month

2: Service Profile

3: Service by Group

4: Provider Profile

5: Patient Profile

6: Patients by Setting

Observations per Month 01/02/17 - 06/08/17 STN: All

	UNIQUE	Number of	Encounters	Cumulative	Encounters
MONTH	SSNs/Mnth	Encounters	% of Total	Encounters	Cumulative %
Apr 2017	11	13	100.00	13	100.00

Note: A missing month indicates no data for that month

Selected Provider(s): DRMPROVIDER.DENTIST

Save All to Excel

Save to Excel

Print All

Print

Close

**Figure 83: Service Reports Screen**

- If more than one report, or all report options are checked on the **Service Reports Selection** screen, the reports appear in tabs on the **Service Reports** screen.
- Save the report to Excel or print.

## Data Warehouse Reports

The **Load Data Warehouse Reports** submenu is required to be selected first; opening an informational screen stating Data Warehouse Reports have been loaded. After closing the informational screen and then selecting the Reports menu again will display the **Data Warehouse Reports** submenu option as active.




Figure 84: Reports Menu

The **Data Warehouse Reports** submenu option will connect to specific reports as maintained by the Office of Dentistry. The Office of Dentistry may change the Warehouse Reports XML hyperlink at any time. The Data Warehouse Report submenu option is dynamically created based on settings of web addresses in the **Administrative Toolbox** settings in DRM Plus.

The DRM Plus Administrator can maintain the web addresses for the **Data Warehouse Reports**. The website addresses should end in '.xml'. These XML addresses/files build the structure for the **Data Warehouse Reports** submenu options.

When selecting the **Load Data Warehouse Reports** submenu there may be a few seconds delay. Once the submenu options have been loaded, there should no longer be a delay. If there is a failure to load the submenu options, the website is down or you have lost connection; the next time you try the **Load Data Warehouse Reports** submenu, it will try to load the submenu options again.

So, in the case of a complete failure and a setting of 0-seconds on the timeout – you may still experience a delay of 20-seconds or so when clicking on that menu. The **Data Warehouse Reports** submenu option would then display 'currently unavailable – try again later' as the submenu. It is suggested that you wait and try later in the day or even possible the next day. If the connection failure persists, please contact the local DRM Plus Administrator.




Figure 85: Data Warehouse Reports Submenu Option Failure

## Adverse Events

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Device Tracking

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Extract History File

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Queued Extract History File

This is an administrative function. For more information, please see the DRM Plus Administrator Manual or speak to a local DRM Plus Administrator.

## Help

Use the **Help** menu to access more information on DRM Plus. There are seven submenus: **Contents**; **Version Release Notes**; **Last Broker Call**; **VA Intranet Website**; **Find your DRM Plus Administrators**; **Have a Question, Comment or Suggestion about DRM?**; and **About...**



Figure 86: Help Menu

## Contents

Use the **Contents** to find information on using DRM Plus.





Figure 87: DRM Plus Contents Screen

## Version Release Notes

Select this submenu to view what was introduced in the current version of DRM Plus.



1. Report Column adjustment corrected.
2. Requirement of D0120 and D0150 PSR added including edentulous exemption.

#### M corrections to the DRM Plus Application

1. Modified the coding check error message for D9630.
2. Changed the RVU value for CPT code 94760 from 5 to 0.
3. Added a system parameter, DENTV PATIENT FLAGS.


Print

Close

Figure 88: Release Notes Screen

## Last Broker Call

Select this submenu to see the last broker call.



**Figure 89: Broker Calls Screen**

Use the up and down arrows next to the **Maximum Calls** text box to adjust how many broker calls are retrieved. Use the previous and next arrows to scroll through the broker calls. Use the **Show All** button to scroll through the list of broker calls.

## VA Intranet Website

Select the submenu to access the VA Intranet website.

**Note:** Clicking this submenu connects the user to the VA Intranet website, where additional dental-related information can be obtained. This includes the latest DRM Plus manuals.

## Find your DRM Plus Administrators

Select the submenu to display a list of all the DRM Plus administrators that have full admin functionality. The list of DRM Plus administrators may be from the local dental clinic or every dental clinic that shares the same VistA server.

## Have a Question, Comment or Suggestion about DRM?

Select the submenu and it will display a VA email screen with the **VHA Dental Software Support** group pre-selected in the 'To' email field. The provider will receive a reply from the **VHA Dental Software Support** group in response from the provider's question, comment or suggestion when requested.

## About

The **About** screen contains information on the DRM Plus application currently in use at the facility, including the version number of the executable.





Figure 90: DRM Plus About Screen

The **System Report** button and the **Reporting Off** button are located on the **About** screen.

Both buttons are useful for **Document Storage Systems** (DSS) in helping to evaluate any issue that may have occurred in DRM Plus. It is strongly suggested NOT to select either button except when working with a DSS employee who will be helping the provider in finding details about the issue occurring in the DRM Plus application.

The **System Report** button will cause the DRM Plus application to freeze. When selecting this button any unfiled data entered this session and NOT saved as unfiled data will be lost. The provider will NOT be able to recover the unsaved unfiled data and would have to re-enter all the data after restarting the DRM Plus application.

Selecting the **System Report** button will give the provider an option NOT to continue if they desire. The following informational screen will inform the provider to select the **No** button if they would prefer NOT to lose any unsaved unfiled data recently entered this session.



**Figure 161: System Report Informational Screen**

Selecting the **Yes** button from the informational screen will freeze the DRM Plus application and give the provider an error system report. All unsaved unfiled data at this point will be lost. The provider may select the **restart application** button to re-launch DRM Plus and continue with their dental charting.

When the second button displays **Reporting Off**, the default setting, then the madExcept application, system report, is inactive. The madExcept application is active and will NOT affect the provider's work process during the session unless they experience an error using the application.

When the second button is changed to **Reporting On**, the madExcept application or system report, is inactive. The only visible change the provider may notice when the **Reporting Off** button is selected will be the changing the button's name to **Reporting On**.

**Note:** The **Reporting On** button's name will change back to the default **Reporting Off** button's name each time EDRM is closed and relaunched.

**Note:** The madExcept system report is explained in detail when reading the **madExcept Application** appendix from the EDRM manual.



# DRM Plus Banner

The DRM Plus banner contains vital information about the patient, providers, also coding standards and alerts.




Figure 91: DRM Plus Banner

## Patient Information

Patient information displays on the far-left side of the DRM Plus banner.




Figure 92: Patient Inquiry Button

The patient information area shows the patient's first and last name, social security number, date of birth and current age. Click the **Patient Inquiry** button of the banner to open the **Patient Inquiry** screen.




Figure 93: Patient Inquiry Screen

The **Patient Inquiry** screen contains more detailed information about the patient; including address, phone number, means test information, status, and admissions information.

## Visit Information

Visit information is displayed on the banner in the second box from the left. It contains information on the current encounter. A scheduled appointment automatically fills the field when the provider enters data into DRM Plus on the same day of the appointment.

When no scheduled appointment automatically imports into banner or if one is NOT manually selected from the banner after opening the dental chart then the following will result. After the first completed procedure code has been selected during the current session results in the launch of the **Provider and Location for Current Activities** screen. This will force the dental user to select a scheduled appointment before the first diagnosis code can be selected for the completed transaction.

### Visit Not Selected

DRMPROVIDER,ADMINDENTIST

**Figure 94: Provider/Visit Button**

To change visit or provider information, click the **Provider/Visit** button of the banner. The **Provider and Location for Current Activities** screen displays.

**Figure 95: Provider and Location for Current Activities Screen**

Select the provider in the top part of the screen and the correct appointment information in the bottom part of the screen. The information in the banner changes to reflect the adjustments made on this screen. For more information on navigating this screen, please see the Dental Encounter Data section in the Using the DRM Plus Drop-Down Menus chapter of this manual.

## Dental Provider Information

This section of the banner displays information on the primary and secondary dental providers. To assign a primary and/or secondary dental provider for a patient requires an administrative parameter option, given by the DRM Plus Administrator. Primary/Secondary providers are for planned care, and show who is responsible for a given patient, and what is upcoming for the patient, regardless of who entered the information.

PDP: DRMPROVIDER,ADMINDE  
SDP:

**Figure 96: Designate Dental Providers Button**

The primary dental provider or secondary dental provider can be set as **Fee Basis** in DRM Plus. The **Fee Basis** provider does NOT exist in VistA. This is a “free-text” entry and therefore, most reports in DRM Plus do NOT recognize the **Fee Basis** provider.

**Note:** This option is utilized for these DES reports: **Recare, Provider Planning, Planned Items List, Planned Non-VA Care** and **Active Patients by Provider**.


## Vitals Button

This section of the banner contains a **Vitals** button to launch the **Enter Vitals** screen which will add specific vitals to the patient’s VistA database. Only the blood pressure, pulse and general pain maybe added using the **Enter Vitals** screen.



**Figure 97: Vitals Button**

After selecting the **Vitals** button from the banner, the **Enter Vitals** screen will display for the provider to enter and save three specific vitals to the VistA database. After entering any of the three possible vitals select the **OK** button to save to the VistA database.

A screenshot of the 'Enter Vitals' dialog box. It has a title bar with 'DRM Plus' and 'Enter Vitals'. Inside, there are three input fields: 'Blood Pressure (mmHg):', 'Pulse (BPM):', and 'General Pain (0-10):'. The 'General Pain' field has a dropdown arrow. At the bottom right are 'OK' and 'Cancel' buttons.

**Figure 98: Enter Vitals Screen**

## Adverse Events Button

The **Adverse Events** (AE) button in the banner when selected will launch a screen where an adverse event may be entered and filed with an encounter for the selected patient's chart.




<b>DRMPATIENT,ONE</b> 000-00-0001 01/01/1960 Age: 57	<b>V: Jun 08, 2017@10:00 DENTAL</b> DRMPROVIDER,DENTIST	PDP: DRMPROVIDER,ADMINDE SDP:	  Dental Class OPC, Class IV
---	--	----------------------------------	---

Figure 99: Adverse Events Button

The **AE** button, displayed as a yellow button in the banner screen above, only results as a yellow icon when at least one **Active** adverse event entry has been filed with the patient's chart. Selecting the **AE** button will open the following **Adverse Events** screen. When only filed **Resolved** entries are listed in the **History** section of the **Adverse Events** screen then the **AE** button in the banner will remain the same background color as the banner.

### To add an adverse event:

1. Select the drop-down arrow on the **Select Category and Provider Details** window
2. Select one of the twelve categories from the list
3. View the blue **Categories with Examples** window with the category listed at the top
4. Add or provide detailed text into the upper text window about the adverse event
5. Select the correct **Status** radio button (**Resolved** or **Remain Active**)
6. Select the **Add** button to temporary save the new adverse event
7. May add multiple adverse events this session if desired
8. When all adverse events have been added select the **Done** button to save as unfiled data
9. Filing an encounter is required when adding any adverse event during this session



Visit	Provider	Category	Status	Details
MAR 17, 20...	DRMPROVIDER.ADM...	Hard Tissue Damage	Active	Bone/mandible fr...
MAR 17, 20...	DRMPROVIDER.ADM...	Allergy/Toxicity/For...	Resolved	Nitrous oxide toxi...

Figure 100: Adverse Events Screen

**Note:** The **Add** button from the **Select Category and Provide Details** section will temporary save each adverse event that may be entered this session. The name on the button will change to **Save Changes** when a filed adverse event is selected from the **History** table and edited by the provider. After all the adverse events have been added/edited/deleted during this session select the **Done** button to save all adverse events as unfiled data.

**Note:** There is one trigger, completed procedure **D9930**, which will automatically open the **Adverse Events** screen when selected this session. The trigger occurs when the provider selects the **Next** button which will automatically open the **Adverse Events** screen. The completed procedure D9930 is only suggested to include an adverse event associated with the procedure however **NOT** required for an adverse event to be entered and filed. Selecting the **Cancel** button will allow the provider to complete the encounter without adding an adverse event.

### To edit a filed Active adverse event:

1. Highlight the filed **Active** adverse event entry from the **History** table
2. May edit or add to the **Provider Details** in the upper text window
3. Or may change the **Status** radio button (to **Resolved**)
4. Select the **Save Changes** button to temporary save the edited adverse event
5. Select the **OK** button on the informational **Adverse event edited** screen
6. When all adverse events have been added select the **Done** button to save the unfiled data
7. Filing an encounter is required when editing any adverse event during this session

The screenshot shows the 'Adverse Events' window. At the top, there's a 'Select Category and Provide Details' section with a dropdown menu set to 'Allergy/Toxicity/Foreign Body Response'. Below it is a text area containing 'Nitrous oxide toxicity during the extraction.' and 'I am adding more details about this adverse event and reviewing the status.' To the right of the text area are 'Status' radio buttons for 'Resolved' (selected) and 'Remain Active', and an 'Examples' button. Below the text area are checkboxes for 'Display Resolved' (checked) and 'Display Deleted' (unchecked), and a 'Save Changes' button. On the right side, there's a 'Categories with Examples' panel listing various categories and their examples, such as 'Allergy/Toxicity/Foreign Body Response examples', 'Aspiration or Ingestion of Foreign Body examples', 'Bleeding examples', 'Hard Tissue Damage examples', 'Infection examples', and 'Nerve Injury examples'. At the bottom, there's a 'History' table with columns for 'Visit', 'Provider', 'Category', 'Status', and 'Details'. The table shows two entries: one for 'MAR 17, 20...' with 'DRMPROVIDER.ADM...' and 'Hard Tissue Damage' (Active), and another for 'MAR 17, 20...' with 'DRMPROVIDER.ADM...' and 'Allergy/Toxicity/For...' (Resolved). At the bottom of the window are 'Delete', 'Done', and 'Cancel' buttons.

Visit	Provider	Category	Status	Details
MAR 17, 20...	DRMPROVIDER.ADM...	Hard Tissue Damage	Active	Bone/mandible fr...
MAR 17, 20...	DRMPROVIDER.ADM...	Allergy/Toxicity/For...	Resolved	Nitrous oxide toxi...

**Figure 101: Adverse Events Screen**

**Note:** Only a filed **Active** adverse event may be edited. When trying to edit a filed **Resolved** adverse event the provider will get a screen informing them, they can **NOT** edit any item that is not **Active**.

**Note:** Only one **Active** adverse event may be edited and filed per encounter. A second **Active** adverse event that should be edited requires to be filed with a different encounter visit.

The **Display Resolved** check box is defaulted as checked when opening the **Adverse Events** screen. The provider may uncheck the **Display Resolved** check box if the provider doesn't want the **History** table to display any resolved adverse events. The provider may select the **Display Deleted** check box to see if there are any deleted adverse events filed for the patient's chart. Any newly entered resolved and deleted adverse event will be required to be filed with an encounter before that adverse event will be permanently viewable with either of these check boxes.

### To delete a filed adverse event:

1. Highlight the adverse event in the **History** table
2. Select the **Delete** button to temporary save the deleted adverse event
3. Select the **OK** button on the informational **Adverse event deleted** screen
4. When all adverse events have been entered select the **Done** button to save the unfiled data
5. Filing an encounter is required when deleting an adverse event this session

The **Delete** button selection for any filed adverse event may only be completed by a full DRM Plus Administrator. The deleted adverse event requires an encounter to be filed to complete the process. When selecting the **Delete** button as a non-admin provider the following screen will display.




Figure 102: Not DRM Plus Administrator Informational Screen

**Note:** All deleted adverse events will display as a line through the entry in the **History** table when the **Display Deleted** check box is selected by any DRM Plus provider that would want to view them.

When the **Cancel** button is selected before the **Done** button while the **Adverse Events** screen is open then the add/edit/delete entries will NOT be saved. The **Cancel** button will undo any action taken with any new entries by the provider when adding a new adverse event, modifying a filed adverse event, or deleting a filed adverse event. All recent adverse event data entered while the **Adverse Events** screen was open will be lost and will need to be re-entered again if the **Cancel** button is selected.

**Note:** When the **Done** button has been selected and any saved adverse events were entered in error during this session follow one of the methods that follow to remove them. The only way to remove the adverse event this session would be by selecting the **Delete** button from the **Adverse Events** screen and then the **Done** button or refresh the patient's chart and don't save the unfiled data that has been entered up to this time.

Each header on the **History** table is selectable for sequence reorder of the data listing in each column. In the next dialog the **Visit** column selected would display the most recent visit at the top with the older visits following in sequence.

History				
Visit	Provider	Category	Status	Details
MAR 02, 2020	DRMPROVIDER.ADM...	Wrong Site/Patient/...	Active	Wrong tooth or sit...
FEB 28, 2020	DRMPROVIDER.ADM...	Bleeding	Resolved	Vascular damage...

**Figure 103: History Section on Adverse Events Screen**

The **Examples** button will open and close the right half of the screen labeled **Categories with Examples** in the header. This blue view-only window allows the provider to view common examples associated with each category that is listed as an adverse event. The common examples in the blue window which may be viewed are NOT an exclusive list.

The next dialog does NOT display the blue **Categories and Examples** view-only window because the **Examples** button was selected. Select the **Examples** button again and the blue view-only window will display.

**Adverse Events** [?] [X]

Select Category and Provide Details:

Status: ☐ Resolved ☐ Remain Active **Examples**

☒ Display Resolved ☐ Display Deleted **Add**


History				
Visit	Provider	Category	Status	Details
MAR 17, 20...	DRMPROVIDER.ADM...	Hard Tissue Damage	Active	Bone/mandible fr...
MAR 17, 20...	DRMPROVIDER.ADM...	Allergy/Toxicity/For...	Resolved	Nitrous oxide toxi...

**Delete** **Done** **Cancel**

**Figure 104: Adverse Events Screen without Categories and Examples**

## Device Tracking Button

The **Device Tracking** (DT) button in the banner when selected will launch a screen allowing device tracking identifiers to be entered and then filed with an encounter for this patient's chart.



**Figure 105: Device Tracking Button on DRM Plus Banner**

Selecting the **Device Tracking** (DT) button, only results as a yellow button in the banner when at least one device completed procedure has been filed without any device identifiers attached, will open the following **Device Tracking** screen. When all filed device completed procedures have device identifiers filed to them in the **Device Tracking** screen then the **DT** button in the banner will remain the same background color as the banner.

[illegible]

**Figure 106: Device Tracking Screen**

From the **Device Tracking** screen, a provider may add device identifier(s) with a completed procedure transaction, a diagnostic finding transaction, or a **No Transaction** by using the **Save Entry** button. The previous dialog displays no current session device procedures entered. There are 26 different completed procedures that may be displayed in this screen if entered this session or they have been filed in the past on the patient's chart.



The device tracking **Transactions** table may display four different views. The other three views occur when selecting the following check boxes; **Display All Findings**, **Display All Completed Care**, or **Display All “No Transactions”**. The check box selections may be individual or accumulative. All four different views may have device identifier(s) filed/attached to a specific procedure transaction or non-transaction.

The following options explain the details about the four views that may result when selecting a check box or combination of check boxes.

- **No check boxes** selected, which is the default selection when opening the **Device Tracking** screen, will display all the unfiled/ filed 26 trigger completed procedures suggested to have device identifiers.
- Selecting the **Display All Findings** check box displays all the unfiled/ filed diagnostic finding transaction(s) for the patient’s chart.
- Selecting the **Display All Completed Care** check box displays all the unfiled/ filed completed care transaction(s) for the patient’s chart.
- Selecting the **Display All “No Transactions”** check box, inactive until the first no transaction device has been added and filed, displays all the filed **No Transactions** for the patient’s chart.
- Any combination of one or three check boxes selected would result in the **Transactions** table displaying the accumulation of transactions/non-transactions from the check boxes that were selected.

### **To add completed procedure device identifiers:**

1. Enter one of the 26 trigger completed procedures from the **Treatment & Exam** screen this session
2. Select the **DT** button from the banner or select the **Next** button when completing the encounter
3. Highlight the completed procedure transaction in the **Transactions** table to add device identifiers
4. Required to enter at least one device identifier with this completed procedure transaction
5. Select the **Save Entry** button to save the new device identifiers to the completed procedure transaction
6. Select the **OK** button on the informational **Device identifiers saved** screen
7. Select the **Done** button to save all the device identifier(s) during this session
8. Filing an encounter is required when adding device identifier(s) this session

**Device Tracking**

**Transactions**

Date	Code	Description	Tooth	Surf
3/22/2017	D6010	ENDOSTEAL IMPLANT BODY PLACE	6	

☐ Display All Findings
 ☐ Display All Completed Care
 ☐ Display All "No Transactions"

No Transaction: 
 Date: 
 Description: 
 Tooth:

**Device Identifiers**

UDI:

Type:

Subtype:

Manufacturer:

Reference:

Lot:

Serial Number:

Details:

☐ Removed

[FDA website on UDI Basics.](#)

**Device Instructions and Examples**

Highlight the transaction that placed the device. Enter the Unique Device Identifiers (UDI) and click the lookup to fill in the details. Then click Save Entry.

The easiest way of recording the UDI is to scan it from the barcode labeled on the packaging. There are several types of barcodes. The DataMatrix (similar to the QRcode) has become more popular recently.

There are 3 issuing agencies for UDI. GSI, HIBCC, and ICCBBA. Sometimes the UDI barcode will be labeled with one of those identifiers.

Lookup values are returned from the U. S. National Library of Medicine. Retrieval may require several lookup attempts depending on network traffic.

If there is not a successful response, check to see that the correct barcode/UDI is being scanned or entered. If these are correct, then it is possible the information is not available from the National Library of Medicine and the information will need to be entered manually.


UDI - unique device identification  
 DI - device identifier fixed portion of UDI  
 PI - production identifier lot, serial, expiration, manufacturer  
 UDI = DI + PI  
 example:  
 (01)00614141999996(17)100101(10)123ABC(21)1234567890

DI = (01)00614141999996  
 PI = (17)100101(10)123ABC(21)1234567890

**Figure 107: Device Tracking Screen with Completed Transaction Entry**

Any unfiled/filed **Diagnostic Finding** transaction or any other unfiled/filed completed procedure transaction, which is **NOT** one of the 26 trigger completed procedures, may have device identifier(s) added and filed following the steps listed previously. Use the two check boxes at the bottom of the **Transactions** table; **Display All Findings** and **Display All Completed Care** to display the transactions to be associated with new **Device Identifiers**.

The easiest way for recording the **UDI** (unique device identification) is to scan it using the barcode labeled on the packaging. There are several types of barcodes. The DataMatrix (like the QRcode) has become more popular recently. There are three issuing agencies for **UDI** details: GSI, HIBCC, and ICCBBA. Sometimes the **UDI** barcode will be labeled with one of those identifiers.

The **UDI** icon  at the end of the **UDI** field in the **Device Identifiers** window allows importing of a few or many **Device Identifiers**. After the UDI has been entered the provider may select the **Enter** key or select the **UDI** icon to import the field values found in the UDI. The imported field values come from the U.S. National Library of Medicine. The retrieval may require several import attempts depending on the network traffic.

**Note:** The drop-down field lists of (**Types:**) and (**Sub-Types:**) in the **Device Identifiers** section is information obtained from the FDA.

When the **Cancel** button is selected before the **Done** button while the **Device Tracking** screen is open then the add/edit/delete entries will **NOT** be saved. The **Cancel** button will undo any action taken with any new entries from the provider when adding new device identifiers, modifying filed device identifiers, or deleting filed device identifiers. Any recent device identifier data entered this session will be entirely lost and will need to be re-entered again.

Under certain circumstances there may be a reason why device identifiers should NOT be filed with a dental transaction. The **No Transaction** section, right below the **Transactions** section, which includes the **Add** button on the **Device Tracking** screen will allow the following.

### To add non-transaction device identifiers:

1. Select the **DT** button in the banner
2. Select the **Add** button from the **No Transaction** section
3. Data may be entered into **Date**, **Description** or **Tooth** fields
4. **Description** field is the only required data needed for entry
5. Required to enter at least one device identifier with the non-transaction entry
6. Select the **Save Entry** button to save the new device identifiers to the non-transaction entry
7. Select the **OK** button on the informational **Device identifiers saved** screen
8. Select the **Done** button to save all the device identifier(s) during this session
9. Filing an encounter is required when adding device identifier(s) this session

When adding a new **No Transaction** entry after selecting the **Add** button; the button will become a **Cancel** button. Selecting the **Cancel** button will NOT save any data entered in **No Transaction** fields or device identifier fields. Any **No Transaction** field data and device identifier data will need to be re-entered.

The next **Device Tracking** screen will display the **No Transaction** filed in the **Display All “No Transactions”** view screen. Highlighting the non-transaction entry displays the device identifiers filed with the non-transaction entry.

The **Transactions** table displays a green (checkmark) icon in the first column which identifies the non-transaction entry having saved unfiled **Device Identifiers**.

The screenshot shows the 'Device Tracking' application window. It features a 'Transactions' table with columns for Date, Code, Description, Tooth, and Surf. A green checkmark icon is visible in the first column of the first row. Below the table, there are checkboxes for 'Display All Findings', 'Display All Completed Care', and 'Display All "No Transactions"'. The 'No Transaction' section includes an 'Add' button and fields for Date, Description, and Tooth. The 'Device Identifiers' section contains fields for UDI, Type, Subtype, Manufacturer, Reference, Lot, and Serial Number, along with a 'Details' text area and 'Delete Entry' and 'Save Entry' buttons. A 'Device Instructions and Examples' panel on the right provides guidance on entering UDIs and includes a link to the FDA website on UDI Basics.

	Date	Code	Description	Tooth	Surf
✓	03/22/2017		Bone graft material on mandibular	30	

☐ Display All Findings   
 ☐ Display All Completed Care   
 ☒ Display All "No Transactions"

**No Transaction**  
 Add    Date: 03/22/2017    Description: Bone graft material on mandibular    Tooth: 30

**Device Identifiers**  
 UDI:    Type: Prosthetic Devices    Subtype: Bone grafting material  
 Manufacturer:    Reference:    Lot:    Serial Number:  
 Details: Bone graft material on mandibular around the tooth #30 region by fee dentist.  
☐ Removed




Delete Entry    Save Entry



Done    Cancel    [FDA website on UDI Basics.](#)


**Device Instructions and Examples**  
 Highlight the transaction that placed the device. Enter the Unique Device Identifiers (UDI) and click the lookup to fill in the details. Then click Save Entry.  
 The easiest way of recording the UDI is to scan it from the barcode labeled on the packaging. There are several types of barcodes. The DataMatrix (similar to the QRcode) has become more popular recently.  
 There are 3 issuing agencies for UDI. GS1, HIBCC, and ICCBBA. Sometimes the UDI barcode will be labeled with one of those identifiers.  
 Lookup values are returned from the U. S. National Library of Medicine. Retrieval may require several lookup attempts depending on network traffic.  
 If there is not a successful response, check to see that the correct barcode/UDI is being scanned or entered. If these are correct, then it is possible the information is not available from the National Library of Medicine and the information will need to be entered manually.  
 UDI - unique device identification  
 DI - device identifier    fixed portion of UDI  
 PI - production identifier    lot, serial, expiration, manufacturer  
 UDI = DI + PI  
 example:  
 (01)00614141999996(17)100101(10)123ABC(21)1234567890  
 DI = (01)00614141999996  
 PI = (17)100101(10)123ABC(21)1234567890


Figure 108: Device Tracking Screen with No Transaction Entry


**Note:** The only way to delete any filed **No Transaction** entry must be done in the **Device Tracking** screen. The **No Transaction** may have been filed by mistake on a patient's chart or filed to the wrong patient's chart. Deleting the encounter on the wrong patient's record will NOT delete the **No Transaction** entry if one was filed with that encounter.

The following explains the meaning of the icons found in the **Device Tracking** screen. The first column icons from any of the multiple Transactions views may be displayed as either a green  icon, red  icon or a blue  icon.




The green  icon represents unfiled/filed device identifier(s) filed with a transaction or non-transaction found in all four views. The green  icon also may represent a filed transaction or non-transaction as a **Removed** device tracking entry.

The red  icon represents any of the 26 trigger completed procedures having no device identifier(s) filed with the completed procedure transaction. The red  icon will only be displayed in the default view or when the **Display All Completed Care** check mark has been selected.

The blue  icon represents all other unfiled/filed completed procedures or diagnostic findings. These completed procedures/findings are left up to the provider if they would prefer to file a device identifier(s) with the transaction. These completed procedure/finding transactions are NOT identified for any device tracking entry by the VHA Office of Dentistry.

**Note:** The red  icon will display for any of the 26 unfiled/filed trigger completed procedures identified by the VHA Office of Dentistry. When one of the 26 trigger completed procedures is entered during the current session it will force the provider to the Device Tracking screen when the provider selects the **Next** button when trying to complete the encounter. Then it will be up to the provider to decide if they want to save any device identifiers to the unfiled completed transaction or NOT. The **Cancel** button on the Device Tracking screen may be selected to complete the encounter by NOT entering any device identifiers to the trigger completed procedure. The **DT** icon button in the banner will also display as a yellow icon informing the provider at least one trigger completed procedure having no device identifiers associated with it is in the **Device Tracking** screen.

**Note:** The 26 unfiled/filed trigger completed procedures identified by the VHA Office of Dentistry are suggested to have device identifiers associated with them however they are NOT required to have device identifiers.

Each header at the top of the **Transactions** table is selectable for sequence order of data listed in that column. In the next dialog the **Icon** column, first column, selected displays the red  icon, green  icon and then blue  icon from top to bottom.

Transactions						
	Date	Code	Description	Tooth	Surf	
	3/22/2017	D3428	BONE GRAFT W/PRAD SURG-FIRST	3		
	3/22/2017	D6010	ENDOSTEAL IMPLANT BODY PLACE	6		
	3/22/2017	D1110	DENTAL PROPHYLAXIS ADULT			
	3/22/2017		Missing	16		
	3/22/2017		Missing	11		
	3/22/2017		Missing	6		
	3/22/2017		Missing	1		
	3/17/2017	D0210	INTRAORAL FULL IMAGE SERIES			

**Figure 109: Transaction Section on Device Tracking Screen**

Here is how to **Remove** filed device identifiers for a transaction or non-transaction. The provider will be required to highlight the filed completed procedure transaction, the filed diagnostic findings transaction, or the filed non-transaction to label the device identifiers as removed.

### **To remove filed device identifiers for a transaction or non-transaction:**


1. Find the filed transaction or non-transaction that should be marked as **Removed**
2. Highlight the filed transaction or non-transaction listed on the **Transactions** table
3. Select the **Removed** check box in lower left corner of the **Device Tracking** screen
4. Select the **Save Entry** button to save the transaction or non-transaction marked as **Removed**
5. Select the **OK** button on the informational **Device identifiers edited** screen
6. Select the **Done** button to save all the **Device Tracking** entries during this session
7. Filing an encounter is required when removal of a device tracking entry this session

Here is how to delete device identifiers that were filed to a transaction or non-transaction. The provider will be required to highlight the filed completed procedure transaction, the filed diagnostic findings transaction, or the filed non-transaction that should be deleted.

### **To delete the filed device identifiers for a transaction or non-transaction:**

1. Find the filed transaction or non-transaction where the device identifiers should be deleted
2. Highlight the filed transaction or non-transaction listed on the **Transactions** table
3. Select the **Delete Entry** button to save the deleted device identifiers for transaction or non-transaction
4. Select the **OK** button on the informational **Device identifiers deleted** screen
5. Select the **Done** button to save all the **Device Tracking** entries during this session
6. Filing an encounter is required when deleting device identifiers this session

Only a full DRM Plus Administrator can delete device identifiers filed to a transaction or non-transaction. When selecting the **Delete** button as a non-admin provider the following informational screen will display.



**Figure 110: Not DRM Plus Administrator Informational Screen**

The **Examples** button will display and hid the right half of the screen labeled **Device Instructions and Examples** as the header. This blue view-only window allows the provider to view device instructions and examples of the **UDI** (unique device identification) information.

The next **Device Tracking** screen doesn't display the blue **Device Instructions and Examples** view-only window because the **Examples** button was selected. Select the **Examples** button again and the blue view-only window will display.

 A screenshot of the "Device Tracking" window. The window has a title bar with "DRM" and "Device Tracking" and standard window controls. The main area is divided into two sections. The top section, titled "Transactions", contains a table with columns: Date, Code, Description, Tooth, and Surf. Below the table are three checkboxes: "Display All Findings", "Display All Completed Care", and "Display All 'No Transactions'", each followed by a "No Transaction" label. Below these is an "Add" button and input fields for "Date:", "Description:", and "Tooth:". The bottom section, titled "Device Identifiers", contains several input fields: "UDI:", "Type:" (with a dropdown arrow), "Subtype:" (with a dropdown arrow), "Manufacturer:", "Reference:", "Lot:", "Serial Number:", and "Details:" (with a scroll bar). To the right of these fields are buttons for "Examples", "Delete Entry", and "Save Entry". At the bottom left of this section is a "Removed" checkbox. At the very bottom of the window are "Done" and "Cancel" buttons.


**Figure 111: Device Tracking Screen without Device Instructions and Examples**

**Note:** The [FDA website on UDI Basics](#) address found on the bottom of the **Device Instructions and Examples** view-only window will launch the U.S. FOOD & DRUG ADMINISTRATION website that will explain the UDI, unique device identifier, in greater detail.

**Note:** If there is not a successful response, check to see that the correct barcode/UDI is being scanned or entered. If these are correct, then it is possible the information is not available from the National Library of Medicine and the information will need to be entered manually.

## Dental Class Information

This section of the banner contains information on the patient's dental class.




**Figure 112: Dental Class Button**

Click the **Dental Class** button of the banner to go to the DRM Plus **Cover Page** tab. However, only a DRM Plus Administrator or an end-user that has this administrative parameter option can change a patient's dental class. Please see the Cover Page chapter of this manual for more information.

## Clean Slate

This section, when present, is located between the **Dental Class** button and the icons on the far right of the banner. **Clean Slate** is only a viewable window that displays the most recent date that a clean slate was performed on the patient's chart. Clean slate may only be performed by a DRM Plus Administrator, or an end-user who has received the administrative parameter option.



**Figure 113: Clean Slate Notification**

The clean slate functionality has been added to clear the graphical portion of the **Treatment & Exam** screens in DRM Plus and delete all planned treatment for the selected patient. The recent clean slate can be restored for this patient at any time until a new encounter has been filed.


## Icons

The icons located on the right side of the banner show patient flags, alerts and provide an easy way for the user to look up general coding standards.



**Figure 114: DRM Plus Banner Icons**

## General Coding Standards


The first icon is the **General Coding Standards**  icon. Clicking this icon takes the provider to the VA Office of Dentistry website where current **General Coding Standards** may be viewed. This hyperlink should be changed by the DRM Plus Administrator if requested to do so by the VA Dental Informatics and Analytics Director.

**Note:** If NOT immediately directed to the website after clicking the icon please notify the Help Desk.

## Patient Flags

There are four possible patient flags that can display in DRM Plus, including: **Clinical Reminders**, **Consult**, **Exam Quality Indicator** and **Fluoride Quality Indicator**.

### Clinical Reminders

The **Clinical Reminders**  icon displays on the right side of the DRM Plus banner when there are clinical reminders due for the selected patient. Providers must still process clinical reminders using CPRS.

Clicking the **Clinical Reminders** icon displays an information screen, stating that the selected patient has clinical reminders due.






Figure 115: Clinical Reminders Due

The **Clinical Reminders** icon should only display if the current end-user is responsible for and may resolve the clinical reminder(s) listed on the CPRS cover sheet. If clinical reminder(s) displaying in the list cannot be resolved by the end-user, contact local IT for assistance.


### Consult

The **Consult**  icon displays when the patient has an incomplete consult in their chart. Click the icon. The **Clinical Record** tab displays with consultations selected. For more information, see the Clinical Record chapter of this manual.

### Exam Quality Indicator

The **Exam Quality Indicator**  icon displays when the patient is due for a follow-up exam. The icon only displays if the patient meets certain class restrictions. Clicking this icon produces a screen, which reveals when the patient last received a qualifying exam.







**Figure 116: Exam Quality Indicator Screen**

Click the **OK** button and the informational screen will close. For more information, see the Cover Page chapter of this manual.

**Note:** The **Exam Quality Indicator** and **Fluoride Quality Indicator** only are active for VA dental care under **Class I, IIA, IIC** or **IV** classification. These dental classes are eligible for any necessary dental care to maintain or restore oral health and masticatory function, including repeat care.

### Fluoride Quality Indicator

The **Fluoride Quality Indicator**  icon displays when the patient is due for a fluoride intervention. Clicking this icon displays an explanation of why the patient is at risk and needs intervention.




**Figure 117: Fluoride Quality Indicator Screen**

Click the **OK** button. The **Cover Page** tab displays. For more information, see the Cover Page chapter of this manual.


**Note:** The **Exam Quality Indicator** and **Fluoride Quality Indicator** only are active for VA dental care under **Class I, IIA, IIC** or **IV** classification. These dental classes are eligible for any necessary dental care to maintain or restore oral health and masticatory function, including repeat care.

## Alerts

The DRM Plus **Alert**  icon shows if the patient has any associated alerts. It can also be used to add alerts to the patient's chart in DRM Plus.


To add an alert:

1. Click the **Alert** icon.
2. The **Dental Alerts** screen displays.




**Figure 118: Dental Alerts Screen**

3. Click the desired **Common Dental Alerts** check box on the left side of the screen.
4. Click the **right arrow** button.
5. The selected alert(s) display on the right side of the screen.
6. Alternatively, type the alert directly into the **Patient's Dental Alerts** screen window.
7. Click the **OK** button.

The **Alert** icon changes when patient alerts are present and displays as a  stop sign.

To view the alerts:

1. Click the **Alert** icon displaying a stop sign.
2. The **Patient's Dental Alerts** screen displays.



**Figure 119: Patient's Dental Alerts Screen**

3. Click the **Common Alerts** button to expand the **Common Dental Alerts** window.
4. The full **Dental Alerts** screen displays. From here, more dental alerts can be added, or the alert can be removed by placing the cursor into the **Patient's Dental Alerts** window and deleting the text.
5. Click the **OK** button to finish.

# Cover Page

The **Cover Page** tab displays important patient information. The tab has 8 major sections: **Dental Eligibility, Demographics, Case Management, Recent Dental Activity, Fluoride Indicator Prescription Date, Dental Alerts, Notes** and **Filed Planned Care**.

**Dental Record Manager Plus in use by: DRMPROVIDER,DENTIST for DRMPATIENT,ONE**

File Edit Dental Encounter Data Treatment & Exam Tools Reports Help

DRMPATIENT, ONE 000-00-0001 01/01/1960 Age: 57 V: Jun 08, 2017@10:00 DENTAL DRMPROVIDER.DENTIST PDP: DRMPROVIDER.ADMINDE SDP: DRMPROVIDER.DENTIST Dental Class: DPC, Class IV

Cover Page Clinical Record Dental History Chart/Treatment Exam

**Dental Eligibility**

Dental Class: T5-DPC, Class IV

Service Connected Teeth/Service Trauma

Other

Adjuvantive Medical Condition(s) Add/Edit AMC

Eligibility Expiration Date Anticipated Rehab Date Save

**Demographics**

Primary Eligibility: NSC

Service Separation Date: JAN 1, 1980 (MARINE CORPS)

**Case Management**

Status: ☒ Active ☐ Inactive ☐ Maintenance

Suggested Recare Date Save

**Recent Dental Activity**

Last Qualifying Exam:

Last Comprehensive Exam: Apr 28, 2017

Last Brief Exam:

Last Periodontal Exam:

Last Panorax Image: May 26, 2017

Last Full Mouth Image: May 31, 2017

Last Bitewing Image:

Last CBCT Image:

Last Prophylaxis:

Last Visit: Jun 08, 2017

Last Provider: DRMPROVIDER.DENTIST

Fluoride Indicator Prescription Date Save

**Dental Alerts**

Antibiotic Premed

**Notes**

Save

**Filed Planned Care**

Treatment Plan: Non-VA Care (D7950) MANDIBLE GRAFT DX: ()


Sequencing Notes: [09:28 AM 05/25/2017 ADMIN] Adding an additional planned note. Select the F5 key first and then enter text data.

**Figure 120: DRM Plus Cover Page Tab**

**Note:** **Dental Eligibility** is a DRM Plus administrative function. Please contact a DRM Plus Administrator for more information about this option.

## **Dental Eligibility**

The patient's dental class, service-connected teeth/service trauma, adjunctive medical conditions, eligibility expiration date, and anticipated rehab date are displayed in this area. **Dental Eligibility** can only be changed by a DRM Plus Administrator or an end-user who has received this administrative parameter option.



**Figure 121: Dental Eligibility**

### **Dental Class**

**Dental Class** can only be changed by a DRM Plus Administrator or a user who has been given this administrative parameter option.

### **Service Connected Teeth/Service Trauma**

**Service Connected Teeth/Service Trauma** can only be changed by a DRM Plus Administrator or user who has been given this administrative parameter option.

### **Adjunctive Medical Condition(s)**

**Adjunctive Medical Condition(s)** can only be changed by a DRM Plus Administrator or user who has been given this administrative parameter option.

### **Eligibility Expiration Date**


**Eligibility Expiration Date** can only be changed by a DRM Plus Administrator or user who has been given this administrative parameter option.

### **Anticipated Rehab Date**

**Anticipated Rehab Date** can only be changed by a DRM Plus Administrator or user who has been given this administrative parameter option.

## Demographics

Patient demographic information is located here and is imported from VistA. The fields cannot be updated or changed in DRM Plus.




The screenshot shows a window titled "Demographics". It contains two fields: "Primary Eligibility:" with the value "NSC" and "Service Separation Date:" with the value "JAN 1, 1980 (MARINE CORPS)". The rest of the window is empty.

Figure 122: Demographics

## Case Management

Use the **Case Management** (Disposition) section to adjust the patient's status and file suggested recare dates. Click the **Save** button to update any new changes.



The screenshot shows a window titled "Case Management". It has a "Status" section with three radio buttons: "Active" (selected), "Inactive", and "Maintenance". Below this is a "Suggested Recare Date" field with an ellipsis button (...). To the right of the field is a "Save" button.

Figure 123: Case Management

### **Status**


Change the **Status** of the patient by clicking the appropriate radio button.

### **Suggested Recare Date**

1. Click the **Ellipsis (...)** button to display a calendar screen.
2. Select the desired date from the menu and click the **OK** button.
3. Click the **Save** button to file the data.

## **Recent Dental Activity**

This section displays dates, if applicable, when specific completed procedures were last performed on the patient. Procedure codes that activate the date in this section may be viewed by hovering the mouse cursor over the recent dental activity description.




Recent Dental Activity	
Last Qualifying Exam:	Apr 03, 2017
Last Comprehensive Exam:	Apr 03, 2017
Last Brief Exam:	
Last Periodontal Exam:	
Last Panorex Image:	Mar 21, 2017
Last Full Mouth Image:	Apr 03, 2017
Last Bitewing Image:	Mar 09, 2017
Last CBCT Image:	Apr 06, 2017
Last Prophylaxis:	Mar 14, 2017
Last Visit:	Apr 06, 2017
Last Provider:	DRMPROVIDER,ADMIND

**Figure 124: Recent Dental Activity**

## **Fluoride Indicator Prescription Date**

Add a fluoride indicator prescription date here.



Fluoride Indicator Prescription Date

...

**Figure 125: Fluoride Indicator Prescription Date**

1. Click the **Ellipsis (...)** button.
2. Select a date from the calendar screen and click the **OK** button.
3. Click the **Save** button to file the date.

To delete the prescription date:


1. Click the **date** to place the cursor in the field.
2. Delete the date.
3. Click the **Save** button to file.

The field is active regardless of whether the patient has a fluoride quality indicator.

**Note:** No future dates are allowed for the fluoride indicator prescription date.

## **Dental Alerts**

The patient's **Dental Alerts**, if any, are listed here.




**Figure 126: Dental Alerts**

Please see the Alerts section in the DRM Plus Banner chapter of this manual for more information.

## **Notes**

Add general notes in this text box. These notes are NOT imported into the TIU progress note.




**Figure 127: Notes**

1. Place the cursor in the text box and begin typing.
2. Click the **Save** button to file the notes.
3. A screen displays showing that the information is saved. Click the **OK** button.



## **Filed Planned Care**

The filed treatment plan for the patient, if applicable, is displayed here. It cannot be edited on this page.



The screenshot shows a window titled "Filed Planned Care". Inside the window, the text is as follows:


Treatment Plan:  
Non-VA Care  
(D6010) ENDOSTEAL IMPLANT BODY PLACE: 29 DX: ()  
second appt  
(D6056) PREFABRICATED ABUTMENT: 29 DX: ()  
(D6059) ABUT SUPPORT CROWN PORCH HNOB: 29 DX: ()  
Phase 1  
(D2160) AMALGAM THREE SURFACES PERMA: 3 (DOL) DX: () (Next-Apppt)  
(D2392) POST 2 SRFC RESINBASED CMPST: 4 (OL) DX: () (Next-Apppt)  
Phase 2  
(D2392) POST 2 SRFC RESINBASED CMPST: 13 (OL) DX: ()  
(D2160) AMALGAM THREE SURFACES PERMA: 14 (DOL) DX: ()  
Sequencing Notes:

**Figure 128: Filed Planned Care**

Please see the Treatment Plan section in the Chart/Treatment-Treatment & Exam chapter of this manual for more information.


# Clinical Record

The **Clinical Record** tab allows access to various areas of the patient's record. From here, view **Problems, Medications, Vital Signs, Radiology, Allergies, Lab Results, Postings, Immunizations, D/C Summaries, Consultations, Health Summary** and **Notes**.




**Figure 129: Clinical Record Tab**

1. Click the **radio button** that corresponds to the type of clinical record to be viewed.



**Figure 130: Clinical Record Radio Buttons**

2. A list of TIU notes, corresponding to the selected **Notes** radio button, displays in the area below the clinical record radio button section.



**Figure 131: Clinical Record Selection Window**

3. Select a TIU note entry. Details of the TIU note displays on the right side of the screen.

LOCAL TITLE: DENTAL NOTE  
 STANDARD TITLE: DENTISTRY NOTE  
 DATE OF NOTE: APR 28, 2017@15:20 ENTRY DATE: APR 28, 2017@15:20:41  
 AUTHOR: DRMPROVIDER,ADMINDE EXP COSIGNER:  
 URGENCY: STATUS: COMPLETED

Patient Name: DRMPATIENT,ONE, DOB: 01/01/1960, Age: 57  
 Visit: 04/28/2017 15:00.  
 Primary PCE Diagnosis: K02.52 (DENTAL CARIES ON PIT AND FISSURE SURFACE PENETRATING INTO DENTIN).  
 Dental Category: 15-OPC, Class IV. Treatment Status: Active.

Completed Care:  
 (D0330) DENTAL PANORAMIC IMAGE.  
 DX: K02.52 Dental Caries on Pit and Fissure Surface Penetrating into Dentin  
 (D0140) LIMIT ORAL EVAL PROBLEM FOCUS.  
 DX: K02.52 Dental Caries on Pit and Fissure Surface Penetrating into Dentin

Presentation/Chief Complaint:  
 Patient presents for limited oral evaluation  
 Patient has no dental complaints

Vital Signs:  
 Vital signs not obtained

Past Medical History and Medications:  
 No significant changes since the last dental visit

Intraoral and Extraoral Screening Exam Findings:  
 04/28/2017 Cancerous condition

Radiographic Findings:  
 Radiographic findings consistent with charted entries

Oral Examination:  
 Oral Health Assessment Findings:  
 Plaque Index: 1 - Slight  
 Xerostomia: 0 - None  
 Caries Risk: 1 - Low  
 Oral Hygiene: 1 - Good

Dental Examination:  
 Missing Teeth: 1, 6, 11, 16.  
 Defective Restoration: 2(DOL) Amalgam, 3(OL) Amalgam, 14(OL) Amalgam, 15(DOL) Amalgam.

Figure 132: Clinical Record TIU Note Window

## Problems

Click the **Problems** radio button to view a list of active problems previously entered.

To inactivate a problem:

1. Select the **problem** from the list that should be inactivated.
2. Right-click on the problem selected.
3. The **Inactive Highlighted Problem List Entry** button displays.




Figure 133: Inactivate Highlighted Problem List Entry Button

4. Select the **Inactive Highlighted Problem List Entry** button to inactivate the problem.
5. Select **Yes** and a screen displays, confirming the inactivation. The problem is removed from the list in DRM Plus and in CPRS the problem is listed as inactive.

**Note:** Use CPRS to reactivate the problem. Please see the File Data Option Screen section in the Completing the Encounter chapter of this manual for more information to add a new problem.

## Consultations

Selecting the **Consultations** radio button displays a complete list of consults. They appear with the abbreviated notation of the consult status included in the listing. Consults can be filtered by status, service or date range.

To filter consultations by status, service or date range:

1. Click the **Consultations** radio button. Consultations, if present, display in the entries area.
2. Right-click in the entries area. The filter menu displays.




Figure 134: Consultations Filter Menu

3. Select the desired filter from the menu.
4. Consults filtered by status, service or date range have a drop-down menu to filter the consults into smaller sub-views.
5. The consultation sub-view list displays results in this consult on the right side of the **Clinical Record** tab screen.

## Notes

Selecting the **Notes** radio button reveals a listing of all completed notes (the default listing). Right clicking the note window, where all notes are listed, brings up an option box.






Figure 135: All TIU Completed Notes List

Select the **TIU Filters** menu to filter the list of TIU progress notes by the listed criteria.



**Figure 136: TIU Filters Submenus**

Open the desired note. Right click the list of notes again to view the functions available for the selected note.




**Figure 137: Note Functions Menu**

The functions available coincide with the **Action** menu options, as seen on the **Notes** tab in CPRS, and works similarly. The note functions can be selected with a left click.

## Adding a New TIU Progress Note


To add a new TIU progress note for informational purposes only, without an ADA procedure code:

1. Select the **New Note** option from the **Note Functions** menu.
2. The **Set Progress Note Title** screen displays.



**Figure 138: Set Progress Note Title Screen**

3. Select the TIU progress note title from the list on the **Set Progress Note Title** screen. Note information displays on the screen.



**Figure 139: New TIU Progress Note without Any Procedures**

4. Use the tools to create a historic visit or select a scheduled visit by using the drop-down menu.
5. Enter notation directly into the note.
6. Enter the provider's electronic signature and click the **Finish** button to complete the note.


**Note:** Historical notes may be entered using this option.

## Adding a New TIU Progress Note Addendum

Creating an addendum to a previously completed TIU progress note to provide additional information, or to clarify any issues, does NOT require entering an ADA/CPT procedure code. This type of addendum can be done from DRM Plus in the **Clinical Record** tab or in the CPRS GUI. An addendum that adds an ADA/CPT procedure to a signed note requires passing information to VistA PCE/DES and requires entering an ADA/CPT code through the **Completed Care** screen entry process.

To record a note addendum without an ADA/CPT procedure code:

1. Select the **note to be appended** from the list of notes.
2. The note displays in the viewer.
3. Right-click the area where the notes are listed to view the **Note Functions** menu.
4. Select **Addendum to a Note** from the menu. Only signed notes can have an addendum.
5. Type the note directly into the note viewer. Right-click on the note viewer to import menu information or cancel the note.



**Figure 140: Import Menu Information Screen**


6. Enter the provider's **Electronic Signature** and click the **Finish** button to add the addendum.
7. A confirmation screen displays. Click the **OK** button.

**Note:** If **Addendum to a Note** is selected in error, right-click the appended note and select the **Cancel Note** submenu.



# Dental History

The **Dental History** tab displays all dental completed care information for each tooth and non-tooth entry filed in DRM Plus.



**Figure 141: Dental History Tab**

## Viewing Dental Information by Tooth

To view dental information by tooth:

1. Click a **tooth** in the tab diagram. Teeth numbered in red with an asterisk have information associated with them. When a tooth is selected, the tooth graphic is colored red in the diagram.




Figure 142: Dental History Teeth Diagram

2. Information about the selected tooth displays on the right side of the screen.
3. To de-select a tooth, click it again. The tooth turns white and the information about the tooth is removed from the right side of the screen.

## Viewing Other Dental History Information

To see dental history that is NOT listed by tooth, select the **Other** check box. The information displays on the right side of the screen.


## Viewing All Dental History Information

To see all dental history associated with this patient, select the **All** check box. The information displays on the right side of the screen.

## Viewing Dental History Information by Episode of Care

To view a patient's dental history by episode of care:

1. Choose whether to view the **Current Episode of Care** or **All Episodes of Care** by selecting the appropriate radio button.



**Figure 143: Episode of Care Radio Buttons**

2. Select the **tooth** or **teeth** to be viewed or click the **Other** or **All** check boxes.
3. The information displays on the right side of the screen.

## Episode in Date Range

To view a patient's dental history information by date:

1. Click the **Start** and **End** drop-down arrows.




**Figure 144: Episode in Date Range**

2. Use the calendar screen to choose the desired start and end dates.
3. Click the **Episode in Date Range** button. The results display on the right side of the screen.
4. Select the **tooth** or **teeth**, or the **Other** or **All** check boxes to view the desired information.

# Chart/Treatment – Treatment & Exam

There are two main sections to the **Chart/Treatment** tab: **Treatment & Exam** and **Periodontal Chart**.



**Figure 145: Chart/Treatment Tab Displaying Treatment & Exam Screen**

The **Treatment & Exam** screen has several important component views. On the upper left side of the screen are the following functions: **Diagnostic Findings**, **Treatment Plan** and **Completed Care** views, with their corresponding buttons. These tools are used to enter information on diagnostic findings, create treatment plans, view previously completed care, or enter dental procedures/diagnoses on today's encounter. The last button **Include "..."** allows the user to view information from a combination of the views on one screen.




**Figure 146: Diagnostic Findings, Treatment Plan, Completed Care and Include Buttons**

The **Seq Plan** (Sequencing), **Chart Hx** (History), **Summary**, **TMJ**, **Habits** (Parafunctional), **H&N**, **PSR**, **OHA** (Oral Health Assessment), **Occl** (Occlusion), and **Social Hx** (History) buttons are on the upper right side of the screen.

* Seq Plan	* H & N
Chart Hx	* PSR
Summary	* OHA
* TMJ	* Occl
Habits	Social Hx

**Figure 147: Treatment & Exam Specialty Buttons**

The center of the screen has a graphic display of all the teeth. Use the visual representation in combination with the **Diagnostic Findings**, **Treatment Plan** and **Completed Care** view buttons and icons to enter information about the patient. There is also a **Key** button on the right side of the display that shows which conditions the various colors, patterns and symbols represent.



**Figure 148: Completed Care Graphic Display**

The lower portion of the screen has buttons to change the view in the graphic display. Additionally, the transaction table, which shows detailed information entered in the **Diagnostic Findings**, **Treatment Plan** and the **Completed Care** view screens. There are tool buttons to the right of the transaction table, which allow the user to **Edit**, **Delete** or **Complete** the planned treatment; the **Next** button moves on to updating the encounter. Toggle between the teeth of the upper and lower arch by clicking the **Upper** and **Lower** buttons to the left of the transaction table on the bottom section of the screen. The **Full** button allows the user to view both arches.

	Visit Date	Stat	Category	Tooth/Quad	Surf/Root	Diag	Code	Description	Prov	
Upper	06/08/2017	C	Restored	18	DOL	K08.530	D2160	AMALGAM THREE SURFACES PERMA	DEN	Edit
	05/31/2017	C	Restored	14	OL	K08.530	D2150	AMALGAM TWO SURFACES PERMANE	ADMIN	
	05/31/2017	C	Restored	15	DOL	K08.530	D2160	AMALGAM THREE SURFACES PERMA	ADMIN	
Lower	05/31/2017	C	Restored	30	OL	K02.52	D2150	AMALGAM TWO SURFACES PERMANE	RES2	Delete
	05/31/2017	C	Restored	31	DOL	K02.52	D2160	AMALGAM THREE SURFACES PERMA	RES2	
	05/31/2017	C	Diagnost			K02.53	D0210	INTRADURAL FULL IMAGE SERIES	ADMIN	
Full	05/26/2017	C	Diagnost			K02.61	D0330	DENTAL PANORAMIC IMAGE	ADMIN	Complete
	05/26/2017	C	Diagnost			K02.62	D0191	BRIEF ASSESSMENT	ADMIN	
	04/28/2017	C	Diagnost			K02.52	D0140	LIMIT ORAL EVAL PROBLM FOCUS	ADMIN	
	04/28/2017	C	Diagnost			K02.52	D0330	DENTAL PANORAMIC IMAGE	ADMIN	Next =>
	03/22/2017	C	ImplPost	6		K08.434	D6010	ENDOSTEAL IMPLANT BODY PLACE	ADMIN	

**Figure 149: Completed Care Transaction Table**

## Diagnostic Findings

1. Click the **Diagnostic Findings** view button on the left side of the screen.
2. Select the desired finding from the icons to the right of the **Diagnostic Findings** view button.




Figure 150: Diagnostic Findings View Button with Icons

3. Click the **tooth/area of the tooth** in the graphic display. Use the **Upper** and **Lower** buttons on the left side of the text display to view the arch and the previous diagnostic findings on the upper and lower arches.




Figure 151: Diagnostic Findings Graphic Chart and Transaction Table

4. The finding appears in both the text and graphic display. Click the **Key** button in the upper right corner of the graphic display to see how various findings are shown in the graphic.
5. Use the **Clear** icon to remove any finding entered during today's encounter only. Click the **Clear** icon and then click the desired finding on the graphic to remove the finding from both the graphic and the transaction table.

**Note:** The **Stats** column in the transaction table displays an 'F' when the transaction is a finding.

## Editing Diagnostic Finding Descriptions

**Diagnostic Finding** descriptions that have been entered, but for which no TIU progress note has yet been filed, can be edited.

To edit a diagnostic finding description:

1. Select the **finding** by highlighting it in the transaction table.
2. Click the **Edit** button to the right of the transaction table.
3. The **Edit Transaction** screen displays.




Figure 152: Edit Transaction Screen

4. Enter the new description in the text box. Note that only the description can be edited. The other information, such as **Visit Date**, **Tooth/Quad** or **Category** cannot be edited with this button.
5. Click the **OK** button. The edited information displays in the transaction table.

## Deleting a Diagnostic Finding

**Diagnostic Findings** that have been entered, but for which no TIU progress note has yet been filed, can be deleted.

To delete a diagnostic finding:

1. Select the **finding** by highlighting it in the transaction table.
2. Click the **Delete** button on the right side of the transaction table.
3. The finding is deleted.

**Note:** If the user attempts to delete a finding that has already been filed with an old encounter, the item is removed from the graphic display, but remains in the transaction table with a line through it. The DRM Plus Administrator may delete the finding from the graphic and the transaction table, unless the transaction has previously been deleted by a DRM Plus non-Administrator.

**Note:** Clicking a transaction table **column heading** sorts the table. Generally, in ascending order depending on the current view. Clicking the **column heading** a second time returns the table to the original descending view. This functionality works the same for all three **Treatment & Exam** transaction table views.

## Treatment Plan

### Entering a Treatment Plan

There are multiple ways to enter a planned procedure code for a patient: by adding the code directly utilizing the **Add** button with text box, or by selecting the icon that corresponds to the planned treatment and choosing the tooth from the graphic display. Use the **ADA Codes** icon, **CPT Codes** icon, **Quick Codes** icon or **Personal Speed Code** icon as additional ways to enter planned treatment for the patient.

Rules for entering a procedure code for a planned item in the **Treatment Plan** view include:

1. Always use a **standard** icon, first four columns and **P&C**, if one is available first.
2. When no standard icon is available, use the **ADA Codes**, **CPT Codes**, **Quick Codes** or a **Speed Code** icon.
3. The **Add** button with text box may be used interchangeably with rule 2.
4. Always enter transactions in the same order that they are performed on the patient.



Figure 153: Treatment Plan View Button with Icons

To enter a planned treatment using the **Add** button and text box:

1. Click the **Treatment Plan** view button. Notice that the **Add** button and text box are active.
2. Type a word or procedure code into the **Add** text box and click the **Add** button.
3. The **Code Details** screen displays if the procedure code is tooth or quadrant related.




Figure 154: Code Details Screen

Enter the **Tooth** number and, if applicable, the **Area** or **Surface** modifiers and click the **OK** button. The graphic and the displays adjust to reflect the addition.

**Note:** When surface modifiers are required for a procedure code, they must be entered using uppercase letters (**M, O, D, F, B, L** and **I**). When root modifiers are required for a procedure code, they must be entered using lowercase letters (**r, b, l, d** and **m**).



To enter a planned treatment using the standard **Treatment Plan** icons:

1. Click the **Treatment Plan** view button.
2. Click the icon to the right that corresponds to the desired planned treatment.
3. Click the appropriate **tooth, area** and/or **surface**.
4. As in diagnostic findings, toggle between the upper and lower arch by clicking the **Upper** or **Lower** button on the left side of the display.
5. The graphic display and the transaction table display the new addition.

To enter a planned treatment using the **ADA Codes** icon or the **CPT Codes** icon:

1. Click the **ADA Codes** icon.




Figure 155: ADA, CPT and Quick Codes Icons

2. The **Procedure Code Selection** screen displays for code selection.

The image shows a software interface for selecting procedure codes. It has a blue background. At the top, it says 'Please select the code:'. Below this is a 'Filter:' text box. Underneath is a 'Codes:' section containing a table with two columns: 'Code' and 'Description'. The table lists various dental codes and their descriptions. Below the table is a 'Full Description:' text box. At the bottom, there are four buttons: 'Category', 'Multi-Add', 'Add', and 'Exit'.

Code	Description
D0100	DIAGNOSTIC
D1100	PREVENTIVE
D2100	RESTORATIVE
D3000	ENDODONTICS
D4200	PERIODONTICS
D5100	REMOVABLE
D5900	MAXILLOFACIAL PROSTH
D6000	IMPLANT SERVICES
D6200	FIXED PROSTHODONTICS
D7100	SURGERY
D8000	ORTHODONTICS
D9000	ADJUNCTIVE
D0100	DIAGNOSTIC

Figure 156: Procedure Code Selection Screen

3. Type the procedure code into the filter or use the scroll bar to search for a code.
4. A description of the highlighted procedure code displays in the **Full Description** text box.
5. Click the **Category** button to return to the list of categories at the top of the scroll sheet. Click the **Multi-Add** button to add multiple procedure codes on other teeth.
6. When the desired code is highlighted, click the **Add** button.
7. If the procedure code needs to be attached to a specific tooth and surface modifier, the **Code Details** screen displays.

	Code	Description	Area	Tooth	Surfaces
1. <input checked="" type="checkbox"/>	D2140	AMALGAM ONE SURFACE PERM			

Ok Cancel

**Figure 157: Code Details Screen**

8. Fill in the fields with the requested information and click the **OK** button.
9. The **Diagnosis Code** screen displays if the parameter is activated.
10. Select the planned treatment procedure on the left side of the screen (default is all selected) and the diagnosis code(s) on the right side. If the correct diagnosis code does NOT appear, use the **Additional Diagnosis Code Search** window to find a different diagnosis code.
11. Click the **OK** button. The information displays in the transaction table on the graphic chart.
12. To undo any graphical entry on today's encounter, use the **Clear** icon as described in the **Diagnostic Findings** section of this chapter.

Entering a planned treatment using the **Quick Codes** icon is like adding a planned treatment with the **ADA Code** icon. To enter a planned treatment using the personal **Speed Code** icons, please see the **Perio Buttons** section later in this chapter.

## Editing a Treatment Plan Description

1. Select the desired planned entry in the transaction table.
2. Click the **Edit** button. The **Edit Transaction** screen displays.

Description:

AMALGAM ONE SURFACE PERMANEN

Ok Cancel

**Figure 158: Edit Transaction Screen**

3. Enter the new description in the text box.
4. Click the **OK** button.
5. The description is changed in the transaction table. Note that the **Code**, **Category** and other information cannot be changed with the **Edit** button.

## Deleting a Treatment Plan

1. Select the desired planned entry in the transaction table.
2. Click the **Delete** button.
3. The planned treatment is removed from the graphic and the transaction table.

Every DRM Plus user can delete any planned item, regardless of who entered it. The planned entry is removed from the graphic and transaction table in the **Treatment Plan** and **Sequencing** screens.

## Completing a Treatment Plan

1. Select the planned procedure to be completed from the transaction table.
2. Click the **Complete** view button.
3. The **Diagnosis Code** screen displays.

**Completed Care Code: D2140, in use by: DRMPROVIDER,DENTIST, for: DRMPATIENT,ONE**

**Includes all adhesives (amalgam bonding agents), liners, bases. If pin(s) used also include code D2951.**

Select your Completed Procedure(s) on the left and their related Diagnosis Code(s) on the right.  
If the suggested Diagnosis Code(s) are not applicable or are blank, use the Code Search to find the appropriate Diagnosis Code(s).

Note: Previously filed Diagnosis Code(s) are listed in the lower right checklist box.

☒ Tooth# / Surface(s): 3 / L

☐ K02.52-DENTAL CARIES ON PIT AND FISSURE SURFACE PENETRATING INTO DENTIN  
☐ K02.62-DENTAL CARIES ON SMOOTH SURFACE PENETRATING INTO DENTIN  
☐ K02.7-DENTAL ROOT CARIES  
☐ K03.0-EXCESSIVE ATTRITION OF TEETH  
☐ K03.1-ABRASION OF TEETH  
☐ K03.81-CRACKED TOOTH  
☐ K08.51-OPEN RESTORATION MARGINS OF TOOTH  
☐ K08.52-UNREPAIRABLE OVERHANGING OF DENTAL RESTORATIVE MATERIALS  
☐ K08.530-FRACTURED DENTAL RESTORATIVE MATERIAL WITHOUT LOSS OF MATERIAL  
☐ K08.531-FRACTURED DENTAL RESTORATIVE MATERIAL WITH LOSS OF MATERIAL  
☐ K08.54-CONTOUR OF EXISTING RESTORATION OF TOOTH BIOLOGICALLY INCOMPATIBLE WITH ORAL HEALTH  
☐ K08.56-POOR AESTHETIC OF EXISTING RESTORATION OF TOOTH  
☐ K02.5XXA-FRACTURE OF TOOTH (TRAUMATIC), INITIAL ENCOUNTER FOR CLOSED FRACTURE  
☐ K02.53-DENTAL CARIES ON PIT AND FISSURE SURFACE PENETRATING INTO PULP  
☐ K02.63-DENTAL CARIES ON SMOOTH SURFACE PENETRATING INTO PULP

Additional Diagnosis Code Search  
 <Type 3 or more characters and press enter.>

☐ K08.530-FRACTURED DENTAL RESTORATIVE MATERIAL WITHOUT LOSS OF MATERIAL  
☐ K02.52-DENTAL CARIES ON PIT AND FISSURE SURFACE PENETRATING INTO DENTIN  
☐ K02.53-DENTAL CARIES ON PIT AND FISSURE SURFACE PENETRATING INTO PULP  
☐ K02.62-DENTAL CARIES ON SMOOTH SURFACE PENETRATING INTO DENTIN  
☐ K02.61-DENTAL CARIES ON SMOOTH SURFACE LIMITED TO ENAMEL

Primary	Selected Procedure(s)	Selected Diagnosis(es)

General Coding Standards ☒ File in PCE ☒ File in Dental History

OK Cancel

**Figure 159: Diagnosis Code Screen**

4. Select the planned treatment procedure on the left side of the screen (default is all selected) and the diagnosis code(s) are mapped on the right side of the screen. If the correct diagnosis code is **NOT** listed, use the **Additional Diagnosis Code Search** to find a different diagnosis code.
5. The diagnosis codes in the lower section comprise the **Diagnosis Codes QuickList** which have been filed with every encounter for this patient's chart during the past 24 months. This should help the provider decrease the entry time for searching and locating past diagnoses that may still be relevant for the patient's visit.
6. Notice the **File in PCE** check box. If the completed treatment is **NOT** to be filed in **PCE**, uncheck this box.
7. Click the **OK** button.
8. The planned procedure is removed from the transaction table and is now part of **Completed Care** transactions for the patient.

**Note:** When no diagnosis code(s) are mapped to (or listed with) a procedure code on the **Diagnosis Code** screen the **Additional Diagnosis Code Search** must be used to find the appropriate diagnosis code(s) to be filed with the selected procedure code.

**Note:** The **PCE** check box option should only be used by advanced users. Changing this option effects, the data being sent to VistA PCE. Do NOT change this check box unless the user does NOT want the data to be sent to VistA PCE.

**Note:** Planned procedure codes can be designated with or without a related diagnosis. The parameter for dental users to acquire this functionality is **Tools** menu → **User Options** submenu → **Treatment System** tab → select the top check box.

## Completed Care

Click the **Completed Care** view button to see all treatments that have been previously completed in the VA for the patient or entering any new completed treatment for today's visit.



**Figure 160: Completed Care View Button Active with Icons**

The rules for entering a procedure code for a completed treatment in the **Completed Care** view include:

1. Always use a **standard** icon, first four columns and P&C, if one is available first.
2. When no standard icon is available use **ADA, CPT, Quick Codes** or a **Speed Code** icon.
3. The **Add** button and corresponding text box may be used interchangeably with rule 2.
4. Always enter transactions in the same order they are performed on a patient.

## **Entering Completed Care**

There are several ways to enter completed treatment. Planned treatments, which are completed (see Completing a Treatment Plan), display in the **Completed Care** transaction table. Completed treatment can also be entered manually:

1. Click the **Completed Care** view button.
2. Select the desired associated **standard** icon.
3. Choose the appropriate tooth/area on the graphic display.
4. Complete the **Diagnosis Code** screen.
5. The entry displays in the graphical/transaction tables.
6. To undo any graphical entry on today's encounter, use the **Clear** icon as described in the **Diagnostic Findings** section of this manual.

Completed treatment can also be entered through the **Add** button and text box, the **ADA**, **CPT** and **Quick Codes** icons, as well as **Speed Code** icons. Please see the **Treatment Plan** section of this manual for further information on these functions.

## Editing Completed Care Description

**Completed Care** description can only be edited if it has NOT yet been made a part of the TIU progress note. To edit a completed care entry:

1. Choose a completed care entry from the transaction table.
2. Click the **Edit** button.
3. The **Edit Transaction** screen displays.




Figure 161: Edit Transaction Screen

4. Type the new description into the text box.
5. Click the **OK** button. Note that only the description can be edited.

## Deleting a Completed Care

1. Click the desired completed entry in the transaction table entered this session.
2. Click the **Delete** button.
3. The completed treatment is removed from the graphic and transaction tables.
4. To finalize the deletion, the encounter with the patient must be completed and filed.


DRM Plus users are NOT allowed to delete any previously filed completed treatment. This can be performed only by DRM Plus Administrators.

## Include “Completed”/Include “Findings and Completed”/Include “Findings”

This button can be used in conjunction with the other buttons to include more than one type of information on the display. When the **Diagnostic Findings** view button is active, and the associated information is displayed in the graphic and transaction tables, clicking this button adds the **Completed Care** information to the display. If the **Treatment Plan** view button is active, and the associated information is displayed, clicking this button adds the **Diagnostic Findings** and the **Completed Care** information to the graphic and transaction tables. Finally, when the **Completed Care** view button is active, and the associated information is displayed, clicking this button adds the **Diagnostic Findings**.

## Perio Buttons Icon

The icon table includes the **Perio Buttons** icon. Click this **Perio Buttons** icon to see the second set of speed codes associated with perio mode, if entered by the user. Please see the Speed Codes section in the Using the DRM Plus Drop-Down Menu chapter of this manual for further information.




**Figure 162: Perio Buttons and Speed Code Icons**

This option is only available while **Treatment Plan** or **Completed Care** are active, but NOT with **Diagnostic Findings**.

Click the desired icon to add it to the planned treatment. If these buttons are clicked when **Completed Care** screen is active, the **Diagnosis Code** screen displays, which allows the selection of diagnosis codes that are mapped to that procedure code to be entered.

To change the location of the speed code icon:

1. Click one of the **blank** icons where the speed code icon is to be moved.
2. The **Configure Button** screen displays.



**Figure 163: Configure Button Screen**

3. Use the drop-down menu to assign the speed code.
4. Use the **Perio Mode** check box to link the speed code icon with the **Perio Buttons** icon.
5. Click the **OK** button. The old speed code icon location is cleared, and the speed code icon is now in the new assigned icon location.

## Seq Plan/Sequencing Button

Use sequencing in combination with the **Treatment Plan** screen to organize when to perform specific planned treatment. There are nine buttons on the **Tx Planning/Sequencing** screen. The nine buttons are **Add Phase**, **Add Sub-phase**, **Add Non-VA Care Phase**, **Copy Non-VA Care Phase to Clipboard**, **Copy All to Clipboard**, **Print**, **Save & Exit**, **Save** and **Cancel**. The **Tx Planning/Sequencing** screen has the nine buttons positioned or grouped on the right side of the screen.

The screenshot shows the 'Tx Planning/Sequencing' window. It features a tree view on the left under 'Planned Treatment' with the following structure:

- Unsequenced Treatment
  - Non-VA Care
    - D7950 - MANDIBLE GRAFT
  - Phase 1
  - Phase 2
    - D2150 - AMALGAM TWO SURFACES PERMANE [14/DOL]
    - D2160 - AMALGAM THREE SURFACES PERMA [15/DOL]

To the right of the tree is a table with three columns: 'Next Appt.', 'Complete', and 'Delete'. Each row in the tree has a corresponding row of checkboxes in this table.

Planned Treatment	Next Appt.	Complete	Delete
Unsequenced Treatment	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Non-VA Care	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
D7950 - MANDIBLE GRAFT	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Phase 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Phase 2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
D2150 - AMALGAM TWO SURFACES PERMANE [14/DOL]	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
D2160 - AMALGAM THREE SURFACES PERMA [15/DOL]	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

On the right side of the window, there are nine buttons: 'Add Phase', 'Add Sub-phase', 'Add Non-VA Care Phase', 'Copy Non-VA Care Phase to Clipboard', 'Copy All to Clipboard', 'Print', 'Save & Exit', 'Save', and 'Cancel'. At the bottom, there is a text area for 'Additional Dental Treatment Plan Notes' with a timestamp and a message: '[09:28 AM 05/25/2017 ADMIN] Adding an additional planned note. Select the F5 key first and then enter text data.'

**Figure 164: Tx Planning/Sequencing Screen**

All planned transactions sequenced into a **Phase 1**, **Phase 2**, etc. and with all sub-phases associated with numbered phases has full functionality of **Next Appt.**, **Complete** and **Delete** check boxes.

The **Add Non-VA Care Phase** button allows only one **Non-VA Care** phase to be added. The **Non-VA Care** phase will be listed below the **Unsequenced Treatment** and above **Phase 1**. Users may add an unlimited number of sub-phases in a **Non-VA Care** phase. When trying to add a second **Non-VA Care** phase a warning message appears stating that only one **Non-VA Care** phase may be added.

The **Next Appt.** and **Complete** check boxes for the **Non-VA Care** phase are inactive; if selected, an informational screen appears stating these check boxes may NOT be used in the **Non-VA Care** phase. The **Delete** check box remains active for the **Non-VA Care** phase and all **Non-VA Care** sub-phases.

The DRM Plus user should delete a planned **Non-VA Care** phase, sub-phases and transactions after receiving verification from the Non-VA Care provider that the procedure(s) has/have been completed.

The two buttons of **Copy Non-VA Care Phase to Clipboard** and the **Copy All to Clipboard** when selected will allow the user to copy the planned transactions listed in the **Tx Planning/Sequencing** screen. The copied data includes filed and unfiled planned treatment from only the **Tx Planning/Sequencing** screen to paste on any word document, text document or any application window if allowed.

The end-user may use the <F5> key now in the **Additional Dental Treatment Plan Notes** to import a date/time stamp identifier before or after the text entry if desired.

## Plan a Treatment Sequence

1. Click the **Seq Plan** (Sequencing) button.
2. The **Tx Planning/Sequencing** screen displays.
3. Information from the **Treatment Plan** transaction table is shown on the screen. Use the **Add Phase**, **Add Sub-phase** or **Add Non-VA Care Phase** buttons to add a new phase and/or sub-phase. Highlight the desired phase listed in the screen to add the sub-phase under this phase.
4. Change the sequence of the planned treatments by dragging and dropping them into the correct phase.
5. If the planned treatment is to be completed at the next appointment, click the corresponding check box. **Non-VA Care** phase is NOT allowed a next appointment.
6. Add Additional Dental Treatment Plan Notes in the text box.
7. Click the **Save & Exit** button if only planned items have been added or sequenced for this patient. This option requires that no new data be entered as completed transactions, Perio, H&N or any other modal at the same time for the option to work.
8. Click the **Save** button to save the progress in sequencing and keep working on this encounter.

**Note:** Always enter planned treatment in the same order they are performed on a patient.

**Note:** The **Save & Exit** button from the sequencing screen files any changes and minimize DRM Plus. Any new planned entries added have the same **Visit** date as the latest TIU progress note filed on this patient. The most recent dental encounter must have an **Active** status for this feature to work.

## Complete a Planned Treatment in the Sequencing Screen

1. Click the **Seq Plan** (Sequencing) button.
2. The **Tx Planning/Sequencing** screen displays.
3. Choose the planned treatment that is to be completed by checking the corresponding check box in the **Complete** column. **Non-VA Care** phase is NOT allowed to be completed in the **Tx Planning/Sequencing** screen.
4. Click the **Save** button.
5. The **Diagnosis Code** screen displays. Please see the Completing a Treatment Plan portion in the Treatment Plan section of this chapter for further information.

## Deleting a Planned Treatment in the Sequencing Screen


1. Click the **Seq Plan** (Sequencing) button.
2. The **Tx Planning/Sequencing** screen displays.
3. Choose the planned treatment that is to be deleted by checking the corresponding check box in the **Delete** column. The DRM Plus user should delete a planned **Non-VA Care** phase, sub-phases and transactions after receiving verification from the Non-VA Care provider that the procedure(s) has/have been completed.



4. Click the **Save** button.
5. The planned treatment is deleted from the sequencing screen, the transaction table and the graphical chart on the **Treatment Plan** screen.

## **Chart Hx (History) Button**

Click the **Chart Hx** button to see a completed care chart of the patient's dental history. The transaction table includes the text details of the **Visit Date**, **Stat**, **Category**, **Tooth/Quadrant**, **Surface/Root** modifiers, **Codes**, **CPT Description** and **Provider** initials.



**Figure 165: Tx History Screen**

Use the **Appointments** drop-down menu to see the patient's history by different appointment dates. View tooth notes on the patient's file by clicking the **Notes** button tied to the note's appointment date.

**Note:** To display any past tooth-specific note, click the **Chart Hx** button. Continue by clicking the drop-down arrow of the **Appointment** field, in the top left corner of the **Tx History** screen and selecting the appropriate date the tooth-specific note was entered. Once the date is selected, the tooth numbers in the graph displays in yellow. Once the desired tooth is found, click the **Notes** button. Then use the drop-down arrow to select the tooth-specific note of interest and view its contents.

## Summary Button

Click the **Summary** button to view a summary of a patient's chart. Periodontal information displays in the summary as well. Use the tools on the screen to view the information by quadrant or by tooth.




Figure 166: Tooth/Quadrant Summary Screen

**Note:** On the **Treatment & Exam** screen, with the following views: **Diagnostic Findings**, **Treatment Plan** or **Completed Care**, clicking the **Summary** button displays the history of the selected primary view. The upper half of the window shows the summary of the primary view that is active when the **Summary** button is selected. Only the activated primary view, with all its entries, displays in the restorative top window. The lower half of the screen displays the periodontal summary, which includes the latest exam in periodontal history with filed **Diagnostic Findings** and **Completed Care**; however, NOT including surfaces or roots.

## H&N Button

Use the functions in the **H&N** (Head and Neck) button to enter and view diagnostic information on the patient's head and neck.

Date	Provider	Type	Description
04/28/2017	ADMIN	INITIAL	Cancerous condition

**Figure 167: Head and Neck Findings Screen**

To enter an **H&N** finding:

1. Select the graphic shape that best represents the finding by using the **Shape** drop-down menu.
2. Select the contrasting color for the finding by clicking the **Color** box. A visual list of possible colors displays.
3. Select the size of the graphic by using the **Size** drop-down menu.
4. The date box defaults to **Today**, which is required for new data entries. After the first filed TIU progress note entry of **H&N**, the date box defaults to **All**. The user may click the drop-down arrow and highlight a previous exam to view the entries on a previous date.
5. Click the graphic to show where the lesion is located on the patient.
6. The **H&N Detail** screen displays.

The screenshot shows a software window titled "H & N Detail". Inside, there is a section labeled "Initial Finding" containing a table with three columns: "Date", "Provider", and "Description". The table has one row with the values "06/09/2017", "DEN", and "Amalgam Tattod". To the right of the table is a button labeled "Common Findings". Below the table section is a large, empty rectangular area labeled "Progress". At the bottom right of the window, there is a checkbox labeled "Resolved" which is currently unchecked, and two buttons labeled "Save" and "Cancel" stacked vertically.

Date	Provider	Description
06/09/2017	DEN	Amalgam Tattod

Common Findings

Progress

☐ Resolved

Save

Cancel

**Figure 168: H&N Detail Screen**

7. Enter the description of the lesion in the **Description** column or click the **Common Findings** button to see a list of commonly appearing lesions and add the description to the finding.
8. Click the **Save** button. The finding displays on the **H&N Findings** transaction table.

To add new details of an **H&N** finding:

1. Click the desired finding in the **H&N Findings** transaction table and click the **Details** button.
2. The **H&N Detail** screen displays with progress information to be entered.

Date	Provider	Description
04/28/2017	ADMIN	Cancerous condition

Common Findings

Date	Provider	Description
06/09/2017	DEN	<<Enter description>>

☐ Resolved

Save

Cancel

**Figure 169: H&N Detail Screen with Progress Description**

3. Enter information in the **Description** column of the **Progress** window.
4. Click the **Save** button to save the entered information and return to the **H&N Findings** screen.


To delete a lesion:

1. Highlight the finding entry on the **H&N Findings** transaction table.
2. Click the **Delete** button.
3. A screen confirming that the entry is to be deleted displays. Click the **OK** button. Note that this screen only displays if the finding entry is made during this session.
4. The entry and the mark on the **H&N Findings** graphic are removed if entered this session.

**Note:** Clicking the **Delete** button deletes the highlighted finding entry(s) during the same session the **H&N Findings** transaction was entered. Deleting any filed **H&N Findings** transaction by any DRM Plus end-user results in a line through the entry, and it remains in the transaction table.

## **PSR Button**

Click the **PSR** button to view the **Periodontal Screening/Recording (PSR)** screen. Use this screen to view previously entered PSR information if present, or to enter new PSR information.



The screenshot shows a software window titled "Periodontal Screening/Recording...". Inside the window, there is a label "Exam Date:" followed by a dropdown menu showing "04/28/2017" and a "New" button. Below this is a large rectangular area containing six text boxes arranged in two rows of three. The top row is labeled "Sextant #1:", "Sextant #2:", and "Sextant #3:". The bottom row is labeled "Sextant #6:", "Sextant #5:", and "Sextant #4:". Each text box contains the number "1". At the bottom of the window are three buttons: "Definitions", "Ok", and "Cancel".

**Figure 170: Periodontal Screening/Recording (PSR) Screen**

Use the **Exam Date** drop-down menu to view past periodontal screening information.

To enter new information:

1. Click the **New** button. All the **Sextant** text boxes will default to a 0-dental value.
2. Enter the desired national dental value in each **Sextant** text box. Entering an “\*” requires a number added with the symbol to be saved (i.e., “3\*”).
3. Click the **Definitions** button to view PSR definitions of national dental values.
4. Click the **OK** button to complete.

**Note:** The **PSR** modal allows two providers to enter a PSR exam on the same day; however, it only displays the last PSR exam that was entered on that day. The first entered PSR exam is only viewable in the TIU progress note of the provider that filed the encounter. The second provider’s filed note has a different header, which includes the word ‘modified’ in the PSR exam of the TIU progress note.

## OHA (Oral Health Assessment) Button

The specialty **Plaque** button was combined with **Xerostomia**, **Caries Risk** and **Oral Hygiene** to create the **OHA** button. This also applies to the **Plaque** button on the **Periodontal Chart** screen.

To enter new data in the **Oral Health Assessment (OHA)** screen, click the **New** button. Today's date is imported into the **Date** field on the screen.

Date	Provider	NFT	PI	X	CR	OH
06/09/2017	DEN		1	0	2	2
05/31/2017	ADMIN		1	0	2	2
04/28/2017	ADMIN		1	0	1	1

**Figure 171: Oral Health Assessment (OHA) Screen**

The **NFT** check box, **Patient has no remaining functional teeth, roots or implants**, may be selected in the event the patient meets these criteria and no findings can be entered in the **Diagnostic Findings** chart.

Checking this box automatically completes the **Diagnostic Findings** element and the **Periodontal Assessment** element when filing any exam/consult code during a dental encounter. It also automatically selects the **0-Edentulous** radio button in the **Caries Risk** section.

The radio buttons default to **4 – Not Recorded** in all four fields. This selection does NOT import any clinical finding into the TIU progress note, nor does it display in the transaction table of the **OHA** screen. The provider has the option of selecting the appropriate radio button (0-3) for each of **Plaque Index**, **Xerostomia**, **Caries Risk** and **Oral Hygiene**, or simply leaving the default setting.

The entry of date, provider's initials and each field value entered is captured in the transaction table at the bottom of the screen. The provider must enter at least one value (between 0-3) in one of the four fields to save and file an oral health assessment.

The **Definitions** button has the American Dental Association definitions for field values when entering **Plaque Index**, **Xerostomia Risk** and **Caries Risk**. The rest of the **Xerostomia** and **Caries Risk** definitions may be viewed using the scroll bar on the right side of the screen.

**DRM Plus - Oral Health Assessment**

Date: 06/09/2017 [New] [Definitions]

**Clinical Findings**

☐ Patient has no remaining functional teeth, roots or implants.

**Plaque Index**

☐ 0 - None  
☒ 1 - Slight  
☐ 2 - Moderate  
☐ 3 - Heavy  
☐ 4 - Not Recorded

**Xerostomia**

☒ 0 - None  
☐ 1 - Slight  
☐ 2 - Moderate  
☐ 3 - Significant  
☐ 4 - Not Recorded

**Caries Risk**

☐ 0 - Edentulous  
☐ 1 - Low  
☒ 2 - Moderate  
☐ 3 - High  
☐ 4 - Not Recorded

**Oral Hygiene**

☐ 0 - Excellent  
☐ 1 - Good  
☒ 2 - Fair  
☐ 3 - Poor  
☐ 4 - Not Recorded

Date	Provider	NFT	PI	X	CR	OH
06/09/2017	DEN		1	0	2	2
05/31/2017	ADMIN		1	0	2	2
04/28/2017	ADMIN		1	0	1	1

[OK] [Cancel]

**Definitions Panel**

**PLAQUE INDEX**

0 - None:  
No Plaque in the gingival area.

1 - Slight:  
A film of plaque adhering to the free gingival margin and adjacent area of the tooth. The plaque may be recognized only by running a probe across the tooth surface.

2 - Moderate:  
Moderate accumulation of soft deposits within the gingival pocket and on the gingival margin and/or adjacent tooth surface that can be seen by the naked eye.

3 - Heavy:  
Abundance of soft matter within the gingival pocket and/or on the gingival margin and adjacent tooth surface.

**XEROSTOMIA**

0 - None:  
None.

1 - Slight:  
Moist mucosa with limited flow and no pooling of saliva. Patient may report subjective dryness and/or stale breath.

2 - Moderate:  
In addition to slight symptoms, initially moist mucosa but dries during oral examination. Saliva may be thick and ropery. Patient may reports difficulty with dry foods, altered taste

[VA Dental Definitions](#)

**Figure 172: OHA Definitions Displaying Plaque Index/Xerostomia**

The Oral Health Assessment (OHA) **Definitions Panel** is displayed in the previous screen. The **Definitions Panel** may be automatically expanded due to the parameter selection that is defaulted when first loaded. The user may change this parameter by navigating to **Tools** menu → **User Options** submenu → **Exam Settings** tab → uncheck the **Requirements**.

At the bottom of the OHA definitions panel the **VA Dental Definitions** hyperlink takes the provider to the VA dental website to view current OHA definitions. This hyperlink should be changed by the DRM Plus Administrator when instructed to do so by the VA Dental Informatics and Analytics Director.

**Note:** If NOT immediately directed to the website when selecting the hyperlink, please notify the Help Desk.

**Note:** The **Plaque Index** definitions have been reprogrammed and only allow whole number entries. Most **Plaque Index** values filed before the loading of the exam template patch retain the decimal value, if entered with one, located in the **PI** column of the **OHA** transaction table.



## **TMJ Button**

The **TMJ** modal functions like the **OHA** modal when entering a new exam. Click the **New** button and today's date imports into the **Date** field.

At least one entry in either the **History** or **Clinical Findings** sections from this **TMJ** screen requires data to be selected in order to save. Selecting the second **History** radio button, **Patient reports symptoms associated with TMJ's**, allows multiple check box selections and at least one is required for this option to save.

The text windows below the **Other** check box in **History** or **Clinical Findings** only opens if the check box has been selected and each requires a text entry. The **Other** text boxes allow an unlimited text field.

DRM Plus - TMJ

Date: 05/31/2017 [New]

**History**

☒ Patient reports no symptoms associated with TMJ.

☐ Patient reports symptoms associated with TMJ's:

☐ History of trauma ☐ Popping/Clicking ☐ Pain upon opening

☐ Spontaneous pain ☐ Crepitus ☐ Pain upon chewing

☐ Limited opening

☐ Other:

**Clinical Findings**

0 (mm) Max Incisal Opening None Popping/Clicking

0 (mm) Left Lateral None Crepitus

0 (mm) Right Lateral - Pain to manipulation

☐ Other: - Deviation upon opening

[OK] [Cancel]


**Figure 173: TMJ Screen**

The **Clinical Findings** section has three numerical fields to enter a millimeter value and four drop-down menu options in selecting popping/clicking, crepitus, pain to manipulation, and deviation upon opening. The **Other** check box allows an unlimited text field for additional text information if selected.

The minimum requirement to enter a new **TMJ** finding is the selection of only one historical or clinical finding from this screen.

## **Occl (Occlusion) Button**


The **Occl** (Occlusion) button functions are different from the other modals when entering a new occlusion finding. Click the **New** button and today's date imports into the **Date** field. When there is previous filed data present then all that filed data imports into the new exam. The user needs to add/delete any new occlusion findings and click the **OK** button to save.



**Figure 174: Occlusion Screen**

The **Clinical Findings** drop-down menu option **Mandibular relationship\*** is the only required (\*) field on this screen. The six other drop-down menu options and the two numerical box selections are optional entries.

The **Definitions Panel** displays the **Angle's Classification** definitions. These angle's classifications are for the selections displayed in the left bottom four drop-down menus. The **Definitions Panel** maybe reduced to display only the **Occlusion** screen by selecting the **Definitions** button.



**Figure 175: Occlusion Angle's Classification Definitions**

The **Definitions Panel** may be expanded due to the parameter selection that is defaulted when first loaded. The user may change this parameter by going to the **Tools** menu → **User Options** submenu → **Exam Settings** tab → uncheck the **Requirements**.

At the bottom of the occlusion **Definitions Panel**, the **VA Dental Definitions** hyperlink takes the provider to the VA dental website to view current occlusion definitions. This hyperlink should be changed by the DRM Plus Administrator when instructed to do so by the VA Dental Informatics and Analytics Director.

**Note:** If NOT immediately directed to the website when selecting the hyperlink, please notify the Help Desk.

## **Habits (Parafunctional) Button**

The **Habits** (Parafunctional) modal functions like the **OHA** modal when entering a new exam. Click the **New** button and today's date imports into the **Date** field.

The **History** and **Clinical Findings** fields each have two radio buttons for selection. When the second radio button is selected in each field, multiple options become active for selection. The **Other** check box allows an unlimited text field for additional text information.

The minimum requirement to enter a new parafunctional habit finding is to select only one historical or clinical finding from the **Parafunctional Habits** screen.

DRM Plus - Parafunctional Habits

Date: 06/12/2017 [New]

History

☐ Patient reports no known parafunctional habits.

☐ Patient reports the following habits:

☐ Bruxing ☐ History of eating disorder(s)

☐ Clenching

☐ Other:

Clinical Findings

☐ No evidence of parafunctional habits.

☐ Parafunctional habits evidenced by:

☐ Attrition ☐ Erosion

☐ Abrasion ☐ Hypertrophy of masticatory muscles

☐ Other:

[OK] [Cancel]

**Figure 176: Parafunctional Habits Screen**

## Social Hx (Social History) Button

The **Social History** modal functions like the **OHA** modal when entering a new exam. Click the **New** button and today's date imports into the **Date** field.

The minimum requirement is the selection of one of the two **History** radio buttons. When selecting the second radio button then at least one check box option is required to save the new historical data.

Any combination of check boxes maybe selected for **Present/Past**. The tobacco and alcohol drop-down options are per day, per week, per month, and per year except for cigarettes. **Cigarettes** have only the drop-down options of pack year history, per day, and per week. The text box with the **Drug Abuse** selection is optional when one of the check boxes is selected. The bottom text box is optional and allows an unlimited text field of formation if selected to enter data about eating disorders, dietary concerns, piercings, etc.

	Present	Past	
Tobacco			
Cigarettes	<input type="checkbox"/>	<input type="checkbox"/>	0 pack year history 0 years
Pipe/Cigar	<input type="checkbox"/>	<input type="checkbox"/>	0 times per day for 0 years
Smokeless	<input type="checkbox"/>	<input type="checkbox"/>	0 times per day for 0 years
Alcohol	<input type="checkbox"/>	<input type="checkbox"/>	0 drinks per day for 0 years
Drug Abuse	<input type="checkbox"/>	<input type="checkbox"/>	

Figure 177: Social Hx (History) Screen

**Note:** All entries made into the **Social History**, **OHA**, **Occlusion**, **Parafunctional Habits** and **TMJ** modals NOT filed are cleared by selecting the **Cancel** button.

## Multiple Filings to Same Modal on Same Day

The following occurs when any modal has multiple filed findings and/or history during the same day on the same patient's chart; subsequent filings may be changed/edited or remain the same as the previous filing.

The last modal filing during one calendar day is the only record present from the date drop-down menu of that modal screen. The TIU progress notes filed that day would have the record of the modal data entered that session/encounter.

The first provider (HYG) may file findings in the **OHA** modal and that data displays in the screen for every other user of DRM Plus to review.

**DRM Plus - Oral Health Assessment**

Date: 05/26/2017 [New] [Definitions]

**Clinical Findings**

☐ Patient has no remaining functional teeth, roots or implants.

**Plaque Index**

☐ 0 - None  
☐ 1 - Slight  
☒ 2 - Moderate  
☐ 3 - Heavy  
☐ 4 - Not Recorded

**Xerostomia**

☐ 0 - None  
☐ 1 - Slight  
☐ 2 - Moderate  
☐ 3 - Significant  
☒ 4 - Not Recorded

**Caries Risk**

☐ 0 - Edentulous  
☐ 1 - Low  
☐ 2 - Moderate  
☐ 3 - High  
☒ 4 - Not Recorded

**Oral Hygiene**

☐ 0 - Excellent  
☒ 1 - Good  
☐ 2 - Fair  
☐ 3 - Poor  
☐ 4 - Not Recorded

Date	Provider	NFT	PI	X	CR	OH
05/26/2017	HYG		2			1
04/28/2017	ADMIN		1	0	1	1

[OK] [Cancel]

[VA Dental Definitions](#)

**Definitions Panel**

**PLAQUE INDEX**

0 - None:  
No Plaque in the gingival area.

1 - Slight:  
A film of plaque adhering to the free gingival margin and adjacent area of the tooth. The plaque may be recognized only by running a probe across the tooth surface.

2 - Moderate:  
Moderate accumulation of soft deposits within the gingival pocket and on the gingival margin and/or adjacent tooth surface that can be seen by the naked eye.

3 - Heavy:  
Abundance of soft matter within the gingival pocket and/or on the gingival margin and adjacent tooth surface.

**XEROSTOMIA**

0 - None:  
None.

1 - Slight:  
Moist mucosa with limited flow and no pooling of saliva. Patient may report subjective dryness and/or stale breath.

2 - Moderate:  
In addition to slight symptoms, initially moist mucosa but dries during oral examination. Saliva may be thick and ropery. Patient may reports difficulty with dry foods, altered taste

Figure 178: Provider HYG Filed OHA Data First Today

If the second provider (ADMIN) changes/edits the filed **OHA** findings from the first provider (HYG) on the same calendar day; modifies the findings in the **OHA** modal. The modified **OHA** findings filed by the second provider will be the data viewed by the local clinical providers in the patient's **OHA** screen.

**DRM Plus - Oral Health Assessment**

Date: 05/26/2017 [New] [Definitions]

**Clinical Findings**

☐ Patient has no remaining functional teeth, roots or implants.

**Plaque Index**

☐ 0 - None  
☐ 1 - Slight  
☒ 2 - Moderate  
☐ 3 - Heavy  
☐ 4 - Not Recorded

**Xerostomia**

☒ 0 - None  
☐ 1 - Slight  
☐ 2 - Moderate  
☐ 3 - Significant  
☐ 4 - Not Recorded

**Caries Risk**

☐ 0 - Edentulous  
☐ 1 - Low  
☒ 2 - Moderate  
☐ 3 - High  
☐ 4 - Not Recorded

**Oral Hygiene**

☐ 0 - Excellent  
☒ 1 - Good  
☐ 2 - Fair  
☐ 3 - Poor  
☐ 4 - Not Recorded

Date	Provider	NFT	PI	X	CR	OH
05/26/2017	ADMIN		2	0	2	1
04/28/2017	ADMIN		1	0	1	1

[OK] [Cancel]

[VA Dental Definitions](#)

**Definitions Panel**

**PLAQUE INDEX**

0 - None:  
No Plaque in the gingival area.

1 - Slight:  
A film of plaque adhering to the free gingival margin and adjacent area of the tooth. The plaque may be recognized only by running a probe across the tooth surface.

2 - Moderate:  
Moderate accumulation of soft deposits within the gingival pocket and on the gingival margin and/or adjacent tooth surface that can be seen by the naked eye.

3 - Heavy:  
Abundance of soft matter within the gingival pocket and/or on the gingival margin and adjacent tooth surface.

**XEROSTOMIA**

0 - None:  
None.

1 - Slight:  
Moist mucosa with limited flow and no pooling of saliva. Patient may report subjective dryness and/or stale breath.

2 - Moderate:  
In addition to slight symptoms, initially moist mucosa but dries during oral examination. Saliva may be thick and ropery. Patient may reports difficulty with dry foods, altered taste

**Figure 179: Provider ADMIN Edited and Filed the second OHA Data Today**

**Note:** The first provider's **OHA** findings are only present in their filed TIU progress note.

**Note:** The **OHA** modal and **H&N** modal are the only screens that displays the date and the provider's initials filed with the findings and/or history data.

# Chart/Treatment – Periodontal Chart

Enter periodontal information on the patient from the **Periodontal Chart** screen, within the **Chart/Treatment** tab.

**Dental Record Manager Plus in use by: DRMPROVIDER,DENTIST for DRMPATIENT,ONE**

File Edit Dental Encounter Data Treatment & Exam Tools Reports Help

DRMPATIENT,ONE 000-00-0001 01/01/1960 Age: 57 Jun 12, 17@09:00 DENTAL DRMPROVIDER,DENTIST PDP: DRMPROVIDER,ADMINDE SDP: Dental Class OPC, Class IV

Cover Page Clinical Record Dental History Chart/Treatment Exam

**Treatment & Exam**

History \*PSR\* Compare Stats Summary \*OHA\* \*H&N\* Notes

Pocket FGM MGJ Bleeding Delayed Supp'n Mobility Furcation Clear

0 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 X

Cal Undo Keyboard Mode [F10] Adv Back =>

Upper Facial

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

	1	2	3	4	5	6	7	8
<b>Pocket</b>		4 4 4	4 3 3	3 2 2	2 3 4	- - -	4 3 2	2 2 2
<b>Bleeding</b>		B b			B		b	
<b>Suppuration</b>								
<b>FGM</b>		2 2 2	1 - -	- - -	- - 1	- - -	1 - -	- - -
<b>MGJ</b>		- - -	- - -	- - -	- - -	- - -	- - -	- - -
<b>Furcation</b>		- - -	- - -	- - -	- - -	- - -	- - -	- - -
<b>Attachment</b>		6 6 6	5 3 3	3 2 2	2 3 5	- - -	5 3 2	2 2 2
<b>Attach Chg</b>		0 0 0	0 0 0	0 0 0	0 0 0	0 0 0	0 0 0	0 0 0
<b>Mobility</b>		0½	-	-	0½	-	0½	-


Next =>

**Figure 180: Periodontal Chart Screen**

The upper left side of the screen shows the **History**, **Compare**, **Summary**, **H&N**, **PSR**, **Stats**, **OHA** and **Notes** buttons. The top center of the screen displays various periodontal condition-specific icons. Use these to mark periodontal findings on the patient's chart.

Use the options on the upper-right side of the screen to adjust the view of the tooth/arch graphic.

The center of the screen features the tooth/arch graphic. Clicking various areas in combination with the condition-specific icons located in the top center of the screen enters information into the patient's chart.


The bottom of the screen shows text (only a quadrant is viewable depending where the cursor shield  is located) of the periodontal information entered using the graphic and the condition-specific icons in the top center of the screen.

**Note:** The **Periodontal Chart** screen allows two providers to enter perio data on the same day. However, it only displays the last perio data that was filed on that day. The first entered perio data is only viewable in the TIU progress note of the provider that filed the first encounter.



## History and Compare Buttons

Clicking the **History** or **Compare** buttons display similar screens. Use the **Compare** function to see two periodontal charts for a patient.



**Figure 181: Periodontal History/Compare Screen – Compare Viewable**

The information is color-coded. The data from the first date is displayed in red, while the data from the second date is displayed in blue. Use the drop-down menus to change the dates.

Click the **History** button to view a graphic of the patient's history.




Figure 182: Periodontal History/Compare Screen - History Viewable

Use the drop-down menu to change the date. Use the **Notes** button to view previously entered notes concerning the patient's periodontal history.

### Summary Button

See the **Summary Button** section in the Chart/Treatment – Treatment & Exam chapter of this manual.

### H&N Button


See the **H&N Button** section in the Chart/Treatment – Treatment & Exam chapter of this manual.

### PSR Button

See the **PSR Button** section in the Chart/Treatment – Treatment & Exam chapter of this manual.

## Stats Button

Select the **Stats** button to view the patient's total number of periodontal warning levels.



The screenshot shows a window titled "Periodontal Statistics" with a close button (X) in the top right corner. The window contains a table with five columns: "Description (Warning Level)", "Teeth < Warning", "Sites < Warning", "Teeth >= Warning", and "Sites >= Warning". The table lists eight conditions: Pocket Depth (4), Free Gingival Margin (4), Mucogingival Junction (3), Furcation (2), Bleeding (0), Delayed Bleeding (0), Suppuration (0), and Mobility (1). The values for "Teeth >= Warning" and "Sites >= Warning" are displayed in red text. An "Ok" button is located at the bottom right of the window.


Description (Warning Level)	Teeth < Warning	Sites < Warning	Teeth >= Warning	Sites >= Warning
Pocket Depth (4)	15	124	12	38
Free Gingival Margin (4)	27	162	0	0
Mucogingival Junction (3)	27	123	0	0
Furcation (2)	10	28	0	0
Bleeding (0)	23	157	4	5
Delayed Bleeding (0)	22	156	5	6
Suppuration (0)	26	161	1	1
Mobility (1)	26		1	


Figure 183: Periodontal Statistics Screen

## OHA (Oral Health Assessment) Button

See the **OHA Button** section in the Chart/Treatment – Treatment & Exam chapter of this manual.

## Notes Button

Select the **Notes** button  to enter general patient notes or tooth specific notes into the chart.



**Figure 184: Notes Screen**


Choose between creating a more general note on the patient by using the tools on the left side of the screen, or a note about a specific tooth on the right side of the screen.

- To create a new note, click the **New Entry** button.
- To view all general patient notes in the patient's chart, click the **Show All** button.
- To make an entry about a specific tooth, use the drop-down menu on the right side of the screen to select a tooth before making a new entry.

To view any tooth-specific note that has been entered during the past year into the patient's chart, click the **History** button. View the note by clicking the selected history's **Notes** button, according to the appropriate visit date. Note that this is the same screen that displays if the **Notes** icon is clicked on the **Treatment & Exam** screen.


## Entering Periodontal Information

The condition-specific icons in the top center of the screen works in combination with the tooth/arch graphic.



**Figure 185: Periodontal Condition-Specific Icons and Periodontal Graphical Chart**

To enter periodontal information:


1. Click the desired condition-specific icon (**Pocket, FGM, Bleeding, Mobility, etc.**).
2. On the graphic, the cursor shield  is the graphical pocket location for perio data entry.
3. Click the desired **number** below the list of condition icons (if applicable). **Bleeding, Delayed Bleeding** and **Suppuration** do NOT require any clicks on a number.
4. Use the buttons on the top right of the screen to view different areas of the tooth/arch, or to change the view.
5. The condition and location display on the graphic and on the transaction table below the graphic. The transaction table only displays a quadrant of the upper/lower arch.


**Note:** Only one periodontal chart exam should be completed per day, per patient. When a second exam is completed during the same day, the periodontal chart history only saves the second exam. TIU progress notes in VistA still have all the data entered.


## **Other Tools**


The last line of buttons in the top center of the screen features several tools.


**Key** button  displays keyboard shortcuts.


**Undo** button  undoes the last action performed on the screen. This action is limited to the last nine actions performed.

**Keyboard Mode** button  switches the program to keyboard mode.

**Advance** button  moves the cursor shield on the graphic forward due to the direction of the arrow.

**Back** button  moves the cursor shield on the graphic back due to the direction of the arrow.

**Direction** button  shows the direction the cursor shield moves on the graphic. The direction of advancement depends on the orientation of this button.


**Auto Advance** button  when the button is active, the cursor shield automatically moves in the way designated for the perio exam entry of the condition-specific icons. It does **NOT** allow **Bleeding**, **Delayed Bleeding** or **Suppuration**, to automatically move.

The auto exam sequence is a parameter adjustment that may be redesigned by each provider. For more information, please see the Treatment & Exam/Show Configuration/Periodontal section of the Using the DRM Plus Drop-Down Menus chapter of this manual.

# Exam




Providers can file required data using a national standard exam style format for each exam/consult code (**D0120, D0140, D0150, D0160, D0170, D0180** and **D9310**) in conjunction with the **Exam** tab in DRM Plus. Mandatory elements for each exam/consult code and requirements for each element are based on the user's exam/consult code selection. Initially, each element is marked with required or optional icon. The **Exam** tab interfaces automatically with existing DRM Plus modals (i.e. Head & Neck) for easy data entry. The **Exam** tab, when activated, generates a TIU progress note associated with a specific visit containing the entire exam or consult's required information, along with other DRM Plus note objects (i.e. dental alerts, etc.).

To proceed to the **Exam** tab, the user is prompted to select one of the six exam or consult procedure codes from the **Treatment & Exam/Completed Care** view screen.




**Figure 186: Exam Tab Message**

Entering the exam/consult procedure activates the **Exam** tab and displays the procedure in the **Exam Type** drop-down menu. One way to change the exam/consult procedure after selecting it and determining it was incorrect is to use the drop-down menu located on the **Exam** tab. Another way to change the exam/consult procedure is by deleting the procedure from the **Completed Care** view screen and entering a new one.

Selecting the exam/consult code from the **Completed Care** view screen triggers all the elements on the **Exam** tab with a required  icon or optional  icon. When these elements are satisfied, a completed  icon displays. Some elements automatically pull data from the modals when entered from the **Treatment & Exam** screen.

**Note:** Users may enter additional optional information in each element, if desired, for the selected exam/consult code.

In the following example, the **D0150**, comprehensive exam, was selected from the **Treatment & Exam/Completed Care** view screen. Upon selecting the **Exam** tab, **D0150 Comprehensive** displays in the **Exam Type** drop-down menu. Twelve of the sixteen elements require data entry by the provider when selecting the D0150, comprehensive exam.




**Figure 187: Exam Tab**

The seven exam/consult codes each have a different set of required and optional elements activated when selecting a specific code. The dialogue found in each element section that follows are those that would display if the user selected D0150, comprehensive exam, during a session.

The **Back** button located on the **Exam** tab screen returns the end-user to the **Treatment & Exam** screen.

The **Next** button located on the **Exam** tab allows the end-user to proceed to the **Filing Options** screen, which is the next screen when completing the encounter. This **Next** button also opens any required element if that element has NOT been completed.

The **Back/Next** buttons located on each element screen only move backward or forward to other element screens. The **Next** button allows the provider to move to the next required element for this exam type depending on the parameter the user selected. Usage of the **Back** button is NOT dependent on the element's completion.

**Note:** The  icon on the right of the element button displays when there was previously filed data on the same day and that data is associated with the element it is displayed on.



## Exam Elements

### Presentation/Chief Complaint Element

The **Presentation/Chief Complaint** element is required for all seven exam/consult procedure codes and defaults open when the **Exam** tab is selected. The presentation of the exam/consult code is automatically imported and displayed at the top of the element. This element requires one of the two radio buttons to be selected. The selection of the second radio button, **Patient presents with dental complaint(s)** opens two text boxes which require a text entry only in the first text box intended for the dental complaint of the patient. The second text box (optional text entry) allows data entry for the history of the patient's present illness.

**Presentation/Chief Complaint**

Presentation/Chief Complaint

☒ Presentation  
Patient presents for comprehensive oral evaluation

☐ Chief Complaint

☐ Patient has no dental complaints

☐ Patient presents with dental complaint(s):

Additional Annotations

Annotations

Presentation/Chief Complaint:  
Patient presents for comprehensive oral evaluation

OK Cancel Next >>

**Element Requirements Panel**

D0150 Comprehensive  
Presentation/Chief Complaint requirements:

This element requires one of the two radio buttons to be selected.

The selection of the second radio button will open two text boxes which require at a minimum a statement entered in the first text box.

**Figure 188: Presentation/Chief Complaint Element Screen**

The following information is the same with all sixteen elements except the **Presentation/Chief Complaint** does NOT have a **Back** button, and **Disposition** does NOT have a **Next** button.

**Additional Annotations** is a free text window, which allows the provider to enter additional information about the patient's chief complaint. It offers right-click functionality for importing .txt files, if desired.


**Annotations** is view-only and captures everything entered using this element. Select the **OK** button to save all required information entered and close the element or click the **Next** button to move to the next required element, depending on the selected parameter.

## Vitals Elements

The **Vitals** element is required for all seven exam/consult procedures. This element requires one of the two radio buttons to be selected. The second radio button, **Vital signs obtained**, defaults to any vitals that have been entered, using the **Enter Vitals** screen, during 24 hours of the visit date. The **Provider/Visit** for this encounter also must be entered in the banner for this feature to work. Entering the **Provider/Visit** should be the first action taken when opening any patient's chart for a new encounter.

Figure 189: Vitals Element Screen

If no vitals have been entered using the **Enter Vitals** screen, the information can be entered manually; however, no date is attached to these entries. **Dental Pain** is the only required vital sign.

To enter vitals, first select the scheduled appointment or visit and then click the **Vitals** button  located on the DRM Plus banner. The user can also click the **Enter Vitals** button in the lower left portion of the **Vitals** element screen.

Specific vital information can be entered into the **Additional Annotations** free text window. Right-click this window to import text files, if desired. All information entered using the **Vitals** element appear in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move to the next element, or the **OK** button to close the element. The **Back** button moves to the previous required element, depending on the selected parameter, and is NOT dependent on the element's completion.

**Note:** The only vital sign, **Dental Pain**, may be saved as unfiled data for this element.

## PMH (Past Medical History) and Medications Element

The **PMH** element is required for all seven exam/consult procedures. This element requires one of the three radio buttons to be selected. The selection of the first radio button opens an optional text box to enter additional information if the patient is new to the clinic. The selection of the third radio button opens a required text box to enter any significant changes noted since the last dental visit.

The eight positive/negative check box conditions, one free text positive condition or the five **Imports** check boxes are optional entries of patient information for this element. The user may select one import such as the patient's medications or use the **Select ALL Imports** button to import all four previously filed medical histories about the patient which are being stored in a VistA database.

Figure 190: PMH Element Screen

**Additional Annotations** is a free text window, which allows the provider to enter additional information about the patient's past medical history. It offers right-click functionality of import text files, if desired.

**Annotations** is view-only and captures everything entered using this element.

Once the element is complete, click the **Next** button to move to the next element, or the **OK** button to close the element. The **Back** button moves to the previous required element, depending on the selected parameter, and is NOT dependent on the element's completion.

## Social History Element

The **Social History** element is required for the D0150 and D0180 exams. This element requires new social history findings entered with the **Social History** screen when completing one of the two required exams. The **Social History** screen may be opened with the specialty button located on the **Treatment & Exam** screen, or by the **Social History** button located in the lower left corner of this screen. The minimum requirement to enter a new social history entry is to select at least one historical finding from the **Social History** screen.

**Social History**

Social History

You have not completed the following section: Social History Findings

You can do this in the Chart/Treatment tab of DRM Plus.

Alternately, you can click on the button below for the specific section.

Additional Annotations

Annotations

**Element Requirements Panel**

D0150 Comprehensive Social History requirements:

The minimal requirement to enter a new Social History entry is to select at least one entry from the Social History screen.

Selection of the second radio button requires at least one selection from the field now activated.

Social History

<< Back OK Cancel Next >>

**Figure 191: Social History Element Screen**

Additional information regarding the patient can be entered using the **Additional Annotations** window. Right-click in this window to import text files, if desired. All information added to this element displays in the read-only **Annotations** window. Once this element is complete, click the **Next** button to move on to the next element.

## H&N (Head and Neck) Findings Element

The **H&N Findings** element is required for the D0120, D0150 and D0180 exams. This element requires a new H&N finding or historical entry using the **H&N Findings** screen. This element imports data entered from the **H&N Findings** screen. The **Screening Negative** button on this element's screen allows a new screening negative entry directly into the element and import it into the **H&N Findings** screen for the patient's permanent record.

**H&N Findings**

H&N Findings

You have not completed the following section: Head and Neck Findings

You can do this in the Chart/Treatment tab of DRM Plus.

Alternately, you can click on the button below for the specific section.

Additional Annotations

Annotations

H&N Findings Screening Negative << Back OK Cancel Next >>

**Element Requirements Panel**

D0150 Comprehensive  
H&N Findings requirements:

This element will require one new H&N finding or historical entry using the H&N Findings screen.

The "Screening Negative" button on this element's screen will allow a new negative screening entry directly into the element and import it into the H&N Findings screen for the patient's permanent record.

Figure 192: H&N Findings Element Screen

The **H&N Findings** screen can be accessed by either clicking the specialty button from the **Treatment & Exam** screen or clicking the **H&N Findings** button in the left bottom corner of this element. Once the H&N findings are added to this element, additional information regarding the patient can be entered using the **Additional Annotations** free text window. Right-click this window to import text files, if desired. All information entered using this element displays in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

## Radiographic Findings Element

The **Radiographic Findings** element is required for the D0150 and D0180 exams. The radiographic element requires at least one selected check box from the top six options. The provider may select any combination of the top six check boxes for the patient's TIU progress note. The fourth check box down the left column requires some data entry in the text box or at least one per-defined statement to satisfy the requirements.

The **Exams Settings** tab located using the **Tools** menu → **User Options** submenu → **Exam Settings** tab → **Configure Radiographs Requirement** parameter that allows each user to require a radiographic finding for all exam/consult codes with mandatory entries in the **Exam** tab.

**Radiographic Findings**

Radiographic Findings

☐ No radiographs obtained this visit

☐ Patient declined radiographs

☐ Radiographic findings consistent with charted entries

☐ Radiographs reviewed, findings noted below:

☐ No radiographic caries noted

☐ No apparent bony pathology noted

Additional Annotations

Annotations

Radiographic Findings:

Element Requirements Panel

D0150 Comprehensive  
Radiographic Findings requirements:

This element will require at least one check box selected from the top six.

The provider may select any combination of the top six check boxes for the patient's record.

The selection of the fourth check box in the left column requires a text description about radiographic findings and/or selection of at least one pre-defined statement.

<< Back OK Cancel Next >>

**Figure 193: Radiographic Findings Element Screen**

Up to twelve pre-defined statements on radiographic findings can be selected from this screen. The check boxes found in the pre-defined statements window have three national radiographic findings statements pre-loaded; however, all twelve statements can be created locally.

The DRM Administrator has priority over all other users and may add or delete all twelve, if desired.

Additional information regarding the patient can be added using the **Additional Annotations** free text window. The user can right-click to import a text file, if desired. All information entered using this element displays in the **Annotations** read-only window. Once the element is complete, click the **Next** button to move on to the next element.

## Diagnostic Findings Element

The **Diagnostic Findings** element is required for all seven exam/consult procedures. All exam/consult codes are required as appropriate for new and updated findings. The second check box only displays after data is entered from at least one **Treatment & Exam** findings screen that satisfies this option. Informational screens advise the user of any missing requirements for a specific exam code.

The **NFT**, no functional teeth, check box when selected in the **OHA** screen, bypasses all requirements in this element for all exam/consult procedures. The **OHA** screen may be opened from the **Treatment & Exam** screen or by selecting the **OHA** button in the lower left area of this screen.

- **D0120:** Requires a **Plaque Index** entry from the **OHA** screen. Also requires the selection of a **Mobility** radio button, as it pertains to the patient.
- **D0150:** Requires at least one entry from the **Diagnostic Findings** view screen or the first check box of no apparent pathology selected. Requires a **Plaque Index** entry from the **OHA** screen. Also requires the selection of a **Mobility** radio button, as it pertains to the patient.
- **D0180:** Requires an **Oral Hygiene** entry from the **OHA** screen.

**Diagnostic Findings**

Diagnostic Findings

Charting

☐ Dentition exhibits no apparent evidence of dental pathology at this visit

Screening Negative

☐ Other:

H&N Findings

Parafunctional Habits

TMJ Findings

Occlusal Findings

Tooth Mobility

☐ Tooth Mobility Not Assessed

☐ No Significant Tooth Mobility Noted

☐ Tooth Mobility Noted (see Periodontal charting)

☐ Tooth Mobility Noted:

Additional Annotations

Annotations

Element Requirements Panel

D0150 Comprehensive Diagnostic Findings requirements:

This element requires complete exam findings charted.

This element requires an entry in Diagnostic Findings section of the Chart or the first statement of no apparent pathology.

The second check box becomes viewable and will only display after new patient findings are entered into any other DRM Plus findings screen.

This element requires the second check box to have a Plaque Index entry from the OHA screen.

This element requires one of the four radio buttons for Tooth Mobility. If "Tooth Mobility Noted" is selected, then details must be entered.

When the OHA check box "Patient has no remaining functional teeth, roots or implants" is selected this element will be satisfied.

Diagnostic Findings Oral Health Assessment << Back OK Cancel Next >>

Figure 194: Diagnostic Findings Element Screen

The user can enter any additional information regarding the patient using the **Additional Annotations** free text window. Right-click this window to import a text file, if desired. All information entered using this element displays in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

## Periodontal Assessment Element

The **Periodontal Assessment** element is required for the D0120, D0150 and D0180 exams. The D0120 and D0150 exams required at least one selection from the **Periodontal General Assessment** section and a **PSR**, Periodontal Screening/Recording, screening. The **Detailed Assessment** button allows the user to enter additional perio data; however, this is optional for the D0120 and D0150 exam codes.

The check box, **Include Last Perio Chart**, defaults as unchecked. If the provider wants to import the last filed periodontal chart into this element the provider will need to select the check box. When the **Periodontal Assessment** element is **NOT** required for an exam and the provider wants to enter some perio data; at least one selection from the **Periodontal General Assessment** section, or entering new periodontal charting, or importing the last periodontal chart would be the minimal requirement.

The **NFT**, no functional teeth, check box when selected in the **OHA** screen, bypasses all requirements for the D0120, D0150 or D0180 exams in the periodontal element. The user may access the **Periodontal Chart** screen, **OHA** screen or **PSR** screen using the buttons found on this **Periodontal Assessment** element screen.

**Periodontal Assessment**

Periodontal General Assessment

☐ Good Periodontal Health

	Acute	Chronic	Generalized	Localized	Slight	Moderate	Severe	Aggressive
<input type="checkbox"/> Gingivitis								
<input type="checkbox"/> Gingival Enlargement								
<input type="checkbox"/> Generalized Periodontitis								
<input type="checkbox"/> Localized Periodontitis								
<input type="checkbox"/> Peri-implantitis								

☐ Include Last Perio Chart

Additional Annotations

Annotations

PSR Oral Health Assessment << Back OK Cancel Next >>

**Element Requirements Panel**

D0150 Comprehensive Periodontal Assessment requirements:

This element requires at least one selection from the Periodontal General Assessment section and completion of a PSR.

The "Detailed Assessment" button allows the user to enter additional periodontal data; however, this is optional for this exam code.

The "Include Last Perio Chart" check box will include the last periodontal chart when selected.

When the OHA check box "Patient has no remaining functional teeth, roots or implants" is selected this element will be satisfied.

Figure 195: Periodontal Assessment Element Screen

The D0180 exam requires one selection from the **Periodontal General Assessment** section as well. The **Periodontal Detailed Assessment** section is optional and has optional text boxes with each selection if more descriptive detail is needed.



The D0180 exam also requires the first four rows in the **Additional Periodontal Detail** section to have at least one selection. The **Additional Periodontal Comments** text box is optional in this section. When the **Other** check box is selected from the **Past Periodontal Tx History** options, data entry is required in the **Additional Periodontal Comments** text box.

**Periodontal Assessment**

Periodontal General Assessment

☐ Good Periodontal Health

☐ Gingivitis

☐ Gingival Enlargement

☐ Generalized Periodontitis

☐ Localized Periodontitis

☐ Peri-implantitis

Periodontal Detail Assessment

☐ Mucogingival defect

☐ Failed implant

Additional Annotations

Annotations

PSR Oral Health Assessment << Back OK Cancel Next >>

**Element Requirements Panel**

D0180 Perio  
Periodontal Assessment requirements:

This element requires one selection from the Periodontal General Assessment section.

The Periodontal Detail Assessment section is optional and can be selected with details completed as needed.

This element requires every row in the Additional Periodontal Details to have at least one selection.

The Additional Periodontal Comments are optional.

When the "Other" check box is selected in the "Past Periodontal Tx History" row a text entry is required in the "Additional Periodontal Comments" text box.

When the OHA check box "Patient has no remaining functional teeth, roots or implants" is selected this element will be satisfied.


**Figure 196: Periodontal Detailed Assessment Button Selected**

**Note:** The D0180 exam code does NOT allow the user to select the **Brief Assessment** button from this screen.

Additional information regarding the patient may be entered using the **Additional Annotations** free text window. The user may right-click in this window to import a text file, if desired. All information entered using this element appears in the read-only **Annotations** window. Once this element is complete, click the **Next** button to move on to the next element.

The provider may open the **OHA** screen or the **PSR** screen using the buttons in the lower left area of this screen or from the **Treatment & Exam** screen to enter an OHA finding or a PSR exam. The **Periodontal Chart** button on the middle right-side of this screen allows access to that chart to enter any new findings for this session.

The following dialog displays the VA Office of Dentistry perio definitions when selecting the **Definitions** button from the **Periodontal Assessment** screen.



**Figure 197: VA Dentistry Periodontal Definitions**

## Parafunctional Habits Element

The **Parafunctional Habits** element is optional for all seven exam/consult types. It imports all data entered in this session from the **Parafunctional Habits** screen, or simply is left blank. The **Parafunctional Habits** screen can be opened by clicking the **Habits** (Parafunctional) button from the **Treatment & Exam** screen, or by clicking the **Parafunctional Habits** button in the bottom left portion of this element.

The minimum requirement to enter a new parafunctional habit finding is to select at least one history or one clinical finding from the **Parafunctional Habits** screen.

The screenshot shows the 'Parafunctional Habits' window. The main area contains a message: 'You have not completed the following section: Parafunctional Habits. You can do this in the Chart/Treatment tab of DRM Plus. Alternately, you can click on the button below for the specific section.' Below this is an 'Additional Annotations' text box and an 'Annotations' list box. On the right, the 'Element Requirements Panel' for 'D0150 Comprehensive Parafunctional Habits requirements' states: 'This element is optional. However if you choose to enter any data for this element, please note: The minimal requirement is to enter at least one historical or one clinical finding from the Parafunctional Habit screen.' The bottom of the window has a 'Parafunctional Habits' button, '<< Back', 'OK', 'Cancel', and 'Next >>' buttons.

**Figure 198: Parafunctional Habits Element Screen**

Additional information regarding the patient may be entered using the **Additional Annotations** free text window. The user can right-click in this window to import a text file, if desired. All information entered using this element displays in the read-only **Annotations** window. Once this element is complete, click the **Next** button to move on to the next element.

## TMJ Findings Element

The **TMJ Findings** element is required for the D0120, D0150 and D0180 exams. It requires new TMJ findings entered from the **TMJ** screen, which can be opened either by clicking the **TMJ** button from the **Treatment & Exam** screen, or the **TMJ Findings** button located in the bottom left portion of this element.

The minimum requirement to enter a new TMJ exam finding is to select at least one historical or clinical finding from the **TMJ** screen. Additional information regarding the patient can be entered using the **Additional Annotations** free text window. The user can right-click in this window to import a text file, if desired.

**Figure 199: TMJ Findings Element Screen**

All information entered using this element appears in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

## Occlusal Findings Element

The **Occlusal Findings** element is required for every D0150 and D0180 exams. Occlusal findings must be entered into and imported from the **Occlusion** screen when completing one of these exams. The **Occlusion** screen can either be opened by clicking the **Occl** (Occlusion) button from the **Treatment & Exam** screen, or by clicking the **Occlusal Findings** button in the bottom left corner of this screen.

From the **Clinical Findings** section, the **Mandibular relationship\*** is the only required field. All data from the last filed occlusion exam imports into the screen with a new exam, which requires the provider to remove and/or add correct data for the new exam.

**Occlusal Findings**

Occlusal Findings

You have not completed the following section: Occlusal Findings

You can do this in the Chart/Treatment tab of DRM Plus.

Alternately, you can click on the button below for the specific section.

Additional Annotations

Annotations

**Element Requirements Panel**

D0150 Comprehensive  
Occlusal Findings requirements:

This element requires Occlusal Findings entered from the Occlusion screen.

In the "Clinical Findings" section, the drop-down menu option 'Mandibular relationship' is the only required field.

All data from the last filed Occlusion exam will import into the screen when a new exam is selected allowing the provider to accept or change findings for the current exam.

Occlusal Findings << Back OK Cancel Next >>

**Figure 200: Occlusal Findings Element Screen**

Enter additional information regarding the patient using the **Additional Annotations** free-text window, which offers right-click functionality to import text files, if desired. All information entered using this element displays in the read-only **Annotations** window. When the element is complete, click the **Next** button to move on to the next element.

## Salivary Flow Element

The **Salivary Flow** element is optional for all seven exam/consult procedures. This element requires one of the two radio buttons to be selected when entering data. The second radio button; **Clinically abnormal salivary quantity and/or quality noted**, requires a statement to be entered in the text box.

The **Xerostomia** value and description imports for viewing on this **Salivary Flow** element screen if entered from the **OHA** screen during this session.

**Salivary Flow**

☐ Clinically normal salivary quantity and quality noted

☐ Clinically abnormal salivary quantity and/or quality noted:

NOTE: Xerostomia may be entered in Oral Health Assessment (OHA) in addition to any comments made here.

Additional Annotations

Annotations

Salivary Flow Findings:

Oral Health Assessment

<< Back OK Cancel Next >>

**Element Requirements Panel**

D0150 Comprehensive Salivary Flow requirements:

This element is optional.

However if you choose to enter any data for this element, please note:

The minimal requirement is to select one of the two radio buttons in the element.

The second radio button option requires a text entry in the text box.

**Figure 201: Salivary Flow Element Screen**

The user can enter additional information regarding the patient using the **Additional Annotations** free text window. This window offers right-click functionality to import text files, if desired. All information entered using the element displays in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

## Removable Protheses Element

The **Removable Protheses** element is optional for all seven exam/consult types. This element requires one of the top three radio buttons to be selected when entering data. If the user selects the third radio button; **Patient presents with removable prothes(es)**, the user must then select one of either the **Maxillary** or **Mandibular** radio buttons: **Partial** or **Complete**. Only one is required and, once selected, it pens two more radio buttons: **Satisfactory** or **Unsatisfactory**.

If the user selects **Unsatisfactory**, four check boxes (**Occlusion**, **Retention**, **Stability** and **Esthetics**) and a text box become active. The user needs to check any box pertaining to the patient and note any additional prostheses in the text box.

The **Other Protheses** check box opens a required text box for any other prostheses that should be added to the TIU progress note for the patient.

**Removable Protheses**

Removable Protheses

☐ Patient has no removable prothes(es)

☐ Patient's protheses not evaluated at this time

☐ Patient presents with removable prothes(es):

**Maxillary**

☐ Partial ☐ Complete

☐ Satisfactory ☐ Unsatisfactory

☐ Occlusion ☐ Retention

☐ Stability ☐ Esthetics

**Mandibular**

☐ Partial ☐ Complete

☐ Satisfactory ☐ Unsatisfactory

☐ Occlusion ☐ Retention

☐ Stability ☐ Esthetics

☐ Other Protheses:

Additional Annotations

Annotations

Removable Protheses Findings:

Element Requirements Panel

D0150 Comprehensive Removable Protheses requirements:

This element is optional.

However if you choose to enter any data for this element, please note:

The minimal requirement is to select one of the top three radio buttons in the element.

The selection of the third radio button will allow the selection of another radio button listed in a maxillary or mandibular column.

Only one selection of a maxillary or mandibular radio button is allowed which is followed by another selection of a satisfactory or unsatisfactory radio button.

When selecting the unsatisfactory option there may be up to four possible descriptive words for selection or a text box to enter something that may not be listed.

The check box Other Protheses opens a required text box for any other prostheses that should be added in the progress note for the patient.

<< Back OK Cancel Next >>

Figure 202: Removable Protheses Element Screen

The user may enter any additional information regarding the patient using the **Additional Annotations** text window. This window offers right-click functionality to import text files, if desired. All information entered using the element displays in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

## Assessment/Plan Element

The **Assessment/Plan** element, comprised of an assessment and planned section, is required for all seven exam/consult procedures. The top **Assessment Summary** section is optional for the completion of this element.

The **Treatment Plan** section requires one of the four check boxes or only one pre-defined statement to be selected to complete the element. The first check box; **Include charted treatment plan**, loads automatically and imports the patients newly entered and/or past planned treatment.

Up to twelve pre-defined statements for the **Assessment Summary** and twelve **Treatment Plan** statements can be selected from this screen. The check boxes found in the pre-defined statement windows have three national assessment and three national planned statements pre-loaded; however, all twelve statements from either can be created locally. The DRM Administrator has priority over all other users and may add or delete all twelve, if desired.

DRM Assessment/Plan (review/edit)

Assessment/Plan (review/edit)

Assessment Summary

☐ No urgent dental needs or acute dental infections noted on examination.  
☐ No contraindications for planned procedure(s).  
☐ Dentition in satisfactory repair and function.  
☐ Clinic wide Assessment Summary pre-defined statement entered.  
☐ First end-user Assessment Summary pre-defined statement entered.  
☐ Second end-user Assessment Summary pre-defined statement entered.

Treatment Plan

☒ Include charted treatment plan ☐ Referred to assigned dentist for further treatment planning  
☐ No treatment required at this time ☐ Final treatment plan to be completed at a later time

☐ Reviewed risks/benefits/alternatives associated with the proposed treatment plan. Patient agrees to  
☐ Adequate dentition for mastication. Replacement of missing teeth is not indicated.  
☐ Patient is not eligible for replacement of teeth through VA.

Additional Annotations

Annotations

Assessment/Plan:  
Planned Procedures:  
Non-VA Care  
(D7950) MANDIBLE GRAFT: . DX: ( ) .  
Phase 2  
(D2150) AMALGAM TWO SURFACES PERMANE: 14(OL) . DX: ( ) . (Next-Appt)  
(D2160) AMALGAM THREE SURFACES PERMA: 15(DOL) . DX: ( ) . (Next-Appt)

Treatment Plan << Back OK Cancel Next >>

Element Requirements Panel

D0150 Comprehensive Assessment/Plan (review/edit) requirements:

In the "Assessment Summary" section, the text box and pre-defined statements are optional. It is highly recommended at a minimum an appropriate pre-defined statement is selected for detail.

The "Treatment Plan" section requires one of the four check boxes, and/or one pre-defined statement to be selected.

Planned data will import automatically and the check box will appear preselected if there is any previously planned data present.

Figure 203: Assessment/Plan Element Screen

The user may enter additional information regarding the patient using the **Additional Annotations** free text window. This window offers right-click functionality to import text files, if desired. All information entered using this element appears in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

Usage of the **Back** button is NOT dependent on the element's completion.



**Note:** The **Assessment/Plan** element imports incomplete when it is saved as unfiled data and reloaded. The provider is required to review/edit this element again at this time.

## Patient Education Element

The **Patient Education** element is optional for all seven exam/consult procedures. This element requires one of the two radio buttons to be selected when entering data. The selection of the second radio button opens a text box that requires specific information regarding the patient, or one pre-defined statement may be selected.

Up to twelve pre-defined statements on patient education can be selected from this screen. The check boxes found in the pre-defined statements window have two national patient education statements pre-loaded; however, all twelve statements can be created locally. The DRM Administrator has priority over all other users and may add or delete all twelve, if desired.

The screenshot shows the 'Patient Education' window. The 'Patient Education' section has two radio buttons: 'No barriers to learning identified' and 'Barriers to learning identified:'. Below the second radio button is a large text area for pre-defined statements. There are two checkboxes: 'Reviewed risks/benefits/alternatives associated with the proposed treatment plan. Patient agrees to' and 'Patient received post operative care instructions and verbally indicates understanding.' Below this is an 'Additional Annotations' text box. At the bottom left is an 'Annotations' text box. On the right is an 'Element Requirements Panel' with a blue background, containing text about the element's requirements. At the bottom right are four buttons: '<< Back', 'OK', 'Cancel', and 'Next >>'.

**Figure 204: Patient Education Element Screen**

The user may enter additional information regarding the patient using the **Additional Annotations** free text window. This window offers right-click functionality to import text files, if desired. All information entered in this element displays in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

## Disposition Element

The **Disposition** element is required for all seven exam/consult procedures. This element requires one of the two radio buttons to be selected when entering data. The selection of the second radio button; **Next visit**, requires at least one of the following: one selection of the eight date ranges, a text description about the next visit typed in the text box, or one selection from the pre-defined statements.

The **Next Appointment** check boxes, if selected in the **Tx Planning/Sequencing** screen, import automatically into the **Annotations** view window of this element.

Up to twelve pre-defined statements on disposition can be selected from this screen. The check boxes found in the pre-defined statements window have four national disposition statements pre-loaded; however, all twelve statements can be created locally. The DRM Administrator has priority over all other users and may add or delete all twelve, if desired.

**Disposition**

Disposition

☐ No follow up appointment indicated

☒ Next visit:

☐ within 1 week ☐ 1-2 weeks ☐ 2-4 weeks ☐ 1-2 months

☐ 2-3 months ☐ 3-4 months ☐ 4-6 months ☐ recall

☐ Patient has no eligibility for VA dental benefits and was recommended to the private sector for routine

☐ Patient to return to dental clinic for continuing care.

☐ Patient to be scheduled for continuing development of treatment plan.

☐ Patient provided instructions for obtaining fee dental care subject to VA authorization of proposed tre

Additional Annotations

Annotations

Disposition:

Next Appointment:

(D2150) AMALGAM TWO SURFACES PERMANE: 14(OL) Amalgam

(D2160) AMALGAM THREE SURFACES PERMA: 15(DOL) Amalgam

<< Back OK Cancel

**Element Requirements Panel**

D0150 Comprehensive Disposition requirements:

This element requires one of the two radio buttons to be selected when entering data.

The selection of the second radio button will require at least one of the following:

- one selection of the eight date ranges or
- a text description about the next visit typed in the text box or
- one selection of the pre-defined statements.

Figure 205: Disposition Element Screen

The user may enter additional information regarding the patient using the **Additional Annotations** free text window. This window offers right-click functionality to import text files, if desired. All information entered using this element displays in the read-only **Annotations** window. Once the element is complete, click the **Next** button to move on to the next element.

Usage of the **Back** button is NOT dependent on the element's completion.

## Import Previously Filed Data Screen

When any provider has filed patient dental data previously that day in any of the following screens: **Social History, OHA, TMJ, Parafunctional Habits, Occlusion, Diagnostic Findings, Head and Neck Findings, PSR or Periodontal Chart**; then that data may be imported by a second provider entering an exam encounter that same day. The second provider, after selecting an exam code and then selecting the **Exam** tab, has the following **Import Previously Filed Data** screen display. This screen allows the provider an option to import the data into his exam template to satisfy some possible requirement from the exam code that they may have selected.

Selecting the check boxes of any or all previously filed data that day imports that data into their present exam template session. There is a **Check ALL the above** check box at the bottom of the screen which allows all that day's filed data to be imported into this new TIU progress note. After selecting the desired check boxes, the provider clicks the **Import** button to incorporate this data into the current exam template. If none of the data should be imported into the current exam template, then select the **Cancel** button.

DRM Plus - Import Previously Filed Data

Previously Filed Today

☐ Diagnostic Finding(s) were filed earlier today

☐ Dental Examination:  
Missing Teeth: 17, 32.

☐ Periodontal Charting was filed earlier today

☐ Periodontal Charting  
(This is textual display of periodontal findings. Please see DRM Plus for perio charting graphic.)

TOOTH	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
F (SUP)																
F (BLD)	B	b														
F (FGM)	222	1														
F (PD)	444	433	322	234												
F (MOB)	0.5															
L (PD)	443	222	222	344												
L (FGM)	211															

TOOTH 32



32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17
F (FGM)	1														
F (PD)	444	332	222	222	222	222	222	222	222	222	222	222	222	223	344
L (PD)	443	322	222	222	222	222	222	222	222	222	222	222	222	233	444
L (FGM)	1														

KEY: B=Bleeding b=Delayed Bleeding S=Suppuration  
 +=Pocket/MSJ > 5 +=FGM coronal to CEJ  
 (Periodontal chart may be a composite of entries from different dates)

☐ Check ALL the above.

Import Cancel

Figure 206: Import Previously Filed Data Screen

The  icon on the right of the element button displays when there was previously filed data that day and that data is associated with the element. The **Diagnostic Findings** and **Periodontal Assessment** elements inform the provider, by displaying the  icon, there is one or any combination of diagnostic findings and periodontal charting filed earlier the same day by another provider.

Below the element icons located on the left side of the **Exam** tab, there is a **Review Data Filed Today** button. This button only displays when data was filed to the VistA database earlier that day by all providers working with this patient. This button allows the provider to open the **DRM Plus - Import Previously Filed Data** screen to make edits or corrections by the present provider on what was filed previously that day.

**Note:** The **yellow asterisk** icon will always display when vitals have been entered for this patient during the same day as the **Visit** date (24-hour range) and taken anywhere in the hospital for this patient.

# Completing the Encounter

Click the **Next** button, located in the bottom right corner of both the **Treatment & Exam** and **Periodontal Chart** screens, to begin completing the encounter for the patient's chart. If there was no provider/visit in the banner, imported or selected, when the patient's chart was opened and NOT one completed transaction selected; then clicking the **Next** button brings up the **Provider and Location for Current Activity** screen to select/enter a visit. After selecting one procedure code the **File Data Option** screen displays, which is explained later in this chapter. In addition, if the system finds possible duplicate procedure codes in the planned and completed transactions, the **Potential Duplicate Transaction** screen displays.

To complete the patient encounter with completed transactions:

1. Click the **Next** button on the bottom right corner of either the **Treatment & Exam** or **Periodontal Chart** screens.
2. The **Filing Options** screen displays.

**Dental Record Manager Plus in use by: DRMPROVIDER,DENTIST for DRMPATIENT,ONE**

File Edit Dental Encounter Data Treatment & Exam Tools Reports Help

DRMPATIENT, ONE Jun 12, 17@09:00 DENTAL PDP: DRMPROVIDER,ADMINDE Dental Class  
000-00-0001 01/01/1960 Age: 57 DRMPROVIDER,DENTIST SDP: OPC, Class IV

Cover Page Clinical Record Dental History Chart/Treatment Exam

**Filing Options**

☒ File Data With a Note  
☐ File Data With a Note Addendum  
☐ File Data Without a Note

**Visit Date/Time**  
06/12/2017 09:00

**Encounter Dental Class**  
15-OPC, Class IV

**Suggested Recare Date**  
...

**Disposition**  
☒ Active  
☐ Inactive  
☐ Maintenance

**Primary PCE Diagnosis**  
Select a primary PCE diagnosis by clicking on a checkbox below.

Select	Diagnosis	Diagnosis Description	Send Dx to CPRS Problem List	Additional Information	Procedure	Procedure Description
<input checked="" type="checkbox"/>	K03.0	EXCESSIVE ATTRITION OF TEETH	<input type="checkbox"/>	Tooth: 5	D2750	CROWN PORCELAIN W/ H NOBLE M
<input type="checkbox"/>	K03.0	EXCESSIVE ATTRITION OF TEETH	<input type="checkbox"/>	Tooth: 7	D2750	CROWN PORCELAIN W/ H NOBLE M

**Service Connection**  
☐ Encounter is associated with service-connected condition  
☐ Southwest Asia Conditions (SWAC) ☐ Ionizing Radiation Treatment ☐ Combat Veteran ☐ Agent Orange Treatment ☐ Head and Neck  
☐ Military Sexual Trauma ☐ PR0J T112/SHAD

**Additional Providers**  
Add ... Remove

**Additional Signers**  
Add ... Remove

**Station**  
☒ 500 ☐ 500B2

Return to Chart <= Back Next >= Finish

**Figure 207: Filing Options Screen**

3. Select the correct **Filing Option**. The information displayed in the **Visit Date/Time** and the **Disposition** defaults can be changed. The **Encounter Dental Class** defaults if saved on the cover page. Select the **Suggested Recare Date** if applicable to the patient.
4. Select the **Primary PCE Diagnosis** for the encounter, unless there is only one no matter how many transactions there are which defaults checked.
5. Provider may select the diagnosis to be sent to the patient's problem list.

6. Enter the locally required information for service connection.
7. Enter the optional information of **Additional Providers** or **Additional Signers** on this screen. Please see the **Additional Providers/Additional Signers** section, found later in this chapter, for more information.
8. Select the appropriate facility (station) by clicking the appropriate radio button.
9. Click the **Next** button.
10. The **Dental Class Discrepancy** screen may display, if the encounter dental class and cover page dental class do NOT match. This does NOT stop the user from completing the encounter.
11. The **Set Progress Note Title** screen displays. Note that this screen does NOT display if the **File Data with a Note Addendum** or **File Data Without a Note** options were selected from the filing options.

**Set Progress Note Title**

Selection of TIU Progress Note or Consult Title for Note  
Select title from list or type the desired title.

Typing the desired title in top box will cause the list of titles to change using the letters typed.

Progress Note Title

DENTAL NOTE  
DENTAL CONSULT

ADDENDUM <TOM TEST ADDENDUM>  
ADDICTION <ASI-ADDICTION SEVERITY INDEX>  
ADMISSION <ADMISSION ASSESSMENT>  
ADMISSION ASSESSMENT  
ADIR <ADVANCE DIRECTIVE>  
ADMISSION <ADMISSION ASSESSMENT ONE>  
ADMISSION <ADMISSION ASSESSMENT TWO>  
ADMISSION ASSESSMENT ONE  
ADMISSION ASSESSMENT TWO  
ADMIT <JOY ADMIT (BP)>  
ADMIT <JOY ADMIT>  
ADVANCE <ADVANCE DIRECTIVE NOTIFICATION AND SCREENING>  
ADVANCE <ADVANCE DIRECTIVE>  
ADVANCE DIRECTIVE

OK Cancel

**Figure 208: Set Progress Note Title Screen**

12. Search for the title using the **Progress Note Title** search box or use the scroll bar to select the desired title from the list. Default TIU progress note titles at the top of the list are set in the **CPRS Tools menu → Options submenu → Notes tab → Document Titles** button.
13. Click the **OK** button.
14. The **Progress Note** screen displays.

**Figure 209: Progress Note Screen**

15. Ensure the information on the screen is correct. Enter the electronic signature and click the **Finish** button. An electronic signature is **NOT** required to file the note as unsigned. Please see the **Electronic Signature** section, found later in this chapter, for more information.
16. The **Change Provider** screen displays.

**Figure 210: Change Provider Screen**

17. Click the **Yes** button to change the provider. The **Search for Provider** screen displays. If **No** is selected, an informational screen displays. Please see step 19.





## Potential Duplicate Transactions Screen

If the system detects that a possible duplicate transaction exists, the **Potential Duplicate Transactions** screen displays, when the **Next** button on the **Treatment & Exam** or **Periodontal Chart** screen is clicked.

These transactions have already been **PLANNED** for the following teeth:

Tooth	Surface	ADA Code	Description	Keep	Delete
7	r	D3310	END THXPY ANTERIOR TOOTH	<input type="checkbox"/>	<input checked="" type="checkbox"/>

These are the **COMPLETED** transactions you just entered:

Tooth	Surface	ADA Code	Description	Keep	Delete
7		D2750	CROWN PORCELAIN W/ H NOBLE M	<input checked="" type="checkbox"/>	<input type="checkbox"/>

☐ Keep Planned      ☐ Let Me Decide  
☒ Keep Completed      ☐ Keep All

Process and Go Back      Cancel      Process and Continue

Figure 213: Potential Duplicate Transactions! Screen

A list of planned transactions for a tooth, along with the procedure code and the description, are listed in the top portion of the screen. Completed transactions taken during this encounter are listed in the bottom portion of the screen, with elements that match other planned transactions for the same tooth.


The radio buttons on the lower left portion of the screen designate which information is kept or discarded. The **Keep Planned** radio button keeps the planned transaction and deletes the conflicting completed transaction entered during the encounter. The **Keep Completed** radio button keeps the completed transaction entered during this encounter and deletes the planned transaction. The **Let Me Decide** radio button clears all check boxes and allows for picking and choosing among the planned and completed transactions. **Keep All** allows all conflicting procedure codes that are planned and completed to be processed.

After clicking the desired radio button and selecting, if necessary, which transactions are kept and/or discarded, click the **Process and Go Back**, **Cancel** or **Process and Continue** buttons to continue completing the encounter. **Process and Go Back** processes the procedure codes and returns to the **Treatment & Exam** or **Periodontal Chart** screen. **Cancel** displays the **Periodontal Chart** or **Treatment & Exam** screen without processing any information. **Process and Continue** processes the procedure codes and continues updating the TIU progress note. The **Filing Options** screen displays. Continue completing the encounter from this point, as outlined previously in this chapter.

**Note:** The system does NOT present a user with the **Potential Duplicate Transactions** screen when the potential duplicate is a tooth-related radiographic procedure.

## **File Data Option Screen**

If no procedure code has been entered as completed treatment on the **Completed Care** screen and the user clicks the **Next** button, the **Provider and Locations for Current Activities** screen displays if no visit is present in the banner. When the visit is present or after one has been selected then the **File Data Option** screen displays.



**Figure 214: File Data Option Screen**

Click the appropriate radio button to designate whether the data is to be filed to PCE/DES with a procedure code, or to file to DES-Only data.

### **File to PCE/DES with Code**

To file data to PCE/DES with a procedure code:

1. The **Provider and Location for Current Activities** screen displays. Select the correct visit.
2. Click the File to PCE/DES with code radio button on the File Data Option screen.
3. The program defaults to **D4999** if entering this screen from the **Periodontal Chart**, and **D9999** if entering from **Treatment & Exam**.
4. If the default procedure code is incorrect, click the code drop-down menu and select the desired procedure code.
5. Click the **OK** button.
6. The **Diagnosis Code** screen displays.
7. Select the procedure code on the left side of the screen (default is all selected) and the diagnosis code(s) on the right side. If the correct diagnosis code does NOT appear, use the **Additional Diagnosis Code Search** to find a different diagnosis code.

## File to DES-Only Data

To file to DES-Only data:

1. The **Provider and Location for Current Activities** screen displays. Select the correct visit.
2. Click the **File to DES-Only Data** radio button.
3. If a procedure code is desired (NOT required), check the **Add code** check box and use the drop-down menu to select a procedure code.
4. The **Diagnosis Code** screen displays.
5. Select the procedure code on the left side of the screen (default is all selected) and the diagnosis code(s) on the right side. If the correct diagnosis code does NOT appear, use the **Additional Diagnosis Code Search** to find a different diagnosis code.

## Filing Options Screen

The **Filing Options** screen is divided into ten main sections: **Filing Options**, **Visit Date/Time**, **Encounter Dental Class**, **Disposition**, **Suggested Recare Date**, **Primary PCE Diagnosis & Send DX to CPRS**, **Problem List**, **Service Connection**, **Additional Providers**, **Additional Signers** and **Station**.

### Filing Options

Use the **Filing Options** radio buttons to choose how the encounter is to be filed.




Figure 215: Filing Options

The options are: **File to Data with a Note**, **File Data With a Note Addendum** and **File Data Without a Note**. The **File Data Without a Note** creates no TIU progress note.

### Visit Date/Time

Adjust the **Visit Date/Time** using this function. The program defaults to whatever was entered on the **Provider and Location for Current Activities** screen.




Figure 216: Visit Date/Time

To change the visit date and time:

1. Click the **Ellipsis (...)** button.
2. The **Provider and Location for Current Activities** screen displays.

- Choose the correct provider and appointment from the **Provider and Location for Current Activities** screen and click the **OK** button. Please see the Dental Encounter Data section in the Using DRM Plus Drop-Down Menus chapter of this manual, for further information.

## Encounter Dental Class

Use the drop-down menu to change the **Encounter Dental Class** for this encounter.



The screenshot shows a window titled "Encounter Dental Class" with a drop-down menu. The selected option is "15-OPC, Class IV".

Figure 217: Encounter Dental Class

## Disposition

Use the radio buttons to change the patient's **Disposition** or case management status.



The screenshot shows a window titled "Disposition" with three radio buttons: "Active" (selected), "Inactive", and "Maintenance".

Figure 218: Disposition

## Suggested Recare Date


Use the **Ellipsis (...)** button to add/change the suggested recare date.



The screenshot shows a window titled "Suggested Recare Date" with a text field and an ellipsis button (...).

Figure 219: Suggested Recare Date

Click the **Ellipsis (...)** button to open the **Select Date** screen.



The screenshot shows a window titled "Select Date" with a calendar for May 2017 and a table of month increments. The calendar shows the date 5/31/2017 as today. The table of month increments is as follows:

Select Number of Months from Date Selected	
+3 Months	+12 Months
+4 Months	+18 Months
+6 Months	+21 Months
+9 Months	+24 Months

Buttons for "OK" and "Cancel" are at the bottom.

Figure 220: Select Date Screen

Select the date from the calendar or select the number of additional months from the date selected first on the calendar and then click the **OK** button.

## Primary PCE Diagnosis & Send Dx to CPRS Problem List

This window features a list of the diagnosis codes with descriptions, selections to send diagnosis to CPRS problem list, procedure codes with descriptions, and additional information.

Select the primary PCE diagnosis by clicking the check box under the **Select** column in the **Primary PCE Diagnosis** window. If there is only one diagnosis, or several of the exact same diagnosis, the selection is checked by default. All encounters require only one primary diagnosis filed to VistA PCE.

Select diagnoses to **Send Dx to CPRS Problem List** by clicking the check box in that column.

Select	Diagnosis	Diagnosis Description	Send Dx to CPRS Problem List	Additional Information	Procedure	Procedure Description
<input checked="" type="checkbox"/>	K03.0	EXCESSIVE ATTRITION OF TEETH	<input type="checkbox"/>	Tooth: 5	D2750	CROWN PORCELAIN W/ H NOBLE M
<input type="checkbox"/>	K03.0	EXCESSIVE ATTRITION OF TEETH	<input type="checkbox"/>	Tooth: 7	D2750	CROWN PORCELAIN W/ H NOBLE M

Figure 221: Primary PCE Diagnosis Window

The **Primary PCE Diagnosis** window allows one completed transaction at a time to be deleted when using the **Delete Highlighted Transaction** button.

First highlight the completed transaction that should be deleted and then select the **Delete Highlighted Transaction** button located in the upper right corner of the window.

Figure 222: Transaction Deletion Screen

The resulting pop-up will display asking are sure you want to proceed. Selecting the **Yes** button will delete that completed transaction.

## Service Connection

Use the check boxes to denote service connection, if applicable.

Figure 223: Service Connection

## Additional Providers/Additional Signers

Add or remove providers from this patient encounter by using the tools in this area.

The image shows two side-by-side panels. The left panel is titled 'Additional Providers' and contains a large empty text box, an 'Add ...' button, and a 'Remove' button. The right panel is titled 'Additional Signers' and contains a similar layout with an empty text box, an 'Add ...' button, and a 'Remove' button.

Figure 224: Additional Providers/Additional Signers

To add a provider or signer:

1. Click the **Add** button.
2. The Search for Provider or the Search for Signer screen displays.

The image shows a window titled 'Search for Provider' with a close button (X) in the top right corner. Inside the window, there is a text box with the placeholder text 'Enter search criteria and press Enter to view results below'. Below the text box is a table with four columns: 'Name', 'Initials', 'Title', and 'Service'. The table has 11 rows, all of which are empty. At the bottom of the window, there are two buttons: 'OK' and 'Cancel'.

Name	Initials	Title	Service

Figure 225: Search for Provider Screen


3. Enter the name into the text box and press **<Enter>**.
4. Select the correct name from the results and click the **OK** button.
5. The provider displays on the **Additional Provider/Signer** window. Repeat as necessary.

To remove a provider or signer:

1. Highlight the provider or signer name.
2. Click the **Remove** button.
3. A confirmation screen displays. Click the **Yes** button.
4. The name is removed from the **Additional Provider/Signer** window.

## Station

Select the appropriate **Station** (facility) by clicking the appropriate radio button, if applicable.




Station

☒ 500      ☐ 500BZ

Figure 226: Station

## Progress Note Screen

The **Progress Note** screen has several major areas and functions. This is the final screen in completing an encounter with the patient, and has different incarnations depending on how the transaction is to be filed: **File Data With a Note**, **File Data With a Note Addendum** or **File Data Without a Note**. The following screen displays when the transaction is filed using the **File Data With a Note** option.



Dental Record Manager Plus in use by: DRMPROVIDER,DENTIST for DRMPATIENT,ONE

File Edit Dental Encounter Data Treatment & Exam Tools Reports Help

DRMPATIENT,ONE 000-00-0001 01/01/1960 Age: 57 Jun 12, 17@09:00 DENTAL DRMPROVIDER,DENTIST PDP: DRMPROVIDER,ADMINDE SDP: Dental Class OPC, Class IV

Cover Page Clinical Record Dental History Chart/Treatment Exam

TIU Note Title: DENTAL NOTE

Note Date/Time: Monday, June 12, 2017 11:53:59 AM

Diagnostic Findings  
Added  
Modified  
Deleted

Planned Items  
Added  
Modified  
Deleted

Sequencing Notes

Periodontal Exam  
PSR Exam  
Modified

Head/Neck Findings  
Added  
Modified  
Deleted

Tooth Notes  
Code Boilerplate  
Note Boilerplate

Completed Items  
Added  
Modified  
Deleted

Next Appointment

Dental Alerts

Social History

OHA Findings

Occlusal Findings

Parafunctional Habits

TMJ Findings

Adverse Events

Device Tracking

Insert Text Marker

View CPRS Templates

CNT Navigator

Electronic Signature Code

Return to Chart <= Back Next >= Finish

Patient Name: DRMPATIENT,ONE, DOB: 01/01/1960, Age: 57  
Visit: 06/12/2017 09:00  
Primary PCE Diagnosis: K03.0 (EXCESSIVE ATTRITION OF TEETH)  
Dental Category: 15-OPC, Class IV. Treatment Status: Active.

Completed Care:  
(D2750) CROWN PORCELAIN W/ H NOBLE M. Tooth: 5.  
DK: K03.0 Excessive Attrition of Teeth  
(D2750) CROWN PORCELAIN W/ H NOBLE M. Tooth: 7.  
DK: K03.0 Excessive Attrition of Teeth

Periodontal Examination:  
Periodontal Charting  
(This is textual display of periodontal findings. Please see DRM Plus for perio charting graphic.)

TOOTH	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
F SUP	B	b			B	b									S	
F BLD	222	1			1					1		1			b	B
F FGM	444	433	322	234			432	222	222	223		322	222	234	444	
F FID	444	433	322	234			432	222	222	223		322	222	234	444	
L FID	443	222	222	344			432	222	222	334		443	222	333	344	
L FGM	211				1										112	
TOOTH	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28
F SUP	222	221	11													
F BLD	444	443	332	222	222	222	222	222	222	222	222	222	222	222	222	222
F FGM	444	432	222	222	222	222	222	222	222	222	222	222	222	222	222	222
L FID	444	432	222	222	222	222	222	222	222	222	222	222	222	222	222	222
L FGM	222	111	1													
L BLD	b	b														

KEY: B=Bleeding b=Delayed Bleeding S=Suppuration  
\*\*Pockets/MDJ > 9 \*\*FGM coronal to CEJ  
(Periodontal chart may be a composite of entries from different dates)

Tooth and Patient Notes:  
Patient Note

Figure 227: Progress Note Screen

Use the **Ellipsis (...)** button to change the TIU progress note title here if needed. Use the **up** and **down arrows** to change the **Note Date/Time** if desired. From this screen, the user can also: view and/or import DRM Plus objects, view and import CPRS templates, launch the **CNT Navigator** (DSS clinical note templates), view the patient's TIU progress note, import VistA medical information, and add an electronic signature.

In addition to the imports that can be made into a TIU progress note, free texting for narrative may also be entered below or above the **Text Marker** (dashed line) prior to filing the note by the user.

**Note:** All information added below the **Text Marker** may be saved as unfiled data when entered. All information edited or added above the **Text Marker** may not be saved as unfiled data when entered.

## Viewing/Importing DRM Object/Progress Note

To view/import DRM Plus objects:

1. Click the **View DRM Plus Objects** button (defaulted) below the tree on the left-side of the screen.
2. Double-click the DRM Plus object and it displays in the viewer on the right-side of the screen. The DRM Plus object imports where the cursor is positioned.
3. Information can be added or deleted by typing directly into the progress note on the right-side of the screen.

## Viewing/Importing CPRS Templates


To view/import CPRS templates:

1. Click the **View CPRS Template** button below the tree on the left-side of the screen.
2. The shared/personal templates tree, if expanded, displays on the left-side of the screen. Most functionality of CPRS templates in DRM Plus is the same as in CPRS.
3. Click the **View DRM Plus Objects** button to return to the DRM Plus object tree.

## Importing VistA Medical Information

To import VistA medical information into the patient's TIU progress note:

1. Right-click in the progress note area.
2. The **Import Menu** displays.



**Figure 228: Import Menu Screen**

3. Choose the information that is to be imported from the top half of the menu.
4. The information displays in the progress note where the cursor was positioned.



## Other Options in the Import Menu

- **Hide Note Objects:** Hides or closes the objects tree on the left side of the screen and enlarges the progress note viewer.
- **View Note Objects:** Undoes the hide note objects.
- **Copy:** Copies selected text in the progress note.
- **Cut:** Cuts selected information in the progress note.
- **Paste:** Pastes information into the progress note.
- **Select All:** Selects all text in the progress note.
- **Print Note:** Prints the progress note.
- **Import Text File to Note:** Navigates to a text file to import into a note. See the User Options section in the Using DRM Plus Drop-Down Menus chapter of this manual for more information on automatically setting the location for this text file.

## Accessing Dental CNTs

Click the **CNT Navigator** button to access dental CNTs. These are DSS product clinical note templates. The dental CNTs may NOT be mapped for DRM Plus at the user's site; it would require IT assistance. Please see the CNT Navigator section, further in this chapter.

## Electronic Signature

Enter the **Electronic Signature Code** and click the **Finish** button to complete the TIU progress note. Clicking the **Finish** button without entering an electronic signature leaves the patient TIU progress note status as unsigned.

## Progress Note Addendum

When the **File Data With Note Addendum** option has been chosen, the **Progress Note** screen displays slightly different. The functions are the same as in the previously shown progress note screen with the additions of selecting the **Note Categories** radio button options and with the parent note displaying in the upper right window.

**Dental Record Manager Plus in use by: DRMPROVIDER,ADMINDENTIST for DRMPATIENT,ONE**

File Edit Dental Encounter Data Treatment & Exam Tools Reports Help

DRMPATIENT,ONE 000-00-0001 01/01/1960 Age: 57 V: May 31, 2017@09:00 DENTAL DRMPROVIDER,ADMINDENTIST PDP: DRMPROVIDER,ADMINDENTIST SDP: Dental Class OPC, Class IV

Cover Page Clinical Record Dental History Chart/Treatment Exam

Note Categories:  
☒ My Signed  
☐ All Completed  
☐ By Title  
☐ By Location  
☐ By Author

05/31/2017 DENTAL NOTE  
 05/31/2017 DENTAL NOTE  
 05/26/2017 DENTAL NOTE  
 05/25/2017 DENTAL NOTE  
 05/05/2017 DENTAL NOTE  
 04/28/2017 DENTAL NOTE  
 03/22/2017 DENTAL NOTE  
 03/17/2017 DENTAL NOTE

LOCAL TITLE: DENTAL NOTE  
 STANDARD TITLE: DENTISTRY NOTE  
 DATE OF NOTE: MAY 31, 2017@09:22 ENTRY DATE: MAY 31, 2017@09:22:27  
 AUTHOR: DRMPROVIDER,ADMINDENTIST EXP COSIGNER:  
 URGENCY: STATUS: COMPLETED

Patient Name: DRMPATIENT,ONE, DOB: 01/01/1960, Age: 57  
 Visit: 05/31/2017 09:00  
 Primary PCE Diagnosis: K08.530 (FRACTURED DENTAL RESTORATIVE MATERIAL WITHOUT LOSS OF MATERIAL).  
 Dental Category: 15-OPC, Class IV. Treatment Status: Active.

Diagnostic Findings  
 Added  
 Modified  
 Deleted

Planned Items  
 Added  
 Modified  
 Deleted

Sequencing Notes  
 Periodontal Exam  
 PSR Exam  
 Modified  
 Head/Neck Findings  
 Added  
 Modified  
 Deleted  
 Tooth Notes  
 Code Boilerplate  
 Note Boilerplate  
 Completed Items  
 Added  
 Modified  
 Deleted  
 Next Appointment  
 Dental Alerts  
 Social History  
 OHA Findings  
 Occlusal Findings  
 Parafunctional Habits  
 TMJ Findings  
 Adverse Events

Patient Name: DRMPATIENT,ONE, DOB: 01/01/1960, Age: 57  
 Visit: V: May 31, 2017@09:00 DENTAL. Primary PCE Diagnosis: K08.530.  
 Dental Category: 15-OPC, Class IV. Treatment Status: Active.

Completed Care:  
 (D2160) AMALGAM THREE SURFACES PERMA. Tooth: 31. Surface(s): DOL.  
 DX: K02.53 Dental Caries on Pit and Fissure Surface Penetrating into Pulp  
 (D2150) AMALGAM TWO SURFACES PERMANE. Tooth: 30. Surface(s): OL.  
 DX: K02.53 Dental Caries on Pit and Fissure Surface Penetrating into Pulp

Dental Alerts:  
 Antibiotic Premed

Electronic Signature: \_\_\_\_\_

Return to Chart <= Back Next >= Finish


**Figure 229: Progress Note Addendum Screen**

The notes associated with the selected **Note Categories** appear in the narrow window below the note categories area. Use the drop-down menu below the note categories to filter the results. Click a parent note to view details in the upper right window of the screen.

Finish the process of filing the encounter with a note addendum, which displays in the lower right window of the screen.


## CNT Navigator

Clicking the **CNT Navigator** button displays up a directory of clinical note templates in use, as a tool, to assist the user in writing a note, or additional information to a note. The **CNT Navigator** button is located on the bottom left side of the **Progress Note** screen. These are Document Storage Systems (DSS) product clinical note templates, and may NOT function unless they are mapped correctly, which requires IT assistance.



**Figure 230: CNT Navigator Button**

Using the mouse, point and click to select pre-determined text in the development of a note. Type the text within the CNT if the pre-determined text or statement does NOT contain the necessary verbiage.



**Figure 231: Clinical Note Template Navigator Screen**

To access a CNT:

1. Click any of the tabs.
2. A listing of CNTs specific to the selected tab displays. Either double-click the desired template, or click once to select it and then click the **Run** button.

## Navigating Within CNTs

To navigate within a CNT:

1. Point and click within the windows, tabs, drop-down arrows, check boxes, and radio buttons. Each navigational method provides the user with a different method of entering or selecting information.

The screenshot shows a software window titled "Dental Hygiene Progress Note" with a menu bar (File, Tools, Help) and standard window controls. The form is divided into several sections:

- Vital Signs:** Contains input fields for Temp., Pulse, Systolic, Height, Resp., Diastolic, BMI, Age, Weight, Pain, SP02%, and BSA. Each field has a date/time label and a small up/down arrow icon.
- Patient Sex:** Radio buttons for Male and Female.
- Chief Complaint:** Includes buttons for "Click for Normal", "Select from Pick List", and "Clear Text". A list box on the right shows a scrollable list of options: Headache, Pain to heat, Pain to cold, Pain to pressure, Toothache, Swelling, Pain with Postural Changes, Teeth Sensitive to Chewing, and Minimal to No Relief from OTC Pain Medication.
- Pre-Medication:** A dropdown menu and a checkbox.
- ASA:** Radio buttons for levels 1, 2, 3, and 4.
- Footer:** A row of buttons: Open, Save, Reminders, Preview, Cancel, Finish, Prev, 1 of 6, and Next.

Figure 232: Navigating Within CNT

2. Preview the note by clicking the **Preview** button.
3. Click the **Return** button to continue writing the note.

**Preview Note (Not Connected To VISTA).** [X]

MEDICATIONS:  
Ibuprofen

ALLERGIES:  
Sulfa

ORAL HYGIENE:  
Fair.

CALCULUS ASSESSMENT:  
Slight.

ORAL ASSESSMENT:  
Caries.

ANESTHETIC:  
1.8 ml. of Lidocaine 2% with epinephrine 1:100,000.

REACTION TO TREATMENT:  
Fairly well.

RECAP PLAN:  
6 months.

TREATMENT NEEDED:  
Oral Hygiene

Accept Note

Return To Form

Integration:

☒ Send to Clipboard.


☐ Send to TIU/PCE.

**Figure 233: Preview Note**

When the note is complete, click the **Finish** button. Then click the **Accept Note** button.

## Consult Notes

The option to complete a consult displays after the **Filing Options** screen, during the process of completing a patient encounter, when the **Set Progress Note Title** screen displays.



The screenshot shows a software window titled "Set Progress Note Title". Inside the window, there is a header section with the text "Selection of TIU Progress Note or Consult Title for Note" and a sub-instruction "Select title from list or type the desired title." Below this, a note states: "Typing the desired title in top box will cause the list of titles to change using the letters typed." The main area of the window is a list box labeled "Progress Note Title" at the top. The list box contains the following items: "DENTAL NOTE", "DENTAL CONSULT", "ADDENDUM <TOM TEST ADDENDUM>", "ADDICTION <ASI-ADDICTION SEVERITY INDEX>", "ADMISSION <ADMISSION ASSESSMENT>", "ADMISSION ASSESSMENT", "ADIR <ADVANCE DIRECTIVE>", "ADMISSION <ADMISSION ASSESSMENT ONE>", "ADMISSION <ADMISSION ASSESSMENT TWO>", "ADMISSION ASSESSMENT ONE", "ADMISSION ASSESSMENT TWO", "ADMIT <JOY ADMIT (BP)>", "ADMIT <JOY ADMIT>", "ADVANCE <ADVANCE DIRECTIVE NOTIFICATION AND SCREENING>", "ADVANCE <ADVANCE DIRECTIVE>", and "ADVANCE DIRECTIVE". At the bottom right of the window are two buttons: "OK" and "Cancel".

**Figure 234: Set Progress Note Title Screen**

To complete the consult:

1. Choose the consult title from the **Set Progress Note Title** screen.
2. The patient's pending consults display on the screen.

**Set Progress Note Title**

Selection of TIU Progress Note or Consult Title for Note  
Select title from list or type the desired title.

Typing the desired title in top box will cause the list of titles to change using the letters typed.

Progress Note Title

DENTAL NOTE  
**DENTAL CONSULT**  
 ADDENDUM <TOM TEST ADDENDUM>  
 ADDICTION <ASI-ADDICTION SEVERITY INDEX>  
 ADMISSION <ADMISSION ASSESSMENT>

The above title must be associated with a consult request. Select one of the following or chose a different title.

Consult Request	Date	Service	Procedure	Status	# Notes
07/17/2014	14.24	DENTAL		PENDING	0

OK Cancel

**Figure 235: Set Progress Note Title with Pending Consult**


3. Select the consult from the list.
4. The consult is added to the progress note on the **Progress Note** screen. Once the electronic signature is entered on the **Progress Note** screen, the consult is complete.

## Resident Filing as Cosigners or Distributed Providers

A 2006 VA Directive stated that residents are users with a VistA **Person Class** of **V030300**, V11550 or V115600. All residents are required to have a distributed provider (attending) to complete the encounter with the patient. Since most sites require residents to enter a cosigner for the note, the cosigner defaults as the distributed provider in PCE. If there is no cosigner required for the resident, or the user filing data to a resident, they must enter the distributed provider before filing.

To add a distributed provider:


1. Click the **Finish** button.
2. The **Change Provider** screen displays.



The 'Change Provider' dialog box has a blue title bar. The main text asks 'DO YOU WANT TO FILE THIS ENCOUNTER UNDER ANOTHER PROVIDERS NAME?'. Below this is a checkbox labeled 'Print New Signed Progress Note' with a subtext '( Only if you are the provider and have a default printer )'. At the bottom are three buttons: 'Yes', 'No', and 'Cancel'.

Figure 236: Change Provider Screen

3. Click the desired response. (Residents click the **No** button).
4. The **Search for Cosigner Provider** screen or the **Search for Distributed (PCE Primary) Provider** screen displays.



The 'Search for Provider to be the DISTRIBUTED (PCE Primary)' dialog box has a blue title bar. It contains a text input field with the prompt 'Enter search criteria and press Enter to view results below'. Below the input field is a table with four columns: 'Name', 'Initials', 'Title', and 'Service'. The table has 11 empty rows. At the bottom right are 'OK' and 'Cancel' buttons.

Name	Initials	Title	Service

Figure 237: Search for Distributed Provider Screen



5. Enter the search terms in the text box and press the <Enter> key.
6. Choose the provider from the search results and click the **OK** button.
7. An information screen displays confirming that the information was filed.

**Note:** VA dental gives credit/RVU time to the resident actually performing the procedures, all of which is filed to DES; to meet the VA requirement that the attending (distributed) gains credit for the encounter, this has to be filed in PCE.

When the encounter is filed to PCE, the resident becomes the secondary provider, and the distributed provider becomes the primary provider for the PCE encounter. All procedures and diagnoses are assigned to the distributed provider in VistA PCE.

**Cache TRM:3816 (CACHE)**

File Edit Help

PCE Update Encounter May 31, 2017@10:31:42 Page: 1 of 2

DRMPATIENT, ONE 000-00-0001 Clinic: DENTAL

Encounter Date 5/31/2017 10:00 Clinic Stop: 180 DENTAL

1 Encounter Date and Time: MAY 31, 2017@10:00

2 Provider: DRMPROVIDER, ADMINDENTIST PRIMARY Dental/Dentist/Oral & Maxill

3 Provider: DRMPROVIDER, RESIDENT TWO Dental/Dental Resident

4 ICD-10 Code or Diagnosis: K02.52 Dental caries on pit and fissure surface penetrating into dentin

Provider Narrative: DENTAL CARIES ON PIT AND FISSURE SURFACE PENETRATING INTO DENTIN

Primary/Secondary ICD-10 Diagnosis for the Encounter: PRIMARY

5 CPT Code: D2160 AMALGAM THREE SURFACES PERMA

Primary ICD-10 Diagnosis: K02.52 Dental caries on pit and fissure surface penetrating into dentin

+ + Next Screen - Prev Screen ?? More Actions

ED Edit an Item TR Treatment DD Display Detail

DE Delete an Item IM Immunization DB Display Brief

EN Encounter PE Patient Ed IN Check Out Interview

PR Provider ST Skin Test CR Contra/Refusal Event

DX Diagnosis (ICD) XA Exam QU Quit

CP CPT (Procedure) HF Health Factors

Select Action: Next Screen//

**Figure 238: PCE Encounter Information in VistA**

## Appendix A – Glossary of VA Terms

ADPAC	Automated Data Processing Applications Coordinator
AICS	Automated Information Collection System: formerly Integrated Billing, the program that manages the definition, scanning and tracking of Encounter Forms.
APPOINTMENT	A scheduled meeting with a provider at a clinic; an appointment can include several encounters involving other providers, tests, procedures, etc.
CBOC	Community Based Outpatient Clinic
CC	Coordinating Committee
CFO	Chief Financial Officer
CHECKOUT PROCESS	Part of the Medical Administration (PIMS) appointment processing. The checkout process documents administrative and clinical data related to the appointment.
CIO	Chief Information Officer
CIR	Corporate Information Repository
CIRN	Clinical Information Resource Network
CLINICIAN	A doctor or other provider in the medical center authorized to provide patient care.
CNT	Clinical Note Template (Used to format TIU Progress Notes)
CPR	Cardiopulmonary Resuscitation
CPRS	Computerized Patient Record System
CPT	Common Procedure Terminology
CQI	Continuous Quality Improvement
DAS	Dental Activity System (also called AMIS)
DES	Dental Encounter System (also called DES)
DHCP	Decentralized Hospital Computer Program (See: VistA)
DNR	Do Not Resuscitate
DOD	Department of Defense
DRG	Diagnostic Related Group
DSS	Decision Support System
DSS	Document Storage Systems, Inc.

DVA	Department of Veterans Affairs
EDI	Electronic Data Interchange
ELC	Executive Leadership Council
EMR	Electronic Medical Record
ENCOUNTER	A contact between a patient and a provider who has responsibility for assessing and treating the patient at a given contact, exercising independent judgment. A patient can have multiple encounters per visit.
ENCOUNTER FORM	A paper form used to display and collect data pertaining to an outpatient encounter, developed by the AICS package.
EPISODE OF CARE	Many encounters for the same problem constitute an episode of care. An outpatient episode of care may be a single encounter or can encompass multiple encounters over a long period of time.
FEMA	Federal Emergency Management Agency
FIM	Federal Independence Measure
FRP	Federal Response Plan
GAF	Global Assessment of Functioning
GPRA	Government Performance and Results Act
GUI	Graphic User Interface
HEALTH SUMMARY	A Health Summary is a clinically oriented, structured report that extracts multiple kinds of data from VistA and displays it in a standard format.
HR IGA	Human Resources Office of Intergovernmental Affairs
INPATIENT VISIT	Inpatient encounters include the admission of a patient to a VAMC and any clinically significant change related to treatment of that patient.
IOM	Institute of Medicine
ISDA	Intensity Severity Admission Discharge (criteria)
IT	Information Technology
JCAHO	Joint Commission on the Accreditation of Healthcare Organizations
LAN	Local Area Network
MVV	Mission Vision Values
MAC	Management Assistance Council
MCCR	Medical Care Cost Recovery
MDS	Minimum Data Set

NHCU	Nursing Home Care Unit
OERR	Order Entry Results Reported
OSHA	Occupational Safety & Health Administration
OUTPATIENT ENCOUNTER	Outpatient encounters include scheduled appointments and walk-in unscheduled visits
OUTPATIENT VISIT	The visit of an outpatient to one or more units or facilities located in or directed by the provider maintaining the outpatient health care services (clinic, physician's office, hospital medical center) within one calendar day.
PACS	Picture Archiving and Communications System
PAI	Patient Assessment Instruction
PCE	Patient Care Encounter
PI	Prevention Index
PM	Performance Management
PROCEDURE	A test or action done for or to a patient that can be coded with the CPT coding process.
PROVIDER	The entity which furnishes health care to a consumer.
PSA	Patient Service Area
PTSD	Post-Traumatic Stress Disorder
RM	Risk Management
RPM	Resource Planning Methodology
SD	Standard Deviation
SMI	Seriously Mentally Ill
SSC	Shared Service Center
TIU	Text Integrated Utility
TQI	Total Quality Management
UM	Utilization Management
UNC	Universal Naming Convention. Used in place of Drive letters.
VA	Department of Veteran Affairs
VAVS	Veterans Administration Voluntary Service
VHA	Veterans Healthcare Administration
VISN	Veterans Integrated Service Network
VISIT	The visit of a patient to one or more units of a facility within one calendar day.

VistA	Veterans Information System Technology Architecture, the new name for DHCP.
VSO	Veterans Service Organization

## Appendix B – Common Application Functions

Microsoft Windows tools are used in DRM Plus. Left-clicking, right-clicking, double-clicking, drop-down arrows, radio buttons, check boxes, text boxes, highlighting and scroll bars are used throughout DRM Plus. Certain buttons and clicking options are common to most screens and are discussed below.

**OK:** Clicking the **OK** button is used to finalize a selection or end a process. The open screen is closed, and the user is moved to another screen.

**Cancel:** Clicking the **Cancel** button cancels the action taken on a screen and returns the user to the previous screen.

**Next:** Clicking the **Next** button moves the user to the next screen.

**Back:** Clicking the **Back** button moves the user to the previous screen.

**Add:** Clicking the **Add** button adds a selected item to a function.

**Edit:** Clicking the **Edit** button allows the user to edit a selection.

**Delete:** Clicking the **Delete** button allows the user to delete a selection.

**Reset:** Clicking the **Reset** button resets changed settings to their original settings during this session.

**Finish:** Same as the **OK** button.

**Browse Buttons:** Clicking the **Browse** button moves the user to a previously programmed selection screen.

**Radio Buttons:** Clicking a **radio button** displays a dot in the button, designating a specific option. Only one radio button is allowed for a section in a group.

**Check boxes:** Clicking a **check box** works the same as a radio button, however multiple selections may be added from one group.

**Text Boxes:** Clicking in a text box allows the user to type text into the box.

**Drop-Down Arrows:** Clicking these **arrows** displays a menu of selections.

**Selection Arrows:** Clicking these **arrows** allows a selected item to be moved from one dialogue box to another.

**Search Boxes:** Typing selection criteria in a search box causes the criteria to be matched to a master file. Matches are displayed, allowing the user to highlight the desired selection for further action. DRM Plus requires the user to press the <Enter> key after entering the criteria.

**Sorting:** Clicking a transaction table column heading sorts the table, usually in ascending order, depending on the current view. Clicking the column heading a second time returns the table to its original view.

**Highlighting:** Clicking an item results in its being highlighted and selected for the next action to be completed.

**Shift Key:** Generally, holding the <Shift> key down allows for selection of multiple consecutive items in Windows applications.

**Control Key:** Generally, holding the <Ctrl> key down allows for selection of multiple items in Windows applications.

**Keyboard Use:** When a letter or a button name is underlined (Add or Speed Code) the keyboard can be used to activate the button. The action required is to press and hold the <Alt> key, then press the underlined letter.

## Appendix C – Hints and Notes

### Save Unfiled Data

If the user has entered any data and attempts to close DRM Plus, or switch to a different patient, DRM Plus displays a screen prompting the user to save the current patient's entries. Clicking the **Yes** button initiates the save unfiled data function.

When the user has saved unfiled data and no longer needs this data for the patient, there are two options to eliminate this saved, unfiled data. The first is to select the **Delete** button when opening the patient's chart; a dialogue screen asks the user if they want to load the saved data. The second option is to delete the saved unfiled data from the **Unfiled Data by Provider** report.

### Dental Class Displayed on Banner

The patient's **Dental Class** displays in the banner area, only if the information was entered on the cover page in the **Dental Eligibility/Dental Class** field by a DRM Plus Administrator. As soon as the **Dental Class** is selected, it displays in both fields: the dental class field on the **Cover Page** tab screen, and the **Dental Class** box in the DRM Plus banner.

### Diagnostic Findings

**Diagnostic Findings** are NOT updated automatically from **Completed Care** entries for any encounter. Any **Completed Care** entries that are filed need to be entered as **Diagnostic Findings** during a patient's future dental examination.

Always mark teeth missing in **Diagnostic Findings** before entering **Partials, Dentures, Implants** or **Bridge** findings. DRM Plus works best when missing teeth or edentulous arch(es) are entered before any other findings. **Dentition** is always entered first on a new patient before any dental data is entered in DRM Plus, if the patient is a juvenile.

**Diagnostic Findings** may be deleted after the encounter has been filed by a user. If this happens after the encounter has been filed, the deleted findings are removed from the graphic; however, the text entry remains in the transaction table with a line through it. DRM Plus Administrators can completely remove any **Diagnostic Finding** entry from the transaction table, unless the entry was already deleted by a user.

### Treatment Plan

For implant procedures entered in the **Treatment Plan** screen, and if a related finding of **Missing** has NOT been entered, DRM Plus does NOT allow this corrective planned treatment procedure to be entered.

It is recommended to use the **Include Findings and Completed** button to temporarily combine screens of the planned treatment with the findings and NOT use the automatic **Include** parameter in DRM Plus. The automatic **Include** parameter is the original default in DRM Plus when the user is viewing the **Treatment Plan** screen. End-users may edit this parameter by accessing the **Treatment & Exam** menu → **Show Configuration** submenu → **Tx & Exam** tab → **Display Defaults** drop-down menu → remove check marks in check boxes.



## **Multi-Add Screen**

Small buttons [<] and [>] have been added to the **Procedure Code** selection screen to enable the user to move the screen to the other half of the graphic chart. This may be necessary to see what is beneath, especially when entering multi-add codes.

Missing teeth display as white text on a blue background. This can help the user visualize the mouth while entering multi-add codes. The missing teeth are still selectable if needed for the procedures partial, denture, implant or bridge.

## **Ranged Codes**

Certain codes, designated as ranged codes, used for multiple teeth procedures, mostly in the prosthetics area, have been restricted and can only be used with hard coded DRM Plus icons. These icons are **Partial, Bridge, Conn Bar** and some **Denture** procedure codes. These icons are only found with the **Treatment Plan** and **Completed Care** screens. When other options are utilized for selecting these procedure codes, when using the **ADA Codes** icon or **Add** box, the ranged codes are disabled and cannot be selected.

## **Speed Codes**

**Speed Codes** do NOT have multi-add or suggestion links functionality. The multi-add and suggestion links functionalities work with the **Quick Codes** icon, **ADA Codes** icon, **CPT Codes** icon, using the **Add** box to select a procedure code. Some of the DRM Plus standard icons allow multi-add functionality. Speed codes are only used with the **Treatment Plan** and **Completed Care** view screens for selecting a procedure code.

**Speed Codes** that include codes violating coding compliance rules need to be edited or deleted.

## **Tx Planning/Sequencing Screen**

There is a maximum of nine phases which may be added to the sequenced planned treatment for one patient. There are an unlimited number of sub-phases possible in each phase. There is only one **Non-VA Care** phase that may be added however with an unlimited number of sub-phases.

Ranged codes, mostly prosthetics, move as a block when one code is highlighted, by left-clicking the item, holding and dragging to the proper phase.

End-users can re-size and move the **Tx Planning/Sequencing** screen in all directions. DRM Plus allows users to drag and drop multiple procedures at one time by pressing the <Shift> and <Ctrl> keys.

Users may add or modify the **Treatment Plan** and **Tx Planning/Sequencing** screen and file the changes without having to click the **Next** button and create a TIU progress note as long as no new completed transactions, perio, PSR or Head & Neck data are entered. In addition, the most recent dental encounter must have an **Active** status for this feature to work. Clicking the **Save & Exit** button from the **Tx Planning/Sequencing** screen files any changes made and minimizes DRM Plus. Any new planned entries added have the same visit date as the latest active TIU progress note filed for the patient.

If the completed care screen is set as the display default in DRM Plus, transactions completed from the sequencing screen may require a refresh for those transactions to display in the **Completed Care** screen. Selecting the **Diagnostic Findings** button or the **Treatment Plan** button and then reselecting the **Completed Care** button refreshes the application.

## **Completed Care**

Completed procedures that have been filed can only be deleted by a DRM Plus Administrator. If an end-user, who does NOT have this administrative parameter option, highlights an entry in the **Completed Care** transaction table and clicks the **Delete** button, a screen stating that the transaction cannot be deleted displays. When a DRM Plus Administrator deletes an entry from DRM Plus, appropriate procedures must be followed to correct any associated entry filed in VistA TIU or possible VistA PCE.

In those rare cases where pre-existing charted care does NOT allow the user to click and select a procedure code via the graphic chart, the user may enter the desired code via the **ADA Codes** icon.

If one or more of the procedures entered with multi-add functionality require an assignment of a different diagnosis code, deselect the procedure code by clicking the check box of the procedure description located on the upper left side of the **Diagnosis Code** screen. The deselected procedure code no longer appears highlighted and is removed from the table at the bottom of the screen, indicating that it does NOT associate with the selected diagnosis code. Click the **OK** button at the bottom of the screen, and the remaining unchecked procedure code(s) recycle and appear highlighted in the **Diagnosis Code** screen for the user to assign a diagnosis code. This process of recycling continues until all procedure codes added with the multi-add functionality have a diagnosis code assigned by the user.

## **Periodontal Chart**

The **History** button maintains graphical entries from previous perio examinations; therefore, any prior perio exam graph may be viewed with the **History** button, selected from the **Periodontal Chart** screen. The **Periodontal History/Compare** screen may be vertically extended to view any data that is NOT visible.

To use the **Furcation** icon, the cursor shield must be positioned at the root location where an entry in the graphics would be appropriate.

The **X** button located at the end of the pre-defined measurement scale results in a null entry in the transaction table for **Pocket, FGM, MGJ, Mobility** and **Furcation** entries. This null entry only works when the specific icon is active for **Pocket, FGM, MGJ, Mobility** and/or **Furcation**. The null entry remains as a “-” mark in the transaction table, and there is no entry in the TIU progress note.

If the error is recognized immediately, the **Undo** button may be clicked. Otherwise, place the cursor shield on the pocket where the incorrect value was entered. If the value is for **Pocket, FGM, MGJ, Mobility** or **Furcation**; the incorrect value can be replaced by entering the correct value.

A zero entry results in no graphical view; however, it results in a zero entry in the transaction table and TIU progress note, since it is a measurement. If no recording should be present for a given icon, a null entry can be created by clicking the **X** button. If the incorrect entry is for **Bleeding, Delayed Bleeding** or **Suppuration**; click the identical icon again to remove the graphic display and transaction table entry.

When perio data is imported into the TIU progress note, it has each tooth displayed with each surface and condition shown in the vertical column, under the tooth number. The key at the bottom of the perio data explains certain symbols. There may be a statement at the bottom of the perio data, which informs other providers that this TIU progress note contains perio data from the current exam, as well as data that has NOT changed from at least one previous exam. If the **Clear** icon was used at the beginning of the perio exam, then only data from this current exam is imported into the note.

Warning levels can be changed by the end-user and then displayed. The pocket depth warning level should be the same for the perio chart graphics, and for the pocket depth warning level listed on the **Statistics** tab. The pocket depth warning level on the **Statistics** tab must be the same as the pocket depth warning level in the **Periodontal** tab. Both tabs are found by utilizing the **Treatment & Exam** menu → **Show Configurations** submenu.

## **Completing the Encounter**

Selecting any **Service Connection** check box sends a flag to PCE for the encounter.

DRM Plus is now aligned with CPRS for patients who are **Combat Veteran** service connected. If appropriate, the Combat Veteran option defaults to a check (**Yes**) in DRM Plus. To remove the check (change to **No**), click the check box to the left of the **Combat Veteran** field. The check is removed.

VistA has co-signature functionality and is checked by CPRS and DRM Plus. Both GUIs also have additional signer functionality. Additional signers are NOT required but may be added to a TIU progress note by a provider. Do NOT confuse additional signers with co-signers. Co-signers are built into VistA by facility management based on business rules. If the software detects that a co-signer is required, a screen displays, requesting a co-signature. A provider may need a co-signer for one or all TIU progress notes.

Since most sites require residents to enter a co-signer, the co-signer defaults as the distributed provider in PCE. If there is no co-signer (i.e., no note exists when using **File Data Without a Note**, or the resident is NOT required to have a co-signer) then the resident, or user filing resident data, must enter the distributed provider prior to filing. When the encounter is filed to VistA PCE, the resident becomes the secondary provider and the distributed (attending) provider becomes the primary provider for the encounter. All procedures and diagnoses are assigned to the distributed provider in VistA PCE.

VA Dental wants to give credit/RVU time, etc. to residents who perform the procedures, and all of this is filed to VistA DES. However, to meet the VA requirement that gives credit for the encounter to the attending (distributed), requires this to be filed in VistA PCE.

DRM Plus users may import a **TIU Note Boilerplate** into a patient's TIU progress note. If the **Note Boilerplate** parameter is set with a check mark, an informational screen displays when selecting a TIU progress note or consult title, if that title has a note boilerplate associated with it. The informational screen allows a user to select **Yes** or **No** to the question of importing the note boilerplate. If this parameter is NOT set with a checkmark, then the note boilerplate imports into the patient's TIU progress note without an opportunity to decline this action. This parameter is located under the **Tools** menu → **User Options** submenu → **General** tab → **Other Parameters** button.

CPRS templates automatically import into a patient's TIU progress note if the TIU note/consult title selected is associated with a CPRS template and there is no option for the user to decline this template import. When the template displays or opens, complete or fill-in the appropriate information on the template and close or finish it. The information entered on the template imports into the patient's TIU progress note. **Please note**, DRM Plus does NOT support **Reminder Dialog** or **COM objects** CPRS templates.

There are generally two types of TIU progress notes created using DRM Plus: 1) using the **Exam** tab or 2) NOT using the **Exam** tab. A TIU progress note created using the **Exam** tab sequences DRM Plus objects in the order designed and approved by the VA Dental Exam Committee. The sequencing of DRM Plus objects in the **Note Objects Sequence** parameter screen is overwritten by the **Exam** tab sequence design when the **Exam** tab is used. TIU progress notes created without using the **Exam** tab sequence DRM Plus objects in the order set by the user in the **Note Objects Sequence** parameter screen. This parameter is located by selecting the following: **Tools** menu → **User Options** submenu → **Progress Note** tab → **Set Note Sequence** button. The **Note Boilerplate** in the list of **Note Objects Sequence** parameter screen includes both TIU **Note Boilerplates** and supported CPRS templates.

All DRM Plus data objects displayed in the left narrow window of the **Progress Note** screen may be imported into the TIU progress note automatically. This depends on whether the parameters are activated by the end-user. However, most of the data objects are NOT allowed to be de-selected by the end-user for the automatic importing process. These parameters are located by selecting **Tools** menu → **User Options** submenu → **Progress Note** tab → **Progress Note Data** button.

DRM Plus **Code Boilerplates** are listed individually with the DRM Plus objects on the **Progress Note** screen. DRM Plus users who have created code boilerplates can now import them into a TIU progress note by clicking the desired **Code Boilerplate** object, listed with DRM Plus objects. Using the cursor, set the object where it is sequenced in the TIU progress note, then double-click to import the object.

Text files may be created and saved in a preferred directory. Right-click in a **Progress Note** window and click **Import Text File** option to navigate to and open the saved file. The file is then placed in the TIU progress note as designated by the cursor placement.

Follow these steps to set up a file for importing:

1. Create a folder in an appropriate directory (usually a server drive).
2. From the **Tools** menu → **User Options** submenu → **Set File Folder** button, navigate to the folder created in Step 1. Set the folder by double-clicking it.
3. Create a text file from the **Tools** menu → **Windows Notepad** submenu and save it as a .txt file in the designated folder from Step 1.

If the user clicks the **Finish** button before the provider's electronic signature is entered, the dental data is filed in VistA (PCE, DES and TIU) as an unsigned progress note, and may be viewed with CPRS. If the electronic signature is entered before clicking the **Finish** button, the dental data is filed as a signed progress note in VistA (PCE, DES and TIU).

When the user clicks the **Finish** button, a prompt may display if there are no planned items, and if the patient status is **Active**. Correctly identifying the patient's status is important for reporting. The user should click the proper radio button (**Active**, **Inactive** or **Maintenance**) and then click the **OK** button.

## **Reports – Non-Clinical Time by Provider**

The **Non-Clinical Time by Provider** report displays an approximate numerical unit of days (1 day = 8 hours). Accumulation of less than 4-hours results in rounding down to the nearest whole number day, and accumulation of 4 or greater hours round up to the nearest whole number day.

## **Code Boilerplates**

Multiple boilerplates may be added for a single code, or multiple codes may be associated with the same boilerplate. The user may establish as many boilerplates and related codes as necessary.

If a DRM Plus Administrator creates a code boilerplate in their **Administrative Toolbox** submenu, then every user may use the code boilerplate by entering the name of that code boilerplate to their parameter. Enter the name precisely as it was entered by the DRM Plus Administrator and click the **OK** button in the end-user's **User Settings** screen. This action imports the administrative code boilerplate to the end-user's code boilerplates.

## **Last Broker Call**

The **Last Broker Call** submenu is used by the IT or ADPAC personnel to document problems. It is NOT usually accessed by providers.

## **Recent Dental Activity**

This section on the **Cover Page** tab displays the most recent date for selected types of procedure codes. Hover the cursor over the heading to display all the ADA procedure codes, which comprise the data for that heading.

- **Last Qualifying Exam** = D0120, D0150 or D0180
- **Last Comprehensive Exam** = D0150 or D0160
- **Last Brief Exam** = D0120, D0140 or D0170
- **Last Periodontal Exam** = D0180
- **Last Panorex Image** = D0330
- **Last Full Mouth Image** = D0210
- **Last Bitewing Image** = D0270, D0272, D0273, D0274 or D0277
- **Last CBCT Image** = D0364, D0365, D0366, D0367, D0368, D0380, D0381, D0382, D0383 or D0384
- **Last Prophylaxis** = D1110, D4341, D4342, D4346, D4355 or D4910
- **Last Visit** = the last dental visit date the encounter was filed on.
- **Last Provider** = the provider that filed the last encounter.

## Appendix D – Icon Definitions

### Diagnostic Findings

The following table explains the actions required to enter a **Diagnostic Finding**:

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Restore	As many as required	As many as required	No	Material or surfaces	Graphic in green
Missing	As many as required	No	As many as required	Selected Roots	Graphic for roots is outlined
Observe	As many as required	No	As many as required	Selected Roots	Tooth outlined in red
Partial	As many as required	No	No	Cancel or Complete	Graphic allows root condition graphics to show. Graphic in purple.
Denture	Any tooth	No	No	None	Graphic allows root condition graphics to show. Graphic in purple.
Implant	As many as required	No	No	None	Graphic in violet
Sealant	No	As many as required	No	None	Graphic in green
Endo	No	No	Yes	Materials and roots	Graphic color denotes material
Apico	No	No	Yes	Select roots	Graphic in red
Retro	No	No	Yes	Selected roots and materials	Requires Apico to be present. Graphic denotes material
Bridge					See special instructions
Conn Bar					See Special Instructions
Hemi	No	No	As many as required	Selected roots	Graphic in dark gray

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Coping	Yes	No	No	None	Graphic in green
P&C	No	No	As many as required	Selected roots and materials	Graphic in green and denotes material
Impact	As many as required	As many as required	As many as required	Selected surfaces and roots	Graphic in light blue, roots in blue-green
Def Rest	Yes	As many as required	Yes	Selected surfaces, roots and materials	Graphic in yellow, denotes material
Caries	Yes	As many as required	As many as required	Selected surfaces and roots	Graphic in red, root caries initiates description box
Drifting	Yes	No	No	Direction	Graphic is yellow arrow to the left of tooth
Tipped	Yes	No	No	Direction	Graphic is light blue arrow to left of tooth
Rotated	Yes	No	No	Direction	Graphic is green arrow to left of tooth
Ret Root	No	No	As many as required	None	Graphic removes crown
UndrCont	No	As many as required	No	Selected surfaces	Graphic is red and yellow
OverCont	No	As many as required	No	Selected surfaces	Graphic is red and yellow
Overhang	No	As many as required	No	Selected surface	Graphic is red and yellow
Lesion	No	No	As many root surfaces as required	Selected roots	Can also click implant. Graphic is a red circle.
Faceted	No	As many as required	No	Selected surfaces	Graphic is red and yellow

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Cracked	No	As many as required	As many as required	Selected surfaces and roots	Graphic is red and yellow
Chipped	No	As many as required	No	Selected surfaces	Graphic is red and yellow
Supr/Sub	Yes	No	No	Direction	Graphic is red arrow
Open Ct	Yes	No	No	None	Graphic is red arrow to right of tooth
Abfract	No	Facial or lingual	No	None	Graphic is blue arrows
Dentition	No	No	No	Dentition box	Converts graphic to juvenile. Must be done before other entries
Perm/Prim	Yes	No	No	Dentition Box	Designate selected tooth as primary or permanent
Endentulous	Any tooth	No	No	None	Graphic removes all teeth and roots in arch

**Note:** Certain **Diagnostic Findings** or **Completed Care** procedures, once entered, display graphically on all screen views. These items are, if entered from the **Diagnostic Findings** screen: missing, implant, impacted, retained root, hemi section, dentition, and observe. If entered from the **Completed Care** screen: extract, hemi section, implant and observe.



## **Treatment Plan**

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Restore (2)	Yes	As many as required	No	Code Selection box	Graphic in Blue
Extract	Yes	No	Yes	Code Selection Box	Graphic in shadow
Observe	As many as required	No	No	None	Tooth outlined in red
Partial (1)	As many as required	No	No	Cancel or Complete	Graphic allows root condition graphics to show
Denture	Any tooth	No	No	Code Selection box	Graphic allows root condition graphics to show.
Implant	Yes	No	No	Code Selection Box	Graphic in violet, Diagnostic Finding must be Missing
Sealant	No	As many as required	Yes	Code Selection Box	Graphic in blue, no root graphic
Endo	No	No	Yes	Code Selection Box	Graphic in blue
Apico	No	No	Yes	Code Selection Box	Graphic in red
Retro	No	No	Yes	Code Selection Box	Requires Apico to be present. Graphic in blue.
Bridge (1)				Code Selection Box	See special instructions.
Conn Bar (1)				Code Selection Box	See special instructions

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Hemi	No	No	As many as required	Code Selection Box	Graphic in gray
Coping	Yes	No	No	Code Selection Box	Graphic in blue
P&C	No	No	As many as required	Code Selection Box	Graphic in blue
Perio But- tons	No	No	No	None	No Graphic. See Speed Code instructions.

Certain codes associated with these icons are defined as ranged codes. Ranged codes can only be entered by clicking the proper icon button. The **ADA Codes** icon, **CPT Codes** icon, **Quick Codes** icon or **Speed Code** icons are NOT allowed for entering ranged codes.

To designate a root restoration, click the **Restore** icon. Click the tooth surface that corresponds to the root surface. Using a tooth note, or the description edit feature in the transaction table, explain the root restoration.

**Note:** Certain **Diagnostic Findings** or **Completed Care** procedures, once entered, appear graphically on all screen views. If these items are entered from the **Diagnostic Findings** screen, they are missing, implant, impacted, retained root, hemi section, dentition, and observe. If entered from the **Completed Care** screen, they are extract, hemi section, implant, and observe.

## **Completed Care**

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Restore (2)	Yes	As many as required	No	Code Selection Box	Graphic in Green
Extract	Yes	No	Yes	Code Selection Box	Tooth disappears. Root graphic in dark gray.
Observe	As many as required	No	No	None	Tooth outlined in red
Partial (1)	As many as required	No	No	Cancel or Complete. Then Code Selection box.	Graphic allows root condition graphics to show. Graphic in purple.
Denture	Any tooth	No	No	Code Selection Box	Graphic allows root condition to show. Graphic in blue-purple.
Implant	Yes	No	No	Code Selection Box	Graphic in violet, Diagnostic Finding must be Missing
Sealant	No	As many as required	Yes	Code Selection Box	Graphic in green. No root graphic.
Endo	No	No	Yes	Code Selection Box	Graphic in pink.
Apico	No	No	Yes	Code Selection Box	Graphic in red.
Retro	No	No	Yes	Code Selection Box	Requires Apico to be present. Graphic in green.
Bridge (1)				Code Selection Box	See Special Instructions

ICON	CLICK TOOTH	CLICK SURFACE	CLICK ROOT	POP-UP SCREEN	ADDITIONAL COMMENTS
Conn Bar (1)				Code Selection Box	See Special Instructions
Hemi	No	No	As many as required	Code Selection Box	Graphic in gray
Coping	Yes	No	No	Code Selection Box	Graphic in green
P&C	No	No	As many as required	Code Selection Box	Graphic in green
Perio But- tons	No	No	No	None	No Graphic. See Speed Code instructions.

Certain codes associated with these icons are defined as ranged codes. Ranged codes can only be entered by clicking the proper icon button. The **ADA Codes** icon, **CPT Codes** icon, **Quick Codes** icon or **Speed Code** icons are NOT allowed for entering ranged codes.

To designate a root restoration, click the **Restore** icon. Click the tooth surface that corresponds to the root surface. Using a tooth note, or the description edit feature in the transaction table, explain the root restoration.

**Note:** Certain **Diagnostic Findings** or **Completed Care** procedures, once entered, appear graphically on all screen views. If these items are entered from the **Diagnostic Findings** screen, they are missing, implant, impacted, retained root, hemi section, dentition, and observe. If entered from the **Completed Care** screen, they are extract, hemi section, implant, and observe.

## Special Descriptions – Bridge Icon

Place the cursor on the first tooth included in the bridge and **drag** the cursor to the last tooth in the range of teeth included in the bridge so that the 'bar' drawn for the bridge is **one continuous line** encompassing all of the areas (teeth, whether missing or not) included in the bridge. The Code Selection screen displays, with the lowest number tooth selected, as shown below.

Code	Description
D6750	RET CROWN PORCLN HIGH NOBLE
D6751	RET CROWN PORCLN BASE METAL
D6752	RET CROWN PORCLN NOBLE METAL
D6780	RET CROWN 3/4 HIGH NOBLE
D6781	RET CROWN 3/4 BASE METAL
D6782	RET CROWN 3/4 NOBLE METAL
D6783	RET CROWN 3/4 PORC/CERAMIC
D6790	RET CROWN FULL HIGH NOBLE

**Figure 239: Procedure Code Selection Screen**

Select the appropriate code and click the **Add** button. The **Procedure Code** selection screen moves to the first pontic tooth. Select the appropriate code and click the **Add** button again. Continue this process until all required teeth and codes have been selected. Click the **Finished** button once the selection process is complete.

**Note:** The [<<] and [>>] buttons may be clicked to move backward or forward to different teeth for code selection. The **Reset** button, when activated by using the [<<] and [>>] buttons, clears all previously entered codes for the selected tooth.

## **Special Descriptions – Conn Bar Icon**

This functionality works for all three view screens of **Diagnostic Findings, Treatment Plan** and **Completed Care**. Place the cursor on the first tooth location and drag it to the final tooth location. When entering a connector bar from the **Treatment Plan** screen or the **Completed Care** screens, the **Procedure Code** selection screen displays, with no tooth selected. DRM Plus defaults to the correct connector procedure code, depending on what conditions the connector bar was entered on. Click the **OK** button to complete the connector bar entry.

## **Special Descriptions – Notes Icon**

Please see the **Notes** section in the Chart/Treatment – Periodontal Chart chapter of this manual. This icon works the same for all three **Treatment & Exam** screen views.

Items to consider on tooth notes:

- Teeth designated as primary show in the tooth drop-down menu with the appropriate letter without a number. A tooth designation for **Supernumerary** teeth is displayed after tooth #32 in the drop-down menu.
- When a tooth-specific note has been entered, the tooth number in the graphical chart on the **Diagnostic Findings, Treatment Plan, Completed Care** and **Periodontal Chart** screens displays in yellow.
- Previously entered notes (tooth/patient note) appear disabled. The note appears disabled if it was saved as unfiled data and DRM Plus is closed and reopened. It also becomes disabled when an assistant saves unfiled data to a provider. When the provider re-opens DRM Plus to complete the encounter, they may edit-delete the note (tooth/patient note) by clicking the **New Entry** button. This activates the disabled entry, so the provider can modify or delete the note before the encounter is finished.
- After entering any new dental data into **Patient Notes** or **Tooth Notes** windows select the **Notes** icon to save and close the **Notes** screen. Closing the **Notes** screen before opening any other screen in DRM Plus is the only action that saves the data entered in the **Notes** screen.

# Appendix E – Using the Keyboard to Enter Periodontal Data

## Overview

The **Periodontal** screen is designed for data entry using the mouse. Data entry using the keyboard is also an option. Clicking the **Keyboard Mode** key, <F10>, initiates keyboard navigation and data entry on this screen.




Figure 240: Keyboard Mode Button on Periodontal Chart Screen

When this button is clicked, all the icons are disabled, and the keyboard is activated. The keyboard functions are described in this section.

## Navigating the Periodontal Screen

Use the following keystrokes to change the screen views.

### Arch Views

The screen moves from the existing view to any of the other views by using the following keys:

**U** = Upper

**L** = Lower

**N** = Lingual

**F** = Facial

**F11** = Full (This screen must be closed by using the mouse.)

### Cursor Movement

There are four options for moving the cursor to select a tooth/surface.

**Enter:** Moves the cursor one surface in the direction of the higher numbered tooth.

**Backspace:** Moves the cursor one surface in the direction of the lower numbered tooth.

**">" with or without the <Shift> key:** Moves the cursor one surface in the direction of the higher numbered tooth.

**"<" with or without the <Shift> key:** Moves the cursor one surface in the direction of the lower numbered tooth.

Press the <A> key to toggle the **Auto Advance** function on or off.

## Entering Data

Entering data from the keyboard requires the cursor to be placed on the desired tooth/surface. The user must then select the desired condition and enter the data values in the appropriate manner.

**Note:** All numeric values must be entered with two digits (1 = 01, 2 = 02 and 10 = 10). Entering FGM and Mobility require the values to have a prefix (see below).

**K = Pocket:** Press the <K> key and with the cursor in the correct position, enter a two-digit value. Then move the cursor to the next surface.

**G = FGM:** Press the <G> key and with the cursor in the correct position, enter a plus sign (+) and a two-digit value, or just a two-digit value. Then move the cursor to the next surface.

**J = MGJ:** Press the <J> key and with the cursor in the correct position, enter a two-digit value. Then move the cursor to the next surface.

**B = Bleeding:** Press the <B> key with the cursor in the correct position.

**D = Delayed Bleeding:** With the cursor in the correct position, press the <D> key.

**S = Suppuration:** With the cursor in the correct position, press the <S> key.

**O = Mobility:** Press the <O> key with the cursor in the correct position and enter a two-digit value. For a value of 1 1/2, 2 1/2, etc., enter the two-digit value preceded by a plus sign “+”.

**I = Furcation:** Press the <I> key with the cursor in the correct position and enter a two-digit value. Move the cursor to the next surface.

**R = Reset:** Pressing the <R> key resets all values to zero (use this functionality with extreme caution).

**Ctrl Z = Undo.**

## Special Buttons

Viewing the special buttons screen requires pressing the following keys:

**H** = History

**C** = Compare

**M** = Summary

**E** = Head & Neck

**P** = PSR

**Q** = Stats

**Note:** Displaying these screens using the keyboard turns off the keyboard function for these screens. The mouse is required to navigate these screens.

## Other Functions

**Z** = Cal and **X** = Lock

For convenience, a tear-out of the **Periodontal Keyboard** shortcut chart is available on the last page of this manual.



## Appendix F – Ranged Codes


Using the graphic icon on the **Treatment Plan** and **Completed Care** screens is not only straightforward, it minimizes the potential for errors. Certain codes, designated as ranged codes, used for multiple teeth procedures, mostly in the prosthetics area, have been restricted to icon use only for procedure entry. These codes are listed below.

Removable Prosthodontics	Maxillofacial Prosthetics	Fixed Prosthodontics	Fixed Prosthodontics
D5130	D5934	D6205	D6782
D5140	D5935	D6210	D6783
D5211		D6211	D6790
D5212		D6212	D6791
D5213		D6214	D6792
D5214		D6240	D6793
D5221		D6241	D6794
D5222	Implant Services	D6242	D6920
D5223	D6055	D6245	
D5223	D6068	D6250	
D5224	D6069	D6251	
D5225	D6070	D6252	
D5226	D6071	D6253	
D5281	D6072	D6545	
D5820	D6073	D6548	
D5821	D6074	D6549	
D5863	D6075	D6600	
D5864	D6076	D6601	
D5865	D6077	D6602	
D5866	D6112	D6603	
	D6113	D6604	
	D6114	D6605	
	D6115	D6606	
	D6116	D6607	
	D6117	D6608	
	D6194	D6609	
		D6610	
		D6611	
		D6612	
		D6613	
		D6614	
		D6615	
		D6624	
		D6634	
		D6710	
		D6720	
		D6721	
		D6722	
		D6740	
		D6750	
		D6751	
		D6752	
		D6780	
		D6781	

## Appendix G – Option to Set Dental Patients to Inactive Status

A new option that IT may run (or may give as a secondary VistA option to DRM Plus users) allows the user to check the system for patient activity and set patients without current encounters to **Inactive** status.

Use the **OPTION NAME: DENTV INACTIVATE PATIENTS** and then select an inactivate date (defaulted to 2 years).

A screenshot of a Telnet window titled "127.0.0.1 NTI - Cache Telnet". The window has a menu bar with "File", "Edit", and "Help". The main text area shows a VistA command prompt "VAH>D ^XUP". The output includes: "Setting up programmer environment", "This is a TEST account.", "Access Code: \*\*\*\*\*", "Terminal Type set to: C-VT320", "You have 352 new messages.", "Select OPTION NAME: DENTV INACTIVATE PATIENTS", "Inactivate Dental Patients", "Inactivate Dental Patients", "Patients not seen in dental within the selected time frame (defaulted to 2 years) will be inactivated!", "Select inactivation date: T-730//", "Successfully queued, task# 3516731", and "VAH>".

```
127.0.0.1 NTI - Cache Telnet
File Edit Help

VAH>D ^XUP

Setting up programmer environment
This is a TEST account.

Access Code: *****
Terminal Type set to: C-VT320

You have 352 new messages.
Select OPTION NAME: DENTV INACTIVATE PATIENTS      Inactivate Dental Patients
Inactivate Dental Patients
Patients not seen in dental within the selected time frame
(defaulted to 2 years) will be inactivated!
Select inactivation date: T-730//
Successfully queued, task# 3516731
VAH>
```


Figure 241: VistA Option to Set an Inactivation Date

## Appendix H – How to Map Dental CNTs

**Note:** These procedures require IT assistance. The dental staff should provide these instructions to IT.

In order to confirm that the CNTs are mapped correctly:

1. Open the **\DOCSTORE\Array\dnarray.txt** file located on the server to confirm that within the dnarray.txt file, the patch for each CNT is correct.



**Figure 242: Map Dental CNTs from dnarray.txt file**

2. Open one of the CNT .ini files within the \DOCSTORE\FORMS directory (i.e., \\vhaserver-name\DOCSTORE\FORMS\160\_DENT\04160001\DRMEval.ini)




Figure 243: DRMEval.ini File

3. Confirm that the VHAservername path matches what was in the \Array\dntarray.txt file and is correct.
4. If the VHAservername within any of the \FORMS\160\_DENT\#####.ini files does NOT match the directory shown above (in the ARRAY folder) then each .ini file within the FORMS directory need to be opened and edited to reflect the correct path.

## Appendix I – Recommendations for Coding of Prosthetic Appliance

Coding for prosthetic appliances is to be done at the time when the prostheses are delivered to the patient and/or home care instructions are provided and documented.

Taking workload credit for the undeliverable prosthetic appliances should occur:

1. After death.
2. After 6 months and 3 attempted telephonic contacts. Attempts to reach patient should occur at a timely interval (i.e., at least one week apart).
3. No response to a final letter to patient's last known address.
4. CPRS documentation to insure no citation for patient abandonment.

To take workload credit with no associated patient visit (i.e., phone contact):

1. Follow local policies for documentation of telephonic contact. Codes for phone contacts in DRM Plus are listed under the CPT codes (99441, 99442 and 99443).

To file workload credit:

1. Enter code(s) in the **Completed Care** screen.
2. Select the correct previously filed **Visit Date/Time** associated with this addendum.
3. Click the **Next** button and select the **File Data With a Note Addendum** radio button.

### Figure 244: Filing Options Screen

- 222

## Appendix J – Business Use of DRM Plus

As important as it is to know how to use DRM Plus, it is equally important to be aware of how to enter codes, fees and similar DRM Plus business considerations.

Entering DRM Plus business information may be viewed in three different components:

- Local policy and practice
- National policy and practice
- National business practice

Adhering to the following guidelines for each is essential and critical to the success of treating our Veterans. Entering valid data can provide important information on the allocation of VA resources on a local and national level, clarify current and future funding issues, and determine how to provide even better care to patients.

Review this information carefully. If the user has any questions about this software or the accompanying business policies, contact the local DRM Plus Subject Matter Expert on site. If this person cannot answer the question, s/he knows who to contact to find the correct solution.

### **Local Policy and Practice**

This includes the workflow process addressing the provider's data entries for the following:

- Diagnostic Findings
- Observations
- Approval for proposed Patient Treatment Plans are to be established by the local dental manager.

### **National Policy and Practice Coding Standards**

All completed procedures on site must be entered into the DES through DRM Plus. This includes the following:

- Procedures by staff
- Fee-basis on site
- Sharing on site
- Contract on site
- Residents
- Without compensation
- Students
- Hygienists

For example, if a surgeon goes to the OR and enters the procedure into PCE through the surgery package, or through an encounter, that procedure must be entered into DES through DRM.

Coding standards should be followed to calibrate providers, if the same encounters are observed at two separate clinics. This ensures that the encounters are coded the same.

## Appendix K – Data Security

As a VA computer user, one of the best and most important ways to contribute to good computer security is to know all data, its level of sensitivity, that it is virus-free, what would happen if it were unavailable, how long it could be done-without, and the effect of another user changing it without approval.

Classifying data involves determining how sensitive and valuable it is, and what protection it needs. Information is classified according to sensitivity, which is based on its need for:

- **Confidentiality:** The information must be kept private as its owner instructs.
- **Integrity:** The information must not be inappropriately changed or destroyed.
- **Availability:** The information must be ready for use, as needed.

The amount of information and the context in which it is found can affect its value. Some information is confidential only at certain times (i.e., contracting or economic forecast information, which is sensitive until its publication or release date, after which it is made public). Current information is generally more valuable than older information.

When protecting data, all employees and contractors have a responsibility to:

- Be familiar with VA security policies, procedures, rules and regulations (i.e., know what to do, how to do it and why).
- The user should ask a supervisor or ISO any questions about these security responsibilities.

The user is responsible for:

- Reporting known or suspected incidents immediately to the ISO.
- Using VA computers only for lawful and authorized purposes.
- Choosing good passwords and changing them every 90 days. Do not write down or share log-in information with anyone, including Help Desk.
- Complying with safeguards, policies and procedures to prevent unauthorized access to VA computer systems.
- Recognizing the accountability assigned to the user's UserID and password. Each user must have a unique ID to access the VA systems. Recognize that UserIDs are used to identify an individual's actions on VA systems and the Internet. Individual user activity is recorded, including sites and files accessed on the Internet (recorded as the files go through the firewall).
- Ensuring that data is backed up, tested and stored safely.
- Not generating or sending offensive or inappropriate email messages, graphical images or sound files. Limit distribution of email to only those who need to receive it. Realize that the user is identified as a user of the VA computer systems when logged on to the Internet.
- Using authorized virus scanning software on the workstation or PC and home computer. Know the source before using discs or downloading files. Scan files for viruses before execution.
- Complying with terms of software licenses and only using VA-licensed and authorized software. Do not install single-license software on shared hard drives or servers.
- Complying with terms of software licenses and using only VA-licensed and authorized software. Do not install single-license software on shared hard drives or servers.



- Knowing data and properly classifying and protecting it, as well as inputs and outputs, according to their sensitivity and value. Label sensitive media use a screen saver with a password, logoff when leaving the work area, and secure that sensitive information is removed from hard disks sent out for maintenance. Do not send sensitive information over the Internet unless it has been encrypted.
- Learning as much as possible about information security to assist the user's ISO. Numbers alone make users the most important security asset. Compared to one ISO for a system, users offer a chance for numerous eyes and ears to remain alert to potential threats to information systems.

## Appendix L – madExcept Application

The **madExcept** application is a tool that has been added to DRM Plus to assist with error reporting from the field and implementing a fix in the application. In order to be prepared to use this tool if an error occurs while using DRM Plus, please review the following directions:

Select all the **OK** buttons from any traditional error screen that may display in DRM Plus. There may be more than one traditional error screen or there may be none, but no matter how many, select the **OK** button on all. If any informational screen displays asking if the user would like to view the last broker call, select the **No** button from that screen and continue through this process.

As soon as the **DENTALMRMTX.EXE** error screen appears, select the ‘**send system report**’ button (first button from the left) as displayed:





Figure 245: DENTALMRMTX.EXE Screen

After selecting the ‘**send system report**’ button, the **Contact Information** screen will appear. The screen asks for contact information so enter provider’s name and VA email address. Select the **Continue** button on the **Contact Information** screen to move on to the next step.

The image shows a "Send Assistant" dialog box with a "Contact Information" tab selected. The tab is labeled with a "1" in a box, and a "2" is visible next to it. The form contains two text input fields: "your name:" with the placeholder text "your name" and "your email:" with the placeholder text "name@va.gov". Below these fields is a checkbox labeled "remember me" which is checked. At the bottom of the dialog, there are three buttons: "Continue", "Skip", and "Cancel".

Figure 246: Contact Information Screen


The **Error Details** screen will immediately appear. Add information in the in the **Error Details** screen providing as much information as possible. If known, list every click that occurred prior to the error message. The more details listed; the easier it will be to reproduce the error and fix it.



**Figure 247: Error Details Screen**

Select the **Continue** button from the **Error Details** screen.

The screen that follows is when the application opens the dental users email screen. The example used is MS Outlook. The email will not be addressed or populated with any individual or group in the email address field. The dental user will need to add the individuals or group where the system report should be sent too. The system report will be attached automatically.



**Figure 238: Outlook Window**

Enter the name and location (include city and state) of the VA dental clinic where the error occurred. For example, if the error occurred at the Daytona Beach Dental Clinic, enter Daytona Beach, Florida and do NOT enter Gainesville, Florida, even though Daytona Beach is a Gainesville satellite clinic. Also include, if known, the person class or provider type/specialty and the phone number that may be used to contact the person reporting the error. Enter all this data in the email address window.


The user may add any other individuals as recipients (To and CC) in the email address fields as appropriate. Additional information may be included by entering that information in the email address window.

Click on the **Send** button from the email screen.

Select the '**restart application**' button (bottom button on the right) from the **DENTALMRMTX.EXE** screen. This will close the **madExcept** tool and restart DRM Plus.

**Note:** After an error occurs while using DRM Plus, please reboot the computer before continuing.

# Periodontal Keyboard Shortcuts Tear-Out



**Figure 249: Periodontal Keyboard Entry Screen**