



**Prosthetics
Remove IFCAP Item Description
From 2319**

Release Notes

Patch RMPR*3*99

Version 3.0

November 2004

Table of Contents

Remove IFCAP Item Description from 2319 Release Notes	1
Patch RMPR*3*99 Overview.....	1
NIF IFCAP Changes in the 2319 (Roll and Scroll)	3
View Prosthetics Item Transactions (Read Only).....	6
GUI 2319 Button on DOR Menu – Appliance Transactions (Tab 4).....	7
GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8).....	9
Appendix A.....	11
Glossary	11

Remove IFCAP Item Description from 2319 Release Notes

Patch RMPR*3*99 Overview

Introduction

Patch **RMPR*3*99 (Remove IFCAP Item Description from 2319)** was developed for Prosthetics to provide an interim solution for the continuation of the National Item File (NIF) data cleansing effort.

On the national level, the IFCAP Item file is pulled from each facility's database and merged into a VISN database. Then the VISN's NIF team reviews the Items and edits, deletes and adds descriptions (which are considered to be the "*data cleansing*" process).

If the IFCAP Item (at the national level) was a Prosthetic Item that was changed, added or deleted, then the Prosthetic Patient Treatment Record (2319) was historically being changed to match the data cleansing changes.

Patch RMPR*3*99 prevents the Prosthetic 2319 from historically being changed when data cleansing happens. This patch protects the **PROSTHETICS ITEM DESCRIPTIONS** in our Prosthetic 2319 from being overridden when the NIF make changes to the **IFCAP DESCRIPTIONS**.

Patch RMPR*3*99 changes

Patch RMPR*3*99 affects multiple **Prosthetics Menu** options. All three types of the Prosthetics 2319 have changes including:

- 1) Roll and Scroll
- 2) Read Only Roll and Scroll
- 3) GUI 2319 from the Delayed Order Report (DOR) application.

This patch will change the Prosthetics 2319 *Appliance Transactions* and *Home Oxygen Items*, switching the **IFCAP Item Description** to the **HCPCS Description** display.

Note to IRM

There is no new .EXE file to install. All the changes are M routines and a new cross reference in File 660. The **IFCAP Item Description** will remain in the database until phase 2 of the NIF implementation.

This is the menu option name that IRM references: **RMPR VIEW 2319 READ ONLY**.

Note: Please read these Release Notes for further information.

Continued on next page

Patch RMPR*3*99 Overview, Continued

No HCPCS entered prior to 1997

Since HCPCS codes were not used prior to 1997, you will continue to see the **IFCAP Item Description** that was selected at that time. This is considered to be protected data and will not change.

Even though the **Appliance Transaction (Screen 4)** column heading displays **HCPCS** (in the Roll and Scroll), the information will continue to display the original **IFCAP Item Description** for those transactions entered prior to 1997.

Even though the GUI version of the **Appliance Transaction (Tab 4)** displays **Item** as the column header, it is actually the **IFCAP Item** for those transactions entered prior to 1997. The column heading was **not** changed to **HCPCS**.

NIF IFCAP Changes in the 2319 (Roll and Scroll)

2319 Roll and Scroll changes

Access the **Prosthetics Official's Menu** option, **Display/Print Patient 2319 (23)** to view the **Appliance Transactions (Screen 4)** or **Home Oxygen (Screen 8)**. Below are the sample screens for the Roll and Scroll 2319 that were affected by the NIF IFCAP changes.

Display/Print Patient 2319 Menu (23)

```
Select Prosthetic Official's Menu Option: DD <Enter> Display/Print

23   Display/Print Patient 2319
TI   Transaction Inquiry
IP   Print All Prosthetic Items
BI   Print Prosthetic Billings for MAS
IH   Item History
SE   Search for Recalled Item
SI   Site Parameter Inquiry
VI   Vendor Inquiry

Select Display/Print Option: 23 <Enter> Display/Print Patient 2319
```

Select the patient

```
Select PROSTHETIC PATIENT:   PROSPATIENT,ONE <Enter>   LOMA
LINDA VAMC      6-10-31
000040707      YES      70%      SC VETERAN      LL/SD/
MODULE III - GIM PCP:
PROSPROVIDER,ONE CX  **ADVANCE DIRECTIVE - NO **

DEVICE: HOME// <Enter> VIRTUAL CONNECTION   Right Margin: 80// <Enter>

Current Disability Codes are:
AO/DIS          ALL OTHER S/C  S/C
```

Select the Appliance Transactions (Screen 4)

```
Select one of the following:

      1      PATIENT DEMOGRAPHICS
      2      CLINIC ENROLLMENTS/CORRESPONDENCE
      3      ENTITLEMENT INFORMATION
      4      APPLIANCE TRANSACTIONS
      5      AUTO ADAPTIVE INFORMATION
      6      CRITICAL COMMENTS
      7      ADD/EDIT DISABILITY CODE
      8      HOME OXYGEN ITEMS

Enter 10-2319 screen to VIEW (1-8), '^' to EXIT, or 'return' to
continue : 4 <Enter> APPLIANCE TRANSACTIONS

HHUU,UXXHUS G          SSN: 000-04-0707   DOB: JUN 10,1931   CLAIM# SS
```

Continued on next page

NIF IFCAP Changes in the 2319 (Roll and Scroll), Continued

HCPCS Column

Notice the **HCPCS** column heading below is highlighted. This is the column of data that has been changed with Patch RMPR*3*99 to the **HCPCS Description** from the **IFCAP Item Description** (which was the **Item** column heading).


New HCPCS Column Title

Date	Qty	HCPCS	Type	Vendor	Sta	Serial	Delivery Date	Tot Cost
1. 05/29/02	1	AUTO BLOOD	I	VENDOR,ONE	605		05/29/02	42.01
2. 07/26/00	1	VASCULAR G	I	VENDOR,ONE	605	019579	07/26/00	99.19

End of Appliance/Repair records for this veteran!
 +=Turned-In *=Historical Data I=Initial X=Repair S=Spare R=Replacement
 Enter 1-2 to show full entry, '^' to exit or `return` to continue. 1 <Enter>

PROSPATIENT,ONE SSN: 000-00-0001 LOMA LINDA VAMC DOB: 99-99 -1931
APPLIANCE/REPAIR LINE ITEM DETAIL <4-1>
 TYPE OF FORM: STOCK ISSUE INITIATOR: PROSPROVIDER,THREE
 DATE: MAY 29, 2002
 DELIVER TO:
 TYPE TRANS: INITIAL ISSUE QTY: 1 SOURCE: COMMERCIAL
 VENDOR: PROSVENDOR,ONE
 VENDOR PHONE: 800-555-5555
 28690 N. Any DRIVE
 LAKE FORREST, ILLINOIS 60045
 DELIVERY DATE: MAY 29, 2002
 TOTAL COST: \$42.01 OBL:
 REMARKS: [If you type a specific note here, it appears on Screen 4 under the transaction with notes...see below.]
 DISABILITY SERVED: SC/OP
ITEM DESCRIPTION: BLOOD PRESSURE UNIT DIGITAL
 APPLIANCE: BLOOD PRESSURE UNIT DIGITAL
 PSAS HCPCS: A4670 AUTO BLOOD PRESSURE MONITOR
 ICD-9 Code: 401.9 HYPERTENSION NOS
 CPT MODIFIER:
 DESCRIPTION: AUTO BLOOD PRESSURE MONITOR
 EXTENDED DESCRIPTION:
 AUTOMATIC BLOOD PRESSURE MONITOR
 Enter RETURN to continue or '^' to exit:

The Item Description is new!!



Remarks

If you enter a specific note in the **Remarks** field, it will be easier to locate it to identify it for future Item needs. The note displays on Screen 4 below the transaction.

Item Description / Appliance prompt

The **IFCAP Item Description** field display above is new. This is the field you should review to determine which **Prosthetic Item** was provided to the patient. You should no longer review the **Appliance** field. When the **Item Description** and the **Appliance** field are identical (as shown above), then the “new note” that is described on the next page will not appear. See next page.

Continued on next page

NIF IFCAP Changes in the 2319 (Roll and Scroll), Continued

Item Description NEW Note

With Patch RMPR*3*99, whenever an **IFCAP Item Description** has been changed, edited or deleted, you will see a new note under this field as follows:

***** See Above for Original Item Description *****

Due to the data cleansing, the **Appliance** field (which was the original Item selected) has a new description and the original description is now the **Item Description** field.

After Patch Install - Item Description with NOTE

=====>
=====>

```
TYPE OF FORM: STOCK ISSUE
INITIATOR: PROSPROVIDER,ONE
DATE: MAY 29, 2002
DELIVER TO:
TYPE TRANS: INITIAL ISSUE      QTY: 1      SOURCE: COMMERCIAL
VENDOR:      PROSVENDOR,ONE
VENDOR PHONE: 800-555-5555
              28690 N. Any DRIVE
              LAKE FORREST,    ILLINOIS      60045
DELIVERY DATE: MAY 29, 2002
TOTAL COST: $42.01              OBL:
REMARKS:
DISABILITY SERVED: SC/OP
ITEM DESCRIPTION: BLOOD PRESSURE UNIT DIGITAL
*** See Above For Original Item Description ***
APPLIANCE: BLOOD PRESSURE CUFF DIGITAL
PSAS HCPCS: A4670      AUTO BLOOD PRESSURE MONITOR
ICD-9 Code: 401.9    HYPERTENSION NOS
CPT MODIFIER:
DESCRIPTION: AUTO BLOOD PRESSURE MONITOR
EXTENDED DESCRIPTION:

AUTOMATIC BLOOD PRESSURE MONITOR

Enter RETURN to continue or '^' to exit:
```

Example above

To interpret the example above, the original **Item Description** was “*Blood Pressure Unit Digital.*” Due to data cleansing, the description was changed to “*Blood Pressure Cuff Digital.*”

The **Appliance Field** was changed to the new description: “*Blood Pressure Cuff Digital*” and the original or saved description appears in the new field: **Item Description**.

View Prosthetics Item Transactions (Read Only)

“Read Only” patch change

The change with Patch RMPR*3*99 that was previously explained also affects the Roll and Scroll version of Prosthetics and provides the following new (Read only) menu option:

View Prosthetics Item Transactions

Below is a sample screen.

View 2319 (Read Only)

```
Select OPTION NAME: VIEW PROSTHETICS ITEM TRANSACT <Enter> RMPR VIEW 2319
READ ONLY

View Prosthetics Item Transactions

SITE: SUPPORT ISC// <Enter>                               578A5

Select PROSTHETIC PATIENT: PROSPATIENT, ONE <Enter>      SUPPORT ISC      4-5-61
0001111111
    YES      SC VETERAN

PATIENT IS DECEASED. DATE OF DEATH WAS OCT 21,1993@08:00

Would you Like to continue Processing this Patient? NO// YES <Enter>

DEVICE: HOME// <Enter>  TELNET      Right Margin: 80// <Enter>

    Select one of the following:

        I      ITEM TRANSACTIONS
        H      HOME OXYGEN ITEMS

Enter DATA screen to VIEW (Item Transactions or Home Oxygen), '^' to EXIT, or
'return' to continue :
```


GUI 2319 Button on DOR Menu – Appliance Transactions (Tab 4)

Delayed Order Reports (DOR)

The **GUI 2319 Appliance Transactions (Tab 4)** window from the Delayed Order Report (DOR) below shows the heading: **Item** column; however this information is now the **HCPCS Description** with Patch RMPR*3*99.

The column heading **Item** will not be changed to **HCPCS**.

Appliance Transactions (Tab 4)

Date	Qty	Item	Typ	Vendor	Station	Serial	HCPC	Total Cost
12/17/03	1	LATERAL TROCHANTERIC PAD	I	PROSVENDOR.TWO	605		L1290	134.50
12/17/03		SHIPPING	X	PROSVENDOR.TWO	605			50.00
08/22/03	1	OTHER ACCESSORIES	I	PROSVENDOR.THRE	605		K0108	35.89
07/22/03	1	OXYGEN CONCENTRATOR	I	PROSVENDOR.FOUR	605	0296AZ1862	E1390	1994.85
01/28/03	1	DME DELIVERY, SET-UP	X	PROSVENDOR.FIVE	605		A9901	72.15
01/27/03	1	STANDARD WHEELCHAIR	I	PROSVENDOR.THRE	605	20502086256	K0001	161.02
01/09/03	1	DME DELIVERY, SET-UP	X	PROSVENDOR.SIX	605		A9901	108.75
01/09/03	1	DME DELIVERY, SET-UP	X	PROSVENDOR.FIVE	605		A9901	72.15
01/09/03	1	BATH/SHOWER CHAIR	I	PROSVENDOR.FIVE	605		E0240	12.06
01/09/03	1	WAIST BELT	I	PROSVENDOR.EIGH	605		L2190	10.50
01/09/03		SHIPPING	X	PROSVENDOR.EIGH	605			9.95
01/08/03	1	PATIENT LIFT SLING OR SEAT	I	PROSVENDOR.SEVE	605		E0621	63.19
01/08/03	1	PATIENT LIFT HYDRAULIC	I	PROSVENDOR.FIVE	605	02L360416	E0630	543.21
01/08/03	1	DME DELIVERY, SET-UP	X	PROSVENDOR.FIVE	605		A9901	72.15
12/30/02	1	MATTRESS INNERSPRING	I	PROSVENDOR.FIVE	605		E0271	270.00
12/27/02	1	LIGHTWEIGHT WHEELCHAIR	I	PROSVENDOR.SEVE	605		K0003	296.03
12/27/02	1	GEL PRESSRE PAD/CUSHION NOMP	I	PROSVENDOR.NINE	605		E0178	54.75
12/27/02	1	DME DELIVERY, SET-UP	X	PROSVENDOR.FIVE	605		A9901	72.15
12/27/02	1	DME DELIVERY, SET-UP	X	PROSVENDOR.FIVE	605		A9901	72.15
12/27/02	1	DME DELIVERY, SET-UP	X	PROSVENDOR.FIVE	605		A9901	72.15
12/27/02	1	HOSP BED SEMI-ELECTR W/ MATT	I	PROSVENDOR.FIVE	605	89002L35456	E0260	575.21

Buttons: View Detail, Close. Total records found -> 33

Next Step

Click the **View Detail** button at the bottom of this window to display the **Appliance Transaction Detail** window and proceed to the next page.

Continued on next page

GUI 2319 Button on DOR Menu – Appliance Transactions (Tab 4), Continued

Appliance Field The **Appliance** field listed in the **Transaction Detail** window below is protected from changes to IFCAP (from the NIF group’s data cleansing).

The example below, “Wheelchair 9000XT 18-inch” in the **Appliance** field will remain the same with Patch RMPR*3*99.

Note: The **Appliance** field below will always display the originally selected **IFCAP Item Description**. The new “**Item Description**: field in Roll and Scroll and the new note “*See Above...*” is not necessary in the GUI 2319. This is because the GUI 2319 **IFCAP Item Description** will always be displayed and protected from changes.

Transaction Detail window

The screenshot shows the 'Transaction Detail' window with the following data:

Appliance Item		Lab		IFCAP Vendor		Historical Data	
Appliance:	WHEELCHAIR 9000XT 18-INCH	Order #:		Name:	PROVIDOR SEVEN	Name:	
Description:	9000XT 18" Recliner Wheelchair	Completed:		Addr:	899 A WY ST	Addr:	
HCPCS:	K0003 Qty: 1	Tech:		City:	ELYRIA	City:	
HCPCS Desc:	LIGHTWEIGHT WHEELCHAIR	Lab Cost:		State:	OHIO Zip: 44036-2125	State:	
ICD-9 Code:	203.00 CPT Modifier: NU Disability: SC/OP	Material Cost:		Phone:	555/555-5555	Phone:	
Invent Point:		Labor Hrs:		Track #:		Item:	
Source:	COMMERCIAL	Labor Cost:				Station:	
Deliver To:	Ink medical Completion Date: JAN 08, 2003	Ortho Code:					
Recv Stat:		Restor Code:					
Return Stat:		Lab Remarks:					
Remarks:							
Extended Description							
w/Swingaway Foot rest, Composite footrest, Std. caster, Black, Adjustable Desk Arm							


Close button To close the **Transaction Detail** window and return to the **View 2319** window, click the **Close** button.

GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8)

Home Oxygen (Tab 8)

The 2319 in the Delayed Order Report (DOR) shown below of the **Home Oxygen (Tab 8)** window shows the **Item** column which is the **HCPCS Description** with Patch RMPR*3*99.

Home Oxygen window



Date	Qty	Item	Typ	Vendor	Station	Serial	HCPC	Total Cost
01/15/04	1	PORTABLE GASEOUS O2	X	PROSVENDOR.TEN	605		E0431	50.00
12/22/03	1	OXYGEN CONCENTRATOR	X	PROSVENDOR.TEN	605		E1390	83.00
12/22/03	1	PORTABLE GASEOUS O2	X	PROSVENDOR.TEN	605		E0431	50.00
11/14/03	1	OXYGEN CONCENTRATOR	X	PROSVENDOR.TEN	605		E1390	83.00
11/14/03	1	PORTABLE GASEOUS O2	X	PROSVENDOR.TEN	605		E0431	50.00
09/03/03	1	OXYGEN CONCENTRATOR	I	PROSVENDOR.TEN	605		E1390	49.00
09/03/03	1	PORTABLE GASEOUS O2	I	PROSVENDOR.TEN	605		E0431	30.00
09/03/03	1	OXYGEN CONCENTRATOR	I	PROSVENDOR.TEN	605		E1390	83.00
09/03/03	1	PORTABLE GASEOUS O2	I	PROSVENDOR.TEN	605		E0431	50.00

Click on Header to SORT

View Detail Close Total records found -> 9

View Detail button

Click the **View Detail** button at the bottom of this window to display the **Home Oxygen Detail** and proceed to the next page.

Continued on next page

GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8), Continued

Home Oxygen (Tab 8) Transaction Detail

Notice the **Appliance** field on the **Transaction Detail** window from the **Home Oxygen (Tab 8)** remains the same, even when changes are made to the NIF IFCAP Item Description.

Note: The **Appliance** field below will always display the originally selected **IFCAP Item Description**. The new “**Item Description:** field in Roll and Scroll and the new note “*See Above...*” is not necessary in the GUI 2319. This is because the GUI 2319 **IFCAP Item Description** will always be displayed and protected from changes.

Transaction Detail

The screenshot shows the 'Transaction Detail' window with the following fields:

- Appliance Item:** Appliance: RENTAL MONTHLY PORTABLE OXYGEN SYSTEMS; Description: [Empty]; HCPCS: E0431 Qty: 1; HCPCS Desc: PORTABLE GASEOUS O2; ICD-9 Code: [Empty] CPT Modifier: RP Disability: SCOP; Invent Point: [Empty]; Source: COMMERCIAL; Deliver To: [Empty] Completion Date: JAN 24, 2004; Recv Site: [Empty]; Return Stat: [Empty]; Remarks: [Empty]
- Lab:** Order #: [Empty]; Completed: [Empty]; Tech: [Empty]; Lab Cost: [Empty]; Material Cost: [Empty]; Labor Hrs: [Empty]; Labor Cost: [Empty]; Ortho Code: [Empty]; Restor Code: [Empty]; Lab Remarks: [Empty]
- Transaction Info:** Date: JAN 15, 2004@; Transaction: P49082; Type Trans: REPAIR; Form Type: VISA; Initiator: PROSPROVIDER.FOUR; Site: LOMA LINDA VAMC; Bank Auth: [Empty]; Total Cost: \$50.00
- IFCAP Vendor:** Name: PROSPROVIDER.TEN; Add: 935 S. Milliken Ave Ste. B; City: Ontario; State: CALIFORNIA; Zip: 91761; Phone: 555-555-5555; Track #: [Empty]
- Historical Data:** Name: [Empty]; Add: [Empty]; City: [Empty]; State: [Empty]; Zip: [Empty]; Phone: [Empty]; Item: [Empty]; Station: [Empty]
- Extended Description:** [Empty]

Appendix A

Glossary

DOR	Delayed Order Report – Prosthetics application.
------------	---

GUI	Graphical User Interface. This is a windows-based application with buttons and tabs instead of prompts in a Roll and Scroll environment for an application. Microsoft Windows is a GUI environment.
------------	---

HCPCS	Healthcare Common Procedure Coding System – formerly Health Care Financing Administration Common Procedures Coding System - is distributed by Centers for Medicare & Medicaid Services (CMS) and the Health and Human Services (HHS) for all other substances, equipment, supplies, or other items used in Health Care Services. These items include, but are not limited to: medical supplies; orthotic and prosthetic devices; and durable medical equipment.
--------------	---

NIF	<p>The National Item File (NIF) database is used to standardize naming practices so that facilities will be able to identify, track and evaluate purchases.</p> <p>The NIF is also a group that works on this database file.</p>
------------	--

VistA	Veterans Health Information Systems and Technology Architecture, formerly known as Decentralized Hospital Computer Program, encompasses the complete information environment at VA medical facilities. It consists of hardware, software packages, and comprehensive support for system-wide and station specific, clinical and administrative automation needs.
--------------	--
