

Department of Veterans Affairs  
Decentralized Hospital Computer Program

# **EQUIPMENT/TURN-IN REQUEST TECHNICAL MANUAL**

Version 1.0  
June 1996

Information Resource Management Field Office  
Washington, DC







# Table of Contents

<b>CHAPTER 1 INTRODUCTION.....</b>	<b>1</b>
1.1 OVERVIEW.....	1
1.2 ELECTRONIC SIGNATURE.....	1
<b>CHAPTER 2 IMPLEMENTATION AND MAINTENANCE.....</b>	<b>3</b>
2.1 OPERATING SPECIFICS.....	3
2.1.1 SIZING INFORMATION.....	3
2.2 GLOBAL JOURNALING.....	3
2.3 TRANSLATION TABLES.....	3
<b>CHAPTER 3 ROUTINE DESCRIPTIONS.....</b>	<b>5</b>
3.1 OVERVIEW.....	5
3.2 LIST OF ROUTINES.....	5
<b>CHAPTER 4 FILE LIST.....</b>	<b>9</b>
4.1 OVERVIEW.....	9
4.2 FILE LIST (BY NUMBER).....	9
4.3 FILE PROTECTION.....	10
4.4 FILES WITH SECURITY ACCESS.....	10
<b>CHAPTER 5 EXPORTED OPTIONS.....</b>	<b>11</b>
5.1 MENU STRUCTURE AND OPTION DEFINITIONS.....	11
5.2 ASSIGNING MENUS AND SECURITY KEYS TO USERS.....	11
<i>Equipment Request Menus.....</i>	<i>11</i>
<i>Turn-In Menus.....</i>	<i>12</i>
5.3 OPTION LIST BY NAME.....	12
5.4 OPTION LIST (BY MENU).....	25
5.4.1 Menu Diagram for PRCN NX MENU.....	25
5.4.2 Menu Diagram for PRCN NX CMR MENU.....	26
5.4.3 Menu Diagram for PRCN NX PPM MENU.....	26
5.4.4 Menu Diagram for PRCN NX EQ COMM MENU.....	27
5.4.5 Menu Diagram for PRCN NX ENG MENU.....	27
5.4.6 Menu Diagram for PRCN NX CONCUR MENU.....	27
5.4.7 Menu Diagram for PRCN TURN WHSE MENU.....	27
<b>CHAPTER 6 CROSS REFERENCES.....</b>	<b>29</b>
<b>CHAPTER 7 ARCHIVING AND PURGING.....</b>	<b>36</b>
7.1 ARCHIVING.....	36
7.2 PURGING.....	36
<b>CHAPTER 8 CALLABLE ROUTINES.....</b>	<b>39</b>
<b>CHAPTER 9 EXTERNAL RELATIONS.....</b>	<b>41</b>
<b>CHAPTER 10 INTERNAL RELATIONS.....</b>	<b>49</b>
<b>CHAPTER 11 PACKAGE-WIDE VARIABLES.....</b>	<b>51</b>
<b>CHAPTER 12 ON-LINE DOCUMENTATION.....</b>	<b>53</b>



# Revision History

Initiated on 12/29/04

Date	Description (Patch # if applic.)	Project Manager	Technical Writer
12/29/04	Updated to comply with SOP 192-352 Displaying Sensitive Data.		Mary Ellen Gray
12/29/04	Pdf file checked for accessibility to readers with disabilities.		Mary Ellen Gray















**PRCNPURG**

This routine is used to purge completed or cancelled requests.

**PRCNREQE**

This routine is used to edit a equipment request.

**PRCNREQN**

This routine is used to enter a new equipment request.

**PRCNREQQ**

This routine is used to process changes in the line item quantity.

**PRCNRNK**

This routine is used to rank the equipment requests by the PPM.

**PRCNRPT**

This routine is used to print some reports.

**PRCNRSMT**

This routine is used to resubmit a request.

**PRCNSERP**

This routines prints the Service Priority report.

**PRCNSPL**

This routine is used to split an equipment request transaction.

**PRCNSTAT**

This routine is used to print the Request Status Report.

**PRCNTICM**

This routine is used by the CMR Responsible Official to review and approve a turn-in request.

**PRCNTICN**

This routine is used to cancel a turn-in request and associated equipment request.

**PRCNTIED**

This routine is used to edit a turn-in request.

**PRCNTIPP**

This routine is used by the Personal Property Manager (PPM) to review and complete a turn-in request.













Name: PRCN NX CMR DISPLAY  
 Type: Run Routine Routine: CMR^PRCNDISP  
 Lock: PRCNCMR  
 Entry Action:  
 Exit Action:  
 Menu Text: Display/Print Equipment Request (CMR)  
 Description: This option allows the CMR to display to screen or print equipment requests.

Name: PRCN NX CMR MENU  
 Type: Menu Routine:  
 Lock: PRCNCMR  
 Entry Action: S PRCNCMF=1  
 Exit Action: K STA,PRCNCMF  
 Header: S STA=3 D EMSG^PRCNUTL S STA=45 D  
 EMSG^PRCNUTL,TMSG^PRCNUTL  
 Menu Text: CMR Official Equipment Request Menu  
 Description: This is the top level master menu for the CMR Official.  
 Menu:

Enter New Equipment Request	PRCN NX NEW
Edit Equipment Request	PRCN NX EDIT
Cancel Equipment Request	PRCN NX CANCEL
Request Status Report	PRCN NX CMR STATUS
Display/Print Equipment Request (CMR)	PRCN NX CMR DISPLAY
Approve Equipment Requests (CMR Official)	PRCN NX CMR APP
Resubmit Request (CMR)	PRCN NX CMR RESUBMIT
Process Equipment Turn-Ins Menu ...	PRCN TURN MENU

Name: PRCN NX CMR RESUBMIT  
 Type: Run Routine Routine: EN^PRCNRSMT  
 Lock: PRCNCMR  
 Enter Action: S PRCNCMF=1  
 Exit Action: K PRCNCMF  
 Menu Text: Resubmit Request (CMR)  
 Description: This option allows the responsible CMR Official to resubmit transactions with a particular status.

Name: PRCN NX CMR STATUS  
 Type: Run Routine Routine: CMR^PRCNSTAT  
 Lock: PRCNCMR  
 Entry Action:  
 Exit Action:  
 Menu Text: Request Status Report  
 Description: This option will print the status of all requests pertaining to this person's CMR responsibilities.



Name: PRCN NX EDIT  
 Type: Run Routine Routine: EN^PRCNREQE  
 Lock:  
 Entry Action:  
 Exit Action:  
 Menu Text: Edit Equipment Request  
 Description: This option allows the user to review and edit their equipment request, as long as the CMR Official has not yet approved it.

Name: PRCN NX ENG  
 Type: Run Routine Routine: PRCNENG  
 Lock: PRCNEN  
 Entry Action:  
 Exit Action:  
 Menu Text: Approve Equipment Request (Engineering)  
 Description: This option should be assigned to those Engineering personnel who may need to complete important information about equipment requests.

Name: PRCN NX ENG MENU  
 Type: Menu Routine:  
 Lock: PRCNEN  
 Entry Action: F STA=8,26,30 D EMSG^PRCNUTL  
 Exit Action:  
 Header:  
 Menu Text: Engineering Equipment Request Menu  
 Description: This menu option is the main menu option for Engineering personnel who may review and process non-expendable equipment requests.  
 Menu:  
     Approve Equipment Request (Engineering) PRCN NX ENG

Name: PRCN NX EQ COMM APPROVE  
 Type: Run Routine Routine: EN^PRCNEQAP  
 Lock: PRCNEQP  
 Entry Action:  
 Exit Action:  
 Menu Text: Process Equip Committee Decisions  
 Description: This menu option is to be assigned to either an Equipment Committee member or to the PPM to process the decisions of the Equipment Committee. After the Equipment Committee has reviewed equipment requests, their decisions must be entered.



Menu Text: Equipment Request Menu (Requestor)  
 Description: This is the main menu for the requestor.  
 Menu:

Enter New Equipment Request	PRCN NX NEW
Edit Equipment Request	PRCN NX EDIT
Cancel Equipment Request	PRCN NX CANCEL
Request Status Report	PRCN NX STATUS
Display/Print Equipment Request	PRCN NX DISPLAY
Split Equipment Request	PRCN NX SPLIT
Resubmit Equipment Requests	PRCN NX RESUBMIT
Process Equipment Turn-Ins Menu	PRCN TURN MENU

Name: PRCN NX NEW  
 Type: Run Routine Routine: EN^PRCNREQN  
 Lock:  
 Entry Action:  
 Exit Action:  
 Menu Text: Enter New Equipment Request  
 Description: This option is to allow a user to enter a new equipment request for processing.

Name: PRCN NX OUTSTAND REPORT  
 Type: Print Routine:  
 Lock:  
 Entry Action:  
 Exit Action: I IOT="TRM" W !,"Press RETURN to continue..." R C:DTIME  
 DIC {DIP}: PRCN(413, L.: 0  
 FLDS: [PRCNBRIE] BY: [PRCNSRPR]  
 DIS(0): I \$P(^PRCN(413,D0,0),U,7)=31!(\$P^(0),U,7)=10)  
 Menu Text: Outstanding Equipment Request Report  
 Description: This report prints a list of requests by service and priority that have not been processed yet.









Name: PRCN NX TURN BY CMR  
 Type: Print Routine: CMR^PRCNRPT  
 Lock:  
 Entry Action:  
 Exit Action: R:IOT="TRM" !!,"Press RETURN to continue...", C:DTIME  
 Menu Text: Turn-In Status by CMR Official Report  
 Description: This report prints the status of turn-in requests by CMR Official.

Name: PRCN NX TURN BY SERVICE  
 Type: Print Routine: SER^PRCNRPT  
 Lock:  
 Entry Action:  
 Exit Action: R:IOT="TRM" !!,"Press RETURN to continue...", C:DTIME  
 Menu Text: Turn-In Status by Service Report  
 Description: This report prints the status of turn-ins by service.

Name: PRCN NX TURN BY STATUS  
 Type: Print Routine:  
 Lock:  
 Entry Action:  
 Exit Action: R:IOT="TRM" !!,"Press RETURN to continue...", C:DTIME  
 DIC {DIP}: PRCN(413.1, L.: 0  
 FLDS: [PRCNTIST] BY: 6  
 FR: , TO: ,  
 Menu Text: Turn-In Status Report  
 Description: This report prints the status of turn-ins sorted by status.

Name: PRCN NX TURN BY TRANS  
 Type: Print Routine:  
 Lock:  
 Entry Action:  
 Exit Action: R:IOT="TRM" !!,"Press RETURN to continue...", C:DTIME  
 DIC {DIP}: PRCN(413.1, L.: 0  
 FLDS: [PRCNTIST] BY: .01  
 FR: , TO: ,  
 Menu Text: Turn-In Status by Transaction Number Report  
 Description: This report prints the status of turn-ins sorted by transaction number.

























Kill: K ^PRCN(413.2,DA(1),1,"B",\$E(X,1,30),DA)  
 Desc: This cross-reference is by the scheduled meeting date.

---

-----  
 ----  
**EQUIPMENT COMMITTEE - File 413.2**

(413.3,.01) NAME 0;1

Xref 1: 413.3^B  
 Set: S ^PRCN(413.3,"B",\$E(X,1,30),DA)=""  
 Kill: K ^PRCN(413.3,"B",\$E(X,1,30),DA)  
 Desc: This cross-reference is by concurring official.

Xref 2:TRIGGER^413.3^1  
 Create: S X=\$P(^VA(200,X,0),U,9)  
 Delete: NO EFFECT  
 Field: #1  
 Desc: This triggers the title from the New Person file.

---

-----  
 ----  
**SPECIAL HANDLING CODES - File 413.4**

(413.4,.01) CODE 0;1

Xref 1: 413.4^B  
 Set: S ^PRCN(413.4,"B",\$E(X,1,30),DA)=""  
 Kill: K ^PRCN(413.4,"B",\$E(X,1,30),DA)  
 Desc: This cross-reference is by the special handling code.

---

-----  
 ----  
**NX STATUS - File 413.5**

(413.5,.01) STATUS 0;1

Xref 1: 413.5^B  
 Set: S ^PRCN(413.5,"B",\$E(X,1,30),DA)=""  
 Kill: K ^PRCN(413.5,"B",\$E(X,1,30),DA)  
 Desc: This cross-reference is by the status.

---

-----  
 ----  
**COUNTER - File 413.7**

(413.7,.01) CODE 0;1













## CHAPTER 8 CALLABLE ROUTINES

Please see Chapter 9 External Relations for information on access to other packages.













**Integration Reference #1523**

NAME: DBIA1520-D	ENTRY: 1523
CUSTODIAL PACKAGE: IFCAP	Washington
SUBSCRIBING PACKAGE: EQUIPMENT/TURN-IN REQUEST	Washington
USAGE: Private	APPROVED: APPROVED
STATUS: Active	EXPIRES:
DURATION: Till Otherwise Agr	VERSION:
	TYPE: Routine

This agreement will allow the NX (Equipment/Turn-In) module to call the IFCAP transaction utility program when creating 2237s which will create the record in file 410 and process all checks on creating a 2237.

ROUTINE: PRCSUT3

**Integration Reference #1524**

NAME: DBIA1520-E	ENTRY: 1524
CUSTODIAL PACKAGE: IFCAP	Washington
SUBSCRIBING PACKAGE: EQUIPMENT/TURN-IN REQUEST	Washington
USAGE: Private	APPROVED: APPROVED
STATUS: Active	EXPIRES:
DURATION: Till Otherwise Agr	VERSION:
FILE: 410	ROOT: PRCS(410, TYPE: File

The NX (Equipment/Turn-In) module requests permission to reference file 410 to create/edit 2237s which are the end product of this module. Includes addition of two templates to file 410, PRCN2237 and PRCN2237E.

**Integration Reference #1525**

NAME: DBIA1520-F	ENTRY: 1525
CUSTODIAL PACKAGE: IFCAP	Washington
SUBSCRIBING PACKAGE: EQUIPMENT/TURN-IN REQUEST	Washington
USAGE: Private	APPROVED: APPROVED
STATUS: Active	EXPIRES:
DURATION: Till Otherwise Agr	VERSION:
FILE: 440	ROOT: PRC(440, TYPE: File

This agreement requests permission for the NX (Equipment/Turn-In) module to point, with read access only, to the Vendor file (440).













# CHAPTER 11 PACKAGE-WIDE VARIABLES

Currently there are no package-wide variables used in the Equipment/Turn-In Request package.









